



ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA
CTO USING THE SEISMIQ™ INTRAVASCULAR
LITHOTRIPSY SYSTEM AND THE ELUVIA™
DRUG-ELUTING VASCULAR STENT SYSTEM

PRESENTED BY:

LIZ GENOVESE, M.D., M.S., F.A.C.S.

- assistant professor of surgery
- director, Penn Advanced Limb Preservation Program at HUP
- Division of Vascular Surgery and Endovascular Therapy
- University of Pennsylvania

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

HISTORY

PATIENT DEMOGRAPHICS

71-year-old patient

PRESENTING SYMPTOMS

Rapidly progressive left lower extremity rest pain and tissue loss

MEDICAL HISTORY

Dialysis-dependent end-stage renal disease

- Diabetes mellitus

Multilevel peripheral arterial disease

High likelihood of severe arterial calcification

NONINVASIVE VASCULAR ASSESSMENT

Markedly diminished arterial waveforms

Toe pressure:

- 20 mmHg
- Flatline contour

CLINICAL INTERPRETATION

Clinical presentation consistent with critical limb threatening ischemia (CLTI)

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

DIAGNOSTIC FINDINGS

CROSS-SECTIONAL IMAGING

CT scan findings:

Common femoral artery: widely patent

Calcification increased upon entering SFA

FEMOROPOPLITEAL SEGMENT

Superficial femoral artery (SFA):

- Dense, nodular, circumferential, boulder-like calcium

Above-the-knee popliteal artery:

- Severe calcification
- Remained patent despite disease burden

DISTAL VESSEL ASSESSMENT

Tibial runoff:

- Significant multilevel disease
- Required additional below-the-knee intervention

OVERALL ASSESSMENT

Long, severely calcified chronic total occlusion (CTO)

Complex anatomy with limited distal targets

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES






DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

DIAGNOSTIC FINDINGS

INITIAL ANGIOGRAM



HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

TREATMENT APPROACH

LESION CHARACTERISTICS

Long segment CTO

Severe, nodular, circumferential calcification

KEY PROCEDURAL RISKS

High risk of:

- Dissection
- Vessel recoil
- Underexpansion
- Restenosis

ANATOMIC CHALLENGES

Dense calcium significantly limited:

- Lesion crossing
- Vessel compliance

SURGICAL CONSIDERATIONS

Open bypass highly morbid and unfavorable due to:

- Renal failure
- Poor distal targets
- **Limitations of alternative endovascular strategies**

“Leave nothing behind” approaches deemed unsuitable due to:

- High risk of recoil
- Residual stenosis
- Restenosis in heavily calcified disease

TREATMENT STRATEGY GOALS

Effective calcium modification

Creation of a durable flow channel

Long-term suppression of neointimal hyperplasia

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

TREATMENT APPROACH | CALCIUM MODIFICATION STRATEGY

SELECTED TECHNOLOGY

SEISMIQ™ Intravascular Lithotripsy (IVL) System

RATIONALE FOR SEISMIQ™ IVL SELECTION

Directionally align emitters to calcified segments

Select specific emitters (on/off) to preserve pulses

Precisely target acoustic energy

Uniform energy delivery with laser-based IVL technology

IVL TREATMENT TECHNIQUE

Delivery from below the knee upward

Focused treatment of:

- Nodular distal SFA lesion
- Circumferentially calcified proximal SFA

RESULT OF CALCIUM MODIFICATION

Improved vessel compliance

Creation of a functional flow channel following dissection plane crossing

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series




CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMOIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

TREATMENT APPROACH | CALCIUM MODIFICATION STRATEGY

FOCUSING ON
NODULAR CALCIUM

TREATING AREAS OF
CIRCUMFERENTIAL CALCIUM

EMITTER SELECTION TO
UTILIZE ALL PULSES AT
DIFFICULT AREAS



HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

TREATMENT APPROACH | DEFINITIVE TREATMENT STRATEGY

POST PREPARATION ASSESSMENT

Persistent high risk lesion morphology

Long lesion length

Presence of dissection

THERAPEUTIC GOAL

In heavily calcified anatomy:

- Maintain durable luminal gain
- Prevent elastic recoil
- Suppress neointimal hyperplasia

DEVICE STRATEGY

Use of a drug-eluting stent scaffold to:

- Maintain vessel expansion
- Address the high likelihood of restenosis with uncoated balloon angioplasty
- Deliver sustained antiproliferative therapy

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

IMPLANT PROCEDURE

STENT DEPLOYMENT

Vessel predilation

Three ELUVIA™ Drug-Eluting Vascular Stents were deployed

Achieved full coverage of long diseased segment

OPTIMIZATION

Postdilatation performed to:

- Optimize stent expansion
- Ensure adequate apposition

ADJUNCTIVE TREATMENT

Additional below-the-knee intervention:

- Addressed outflow disease
- Supported procedural durability

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES




DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

IMPLANT PROCEDURE

ANGIO AFTER
IVL ALONE



PTA WITH NC BALLOON



ELUVIA™
STENTING



HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

TREATMENT RESULTS & OUTCOMES

COMPLETION ANGIOGRAPHIC RESULTS

Excellent luminal gain

Complete resolution of:

- Balloon “wasting” observed prior to IVL

No significant residual stenosis


RUNOFF ASSESSMENT

Strong two-vessel runoff preserved


TECHNICAL OUTCOME

Durable expansion of stented SFA/popliteal segment

COMPLETION ANGIO



SFA/POP



HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES

*Results from a single case. Individual outcomes may vary.



TREATMENT RESULTS & OUTCOMES

CLINICAL FOLLOW UP

Six-week follow-up demonstrated:

- Positive tissue healing trajectory
- Early confirmation of treatment effectiveness



CLINICAL SIGNIFICANCE

Highlights the critical role of the ELUVIA™ Drug-Eluting Vascular Stent System in:

- Long, heavily calcified femoropopliteal lesions
- CLTI patients with limited surgical options

INTEGRATED TREATMENT STRATEGY

SEISMIQ™ IVL:

- Enabled precise and effective vessel preparation

ELUVIA™ Drug-Eluting Vascular Stents:

- Provided sustained, controlled drug elution
- Maintained long-term luminal expansion
- Suppressed neointimal hyperplasia

PATIENT OUTCOME

- Durable revascularization
- Favorable early healing in a highly comorbid patient
- Successful outcome in a complex, calcified CTO

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES

*Results from a single case. Individual outcomes may vary.



DE Clinical Spotlight Series

CASE 2 | TREATING A HIGHLY CALCIFIED LONG SFA CTO USING THE SEISMIQ™ INTRAVASCULAR LITHOTRIPSY SYSTEM AND THE ELUVIA™ DRUG-ELUTING VASCULAR STENT SYSTEM

**Boston
Scientific**

Advancing science for life™

HISTORY

DIAGNOSTIC FINDINGS

TREATMENT APPROACH

CALCIUM MODIFICATION STRATEGY

DEFINITIVE TREATMENT STRATEGY

IMPLANT PROCEDURE

TREATMENT RESULTS & OUTCOMES

**Boston
Scientific**
Advancing science for life™

**Interventional Oncology
& Embolization**

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2026 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-2492814-AA