

# IMPERIAL RCT SUMMARY

The world's first head-to-head DES SFA Trial, evaluating Boston Scientific Corporation's Eluvia™ Drug-Eluting Vascular Stent System and Cook Medical's Zilver™ PTX™ Stent




## Sustaining strong results through five years

The results of the IMPERIAL RCT show that Eluvia Drug-Eluting Stent (DES) is clinically effective and safe in treating patients with symptomatic SFA disease both in the short-term during the height of restenosis risk, and long-term out to five years.

Eluvia DES demonstrated **superiority over Zilver PTX<sup>1</sup>** with a statistically significant primary patency through 1-Year


### IMPERIAL RCT 1-YEAR RESULTS

K-M Primary Patency<sup>2</sup>



Eluvia DES showed lower revascularization rates than Zilver PTX through 5 years<sup>3</sup> with **statistical significance<sup>4</sup>** at 2-Years

### FREEDOM FROM CD-TLR RATES



### IMPERIAL TRIAL 2-YEAR CLINICAL RESULTS

Excellent Patient Follow-up at 24-Months (~90%)

**Durable, consistent outcomes in long and complex lesions**

|                                       | IMPERIAL RCT <sup>5</sup><br>(n = 309) | IMPERIAL Long Lesions <sup>6,7</sup><br>(n = 56) | Diabetic Subgroup <sup>2,8</sup><br>(n = 116) | Severe/Moderate Calcium Subgroup <sup>2</sup><br>(n = 193) | CTO Subgroup <sup>2</sup><br>(n = 96) |
|---------------------------------------|--|--|---|--|---------------------------------------|
| Study Design                          | RCT, global multicenter                | Single arm, global multicenter                   | RCT, global multicenter                       | RCT, global multicenter                                    | RCT, global multicenter               |
| <b>24-month primary patency rate*</b> | <b>83.0%</b>                           | <b>77.2%</b>                                     | <b>85.7%</b>                                  | <b>85.0%</b>   | <b>76.4%</b>                          |
| Lesion length (mm)                    | 86.5                                   | 162.8  | 87.0  | 89.9   | 94.4                                  |
| Severe calcification                  | 40%                                    | 28%  | 46%   | n/a  | n/a                                   |
| Total occlusions                      | 31%                                    | 32%  | 25%   | n/a  | 100%                                  |

Highest primary patency ever reported at 2 years\*\*

Highly durable outcomes in ~16cm lesions at 2 years

TLR (12%) in line with overall cohort and low stent thrombosis rate (0.9%)

Remarkable primary patency and <10% TLR in heavy calcium

Highly durable outcomes in CTOs at 2 years

**Eluvia DES patients  
on average avoided  
reintervention 6 months  
longer than Zilver PTX  
patients at 3-Years<sup>3†</sup>**

**Zilver PTX**

**13 months**

p=0.0058

**Eluvia DES**

**19 months**

months ..... 6 ..... 12 ..... 18 ..... // ..... 36

## IMPERIAL TRIAL OBJECTIVE

Evaluate the safety and effectiveness of the Boston Scientific Corporation Eluvia™ Drug-Eluting Vascular Stent System for treating Superficial Femoral Artery (SFA) and/or Proximal Popliteal Artery (PPA) lesions up to 140mm in length.

## IMPERIAL TRIAL DESIGN

Global multi-center, 2:1 randomization against Cook Medical's Zilver™ PTX™ Stent, controlled, single-blind, non-inferiority trial; core lab adjudicated.

- 465 (RCT) patients across 64 sites
- 5-year follow-up
- Degree of stenosis ≥ 70% (visual angiographic assessment)
- Vessel diameter ≥ 4mm and ≤ 6mm
- Total lesion length ≥ 30mm and ≤ 140mm

## BASELINE CHARACTERISTICS

| Patient Demographics | Eluvia (n=309) | Zilver PTX (n=156) |
|----------------------|----------------|--------------------|
| Age (Years)          | 68.5±9.5       | 67.8±9.4           |
| Male Gender          | 66.0%          | 66.7%              |
| Diabetes Mellitus    | 41.7%          | 43.6%              |
| History of Smoking   | 86.1%          | 84.0%              |

| Lesion Characteristics    | Eluvia (n=309) | Zilver PTX (n=156) |
|---------------------------|----------------|--------------------|
| Target Lesion Length (mm) | 86.5±36.9      | 81.8±37.3          |
| Severely Calcified        | 40.1%          | 32.3%              |
| Total Occlusions          | 31.2%          | 30.3%              |
| Extending into Distal SFA | 66.3%          | 65.4%              |

\* Intention to treat. Kaplan-Meier estimate utilizing time-to-event of clinically-driven TLR up to 730 days and Duplex Ultrasound data at 24 months. Primary patency defined as duplex ultrasound PSVR ≤2.4, in the absence of clinically-driven target lesion revascularization or bypass of the target lesion, as assessed by the DUS core lab.

\*\*Highest-two-year primary patency based on 24-month Kaplan-Meier estimates reported for IMPERIAL, IN.PACT SFA, ILLUMINATE, LEVANT II and Primary Randomization for Zilver PTX RCT.

† Among patients who underwent a CD-TLR within 3 years of the index procedure

1. IMPERIAL Trial: A global randomized controlled multi-center trial with 2:1 randomization of the Eluvia™ Drug-Eluting Stent against Cook Medical's Zilver™ PTX™ Stent, single-blind, non-inferiority design; independent core lab adjudication. Superiority determined in a post hoc analysis that was specified prior to unblinding. 12-Month Primary Patency rate of 86.8% in the Eluvia arm vs. 77.5% in the Zilver PTX arm (p-value = 0.0144). Gray WA, Lancet. 2018 Sep 24. pii: S0140-6736(18)32262-1.

2. Gray, W. 2 year Outcomes from the IMPERIAL Randomized Head to Head Study of Eluvia DES and ZilverPTX. LINC 2020.

3. Gray W. 5-year Results from the IMPERIAL Randomized Study of Eluvia and Zilver PTX Drug-eluting Stents and Long Lesion Substudy for Femoropopliteal Artery Disease; CRT 2023, Washington DC Feb 27, 2023.

4. Müller-Hülsbeck S, et al. Two-Year Efficacy and Safety Results from the IMPERIAL Randomized Study of the Eluvia Polymer-Coated Drug-Eluting Stent and the Zilver PTX Polymer-free Drug-Coated Stent. Cardiovasc Intervent Radiol. 2021;44(3):368-375.

5. In IMPERIAL RCT, Eluvia K-M Primary Patency was 83% vs. 77.1% for Zilver PTX at 24 months, p=0.1008.

6. Golzar, J. et al, et al. Effectiveness and Safety of a Paclitaxel-Eluting Stent for Superficial Femoral Artery Lesions up to 190 mm: One-Year Outcomes of the Single-Arm IMPERIAL Long Lesion Substudy of the Eluvia Drug-Eluting Stent. J Endovasc Ther. 2020;27(2):296-303.

7. Vermassen, F. VIVA Late-Breaking

8. Müller-Hülsbeck S. LINC 2019

**Boston  
Scientific**

Advancing science for life™

### Peripheral Interventions

300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[bostonscientific.com](http://bostonscientific.com)

To order product or for more  
information contact customer  
service at 1.888.272.1001

© 2025 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-1528409-AD



**ELUVIA™ Drug-Eluting Stent  
Indications, Safety, and Warnings**