



## Prospective, Multicenter Evaluation of Transcarotid Artery Revascularization (TCAR) in Standard Risk Patients: 30-Day Outcomes of the ROADSTER 3 Study

*Dermody, M. (2024). Prospective, Multicenter Evaluation of Transcarotid Artery Revascularization (TCAR) in Standard Risk Patients: 30-Day Outcomes of the ROADSTER 3 Study. Presented at VIVA Conference 2024.*

### Objective

First-ever prospective, multicenter, trial evaluating real world usage of TCAR using the ENROUTE® Transcarotid Stent System (TSS) in conjunction with the ENROUTE® Transcarotid Neuroprotection System (NPS) for the treatment of carotid stenosis in standard surgical risk patients.

### Methods

Between 2022–2024, 344 intent-to-treat (ITT) patients were enrolled (320 treated per-protocol (PP)) at 53 US sites. The primary endpoint for this single-arm, post-approval study is a composite of major adverse events (stroke, death, or myocardial infarction (MI)) through 30 days post-procedure, plus ipsilateral stroke from day 31 to 365 post-procedure. The incidence of cranial nerve injury (CNI) within 30 days post-procedure is a key secondary endpoint. Independent neurological assessments are performed for all patients before the procedure, within 24 hours, at 30-days, and at 1-year after TCAR. Events were adjudicated by an independent clinical events committee.

### Findings

In the ITT population, 75.3% were less than 75 years of age, 42.7% were female, and 16.3% were symptomatic. Among symptomatic patients, 25.0% experienced a neurologic event within 2 weeks preceding the TCAR procedure. The mean lesion length was 23.3mm, 47.4% had a Type II or Type III aortic arch, and 17.2% of lesions had severe calcification.

In the ITT population, the rate of stroke/death/MI at 30-days was 0.9% (0.6% PP) with a 30-day stroke rate of 0.9% (0.6% PP, n=2). There were no deaths or MIs through 30-day follow-up. The incidence of CNI within 30 days was 0.6% (0.6% PP); both resolved within 6 months.

### Conclusions

These 30-day results of the ROADSTER 3 study demonstrate that TCAR, using the ENROUTE TSS in conjunction with the ENROUTE NPS, is safe and effective in patients at standard risk for adverse events from carotid endarterectomy.

	ITT n=344	PP n=320
Stroke	0.9% (3)	0.6% (2)
Death	None	
MI	None	
Stroke/Death/MI	0.9% (3)	0.6% (2)



