



THERASPHERE™ Y-90 Glass Microspheres | PUBLICATION SUMMARY

Radiation Segmentectomy: Potential Curative Therapy for Early Hepatocellular Carcinoma

Radiology, 2018 Jun; 287(3); <https://pubs.rsna.org/doi/10.1148/radiol.2018171768> Lewandowski RJ, Gabr A, Abouchaleh N, Ali R, Al Asadi A, Mora RA, Kulik L, Ganger D, Desai K, Thornburg B, Mouli S, Hickey R, Caicedo JC, Abecassis M, Riaz A and Salem R

OVERVIEW

Curative treatment options for early-stage HCC (BCLC 0 or A) include transplantation, surgical resection and RF ablation with good survival outcomes ranging between 60-80%^{*1}, however many patients are not candidates for these therapies

- Radiation segmentectomy is the application of selective ablative radiation doses of Yttrium-90 (Y-90) to tumors, usually delivered to no more than two hepatic segments²
- The threshold dose of 190 Gy has been confirmed³ to maximize cytotoxicity and selective delivery minimizes risk of damage to surrounding parenchyma⁴

OBJECTIVES

- To report one center's long term outcomes of patients with HCC $\leq$ 5 cm, not amenable to transplantation, resection or RF ablation, who underwent radiation segmentectomy
- The authors hypothesized radiation segmentectomy could be considered potentially curative based on the same rationale as transplantation, resection and RF ablation

METHODS

- Retrospective, single center study looked at 70 patients with solitary HCC $\leq$ 5 cm, preserved liver function (Child Pugh A) and no vascular invasion or extrahepatic metastases who underwent radiation segmentectomy with Y-90 glass microspheres (target dose >190 Gy)
- Patients who had surgical resection or transplant after a radioembolization procedure were excluded
- A sub-analysis of patients with HCC $\leq$ 3 cm was also performed (cohort comparable to RF ablation)
- All patients underwent long-term imaging (contrast material-enhanced magnetic resonance [MR] imaging or computed tomography [CT]) and clinical follow-up (toxicity assessment at 1 and 3 months and response assessment in clinic 1 month post and subsequently at 3-month intervals)

KEY RESULTS


RADIOGRAPHIC RESPONSE	TIME TO PROGRESSION (TTP)	LOCAL TUMOR CONTROL	OVERALL SURVIVAL (OS)
63 out of 70 patients (90%) responded according to EASL criteria, of whom 41 (59%) showed complete response (CR).	Median TTP was 2.4 years, or 29 months.	72% of patients had no target lesion progression at 5 years. Local recurrence in complete responders occurred in 4 patients (9.8%).	Median OS (n=70) was 6.7 years, or 80 months. 1-, 3-, and 5-year survival probabilities were 98%, 66% and 57%, respectively. A sub-analysis of patients with tumor size $\leq$ 3 cm (n = 45) resulted in 1-, 3-, and 5-year survival probabilities of 100%, 82% and 75%, respectively.

* RF ablation shows similar reported survival outcomes to resection and transplantation for HCC $\leq$ to 3 cm

BCLC = Barcelona Clinic Liver Cancer; HCC= hepatocellular carcinoma; RF= radiofrequency; EASL= European Association for the Study of the Liver

RESULTS


Overall Survival for All Patients



OVERALL SURVIVAL RATE	Tumor Sizes		
	≤ 3cm (n=45)	> 3cm (n=25)	All ≤ 5cm (n=70)
1-year	100%	96%	98%
3-year	82%	46%	66%
5-year	75%	37%	57%



Contrast material-enhanced CT scan before Y-90 of an 87-year-old man with 4 cm hepatocellular carcinoma in right lobe.



Contrast-enhanced MR image at subsequent 9-year follow-up (now aged 96 years) shows complete necrosis.

CONCLUSION

- Liver transplantation, surgical resection and RF ablation are considered curative treatment options based on phase II randomized studies with limited data demonstrating improved survival
- Radiation segmentectomy data from the present publication showed similar outcomes as therapies considered curative in patients with unablatable BCLC stage 0 or A lesions ≤ 5cm with preserved liver function (Child Pugh A):
 - Provided local tumor control
 - Prolonged time to progression
 - Overall survival outcomes comparable to RF ablation, resection, and transplantation for patients with BCLC stage 0 or A HCC
- Additionally, radiation segmentectomy is an outpatient, minimally invasive intra-arterial therapy with a low toxicity profile that may be a convenient treatment option for patients
- Study Strengths: homogeneity of patient cohort, >10 years of follow-up, strict patient selection
- Study Limitations: retrospective and nonrandomized analysis, selection bias, comparisons to published literature versus an internal control group

1. Bruix J, Reig M, Sherman M. Evidence-based diagnosis, staging, and treatment of patients with hepatocellular carcinoma. *Gastroenterology* 2016;150(4):835-853. 2. Riaz A, Gates VL, Atassi B et al. Radiation segmentectomy: a novel approach to increase safety and efficacy of radioembolization. *Int J Radiat Oncol Biol Phys* 2011;79(1):163-171. 3. Vouche M, Habib A, Ward TJ et al. Unresectable solitary hepatocellular carcinoma not amenable to radiofrequency ablation: multicenter radiology-pathology correlation and survival of radiation segmentectomy. *Hepatology* 2014;60(1):192-201. 4. Vouche M, Lewandowski RJ, Atassi R et al. Radiation lobectomy: time-dependent analysis of future liver remnant volume in unresectable liver cancer as a bridge to resection. *J Hepatol* 2013;59(5):1029-1036

Boston Scientific
Advancing science for life™

Interventional Oncology and Embolization
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2026 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-789805-AB



TheraSphere™ Y-90 Glass Microspheres
Indications, Safety, and Warnings

<https://www.bostonscientific.com/therasphere-indications>