INSTRUCTIONS FOR USE

EMBLEM™ S-ICD

Automated Screening Tool (AST)

EMBLEM™ S-ICD

Automatisches screening tool (AST)

EMBLEM™ S-ICD

Outil de screening automatisé (AST)

EMBLEM™ S-ICD

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Automated Screening Tool (AST)

Instructions For Use

This manual contains instructions for use of Model 2889 EMBLEM S-ICD Automated screening tool (AST).

The following is a trademark of Boston Scientific Corporation or its affiliates: EMBLEM, LATITUDE, ZOOM.

Table of Contents

Description	5
Intended Use	5
Intended Audience	5
Indications for Use / Contraindications / Warnings and Precautions	5
Equipment Needed	6
Collecting Surface ECG	6
Navigating to the AST Application	6
Establishing the ECG Connection	6
Preparing skin prior to placing surface electrodes	6
Positioning the surface ECG electrodes	6
Verifying a Clean ECG	8
Observing the ECG on the screen	8
Using the AST Application	8
Information tab	8
Screening tab	9
Using the Screening Report to Determine an Acceptable Sense Vector	10
Report Description	10
Saving, printing and exporting the report	10
Determining an Acceptable Sense Vector	11
Symbols Table	12

Description

The Model 2889 EMBLEM S-ICD Automated Screening Tool (AST) is software used on the Model 3120 ZOOM LATITUDE Programmer to screen patients to assess them for implant of the EMBLEM Subcutaneous Implantable Cardioverter Defibrillation (S-ICD) system. These instructions describe how to use the AST.

Intended Use

The AST application is solely intended to screen patients for S-ICD implant and is not intended to be used for any cardiac diagnostic purpose. The AST is an alternative to the Model 4744 Patient Screening Tool. The two screening tools serve the same purpose, and may be used independently or together.

Intended Audience

This AST is intended for use by professionals trained or experienced in device implant and/or follow-up procedures.

Indications for Use / Contraindications / Warnings and Precautions

There are no indications, contraindications, warnings or precautions specific to the Model 2889 EMBLEM S-ICD Automated Screening Tool (AST). Refer to the appropriate S-ICD System Pulse Generator (PG) User Manual for a list of indications for use, allowed concomitant device settings, contraindications, and warnings and precautions for the S-ICD system.

Equipment Needed

- Model 3120 Programmer with AST software application
- External ECG cables compatible with the Model 3120 Programmer and associated accessories (e.g. electrodes and skin prep material)
- Measuring tool for measuring 14 cm (Model 4744 Patient Screening Tool may be used)
- External printer or pen drive connected to the Model 3120 Programmer

NOTE: The AST screening report cannot be printed using the internal Model 3120 Programmer strip chart printer.

Collecting Surface ECG

In order to perform the patient screening process, a surface equivalent of the subcutaneous sensing vectors used by the implanted S-ICD system must be obtained. The Model 3120 Programmer is used to collect the surface ECG and the AST software application evaluates the QRS complexes for each patient posture tested.

Navigating to the AST Application

1. Power on the Model 3120 Programmer.
2. Click **EMBLEM S-ICD Automated Screening Tool** from the home screen.

Establishing the ECG Connection

Preparing skin prior to placing surface electrodes


1. Remove hair.
2. Clean skin using a non-alcohol wipe and/or skin prep gel.

Positioning the surface ECG electrodes

It is important to collect the surface ECG in the location that represents the intended position of the implanted S-ICD System. When placing the S-ICD System in the typical implant location, the surface ECG electrodes should be positioned as follows (**Figure 1 Typical Placement of Surface ECG Electrodes for Patient Screening**). If a non-standard S-ICD System subcutaneous electrode or pulse generator placement is desired, the surface ECG electrode locations should be modified accordingly.

For typical implantation:

- **ECG Electrode LL** should be placed in a lateral location, at the 5th intercostal space along the midaxillary line to represent the intended location of the implanted pulse generator.
- **ECG Electrode LA** should be placed 1 cm left lateral of the xiphoid midline to represent the intended location of the proximal sensing electrode of the implanted subcutaneous electrode.
- **ECG Electrode RA** should be placed 14 cm superior to the ECG Electrode LA, to represent the intended position of the distal sensing tip of the implanted subcutaneous electrode.
 - Measure 14 cm using Model 4744 Patient Screening Tool or other measuring tool.
- **ECG Electrode RL** is recommended at the desired location to serve as the patient reference electrode.

1 - Alternate (LEAD I); 2 - Secondary (LEAD II); 3 - Primary (LEAD III)
RA - Right Arm; LA - Left Arm; LL - Left Leg

Figure 1 - Typical Placement of Surface ECG Electrodes for Patient Screening

The following best practices are recommended to obtain an ECG with a stable baseline, free of noise and motion artifacts:

- Ensure surface electrodes are new. Avoid using electrodes from packaging that was previously opened. Electrodes dry out quickly after the package is opened.
- Ensure enough electrode gel is on the electrode contact point.

Verifying a Clean ECG

Observing the ECG on the screen

1. Observe the ECG signal on the Model 3120 Programmer screen. Each ECG signal should have a **stable baseline** and be **free of noise** and **motion artifacts**.
2. Enter the full screen view to zoom in on the signal by clicking the magnifying glass  button in the upper right-hand corner. Gain may be adjusted to enhance view. Manual gain adjustments will not affect AST results or printed reports.

Consider the following troubleshooting tips if noise or a wandering baseline is observed:

- Replace ECG cables used to collect surface ECGs.
- In particularly difficult cases, re-prepare the skin with some light abrasion at the electrode sites (using a gauze pad, for example).
- Ask patient to inhale then exhale and hold for 10 seconds to manage respiratory artifact.

Using the AST Application

The AST application is layed out in tabs. Complete each tab to run the screening tests.

Information tab

Complete the patient information (**Patient Name or ID** and **Date of Birth** are required fields).

NOTE: This information can be edited at any point in the screening process.

Screening tab

Run the ECG screening tests on the **Screening** tab.

1. Select the **Sternal Lead Position** by clicking the pull-down menu (defaults to the most common placement – **Left Sternal Margin**; other options are **Right Sternal Margin** and **Medial**). The **Sternal Lead Position** represents the planned location of the vertical portion of the S-ICD implanted electrode.
2. Run **Supine** posture (required) - Put the patient into a supine posture and click the corresponding Run  button to run the ECG screening test. Results are displayed when the test is complete.
3. Run **Standing/Sitting** posture (required) - Have the patient assume either a standing or sitting posture and click the corresponding Run  button to run the ECG screening test.
4. Run **Other** postures or conditions (optional) - Click the Run  button associated with an **Other** column to add patient-specific postures or conditions, as needed. Click the text field to add posture or condition descriptions.
5. Click the magnifying glass  button below each column to access the **Snapshot** screen and review the details of each test. Each ECG signal should have a **stable baseline** and be **free of noise** and **motion artifacts**.

*NOTE: It is possible to repeat a test. A clinician might want to repeat a test if the test result is **FAIL** but a review of the ECGs reveals a disruption in the signal from a noise source, patient movement, or intermittent issue with the connection to the surface leads. To repeat a test, simply click the Run  button for the desired posture again, then observe the results.*

NOTE: If testing more than 6 postures or conditions, save the report and start a new session. The user interface only provides 6 columns for postures or conditions for a single session.

*NOTE: The **Advanced View** option on the **Snapshot** screen adds the filtered surface ECG data used by the AST. Your local Boston Scientific representative or Technical Services may assist in using this view to further understand the AST results. Printing is disabled in this view.*

Using the Screening Report to Determine an Acceptable Sense Vector

Report Description

The EMBLEM S-ICD Automated Screening Report PDF summarizes screening results and instructs users to check QRS complex morphology across postures to identify acceptable leads. The ECG data used by the AST is provided for each lead/posture combination tested.

NOTE: The gain of the ECG signal for each vector is automatically adjusted prior to printing. This allows the printed report ECGs to be assessed with the Model 4744 Patient Screening Tool, if desired.

Saving, printing and exporting the report

1. To save the report, click the **Save Report** button on the lower left-hand corner of the **Screening** tab - a PDF of the report is saved to the Model 3120 Programmer hard drive.

NOTE: It may take a couple of minutes to generate the PDF report, depending on the number of postures or conditions run.

NOTE: A maximum of 100 reports may be saved to the hard drive. A capacity warning appears when the maximum has been reached. The user has the option to proceed with the oldest session being deleted or abort the session.

2. To print the report, click the **Print Report** button to send the report to a connected external printer.
3. To export the data to a pen drive:
 - Click **End Screening Session** to return to the programmer home screen.
 - Insert pen drive into the USB port of the Model 3120 programmer.
 - Click **Patient Data Management**, then select the **Export Screening Data** tab.
 - Click **Export** or **Export with Password Protect** and follow the prompts (all saved reports are exported).

Determining an Acceptable Sense Vector

The physician must review the report and determine whether the patient is suitable for implant of the S-ICD System. Confirm that the following conditions are met:


- At least one common ECG lead must be deemed acceptable for all tested postures. At a minimum, **Supine** and **Standing/Sitting** postures must be tested.
- The morphology of the intrinsic/paced QRS complex is stable across postures (similar positive/negative peak amplitudes and QRS widths). No significant change to the QRS complex is noted as a result of postural changes. For notched signals, ensure that the location of the larger peak is consistent in relation to the smaller peak.

NOTE: Special circumstances may present in which the physician elects to proceed with the implantation of the S-ICD System despite failing the screening process. In this case, careful attention should be applied to the device setup process of the S-ICD System as the risk of poor sensing and/or inappropriate shock is increased.

Symbols Table

The following symbols may be used on packaging and labeling.

Symbol	Meaning
	Authorized Representative in the European Community
	Manufacturer
	CE mark of conformity with the identification of the notified body authorizing use of the mark
	Reference Number
	Lot Number
	Consult instructions for use on this website: www.bostonscientific-elabeling.com
	Date of manufacture
	Australian Sponsor Address

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or its affiliates. All rights reserved.
Co. Ltd.
© 2016 Boston Scientific Corporation