INSTRUCTIONS FOR USE / BRUGSANVISNING / GEBRAUCHSANWEISUNG / KASUTUSJUHISED / KÄYTTÖOHJEET / INSTRUCTIONS D'UTILISATION / ISTRUZIONI PER L'USO / 한국어 / NAUDOJIMO INSTRUKCIJOS / LIETOŠANAS INSTRUKCIJA / GEBRUIKSAANWIJZING / BRUKSANVISNING / BRUKSANVISNING

VEIN PICK / VENELØFTER / VENENHAKEN / VEENITÖSTJA / LASKIMOASETIN / CHAUSSE-VEINE / DILATA-VENA / 정맥 후비개 / VENOS PARINKIMAS / VĖNAS IZVĒLE / VENEHAAKJE / VENEHOLDER / VENHAKE

REF 6541

Contains: / Indeholder: / Inhalt: / Sisaldab: / Sisältö: / Contient: / Contiene: / 내용물: / Komplektą sudaro: / Satur: / Inhoud: / Inneholder: / Innehåller: (1) Vein Pick / Veneløfter / Venenhaken / Veenitöstja / Laskimoasetin / Chausse-veine / Dilata-vena / 정맥 후비개 / Venos parinkimas / Vėnas izvēle / Venehaakje / Veneholder / Venhake


Figure 1. Using vein pick / Figur 1. Brug af veneløfteren / Abbildung 1. Verwendung des Venenhakens / Joonis 1. Veenitöstja kasutamine / Kuva 1. Laskimoasettimen käyttö / Figure 1. Utilisation du chausse-veine / Figura 1. Utilizzo del dilata-vena / 그림 1. 정맥 후비개 사용 / 1 pav. Venos pakeliklio naudojimas / 1. attēls. Vēnas pletēja izmantošana / Figuur 1. Gebruik van venehaakje / Figur 1. Bruk av veneholder / Figur 1. Användning av venhake

ENGLISH

This literature is intended for use by professionals trained or experienced in device implant and/or follow-up procedures.

**WARNING:** For single patient/single procedure use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

**CAUTION:** The accessory is sterilized with ethylene oxide (EO) gas before final packaging. When the accessory is received, it is sterile provided the container is intact. If the container is wet, punctured, opened, or otherwise damaged, return the accessory to Boston Scientific.

**CAUTION:** Do not use vein pick to puncture vein or dissect tissue during cutdown procedures.

**CAUTION:** Be sure that the vein pick does not puncture the insulation of the lead. This could prevent proper lead function.

STERILE

Sterilized using ethylene oxide.

INTENDED USE

Use to open vein incision before placing lead tip into vein.

The overall clinical benefit for this accessory is derived from the device system with which it is intended to be used.

DESCRIPTION

The vein pick is made from Hoechst Celanese Acetal.

INSTRUCTIONS FOR USE

Open the vein and place point of vein pick into lumen of vein. With point of vein pick facing in direction of desired lead passage, gently raise and lift the pick from where the lead will enter the vein. Pass lead under vein pick and into vein. (Figure 1)

For additional information on how and when to use this accessory (including any warnings, precautions, and other safety information related to its use for the specific procedure being performed), see the literature accompanying the lead to be implanted.

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE** For customers in the European Union, use the accessory name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website:

https://ec.europa.eu/tools/eudamed

DISPOSAL

**CAUTION:** Dispose of all items used during a procedure using standard biohazard handling techniques since they are considered biohazardous.

All items used during a procedure, such as accessories, consumables, and packaging, may be contaminated with infectious substances. Consider the following to minimize the risk of infection, microbial hazards, or physical harm:

- Biohazardous waste should be disposed in a biohazard container that is labeled with the biological hazard symbol and taken to a designated facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and/or local government policy.
- Biohazardous waste should be treated with an appropriate thermal or chemical process.

**NOTE:** *Untreated biohazardous substances should not be disposed of in the municipal waste system.*

WARRANTY DISCLAIMER

Boston Scientific disclaims all expressed and implied warranties for this product, including without limitation any implied warranties of merchantability or fitness for a particular purpose. Buyer assumes all risk of loss or damages arising from use of this product.

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© 2021 Boston Scientific Corporation or its affiliates. All rights reserved.