

Product Information for Patients


RELIANCE 4-FRONT™

Cardiac Lead

REF 0652, 0653, 0657, 0658, 0672, 0673, 0675, 0676, 0692, 0693, 0695, 0696
0636, 0650, 0651, 0654, 0655, 0662, 0663, 0665, 0682, 0683, 0685, 0686

Device Information

A lead is an insulated wire that is passed into your heart through a vein. One end is implanted in the heart wall and the other end is connected to a pulse generator, which contains electronics and a battery. The lead senses your heartbeat and delivers pacing pulses and/or shocks from the pulse generator to the heart.



Information on Safe Use

Before implant, tell your healthcare provider if you are sensitive or allergic to any of the materials or substances listed in the "Patient Contacting Materials and Substances" on page 2.

Carry your Implant Card at all times and present it before entering protected environments, such as for an MRI scan.

Tell your healthcare providers, such as your doctor, dentist, or technician, that you have an implanted device, and show them your Implant Card(s) so that they can take any necessary precautions.

Information about your implanted device may be found using the websites shown on the back of your Implant Card. Device information available on the websites varies by region.

Contact your healthcare provider that follows your device:

- If you have an adverse event (in other words, an undesirable side effect) or unusual/unexpected symptoms. These include new symptoms or symptoms like the ones felt before implant of your device.
- After any medical procedure and/or surgery in order to have your device checked.

Warnings and/or Precautions Related to Reciprocal Interference

Some procedures use equipment that can potentially interfere with the operation of your device and vice versa (i.e., reciprocal interference). Talk with your healthcare provider that follows your device before undergoing these types of procedures. Examples of these types of procedures are listed below.

- **Magnetic resonance imaging (MRI):** This is a diagnostic test that uses a strong electromagnetic field. Some devices have been evaluated to allow the patient to undergo MRI scans under specific conditions. Talk to your healthcare provider that follows your device about the eligibility of your device. If your device is not eligible to be scanned (i.e., not part of an ImageReady™ system), or if the required conditions are not met, MRI scans can severely damage your device and should not be performed. Your healthcare provider that follows your device must always confirm that both you and your device are eligible and ready for an MRI scan in order for you to undergo this procedure. Hospitals/healthcare centers keep MRI equipment in rooms marked with signs that indicate magnets are inside. Do not go inside these rooms unless your healthcare provider that follows your device has confirmed that your device is eligible and you meet the requirements for an MRI scan.
For information about MRI scanning, contact your healthcare provider that follows your device and/or visit the following website:
www.bostonscientific.com/imageready
- **Diathermy:** This uses an electrical field to apply heat to tissues in the body and could damage your device or injure you. Diathermy should not be performed.

Expected Lifetime and Follow-up

Based on results of testing, the expected lifetime of this lead is typically a minimum of 10 years. Your healthcare provider will monitor the lead and will decide if and when it may need to be replaced.

Discuss a follow-up plan with your healthcare provider that follows your device, including the frequency and type of follow up.

Report any serious incident that involves this lead (i.e., an event requiring medical attention) to your healthcare provider that follows your device, as well as to Boston Scientific using the information on the last page and to the relevant local regulatory authority in your country.

Patient Contacting Materials and Substances

The lead contains the following materials and substances that come in contact with the body:

Material/Substance	Percentage (%) of patient-contacting surface area ^a
Models 0652, 0653, 0657, 0658, 0672, 0673, 0675, 0676, 0692, 0693, 0695, 0696	
Silicone	70% – 90%
Platinum, ePTFE (expanded polytetrafluoroethylene) ^b	5% – 25%
IROX™ (iridium oxide), polyurethane, TiO ₂ (titanium dioxide), BaSO ₄ (barium sulfate), dexamethasone acetate (drug) ^c , PEEK (polyetheretherketone), titanium (grade 2), platinum/iridium	Additive, trace amount, and/or <5%











- Typical total surface area of lead = 45-60 cm² (depends on model).
- GORE™ models (0657, 0658, 0692, 0693, 0695, 0696) only: The platinum shocking coils are covered with an ePTFE coating. This coating is the primary patient-contacting surface. Both the platinum and ePTFE are included in the surface area percentage. GORE is a trademark of W.L. Gore and Associates.
- 0.96 mg for models 0652, 0653, 0657, 0658, 0672, 0673, 0675, 0676, 0692, 0693, 0695, 0696

Material/Substance	Percentage (%) of patient-contacting surface area^a
Models 0636, 0650, 0651, 0654, 0655, 0662, 0663, 0665, 0682, 0683, 0685, 0686	
Silicone	75% – 95%
Platinum, ePTFE (expanded polytetrafluoroethylene) ^b	5% – 25%
IROX™ (iridium oxide), polyurethane, TiO ₂ (titanium dioxide), BaSO ₄ (barium sulfate), dexamethasone acetate (drug) ^c , titanium (grade 2), platinum/iridium, PEG (polyethylene glycol)	Additive, trace amount, and/or <5%

- a. Typical total surface area of lead = 45-60 cm² (depends on model).
- b. GORE™ models (0654, 0655, 0682, 0683, 0685, 0686) only: The platinum shocking coils are covered with an ePTFE coating. This coating is the primary patient-contacting surface. Both the platinum and ePTFE are included in the surface area percentage.
- c. 0.97 mg for models 0636, 0650, 0651, 0654, 0655, 0662, 0663, 0665, 0682, 0683, 0685, 0686

Symbols Table

The following symbols are used for patient information:

Symbol and Definition	Patient Information
 Date	Date of device implant
 Health Care Center or Doctor	Name and contact information of healthcare institution where device implanted or clinician who performed implant
 Person Identification	Name of patient receiving implanted device
 Reference Number	Model number of implanted device
 Manufacturer	Company that made implanted device
 Authorized Representative in the European Community	European contact information for company that made implanted device
 Serial Number	Serial number of implanted device
 Use by date	Date the device must be implanted by
 Date of Manufacture	Date the device was made
 Unique Device Identifier	Barcode with information about the device

Остаряла версия. Да не се използва.
Zastaralá verze. Nepoužívat.
Forældet version. Må ikke anvendes.
Version überholt. Nicht kasutage.
Aegunud versioon. Mην την χρησιμοποιείτε.
Παλιά έκδοση. Μην την χρησιμοποιείτε.
Outdated version. Do not use.
Version obsolete. Ne utilizar.
Zastarjela verzija. Nemojte upotrebljavati.
Version périmée. Ne pas utiliser.
Zastarjela verzija. Nemojte upotrebljavati.
Úreлт útgáfa. Notið ekki.
Versione obsoleta. Non utilizzare.
Pasenusi versija. Nenaudokite.
Elavult verzió. Ne használja!
Utdatert versjon. Skal ikke brukes.
Wersja przeterminowana. Nie używać.
Versão obsoleta. Não utilize.
Versiune expirată. A nu se utiliza.
Zastarana verzija. Nepoužívať.
Vanhentunut versio. Älä käytä.
Föråldrad version. Använd ej.
Güncel olmayan sürüm. Kullanmayın.

EC REP

Guidant Europe NV/SA
Boston Scientific
Green Square,
Lambroekstraat 5D
1831 Diegem, Belgium



Boston Scientific Corporation
4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA

1.800.CARDIAC (227.3422)
Worldwide: +1.651.582.4000

www.bostonscientific.com
www.bostonscientific.com/patientlabeling

CE 2797

© 2021 Boston Scientific Corporation or its affiliates.
All rights reserved.



51070663-001 en Europe 2021-12