Agile™ Esophageal

PARTIALLY COVERED

Stent System

Directions for Use	2
Instrucciones de uso	8
Mode d'emploi	14
Gebrauchsanweisung	20
Istruzioni per l'uso	27
Instructies voor gebruik	33
Instruções de Utilização	39

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION	3
Contents	3
Figure 1. Agile™ Esophageal Partially Covered Stent.....	3
User Information	3
MRI Safety Information.....	3
INTENDED USE/INDICATIONS FOR USE	4
CONTRAINDICATIONS	4
WARNINGS AND CAUTIONS	4
ADVERSE EVENTS	4
HOW SUPPLIED	5
Handling and Storage	5
INSTRUCTIONS FOR USE	5
Required Equipment	5
Initial preparation of delivery system	5
OPERATIONAL INSTRUCTIONS	5
Start of initial stent placement procedure.	5
1. Locate the Stricture	5
2. Examine the Stricture (Endoscopically and/or Fluoroscopically)	5
3. Choose the Stent Size	5
4. Insert Guidewire and Place Through Stricture	5
5. Advance the Delivery System over the Guidewire and Position Stent	6
Figure 2. Delivery System and Radiopaque (RO) Markers	6
6. Deploy Stent	6
Figure 3. Delivery System, Visual Markers and Handles	6
Reconstriction Technique	6
7. Assess Deployed Stent Position and Remove the Delivery System	7
8. Remove Gastroscope	7
POST-PROCEDURE	7
Implantable Device Patient Information	7
Implant Card Instructions	7
WARRANTY	7

Agile™ Esophageal

PARTIALLY COVERED

Stent System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

Contents

- (1) Agile Esophageal Partially Covered Stent System

The Agile Esophageal Partially Covered Stent System is comprised of a metallic implantable stent pre-loaded inside of a flexible delivery system. The system is compatible with gastroscopes with a minimum 3.7 mm working channel. The stent is made from braided Nitinol wires which form a self-expanding, radiopaque (RO) cylindrical mesh. The stent has flares at each end to aid in minimizing migration after the stent is placed in the esophagus. The flares are a wider diameter than the stent body. The wire ends are looped at the end of the stent. The proximal and distal stent ends each have a continuous suture threaded around their circumferences. The suture is intended to aid in removal or repositioning during the initial stent placement procedure, to be used in the event of incorrect placement.

The stent is partially covered with a silicone polymer to restrict tumor ingrowth through the wire mesh and to occlude concurrent esophageal fistulas. (Figure 1)

The delivery system is a co-axial tube design. The exterior tube is used to constrain the stent before deployment and to reconstrain the stent after partial deployment. The exterior tube has clear distal section so that the constrained stent is visible. A yellow transition zone on the inner tube of the delivery system is visible between the stent and the blue outer sheath. The system has RO and visual markers to aid in accurate stent placement.

There are three visual markers on the delivery system handle to aid in stent placement. (Figure 3) The most distal visual marker indicates the stent is fully

constrained on the delivery system, (Figure 3 #6) one visual marker indicates that the stent is 50% deployed (mid-point) (Figure 3 #7) and the most proximal visual marker indicates the point at which stent reconstraint is no longer possible. (Figure 3 #8)

There are five radiopaque (RO) markers to aid in the deployment of the stent while using fluoroscopy (Figure 2). There are two RO markers on the inner tube of the delivery system identifying the ends of the constrained stent (Figure 2, #1 and #4). Between these RO markers are two additional RO markers. One RO marker indicates the mid-point of the stent (Figure 2 #2). One RO marker on the interior tube indicates the point at which stent reconstraint is no longer possible (Figure 2 #3). The fifth RO marker at the leading end of the exterior tube indicates how far the stent has been deployed (Figure 2, #5).

The system has a single central lumen to accommodate a 0.035 in (0.89 mm) guidewire.


Figure 1. Agile Esophageal Partially Covered Stent

Caution: Read the entire Instructions for Use thoroughly before using the Agile Esophageal Partially Covered Stent System.

User Information

The Agile Esophageal Partially Covered Stent System should only be used by or under the supervision of physicians thoroughly trained in esophageal prosthesis placement. A thorough understanding of the technical principles, clinical applications, and risks associated with this procedure is necessary before using this device.

MRI Safety Information

A person with the Agile Esophageal Partially Covered Stent may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Name	Agile Esophageal Partially Covered Stent
Static magnetic Field Strength (B₀)	1.5T or 3.0T

Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Cylindrical Whole-body Coil Cylindrical Head Coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
Scan Duration	Under scan conditions defined above, patient can be scanned for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR Image Artifact	Image artifact can be produced

INTENDED USE/INDICATIONS FOR USE

The Agile™ Esophageal Partially Covered Stent System is intended for maintaining esophageal luminal patency in esophageal strictures caused by intrinsic and/or extrinsic malignant tumors, and occlusion of concurrent esophageal fistulas.

CONTRAINDICATIONS

The Agile Esophageal Partially Covered Stent System is contraindicated for:

- Placement in esophageal strictures caused by benign tumors, as the long-term effects of the stent in the esophagus are unknown.
- Placement in strictures that cannot be dilated enough to pass the gastroscopie or the delivery system.
- Placement of the proximal end of stent within 2 cm of the cricopharyngeal muscle.
- Placement in an esophago-jejunostomy (following gastrectomy), as peristalsis and altered anatomy may displace stent.
- Placement in necrotic chronically bleeding tumors, if bleeding is active at the time of placement.
- Placement in polypoid lesions.
- Those patients for whom endoscopic techniques are contraindicated.
- Any use other than those specifically outlined under indications for use.
- Placement in patients who have an underlying bleeding diathesis.

WARNINGS AND CAUTIONS

- The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.
- As perforation is a known risk, the stent should be used with caution and only after careful consideration in patients who are:
 - undergoing radiation therapy and/or chemotherapy
 - in advanced stages of cancer
- The Agile Esophageal Partially Covered Stent System should be used with caution and only after careful consideration in patients with:
 - Strictures exceeding 12 cm in length
 - Significant preexisting pulmonary or cardiac disease

- This device contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

Warning: Stent is considered to be a permanent device. Once stent placement is permanently achieved, stent removal or repositioning is not recommended.

Warning: Visually inspect the system for any signs of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this warning may result in patient injury.

ADVERSE EVENTS

The potential adverse effects associated with esophageal stent placement may include:

- Bleeding
- Perforation
- Pain
- Aspiration
- Stent migration
- Tumor ingrowth through uncovered portion of stent
- Tumor overgrowth around stent ends
- Foreign body sensation
- Food bolus impaction
- Reflux
- Esophagitis
- Edema
- Ulceration
- Fever
- Infection
- Sepsis
- Septicemia
- Recurrent dysphagia
- Fistula formation
- Tracheal compression/obstruction (or acute airway compression)
- Hematemesis
- Death (other than that due to normal disease progression)
- Stent fracture

Possible Post Stent Complications

- Sensitivity to the metal component of the stent
- Mediastinitis
- Aspiration
- Intestinal obstruction (secondary to stent migration)
- Granulation tissue around stent ends
- Aorto and arterio-esophageal fistula
- Erosion or perforation of stent into adjacent vascular strictures

HOW SUPPLIED

The device is supplied sterile and intended for single use only. The packaging and device should be inspected prior to use. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place. See product label for expiration date.

INSTRUCTIONS FOR USE

Required Equipment

- Gastroscope with a minimum 3.7 mm working channel
- Rat Tooth forceps
- Fluoroscopic capability for pre-stent placement and stent placement confirmation
- 0.035 in (0.89 mm), 450 cm guidewire
- Agile™ Esophageal Partially Covered Stent System containing a stent of the appropriate length and diameter

Pre-Procedure

Radiography of the esophagus performed no more than 10 days before the procedure should be available.

Prepare for the procedure as you would an upper endoscopy. Administer a mild sedative, if necessary and a topical anesthesia to the throat to be repeated during the procedure as patient comfort warrants.

Initial preparation of delivery system

- Carefully remove the delivery system from the protective packaging.
- Visually inspect the device for damage or defects.

OPERATIONAL INSTRUCTIONS

Start of initial stent placement procedure.

1. Locate the Stricture

Intubate the patient using a standard gastroscopie per standard technique. Access the stricture location upon direct visualization. Fluoroscopy can also be used to locate the stricture with the aid of a contrast medium.

2. Examine the Stricture (Endoscopically and/or Fluoroscopically)

A. Examine Stricture Endoscopically

Endoscopically examine both the proximal and distal segments of the stricture. Using the external ruler on the gastroscopie, measure the distance between the distal margin of the stricture to the patient's incisors. Withdraw the gastroscopie to the proximal margin of the stricture and measure the distance to the patient's incisors. The stricture length is calculated as the difference between those two distances. To minimize the potential of stent migration, dilate the stricture ONLY if passage of the gastroscopie or the delivery system through the stricture lumen is not possible.

Caution: In some patients, tumor encroachment will make dilation of the stricture challenging. Physicians should use judgment based on experience in dilating esophageal strictures. Perforation or bleeding of an esophageal tumor is a risk during a tumor dilation procedure.

Warning: Placement of the Agile Esophageal Partially Covered Stent should not be attempted in patients with esophageal strictures that cannot be dilated wide enough for passage of the gastroscopie or delivery system.

B. Examine the Stricture Fluoroscopically

The stricture may also be examined fluoroscopically. Leaving the gastroscopie in place, observe both the proximal and distal margins of the tumor fluoroscopically. Mark the locations with either radiopaque markers or use anatomical landmarks such as ribs or vertebrae. It is recommended to re-measure the stricture length by measuring the distance between the radiopaque markers.

Warning: Physicians should use judgement based on experience when dilating esophageal strictures. Perforation or bleeding of the esophageal tumor is a risk during any dilation of the tumor.

3. Choose the Stent Size

The size of the stricture must be accurately calculated to ensure the ideal stent size is used. The Agile Esophageal Partially Covered Stent should bridge the tumor and/or the fistula and should extend >1 cm above and below the stricture or fistula. For stent use with a fistula, it is critical to ensure that the covered portion of the stent completely covers the fistula to avoid leakage and facilitate healing. If stent length choice is questionable, always use the longer stent. A second stent of the same diameter may be placed if the first stent does not cover the entire stricture length. The second stent should be placed to ensure complete tumor coverage and a smooth transition between the stents. It is recommended that the proximal stent be placed first followed by the distal stent in order to maximize luminal diameter of the interlocked stents. Care should be used when passing the delivery system through the first stent. As the stent is deployed, the stent will foreshorten. Foreshortening is defined as the percentage decrease from constrained stent length within the delivery system to deployed stent length. Benchtop testing has demonstrated the Agile Esophageal Stent will foreshorten no more than 50% from its constrained length on the delivery system. However, actual foreshortening depends on the anatomy of the lumen and stricture.

Warning: Passing the scope through a just deployed stent is not recommended and could cause the stent to dislodge.

Warning: Do not use in combination with stents from other manufacturers.

4. Insert Guidewire and Place Through Stricture

Pass a guidewire through the working channel of the gastroscopie and then through the stricture and into the stomach. A floppy tip guidewire is recommended in order to reduce potential trauma from the tip of the wire. Endoscopic and fluoroscopic placement of the guidewire is also recommended to assure proper passage through the stricture and proper placement in the stomach. Maintain guidewire position throughout procedure.

Caution: A 0.035 in (0.89 mm) guidewire with a floppy tip is recommended to facilitate passage through tortuous anatomy. Dreamwire™ 0.035 in (0.89 mm) Standard M00556141, Dreamwire 0.035 in (0.89 mm) Stiff M00556161, Jagwire™ 0.035 in (0.89 mm) Standard M00556581 or Jagwire 0.035 in (0.89 mm) Stiff M00556601 is recommended.

5. Advance the Delivery System over the Guidewire and Position Stent


Figure 2. Delivery System and Radiopaque (RO) Markers

There are three visual markers on the delivery system handle to aid in stent placement while using endoscopic visualization (Figure 3). The most distal visual marker indicates the stent is fully constrained on the delivery system (Figure 3 #6), one visual marker indicates that the stent is 50% deployed (mid-point) (Figure 3 #7) and the most proximal visual marker indicates the point at which stent reconstraint is no longer possible (Figure 3 #8).

There are five radiopaque (RO) markers to aid in the deployment of the stent while using fluoroscopy (Figure 2). There are two RO markers on the inner tube of the delivery system identifying the ends of the constrained stent (Figure 2, #1 and #4). Between these RO markers are two additional RO markers. One RO marker indicates the mid-point of the stent (Figure 2 #2). One RO marker on the interior tube indicates the point at which stent reconstraint is no longer possible (Figure 2 #3). The fifth RO marker at the leading end of the exterior tube indicates how far the stent has been deployed (Figure 2, #5).

The Agile™ Esophageal Stent System is passed over a guidewire and through the working channel of the gastroscope.

Under direct endoscopic visualization, position the stent keeping direct endoscopic visualization of the yellow transition zone.

If using fluoroscopic guidance, position the stent so marker 2 (figure 2) is in the center of the tumor or fistula. This ensures that the stent will properly bridge the tumor.

If it is not necessary to cross the Lower Esophageal Sphincter (LES), the distal end of the stent should remain above the LES in order to leave the LES functional and reduce gastric reflux. The stent can cross the LES if necessary due to tumor involvement and stricture.

6. Deploy Stent


Figure 3. Delivery System, Visual Markers and Handles

Caution: Do not twist the delivery system or use a boring motion during stent deployment as this may affect stent positioning and ultimately, the stent function.

Begin stent deployment by holding the distal handle (farthest from the operator, Figure 3 #9) of the delivery system with one hand, and using the other hand, grasp the proximal handle (closest to the operator, Figure 3 #10) and hold this handle stationary. Between the handles is a segment with visual markers. These visual markers are intended to aid in stent deployment. Prior to deploying the stent, you can visualize a marker at the fully constrained/undeployed position (Figure 3 #5).

In order to deploy the stent, hold the distal handle (farthest from the operator, (Figure 3 #9) with one hand and the proximal handle (handle closest to the operator, Figure 3 #10) with the other hand. To deploy the stent, slowly pull the distal handle toward the proximal handle while holding the proximal handle stationary. Monitor the stent release fluoroscopically and/or endoscopically, keeping the markers on the delivery system between the identified stricture margins. If necessary, it is possible to stop deployment and adjust the stent position proximally without reconstraining the stent prior to passing the reconstraint marker (Figure 3 #8). If satisfied with stent placement, proceed with full deployment.

Reconstraint Technique

- The stent can be reconstrained at any point up to the reconstraint markers (radiopaque Figure 2, #3 and visual Figure 3, #8).

Note: Once the visual reconstraint marker cannot be seen during deployment the stent cannot be reconstrained.

- Reconstraint is done by reversing the direction of deployment, by holding the proximal handle (closest to the operator) steady while pushing the distal handle (farthest from the operator) away.
- The stent has been designed to be reconstrained no more than two times.
- Prior to full deployment, if repositioning is desired, the stent can be pulled proximally by slowly pulling back on the delivery system. The ability to pull proximally will be limited by the amount of stent deployed and the tightness of the stricture. Full reconstraint, when possible, is always preferred and recommended over pulling the device proximally.

Note: The stent is fully constrained if visual marker (Figure 3, #6) is fully visible.

Caution: Pulling proximally when partially deployed could further deploy the stent if there is resistance on the stent.

Caution: Do not push the delivery system forward once deployment has begun. The delivery system can be pulled proximally if necessary. The ability to pull proximally will be limited by the amount of stent deployed and the tightness of the stricture.

Warning: Stent is considered to be a permanent device. Once stent placement is permanently achieved, stent removal or repositioning is not recommended.

If the positioning of an Agile Esophageal Partially Covered Stent is not correct and one of the following has occurred, continue to fully deploy the stent:

A. The stent has already been deployed past the reconstraint limit OR

B. The stent has already been reconstrained twice

Then in either case, using rat tooth forceps, grasp the suture on the proximal or distal end of the stent (Figure 1). Gently pull the stent back with the scope to remove the stent during the initial stent placement procedure.

Caution: Fully grasp around the suture when repositioning or removing the Agile Esophageal Partially Covered Stent.

Caution: Do not remove the stent by pulling the stent through the scope. After grasping the stent remove both the gastroscope and stent together.

7. Assess Deployed Stent Position and Remove the Delivery System

Following stent deployment, view the stent endoscopically and/or fluoroscopically to confirm stent expansion as tumor impingement may prevent the stent from achieving its maximum diameter immediately.

Carefully remove the delivery system and the guidewire.

Note: A stent may require 24 hours to expand fully.

Warning: Once stent is in desired location, passing the scope through a just deployed stent is not recommended and could cause the stent to dislodge.

Caution: Never use a rigid-type dilator for post stent placement dilation because the axial force may dislodge the stent. Physicians should use judgement based on experience when dilating.

Caution: An attempt to remove the delivery system and guidewire prior to stent expansion or when a stent is partially deployed may dislodge the stent.

If excessive resistance is felt during delivery system removal due to the stent being partially deployed, then proceed with the following steps:

- A. Wait 3-5 minutes to allow further stent expansion
- B. If the proximal end of the stent is pursed onto the delivery system, use the endoscope to manipulate the delivery system in a circular motion to open the proximal end of the stent.
- C. Re-sheath the exterior tube of the delivery system by pushing the distal handle (Figure 3 #9) away from the operator. Slowly withdraw the delivery system and guidewire.
- D. If removal is still not possible, use a balloon dilation catheter to dilate the stent. It should not be necessary for the balloon diameter/size to be equal to the stent diameter. Judgement should be used when selecting the balloon size. Carefully position the balloon catheter within the stent. Inflate the balloon to its recommended pressure.
- E. Deflate the balloon catheter and withdraw into the gastroscope. Slowly withdraw the delivery system and guidewire.

8. Remove Gastroscope

Withdraw gastroscope from the patient.

This completes the initial stent placement procedure.

Stent placement is considered permanent upon completion of initial stent placement procedure.

POST-PROCEDURE

Patients should have P-A (posteroanterior) and lateral chest films taken as a permanent record of the stent position. Observe the patient for development of any complications of endoscopy, esophageal dilation and stent placement. Vital signs should be monitored and clear liquids given in an upright position during the first 24 hours post-stent placement. Patients being treated for fistula should have no fluids or foods orally until successful sealing of the fistula has been confirmed. Post 24 hours, the patient should be instructed to eat only in an upright sitting position, to chew food thoroughly, to avoid certain foods (such as meats, raw vegetables, and breads), and to drink fluids during and after meals. Patients with stents placed in the distal esophagus or across the LES should be instructed to elevate the head of the bed, and should be prescribed acid suppression therapy to minimize gastric reflux into the stent. Subsequent follow-up at 1 week and at 3

month intervals, thereafter, or for symptomatic dysphagia may be performed to verify patency and placement.

Note: Recurrence or worsening of dysphagia may occur after stent placement due to tumor in-growth or overgrowth over time, severe hyperplasia reaction or stent migration. Repeat endoscopy may be required.

Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

Implant Card Instructions

- Apply the peel-off label from the product onto the supplied patient implant card.
- Fill in the implant date, patient name, healthcare institution and/or physician information.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or sterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2022 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd