Boston
Scientific

RX Locking Device

OLYMPUS® / FUJINON® COMPATIBLE

**Locking Device For Use with
RX Biliary System™**

Directions for Use 2

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Mode d'emploi 11

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51088990-01

2020-03

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
INTENDED USE/INDICATIONS FOR USE	3
Clinical Benefit Statement.....	3
CONTRAINDICATIONS.....	3
OTHER PRECAUTIONS.....	3
HOW SUPPLIED.....	3
Handling and Storage.....	3
OPERATIONAL INSTRUCTIONS	3
Figure 1.....	4
Figure 2.....	4
Figure 3.....	5
Figure 4	5
DEVICE REMOVAL	5
Disposal	5
Post-Procedure	5
WARRANTY	5

RX Locking Device

OLYMPUS® / FUJINON® COMPATIBLE

Locking Device For Use with RX Biliary System™

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

These devices are accessories used with RX Biliary System. Read appropriate Directions for Use.

INTENDED USE/INDICATIONS FOR USE

The RX Locking Device is intended to facilitate placement of the guidewire and locking it into place during ERCP.

Clinical Benefit Statement

To aid in placement of the guidewire and locking it into place during ERCP procedures.

CONTRAINDICATIONS

None known.

OTHER PRECAUTIONS

- The RX Locking Device should only be used by or under the supervision of physicians trained in Endoscopic Retrograde Cholangiopancreatography (ERCP) or Endoscopic Sphincterotomy (ES). A thorough understanding of the technical principles, clinical applications and risks associated with ERCP/ES is necessary before using this device.
- The RX Locking Device is compatible with 0.035 in (0.90 mm) and 0.025 in (0.63 mm) guidewires.
- Any use of this device, other than those indicated in these instructions, is not recommended.

HOW SUPPLIED

RX Locking device is provided non-sterile. Do not use if package is damaged or unintentionally opened before use. Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

The RX Locking Device should be attached to the endoscope prior to the ERCP procedure.

1. Carefully attach the RX Locking Device onto the scope barrel.
2. Rotate so that the body sits flush against the side of the working channel port (See Figure 1).


Figure 1

3. Slide the RX Locking Device up or down the scope barrel and adjust so that, either the "J" slot appears just above the grommet of the working channel. Or, if the device has a locking arm, the bottom of the angle head sits just over the top of the Biopsy Cap (See Figure 2).


Figure 2

4. Place strap around scope in a clockwise direction and connect it to the pin. Make sure the strap is wrapped tightly around the scope to prevent any movement.
5. To insert devices through the Biopsy Cap, align the device with the scope inlet and insert slowly in a straight manner. If not properly done, this may cause more friction on the system and may damage the rubber seal or the scope's working channel. Or to facilitate guidewire introduction into a RX Biliary cannula using a RX Locking Device without a locking arm, place cannula into the slot located on the right side of the RX Locking Device. This will help to align the guidewire in the channel of the cannula.
6. Once the RX system with locking arm is in place, you may lock the guidewire by pushing behind and around the "Locking Arm" and into either side of the "Locking Anchor" (See Figure 3). Or without locking arm, place the guidewire into the "J" slot on the left hand side of the RX Locking Device and slide it to the left side of the "J" slot (See Figure 4).


Figure 3


Figure 4

7. To unlock the guidewire follow the above steps in reverse order.

DEVICE REMOVAL

The RX Locking Device is removed by unlocking the strap and pulling the device off the barrel of the endoscope. The RX Locking Device should be removed and discarded after each patient use.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not to be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

RX Biliary System is a registered trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

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All rights reserved.
© 2020 Boston Scientific Corporation or its affiliates.