IntellaTip MiFi™ Open-Irrigated

Ablation Catheter

Instructions for Use	3
Instrucciones de uso	10
Mode d'emploi	18
Gebrauchsanweisung	25
Istruzioni per l'uso	33
Instructies voor gebruik	41
Instruções de Utilização	49

TABLE OF CONTENTS

DEVICE DESCRIPTION	3
Contents	3
Operating Principle	3
Materials	3
Non-pyrogenic	3
Figure 1. IntellaTip MiFi OI Catheter	3
User Information	3
INTENDED USE / INDICATIONS FOR USE	3
Clinical Benefit Statement	3
Summary of Safety and Clinical Performance	3
CONTRAINDICATIONS	3
PRECAUTIONS	5
ADVERSE EVENTS	5
PATIENT INFORMATION	6
HOW SUPPLIED	6
Device Details	6
Handling and Storage	6
Operating Environment	6
Transport Environment	6
Storage Environment	6
OPERATIONAL INSTRUCTIONS	6
Additional Required Items	6
Accessories	6
Optional Additional Equipment	6
Procedure	6
End of Procedure	7
Disposal	7
Information to Brief the Patient	7
TROUBLESHOOTING	8
WARRANTY	8

IntellaTip MiFi™

Open-Irrigated

Ablation Catheter

ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully read all ancillary device instructions prior to use.

Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

DEVICE DESCRIPTION

The IntellaTip MiFi Open-Irrigated (OI) Ablation Catheter (henceforth, referred to as the IntellaTip MiFi OI Catheter) is a 7.5 F (2.5 mm) quadripolar open-irrigated ablation catheter, designed to deliver Radiofrequency (RF) energy to the 4.5 mm catheter tip electrode for cardiac ablation. The device includes three diagnostic mini electrodes embedded in the tip electrode that are designed to provide additional high resolution and localized electrogram information.

Figure 1 illustrates the IntellaTip MiFi OI Catheter.

The IntellaTip MiFi OI Catheter is designed to be used with a commercially available RF Controller, an Irrigation Pump and Irrigation Tubing Set that meets the catheter flow rate requirements, and a commercially available Mapping System.

Contents

One (1) sterile IntellaTip MiFi OI Catheter

Operating Principle

The IntellaTip MiFi OI Catheter incorporates an open-irrigated cooling mechanism through a tip that is partitioned into two chambers. The proximal chamber circulates heparinized 0.9% normal saline within the tip to cool the proximal end of the tip electrode and mitigate overheating while the distal chamber allows the fluid to flow through six irrigation holes into the patient's vasculature, thereby cooling the tip/tissue interface. A luer connection at the proximal end of the handle connects the catheter to the Irrigation Tubing Set, allowing the Irrigation Pump to generate the flow of saline to the catheter.

The electrode segment is comprised of a tip electrode, three ring electrodes, and includes three diagnostic mini electrodes embedded in the tip. The tip electrode has an embedded temperature sensor and delivers RF energy for cardiac ablation. The ring electrodes record intracardiac Electrogram (EGM) signals for mapping and deliver stimulus for pacing. The three diagnostic mini electrodes are designed to provide additional high resolution, localized electrogram information.

The IntellaTip MiFi OI Catheter interfaces with standard recording equipment. The handle includes the electrical connector for the cable connection to the RF Controller and one luer fitting used to connect the catheter to the Irrigation Tubing Set.

Materials

This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

Non-pyrogenic

This device meets pyrogen limit specifications.


Figure 1. IntellaTip MiFi OI Catheter

User Information

The IntellaTip MiFi OI Catheter is to be used by physicians thoroughly trained in invasive cardiology and in the open-irrigated techniques of RF mapping and ablation, in the specific approach to be used, and in a fully-equipped electrophysiology lab. Assistance to prepare and run the system may only be provided by appropriately trained personnel.

INTENDED USE / INDICATIONS FOR USE

The IntellaTip MiFi OI Catheter is indicated for use in patients who require catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a RF Generator, for cardiac ablation.

Clinical Benefit Statement

When operated by the intended user, the clinical benefits of using cardiac radiofrequency (RF) open irrigated (OI) ablation catheters are the identification and correction of cardiac arrhythmia(s), while decreasing the likelihood of fibrin, thrombus, coagulum and char formation at the catheter tip. Open-irrigated ablation allows for the controlled delivery of higher levels of energy to the heart tissue, changing the ability of the cardiac tissue to conduct abnormal heart rhythms. Greater lesion depth can be achieved with OI ablation, targeting arrhythmic foci effectively through the promotion of adequate lesion penetration. In addition, OI ablation therapy is less invasive compared to open-chest surgical interventions.

The IntellaTip MiFi Open-Irrigated Catheter uses Boston Scientific's Blazer Open-Irrigated RF ablation catheter platform and includes the additional capabilities of the mini-electrodes (ME). The inclusion of the mini-electrodes to the catheter tip provides three localized electrogram (EGM) recordings that improves the ability to resolve three extra intracardiac electrograms during catheter navigation and RF delivery. The ability to capture myocardial electrogram recordings is important in the diagnosis and treatment of arrhythmias. With the addition of mini-electrodes to the catheter, the IntellaTip MiFi Open-Irrigated Catheter provides additional information and context to aid the physician in positioning the catheter for continuous RF ablation lesion lines at target locations.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

The IntellaTip MiFi OI Catheter is contraindicated for use:

- in patients with active systemic infection;
- in patients with a mechanical prosthetic heart valve through which the catheter must pass;
- in patients with conditions where insertion into or manipulation in the cardiac chambers is unsafe as this may increase the risk of perforation or systemic embolic event, such as but not limited to, evidence of myxoma, Left Atrial (LA) thrombus, or intracardiac mural thrombus;
- in patients who are unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation;
- in patients who have vena cava embolic protection filter devices and/or known femoral thrombus who require catheter insertion from the femoral approach;
- in patients who are hemodynamically unstable;
- in patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe, such as but not limited to, a recent previous cardiac surgery: e.g., ventriculotomy or atriotomy, Coronary Artery Bypass Graft (CABG), PTCA/PCI/coronary stent procedure/unstable angina, or in patients with congenital heart disease where the underlying abnormality increases the risk of the ablation, severe rotational anomalies of the heart or great vessels;
- via transseptal approach in patients with an intra-atrial baffle or a foramen ovale patch;
- via retrograde transaortic approach in patients who have had aortic valve replacement.

Do not use this device:

- with a long sheath or a short introducer < 8 F
- in the coronary vasculature

WARNINGS

- If the visibility of the EP catheters are compromised for any reason, the user should stop and not resume ablation therapy until catheter visibility is established in order to prevent patient injuries such as perforation, heart block and injury to adjacent structures.
- Cardiac mapping and ablation procedures should be performed only by physicians thoroughly trained in invasive cardiology and in the techniques of open-irrigated, RF-powered catheter mapping and ablation, and in the specific approach to be used, in a fully-equipped electrophysiology lab.
- Administer appropriate levels of peri-procedural anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedures. There is an increased risk of thromboemboli if appropriate anticoagulation levels are not maintained while the transseptal sheath and/or catheter is in the left side of the heart. Administer anticoagulation therapy during and post-procedure according to the institution's standards to minimize bleeding and thrombotic complications.
- Carefully read all equipment and ancillary device instructions required for the procedure prior to use. Observe all contraindications, warnings, and precautions noted in these directions. Failure to do so may result in patient complications.
- Before using, inspect the IntellaTip MiFi OI Catheter for any defects or physical damage, including electrical insulation on the cables and the catheter shaft that, if used, may cause patient and/or user injury. Do not use defective or damaged devices. Replace damaged equipment if necessary. No modification of this equipment is allowed.
- Contents are supplied STERILE using an Ethylene Oxide (EO) process and should be used by the "Use By" date on the device package. Do not use the device if past the "Use By" date. Do not use if sterile barrier is damaged as use of non-sterile devices may result in patient injury. If damage is found, call your Boston Scientific representative.
- Using the IntellaTip MiFi OI Catheter at lower than prescribed flow rate may increase the potential for thrombus, coagulum, and char that may result in embolism.
- Electromagnetic Interference (EMI) from any source during normal operation may adversely affect the visualization and tracking of the catheter during the procedure, which can cause patient injuries such as perforation, heart block and injury to adjacent structures.
- Start each RF application at low power and carefully follow the power titration and the correlating flow rate procedures as specified in the Operational Instructions section. Too rapid an increase in power during ablation, ablating at high power (> 30 watts) or insufficient flow rate may lead to perforation caused by steam pop, arrhythmias, damage to adjacent structures, and/or embolism.
- Collateral tissue damage is a possibility when using the ablation catheter at the upper power setting (50 watts) with a duration longer than 60 seconds or with a decrease in impedance without moving the tip of the ablation catheter.
- Patients who have had a prior atrial flutter ablation procedure may be at greater risk for perforation and/or pericardial effusion with the use of this catheter system.
- Patients undergoing septal accessory pathway, Atrioventricular (AV) node reentry tachycardia, and/or atrial flutter ablation are at risk for complete AV block which requires the implantation of a temporary and/or permanent pacemaker.
- Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter that may result in embolism.
- During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces to minimize the potential for electrical shock.
- Ensure that the cable/catheter connection remains dry throughout the procedure in order to prevent electric shock or other patient injuries as well as to prevent loss of device function.
- Electrodes and stimulating devices can provide paths of high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes as far away as possible from the ablation site and the Dispersive Pad. Protective impedances may reduce the risk of burns and permit continuous monitoring of the Electrocardiogram (ECG) during energy delivery.
- Before use, ensure irrigation ports are patent and jetting by infusing heparinized normal saline through the catheter tubing. Patency of irrigation ports is important to maintain cooling function and minimize risks of coagulum, fibrin, thrombus and char that may result in embolism as well as perforation caused by steam pop.
- Fibrin may accumulate in or on the sheath/catheter assembly during the procedure. Aspirate when removing the dilator or catheter.
- Do not continue using the catheter if the irrigation ports are occluded or the catheter is not functioning properly.
- To avoid systemic thromboemboli, intravenous heparin or an acceptable alternative must be used when entering the left heart during ablation.
- Patients undergoing a long irrigated ablation procedure have the potential for greater anticoagulation and therefore Activated Coagulation Time (ACT) should be monitored closely due to the increased risk for bleeding/ hemorrhage and/or embolism.
- In the presence of anticoagulation, there may be an increased risk of bleeding from all causes.
- Electrical recording or stimulation equipment must be isolated. Current leakage from any electrical equipment that is connected to the patient must not exceed 10 microamps for intracardiac electrodes.
- Care must be taken to ensure that any equipment used in connection with the BSC catheters be type CF, be defibrillation proof, meet IEC 60601-1 electrical safety requirements, and comply with all local regulatory requirements for specified intended use to reduce the potential risk of inadvertent electrical shock.
- Do not insert or withdraw the IntellaTip MiFi OI ablation catheter without straightening the catheter tip (returning the steering lever to neutral position) in order to prevent entanglement/ entrapment within the valve and/or other device that may result in myocardial trauma and/or may require additional medical/surgical intervention.
- Stimulation of cardiac tissues caused by pacing stimulus and/or RF energy may lead to inadvertent induction of arrhythmias. These arrhythmias may require defibrillation that could also result in skin burns.
- Maximum Catheter Rated Voltage: 150 Vrms (212 Vpk).
- Warnings for patients with implantable pacemakers and Implantable Cardioverter/ Defibrillators (ICDs):
 - Pacemakers, implantable cardioverter/defibrillators, and leads can be adversely affected by radiofrequency energy. It is important to refer to the device manufacturer's instruction for use prior to performing ablation procedures.
 - Do not apply RF energy directly to a lead or to tissue immediately in contact with a lead because it could potentially damage the lead or lead function.
 - Imaging guidance and care must be taken during catheter advancement, manipulation, and withdrawal to avoid lead dislodgement.
 - Temporarily reprogram pacemaker per the manufacturer guidelines during RF ablation to a non-tracking pacing mode if pacing is likely to be required during the ablation. The pacemaker could be damaged by the ablation procedure. Interrogate the device fully after the ablation per the manufacturer guidelines and reprogram to preoperative sensing and pacing parameters.
 - Deactivate ICDs as they could discharge and injure the patient or be damaged by the ablation procedure.
 - Have temporary external sources of pacing and defibrillation available.
 - Perform a complete analysis of the implanted device function after ablation.
 - Fluoroscopic guidance or appropriate imaging and care must be taken during catheter advancement, manipulation, and withdrawal to avoid lead dislodgement.
 - Monitor pre- and post-measurements for sensing and pacing thresholds and impedances to determine the integrity of the lead-patient function.
 - Remember to reactivate the pulse generator after turning off the RF ablation equipment.
 - Do not ablate from within the coronary artery as the resulting myocardial injury can be fatal. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
 - During RF ablation, care must be taken not to deliver RF energy on or near the coronary artery even on the right side of the heart, as the resulting myocardial injury can be fatal.
 - Ablation in contact with any other electrodes alters the function of the catheter and can lead to thrombus, coagulum, or char formation that may result in embolism.
 - At no time should an IntellaTip MiFi OI Catheter be advanced or withdrawn when resistance is felt, without determining the cause. Valve damage, vascular and/or cardiac perforation is a risk with any intracardiac catheter.
 - Catheter entrapment within the heart or blood vessels is a possible complication of cardiac ablation procedures. The potential for catheter entrapment may be increased when the catheter is over torqued and/or positioned in the chordae tendineae. The occurrence of this complication may necessitate surgical intervention and/or repair of injured tissue and/or valve damage.
 - Do not use the IntellaTip MiFi OI Ablation System in the proximity of Magnetic Resonance Imaging (MRI) equipment because the MRI equipment may adversely impact the function of an RF Generator and the ablation system may adversely impact the image quality. This can also lead to loss of visibility during ablation which can cause patient injuries such as perforation, heart block and injury to adjacent structures.
 - Catheter ablation procedures present the potential for significant radiation exposure, which can result in acute radiation injury as well as an increased risk for somatic and genetic effects, to both patients and laboratory staff due to the radiation beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for this use of the device in pregnant women. The long-term risk of protracted fluoroscopy has not been established. Therefore, careful consideration must be given for the use of the device in prepubescent children.
 - There are no data to support the safety and effectiveness of this device in the pediatric population.

- In the event of a suspected failure of the integrity of fluid flow through the catheter or the Irrigation Tubing Set or if there is a rapid temperature rise of $> 15^{\circ}\text{C}$ noted on the RF Controller, the procedure should be stopped, and the catheter withdrawn to reduce the risk of steam pop that could result in adverse events including perforation, embolism or injury to adjacent structures. Both the catheter and the Irrigation Tubing Set should be replaced. The replacement catheter and tubing set must be primed outside the body prior to insertion to reduce risk of air embolism.
- Prior to the procedure, always identify the patient's risk of volume overload. Monitor the patient's fluid balance throughout the procedure and after the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible.
- Excessive curves or kinking of the catheter may damage internal wires and components, including the cooling lumen. This damage may affect steering performance and may cause patient injury.
- Manual bending and/or twisting of the distal curve can damage the steering mechanism and cooling lumens and may cause catheter failure and patient injury.
- Do not scrub the tip electrode, doing so may result in irrigation port(s) occlusion and may lead to catheter failure and patient injury.
- Use both fluoroscopy, or other visualization technique such as echocardiography, and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid conduction pathway injury, cardiac perforation or tamponade.
- Do not deliver RF energy with the catheter outside the target site. RF Generators can deliver significant electrical energy and may cause patient injury.
- In the event of a generator cut-off (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before RF energy is reapplied. Ensure that all of the irrigation holes are patent prior to reuse to reduce the risk of embolism and/or perforation.
- Verify effective contact between the patient and the Dispersive Pad whenever the patient is repositioned as patient movement may disrupt Dispersive Pad contact resulting in patient injury and/or extended procedure times.
- Always verify that the tubing set, catheter and all connections have been properly cleared of air prior to inserting the catheter into the vasculature. Air entrapped in the tubing and catheter can cause potential injury or cardiac arrest. The operator is responsible for removing all air from the system.
- Patients undergoing left-sided ablation procedures should be closely monitored during and post procedure for clinical manifestations of infarction, pulmonary vein injury, nerve damage, embolism and/or atrial esophageal fistula.
- Patients with hemodynamic instability or cardiogenic shock are at increased risk for life-threatening adverse events and ablation must be done with extreme caution.
- The IntellaTip MiFi OI Catheter is not intended to be used for internal cardioversion. Doing so may result in perforation, arrhythmias, embolism, thrombus and/or patient death.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- Inspect irrigation saline for air bubbles and remove any air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause embolism.
- If there is uncertainty regarding the patient's anticoagulation status or rhythm prior to the procedure, there should be a low threshold to perform a Transesophageal Echocardiogram (TEE) prior to the procedure to confirm absence of mural thrombus and/or thrombus in the left atrial appendage.
- Guiding catheters and/or long introducer sheaths present the potential for thromboembolic events. Pre-flush and maintain lumen patency with heparinized intravenous infusion.
- Do not wipe this catheter with organic solvents such as alcohol, or immerse the handle cable connector in fluids. This may result in electrical or mechanical catheter failures. It may also result in an allergic reaction from the patient.
- Irrigation flow during RF ablation may distort distal tip electrogram recordings due to the signal conductivity of the external cooling solution. Careful monitoring of additional intracardiac electrograms during RF application is recommended to reduce the possibility of inadvertent injury to adjacent structures if appropriate. Higher power coupled with higher flow rates may exacerbate the distortion of the EGM signal recordings.
- Pre-procedural anticoagulation therapy is at the discretion of the physician. However, patients with a history of thromboembolic events may require therapeutic anticoagulation therapy, pre-, during and post-ablation to reduce the incidence of major complications. Peri-procedural anticoagulation therapy is recommended for patients undergoing left-sided and transseptal cardiac procedures and should be considered for selected patients undergoing right-sided procedures.
- The safety and/or efficacy of epicardial use of the IntellaTip MiFi OI Catheter has not been evaluated in a clinical trial.
- The Transseptal Procedure (TSP) presents a potential risk for perforation/ tamponade; echocardiography and/or fluoroscopic images should be used to guide the transseptal puncture and a real-time arterial blood pressure monitor should be applied. TSP may induce an air embolus; use proper aspiration and flushing techniques to minimize air embolus.
- Care should be used during multiple sheath/catheter exchanges through the transseptal puncture to avoid causing a residual atrial septal defect that would require repair.
- To avoid patient injury, manipulate the sheath carefully when performing the transseptal puncture especially if the patient has any of the following conditions:
 - Enlarged aortic root
 - Marked right atrial enlargement
 - Small left atrium
 - Marked skeletal deformity or distortion of the thoracic configuration (e.g., scoliosis)
- Any serious incident that has occurred in relation to the device should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS

- The IntellaTip MiFi OI Catheter is designed for use with a compatible RF Controller, Irrigation Pump and Irrigation Tubing Set that meets the catheter flow rate requirements, a compatible Mapping System, and compatible Connection Box.
- Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the electrode will not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode.
- Use only sterile saline and gauze pad to clean the tip.
- The IntellaTip MiFi OI Catheter is not intended to be used with a RF Generator output setting exceeding 50 W or 212 V peak.
- Verify the RF Controller is in the control mode which will deliver the amount of power specified by the power setting unless the measured temperature exceeds the temperature setting. Temperature controlled RF delivery may be affected by the cooling effects of the saline irrigation of the electrode. For example, the Maestro RF Cardiac Ablation Controller has these settings in the power control mode.
- Equipment/accessories carrying high frequency alternating current may cause direct coupled interference and therefore, may disrupt the operation of the RF Generator. It may be necessary to take risk control measures, such as re-orienting, relocating, or shielding the interfering equipment/accessories.
- Use only Dispersive Pads that meet or exceed IEC 60601-2-2 requirements and follow the Dispersive Pad manufacturer's instructions for use. The use of Dispersive Pads which meet ANSI/AAMI requirements (HF18) is recommended.
- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the Dispersive Pad or failure of an electrical lead.
- The IntellaTip MiFi OI Catheter is highly torqueable. Avoid overtorquing. Over-rotating the handle and catheter shaft may cause damage to the distal tip or catheter assembly. Do not rotate the handle and catheter shaft more than one and one-half times the full rotation (540°). If the desired catheter tip position is not achieved, adjust the catheter's curve to disengage the catheter tip from the heart wall before resuming rotation of the handle and catheter shaft.
- Electrophysiology catheters and systems are intended for use only in radiation shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.
- The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
- After use, handle and dispose of product and packaging in accordance with hospital biohazard procedure, administrative and/or local government policy.

ADVERSE EVENTS

Potential adverse events associated with use of the IntellaTip MiFi OI catheter include, but are not limited to:

- Pain or discomfort, for example:
 - Angina
 - Chest pain
 - Non-cardiovascular pain
- Cardiac arrest
- Death
- Hypertension
- Hypotension
- Infection/inflammation/exposure to biohazardous material
- Edema/heart failure/pleural effusion
- Procedural related side effects, for example:
 - Allergic reaction (including anaphylaxis)

- Genitourinary complication
- Side effects related to medication or anesthesia
- Radiation injury/tissue burn
- Renal failure/insufficiency
- Vasovagal response
- Respiratory distress/insufficiency/dyspnea
- Arrhythmia (new or exacerbated)
 - Conduction pathway injury (heart block, nodal injury, etc.)
- Nerve injury, for example:
 - Phrenic nerve injury
 - Vagal nerve injury
- Gastrointestinal disorders
- Vessel trauma, including:
 - Perforation
 - Dissection
 - Coronary artery injury
 - Vasospasm
 - Occlusion
 - Hemothorax
- Cardiac trauma, for example:
 - Cardiac perforation/cardiac tamponade/pericardial effusion
 - Valvular damage
 - Stiff left atrial syndrome
- Injury related to tissue damage and/or adjacent structures, for example:
 - Esophageal injury
 - Pulmonary injury
 - Catheter entrapment
 - Physical trauma
 - Ablation energy
- Fistula, for example:
 - Atrio-esophageal fistula
 - Bronchopericardial fistula
- PV stenosis and its symptoms, for example:
 - Cough
 - Shortness of breath, fatigue
 - Hemoptysis
- Surgical and access complications, for example:
 - Hematoma/seroma
 - AV fistula
 - Bleeding
 - Pseudoaneurysm
 - Pneumothorax
 - Residual atrial septal defect
- Injury due to embolism/thromboembolism/air embolism/foreign body embolism
 - Cerebrovascular Accident (CVA)/stroke
 - Transient Ischemia Attack (TIA)
 - Myocardial infarction
 - Neurological impairment and its symptoms, for example:
 - Cognitive changes, visual disturbances, headache, motor impairment, sensory impairment, and speech impairment
 - Pulmonary embolism
 - Asymptomatic cerebral embolism

The potential adverse events may be related to the ablation catheter(s) and/or the interventional procedure. The severity and/or the frequency of these potential adverse events may vary and may result in prolonged procedure time and/or additional medical and/or surgical intervention, implantation of a permanent device such as a pacemaker, and in rare cases, may result in death.

PATIENT INFORMATION

To screen patients for left ventricular or atrial thrombus or myxoma, it is recommended that patients undergo surface or transesophageal echocardiography or a comparable cardiac imaging study prior to the ablation procedure.

A patient information brochure titled “Understanding Arrhythmias” is available through your BSC representative.

HOW SUPPLIED

One (1) IntellaTip MiFi OI Catheter is supplied sterile using an Ethylene Oxide (EO) process. In addition to the IntellaTip MiFi OI Catheter, please refer to the Additional Required Items section below for a detailed list of other materials typically required in an Electrophysiology (EP) procedure.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Do not use the device if past the “Use By” date.

Device Details

Report any serious incident that occurs in relation to this device to Boston Scientific and to the relevant local regulatory authority for medical devices in your country.

Handling and Storage

Operating Environment

Ambient Temperature: 10 °C to 40 °C

Relative Humidity: 30 % to 75 %

Atmospheric Pressure: 70 kPa to 106 kPa

Transport Environment

Temperature: -29 °C to 60 °C

Relative Humidity: 30 % to 85 %

Atmospheric Pressure: Uncontrolled

Storage Environment

Ambient Temperature: 15 °C to 25 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

OPERATIONAL INSTRUCTIONS

Additional Required Items

Intracardiac electrophysiology and cardiac ablation procedures should be performed in a fully equipped clinical electrophysiology lab.

In addition to the IntellaTip MiFi OI Catheter the following devices and materials are required (refer to corresponding manufacturer’s operator’s manuals for specific material information):

Note: The IntellaTip MiFi OI Catheter is not designed to be compatible with the Maestro 3000 Cardiac Ablation system.

- Compatible RF Controller and accessories
- Compatible Irrigation Pump and accessories
- Irrigation Tubing Set

Note: The RF Generator and pump/tubing system must be capable of delivering flow rates of 2 mL/min (Standby), 17 mL/min (Low Ablation Flow Rate at ≤ 30 W), 30 mL/min (High Ablation Flow Rate at 31 W to 50 W).

- Cable, IntellaTip MiFi OI to RF Controller

Accessories

- Commercially available disposable Dispersive Pads which are in compliance with IEC 60601-2-2
- Sterile, normal (0.9%), heparinized (1 u heparin/mL) saline (commercially available)
- 8F (2.67 mm) Venous Introducer Sheath

Optional Additional Equipment

- Compatible Mapping System and accessories
- Compatible Connection Box
- Filter Module

Procedure

CAUTION: Before use, inspect the packaging for any violation of the sterile barrier and inspect the IntellaTip MiFi OI Catheter for any defects. Do not use potentially contaminated or defective equipment.

Please refer to the operator’s manuals and Instructions for Use for the Irrigation Pump, RF Controller, Mapping System, Connection Box, Filter Module, and the Irrigation Tubing Set for instructions on connecting and operating these systems in conjunction with the IntellaTip MiFi OI Catheter. Use appropriate accessory cables to connect the catheter to the appropriate accessory equipment.

1. Attach the Dispersive Pad to the patient and generator per the manufacturer’s instructions for use.
2. Connect the patient to an ECG recording system to facilitate arrhythmia monitoring per the standard operating procedure of the electrophysiology lab or manufacturer’s operator’s manual.

Note: This should be done prior to introducing any intracardiac catheters.

- Open the IntellaTip MiFi OI Catheter and Cable packages and the Irrigation Tubing Set package. Carefully transfer the package contents into the sterile field, maintaining sterile technique.
- Obtain vascular access via a vein (e.g., a femoral vein) under aseptic conditions. Then place an introducer sheath into the vein using a standard percutaneous technique.
- Connect the RF Controller to the recording system (and the Mapping System if desired) with the appropriate interface cables according to the operator's manuals and/or IFUs.
- Connect the IntellaTip MiFi OI Catheter to the RF Controller. Use the Filter Module and Connection Box as needed per the operator's manuals and IFUs.

Note: If a three-dimensional (3-D) catheter navigation and mapping system is going to be used, please follow the standard operating procedure of the electrophysiology lab or the directions for use contained in the manufacturer's operator's manual.

- Turn ON the power to the RF Controller.
- Set the RF Controller to Power Control mode. The RF Controller's temperature limit must not exceed 50 °C, but can be set lower at physician discretion.
- Turn on the Irrigation Pump.
- Make sure that the Irrigation Pump has the following flow rates: 2 mL/min (Standby); 17 mL/min (Low Ablation Flow—30 W or less); 30 mL/min (High Ablation Flow—above 30 W). Refer to the Irrigation Pump operator's manual for instructions on how to adjust the pump settings if required.
- Refer to either the Irrigation Tubing Set or Irrigation Pump IFU for instructions to connect the Irrigation Tubing Set to irrigation fluid and Irrigation Pump.
- Connect the IntellaTip MiFi OI Catheter to the Irrigation Tubing Set via the luer fitting at the proximal end of the catheter handle. Care must be taken to ensure all luer fittings are secure to prevent leaking.
- Purge the IntellaTip MiFi OI Catheter and Irrigation Tubing Set. Fluid should exit all six (6) irrigation ports during the flushing process. Assure that no air remains within the Irrigation Tubing Set or lumen and all irrigation ports are patent.
- Set the pre-RF delay and post-RF delay on the Irrigation Pump. Reference the Irrigation Pump operator's manual for instructions on how to change the pre-RF delay and post-RF delay.
- Check the catheter steering by articulating the steering knob prior to inserting the catheter in the sheath.
- Before placing the IntellaTip MiFi OI Catheter in the sheath, begin continuous irrigation at a flow rate of 2 mL/min, i.e., standby flow. Check for any leaks at the tip of the catheter (other than normal saline flowing out of the distal ports), at the catheter handle, and at the luer connections and tubing joints.
- Under fluoroscopic guidance, insert the IntellaTip MiFi OI Catheter into the sheath and advance through the vasculature into the heart.

Note: The degree of tip deflection of the IntellaTip MiFi OI Catheter is controlled by the Steering Knob on the catheter handle (see Figure 1). If the Steering Knob is turned in a clockwise direction from its neutral position, the tip will curve proportionately up to a maximum of 270 degrees in one direction depending upon the curve option selected. Turning the Steering Knob in the counter clockwise direction will cause the tip to deflect in the opposite direction. To prevent overstressing the tip, the Steering Knob movement is limited by the handle design. The tension adjust knob may be used when the desired catheter placement is achieved.

- Determine the area of interest for ablation.
- Set the initial power level from 15 W to 20 W.
- Start the procedure from 15 W to 20 W. Power may be increased by 5 W to 10 W increments as needed to create a transmural lesion. A greater than 80 % reduction in unipolar electrogram amplitude or emergence of double potentials of equal and low amplitude may be indicators of a transmural lesion.

Note: Confirm the increased irrigation flow rate prior to onset of RF energy by observation of a decrease in tip electrode temperature of at least a 2 °C. If it is necessary to ablate with power levels of 31 W to 50 W, increase the irrigation flow rate to 30 mL/min before the onset of RF energy delivery. Then return the flow rate to 2 mL/min post-RF energy delivery.

WARNING: Too rapid an increase in power during ablation, ablating at high power (> 30 W) or insufficient flow rate may lead to perforation caused by steam pop, arrhythmias, damage to adjacent structures, and/or embolism.

WARNING: Using the IntellaTip MiFi OI Catheter at a lower than prescribed flow rate may increase the potential for thrombus, coagulum and char that may result in embolism.

- Do not ablate for greater than 60 seconds in duration without moving the tip of the IntellaTip MiFi OI Catheter.
- RF current may be reapplied to the same or alternate sites using the same catheter.

Note: The mini electrodes on the IntellaTip MiFi OI Catheter may also be used to provide additional electrogram information.

End of Procedure

- Prior to removing the catheter, straighten the distal end of the catheter completely.
- Withdraw the catheter when the procedure is finished.
- Turn off RF Generator and pump.
- Carefully monitor patient while in recovery to ensure hemostasis is achieved and any complications are immediately treated.

Disposal

To minimize risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, the catheter may contain biohazardous substances. The catheter and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

In the event that a serious incident occurred in relation to the device, including all patient deaths for procedures where the BSC product was used, the event should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established. Return any catheter related to a complaint, patient harm, injury, or death to Boston Scientific using a BSC Returned Product Kit.

- Returning products for analysis and providing product performance observations helps drive reliability higher on an ongoing basis.
- Be sure to follow the instructions with regard to packaging and shipping any biohazard devices taking care not to expose the handle or connector to fluid that may compromise the catheter and limit analysis.

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of the IntellaTip MiFi OI Catheter in association with the Electrophysiological cardiac interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document.
- Discuss post procedure instructions, including any lifestyle changes, medications, when to call the Healthcare Provider (HCP) and any post procedure follow-up that might be needed.

TROUBLESHOOTING

Problems	Possible Cause	Corrective Action Procedure
Temperature not displayed	Poor catheter/cable connections	1. Verify that the Cable is plugged into the IntellaTip MiFi OI Catheter and/or the RF Controller, Filter Module, and Connection Box. 2. Verify that the Connection Box (if using a Mapping System) is connected to the RF Controller. 3. Replace cable and/or catheter. 4. If the RF Controller still does not display temperature, there may be a malfunction in the temperature sensing system. 5. Consult the operator's manual and correct this malfunction prior to reapplying RF energy.
• Impedance cutoff • Temperature cutoff	Char/coagulum on tip electrode	1. Discontinue RF delivery. 2. Straighten the distal end and withdraw IntellaTip MiFi OI Catheter. 3. Inspect tip electrode for any char/coagulum. 4. If present, gently wipe the tip section with a sterile gauze dampened with sterile saline (do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode). 5. Prior to reinsertion, ensure the irrigations ports are patent. If irrigation port occlusion occurs: a. Ensure IntellaTip MiFi OI Catheter is removed from the patient. b. Fill a 1 mL or 2 mL syringe with sterile saline and attach to the stop-cock sidearm of the IntellaTip MiFi OI Catheter. c. Carefully inject the saline from the syringe into the IntellaTip MiFi OI Catheter. Fluid should exit all six (6) irrigation ports during the flushing process. d. Repeat steps b and c, if necessary. e. If the irrigation ports are cleared, the IntellaTip MiFi OI Catheter can be reintroduced into the patient. WARNING: Do not continue to use the IntellaTip MiFi OI Catheter if the irrigation ports are occluded or the catheter is not functioning properly.
Suspected failure of fluid flow integrity	• Leak in catheter and/or Irrigation Tubing Set • Irrigation Pump out of calibration	1. Discontinue RF delivery. 2. Straightening the distal end and withdraw catheter. 3. Replace IntellaTip MiFi OI Catheter and Irrigation tubing set, prime outside of the patient. 4. Replace IntellaTip MiFi OI Catheter and/or Irrigation tubing set if parameters do not appear normal or if there is any abnormality of the integrity of fluid flow. 5. Refer to the Irrigation Pump operator's manual to verify fluid flow is accurate. 6. Contact BSC representative to replace Irrigation pump.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

IntellaTip and Maestro are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

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All rights reserved.
© 2022 Boston Scientific Corporation or its affiliates.