IntellaNav StablePoint™

Ablation Catheter

Instructions for Use	2
Instrucciones de uso	11
Mode d'emploi	20
Gebrauchsanweisung	29
Istruzioni per l'uso	39
Instructies voor gebruik	48
Instruções de utilização	57

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Operating Principle	3
Figure 1. IntellaNav StablePoint Ablation Catheter	3
Local Impedance During Ablation - DIRECTSENSE	3
Figure 2. Local impedance elements during RF ablation	3
FORCE VISUALIZATION ON RHYTHMIA	4
Figure 3. Visualization of Force on Compatible Rhythmia Mapping System	4
Contents	4
Materials	4
Non-Pyrogenic	4
User Information	4
INTENDED USE / INDICATIONS FOR USE	4
CLINICAL BENEFIT STATEMENT	4
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	4
CONTRAINDICATIONS	4
WARNINGS	4
PRECAUTIONS	6
ADVERSE EVENTS	7
CLINICAL STUDIES	7
HOW SUPPLIED	7
Device Details	7
Handling and Storage	7
Operating Environment	7
Transport Environment	7
Storage Environment	7
OPERATIONAL INSTRUCTIONS	7
Additional Items for Safe Use	7
Accessories	8
Optional Additional Equipment:	8
Preparation	8
Procedure	8
Disposal	9
Post-Procedure	9
INFORMATION TO BRIEF THE PATIENT	9
TROUBLESHOOTING	10
WARRANTY	10

IntellaNav StablePoint™ Ablation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific (BSC) representative.

For single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully read all ancillary device instructions prior to use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to do so may result in patient complications.

DEVICE DESCRIPTION

The IntellaNav StablePoint Ablation Catheter (henceforth, referred to as the IntellaNav StablePoint) is a steerable quadripolar open-irrigated ablation catheter designed to deliver radiofrequency (RF) energy to the 4 mm catheter tip electrode for cardiac ablation. The catheter shaft is 7.5F with 8F ring electrodes. The IntellaNav StablePoint is compatible with introducers or sheaths with a minimum diameter of 8.5F.

The IntellaNav StablePoint incorporates a position sensor for magnetic tracking and navigation, when used with a compatible Rhythmia™ Mapping System. Additionally, the catheter has an embedded force sensor in the distal tip to transmit real-time feedback on the mechanical interaction between the RF tip electrode and myocardial tissue and the catheter is enabled to measure changes in local dielectric properties in the proximity of the tip electrode.

Note: See corresponding Instructions for Use for compatible Rhythmia Mapping System for additional information on the contact force visualization and the local impedance measurement, hereby referred to as DIRECTSENSE™ on IntellaNav StablePoint Catheter.

For ablation, the IntellaNav StablePoint is designed to be used with a commercially available RF Controller (also called RF Generator), an Irrigation Pump and Irrigation Tubing Set that meets the catheter flow rate requirements¹, a commercially available Connection Box, and a dispersive pad (indifferent electrode). For mapping, navigation, and visualization of force and DIRECTSENSE information, the IntellaNav StablePoint is designed to be used with a compatible Rhythmia Mapping System.

Operating Principle

The IntellaNav StablePoint incorporates an open-irrigated cooling mechanism through a tip that is partitioned into two chambers. The proximal chamber circulates heparinized normal saline (0.9%) within the tip to cool the proximal end of the tip electrode and mitigate overheating while the distal chamber allows the fluid to flow through six irrigation holes, thereby cooling the tip/tissue interface. A luer connection at the proximal end of the handle connects the catheter to the Irrigation Tubing Set, allowing the Irrigation Pump to generate the flow of saline to the catheter. A thermocouple temperature sensor embedded in the tip provides feedback on the tip cooling.

The electrode segment is comprised of a tip electrode and three ring electrodes. All the electrodes can be used for recording intracardiac electrograms (EGM) or delivering pacing stimuli from external systems. The tip electrode transmits RF energy for cardiac ablation. The IntellaNav StablePoint interfaces with standard RF Controllers and recording equipment through the Connection Box. The handle includes the electrical connector for the cable connection to the mapping system. The IntellaNav StablePoint is shown in Figure 1.


Figure 1. IntellaNav StablePoint Ablation Catheter

Local Impedance During Ablation - DIRECTSENSE

The DIRECTSENSE on IntellaNav StablePoint feature provides a measure of the dielectric properties closest to the catheter tip electrode. This diagnostic metric can be used in conjunction with other diagnostic elements (e.g., electrogram amplitude, fluoroscopy, intracardiac echocardiography, magnetic and impedance navigation, contact force and tactile feedback) to inform the user on stability and proximity of the catheter electrodes to the endocardial surface before RF energy is applied.

Note: Local impedance is not an indicator of contact force.

During the application of RF energy, the local impedance measure provides additional real-time feedback on dielectric changes near the RF electrode resulting from resistive tissue heating relative to the pre-ablation value. Due to the change in conductivity during tissue heating, local impedance may not represent catheter proximity or relative position of the catheter tip-to-tissue during RF, as it does when RF is OFF.

Note: Local impedance during ablation is only supported by a compatible Rhythmia Mapping System.

Note: RF ablation effectiveness is influenced by a number of factors including: tissue thickness, contact, force, stability, catheter orientation, power, duration, and irrigation flow. Users should consider these factors when performing ablation and confirm effectiveness through functional endpoints such as arrhythmia termination or establishing conduction block.

During the application of RF energy, a compatible Rhythmia Mapping System provides a measure of local impedance as a numerical value widget, a power bar graphic, a tip graphic in orange, and a real-time graph (Figure 2). The value widget is updated to display the change in local impedance from the onset of ablation. The change in local impedance is displayed in orange to match the other ablation color indicators and can be displayed as an absolute, relative or percent value. Once RF energy is terminated, the real-time graph will continue to display an orange overlay to temporally indicate that ablation has occurred at previous data epochs.


Figure 2. Local impedance elements during RF ablation.

A) Generator parameters widget. B) Numerical Value (average impedance) widget during RF. C) Catheter tip graphic during RF. D) Local impedance vs. time trace during RF.

Changes in the local impedance during RF delivery require stable catheter position. Use fluoroscopy or other visualization techniques such as echocardiography to verify catheter location during RF. Incorrect catheter localization may lead to misinterpretation of the impedance measure and an incorrect clinical conclusion or patient injury.

¹Must be capable of delivering flow rates of 2mL/min (Standby), 17 mL/min (Low Power Ablation Flow Rate at 30 W), and 30 mL/min (High Power Ablation Flow Rate at 31 W - 50 W).

FORCE VISUALIZATION ON RHYTHMIA

A compatible Rhythmia Mapping System provides the visualization of the force information in a similar set of widgets comprised of a catheter tip visualization, a force value widget, a force angle indicator, and a real-time force graph (Figure 3). The real time force graph will also have an orange overlay during RF delivery that will persist over the segment that represents the ablation while that data is in the field of view. The average contact force and the variability of contact force can be used within the context of other additional parameters (e.g. catheter position, fluoroscopy, intracardiac electrograms, and tactile feedback) to aid in the positioning of the catheter prior to and during RF delivery.


Figure 3. Visualization of Force on Compatible Rhythmia Mapping System

A) Catheter Tip Visualization, B) Force Value Widget and Force Angle Indicator, and C) Force Real Time Graph with user defined range of interest and user defined high force threshold.

Contents

One (1) Sterile IntellaNav StablePoint Catheter

Materials

This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

Non-Pyrogenic

This device meets pyrogen limit specifications.

User Information

The IntellaNav StablePoint Catheter is to be used by physicians thoroughly trained in invasive cardiology and in the open-irrigated techniques of RF powered catheter mapping and ablation, in the specific approach to be used, and in a fully-equipped electrophysiology lab. Assistance to prepare and run the system may only be provided by appropriately trained personnel.

INTENDED USE / INDICATIONS FOR USE

The IntellaNav StablePoint Catheter is indicated for use in patients who require catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with an RF Controller, for cardiac ablation.

CLINICAL BENEFIT STATEMENT

When operated by the intended user, the clinical benefits of using cardiac radiofrequency (RF) open irrigated (OI) ablation catheters are the identification and correction of cardiac arrhythmia(s), while decreasing the likelihood of fibrin, thrombus, coagulum and char formation at the catheter tip. Open-irrigated ablation allows for the controlled delivery of higher levels of energy to the heart tissue, changing the ability of the cardiac tissue to conduct abnormal heart rhythms. Greater lesion depth can be achieved with OI ablation, targeting arrhythmic foci effectively through the promotion of adequate lesion penetration. In addition, OI ablation therapy is less invasive compared to open-chest surgical interventions.

The IntellaNav StablePoint catheter uses Boston Scientific's Blazer Open Irrigated RF ablation catheter platform and includes the additional state of the art capabilities of a magnetically tracked navigation sensor and contact force technology. The inclusion of a magnetic sensor to the ablation catheter allows the intended user to track the catheter's 3D position when used with the Rhythmia Mapping System. Localization and navigation of intracardiac catheters is important in the diagnosis and treatment of arrhythmias. With the addition of local impedance (through DIRECTSENSE) to the contact force enabled catheter, the IntellaNav StablePoint Catheter provides additional information and context to aid the physician in balancing the formation of transmural lesions with excessive temperatures that may lead to unintended events.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website: <https://ec.europa.eu/tools/eudamed>

CONTRAINDICATIONS

The IntellaNav StablePoint Catheter is contraindicated for use:

- in patients with active systemic infection;
- in patients with a mechanical prosthetic heart valve through which the catheter must pass;
- in patients with conditions where insertion into or manipulation in the cardiac chambers is unsafe as this may increase the risk of perforation or systemic embolic event, such as but not limited to, evidence of myxoma, Left Atrial (LA) thrombus, or intracardiac mural thrombus;
- in patients who are unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation;
- in patients who have vena cava embolic protection filter devices and/or known femoral thrombus who require catheter insertion from the femoral approach;
- in patients who are hemodynamically unstable;
- in patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe, such as but not limited to, a recent previous cardiac surgery (e.g., ventriculotomy or atriotomy, Coronary Artery Bypass Graft [CABG], PTCA/PCI/coronary stent procedure/unstable angina) and/or in patients with congenital heart disease where the underlying abnormality increases the risk of the ablation procedure (e.g., severe rotational anomalies of the heart or great vessels);
- via transseptal approach in patients with an intra-atrial baffle or a foramen ovale patch;
- via retrograde transaortic approach in patients who have had aortic valve replacement.

Do not use this device:

- With a long sheath or a short introducer < 8.5F
- In the coronary vasculature

WARNINGS

- If the visibility of the EP catheters are compromised for any reason, the user should stop and not resume ablation therapy until catheter visibility is established in order to prevent patient injuries such as perforation, heart block and injury to adjacent structures.
- Cardiac mapping and ablation procedures should be performed only by physicians thoroughly trained in invasive cardiology and in the techniques of open-irrigated RF powered catheter mapping and ablation, and in the specific approach to be used, in a fully-equipped electrophysiology lab.
- Administer appropriate levels of peri-procedural anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedures. There is an increased risk of thromboemboli if appropriate anticoagulation levels are not maintained while the transseptal sheath and/or catheter is in the left side of the heart. Administer anticoagulation therapy during and post-procedure according to the institution's standards to minimize bleeding and thrombotic complications.
- Carefully read all equipment and ancillary device instructions required for the procedure prior to use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to do so may result in patient complications.
- Before using, inspect the IntellaNav StablePoint for any defects or physical damage, including electrical insulation on the cables and the catheter shaft that, if used, may cause patient and/or user injury. Do not use defective or damaged devices. Replace damaged equipment if necessary. No modification of this equipment is allowed.
- Contents are supplied STERILE using an Ethylene Oxide (EO) process and should be used by the "Use By" date on the device package. Do not use the device if past the "Use By" date. Do not use if sterile barrier is damaged as use of non-sterile devices may result in patient injury. If damage is found, call your Boston Scientific representative.
- Using the IntellaNav StablePoint at lower than prescribed flow rate may increase the potential for thrombus, coagulum, and char that may result in embolism.
- Electromagnetic Interference (EMI) from any source during normal operation may adversely affect the visualization and tracking of the catheter during the procedure, which can cause patient injuries such as perforation, heart block and injury to adjacent structures.

Note: The IntellaNav StablePoint is not designed to be compatible with the Maestro 3000™ Cardiac Ablation System.

- Carefully follow the power and the correlating flow rate procedures as specified in the Operational Instructions. A drop in impedance may be an indicator of lesion creation. Too rapid an increase in power during ablation, increasing power with a decrease in impedance, ablating at high power (> 30W) or insufficient flow rate may lead to perforation caused by steam pop, arrhythmias, damage to adjacent structures, and/or embolism.
- Collateral tissue damage is a possibility when using the ablation catheter at the upper power setting (50W) with a duration longer than 60 seconds or with a decrease in impedance without moving the tip of the ablation catheter. Power should be increased to > 30W only if lower energies do not achieve the intended result.
- Patients who have had a prior atrial flutter ablation procedure may be at greater risk for perforation and/or pericardial effusion with the use of this catheter system.
- Patients undergoing septal accessory pathway, Atrioventricular (AV) node reentry tachycardia, and/or atrial flutter ablation are at risk for complete AV block which requires the implantation of a temporary and/or permanent pacemaker.
- Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter that may result in embolism.
- During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces to minimize the potential for electrical shock.
- Ensure that the cable/catheter connection remains dry throughout the procedure in order to prevent electric shock or other patient injuries as well as to prevent loss of device function.
- Electrodes and stimulating devices can provide paths of high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes as far away as possible from the ablation site and the dispersive pad. Protective impedances may reduce the risk of burns and permit continuous monitoring of the electrocardiogram (ECG) during energy delivery.
- Before use, ensure irrigation ports are patent and jetting by infusing heparinized normal saline through the catheter tubing. Patency of irrigation ports is important to maintain cooling function and minimize risks of coagulum, fibrin, thrombus, and char that may result in embolism as well as perforation caused by steam pop.
- Fibrin may accumulate in or on the sheath/catheter assembly during the procedure. Aspirate when removing the dilator or catheter.
- Do not continue using the catheter if the irrigation ports are occluded or the catheter is not functioning properly.
- Due to the design of the IntellaNav StablePoint tip, the velocity of fluid exiting the irrigation ports may change based on rate and pressure of flushing. As long as there is fluid jetting out of each port, regardless of the velocity, the catheter is functioning as designed and may be used. However, if any irrigation port has no flow (or extremely low flow compared to adjacent ports) despite attempts to flush the irrigation port, do not insert the catheter in the patient as there may be potential risk of embolism.
- To avoid systemic thromboemboli, intravenous heparin or an acceptable alternative must be used when entering the left heart during ablation.
- Patients undergoing a long irrigated ablation procedure have the potential for greater anticoagulation and therefore Activated Coagulation Time (ACT) should be monitored closely due to the increased risk for bleeding/hemorrhage and/or embolism.
- In the presence of anticoagulation, there may be an increased risk of bleeding from all causes.
- Electrical recording or stimulation equipment must be isolated. Current leakage from any electrical equipment that is connected to the patient must not exceed 10 microamps for intracardiac electrodes.
- Care must be taken to ensure that any equipment used in connection with the BSC catheters be type CF, be defibrillation proof, meet IEC 60601-1 electrical safety requirements, and comply with all local regulatory requirements for specified intended use to reduce the potential risk of inadvertent electrical shock.
- Do not insert or withdraw the IntellaNav StablePoint Catheter without straightening the catheter tip (returning the steering lever to neutral position) in order to prevent entanglement/ entrapment within the valve and/or other device that may result in myocardial trauma and/or may require additional medical/surgical intervention.
- Stimulation of cardiac tissues caused by pacing stimulus and/or RF energy may lead to inadvertent induction of arrhythmias. These arrhythmias may require defibrillation that could also result in skin burns.
- Maximum Catheter Rated Voltage: 122 Vrms (173 Vpk).
- Warnings for patients with implantable pacemakers and implantable cardioverter/defibrillators (ICDs):
 - Pacemakers, implantable cardioverter/defibrillators, and leads can be adversely affected by radiofrequency energy. It is important to refer to the device manufacturer's instruction for use prior to performing ablation procedures.
 - Do not apply RF energy directly to a lead or to tissue immediately in contact with a lead because it could potentially damage the lead or lead function.
- Temporarily reprogram pacemaker per the manufacturer guidelines during RF ablation to a non-tracking pacing mode if pacing is likely to be required during the ablation. The pacemaker could be damaged by the ablation procedure. Interrogate the device fully after the ablation per the manufacturer guidelines and reprogram to preoperative sensing and pacing parameters.
- Deactivate ICDs as they could discharge and injure the patient or be damaged by the ablation procedure.
- Have temporary external sources of pacing and defibrillation available.
- Perform a complete analysis of the implanted device function after ablation.
- Fluoroscopic or appropriate imaging guidance and care must be taken during catheter advancement, manipulation, and withdrawal to avoid lead dislodgement.
- Monitor pre- and post-measurements for sensing and pacing thresholds and impedances to determine the integrity of the lead-patient function.
- Remember to reactivate the pulse generator after turning off the RF ablation equipment.
- Do not ablate from within the coronary artery as the resulting myocardial injury can be fatal. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
- During RF ablation, care must be taken not to deliver RF energy on or near the coronary artery even on the right side of the heart, as the resulting myocardial injury can be fatal.
- Ablation in contact with any other electrodes alters the function of the catheter and can lead to thrombus, coagulum, or char formation that may result in embolism.
- At no time should an IntellaNav StablePoint Catheter be advanced or withdrawn when resistance is felt without determining the cause. Valve damage, vascular and/or cardiac perforation is a risk with any intracardiac catheter.
- Catheter entrapment within the heart or blood vessels is a possible complication of cardiac ablation procedures. The potential for catheter entrapment may be increased when the catheter is over torqued and/or positioned in the chordae tendineae. The occurrence of this complication may necessitate surgical intervention and/or repair of injured tissue and/or valve damage.
- Do not use the IntellaNav StablePoint Catheter in the proximity of magnetic resonance imaging (MRI) equipment because the MRI equipment may adversely impact the function of an RF Controller and the ablation system may adversely impact the image quality. This can also lead to loss of visibility during ablation which can cause patient injuries such as perforation, heart block, and injury to adjacent structures.
- The StablePoint Catheter is not intended to be used with a RF Controller output setting exceeding 50 W or 173 Vpk. The safety and performance of the StablePoint Catheter at powers exceeding 50W has not been evaluated in a clinical trial. Exceeding recommended power settings or using flow rates lower than recommended settings may increase the risk of patient injury.
- Catheter ablation procedures present the potential for significant radiation exposure, which can result in acute radiation injury as well as an increased risk for somatic and genetic effects, to both patients and laboratory staff due to the radiation beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for this use of the device in pregnant women. The long-term risk of protracted fluoroscopy has not been established. Therefore, careful consideration must be given for the use of the device in prepubescent children.
- There are no data to support the safety and effectiveness of this device in the pediatric population.
- In the event of a suspected failure of the integrity of fluid flow through the catheter or the tubing set or if there is a rapid temperature rise of >15 °C noted on the RF Controller, the procedure should be stopped, and the catheter withdrawn to reduce the risk of steam pop that could result in adverse events, including perforation, embolism or injury to adjacent structures. Both the catheter and the Irrigation Tubing Set should be replaced. The replacement catheter and tubing set must be primed outside the body prior to insertion to reduce risk of air embolism.
- Prior to the procedure, always identify the patient's risk of volume overload. Monitor the patient's fluid balance throughout the procedure and after the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible.
- Excessive curves or kinking of the catheter may damage internal wires and components, including the cooling lumen. This damage may affect steering performance and may cause patient injury.
- Manual bending and/or twisting of the distal curve can damage the steering mechanism and cooling lumens and may cause catheter failure and patient injury.
- Do not scrub the tip electrode as this may result in irrigation port(s) occlusion and may lead to catheter failure and patient injury.

- Use both fluoroscopy, or other visualization technique such as echocardiography, and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid conduction pathway injury, cardiac perforation or tamponade.
- Do not use excessive force to advance or withdraw the catheter. The firmness of the tip dictates that care shall be taken to prevent perforation of the heart during catheter manipulation. If the force-sensing feature is active, evaluate applied force to avoid applying excessive loads.
- Do not deliver RF energy with the catheter outside the target site. RF Controllers can deliver significant electrical energy and may cause patient injury.
- In the event of a generator cut-off (impedance or temperature), the IntellaNav StablePoint Catheter must be withdrawn and the tip electrode cleaned of coagulum before RF energy is reapplied. Ensure that all of the irrigation holes are patent prior to reuse to reduce the risk of embolism and/or perforation.
- Verify effective contact between the patient and the Dispersive Pad whenever the patient is repositioned as patient movement may disrupt Dispersive Pad contact resulting in patient injury and/or extended procedure times.
- Always verify that the tubing set, catheter, and all connections have been properly cleared of air prior to inserting the catheter into the vasculature. Air entrapped in the tubing and catheter can cause potential injury or cardiac arrest. The operator is responsible for removing all air from the system.
- Patients undergoing left sided ablation procedures should be closely monitored during and post procedure for clinical manifestations of infarction, pulmonary vein injury, nerve damage, embolism and/or atrial esophageal fistula.
- Patients with hemodynamic instability or cardiogenic shock are at increased risk for life-threatening adverse events and ablation must be done with extreme caution.
- The IntellaNav StablePoint Catheter is not intended to be used for internal cardioversion. Doing so may result in perforation, arrhythmias, embolism, thrombus and/or patient death.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- Inspect irrigation saline for air bubbles and remove any air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause embolism.
- If there is uncertainty regarding the patient's anticoagulation status or rhythm prior to the procedure, there should be a low threshold to perform a Transesophageal Echocardiogram (TEE) prior to the procedure to confirm absence of mural thrombus and/or thrombus in the left atrial appendage.
- Guiding catheters and/or long introducer sheaths present the potential for thromboembolic events. Pre-flush and maintain lumen patency with heparinized intravenous infusion.
- Do not wipe this catheter with organic solvents such as alcohol or immerse the handle cable connector in fluids. This may result in electrical or mechanical catheter failures. It may also result in an allergic reaction from the patient.
- Irrigation flow during RF ablation may distort distal tip electrogram recordings due to the signal conductivity of the external cooling solution. Careful monitoring of additional intracardiac electrograms during RF application is recommended to reduce the possibility of inadvertent injury to adjacent structures if appropriate. Higher power coupled with higher flow rates may exacerbate the distortion of the EGM signal recordings.
- Pre-procedural anticoagulation therapy is at the discretion of the physician. However, patients with a history of thromboembolic events may require therapeutic anticoagulation therapy, pre-, during and post-ablation to reduce the incidence of major complications. Peri-procedural anticoagulation therapy is recommended for patients undergoing left-sided and transseptal cardiac procedures and should be considered for selected patients undergoing right-sided procedures.
- The safety and/or efficacy of epicardial use of the IntellaNav StablePoint Catheter has not been evaluated in a clinical trial.
- The Transseptal Puncture (TSP) presents a potential risk for perforation/ tamponade; echocardiography and/or fluoroscopic images should be used to guide the transseptal puncture and a real-time arterial blood pressure monitor should be applied. TSP may induce an air embolus; use proper aspiration and flushing techniques to minimize air embolus.
- Care should be used during multiple sheath/catheter exchanges through the transseptal puncture to avoid causing a residual atrial septal defect that would require repair.
- To avoid patient injury, manipulate the sheath carefully when performing the transseptal puncture especially if the patient has any of the following conditions:
 - Enlarged aortic root
 - Marked right atrial enlargement
 - Small left atrium
 - Marked skeletal deformity or distortion of the thoracic configuration (e.g., scoliosis)

- Any serious incident that has occurred in relation to the device should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS

- The contact force reading is for information only and is not intended to replace standard handling precautions.
- The local impedance reading is for information only and is not intended to replace standard handling precautions.
- The catheter must be warmed up prior to use. If the catheter has not reached a steady state condition prior to use, there is a potential for measurement drift to occur, which could result in an inaccurate force reading.

Note: Refer to the Instructions for Use on compatible Rhythmia Mapping System for instructions on how to perform warm-up.

- Ensure catheter tip is not in contact with myocardial wall or other device when zeroing the contact force reading. Failure to do so could result in inaccurate force reading.

Note: Refer to the Instructions for Use on compatible Rhythmia Mapping System for instructions on how to zero the catheter.

- Always zero the contact force reading following insertion into the patient or moving the catheter from one heart chamber to another.
- To ensure proper function of the contact force reading, ensure the tip electrode and distal ring electrode are outside of a sheath.

Note: Refer to the Instructions for Use on compatible Rhythmia Mapping System for list of sheaths that are compatible with the sheath detect software feature.

- To ensure proper function of the local impedance reading, ensure the tip electrode and two ring electrodes are outside of a sheath.

Note: Refer to the error messages and notifications section of the compatible Rhythmia Mapping System instructions for use for system-related messages and indications related to DIRECTSENSE.

- When applying high force during mapping and RF application, the user should monitor the contact force visualization on the Rhythmia screen to ensure the force does not exceed the operating range.

Note: Refer to the error messages and notifications section of the compatible Rhythmia Mapping System instructions for use for system-related messages and indications related to force range and inaccurate force reading.

- Do not place the distal end of the catheter near magnets. Magnetization of the catheter may result in degradation of magnetic tracking precision. Such degradation may be manifested by an unstable or complete loss of rendering of the position and/or orientation of the catheter by a magnetic tracking system. If this occurs, the catheter should be replaced.
- The IntellaNav StablePoint Catheter is designed for use with a compatible RF Controller, Irrigation Pump and Irrigation Tubing Set that meets the catheter flow rate requirements, a compatible Mapping System, and compatible Connection Box.
- Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the electrode will not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode.
- Use only sterile saline and gauze pad to clean the tip.
- Verify the RF Controller is in the control mode which will deliver the amount of power specified by the power setting unless the measured temperature exceeds the temperature setting. Temperature controlled RF delivery may be affected by the cooling effects of the saline irrigation of the electrode. For example, the Maestro 4000™ RF Cardiac Ablation Controller has these settings in the power control mode.
- Equipment/accessories carrying high frequency alternating current, may cause direct coupled interference and therefore, may disrupt the operation of the RF Controller. It may be necessary to take risk control measures, such as re-orienting, relocating, or shielding the interfering equipment/accessories.
- Use only dispersive pads that meet or exceed IEC 60601-2-2 requirements and follow the dispersive pad manufacturer's instructions for use. The use of dispersive pads which meet ANSI/AAMI requirements (HF18) is recommended.
- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the dispersive pad or failure of an electrical lead.

- The IntellaNav StablePoint Catheter is highly torqueable. Avoid overtorquing. Over-rotating the handle and catheter shaft may cause damage to the distal tip or catheter assembly. Do not rotate the handle and catheter shaft more than 1½ full rotations (540°). If the desired catheter tip position is not achieved, adjust the catheter's curve to disengage the catheter tip from the heart wall before resuming rotation of the handle and catheter shaft.
- Electrophysiology catheters and systems are intended for use only in radiation shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.
- The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
- After use, handle and dispose of product and packaging in accordance with hospital biohazard procedure, administrative and/or local government policy.

ADVERSE EVENTS

Potential adverse events associated with the use of the IntellaNav StablePoint include, but are not limited to:

- Pain or discomfort, for example:
 - Angina
 - Chest pain
 - Non-cardiovascular pain
- Cardiac arrest
- Death
- Hypertension
- Hypotension
- Infection/inflammation/exposure to biohazardous material
- Edema/heart failure/pleural effusion
- Procedural related side effects, for example:
 - Allergic reaction (including anaphylaxis)
 - Genitourinary complication
 - Side effects related to medication or anesthesia
 - Radiation injury/tissue burn
 - Renal failure/insufficiency
 - Vasovagal response
- Respiratory distress/insufficiency/failure/dyspnea
- Arrhythmia (new or exacerbated)
- Conduction pathway injury (heart block, nodal injury, etc.)
- Nerve injury, for example:
 - Phrenic nerve injury
 - Vagal nerve injury
- Gastrointestinal disorders
- Vessel trauma, including:
 - Perforation
 - Dissection
 - Coronary artery injury
 - Vasospasm
 - Occlusion
 - Hemothorax
- Cardiac trauma, for example:
 - Cardiac perforation/cardiac tamponade/pericardial effusion
 - Valvular damage
 - Stiff left atrial syndrome
- Injury related to tissue damage and/or adjacent structures, for example:
 - Esophageal injury
 - Pulmonary injury
 - Catheter entrapment
 - Physical trauma
- Fistula, for example:
 - Atrio-esophageal fistula
 - Bronchopercardial fistula
- PV stenosis and its symptoms, for example:
 - Cough
 - Shortness of breath, fatigue
 - Hemoptysis

- Surgical and access complications, for example:
 - Hematoma/seroma
 - AV fistula
 - Bleeding
 - Pseudoaneurysm
 - Pneumothorax
 - Residual atrial septal defect
 - Thrombosis
- Injury due to embolism/thromboembolism/air embolism/foreign body embolism
 - Cerebrovascular Accident (CVA)/stroke
 - Transient Ischemia Attack (TIA)
 - Myocardial infarction
- Neurological impairment and its symptoms, for example:
 - Cognitive changes, visual disturbances, headache, motor impairment, sensory impairment, and speech impairment
- Pulmonary embolism
- Asymptomatic cerebral embolism

The potential adverse events may be related to the ablation catheter(s) and/or the interventional procedure. The severity and/or the frequency of these potential adverse events may vary and may result in prolonged procedure time and/or additional medical and/or surgical intervention, implantation of a permanent device such as a pacemaker, and in rare cases, may result in death.

CLINICAL STUDIES

The IntellaNav StablePoint Catheter was developed based on the Blazer™ OI Ablation Catheter platform. The IntellaNav StablePoint Catheter has not been specifically evaluated in a clinical trial. The inclusion of the magnetic sensor, force sensing, local impedance, and associated functionality does not affect the ablation therapy delivery or catheter operation because the use of fluoroscopy, or other visualization techniques such as echocardiography, is still required to confirm device positioning prior to delivery of therapy. Clinical data on the IntellaNav StablePoint Catheter will be available upon receipt.

HOW SUPPLIED

Device Details

The IntellaNav StablePoint Catheter is supplied sterile using an Ethylene Oxide (EO) process. Peel-off labels for device and accessories can be used for device traceability. In addition to the IntellaNav StablePoint Catheter, please refer to the Additional Items for Safe Use section for a detailed list of other materials typically required in an Electrophysiology (EP) procedure.

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.
- Do not use the device if past the "Use By" date.

Handling and Storage

Operating Environment

- Ambient Temperature: 10 °C to 40 °C
- Relative Humidity: 30 % to 75 %
- Atmospheric Pressure: 70 kPa to 106 kPa

Transport Environment

- Temperature: -29 °C to 60 °C
- Relative Humidity: 30 % to 85 %
- Atmospheric Pressure: Uncontrolled

Storage Environment

- Temperature: 15 °C to 25 °C
- Relative Humidity: Uncontrolled
- Atmospheric Pressure: Uncontrolled

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

Intracardiac electrophysiology and cardiac ablation procedures should be performed in a fully equipped clinical electrophysiology lab.

In addition to the IntellaNav StablePoint Catheter, the following materials, devices and equipment will be required:

- Compatible Connection Box
- Compatible Catheter Cable
- Compatible RF Controller and accessories
- Irrigation Pump and accessories
- Irrigation Tubing Set

Note: The RF Controller and pump/tubing system must be capable of delivering flow rates of 2 mL/min (Standby), 17 mL/min (Low Ablation Flow Rate at ≤ 30 W), 30 mL/min (High Ablation Flow Rate at 31 W to 50 W).

Accessories

- Commercially available disposable Dispersive Pads that meet or exceed IEC 60601-1/ IEC 60601-2-2 requirements

Optional Additional Equipment:

- Compatible Mapping System and accessories

Preparation

Caution: Before use, inspect the packaging for any violation of the sterile barrier and inspect the IntellaNav StablePoint Catheter for any defects. Do not use potentially contaminated or defective equipment.

Please refer to the operator's manuals and Instructions for Use for the Irrigation Pump, RF Controller, Mapping System, Connection Box, and the Irrigation Tubing Set for instructions on connecting and operating these systems in conjunction with the IntellaNav StablePoint Catheter. Use the appropriate accessory cables to connect the IntellaNav StablePoint Catheter to accessory equipment.

- Attach the Dispersive Pad to the patient and RF Controller.
- Attach the Location Reference Patch Kit to the patient per the patch IFU.
- Connect the patient to an ECG recording system to facilitate arrhythmia monitoring per the standard operating procedure of the electrophysiology lab or manufacturer's operator's manual.

Note: This should be done prior to introducing any intracardiac catheters.

- Open the IntellaNav StablePoint Catheter, IntellaNav™ Cable packages, and the Irrigation Tubing Set package. Carefully transfer the package contents into the sterile field, maintaining sterile technique.
- Obtain vascular access in a large central vessel (e.g. a femoral vein or artery) under aseptic conditions. Then place an introducer sheath into the vessel using a standard percutaneous technique.
- Connect the Connection Box to the RF Controller (and mapping system if desired) according to the operator's manuals, and/or IFUs.
- Connect the RF Controller to the recording system (and the Mapping System if desired) with the appropriate interface cables according to the operator's manuals, and/or IFUs.
- Connect the IntellaNav StablePoint Catheter to the Connection Box using the compatible Ablation Catheter Cable. Ensure that the cable/catheter connection remains dry throughout the procedure. For connection information, refer to the IFU for additional connection instructions.
- Turn on the power to the RF Controller.
- Set the RF Controller to Power Control mode. The RF Controller's temperature limit must not exceed 50 °C, but can be set lower at physician discretion.
- Turn on the Irrigation Pump.
- Make sure that the Irrigation Pump has the following flow rates: 2mL/min (Standby), 17 mL/min (Low Ablation Flow – 30W or less), 30 mL/min (High Ablation Flow – above 30W). Refer to the Irrigation Pump operator's manual for instructions on how to adjust the pump settings if required.
- Refer to either the Irrigation Tubing Set or the Irrigation Pump IFU for instructions to connect the Irrigation Tubing Set to irrigation fluid and install into the Irrigation Pump.
- Connect the IntellaNav StablePoint Catheter to the Irrigation Tubing Set via the luer fitting at the proximal end of the catheter handle. Care must be taken to ensure all luer fittings are secure to prevent leaking.
- Purge the IntellaNav StablePoint Catheter and Irrigation Tubing Set. Fluid should exit all six (6) irrigation ports during the flushing process. Assure that no air remains within the Irrigation Tubing Set or lumen and all irrigations ports are patent. Pre-condition the tip electrode if needed using programmed stimulation.
- Set the pre-RF delay and post-RF delay on the Irrigation Pump (default is 2 seconds). Reference the Irrigation Pump Operator's Manual for instructions on how to change the pre-RF delay and post-RF delay (default is 2 seconds).
- Check the catheter steering by articulating the steering knob prior to inserting the catheter in the sheath.
- Before placing the IntellaNav StablePoint Catheter in the sheath, begin continuous irrigation at a flow rate of 2 ml/min, i.e. standby flow. Check for any leaks at the tip of the IntellaNav StablePoint Catheter (other than normal saline flowing out of the distal ports), at the IntellaNav StablePoint Catheter handle, and at the luer connections and tubing joints.
- To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion. Any sheath <8.5F is contraindicated.

Warning: The use of excessive force to advance or withdraw the catheter may lead to catheter damage that may result in inaccurate force readings.

Procedure

- Under standard imaging guidance (e.g. fluoroscopy or echocardiography), insert the IntellaNav StablePoint Catheter into the sheath and advance through the vasculature into the heart.

Note: The degree of tip deflection of the IntellaNav StablePoint Catheter is controlled by the Steering Knob on the IntellaNav StablePoint Catheter handle (See Figure 1). If the Steering Knob is turned in a clockwise direction from its neutral position, the tip will curve proportionately in one direction depending upon the curve option selected. Turning the Steering Knob in the counter-clockwise direction will cause the tip to deflect in the opposite direction. To prevent overstressing the tip, the Steering Knob movement is limited by the handle design. The tension adjust knob may be used when the desired catheter placement is achieved.

- In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up after inserting into the patient and prior to the use of the force feature. Refer to instructions for use on compatible Rhythmia Mapping System for details on how to initiate warm-up.
- If using the force feature, zero the contact force reading when the catheter is not in contact with the myocardial wall or other intracardiac devices. Ensure the tip and distal ring are outside the sheath prior to zero-ing. Ensure the catheter tip is not in contact by evaluating the location using standard imaging techniques and on Rhythmia. Refer to instructions for use on compatible Rhythmia Mapping System for details on how to re-zero catheter.
- If using the force feature, confirm the non-contact state and re-zero as needed throughout procedure. Refer to Mapping System user manual for error messaging related to force accuracy.
- If using the force feature, re-zero the force reading when moving the catheter from one chamber to another and upon re-insertion. Warm up is only needed upon first insertion to patient.
- Determine the area of interest for ablation. Use both standard imaging techniques and intracardiac electrograms to aid in proper positioning. Ensure a stable position is achieved before delivering RF.
Note: If using the force and/or local impedance feature(s) to aid in positioning, be sure to use in the context of all available information including intracardiac electrograms, electroanatomic maps, generator information, fluoroscopy, etc.
Note: To ensure proper function of the force feature, ensure the tip electrode and distal ring electrode are outside of a sheath, if using a long sheath.
Note: To ensure correct use of electrograms, ensure the tip electrode and all ring electrodes are outside of a sheath, if using a long sheath.
Note: To ensure proper function of the local impedance feature, ensure the tip electrode and distal two ring electrode are outside of a sheath, if using a long sheath.
- Set the initial power and irrigation flow rate based on the recommended ablation parameters in the table below.

RF Power	Minimum Irrigation Flow Rate
≤30W	17 mL/min
31-50W	30 mL/min

Power levels exceeding 30W may be used when transmural lesions cannot be achieved at lower energy levels but should be accompanied by an increase in irrigation flow rate to 30 mL/min.

Note: Confirm the increased irrigation flow rate prior to onset of RF energy by observation of a decrease in tip electrode temperature of at least 2°C (default).

Warning: Using the IntellaNav StablePoint Catheter at lower than prescribed flow rates may increase the potential for thrombus, coagulum, and char that may result in embolism.

Warning: If using the force feature during RF delivery, contact forces greater than 70g may increase the incidence of steam pop occurrence, especially when used at RF powers greater than 30W.

Note: In some orientations a force >50grams is only displayed as "High". Refer to Instructions for Use on compatible Rhythmia Mapping System for details.

- Start RF energy delivery. Monitor the catheter tip temperature, generator impedance, fluoroscopy, catheter visualization, positional stability, intracardiac electrograms, force and local impedance changes to determine appropriate RF duration.

Note: Tip temperature should be used only as an indicator of adequate irrigation and not as a measure of tissue temperature.

Note: Too rapid an increase in power during ablation, ablating at high power (>30W) or insufficient flow rate may lead to perforation caused by steam pop, arrhythmias, damage to adjacent structures, and/or embolism.

Note: Proximity to the esophagus should be considered when ablating on the posterior atrial wall. Using high power, high contact, or long durations may increase the risk of esophageal injury.

Warning: If using local impedance feature during RF delivery, a decrease in impedance greater than 65Ω from the pre-ablation baseline may not improve characteristics of lesion formation and may increase the risk of steam pop or perforation.

Note: If using the local impedance feature and a sudden rise in local impedance is noted during RF delivery, manually discontinue RF delivery, as this may be an indication of char or coagulum formation.

Note: If using the local impedance feature during RF delivery, a point by point workflow is recommended as the calculation of the impedance drop does not reset with changes in positional stability.

28. Do not ablate for greater than 60 seconds in duration without moving the tip of the IntellaNav StablePoint Catheter.

29. RF energy may be reapplied to the same or alternate sites using the same catheter.

30. To confirm effectiveness of ablation, user shall use functional endpoints such as arrhythmia termination or demonstrated conduction block.

31. Prior to removing the IntellaNav StablePoint Catheter, completely straighten the distal end of the catheter.

32. Withdraw the IntellaNav StablePoint Catheter when the procedure is finished.

33. Turn off the RF Controller and Irrigation Pump.

Disposal

To minimize risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, the catheter may contain biohazardous substances. The catheter and packaging should be treated and disposed of as biohazardous waste in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

In the event that a serious incident occurred in relation to the device, including all patient deaths for procedures where the BSC product was used, the event should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established. Return any catheter related to a complaint, patient harm, injury, or death to Boston Scientific using a BSC Returned Product Kit.

Returning products for analysis and providing product performance observations helps drive higher reliability on an ongoing basis.

Be sure to follow the instructions with regard to packaging and shipping any biohazard devices taking care not to expose the handle or connector to fluid that may compromise the catheter and limit analysis.

Post-Procedure

Carefully monitor patient while in recovery to ensure hemostasis is achieved and any complications are immediately treated.

Note: Any serious incident that has occurred in relation to the device should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of the IntellaNav StablePoint Catheter in association with the Electrophysiological cardiac interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document.
- Discuss post procedure instructions, including any lifestyle changes, medications, when to call the Healthcare Provider (HCP) and any post procedure follow-up that might be needed.

TROUBLESHOOTING

See instructions for Use for compatible Rhythmia Mapping System for troubleshooting related to Force and DIRECTSENSE Error Messaging. If force inaccuracy is suspected, proceed through the following corrective actions: 1) Check the position of the sheath to ensure the tip and distal ring electrode are fully outside the sheath, 2) Reassess the unloaded state in the blood pool and re-zero if necessary, 3) Verify the connections between the catheter, cable, and connection box, or 4) Replace the IntellaNav StablePoint Catheter.

Problems	Possible Cause	Corrective Action Procedure
Temperature not displayed	Poor catheter/cable connections	1. Verify that the Cable is plugged into the Connection Box, and the IntellaNav StablePoint Catheter. 2. Verify that the Connection Box is connected to the RF Controller. 3. Replace cable and/or catheter. 4. If the RF Controller still does not display temperature, there may be a malfunction in the temperature sensing system. 5. Consult the operator's manual and correct this malfunction prior to reapplying RF energy.
• Impedance cutoff • Temperature cutoff	Char/coagulum on tip electrode	1. Discontinue RF delivery. 2. Straighten the distal end and withdraw the IntellaNav StablePoint Catheter. 3. Inspect tip electrode for any char/coagulum. 4. If present, gently wipe the tip section with a sterile gauze dampened with sterile saline (do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode). 5. Prior to reinsertion, ensure the irrigation ports are patent. If irrigation port occlusion occurs: a. Ensure IntellaNav StablePoint Catheter is removed from the patient. b. Fill a 1 mL or 2 mL syringe with sterile saline and attach to the stop-cock sidearm of the IntellaNav StablePoint Catheter. c. Carefully inject the saline from the syringe into the IntellaNav StablePoint Catheter. Fluid should exit all six (6) irrigation ports during the flushing process. d. Repeat steps b and c, if necessary. e. If the irrigation ports are cleared, the IntellaNav StablePoint Catheter can be reintroduced into the patient. WARNING: Do not continue using the IntellaNav StablePoint Catheter if the irrigation ports are occluded or the catheter is not functioning properly.
Suspected failure of fluid flow integrity	• Leak in catheter and/or Irrigation Tubing Set • Irrigation Pump out of calibration	1. Discontinue RF delivery. 2. Straighten the distal end and withdraw catheter. 3. Replace IntellaNav StablePoint Catheter and Irrigation tubing set and prime outside of the patient. 4. Replace IntellaNav StablePoint Catheter and/or Irrigation tubing set if parameters do not appear normal or if there is any abnormality of the integrity of fluid flow. 5. Refer to the Irrigation Pump operator's manual to verify fluid flow is accurate. 6. Contact BSC representative to replace Irrigation pump.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

IntellaNav StablePoint, Rhythmia, DIRECTSENSE, Maestro 3000, Maestro 4000, Blazer, and IntellaNav are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

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All rights reserved.
© 2023 Boston Scientific Corporation or its affiliates.