VIKING™
Fixed Curve Diagnostic Catheter

Instructions for Use	2
Instrucciones de uso	10
Mode d'emploi	18
Gebrauchsanweisung	26
Istruzioni per l'uso	35
Instructies voor gebruik	43
Instruções de Utilização	51

51494765-01

2024-12

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Table 1: Cable Selection3

 Contents.....3

 Non-pyrogenic.....3

 User Information.....3

INTENDED USE4

INDICATIONS FOR USE.....4

 Clinical Benefit Statement.....4

 Summary of Safety and Clinical Performance4

CONTRAINDICATIONS.....4

WARNINGS.....4

PRECAUTIONS6

ADVERSE EVENTS7

HOW SUPPLIED.....8

 Device Details.....8

 Handling and Storage.....8

OPERATIONAL INSTRUCTIONS8

 Preparation8

 Procedure.....8

 Disposal8

 Post-Procedure8

PATIENT COUNSELING INFORMATION9

WARRANTY9

SYMBOL DEFINITIONS.....9

VIKING™

Fixed Curve Diagnostic Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The VIKING Fixed Curve Diagnostic Catheter has an insulated polymer shaft with platinum electrodes located along the distal section of the shaft.

Use appropriate accessory cables to connect the Diagnostic Catheter to the appropriate accessory equipment. Refer to the cable Instructions for Use (IFU) for connection instructions.

Table 1: Cable Selection

Catheter	Cable to Use	Cable Length
VIKING	4-Pin SureLink™ Cable	125 cm
VIKING	4-Pin SureLink Cable (Blue)	125 cm
VIKING	4-Pin SureLink Cable (Red)	125 cm
VIKING	4-Pin SureLink Cable (Yellow)	125 cm
VIKING	4-Pin SureLink Cable	210 cm
VIKING	10-Pin SureLink Cable	125 cm
VIKING	10-Pin SureLink Cable	210 cm

Contents

One (1) VIKING Fixed Curve Diagnostic Catheter

Non-pyrogenic

This device meets pyrogen limit specifications for all patient-contacting parts.

User Information

Cardiac mapping procedures should be performed only by physicians thoroughly trained in the techniques of intracardiac electrophysiology, catheter mapping and pacing, and in the specific approach to be used, in a fully-equipped electrophysiology lab.

INTENDED USE

The VIKING Fixed Curve Diagnostic Catheter is intended for temporary intracardiac sensing, recording, stimulation and temporary pacing during the evaluation of cardiac arrhythmias. The device is intended for use in adult (not pediatric) patients, with the exclusion of pregnant and/or nursing patients.

INDICATIONS FOR USE

The VIKING Fixed Curve Diagnostic Catheter is indicated for use to diagnose cardiac arrhythmias.

Clinical Benefit Statement

The VIKING Fixed Curve Diagnostic Catheter is a basic and essential tool for EP procedures. When operated by the intended user together with the EP recording system, the VIKING Fixed Curve Catheter enables navigation through the cardiovascular anatomy and respective electrical mapping in minimally invasive procedures and assist in identifying the problematic area within the heart chamber. Therefore, the clinical benefit of the Fixed Curve Diagnostic Catheter is to enable the acquisition of electrocardiographic signals and delivery of pacing current for the diagnosis and treatment of health impacting and/or life-threatening cardiac arrhythmias. Indeed, failure to treat a cardiac arrhythmia could result in symptoms such as shortness of breath, palpitations, dizziness, syncope, chest pain or sudden cardiac death.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

The VIKING Fixed Curve Diagnostic Catheter is contraindicated for use:

- in patients with a mechanical prosthetic heart valve through which the catheter must pass;
- in patients with active systemic infection;
- in patients with conditions where insertion into or manipulation in the cardiac chambers is unsafe as this may increase the risk of perforation or systemic embolic event, such as but not limited to, evidence of myxoma, Left Atrial (LA) thrombus, intracardiac mural thrombus, or interatrial baffle or patch;
- from the femoral approach in patients who have vena cava embolic protection filter devices or known femoral thrombus;
- in patients who are hemodynamically unstable;
- in patients with a contraindication to an invasive electrophysiology procedure where insertion or manipulation of a catheter in the cardiac chambers is deemed unsafe, such as but not limited to, a recent previous cardiac surgery (e.g., ventriculotomy or atriotomy, Coronary Artery Bypass Graft [CABG], PTCA/PCI/coronary stent procedure/unstable angina) and/or in patients with congenital heart disease where the underlying abnormality increases the risk of the ablation (e.g. severe rotational anomalies of the heart or great vessels).

WARNINGS

- Cardiac catheterization procedures should be performed only by physicians thoroughly trained in invasive cardiology and in the specific approach to be used, in a fully-equipped electrophysiology lab.
- All equipment used in connection with the BSC Catheter must meet IEC 60601-1 requirements, including those for use with Type CF Defibrillation-Proof applied parts, as well as any other relevant regulatory and safety standards. The use of catheter with devices that do not comply with relevant standards may cause equipment damage or system malfunction, or harm to the patient or user.

- Before using, inspect the VIKING Fixed Curve Diagnostic Catheter for any defects or physical damage, including electrical insulation on the cables and the catheter shaft that, if used, may cause patient and/or user injury. Do not use defective or damaged devices. Replace damaged equipment if necessary. No modification of this equipment is allowed.
- The use of catheters or cables with unprotected male pin connectors presents a risk of electrical hazard. Inadvertent attachment of pin connectors to power supply sockets of connectors could result in electrocution of the patient or operator.
- Do not use the VIKING Fixed Curve Diagnostic Catheter as an internal defibrillation catheter. Doing so may result in perforation, arrhythmias, embolism, thrombus, and/or patient death.
- Cardiac catheterization procedures present the potential for significant radiation exposure, which can result in acute radiation injury as well as an increased risk for somatic and genetic effects, to both patients and laboratory staff due to the radiation beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Due to radiation exposure during cardiac catheterization, the safety and effectiveness of this device has not yet been established in pregnant and/or nursing women and pediatric patients.
- Careful catheter manipulation must be performed to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be performed under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Administer appropriate levels of peri-procedural anticoagulation therapy for selected patients undergoing right-sided procedures. Administer anticoagulation therapy during and post-procedure according to institution standards to minimize bleeding and thrombotic complications.
- Stimulation of cardiac tissues caused by pacing stimuli may lead to inadvertent induction of arrhythmias. These arrhythmias may require defibrillation that could also result in skin burns.
- Do not use if package is opened or damaged.
- Carefully read all equipment and ancillary device instructions required for the procedure prior to use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to do so may result in patient complications.
- This catheter is not indicated for use in Cardiac Ablation or Coronary Artery Mapping.
- Ensure that the cable/catheter connection remains dry throughout the procedure in order to prevent electric shock or other patient injuries as well as to prevent loss of device function.
- Fibrin may accumulate in or on the sheath/catheter assembly during the procedure. Aspirate when removing the dilator or catheter.
- In the presence of anticoagulation, there may be an increased risk of bleeding from all causes.
- Warnings for patients with implantable pacemakers and implantable cardioverter/ defibrillators (ICDs):
 - Have temporary external sources of pacing and defibrillation available.
 - Fluoroscopic or appropriate imaging guidance and care must be taken during catheter advancement, manipulation, and withdrawal to avoid lead dislodgement.
 - Monitor pre- and post-measurements for sensing and

pacing thresholds and impedances to determine the integrity of the lead-patient function.

- Temporarily Deactivate ICDs for the duration of the procedure as they could discharge and injure the patient or be damaged by the procedure.
- Perform a complete analysis of the implanted device function after any electrophysiology procedure.
- Remember to reactivate the pulse generator after the procedure.
- At no time should the VIKING Fixed Curve Diagnostic Catheter be advanced or withdrawn when resistance is felt without determining the cause. Valve damage, vascular and/or cardiac perforation is a risk with any intracardiac catheter.
- Catheter entrapment within the heart or blood vessels is a possible complication of cardiac catheterization procedures. The potential for catheter entrapment may be increased when the catheter is over torqued and/or positioned in the chordae tendineae. The occurrence of this complication may necessitate surgical intervention and/or repair of injured tissue and/or valve damage.
- There are no data to support the safety and effectiveness of this device in the pediatric population.
- Excessive bending or kinking of the catheter shaft may damage internal wires. Manual prebending of the distal curve can damage the electrical wires and may cause patient injury.
- Use both fluoroscopy or other visualization technique, such as echocardiography, and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid conduction pathway injury, cardiac perforation or tamponade.
- Guiding catheters and/or long introducer sheaths present the potential for thromboembolic events. Pre-flush and maintain lumen patency with heparinized intravenous infusion.
- Do not wipe this catheter with organic solvents such as alcohol or immerse the handle cable connector in fluids. This may result in electrical or mechanical catheter failures. It may also result in an allergic reaction from the patient.
- This catheter has not been validated in an MRI environment. Do not use the VIKING Fixed Curve Diagnostic Catheter near MRI equipment. Doing so may result in movement or heating of the catheter, which may cause distortion/loss of the catheter image on electroanatomic mapping. This could lead to patient injuries such as cardiac trauma.
- Portable radiofrequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the VIKING Fixed Curve Diagnostic Catheter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Any serious incident that has occurred in relation to the device should be reported to BSC and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS

- Before using, check shelf life. Do not use catheter after expiration date.
- Before using, inspect for physical damage, including electrical insulation on the cables and the catheter shaft. Replace damaged equipment.
- Do not use the VIKING Fixed Curve Diagnostic Catheter closer than 30 cm (12 inches) to any Wireless Power Transfer (WPT) and 5G cellular devices, otherwise electromagnetic interference from those devices could result in degradation of the performance of this equipment.

ADVERSE EVENTS

Potential adverse events associated with the use of the VIKING Fixed Curve Diagnostic Catheter include, but are not limited to:

- Pain or discomfort, for example:
 - Angina
 - Chest pain
 - Non-cardiovascular pain
- Cardiac arrest
- Death
- Hypertension
- Hypotension
- Infection/inflammation/exposure to biohazardous material
- Edema/heart failure/pleural effusion
- Procedural related side effects, for example:
 - Allergic reaction (including anaphylaxis)
 - Genitourinary complication
 - Side effects related to medication or anesthesia
 - Radiation injury
 - Renal failure/insufficiency
 - Gastrointestinal disorders
 - Vasovagal response
- Respiratory distress/insufficiency/failure/dyspnea
- Arrhythmia (new or exacerbated)
 - Conduction pathway injury (heart block, nodal injury, etc.)
- Vessel trauma, including:
 - Perforation
 - Dissection
 - Coronary artery injury
 - Vasospasm
 - Occlusion
 - Hemothorax
- Cardiac trauma, for example:
 - Cardiac perforation/cardiac tamponade/pericardial effusion
 - Valvular damage
- Injury related to tissue damage and/or adjacent structures, for example:
 - Catheter entrapment
 - Physical trauma
- Surgical and access complications, for example:
 - Hematoma/seroma
 - AV fistula
 - Bleeding
 - Pseudoaneurysm
 - Pneumothorax
- Thrombosis
- Injury due to embolism/thromboembolism/air embolism/foreign body embolism:
 - Cerebrovascular Accident (CVA)/stroke
 - Transient Ischemic Attack (TIA)
 - Myocardial infarction
 - Pulmonary embolism
 - Asymptomatic cerebral embolism
 - Neurological impairment and its symptoms, for example:
 - Cognitive changes, visual disturbances, headache, motor impairment, sensory impairment, and speech impairment

The potential adverse events may be related to the catheter(s) and/

or the interventional procedure. The severity and/or the frequency of these potential adverse events may vary and may result in prolonged procedure time and/or additional medical and/or surgical intervention, implantation of a permanent device such as a pacemaker, and in rare cases, may result in death.

HOW SUPPLIED

Device Details

The VIKING Fixed Curve Diagnostic Catheter is supplied sterile using an Ethylene Oxide (EO) process and intended for single use only.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Catheter does not require any special storage and/or handling conditions or special operating conditions that are beyond what would normally be expected by the end user.

OPERATIONAL INSTRUCTIONS

Preparation

1. Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, do not use and contact your Boston Scientific representative.
2. Open the catheter package and carefully transfer the package contents onto the sterile field, maintaining sterile technique. If the catheter is packaged with a curve retainer, remove the retainer.
3. Before using, inspect for physical damage, including electrical insulation on the cables and the catheter shaft. Replace damaged equipment.

Procedure

The following are general guidelines for use:

1. Insert the catheter by using standard percutaneous catheter introducers. The electrode catheter should be maneuvered under fluoroscopic guidance.
2. Connect the catheter to the appropriate Boston Scientific cable according to Table 1 (see the cable instructions for further details).
3. For stimulation and recording electrograms, connect the lead pins to an amplifier. For temporary pacing, connect the lead pins to an external pulse generator.

Disposal

To minimize risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, the catheter may contain biohazardous substances. The catheter and packaging should be treated and disposed of as biohazardous waste in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

In the event that a serious incident occurred in relation to the device, including all patient deaths for procedures where the BSC product was used, the event should be reported to BSC and the competent authority of the Member State or other relevant local authority in

which the user and/or patient is established. Return any catheter related to a complaint, patient harm, injury, or death to Boston Scientific using a BSC Returned Product Kit.

Returning products for analysis and providing product performance observations helps drive reliability higher on an ongoing basis.

Be sure to follow the instructions with regard to packaging and shipping any biohazard devices taking care not to expose the handle or connector to fluid that may compromise the catheter and limit analysis.

PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling patients on the use of the VIKING Catheter in association with the Electrophysiological cardiac interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document.
- Discuss post procedure instructions, including any lifestyle changes, medications, when to call the Healthcare Provider (HCP) and any post procedure follow-up that might be needed.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

The following are trademarks of Boston Scientific Corporation or its affiliates: VIKING, SureLink.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
©2024 Boston Scientific Corporation