Guidezilla™ II

Guide Extension Catheter

Instructions for Use	2
Instrucciones de uso	5
Mode d'emploi	8
Gebrauchsanweisung	11
Istruzioni per l'uso	14
Instructies voor gebruik	17
Instruções de Utilização	20

Guidezilla™ II

Guide Extension Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. There has been no testing to validate this device beyond a single use. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Guidezilla II Guide Extension Catheter is a single-lumen rapid exchange catheter designed to act as an extension to a traditional guide catheter and facilitate the delivery of interventional devices into the vasculature. The 150 cm device consists of a hydrophilic-coated, single-lumen distal guide segment connected to a stainless steel proximal hypotube. The distal guide segment contains a platinum-iridium radiopaque marker band located 2 mm from the distal tip and a proximal platinum-iridium radiopaque collar. The proximal hypotube contains two positioning marks located at 90 cm and 100 cm from the distal tip. A tab on the proximal end of the hypotube is used for device identification.

Table 1. Compatibility Information

Model	Distal Guide Segment Length	Compatible Guide Catheter	Guidezilla II ID
6F	25 cm	≥ 6F / ≥ 0.070 in ID (1.78 mm)	0.057 in (1.45 mm)
7F	25 cm	≥ 7F / ≥ 0.078 in ID (1.98 mm)	0.063 in (1.60 mm)
8F	25 cm	≥ 8F / ≥ 0.088 in ID (2.24 mm)	0.072 in (1.83 mm)
6F LONG	40 cm	≥ 6F / ≥ 0.070 in ID (1.78 mm)	0.057 in (1.45 mm)

Contents

Qty **Material**

One (1) Single-use Guidezilla II Guide Extension Catheter

Operating Principle

Guidezilla II Guide Extension Catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary vasculature, and to facilitate placement of interventional devices. The Guidezilla II Guide Extension Catheter facilitates back-up support and device delivery, extending the scope of coronary interventions.

Materials

The Guidezilla II Guide Extension Catheter consists of an ABS handle, a stainless steel shaft, a PTFE lined stainless steel reinforced elastomer catheter with a platinum-iridium radiopaque marker and a hydrophilic coating.

Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1 % weight by weight.

Note: This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

Non-pyrogenic

The Guidezilla II Guide Extension Catheter is non-pyrogenic.

User Information

Physicians thoroughly trained in interventional cardiology procedures.

INTENDED USE / INDICATIONS FOR USE

The Guidezilla II Guide Extension Catheter is intended to be used in conjunction with guide catheters to access discrete regions of the coronary vasculature, and to facilitate placement of interventional devices.

Clinical Benefit Statement

The clinical benefit of the Guidezilla II Guide Extension Catheter is to facilitate improvement of myocardial perfusion in patients who are candidates for PCI procedures by providing access to discrete regions of the coronary vasculature.

Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

- Vessels less than 2.5 mm in diameter.
- Vessels in the neurovasculature and the venous system.

WARNINGS

- Use prior to the "use by" date as indicated on the label.
- Never advance the Guidezilla II Guide Extension Catheter into a vessel without a leading guidewire or without confirming location using fluoroscopic guidance. Vessel dissection or perforation may result.
- Due to the size and non-tapered tip of the Guidezilla II Guide Extension Catheter, extreme care must be taken to avoid vessel ischemia or vascular damage.
- Where there is limited clearance between interventional devices and the distal guide segment lumen, those devices must be advanced and withdrawn slowly with the hemostasis valve open to reduce the risk of embolism.
- This is a non-torqueable device. Torqueing the device may result in wire wrap or damage to the device or vessel.

- Never advance the Guidezilla II Guide Extension Catheter into a vessel with an effective diameter less than 2.5 mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guidezilla II Guide Extension Catheter, withdraw the Guidezilla II Guide Extension Catheter until the pressure returns to normal.

PRECAUTIONS

- Inspect the device prior to use for any bends or kinks. Any device damage may decrease the desired performance characteristics.
- The device lumen should be thoroughly flushed with heparinized saline prior to use.
- This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.
- Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.
- Exercise care in handling of the device during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- When the device is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the device without observing the resultant tip response.
- Never advance the Guidezilla II Guide Extension Catheter more than 15 cm beyond the tip of the guide catheter. Further distal advancement of the Guidezilla II Guide Extension Catheter could cause the entire guide segment to track outside of the guide catheter and impede withdrawal of the device.
- If strong resistance is encountered during manipulation of the devices, do not force passage. Determine the cause of the resistance before proceeding. If the cause cannot be removed, withdraw all devices simultaneously.

ADVERSE EVENTS

The adverse events include, but are not limited to:

- Abrupt closure, slow flow/no reflow
- Allergic reaction
- Arrhythmia
- Arterial spasm
- Bleeding, hemorrhage or hematoma
- Death
- Embolism (air, tissue, device fragments, plaque)
- Hypotension/Hypertension
- Infection
- Myocardial ischemia/infarction
- Organ insufficiency/failure (heart, lung, kidney)
- Prolonged procedure time including additional intervention or surgery
- Radiation Injury
- Thrombosis
- Vessel trauma (e.g., perforation, dissection, spasm, intimal disruption, pseudoaneurysm, AV fistula)

HOW SUPPLIED

Device Details

Contents supplied STERILE using an ethylene oxide (EO) process.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

Other items required but not provided:

- Guide catheter with an inner diameter large enough to accommodate the Guidezilla II Guide Extension Catheter
- Y-adaptor with hemostasis valve
- Guidewire
- Sterile syringe (for system flushing)
- Sterile heparinized saline (for system flushing)

Preparation

1. Carefully inspect the Guidezilla II Guide Extension Catheter prior to use to verify that it has not been damaged during handling.
2. Open the pouch using standard sterile handling procedures and place the packaging coil into the sterile field.
3. Carefully remove the Guidezilla II Guide Extension Catheter from its protective tubing. Do not bend or kink the device during removal.
4. Immerse the distal guide segment in heparinized saline solution.

Procedure

Instructions for Use

Deliver the Guidezilla II Guide Extension Catheter according to the following steps:

1. Backload the Guidezilla II Guide Extension Catheter onto the proximal portion of the guidewire while maintaining guidewire position.
2. Fully open the hemostatic valve to allow for passage of the Guidezilla II Guide Extension Catheter. Advance the Guidezilla II Guide Extension Catheter into the hemostatic valve attached to a guide catheter. Advance the Guidezilla II Guide Extension Catheter through the guide catheter.
3. Advance the Guidezilla II Guide Extension Catheter over the guidewire under direct fluoroscopic visualization and position the tip of the Guidezilla II Guide Extension Catheter no more than 15 cm distal to the guide catheter tip. Use the radiopaque marker band and collar as reference points. If the position of the Guidezilla II Guide Extension Catheter is not optimal, it should be carefully repositioned or removed.

Warning: Due to the size and non-tapered tip of the Guidezilla II Guide Extension Catheter, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Warning: Never advance the Guidezilla II Guide Extension Catheter into a vessel with an effective diameter less than 2.5 mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guidezilla II Guide Extension Catheter, withdraw the Guidezilla II Guide Extension Catheter until the pressure returns to normal.

4. Insert the interventional device into the guide catheter and Guidezilla II Guide Extension Catheter through the Y-adapter valve.

Note: Use caution when advancing the interventional device into the distal guide segment.

5. Adequately tighten the hemostatic valve. Perform the interventional procedure according to the directions provided by the manufacturer of the interventional device.
6. Carefully remove all interventional devices. Carefully retract the Guidezilla II Guide Extension Catheter up to the hemostatic valve. Remove the Guidezilla II Guide Extension Catheter from the patient, ensuring that the hemostatic valve is fully open prior to Guidezilla II Guide Extension Catheter removal.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, device may contain biohazardous substances. Devices that contain biohazardous substances are considered biohazardous waste and should be stored in a biohazard container that is labeled with the biological hazard symbol. Untreated biohazardous waste should not be disposed of in the municipal waste system. Biohazardous substances should be removed via incineration prior to disposal. Alternatively, biohazardous waste may be disposed of utilizing a certified facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

INFORMATION TO BRIEF THE PATIENT

The following sections/clauses in the IFU contain information that may be useful for briefing the patient prior to use of this device:

- Clinical Benefit Statement.
- Patients should be informed on post-procedure medications.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Guidezilla is a trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.