# **OptiCross™**

## 40 MHz Coronary Imaging Catheter

<b>en</b>	<b>Instructions for Use</b>	<b>2</b>
<b>es</b>	<b>Instrucciones de uso</b>	<b>8</b>
<b>fr</b>	<b>Mode d'emploi</b>	<b>14</b>
<b>de</b>	<b>Gebrauchsanweisung</b>	<b>20</b>
<b>it</b>	<b>Istruzioni per l'uso</b>	<b>26</b>
<b>nl</b>	<b>Instructies voor gebruik</b>	<b>32</b>
<b>pt-EU</b>	<b>Instruções de Utilização</b>	<b>38</b>

**TABLE OF CONTENTS**

**REUSE WARNING STATEMENT** ..... 3

**DEVICE DESCRIPTION** ..... 3

    Contents ..... 3

    Operating Principle..... 3

    Materials..... 3

    Non-pyrogenic..... 3

    User Information..... 3

    Intended Patient Population ..... 3

**INTENDED USE/INDICATIONS FOR USE** ..... 3

    Clinical Benefit Statement ..... 3

    Summary of Safety and Clinical Performance..... 3

**CONTRAINDICATIONS**..... 3

**WARNINGS**..... 3

**PRECAUTIONS**..... 4

**ADVERSE EVENTS** ..... 4

**CONFORMANCE TO STANDARDS** ..... 4

    Acoustic Output Precaution..... 4

    Imaging Catheter Information..... 4

    Calculation of Estimated In Situ Intensities ..... 4

**ACOUSTIC OUTPUT REPORTING TABLE FOR TRACK 1** ..... 4

    Auto-Scanning Mode..... 4

        Table 1. Acoustic Reporting for Track 1 ..... 5

        Table 2. Terminology..... 5

**ACOUSTIC OUTPUT – COMPLIANCE TO IEC 60601-2-37:2008**..... 5

    Prudent-Use Statement..... 5

    Imaging Catheter Information..... 5

        Table 3. Acoustic Output Reporting Table..... 5

        Table 4. Terminology ..... 5

**HOW SUPPLIED**..... 5

    Device Details ..... 5

    Handling and Storage..... 5

    Operating Environment..... 5

    Transport Environment..... 6

    Storage Environment..... 6

**OPERATIONAL INSTRUCTIONS** ..... 6

    Additional Items for Safe Use ..... 6

    Preparation..... 6

    A. Inspection Prior to Use ..... 6

    B. Preparation for Catheter Use ..... 6

    C. Preparation for MDU5 PLUS Sterile Bag use ..... 6

    Procedure ..... 6

    D. Place Guide Catheter ..... 6

    E. Introduce Imaging Catheter Into Guide Catheter..... 6

    F. Catheter Placement and Imaging..... 6

    G. Troubleshooting ..... 7

    Disposal ..... 7

    Post-Procedure ..... 7

**INFORMATION TO BRIEF THE PATIENT**..... 7

**REFERENCES**..... 7

**WARRANTY**..... 7

    Trademark..... 7

**SYMBOL DEFINITIONS**..... 7

    Figure 1. OptiCross 40 MHz Coronary Imaging Catheter ..... 44

    Figure 2. Normal test image ..... 44

    Table 5. Symbol Definitions ..... 45

# OptiCross™

## 40 MHz Coronary Imaging Catheter

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

#### REUSE WARNING STATEMENT

Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully read all instructions prior to use. Observe all contraindications, warnings, precautions and adverse events noted in these directions. Failure to do so may result in patient complications. Boston Scientific relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.

**Note:** Throughout this IFU, the iLab Guidance System may be referenced, and its function is synonymous with the function of iLab Polaris MultiModality Guidance System and AVVIGO Guidance Systems.

#### DEVICE DESCRIPTION

The OptiCross (40 MHz Coronary Imaging Catheter) is a sterile, short rail Imaging Catheter.

It consists of two main assemblies:

1. Imaging Core
2. Catheter Body

The imaging core is composed of a hi-torque, flexible, rotating drive cable with a radial looking 40 MHz ultrasonic transducer at the distal tip. An electro-mechanical connector interface at the proximal end of the catheter makes the connection to the Motordrive Unit (MDU5 PLUS) Instrument.

The catheter body is comprised of three sections:

1. Distal Imaging Window Lumen
2. Proximal Shaft Lumen
3. Telescope Section

The distal imaging window lumen and proximal shaft lumen sections comprise the "working length" of the catheter, and the telescope section remains outside of the guiding catheter.

The catheter body has a distal imaging window lumen with proximal exit at 1.6 cm from the distal end (Figure 1). A radiopaque (RO) marker is embedded in the catheter body at 0.5 cm from the distal tip. In addition, two insertion depth markers are located on the proximal shaft lumen at 90 cm and 100 cm from the distal tip to aid in estimating catheter position relative to the distal guide catheter tip. The proximal shaft lumen is attached to the telescope section via a strain relief connection.

The telescope section allows the imaging core to be advanced and retracted for 15 cm of linear movement. The corresponding movement of the transducer occurs from the proximal end of the guidewire exit port to the proximal end of the distal imaging window lumen. The telescope section has proximal markers for lesion length assessment, consisting of a series of marks spaced 1 cm apart on the telescope body.

A flush port with a one-way check valve (Figure 1) is used to flush the interior of the catheter body and maintain a flushed condition. The catheter must be flushed with heparinized saline prior to use, as this provides the acoustic coupling media required for ultrasonic imaging. The one-way check valve helps retain saline in the catheter during use.

#### Contents

OptiCross 40 MHz Coronary Imaging Catheter

17.78 cm (7 in) extension tube

3 cm<sup>3</sup> (3 cc) and 10 cm<sup>3</sup> (10 cc) syringes

4-way stopcock

#### Operating Principle

The catheter is introduced into the body via standard percutaneous transluminal technique. It is then advanced and placed at the region of interest with the aid of fluoroscopic visualization. At the site of examination, the transducer is used for ultrasonic imaging as instructed by the Instructions for Use; results are displayed on the compatible BSC IVUS multi-modality guidance console systems branded iLab, iLab Polaris, and AVVIGO. Upon imaging completion, the catheter is withdrawn, and the introduction site is closed. The drive cable and PZT transducer rotate independently of the sheath assembly to provide 360° image resolution. The transducer converts electrical impulses sent by the motor drive into transmittable acoustic energy. Reflected ultrasound signals are converted back to electrical impulses, returned to the motor drive unit, and are ultimately processed by the compatible BSC guidance console system equipment for visualization.

#### Materials

 Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1 % weight by weight.

**Note:** This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices, do not cause an increased risk of cancer or adverse reproductive effects.

#### Non-pyrogenic

This device meets pyrogen limit specifications for all patient-contacting parts.

#### User Information

Intravascular ultrasound examination of coronary anatomy should be performed only by physicians experienced in interventional cardiology or interventional radiology and in the techniques of intravascular ultrasound, and in the specific approach to be used, in a fully-equipped cardiac catheterization lab.

#### Intended Patient Population

Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal coronary interventional procedures.

#### INTENDED USE/INDICATIONS FOR USE

This catheter is intended for ultrasound examination of coronary intravascular pathology only. Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal coronary interventional procedures.

#### Clinical Benefit Statement

The OptiCross 40 MHz Coronary Imaging Catheter provides visualization of coronary vessels to aid decision making for interventions in the coronary vasculature.

#### Summary of Safety and Clinical Performance

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

#### CONTRAINDICATIONS

Use of this Imaging Catheter is contraindicated where introduction of any catheter would constitute a threat to patient safety. The contraindications also include the following patient characteristics:

- Bacteremia or sepsis
- Major coagulation system abnormalities
- Patients diagnosed with coronary artery spasm
- Patients disqualified for CABG surgery
- Patients disqualified for PTCA
- Severe hemodynamic instability or shock
- Use of the imaging catheter to cross a total occlusion

#### WARNINGS

- Do not use device after indicated 'Use By' date. Use of an expired device could result in patient injury due to device degradation.
- The catheter has no user serviceable parts. Do not attempt to repair or to alter any component of the catheter assembly as provided. Using an altered catheter can result in poor image quality or patient complications.
- No modification of this equipment is allowed.
- Air entrapped in the catheter and flushing accessories can cause potential injury. Always verify that the catheter and flushing accessories have been properly cleared of air prior to inserting the catheter into the vasculature.
- Do not pinch, crush, kink or sharply bend the catheter at any time. This can cause poor catheter performance, vessel trauma or patient complications. An insertion angle greater than 45° is considered excessive.

- Do not manipulate, advance and/or withdraw the coated device through a metal cannula or needle. Manipulation, advancement and/or withdrawal through such a metal device may result in destruction and/or separation of the outer hydrophilic coating, resulting in coating material remaining in the vasculature, which may cause adverse events and require additional intervention.
- Never advance or withdraw the Imaging Catheter without fluoroscopic visualization because it may cause vessel trauma or patient complications.
- Do not advance the catheter if resistance is encountered. The catheter should never be forcibly inserted into lumens narrower than the catheter body or forced through a tight stenosis. A catheter that is forcibly advanced may cause catheter damage resulting in vessel trauma or patient complications.
- When advancing the catheter through a stented vessel, catheters that do not completely encapsulate the guidewire may engage the stent between the junction of the catheter and guidewire, resulting in entrapment of catheter/guidewire, catheter tip separation, and/or stent dislocation.
- If resistance is met upon withdrawal of the catheter, verify resistance using fluoroscopy, then remove the entire system simultaneously. A catheter that is forcibly removed may cause vessel trauma or patient complications.
- When readvancing a catheter after deployment of stent(s), at no time should a catheter be advanced across a guidewire that may be passing between one or more stent struts. A guidewire may exit between one or more stent struts when recrossing stent(s). Subsequent advancement of the catheter could cause entanglement between the catheter and the stent(s), resulting in entrapment of catheter/guidewire, catheter tip separation and/or stent dislocation. Use caution when removing the catheter from a stented vessel.
- Inadequately apposed stents, overlapping stents, and/or small stented vessels with distal angulation may lead to entrapment of the catheter with the stent upon retraction. When retracting the catheter, separation of a guidewire from an Imaging Catheter or bending of the guidewire may result in kinking of the guidewire, damage to the catheter distal tip, and/or vessel trauma. The looped guidewire or damaged tip may catch on the stent strut resulting in entrapment.
- If multiple insertions are required, the catheter should not be disconnected from the MDU5 PLUS Motordrive Unit to avoid possible breach of catheter sterility.

## PRECAUTIONS

- Do not attempt to connect the catheter to electronic equipment other than the designated systems because the catheter may not function properly.
- Never attempt to attach or detach the catheter while the motor is running. To do so may damage the connector.
- If difficulty is encountered when backloading the guidewire into the distal end of the catheter, inspect the guidewire exit port for damage before inserting the catheter into the vasculature. The use of a damaged guidewire exit port could increase the resistance of catheter advancement or withdrawal.
- Never advance the Imaging Catheter without guidewire support because it can cause difficulty in reaching the intended region of interest or can cause the distal catheter tip to kink.
- Never advance the distal tip of the Imaging Catheter near the very floppy end of the guidewire. This part of the guidewire will not adequately support the catheter. A catheter advanced to this position may not follow the guidewire when it is retracted and cause the guidewire to buckle into a loop which the catheter may drag along the inside of the vessel and catch on the guide catheter tip. If this occurs, it may be necessary to remove the catheter assembly, guidewire and the guide catheter together. If the catheter is advanced too near the end of the guidewire, advance the guidewire while holding the Imaging Catheter steady. If this fails, withdraw the catheter and guidewire together.
- Never advance or withdraw the Imaging Catheter without the imaging core assembly being positioned at the most distal position of the imaging window because it may cause the catheter to kink.
- During and after the procedure, inspect the catheter carefully for any damage which may have occurred during use. Multiple insertions may lead to catheter exit port dimension change/distortion which could increase the chance of the catheter catching on the stent. Care should be taken when re-inserting and/or retracting catheter to prevent exit port damage.
- Always turn the MDU5 PLUS Motordrive Unit "OFF" before withdrawing the imaging catheter because it could cause the motor drive to overload.
- Cobalt may be present in the device but is not considered to be hazardous.

## ADVERSE EVENTS

Potential adverse events which may be associated with vascular imaging include but are not limited to:

- Allergic reaction
- Angina
- Cardiac arrest

- Cardiac arrhythmias including, but not limited to ventricular tachycardia, atrial/ventricular fibrillation and complete heart block
- Cardiac tamponade/Pericardial effusion
- Death
- Device entrapment requiring surgical intervention
- Embolism (air, foreign body, tissue or thrombus)
- Hemorrhage/Hematoma
- Hypotension
- Infection
- Myocardial infarction
- Myocardial ischemia
- Need for additional intervention or surgery
- Radiation injury
- Stroke and Transient Ischemic Attack
- Thrombosis
- Vessel occlusion and abrupt closure
- Vessel trauma including, but not limited to dissection and perforation

## CONFORMANCE TO STANDARDS

### Acoustic Output Precaution

In principle, it is always beneficial to minimize acoustic output to the patient. There is one scan parameter that can be varied which can cause a change in the radiated ultrasound field. The motor speed (frame rate) can vary downwards from its preset value of 30 frames per second. The maximum in situ intensities will be generated when the motor speed is 30 frames per second. It should also be noticed that the gain setting cannot change the in situ intensity.

Additional acoustic output information can be found in the operator's manual or user's guide.

### Imaging Catheter Information

Acoustic power output varies between different models of imaging catheters. Each imaging catheter supplied by Boston Scientific Corporation is supplied with Instructions for Use (IFUs) that include statements and tables specifying their acoustic power outputs.

US FDA guidelines for measurements and definitions of terms may be found in FDA publication: *Information for Manufacturers Seeking Market Clearance of Diagnostic Ultrasound Systems and Transducers* (September 9, 2008). Where Mechanical and/or Thermal Indices are reported, the MI/TI is displayed on the lower right of the IVUS image next to the imaging catheter identification.

### Calculation of Estimated In Situ Intensities

The estimated spatial peak in situ intensities are calculated from the spatial peak water values using the following equation:

$$I_{in\ situ} = I_{water} \exp(-0.069f_c z_{sp})$$

where  $I_{in\ situ}$  is the estimated in situ intensity,  $I_{water}$  is the measured intensity in water,  $f_c$  is the center frequency of ultrasound in MHz, and  $z_{sp}$  is the distance from the catheter surface to the measurement point in centimeters, 0.075 cm in this case. It should be noted that because of the complex acoustic properties of living tissue, the estimated in situ intensity may not be the same as the actual in situ intensity, and therefore, it should not be interpreted as such.

## ACOUSTIC OUTPUT REPORTING TABLE FOR TRACK 1

### Auto-Scanning Mode

Transducer Model: OptiCross 40 MHz Coronary Imaging Catheter

Operating Mode: B

System Model: iLab Polaris MultiModality Guidance System, AVVIGO Guidance Systems and Motor Drive Units

Application(s): Fetal Imaging & Other

**Note:** The US FDA guidance document "Information for Manufacturers Seeking Market Clearance of Diagnostic Ultrasound Systems and Transducers" dated September 9, 2008 classifies intravascular ultrasound within the application(s) "Fetal Imaging & Other" to determine the maximum allowable acoustic output energy. **The catheter is not intended for fetal imaging.**

**Table 1. Acoustic Reporting for Track 1**

ACOUSTIC OUTPUT			MI	I <sub>SPTA,3</sub> (mW/cm <sup>2</sup> )	I <sub>SPPA,3</sub> (W/cm <sup>2</sup> )
<b>Global Maximum Value</b>			0.302	5.89	100.01
<b>Associated Acoustic Parameters</b>	p <sub>r,3</sub>		(MPa)	1.37	
	W <sub>0</sub>		(mW)		0.026
	f <sub>c</sub>		(MHz)	40.30	40.30
	z <sub>sp</sub>		(cm)	0.075	0.075
	Beam Dimensions	x-6	(cm)		0.039
		y-6	(cm)		0.04
	PD		(μsec)	0.054	0.054
	PRF		(Hz)	7 680	7 680
EDS	Az.	(cm)		0.267	
	Elev.	(cm)		0.053	
<b>Operating Control Conditions</b>	No operator controls affecting acoustic output				

All intensities and total power have uncertainty of +28.7 % to -23.4 %.  
 All pressure values have uncertainty of +17.9 % to -16.0 %.  
 All center frequency values have uncertainty of +6.7 % to -3.7 %.

**Table 2. Terminology**

TERMINOLOGY		
Term	Definition	Units
<b>MI</b>	Mechanical Index, defined as $MI = p_{r,3} / (f_c^{1/2})$	n/a
I <sub>SPTA,3</sub>	Derated Intensity, Spatial Peak Temporal Average	mW/cm <sup>2</sup>
I <sub>SPPA,3</sub>	Derated Intensity, Spatial Peak Pulse Average	W/cm <sup>2</sup>
p <sub>r,3</sub>	Derated Peak Negative Pressure at a location of the maximum derated pulse intensity integral	MPa
W <sub>0</sub>	Total Power	mW
f <sub>c</sub>	Center frequency	MHz
z <sub>sp</sub>	Distance in the z axis direction where the measurements were taken	cm
<b>x-6 and y-6</b>	-6 dB dimensions for In Plane (azimuth) and Out of Plane (elevation) at the x-y plane where z <sub>sp</sub> is obtained	cm
<b>PD</b>	Pulse duration	μs
<b>PRF</b>	Pulse repetition frequency	Hz
<b>EDS</b>	Entrance dimensions of scanning for azimuth and elevation to a plane	cm

**Note:** Since the iLab System are identical with respect to the ultrasound generator, the acoustic output values provided above also apply to the iLab System with MDUS PLUS Motordrive Unit.

**ACOUSTIC OUTPUT – COMPLIANCE TO IEC 60601-2-37:2008**

**Prudent-Use Statement**

It is the responsibility of the system operator to understand the risk of the acoustic outputs generated by the imaging system and its associated imaging catheters. It is also the operator's responsibility to act appropriately to mitigate such risks. To that end, Boston Scientific Corporation has reported Mechanical and/or Thermal Indices that may exceed the requirements of IEC 60601-2-37:2008.

Please note that the Mechanical Index (MI) displayed on the system's screen has not been corrected for finite amplitude effects.

**Imaging Catheter Information**

Acoustic power output varies between different models of imaging catheters. Each imaging catheter supplied by Boston Scientific Corporation is supplied with Instructions for Use (IFUs) that include statements and tables specifying their acoustic power outputs.

IEC requirements for measurements and definition of terms may be found in IEC 60601-2-37: 2008 – "Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment."

Where Mechanical and/or Thermal Indices are reported, the MI/TI is displayed on the lower right of the IVUS image next to the imaging catheter identification.

**Table 3. Acoustic Output Reporting Table**

Table 201.103 ACOUSTIC OUTPUT REPORTING TABLE (in accordance with IEC 60601-2-37)			
Index Label		MI	
Maximum index value		0.302	
Index Component Value			
Associated Acoustic Parameters	p <sub>r,α</sub> at z <sub>MI</sub>	(MPa)	1.37
	P	(mW)	
	P <sub>1x1</sub>	(mW)	
	z <sub>s</sub>	(cm)	
	z <sub>b</sub>	(cm)	
	z <sub>MI</sub> (Depth for Mechanical Index)	(cm)	0.075
	z <sub>pII,α</sub> (Depth for peak attenuated pulse intensity integral)	(cm)	0.075
	f <sub>awf</sub>	(MHz)	40.30
Other Information	t <sub>d</sub>	(μsec)	0.054
	prf	(Hz)	7680
	srr	(Hz)	30
	n <sub>pps</sub>		1
	I <sub>pa,α</sub> at z <sub>pII,α</sub>	(W/cm <sup>2</sup> )	100.01
	I <sub>spta,α</sub> at z <sub>pII,α</sub> or z <sub>sII,α</sub>	(mW/cm <sup>2</sup> )	5.89
	I <sub>spta</sub> at z <sub>pII</sub> or z <sub>sII</sub>	(mW/cm <sup>2</sup> )	7.80
	p <sub>r</sub> at z <sub>pII</sub>	(MPa)	1.53
<b>NOTE 1:</b> There are no user controls that affect the catheter values provided in this table.			
<b>NOTE 2:</b> The Index Component Value is not specified as there is only one operating mode for this IVUS System.			

**Table 4. Terminology**

TERMINOLOGY		
Term	Definition	Units
<b>MI</b>	Mechanical Index, defined as $MI = \frac{p_{r,\alpha} f_{awf}^{-1/2}}{C_{MI}}$	n/a
<b>C<sub>MI</sub></b>	1 MPa MHz <sup>-1/2</sup>	n/a
<b>p<sub>r,α</sub></b>	Attenuated peak-rarefactional acoustic pressure	MPa
<b>f<sub>awf</sub></b>	Acoustic working frequency	MHz
<b>P</b>	Output power	mW
<b>TIS-Scan</b>	Soft tissue thermal index	n/a
<b>z</b>	Distance from the source to a specified point	cm
<b>A<sub>aprt</sub></b>	-12 dB output beam area	cm <sup>2</sup>
<b>t<sub>d</sub></b>	Pulse duration	μs
<b>prf</b>	Pulse repetition rate	Hz
<b>p<sub>r</sub></b>	Peak-rarefactional acoustic pressure	MPa
<b>I<sub>pa,α</sub></b>	Attenuated pulse-average intensity	W/cm <sup>2</sup>
<b>I<sub>pII</sub></b>	Pulse-intensity integral	J/m <sup>2</sup>

**HOW SUPPLIED**

Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

**Device Details**

Do not use if package is damaged or unintentionally opened before use.  
 Do not use if labeling is incomplete or illegible.

**Handling and Storage**

**Operating Environment**

Ambient Temperature: 10 °C to 40 °C  
 Relative Humidity: 30 % to 75 %  
 Atmospheric Pressure: 70 kPa to 106 kPa

## Transport Environment

Temperature: -29 °C to 60 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

## Storage Environment

Ambient Temperature: 15 °C to 30 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

## OPERATIONAL INSTRUCTIONS

### Additional Items for Safe Use

- 0.014 in (0.36 mm) Guidewire of appropriate length
- Introducer/Guide catheter of appropriate size and length
- Heparinized saline
- BSC iLab imaging console
- BSC MDU5 PLUS Motordrive Unit
- BSC Sterile Sled Bag
- BSC Automatic Pullback Sled

### Preparation

**Note:** Medical Electrical Equipment requires special precautions regarding Electromagnetic Compatibility (EMC). This equipment (device) needs to be installed and put into service according to the EMC information contained within the documents accompanying the system.

**Note:** Portable and mobile Radiofrequency (RF) communications equipment can affect Medical Electrical Equipment.

**Note:** This catheter has been tested to IPX7 on its distal end.

### A. Inspection Prior to Use

Before use, inspect the packaging for any violation of the sterile barrier and inspect the catheter and accessories for any damage or defects. Do not use potentially contaminated or defective equipment. If the sterile barrier integrity is compromised or the contents damaged, contact your Boston Scientific representative. Prior to imaging, all equipment to be used during the procedure should be carefully examined to ensure proper performance. If device appears to be compromised, contact your Boston Scientific representative.

Prior to use, verify product is within labeled shelf life. Do not use product if the 'Use By' date has been exceeded. Dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

### B. Preparation for Catheter Use

1. Refer to the setup section of the User's Guide(s) or Instructions for Use for iLab Imaging System, MDU5 PLUS Motordrive Unit, MDU5 PLUS Sterile Bag and Automatic Pullback Sled (if used).

**Note:** All of the following steps have to be performed using Sterile Technique.

2. Using sterile technique, remove the catheter and accessories from tray. Retract the movable imaging core completely to the proximal position via the telescoping shaft. DO NOT pull too hard while retracting the imaging core.
3. Fill the 3 cm<sup>3</sup> (3 cc) and 10 cm<sup>3</sup> (10 cc) syringes with heparinized saline. Connect the 3 cm<sup>3</sup> (3 cc) and 10 cm<sup>3</sup> (10 cc) syringes to the 4-way stopcock, with the 3 cm<sup>3</sup> (3 cc) syringe connected to the inline port of the 4-way stopcock. Then connect the assembly to the extension tube. Ensure that all air is expelled from the assembly (syringes, stopcock, and extension tubing) by flushing, using heparinized saline in syringes. Connect the extension tube to the one way valve on the catheter hub. The 10 cm<sup>3</sup> (10 cc) syringe is to be used as a reservoir for refilling the 3 cm<sup>3</sup> (3 cc) flushing syringe.
4. Flush the imaging catheter TWICE on the prep table continuously with 3 cm<sup>3</sup> (3 cc) volume each time. Do not use excessive pressure. Move the imaging catheter over to the procedure table. Ensure that all air is expelled from the system.
5. Connect the imaging catheter to the MDU5 PLUS Motordrive Unit by aligning the catheter hub and the MDU5 PLUS Motordrive Unit. Push the catheter hub and MDU5 PLUS Motordrive Unit together until the hub clicks into place. To ensure that the hub is fully seated in the MDU5 PLUS Motordrive Unit, gently tug on the hub. If catheter identification is incorrect or missing: See Section Troubleshooting.
6. Remove catheter carefully from sterile hoop dispenser coil. Confirm the imaging core is in the fully retracted position and the catheter is not tightly coiled. Turn on the MDU5 PLUS Motordrive Unit and confirm proper function of the catheter by observing a pattern of partial bright concentric rings on the monitor (Figure 2).
7. While imaging with MDU5 PLUS Motordrive Unit, advance the imaging core to the fully distal position, via the telescoping shaft.

**Note:** Always turn the MDU5 PLUS Motordrive Unit "ON" before using the telescope to advance the imaging core within the catheter.

8. Turn off the MDU5 PLUS Motordrive Unit. The MDU5 PLUS Motordrive Unit should remain off from this point until the catheter is positioned in situ.
9. Refill the 10 cm<sup>3</sup> (10 cc) syringe as needed and reattach to the stopcock without introducing air into the line.
10. To prevent air from being introduced into the catheter lumen, DO NOT retract the imaging core prior to catheter placement. Any amount of retraction of the imaging core prior to catheter placement will require additional flushing.

**Note:** Where pullback device use is desired, flush the catheter one more time while the imaging core is in the full distal position with the catheter installed on the pullback device.

**Note:** If flushing is difficult with the imaging core in the full distal position, then manually retract the imaging core 3 mm - 5 mm and re-flush. Then, manually advance the imaging core to the original full distal position.

**Note:** Exercise care not to kink catheter while handling.

### C. Preparation for MDU5 PLUS Sterile Bag use

Refer to MDU5 PLUS Sterile Bag Instructions for Use.

### Procedure

#### D. Place Guide Catheter

1. Prepare the entry site with a sheath introducer according to standard practice.
2. Before insertion of the imaging catheter, ensure the patient has been prepared using standard procedure for interventional treatment.
3. Place the guide catheter and Y-adapter. Introduce the guidewire and advance it to the region of interest.

#### E. Introduce Imaging Catheter Into Guide Catheter

1. Wet the distal section (approximately 23 cm) of the imaging catheter sheath with heparinized saline to activate the lubricious coating. Always wipe down the guidewire with heparinized saline prior to loading the catheter onto the guidewire.
2. Backload the guidewire into the distal end of the imaging catheter (Figure 1). Advance the guidewire into the imaging catheter until the guidewire exits from the guidewire exit port.

**Note:** Guidewires that supply more stiffness near the distal tips are recommended.

3. Continue to advance the imaging catheter into the guide catheter up to the exit point using the appropriate proximal marker as a reference if needed. Tighten the hemostasis valve on the guide catheter's Y-adapter. Tighten only enough to prevent fluid/blood leakage. An excessively tightened hemostasis valve may distort the image due to binding of the rotating drive cable.
4. Turn on the MDU5 PLUS and check to see that the catheter produces an image. If the image is flickering, some air may still be present in the catheter. Flush the catheter again with the MDU5 PLUS "ON" (Imaging). DO NOT USE EXCESSIVE PRESSURE. The image should appear as a single bright concentric ring. After confirmation of stable image, press imaging button on MDU5 PLUS Motordrive Unit to stop imaging.

#### F. Catheter Placement and Imaging

1. With MDU5 PLUS "OFF" and using fluoroscopy, advance the imaging catheter over the guidewire until the distal marker crosses a minimum of 3 cm beyond the region of interest in the vessel/lesion.
2. Keeping the catheter body and guidewire fixed, turn MDU5 PLUS Motordrive Unit "ON" and retract the imaging core slowly along its length of travel (up to maximum of 15 cm) either manually or using the optional Automatic Pullback Sled to image the region of interest. Retract and advance as desired.

**Note:** Always turn the MDU5 PLUS Motordrive Unit "ON" before advancing or retracting the imaging core within the catheter.

**Note:** If the image fades: See Section Troubleshooting.

3. When done imaging, fully advance the imaging core and turn "OFF" the MDU5 PLUS Motordrive Unit. Maintain the position of the guidewire and remove the imaging catheter.
4. If the imaging catheter is to be re-inserted, flush once with 3 cm<sup>3</sup> (3 cc) volume syringe and coil the catheter and set aside the MDU5 PLUS Motordrive Unit and Automatic Pullback Sled, if used.

**Note:** If multiple insertions are required, the catheter should not be disconnected from the MDU5 PLUS Motordrive Unit to avoid possible breach of catheter sterility.

5. When ready to reinsert the imaging catheter, flush one more time with 3 cm<sup>3</sup> (3 cc) volume syringe.
6. Inspect guidewire exit port prior to re-inserting to verify no damage occurred during withdrawal.

## G. Troubleshooting

1. If there is vibration from the catheter telescope section during advancement of the imaging core, stop imaging. Adjust the location of the sheath such that the imaging window is as straight as possible. Reactivate MDU5 PLUS Motordrive Unit and re-advance the imaging core. Be aware of potential imaging core entanglement if catheter is utilized subsequently.
2. If imaging catheter is not recognized by the system, contact your Boston Scientific representative before proceeding.
3. If the image fades during use, or if shadowed areas persist after flushing in situ, the distal imaging window lumen may contain air bubbles. Remove catheter and repeat flushing procedure in Section, Preparation for Catheter Use, step 4.
4. If the image cannot be recovered as a result of flushing, a drive cable failure or MDU5 PLUS Motordrive Unit disconnection may have occurred. Stop imaging and verify that the hub is fully seated in the MDU5 PLUS Motordrive Unit. If the hub is fully seated and the condition persists, withdraw the catheter. Restart the MDU5 PLUS Motordrive Unit and visually inspect for rotation of the imaging core. If it is not rotating, return the catheter to your Boston Scientific representative for analysis.

## Disposal

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

## Post-Procedure

Assess patient for hematoma and/or other signs of bleeding at the puncture site.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of the imaging catheter in association with the interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document, both for the imaging catheter and for interventional treatments likely to be employed.
- Discuss patient allergies in particular the risk for patients who may be allergic to contrast.
- Discuss the post procedure instructions, including any follow-up requirements, lifestyle changes, medications, and home-care or rehabilitation guidelines.

## REFERENCES

The physician should consult recent literature on current medical practice on imaging catheters.

## WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

## Trademark

The following are trademarks of Boston Scientific Corporation or its affiliates: OptiCross, MDU5 PLUS, iLab and AVVIGO.

All other trademarks are the property of their respective owners.

## SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at [www.bostonscientific.com/SymbolsGlossary](http://www.bostonscientific.com/SymbolsGlossary).

Additional symbols are defined at the end of this document.

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2024 Boston Scientific Corporation or its affiliates.