Vercise™ Neural Navigator 5 Software Programming Manual

Instructions for Use

Instructions for Use	en
Mode d'emploi	fr
Gebrauchsanleitung	de
Instructies voor gebruik	nl
Instrucciones de uso	es
Istruzioni per l'uso	it
Instruções de Utilização	pt-br
Οδηγίες χρήσης	el
Bruksanvisning	no
Brugsanvisning	da
Käyttöohjeet	fi
Bruksanvisning	sv
Pokyny k použití	cs
Návod na použitie	sk
Használati útmutató	hu
Instrukcja obsługi	pl
Инструкции за употреба	bg
Upute za upotrebu	hr
Instrucţiuni de utilizare	ro
Navodila za uporabo	sl
Kasutusjuhised	et
Lietošanas norādījumi	lv
Naudojimo instrukcijos	lt

How to Use This Manual

This manual describes the usage of the Boston Scientific Vercise™ Neural Navigator Software. Throughout this manual, the name “Boston Scientific DBS System” refers to the following: Vercise Genus™, Vercise Gevia™, and Vercise™ PC Deep Brain Stimulation (DBS) Systems.

The products covered in this manual are components of a DBS System. Read all instructions carefully before using the DBS System. For indications for use, contraindications, warnings, precautions, adverse events, storage and handling, sterilization, device disposal, and post-procedure information, refer to the *DBS Information for Prescribers IFU*. For other device-specific information not included in this manual, refer to the appropriate Instructions for Use (IFU) for your Boston Scientific DBS System as listed in your *DBS Reference Guide*. Advise the patient that additional information may be available to them on the Boston Scientific website www.bostonscientific.com/patientlabeling.

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

Trademarks

Vercise™, Vercise Gevia™, Vercise Genus™, Cartesia™, Beamsearch™, Guide™, Illumina 3D™, and STIMVIEW™ are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

The **Bluetooth®** word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Boston Scientific Corporation is under license. See also, Declaration ID D035363.

Warranty

For device warranty information, visit www.bostonscientific.com/warranty.

Summary of Safety and Clinical Performance

For customers in the European Union: visit <https://ec.europa.eu/tools/eudamed> to view a summary of this device's safety and clinical performance information on the EUDAMED website.

Technical Support

There are no user serviceable parts. If you have a specific question or issue, contact your sales representative. To contact Boston Scientific for any other reason, use the contact information provided for your locality via www.bostonscientific.com.

Product Model Numbers

Model Number	Description	Compatible System
*DB-7164 and NM-7164	Clinician Programmer	- Vercise Genus™ Deep Brain Stimulation System - Vercise Gevia™ Deep Brain Stimulation System - Vercise™ PC Deep Brain Stimulation System
*DB-7164-R and NM-7164-R	Clinician Programmer (Refurbished)	
NM-7165	Keyboard	
NM-7171	Programmer Power Supply	
NM-6316	Power Adapter	
DB-7105-N5A and DB-7105-N501A	Vercise™ Neural Navigator 5	- Vercise Gevia™ Deep Brain Stimulation System - Vercise™ PC Deep Brain Stimulation System
DB-7190 and NM-7190	Programming Wand	
DB-5270	Vercise™ Remote Control, RC 4	Vercise Genus™ Deep Brain Stimulation System
DB-6386	Pairing Magnet	

*Applicable after installation of Vercise Neural Navigator 5.0 Software (version 9028429-500 and later).

Table of Contents

Introduction.....	1
Safety Information	1
Intended Use	1
Setting Up.....	2
Place the Implanted IPG into CP Mode (Pairing Mode) (Vercise Genus IPGs Only)	3
Place ETS 3 into CP Mode (Pairing Mode).....	4
Magnet Pairing (Vercise Genus IPGs Only).....	4
Connect the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)	4
Starting a Programming Session.....	5
Connect a Stimulator to Vercise Neural Navigator Software	5
Configuring the Stimulator	7
Link to Brainlab Elements	7
Update Patient Information	7
Configure Leads	8
Check Impedances	8
Viewing the Dashboard.....	9
Program Tile	10
Patient Tile	10
Device Tile	11
Report Tile	11
Programming the Stimulator	12
Conduct a Programming Session	12
Create, Copy, or Delete a Program.....	13
Select an Area	13
Select a Programming Mode	14
Steering Mode	14
Cathodic Steering	14
Anodic Steering	14
Program in Steering Mode.....	15
Program a Linear Lead in Steering Mode:.....	15
Program a Directional Lead in Steering Mode.....	15
Program a Directional Lead using the Beamsearch Methodology.....	16
Custom Mode	17
Program in Custom Mode.....	17
Modify Amplitude.....	18
Amplitude Ramp	18
Modify Pulse Width	19
Modify Rate.....	19
Programming Multiple Areas with Different Rates	19
Modify Program Cycle or Ramp Settings.....	19
Turn Stimulation Off	20
Enable, Disable, or Modify Patient Control	20
Image Guided Programming.....	21
STIMVIEW™ Technology	21
STIMVIEW™ XT Technology.....	22
Load Simulation from GUIDE™ XT Element.....	23
Illumina 3D™ Programming Tool.....	24
To calculate an Illumina 3D setting:	24
To apply a proposed Illumina 3D setting:.....	26
Record and Map Clinical Effects.....	27
Document the Therapeutic Window.....	27

Document Clinical Effect Details.....	28
Visualize the Therapeutic Window.....	29
Map Clinical Effects	30
Stimulator Battery	31
Non-Rechargeable Stimulators.....	31
Energy Use Index	31
Vercise Genus P8, P16, and P32 Non-Rechargeable Stimulators	31
Vercise PC Non-Rechargeable Stimulator.....	31
Elective Replacement Indicator (ERI) Message	32
End of Service (EOS) Message.....	32
Rechargeable Stimulators	32
Estimated Charge Time	32
Charging History.....	32
Voltage.....	32
End a Programming Session	33
Reporting.....	34
Access Report During a Programming Session.....	34
Customize or De-Identify Reports.....	34
Access, Print, or Export Patient Records.....	34
Tools	35
Updates	35
Delete Clinical Effects Data	35
Delete Patient Records	35
Additional Information	36
Stimulator Programmable Characteristics	36
Charge Density	37
Programming Software Information	38
Software Information	38
Label Information.....	38
Programming Wand Information.....	38
Software Licenses	39
Prism.....	39
HIDSharp.....	39
Newtonsoft JSON.NET	41
SharpDX	41
Security.Cryptography	41
Polly.....	42
Combined .NET Foundation Open Source License.....	42
XAML Behaviors.....	43
NLOpt	43
SignalRSwaggerGen	44
Swashbuckle.AspNetCore	44
Clinician Programmer Management	45
Adjust CP Date and Time.....	45
Reset ClinicUser Password	45
Update Password	45
Lost/Forgotten Password.....	46
Software Installation, Uninstallation, and Removal.....	46

Introduction

The Boston Scientific Vercise Neural Navigator Software is used to program the Vercise PC, Vercise Gevia, and Vercise Genus DBS Systems. A programming session may include the following activities:

1. Setting up
2. Starting Vercise Neural Navigator Software
3. Connecting to the Stimulator
4. Configuring the Stimulator and Leads
5. Testing different Stimulation Settings

This manual provides instructions on how to accomplish these steps and perform additional functions, such as exporting reports and backing up data.

Boston Scientific DBS Systems utilize MICC¹ technology which is designed to adapt for impedance changes and maintain consistent therapy over time. MICC technology enables current steering across Lead Contacts which is intended to provide precise positioning of stimulation.

If you have any issues, contact Boston Scientific Technical Support.

Note: Screens depicted in this manual may differ slightly from the screens on your Vercise Neural Navigator Software.

Safety Information

Intended Use

The intended use of the Neural Navigator programming software on a Clinician Programmer platform is to set and adjust stimulation parameters for the compatible Boston Scientific Stimulator.

The Pairing Magnet is intended to initiate pairing of a compatible Boston Scientific Implantable Pulse Generator with a Clinician Programmer or Remote Control.

¹ Multiple independent current control

Setting Up

The Clinician Programmer (CP) communicates with the Stimulator via wireless telemetry. Wireless programming allows patients to move about, if instructed by the clinician, while the clinician adjusts parameters. Vercise Genus devices utilize Bluetooth technology for direct wireless communication between the CP and the Stimulator. Vercise PC and Vercise Gevia devices communicate via a Programming Wand. The Programming Wand uses an inductive telemetry radio frequency (RF) link to communicate with the Stimulator.

Caution: Use only Vercise PC, Vercise Gevia, or Vercise Genus DBS System components with Vercise Neural Navigator Software. Failure to do so may result in the inability to program the Stimulator.

Caution: The CP is not considered as equipment for the patient environment as defined by IEC 60601-1. The CP and the person using the CP should not be in physical contact with the patient while programming.

1. Plug the CP into a power source.
2. Power ON the CP.
3. Login as ClinicUser. If this is the first login, enter password "bsn". You will be prompted to setup a new ClinicUser password when logging into the CP for the first time. The new ClinicUser password must be 10 or more characters in length.

Note: In regions where software is installed by a third party, follow region-specific guidance for the ClinicUser password.

Note: For guidance on resetting your ClinicUser password, see the "Reset ClinicUser Password" section of this manual.

4. Determine the type of Stimulator to be programmed and follow the appropriate sections as listed in Table 1, then continue on to the "Starting a Programming Session" section of this manual.

Table 1: Stimulator Programming Setup Reference Table

Stimulator Type	Model Numbers	Reference Section
Vercise Genus IPG(s)	DB-1408, DB-1416, DB-1432, DB-1216, DB-1232	See "Place the Implanted IPG into CP Mode (Pairing Mode) (Vercise Genus IPGs Only)" on page 3 of this manual.
ETS 3	DB-5170	See "Place ETS 3 into CP Mode (Pairing Mode)" on page 4 of this manual.
Vercise Gevia IPG	DB-1200	See "Connect the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)" on page 4 of this manual.
Vercise PC IPG	DB-1140	See "Connect the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)" on page 4 of this manual.
ETS 2	DB-5132	See "Connect the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)" on page 4 of this manual.

Place the Implanted IPG into CP Mode (Pairing Mode) (Vercise Genus IPGs Only)

1. Use the patient's Remote Control to enter the Stimulator into CP Mode (Pairing Mode). The Remote Control and Stimulator must be linked in order to place the Stimulator into CP Mode (Pairing Mode).
 - a. From the **Lock** screen, press and hold the **P** button until the screen displays "Stimulator Ready for CP Connection" (Figure 1 and Figure 2).


Figure 1. Programs Button on Remote Control


Figure 2. Stimulator Ready for CP Connection

- b. Alternatively, navigate to the  **Clinician Menu** on the patient Remote Control and select  **CP Mode: Stimulator** (Figure 3).

Figure 3. CP Mode: Stimulator

Note: If the Stimulator does not have a linked Remote Control, refer to the appropriate Remote Control IFU as listed in your DBS Reference Guide for instructions on Linking a Remote Control to the Patient's Stimulator.

Note: If a linked Remote Control is not available, see the "Magnet Pairing (Vercise Genus IPGs Only)" section of this manual.

Note: The Stimulator will automatically exit CP Mode (Pairing Mode) after two minutes if a CP connection is not established.

Place ETS 3 into CP Mode (Pairing Mode)

1. Press and hold the  **Stimulation Button** on the side of the ETS until the lights flash (Figure 4). The ETS battery indicator light will alternate blinking green and yellow, indicating that the ETS is in CP Mode (Pairing Mode) and is available for connection.

Note: The Stimulator will automatically exit CP Mode (Pairing Mode) after two minutes if a CP connection is not established.

Figure 4. ETS Battery Indicator Light

Magnet Pairing (Vercise Genus IPGs Only)

If a Remote Control is not available to place the implanted Stimulator into CP Mode (Pairing Mode), you may utilize a magnet. Place the Pairing Magnet directly over the Stimulator for at least 4 seconds, but no more than 10 seconds, and then remove it.

Note: Once the Pairing Magnet is removed, the Stimulator may automatically initiate a device reset cycle. This device reset cycle will take several seconds to complete. During this reset cycle, the Stimulator will not provide Stimulation and cannot communicate. Once the reset is complete, the Stimulator will return to normal operation.

Warning: The Pairing Magnet has a strong magnetic field that may cause interference with devices that are sensitive to magnetic fields. For appropriate guidance, consult the instructions for use for these devices.

Note: If you lose the Pairing Magnet or if the plastic coating surrounding the magnet cracks, contact Boston Scientific for a replacement.

For instructions on using the Remote Control to place the implanted IPG into CP Mode (Pairing Mode), see the “Place the Implanted IPG into CP Mode (Pairing Mode) (Vercise Genus IPGs Only)” section of this manual.

Connect the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)

1. Connect the Programming Wand to the CP using the USB Cable provided with the Programming Wand (Figure 5).
 - a. Plug the Mini USB end of the USB Cable into the USB Port on the side of the Programming Wand.
 - b. Plug the Standard USB end of the USB Cable into the USB Port on the CP.

Figure 5. Vercise PC and Vercise Gevia: Clinician Programmer (CP) and Programming Wand

2. Wait for the Wand to perform a self-test. At the end of the self-test, the Wand will beep.
3. If the Power Light on the Wand is green, place the Wand over the Stimulator.
 - a. If the Power Light on the Wand remains red, contact Technical Support.

Starting a Programming Session

Connect a Stimulator to Vercise Neural Navigator Software

1. Select the Vercise Launcher icon  on the desktop.

2. Select  to start Vercise Neural Navigator Software (Figure 6).

Note: If Brainlab Elements is present on the CP, Vercise Neural Navigator Software may be launched from within Elements .

Note: Multiple software should not be run simultaneously on the same CP (except when launching a software application from within another application).

Note: Vercise Neural Navigator Software may also be launched in Demo Mode using the Vercise Launcher. Demo Mode is used for demonstration purposes only (Figure 6).

Figure 6. Launch Screen with DEMO Option Mode

3. Upon starting Vercise Neural Navigator Software, the interface will show the **Connect Screen** and the software will automatically search for a Stimulator (Figure 7).

a. **Vercise Genus:** A Stimulator that is in CP Mode (Pairing Mode) and is within range will automatically appear on the **Connect Screen** (Figure 7). If no Stimulator is found and displayed on the **Connect Screen**, confirm that the Stimulator is in range and is in CP Mode (Pairing Mode) by viewing the Patient Remote Control Screen. Then, select the **Refresh** button on the CP (Figure 7).

Note: A Vercise Genus Stimulator must be in CP Mode (Pairing Mode) in order to establish a connection with the CP.

b. **Vercise PC and Vercise Gevia:** Press the **Scan Wand** button (Figure 7). The CP will then scan for devices using the Programming Wand. If only one Stimulator is found in range, the CP will automatically connect to that Stimulator. If no Stimulator is found, move the Programming Wand closer to the Stimulator that you are trying to connect to and select the **Scan Wand** button.

Note: If the IPG is in MRI Mode, you will be prompted to exit MRI Mode before connecting the IPG to the Neural Navigator Software (Figure 8). Exit MRI Mode to allow the IPG to connect to the Neural Navigator Software.


Figure 7. Connect Screen


Figure 8. Exit MRI Mode Prompt

4. Press the **Connect** button associated with the Stimulator (Figure 7).

Note: Select the  button at the bottom left of the Connect Screen for Connection Help (Figure 7).

- a. If the Stimulator has not yet been configured, the **Configure Screen** will appear (Figure 9).
- b. If the Stimulator has been previously configured, the **Dashboard** will appear (Figure 11).

Configuring the Stimulator

If a Stimulator is not configured with Leads, the **Configure Screen** will appear automatically upon connection (Figure 9).

Figure 9. Configure Screen

Link to Brainlab Elements

When Brainlab Elements is installed on the CP and Neural Navigator Software is launched from within Elements, you will be prompted during connection to choose whether to link the Stimulator to Elements (Figure 10). This will enable you to import data from Elements to support programming the Stimulator.

⚠ Link with Elements

Would you like to link this Stimulator with Elements data?

No

Yes

Figure 10. Link to Brainlab Elements Prompt

Update Patient Information

The Patient ID automatically defaults to the Stimulator serial number (Figure 9). The Patient ID may be edited by typing into the Patient ID field. The Implant Date displays the date on which a CP first connects to a new Stimulator (Figure 9). The Implant Date may be adjusted by selecting the Implant Date field and choosing a new date.

Configure Leads

During an initial programming session, Lead Configuration must be completed prior to navigating to the **Dashboard** or **Program Screen**. The following Lead types can be programmed by the Neural Navigator Software (Table 2):

Table 2: Programmable Lead Types

Manufacturer	Lead Model Number	VNN Lead Description	Hardware Adapter Required	Linear or Directional	Contact Spacing (Axial)	Number of Contacts	Stimulation Field Model Availability
Boston Scientific (BSC)	DB-2201	DB-2201 Standard	N/A	Linear	0.5 mm	8	Available
Boston Scientific (BSC)	DB-2202	DB-2202 Directional	N/A	Directional	0.5 mm	8	Available
Medtronic (MDT)	3387	MDT 3387 ²	M8	Linear	1.5 mm	4	Unavailable
Medtronic (MDT)	3389	MDT 3389	M8	Linear	0.5 mm	4	Unavailable

For each Lead, select the Lead Type, Stimulator Port to which the Lead is connected, and brain Hemisphere (Figure 9). Enter the Target for each Lead. For Directional Leads, select the Directional Marker Orientation.

Note: When Lead information has been imported from Elements, changes made on the Configure Screen or Device Screen rather than in Elements may cause the Stimulator to be unlinked from Elements.

Before proceeding with programming, verify the Port assignment within the **Configure Screen** accurately reflects the physical connection of the Lead to Stimulator Port(s) (Figure 9). Once a Stimulator has been initially configured, the **Dashboard** will appear automatically upon connection. You may update a Stimulator's Lead configuration at any time from the **Device Screen**.

Warning: Incorrect Lead configuration may result in unintended stimulation and/or tissue damage.

Note: If Brainlab Elements is available on the CP, you may import Lead information and objects into the Neural Navigator Software from Elements. To import or remove Elements data after initial configuration of the Stimulator, Link or Unlink from the Device Screen.

During the implantation procedure of a 32 Contact (4 Port) Stimulator, no more than two Leads should be connected to the Stimulator. Therefore, only two of the four ports are to be configured for programming.

Check Impedances

The **Configure Screen** displays the date, time, and results of the last impedance check (Figure 9). Impedances may be measured using the  button on the **Configure**, **Dashboard**, **Device**, or **Program Screens**. The impedances of each Contact may be used to verify electrical integrity of the system. When an impedance measurement is taken, impedances are assessed between a Contact and the Stimulator Case (monopolar), and between pairs of Contacts (bipolar). Impedances over 8000 Ω between a Contact and the Stimulator Case may be the result of open or unconnected wires and are displayed in yellow on the **Impedance Measurement** window. Impedances less than 200 Ω between pairs of Contacts may be the result of short circuits and are displayed in orange. Contacts that have impedances outside the expected range are marked with a  symbol on the **Dashboard** and the **Program Screen**. The most recent set of impedance measurements are included in a report that may be accessed from the **Dashboard** or the **Program Screen**. Reports may be printed or exported from the **Report Screen**. The **Report Screen** may be accessed by using the VNN menu  in the top right corner of the screen.

Caution: Programming of Contacts that have impedances outside the expected range may result in early battery depletion and/or unintended stimulation.

² When vertically steering stimulation along 1.5 mm spaced Leads, the stimulation field may not be continuous for Contact configurations between Levels.

Viewing the Dashboard

Once the Leads have been configured, select **Continue** to navigate to the **Dashboard**. If a Stimulator is already configured with Leads at the time of connection, the **Dashboard** will appear upon connection. The **Dashboard** provides an overview of relevant patient information including programmed settings, device and charging information, and historical Program usage data. The **Dashboard** displays this information across four **Tiles**: **Program**, **Patient**, **Device**, and **Report** (Figure 11). To make changes to any of the data on the **Dashboard**, select the > button in the top right corner of the corresponding Tile. The Dashboard displays the date, time, and results of the last impedance check. Impedances may be measured using the  button on the **Configure**, **Dashboard**, **Device**, or **Program Screens**.


Figure 11. Dashboard

Program Tile

The **Program Tile** displays an overview of the available Programs on the Stimulator (Figure 11). If a Program is active on the Stimulator, that Program will be displayed on the **Dashboard** by default. You may view any other Programs available on the Stimulator from the **Dashboard** by selecting the Program number (e.g. Program 2) at the top right of the **Program Tile** (Figure 11). Start programming a new Stimulator or create an additional Program by navigating to the desired Program number and selecting the  button in the top right corner of the **Program Tile** (Figure 11) or by selecting **Start Programming** (Figure 12). Make adjustments to an existing Program by selecting the Area you want to change (Figure 11). Changes may be made on the **Program Screen** (Figure 16).


Figure 12. Start Programming

Patient Tile

The **Patient Tile** displays the Patient ID and Implant Date (Figure 13). To make changes to the Patient ID or Implant Date, select the  button in the top right corner of the **Patient Tile**. Changes may be made on the **Patient Screen**.

Figure 13. Patient Tile

Device Tile

The **Device Tile** displays Stimulator information including the type of Stimulator, the battery voltage of the Stimulator, and the results of the most recent impedance check (Figure 14). When Neural Navigator is launched from within Elements, the **Device Tile** also displays the status of link to Elements data. To check impedances or view monopolar and bipolar impedance values, select the Ω Check button. To make changes to the Lead Configuration for the Stimulator or Link/Unlink Elements Data, select the > button in the top right corner of the **Device Tile**.

Note: Changes to the Lead Configuration after Programs have been created may result in deletion of existing Program data.


Figure 14. Device Tile

Report Tile

The **Report Tile** displays the Stimulator's Program Usage over the past 28 days (Figure 15). Any time spent with all stimulation off is represented in grey on the chart. For rechargeable Stimulators, the **Device Tile** also displays estimated daily and weekly charge times for the currently active Program. For non-rechargeable Stimulators, the Energy Use Index for the currently active Program is displayed. See the "Energy Use Index" section of this manual for more information.

To view, print, or export the complete Report, select the > button in the top right corner of the **Report Tile** to navigate to the **Report Screen**. The individual database record for the Stimulator may also be exported, deleted, or anonymized from the **Report Screen**.


Figure 15. Report Tile

Programming the Stimulator

From the **Dashboard**, start programming a new Stimulator or create an additional Program by navigating to the Program number and selecting the > button in the top right corner of the **Program Tile** (Figure 11) or by selecting **Start Programming** (Figure 12). From the **Dashboard**, make adjustments to an existing Program by selecting the Area you want to change (Figure 11). Changes may be made on the **Program Screen** (Figure 16). Impedances may be measured using the Ω Check button on the **Configure, Dashboard, Device, or Program Screens**.

Conduct a Programming Session


Figure 16. Program Parameters

Program parameters may be adjusted as follows:

1. **Choose Program:** Select the **Program Menu** (1, Figure 16) and choose the Program you want to modify from the drop-down.
2. **Select Area:** Select the **Area** (2, Figure 16) within the Program that you want to set up or adjust.
 - See the “Select an Area” section of this manual for more information.
3. **Select Programming Mode:** Select **Custom** or **Steering** from the **Programming Mode Drop-down** (3, Figure 16).
 - Steering is the default mode and allows you to steer monopolar current along the Lead, thus eliminating the need to turn individual Contacts ON and OFF.
 - Custom mode enables programming bipolar configurations and/or configurations using non-contiguous Contacts.
 - See the “Select a Programming Mode” section of this manual for more information.

4. **Select Contacts to Program:** Select the Contact(s) you want to program using the **Steering Controls** (4, Figure 16) or by selecting the Contact(s) directly on the **Lead Diagram** (5, Figure 16).
 - See the “Program in Steering Mode” and “Program in Custom Mode” sections of this manual for more information.
5. Modify **Amplitude** (6, Figure 16), **Pulse Width** (7, Figure 16), and/or **Rate** (8, Figure 16):
 - **Amplitude** can be adjusted in 0.1 mA or 0.5 mA increments using the +/- buttons (9, Figure 16) or set to a Target Amplitude using the **Set** button (10, Figure 16).
 - To toggle to the 0.5 mA step size option, select the **Step Size** button (11, Figure 16).
 - To turn off stimulation for the active Area, select the **Stimulation On/Off Toggle** (12, Figure 16). To turn off stimulation for all active Areas, select **STOP ALL** (13, Figure 16).
 - To give patients Amplitude Control, select the **Patient Control** button (14, Figure 16) and select the toggle to enable Patient Control on all patient Programs.
 - See the “Modify Amplitude”, “Turn Stimulation Off”, and “Enable, Disable, or Modify Patient Control” sections of this manual for more information.
 - **Pulse Width** may be adjusted by selecting the **Pulse Width** button (7, Figure 16) and choosing a new Pulse Width from the table of available options.
 - See the “Modify Pulse Width” section of this manual for more information.
 - **Rate** may be adjusted by selecting the **Rate** (8, Figure 16) and choosing a new rate from the table of available options.
 - See the “Modify Rate” section of this manual for more information.

Caution: High energy stimulations settings may result in early battery depletion.

Create, Copy, or Delete a Program

To create a new Program, select the **Program Menu** and choose one of the four Programs from the drop-down (Figure 17). The system allows you to configure up to four Programs on a Stimulator.

To copy the active Program settings to another Program, select **Copy To** within the **Program Menu**, then choose the destination Program for the copied settings (Figure 17). Program slots with existing settings are represented in bold.

To delete a Program, select **Delete** within the **Program Menu**, then choose the Program you want to delete (Figure 17).


Figure 17. Program Menu

Select an Area

When a new Program is created, an Area is automatically assigned to each brain hemisphere and named based on the defined Target and side of the brain. You may add an additional Area by selecting an empty Area (+) and choosing a configuration (e.g. Left STN). For a given Program, you may configure up to four Areas. You may delete or reassign an Area by selecting the > button on that Area and then selecting **Clear Area**. If desired, you may then select a different Hemisphere configuration.

Select a Programming Mode

Steering Mode

Steering Mode allows you to steer monopolar current along the Lead, thus eliminating the need to turn individual Contacts ON and OFF. Steering Mode incrementally shifts a percentage of the total current to the adjacent Contact(s) using current steering technology to create smooth transitions between Contacts. Steering Mode is limited to a monopolar configuration of either a single Contact or adjacent Contacts with uniform polarity.

While steering, a Virtual Contact will illustrate the axial and rotational location of stimulation on the Lead Diagram. The Virtual Contact is represented by a dotted ring illustrating the axial location of stimulation along the Lead and, for Directional Leads, an arrow illustrating the rotational orientation of stimulation around the Lead (5, Figure 16). The dotted ring and arrow indicator together form the Virtual Contact for a Directional Lead.

Note: When vertically steering stimulation along 1.5 mm spaced Medtronic Leads, the stimulation field may not be continuous for Contact configurations between Levels.

Cathodic Steering

Steering Mode will default to Cathodic Steering Mode, where the Contact(s) act as the cathode(s) and the Stimulator Case acts as the anode.

Anodic Steering

In Anodic Steering Mode, the Contact(s) act as the anode(s) and the Stimulator Case acts as the cathode. To enable Anodic Steering Mode, select the  button, then select "Anodic Steering" within the menu to toggle to the ON state (Figure 18). Once Anodic Steering has been enabled, it will be enabled for the duration of the programming session with that Stimulator unless disabled from the Settings menu. Anodic Steering is disabled by default at the start of each programming session unless the Stimulator has an existing program with settings in an Anodic Steering state.


Figure 18. Anodic Steering Toggle

Select the Area you want to program in Anodic Steering Mode, then select the **+ Anodic** button from the **Steering Drop-down** (Figure 19).

Note: When the polarity of the steering state is switched, the current distribution across Contacts will remain intact, polarities of all active Contacts will be switched, and Amplitude for the selected Area will be set to 0 mA.

Note: Stimulation with only Lead Contacts assigned as anode(s) and the IPG Case assigned as cathode is referred to as Monopolar Anodic Stimulation (MAS).


Figure 19. Anodic Steering Option

Program in Steering Mode

Select the Contact(s) you want to program using the Steering Controls (4, Figure 16) or by selecting the Contact(s) directly on the Lead Diagram (5, Figure 16). The Lead Diagram displays the percentage of anodic (+) or cathodic (-) current assigned to the Lead Contacts and Stimulator Case for a given Area. To select the Units in which Amplitude is displayed on the Contacts and Stimulator Case, select the  button and select “Display Units” within the menu to toggle the available Units. Steering Controls may be used to adjust stimulation location in increments or you may select an option from the drop-down in the center of each set of controls.

Program a Linear Lead in Steering Mode:

1. Select **Steering Mode** (3, Figure 16).
2. Select a Contact to assign it as 100%.
3. Use the  and  buttons to incrementally steer the stimulation in 10% increments along the length of the Lead. Alternatively, you may use the **Level Drop-down** to adjust in preset half-Level increments (Figure 20), or you may select the Contact directly on the Lead Diagram (5, Figure 16).

Note: The Amplitude for the selected Area will be set to 0 mA when another non-adjacent Level is selected, but not when steering in 10% increments.


Figure 20. Level Drop-down

Program a Directional Lead in Steering Mode

1. Select **Steering Mode** (3, Figure 16).
2. Select a Contact to assign it as 100%. You may create an equal spread of current across a Level of Contacts (“ring mode”) by selecting anywhere within that Level, then selecting the center button. To select a single directional segment, select anywhere within that Level, then select the corresponding Contact segment (Figure 21).


Figure 21. Directional Contact Selector

3. Use the  and  buttons to incrementally steer the stimulation in 10% increments along the length of the Lead. You may also use the **Level Drop-down** to adjust in preset half-Level increments (Figure 20). Alternatively, you may select the center button in the Directional Contact Selector to select a new Level (Figure 21).

Note: The Amplitude for the selected Area will be set to 0 mA when another non-adjacent Level is selected, but not when steering in 10% increments.

4. Use the  and  buttons to rotate and steer the stimulation circumferentially around the Lead. Each rotation occurs in 15° increments without turning stimulation off. Alternatively, you may select a single directional Contact using the Directional Contact Selector (Figure 21), or you may select one of nine preset directions for the stimulation field using the **Direction Drop-down** (Figure 22). The directional presets will steer the fully focused stimulation field in one of eight anatomical directions or put the stimulation field into “ring mode.” Ring Mode generates, from a segmented Contact Level, stimulation fields equivalent to those generated by a standard “ring” or cylindrical Contact.

Note: The Amplitude for the selected Area will be set to 0 mA when a non-adjacent direction is selected using the Directional Contact Selector or using the Direction Drop-down, but not when rotating stimulation in 15° increments.

Note: The direction displayed in the Direction Drop-down is the anatomical direction of stimulation based on the Hemisphere and Orientation of the Directional Marker on the Lead selected.


Figure 22. Direction Drop-down

5. Use the  and  buttons to circumferentially spread or shrink the focus of the stimulation field in 10% increments without turning stimulation off. Alternatively, you may select a percentage from the **Focus Drop-down** (Figure 23).

Note: The Amplitude for the selected Area will be set to 0 mA when a non-incremental focus percentage is selected, but not when focusing or spreading stimulation in 10% increments.


Figure 23. Focus Drop-down

6. To choose another starting point or to steer on another Contact, select another Contact. To select a segmented Directional Contact, select the Level, then select one of the three labeled segmented Contacts around the circumference of the center button on the Directional Contact Selector (Figure 21).

Note: The total Amplitude for the selected Area will be set to 0 mA when another Contact is selected.

Note: Stimulation using multiple independent current control (MICC) with a Directional Lead is referred to as Cartesia 3D.

Program a Directional Lead using the Beamsearch Methodology

The Beamsearch Methodology keeps stimulation on while performing the steps. Use the steering controls described in the “Steering Mode” section of this manual to perform the following steps.

Note: Increasing the Pulse Width by more than 10 microseconds (μ s) during a Beamsearch session will reset the Amplitude to 0 mA.

- Find Optimal Level:** Evaluate changes in the type and magnitude of clinical response while steering current along the length of the Lead. Use the  and  buttons.
- Find Optimal Direction:** Evaluate changes in the type and magnitude of clinical response while steering current around the Lead, at the previously defined optimal Level. Use the  and  buttons once a directional Contact or directional preset has been selected. Alternatively, use the  and  buttons to transition into a directional setting while keeping stimulation on.
- Find Optimal Amplitude Setting:** Assess therapeutic window at this best level and direction determined in Steps 1 and 2, and select an Amplitude that provides the best therapeutic benefit while minimizing side effects.

Custom Mode

You may manually assign a percentage of anodic or cathodic current to individual Contacts and the Stimulator Case in Custom Mode. To provide bipolar or multipolar stimulation, you must use Custom Mode. You may assign the Stimulator Case and all the Contacts as anode or cathode individually in Custom Mode. The External Trial Stimulator (ETS) is limited to Custom Mode since the Stimulator Case cannot be assigned as a cathode or anode.

Note: Switching from Custom Mode to Steering Mode may clear the Contact and Stimulator Case assignments and reset the Amplitude to 0 mA.

Program in Custom Mode

1. Select **Custom Mode** (3, Figure 16).
2. Select the IPG Case or Contact that you want to adjust. If it was selected, one tap will assign it as an anode (+). Another tap will reassign it as a cathode (-). Another tap will reassign it as OFF (blank). Tapping on a Contact will first select it without changing the polarity.

Note: Changing the Contact polarities will reset the Amplitude to 0 mA.
3. Select the + and - buttons within the Contact Control section to adjust the percentage of anodic or cathodic current assigned to the selected Contact (Figure 24).
 - The percentage of anodic or cathodic current can be adjusted in 1% or 10% increments.
 - To toggle to the other step size option, select the % **Step** button (Figure 24).
4. To equalize the percentage of current on all Contacts of the same polarity, select a Contact with the polarity you want to equalize, then select **Equalize** (Figure 24).

Note: Equalizing electrodes of the same polarity will reset the Amplitude to 0 mA.

5. **Clear All** will clear the Contact configuration from the Case and all active Contacts and reset the Amplitude to 0 mA (Figure 24).

Note: When programming in Custom Mode, Clinical Effects data may be viewed for the active setting when 3D Split View or 3D Overview is selected in the Display Panel.

Note: When using the ETS, monopolar configurations are not possible since the ETS “case” cannot be assigned as a cathode or anode.

Note: When using the ETS, Clinical Effects data is recorded but not plotted on the Clinical Effects Map (CEM).


Figure 24. Custom Mode Contact Control

Modify Amplitude

Amplitude can be adjusted in 0.1 mA or 0.5 mA increments using the +/- buttons (9, Figure 16) or set to a Target Amplitude using the **Set** button (10, Figure 16). To toggle to the 0.5 mA step size option, select the **Step Size** button (11, Figure 16).

en

Amplitude Ramp

To ramp the Amplitude to a selected value, select the **Set** button (10, Figure 16), then choose a Target Amplitude value (Figure 25). Once you have chosen a Target Amplitude, select your preferred Ramp Speed and apply the new Amplitude (Figure 25). To pause the ramp and keep stimulation on, select the **Pause** button (Figure 26). To pause the ramp and turn stimulation off, use the **Stimulation On/Off Toggle**. Pausing will keep stimulation on and prevent the stimulation from increasing further. From the paused state, you may resume the ramp by selecting the **Resume** button, or you may cancel the ramp by selecting the **Cancel** button (Figure 27). Canceling will keep stimulation on and cancel any further Amplitude increases associated with that ramp. To cancel the ramp and turn stimulation off, select **STOP ALL** (13, Figure 16).

Caution: Applying the Instant Ramp Ramp Speed will instantaneously set the Area Amplitude to Target Amplitude.


Figure 25. Amplitude Ramp


Figure 26. Amplitude Pause


Figure 27. Amplitude Resume or Cancel

Table 3: Amplitude - Measured in milliamperes (mA)				
Definition	Default	Range	Default Step Size	Single Contact Maximum
The Amplitude of the stimulation is the magnitude of current delivered per pulse.	0 mA	0 mA to 20 mA	0.1 mA	12.7 mA

Warning: High stimulation levels can cause permanent tissue damage. If you attempt to exceed a stimulation limit, a message will pop up and notify you that settings that exceed this limit are not allowed.

Modify Pulse Width

Pulse Width may be adjusted by selecting the **Pulse Width** button (7, Figure 16) and choosing a new Pulse Width from the table of available options. Pulse Widths that exceed stimulation limits are greyed out.

Table 4: Pulse Width - Measured in microseconds (µs)		
Definition	Default	Range
The Pulse Width of the stimulation is the length of time that current is applied per pulse.	60 µs	20 µs to 450 µs

Warning: High stimulation levels can cause permanent tissue damage.

Note: Increasing the Pulse Width by more than 10 µs at a time will cause the total Amplitude to be reset to 0 mA.

Modify Rate

Rate may be adjusted by selecting the **Rate** (8, Figure 16) and choosing a new rate from the table of available options. Restricted Rates are greyed out.

Table 5: Rate - Measured in Hertz (Hz) or pulses per second (pps)		
Definition	Default	Range
The pulse rate of the stimulation, often called the Rate or Frequency, dictates how many stimulation pulses are delivered in a second.	130 Hz	2 Hz to 255 Hz

Note: Areas that are assigned to the same Lead Port are restricted from having Rates that sum up to greater than 255 Hz.

Note: Modifying the Rate of an Area will alter the Rates of all configured Areas.

Programming Multiple Areas with Different Rates

Areas are allowed to be programmed with different Rates. By default, the Multiple Rate option is disabled. When you enable Multiple Rates, only the Rates that are compatible with the Rates and Pulse Widths from the other active Areas are available.

Note: If you disable Multiple Rates, the Rate for all Areas will be set to the Rate selected for the active Area.

Note: Modifying the Rate of an Area when Multiple Rates is enabled will alter the available Rates for the other Areas.

Modify Program Cycle or Ramp Settings

To adjust Cycle or Ramp settings for a given Program, select the  button, then select the ON/OFF status display for the setting you want to change.

- **Cycle:** The Cycle setting controls the duration of stimulation that is cycled ON (Cycle On Time) and OFF (Cycle Off Time).

Table 6: Cycle Setting	
Default	Range
Disable	1 s to 90 min

- **Ramp:** The Ramp setting dictates the time to gradually increase the stimulation from zero to the programmed Amplitude when changing Programs or when stimulation is turned ON from an OFF state using the Clinician Programmer or Remote Control.

Table 7: Ramp Setting	
Default	Range
Enabled	1 s to 10 s

Turn Stimulation Off

To turn off stimulation for the active Area, select the **Stimulation On/Off Toggle** (12, Figure 16). To turn off stimulation for all active Areas, select **STOP ALL** (13, Figure 16). This function is only meant for turning all Stimulation OFF, it does not control the stimulation of an individual Area. To turn Stimulation ON, select each Area that you want to turn ON and select the Stimulation ON/OFF switch.

Note: When the Amplitude is at 0 mA, increase the Amplitude to turn ON stimulation.

Enable, Disable, or Modify Patient Control

By default, patients do not have the ability to adjust the Amplitude of their stimulation. However, in some cases, you may want to give a patient the ability to adjust the Amplitude of their stimulation using the Remote Control. To give patients Amplitude Control, select the **Patient Control** button (14, Figure 16) and select the toggle to enable Patient Control on all patient Programs (Figure 28). Once Patient Control is enabled, you may set the allowable Amplitude range by setting a Minimum and Maximum for each area. The step size for the Minimum and Maximum Amplitude is controlled by the Total Amplitude step size.

Note: The Patient Control button displays the current status of the patient for that control area. If Patient Control is enabled, but the Minimum and Maximum are set to the same values, the Patient Control button will display “Range Undefined.” If the Minimum and Maximum are defined, the range will be displayed on the Patient Control button.

Caution: All DBS programming decisions remain in the programming clinician’s judgment. All programmed Patient Amplitude minimum and maximum settings should be screened and selected based on patient physiological response. Untested Patient Amplitude minimum or maximum settings in place of clinical judgment could result in excessive stimulation, inadequate stimulation, or stimulation to an unintended target.

Warning: High stimulation levels can cause permanent tissue damage. If you attempt to exceed a stimulation limit, a message will pop up and notify you that settings that exceed this limit are not allowed.

Warning: Patients may have the ability to change the Amplitude with the Remote Control. The clinician should set and verify the maximum and minimum Amplitude levels allowed by the Remote Control to ensure that stimulation levels remain safe.


Figure 28. Modify Patient Control

Image Guided Programming


Figure 29. Stimulation Field Model (SFM)

STIMVIEW™ Technology

Utilize the **Display Drop-down** (1, Figure 29) to switch between the 3D Split View and 3D Overview of the stimulation field and the Clinical Effects Map (CEM). The 3D Split View and 3D Overview options enable STIMVIEW for Boston Scientific Leads. STIMVIEW, or the Stimulation Field Model (SFM), is a visual representation of the estimated stimulation field for the currently programmed stimulation parameters. 3D Overview presents a three-dimensional view. The 3D Split View provides a dual pane view centered on the Lead. The top pane is in-line with the Lead and the bottom pane is on an axis perpendicular to the Lead. If objects have been imported from Elements, you may show or hide individual objects using the **Object Visibility Drop-down** (3, Figure 29).

Adjust the view of the SFM using the STIMVIEW View Controls (2, Figure 29). Utilize to Zoom, to Rotate, to Pan, or to Reset the original view. When in 3D Split View, both the Lateral and Axial views of the SFM will adjust in unison using these controls. These controls will not affect or adjust any programming parameters.

The orange directional indicator (4, Figure 29), present in both 3D Split View and 3D Overview, is a visual indicator of the orientation of the radiopaque Directional Marker band on the Directional Lead. The orange line and dot (5, Figure 29) correlate to the center of the radiopaque Directional Marker. Note that this only applies to Directional Leads.

The Reference Head (6, Figure 29), located in the bottom right corner of the Display Panel in 3D Split View and 3D Overview, demonstrates the relationship of the Lead currently being programmed to the position of the patient head.

STIMVIEW™ XT Technology

If Brainlab Elements is available on the CP, imported objects may be viewed in STIMVIEW, referred to as STIMVIEW XT when imported objects are displayed.

STIMVIEW XT uses patient-specific anatomy, Lead location and Lead orientation, imported from Elements, along with the Stimulation Field Model (SFM) to inform programming. The SFM provides a modeled, visual representation of the volume of neuronal tissue activated by the currently programmed settings. The use of the SFM does not replace observation of the patient for therapeutic effects or side effects. The SFM is the approximated stimulation field shown in the color red (Figure 29). As programming parameters are adjusted and the stimulation is steered along the Lead, the SFM will adjust accordingly.

Note: *If Brainlab Elements is not available, anatomy and Lead objects cannot be imported. However, a visual representation of the Lead will remain visible for all Leads. The SFM will only remain visible for Leads manufactured by Boston Scientific.*

Caution: *All DBS programming decisions remain in the programming clinician's judgment. All programmed settings should be screened and selected based on patient physiological response. Reliance on anatomy data in place of clinical judgment could result in excessive stimulation, inadequate stimulation, or stimulation to an unintended target.*

Note: *The final stimulation settings should be determined by the clinician when programming the patient.*

You may switch between two different views by selecting either 3D Overview or 3D Split View from the **Display Drop-down** (1, Figure 29).

Load Simulation from GUIDE™ XT Element

If Neural Navigator is linked to Brainlab data with a GUIDE XT Simulation available, you may import simulated Stimulation Field Model (SFM) parameters, referred to as a Simulation. A Simulation may be created in GUIDE XT and imported along with anatomical objects when linking Brainlab data to a Stimulator.

To import a Simulation created in GUIDE XT:

1. Select the Area into which you want to import the Simulation (Figure 30).
2. Select the  button for that Area to open the **Area Options Menu** (Figure 30).
3. Select **Load Simulation** (Figure 30).


Figure 30. Area Options Menu

4. Load your preferred Simulation.
5. Review the contact configuration and setting parameters on the **Load Simulation Settings Screen** (Figure 31). If you do not want to apply the setting as displayed, select  in the top right corner of the **Load Simulation Settings Screen** to cancel (Figure 31).
6. Select a Ramp Speed from the **Ramp Speed Drop-down** (Figure 31).
7. Select **Accept** to apply the setting to the Stimulator (Figure 31).

Note: Stimulation Amplitude will be set to 0 mA, then the Amplitude of the Simulation will increase to the Target Amplitude at the selected Ramp Speed. For more details on the Amplitude Ramp, see the “Modify Amplitude” section of this manual.


Figure 31. Load Simulation Settings

Note: Importing data from Brainlab Elements will revert Program Cycle and Program Ramp settings to the default values. If Program Cycle and/or Program Ramp settings have been modified prior to the import, they must be readjusted by selecting the Program Cycle and/or Program Ramp options from the  menu. See the “Modify Program Cycle or Ramp Settings” section of this manual for more information.

Illumina 3D™ Programming Tool

Illumina 3D is a programming tool within Neural Navigator that uses an algorithm to propose a starting point for programming based on a patient's Lead placement and anatomy. Neural Navigator must be linked to Brainlab data in order to use Illumina 3D. Illumina 3D may be utilized when programming Linear or Directional Boston Scientific Leads.

Note: The Illumina 3D programming tool is not available for Leads manufactured by Medtronic.

Illumina 3D must be enabled from the **Tools Screen** in order to utilize the feature. See the “Updates” section of this manual for instructions on enabling features. Once the feature has been enabled, select the  button on the **Program Screen**, then select **Illumina 3D** within the menu to toggle Illumina 3D to the ON state (Figure 32). Once Illumina 3D has been enabled, it will remain enabled for the duration of the programming session with that Stimulator unless disabled from the Settings menu. Illumina 3D is disabled by default at the start of each programming session unless the Stimulator you are programming has an existing program with Illumina 3D settings.


Figure 32. Illumina 3D Toggle

To calculate an Illumina 3D setting:

1. Select the Area you want to program.
2. Select the desired Pulse Width and Rate.
3. Select the **Illumina 3D** option from the **Steering Drop-down** (Figure 33).


Figure 33. Illumina 3D Steering Option

4. Select an object to target with stimulation and the object(s) you want to avoid stimulating (Figure 34).
5. Use the sliders to select settings for the calculation (Figure 34).
6. Select **Calculate** (Figure 34). You will be provided with an Amplitude estimate for the proposed setting and the corresponding Stimulation Field Model (SFM) for the proposed setting will appear in the **Display Panel** (Figure 34).
7. Review the setting. You may modify the selected objects and/or slider positions as desired. Select **Calculate** again to review the new proposed setting after making modifications (Figure 34). This process may be repeated to refine the setting.
8. Once you finish making changes to selected objects and slider positions, select **Load Settings** to review the contact configuration and parameters of the setting (Figure 34).

Figure 34. Illumina 3D Simulation

To apply a proposed Illumina 3D setting:

1. Review the contact configuration and setting parameters on the **Load Simulation Settings Screen** (Figure 35). If you do not want to apply the setting as displayed, select **✕** to return to Illumina 3D.
2. Select a Ramp Speed from the **Ramp Speed Drop-down** (Figure 35).
3. Select **Accept** to apply the setting to the Stimulator (Figure 35).

Note: Stimulation Amplitude will be set to 0 mA, then the Amplitude of the Simulation will increase to the Target Amplitude at the selected Ramp Speed. For more details on the Amplitude Ramp, see the “Amplitude Ramp” section of this manual.

Caution: All DBS programming decisions remain in the programming clinician’s judgment. All programmed settings should be screened and selected based on patient physiological response. Reliance on anatomy data in place of clinical judgment could result in excessive stimulation, inadequate stimulation, or stimulation to an unintended target.


Figure 35. Load Illumina 3D Simulation Settings

4. Once a setting has been applied to the Stimulator, the Amplitude and anatomical position of the setting may be modified using the Amplitude and Steering Controls, respectively. Adjusting the anatomical position of a setting generated by Illumina 3D is called Illumina 3D Steering.
5. To reopen Illumina 3D, select the **Solution** button (Figure 36) below the Vertical and Rotational Steering Controls.

Note: The Focus Steering Controls are not available when modifying an Illumina 3D setting.


Figure 36. Reopen Illumina 3D

Record and Map Clinical Effects


Figure 37. Clinical Effects Map

Document the Therapeutic Window

While programming, you may make note of a stimulation benefit or side effect for the active stimulation setting by selecting **Benefit** or **Side Effect**, respectively, in the bottom right corner of the **Program Screen** (Figure 37). Document the Therapeutic Window for the active settings (Contact configuration, Pulse Width, and Rate) by selecting **Benefit** at the Amplitude at which you begin to see a clinical impact of stimulation and by selecting **Side Effect** at the Amplitude at which the first side effect occurs. To undo the Benefit or Side Effect selection, select the same button a second time.

Document Clinical Effect Details

You may select **Details** in the bottom right corner of the **Program Screen** to denote more detailed clinical effects for a given stimulation setting. Within the **Clinical Effects Details Panel**, select a 0 to 4 rating for each therapeutic benefit and/or a 0 to 4 rating for each side effect you want to document by selecting the numeric value for that clinical effect (Figure 38). If selection of a numeric rating is not desired, select the name of the Therapeutic Benefit and/or Side Effect without selecting a numeric value. Each clinical effect or numeric score that is selected is captured as data associated with the active stimulation setting for that patient.

Within the **Clinical Effects Details Panel**, select **More Benefits** or **More Side Effects** to view the full list of Benefits or Side Effects or to denote body side, body location, or transience (Figure 38). You may also enter and save up to 250 characters of text associated with each Lead in the Text Note field within the Clinical Effects Details Panel (Figure 38). To remove your Benefit and/or Side Effect selection, select the name or numeric score of the Therapeutic Benefit or Side Effect that you want to remove.


Figure 38. Clinical Effects Details Panel

Visualize the Therapeutic Window

When 3D Split View or 3D Overview is selected in the **Display Panel**, Clinical Effects for the active configuration (Contact configuration, Pulse Width, and Rate) are displayed in the Clinical Effects Map (CEM) Preview at the bottom of the Display Panel (Figure 39). Upon selection of a clinical effect or numeric score for a clinical effect, a dot or ring is plotted on the CEM and the CEM Preview to document the Amplitude. The rating scale of the therapeutic benefit determines the color saturation of the dot. If a side effect is selected, an orange ring is displayed. The active stimulation Amplitude is denoted by a white ring. The CEM Preview may be used to visualize the Therapeutic Window for an active configuration. The same CEM Preview is displayed at the bottom of each Area Panel on the **Dashboard**.

A new CEM Preview is created for programming settings with different Contact configurations or with different Pulse Width and Frequency combinations.

Figure 39. Clinical Effects Map Preview

Map Clinical Effects

The Clinical Effects Map is a graphical summary of assigned Therapeutic Benefits and/or Side Effects at a given position along the DBS Lead array and a stimulation Amplitude (Figure 37). You may utilize the **Display Drop-down** to switch between 3D Split View, 3D Overview, and the Clinical Effects Map View. When the Clinical Effects Map View is selected in the Display Panel, a dot is plotted on the CEM at the Axial Lead position and Amplitude for Linear Leads. When programming directionally, the CEM switches to a polar grid. A new CEM is created for programming settings at different Levels (axial positions along the Lead), as well as for different Pulse Width and Frequency combinations.

The rating scale of the therapeutic benefit determines the color saturation of the dot. Less saturated dots are indicative of greater therapeutic benefit. If a side effect is selected, an orange ring is displayed. The active stimulation position is denoted by a white ring. A visual key indicating the icons representing the active stimulation position, side effect position, and the color saturation for a benefit score appears at the bottom of the CEM. Selecting a benefit dot or side effect ring within the Display Panel will display a pop-up window containing the date and time at which the dot was captured along with the stimulation setting and effects details (Figure 37).

All of this data is saved on the Stimulator and is available for export in the Report.

Note: Clinical effects data is captured in the Clinical Effects Map Preview and in Reports. However, for configurations not possible in Steering Mode and for Directional Lead settings that are not 100% focused or spread, the clinical effects data is not plotted on the CEM.

Note: The Reference Head, located in the bottom right corner of the Display Panel when viewing the CEM, demonstrates the relationship of the Lead currently being programmed to the Hemisphere in which the Lead is implanted and the position of the patient head.

Stimulator Battery

The impact of active Program settings on the Stimulator battery may be viewed on the **Report Tile** on the **Dashboard** or on the **Program Screen** by selecting the **Program Menu** button  and then **Battery** (Figure 17). For a non-rechargeable Stimulator, this information is represented by the Energy Use Index. For a rechargeable stimulator, this information is represented by the Estimated Charge Time.

Non-Rechargeable Stimulators

Energy Use Index

The Energy Use Index applies to non-rechargeable Stimulators only. The Energy Use Index gives you an estimate of longevity of the battery life for the selected Program. After the optimal settings have been identified for a Program, from the **Program Screen**, select the **Program Menu** button  and select **Battery** to obtain the Energy Use Index. The Energy Use Index for the currently active Program may also be viewed from the **Dashboard** on the **Report Tile**.

Use Figure 40 and Figure 41 to identify the longevity that corresponds to this Energy Use Index. The figures take into account nominal non-therapy power consumption, including shelf-life and patient Remote Control use. If the estimate for longevity obtained is below 12 months, consider evaluating a Boston Scientific rechargeable system.

Vercise Genus P8, P16, and P32 Non-Rechargeable Stimulators


Figure 40. Longevity Estimates Based on 24 Hour Per Day Usage

Vercise PC Non-Rechargeable Stimulator


Figure 41. Longevity Estimates Based on 24 Hour Per Day Usage

Elective Replacement Indicator (ERI) Message

You will not be able to connect to a non-rechargeable Stimulator that is nearing the end of its battery life. If you attempt to connect a non-rechargeable Stimulator that is in an ERI state to the Neural Navigator Software, the CP will display an ERI message and the Stimulator battery voltage as shown in Figure 42. During the ERI period, the Stimulator will continue to provide stimulation; however, no changes may be made to the Stimulator settings.

Note: The ERI Message applies to non-rechargeable Stimulators only.

Note: The battery voltage shown in Figure 42 is for demonstration purposes only. The ERI indicator varies per Stimulator type.


Figure 42. ERI Message Displayed on Connect Screen

End of Service (EOS) Message

When the Stimulator has reached its end of service, you will not be able to connect Neural Navigator Software to it and stimulation will no longer be provided. If you attempt to connect a non-rechargeable Stimulator that is in an EOS state to the Neural Navigator Software, the CP will display an EOS message as shown in Figure 43.

Note: The EOS Message applies to non-rechargeable Stimulators only.


Figure 43. EOS Message Displayed on Connect Screen

Rechargeable Stimulators

Estimated Charge Time

The Estimated Charge Time applies to rechargeable Stimulators only. The Estimated Charge Time provides an estimate of the duration and frequency of charging necessary to maintain stimulation for the selected Program. After the optimal settings have been identified for a Program, from the **Program Screen**, select the **Program Menu** button > and select **Battery** to obtain the Estimated Charge Time. Estimated Charge Time may also be viewed from the **Dashboard** on the **Report Tile** or from within the Report.

Charging History

For rechargeable Stimulators, you may view the charging history of the Stimulator from the **Device Screen** by selecting the **Usage** button. A 60 day charge history is displayed by default. If desired, you may adjust the chronological span of the charge history graph.

Voltage

Battery voltage at the time of connection to the Stimulator is displayed on the **Device Screen** and on the **Dashboard** within the **Device Tile**. On the **Device Tile**, the battery icon and associated voltage value is displayed in blue if the battery has sufficient remaining voltage. If the battery voltage is low, the battery icon and associated voltage value is displayed in red.

End a Programming Session

To end a Programming Session on the CP:

1. Select the VNN menu  in the top right corner of the screen.
 - a. Select **Exit Application**, then follow the prompt to end the programming session and close the application.
 - b. Or, select **Disconnect From Stimulator** from the same drop-down, then follow the prompt to end the programming session and disconnect from the patient's Stimulator. This will take you back to the **Connect Screen**.
2. Alternatively, from the **Dashboard**, you may also select the  button to end the programming session and disconnect from the patient's Stimulator. This will take you back to the **Connect Screen**.

All Programs and programming data are automatically saved in real-time during the programming session. No step to actively "save" is required. The Patient's Remote Control automatically syncs with the Stimulator to which it has been linked.

Caution: Do not place a magnet over a Stimulator within 60 seconds of ending a programming session. If you need to reconnect to a Stimulator, either utilize the Remote Control to initiate CP Mode (Pairing Mode), or wait 60 seconds before utilizing a magnet to place the IPG into CP Mode (Pairing Mode).

Reporting

Access Report During a Programming Session

Select the  button at the top right corner of the **Report Tile** on the **Dashboard** to view the Report for the active programming session. Alternatively, you may select **View Report** from the VNN menu . A Report may be printed and exported as a PDF or Excel file.

Customize or De-Identify Reports

Select **Customize** on the **Report Screen** and select the information you want to include in the report by checking on any of the following check boxes:

- Programs
- Configuration
- Clinical Effects Maps
- Clinical Effects Details
- De-Identify Patient Data

Access, Print, or Export Patient Records

To access reports for patients that have been previously programmed using the same Clinician Programmer, select **Patient Records** on the **Connect Screen** when the CP is not connected to a Stimulator. To view an individual patient report, select the **View** button associated with the record you want to view. To view, print, or export a single report or to export a database record for a single patient, select the  button associated with that record, then select the desired action (Figure 44). You may also view, print, or export multiple reports or export the database for multiple reports by selecting the checkbox(es) for each record you want to include and then selecting the desired action from the options at the top of the screen (Figure 44).


Figure 44. Access, Print, or Export Patient Records

Tools

The **Tools** option in the VNN menu  in the top right corner of the screen allows you to activate Leads and features, delete clinical effects data, enable patient data deletion, or change the Neural Navigator software user interface language.

Updates

The **Options Tab** within Tools allows you to enable supported Leads and features.

Note: You can only perform updates when disconnected from the Stimulator.

Note: Only features available in your region will be provided.

To enable a new Lead or feature:

1. Disconnect from the patient's Stimulator by selecting **Disconnect From Stimulator** from the VNN menu  in the top right corner of the screen.
2. Select the **Tools** option from the VNN menu  in the top right corner of the screen.
3. Select the **Options Tab**.
4. Select the grey text box under the ENTER KEY(S) column for the Lead or feature you want to activate.
5. Enter the key provided by your local representative.
6. Select **Verify key(s)**.

Delete Clinical Effects Data

All the Clinical Effects Data for a patient may be deleted from the **Tools Screen** under the **Settings Tab**.

Note: This feature is available only when the CP is linked to a Patient's Stimulator.

To delete the Clinical Effects Data:

1. Select the **Tools** option from the VNN menu  in the top right corner of the screen while connected to a patient's Stimulator.
2. Select the **Delete** button under the CLINICAL EFFECTS DATA section of the screen.
3. Select **Continue**.

Delete Patient Records

Patient records may be deleted after data deletion is enabled within the **Settings Tab** on the **Tools Screen**.

To delete Patient Records for one or all patients:

1. Disconnect from the patient's Stimulator by selecting **Disconnect From Stimulator** from the VNN menu  in the top right corner of the screen.
2. Select the **Tools** option from the VNN menu  in the top right corner of the screen.
3. Select the **Allow data deletion** checkbox.
4. Select the  button in the top left corner of the **Tools Screen** to navigate to the **Connect Screen**.
5. Select **Patient Records**.
6. Select the checkbox(es) within each patient record(s) you want to include.
7. Select **Delete** from the actions at the top of the screen.

Note: You can only enable data deletion when disconnected from a Stimulator.

Additional Information

Stimulator Programmable Characteristics

The stimulation parameters are independent for each DBS Lead such that stimulation of different brain targets may have different Amplitudes, Pulse Widths, Stimulation Rates, and Contact configurations. It is possible to configure one Lead as monopolar, and one as multipolar. It is also possible to configure a single Lead with both monopolar and multipolar Areas.

The programmable parameter ranges for the Stimulator are shown below.

Table 8: Programmable Parameter Ranges		
#	Parameter	Parameter Range
1	Waveform	Charge balanced, asymmetric biphasic
2	Pulse Shape	Rectangular
3	Current or Voltage Regulated	Current
4	Amplitude ³	0 mA to 20 mA
5	Rate ⁴	2 Hz to 255 Hz
6	Pulse Width ⁵	20 μ s to 450 μ s
7	Cycle On/Off	1 s to 90 min
8	Ramp On	1 s to 10 s
9	Contact Connections ⁶	Up to 32
10	Independent Areas of Stimulation (4 Programs with 4 Areas per Program)	16
11	Current Path Options	Monopolar, Bipolar, Multipolar

³ The programmable Amplitude for each individual Contact is limited to 12.7 mA. A programming interlock is enforced to limit the total output current to 20 mA or less per Area. For example, a maximum current output of 12.7 mA on one Contact would limit the total summed current output on the remaining Contacts to 7.3 mA within one Area.

⁴ The Rate is limited to 255 Hz for a given Area. The global Rate limit for each Lead is also 255 Hz.

⁵ The use of lower Pulse Widths than established (60 μ s to 450 μ s) is the sole responsibility of the user.

⁶ During the implantation procedure of a 32 Contact (4 Port) Stimulator, no more than two Leads should be connected to the Stimulator. Therefore, only two of the four ports are to be configured for programming.

Charge Density


Figure 45. Charge Density Limits for Supported DBS Leads

Figure 45 displays the recommended maximum charge density for different combinations of Amplitude (mA) and Pulse Width (μs). The solid black line (Limit: 6 mm²) refers to all Contacts on all cylindrical (ring or tip) Contacts. The dashed black line (Limit: 1.5 mm²) refers to the small directional Contacts of all Directional Leads. These estimates of charge density are applicable for all Leads identified in Table 2.

Warning: Patients may have the ability to change the Amplitude with the Remote Control. The clinician should set and verify the maximum and minimum Amplitude levels allowed by the Remote Control to ensure that stimulation levels remain safe.

Programming Software Information

To access software information, including information regarding the CP, software or Wand, select the VNN menu  in the top right corner of the screen and select either **Clinician Programmer Info** or **Wand Info**.

en

Software Information

To access installation, versioning, and other software information, select **Clinician Programmer Info** within the VNN menu . The following information displays:

- CP Model
- CP Name (Serial)
- CP Platform Version
- Neural Navigator Installation Confirmation Key
- Neural Navigator Install Date and Time
- Neural Navigator Software Version
- Neural Navigator Software Part Number
- Programming Server Installation Confirmation Key
- Programming Server Install Date and Time
- Programming Server Software Version
- Programming Server Part Number

Label Information

To access label information, select **Clinician Programmer Info** within the VNN menu . Select the **Label** button to access the Software Label (Figure 46).


Figure 46. Label Button

Programming Wand Information

To access Programming Wand information, including the serial number, model number, and firmware version, select the **Wand Info** button within the VNN menu  when the Wand is connected to the CP.

Software Licenses

Prism

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Copyright (c) Prism Library

All rights reserved. Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the “Software”), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED “AS IS”, WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

HIDSharp

HIDSharp

Copyright 2010-2019 James F. Bellinger <<http://www.zer7.com/software/hidsharp>>

Licensed under the Apache License, Version 2.0 (the “License”); you may not use this file except in compliance with the License. You may obtain a copy of the License at <http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an “AS IS” BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

Apache License
Version 2.0, January 2004
<http://www.apache.org/licenses/>

TERMS AND CONDITIONS FOR USE, REPRODUCTION, AND DISTRIBUTION

1. Definitions.

“License” shall mean the terms and conditions for use, reproduction, and distribution as defined by Sections 1 through 9 of this document.

“Licensor” shall mean the copyright owner or entity authorized by the copyright owner that is granting the License.

“Legal Entity” shall mean the union of the acting entity and all other entities that control, are controlled by, or are under common control with that entity. For the purposes of this definition, “control” means (i) the power, direct or indirect, to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50%) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

“You” (or “Your”) shall mean an individual or Legal Entity exercising permissions granted by this License.

“Source” form shall mean the preferred form for making modifications, including but not limited to software source code, documentation source, and configuration files.

“Object” form shall mean any form resulting from mechanical transformation or translation of a Source form, including but not limited to compiled object code, generated documentation, and conversions to other media types.

“Work” shall mean the work of authorship, whether in Source or Object form, made available under the License, as indicated by a copyright notice that is included in or attached to the work (an example is provided in the Appendix below).

“Derivative Works” shall mean any work, whether in Source or Object form, that is based on (or derived from) the Work and for which the editorial revisions, annotations, elaborations, or other modifications represent, as a whole, an original work of authorship. For the purposes of this License, Derivative Works shall not include works that remain separable from, or merely link (or bind by name) to the interfaces of, the Work and Derivative Works thereof.

“Contribution” shall mean any work of authorship, including the original version of the Work and any modifications or additions to that Work or Derivative Works thereof, that is intentionally submitted to Licensor for inclusion in the Work by the copyright owner or by an individual or Legal Entity authorized to submit on behalf of the copyright owner. For the purposes of this definition, “submitted” means any form of electronic, verbal, or written communication sent to the Licensor or its representatives, including but not limited to communication on electronic mailing lists, source code control systems, and issue tracking systems that are managed by, or on behalf of, the Licensor for the purpose of discussing and improving the Work, but excluding communication that is conspicuously marked or otherwise designated in writing by the copyright owner as “Not a Contribution.”

“Contributor” shall mean Licensor and any individual or Legal Entity on behalf of whom a Contribution has been received by Licensor and subsequently incorporated within the Work.

2. **Grant of Copyright License.** Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable copyright license to reproduce, prepare Derivative Works of, publicly display, publicly perform, sublicense, and distribute the Work and such Derivative Works in Source or Object form.
3. **Grant of Patent License.** Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in this section) patent license to make, have made, use, offer to sell, sell, import, and otherwise transfer the Work, where such license applies only to those patent claims licensable by such Contributor that are necessarily infringed by their Contribution(s) alone or by combination of their Contribution(s) with the Work to which such Contribution(s) was submitted. If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Work or a Contribution incorporated within the Work constitutes direct or contributory patent infringement, then any patent licenses granted to You under this License for that Work shall terminate as of the date such litigation is filed.
4. **Redistribution.** You may reproduce and distribute copies of the Work or Derivative Works thereof in any medium, with or without modifications, and in Source or Object form, provided that You meet the following conditions:
 - (a). You must give any other recipients of the Work or Derivative Works a copy of this License; and
 - (b). You must cause any modified files to carry prominent notices stating that You changed the files; and
 - (c). You must retain, in the Source form of any Derivative Works that You distribute, all copyright, patent, trademark, and attribution notices from the Source form of the Work, excluding those notices that do not pertain to any part of the Derivative Works; and
 - (d). If the Work includes a “NOTICE” text file as part of its distribution, then any Derivative Works that You distribute must include a readable copy of the attribution notices contained within such NOTICE file, excluding those notices that do not pertain to any part of the Derivative Works, in at least one of the following places: within a NOTICE text file distributed as part of the Derivative Works; within the Source form or documentation, if provided along with the Derivative Works; or, within a display generated by the Derivative Works, if and wherever such third-party notices normally appear. The contents of the NOTICE file are for informational purposes only and do not modify the License. You may add Your own attribution notices within Derivative Works that You distribute, alongside or as an addendum to the NOTICE text from the Work, provided that such additional attribution notices cannot be construed as modifying the License.

You may add Your own copyright statement to Your modifications and may provide additional or different license terms and conditions for use, reproduction, or distribution of Your modifications, or for any such Derivative Works as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.

5. **Submission of Contributions.** Unless You explicitly state otherwise, any Contribution intentionally submitted for inclusion in the Work by You to the Licensor shall be under the terms and conditions of this License, without any additional terms or conditions. Notwithstanding the above, nothing herein shall supersede or modify the terms of any separate license agreement you may have executed with Licensor regarding such Contributions.
6. **Trademarks.** This License does not grant permission to use the trade names, trademarks, service marks, or product names of the Licensor, except as required for reasonable and customary use in describing the origin of the Work and reproducing the content of the NOTICE file.
7. **Disclaimer of Warranty.** Unless required by applicable law or agreed to in writing, Licensor provides the Work (and each Contributor provides its Contributions) on an “AS IS” BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied, including, without limitation, any warranties or conditions of TITLE, NON-INFRINGEMENT, MERCHANTABILITY, or FITNESS FOR A PARTICULAR PURPOSE. You are solely responsible for determining the appropriateness of using or redistributing the Work and assume any risks associated with Your exercise of permissions under this License.
8. **Limitation of Liability.** In no event and under no legal theory, whether in tort (including negligence), contract, or otherwise, unless required by applicable law (such as deliberate and grossly negligent acts) or agreed to in writing, shall any Contributor be liable to You for damages, including any direct, indirect, special, incidental, or consequential damages of any character arising as a result of this License or out of the use or inability to use the Work (including but not limited to damages for loss of goodwill, work stoppage, computer failure or malfunction, or any and all other commercial damages or losses), even if such Contributor has been advised of the possibility of such damages.
9. **Accepting Warranty or Additional Liability.** While redistributing the Work or Derivative Works thereof, You may choose to offer, and charge a fee for, acceptance of support, warranty, indemnity, or other liability obligations and/or rights consistent with this License. However, in accepting such obligations, You may act only on Your own behalf and on Your sole responsibility, not on behalf of any other Contributor, and only if You agree to indemnify, defend, and hold each Contributor harmless for any liability incurred by, or claims asserted against, such Contributor by reason of your accepting any such warranty or additional liability.

END OF TERMS AND CONDITIONS**Newtonsoft JSON.NET**

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Copyright (c) 2007 James Newton-King

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

SharpDX

Copyright (c) 2010-2014 SharpDX - Alexandre Mutel

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Security.Cryptography

MIT License

Copyright (c) 2017 Microsoft

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Polly

New BSD License

=

Copyright (c) 2015-2020, App vNext

All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of App vNext nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL <COPYRIGHT HOLDER> BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Combined .NET Foundation Open Source License

Applicable to the following products:

- SignalR (ASP.NET Core)
- Microsoft.Extensions.Hosting (.NET Runtime)
- System.Runtime.CompilerServices.Unsafe (.NET Runtime)
- Microsoft.AspNetCore.Authentication.JwtBearer (ASP .NET Core)
- System.Configuration.ConfigurationManager (.NET Runtime)
- System.Diagnostics.EventLog (.NET Runtime)
- Microsoft.NETCore.App
- System.Configuration.ConfigurationManager (.NET Runtime)
- EntityFramework

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Copyright (c) .NET Foundation and Contributors

All rights reserved.

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

XAML Behaviors

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Copyright (c) 2015 Microsoft

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

NLOpt

Copyright (c) 1988-2019 the following:

- Association for Computing Machinery
- Carlos Henrique da Silva Santos
- Dieter Kraft
- Enthought, Inc.
- Jean-Sebastien Roy
- K. Madsen
- Massachusetts Institute of Technology
- Michael J. D. Powell
- NLOpt Developers
- North Carolina State University
- S. Zertchaninov
- SciPy Developers
- Sovrasov V.
- Steven G. Johnson

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- (a). Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- (b). Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- (c). Neither the name of the copyright holders nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE REGENTS OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

SignalRSwaggerGen

MIT License

Copyright (c) 2021 Dorin Mocan

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Swashbuckle.AspNetCore

Includes:

- Swashbuckle.AspNetCore.Newtonsoft
- Swashbuckle.AspNetCore.SwaggerUI

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Copyright (c) 2016 Richard Morris

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Clinician Programmer Management

Adjust CP Date and Time

If the date and/or time are incorrect, select the **Adjust Date and Time** icon  on the desktop to launch the Date and Time adjustment window. Select the **Change date and time** button to modify the date and/or time as desired and select **OK** to confirm changes, as shown in Figure 47.


Figure 47. Date and Time Screens

Reset ClinicUser Password

Update Password

If you know your ClinicUser password and would like to change it to a new password, follow the steps below:

1. Press and hold **CTRL + ALT + DEL** (control + alt + delete) on the keyboard and select **Change a password**.
2. Enter your old password in the “Old password” field, and enter your new password in the “New password” and “Confirm password” fields. The ClinicUser password must be 10 or more characters in length.
3. Ensure your new password is the same in each field. Then press **Enter**.

Lost/Forgotten Password

If you are unable to log in to the ClinicUser profile on the CP due to a lost or forgotten password, follow the steps below to reset your ClinicUser password:

1. From the CP login screen, select the PasswordReset user profile (Figure 48).


Figure 48. Login Screen with PasswordReset Profile

2. Refer to the “*Technical Support*” section of this manual and call your local Technical Support team. Technical Support will provide a password that will allow you to log in to the PasswordReset user profile.
3. Enter the password provided by Technical Support to log into the PasswordReset user profile.
4. When prompted, enter a new ClinicUser password in both the “New Password” and “Confirm Password” fields. The new ClinicUser password must be 10 or more characters in length. Ensure the ClinicUser password is the same in each field. Then select **Change Password** (Figure 49).

Note: If desired, you may exit the PasswordReset user profile without changing your ClinicUser password by selecting the  icon at the top right of the screen.


Figure 49. Update Password

5. Once the ClinicUser password has been changed successfully, you will see a message indicating the change was successful. Select **OK** to return to the login screen.
6. Using your new password, log in to the ClinicUser profile.

Software Installation, Uninstallation, and Removal

For instructions on how to install, uninstall, or remove programming software, refer to the *Software Installation Guide* for your Boston Scientific DBS System as listed in your *DBS Reference Guide*.

This page intentionally left blank.

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AREAS
Left STN > 1.9 mA
Right STN > 1.7 mA
+ Area
+ Area

CLINICAL EFFECTS DETAILS

BENEFITS	LEAST	BEST
Benefit	4 3 2 1 0	
Rigidity (1.0)	4 3 2 1 0	1
Tremor (0.0)	4 3 2 1 0	0
Bradykinesia	4 3 2 1 0	
More Benefits >		
SIDE EFFECTS	MILD	SEVERE
Side Effect	0 1 2 3 4	
Speech (Mod)	0 1 2 3 4	1
Paresthesia	0 1 2 3 4	
Muscle Pulling	0 1 2 3 4	
Discomfort	0 1 2 3 4	
Ocular Effect	0 1 2 3 4	
More Side Effects >		
TEXT NOTE		

PROGRAMMING MODE
Steering > **Stimulation** ON 0.1 Step
1.9 mA
60 µs 130 Hz
STEERING Cathodic
LEVEL L2.9
DIRECTION Ant-Lat
FOCUS 40%
PATIENT CONTROL Disabled
CLINICAL EFFECT
Benefit Side Effect
Details

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The screenshot shows the 'CONFIGURE' screen with a dark blue header. On the left, there are four sections: 'PROFILE' with Patient ID (IPG500238) and Implant Date (07 Dec 2022); 'STIMULATOR' with Type (Vercise Genus™ R16), Model (DB1216), and Serial (500238); 'IMPEDANCE' with a 'Check' button and status 'In Range'; and 'ELEMENTS INFO' with a 'Link' status and patient details. On the right, there is a diagram of the stimulator with two leads (L and R) connected to a brain model. Below the diagram are two columns of settings for the Left and Right leads, including Lead Type (DB2202-Direc), Port(s) (L and R), Hemisphere (Left Brain and Right Brain), Target (STN), and Orientation (Anterior-Med and Anterior-Later). A 'Continue' button is at the bottom right.

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AREAS
Left STN 1.9 mA
Right STN 1.7 mA
+ Area
+ Area

CLINICAL EFFECTS DETAILS

BENEFITS	LEAST	BEST
Benefit	4 3 2 1 0	
Rigidity (1.0)	4 3 2 1 0	1
Tremor (0.0)	4 3 2 1 0	0
Bradykinesia	4 3 2 1 0	
More Benefits		

SIDE EFFECTS	MILD	SEVERE
Side Effect	0 1 2 3 4	
Speech (Mod)	0 1 2 3 4	1
Paresthesia	0 1 2 3 4	
Muscle Pulling	0 1 2 3 4	
Discomfort	0 1 2 3 4	
Ocular Effect	0 1 2 3 4	
More Side Effects		

PROGRAMMING MODE
Steering
STIMULATION ON 0.1 Step
1.9 mA
60 µs 130 Hz
Cathodic
LEVEL L2.9
DIRECTION Ant-Lat
FOCUS 40%
PATIENT CONTROL Disabled

CLINICAL EFFECT
Benefit Side Effect
Details

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Copyright© 2007 James Newton-King
Copyright© 2010-2014 SharpDX - Alexandre Mutel
Copyright© 2017 Microsoft
Copyright© 2015 Microsoft
Per "Opera" si intende l'opera dell'autore, in formato Sorgente o Oggetto, resa disponibile in base a questa Licenza, come indicato nelle informazioni sul copyright incluse o allegate all'opera (un esempio è fornito nell'Appendice riportata di seguito).
or its affiliates. All rights reserved.