Boston Scientific

en

es

fr

de

it

nl

sv

fi

no

da

pt

pt BR

Precision™ Spinal Cord Stimulator System 1.5 Tesla MRI Guidelines

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Rx ONLY

CAUTION:

Federal law restricts this device to sale, distribution and use by or on the order of a physician.

ENGLISH

This document is a supplement to the Precision System product manuals and relates specifically to the use of a transmit/receive radio frequency (RF) head coil of a 1.5 Tesla horizontal bore MRI system for patients implanted with the Boston Scientific Precision Spinal Cord Stimulator (SCS) System. The implanted portions of the Precision SCS System can include Boston Scientific percutaneous leads, lead extensions, surgical paddle leads, surgical accessories, and implantable pulse generator (IPG).

Important Note: MRI procedures should be performed using only a transmit/receive RF head coil in a 1.5-Tesla horizontal bore MRI system, ONLY. Do not use open-sided or MRI systems operating at other static magnetic field strengths. The risks of using MRI systems operating at other static magnetic field strengths or using a transmit or transmit/receive RF body coil have not been determined.

The Precision SCS System is “MR Conditional”¹ when exposed to the magnetic resonance imaging (MRI) environment under the specific conditions defined in this document. These include the strength of the static magnetic field, spatial gradient of the static magnetic field, time rate of change of the magnetic field (dB/dt), use of a radio frequency (RF) transmit/receive head coil, head averaged specific absorption

rate (SAR), and Normal Mode MRI scanner operation. Additional conditions, including specific configurations of the Precision SCS System, are required as specified throughout this document.

It is important to read this information in its entirety before conducting or recommending an MRI examination on a patient with a Precision SCS System. These instructions do not apply to other implantable products or other devices, products, or items. Contact Boston Scientific at the address or phone number listed at the end of this document if you have any questions. These instructions also appear on the Boston Scientific website (www.bostonscientific.com).

Definition of terms

ASTM International – the American Society for Testing and Materials, an international standards forming body.

dB/dt – the change of magnetic field strength per unit of time, measured in units of Tesla per second.

Hertz (Hz) – a unit of frequency in Hertz or cycles per second.

MR Conditional¹ – an item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use. Field conditions that define the specified MR environment include field strength, spatial gradient, dB/dt (time rate of change of the magnetic field), radio frequency (RF) fields, and specific absorption rate (SAR). Additional conditions, including

¹ ASTM F 2503-08, “Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment”

specific configurations of the item, may be required.

MRI – Magnetic Resonance Imaging.

Radio Frequency (RF) – electrical fields whose frequencies are in the range of 10,000 Hz and above.

Specific Absorption Rate (SAR)¹ – the mass normalized rate at which RF energy is deposited in biological tissue. SAR is typically indicated in W/kg.

Tesla (T) – the unit of measure of magnetic field strength. One T is equal to 10,000 Gauss.

MRI Transmit/Receive RF Body coil – a coil used to transmit and to receive RF energy that encompasses the entire body region within the MR system bore.

MRI Transmit/Receive RF Head coil – a coil used to transmit and to receive RF energy that is constrained to the head region, only.

W/kg – Watts per kilogram, a measure of the power that is absorbed per kilogram of tissue.

INTRODUCTION

Magnetic Resonance Imaging (MRI) uses a powerful static magnetic field, gradient magnetic fields, and RF energy to create images useful for the diagnosis of various diseases and conditions. These electromagnetic fields may interact with implants and cause movement and/or excessive heating of or around the implanted components, damage to the device electronics and/or voltage induction through the leads and Stimulator causing an uncomfortable “jolting” sensation. In addition, the metallic components of the

neurostimulation systems can affect the MRI image, potentially impacting the diagnostic use of this modality. Therefore, because of these potential deleterious effects, it is imperative to thoroughly read, understand, and comply with the instructions, cautions, and warnings in this document to prevent harm to the patient undergoing MRI examination.

MRI procedures involving the use of transmit or transmit/receive RF body coils are strictly prohibited due to the high likelihood of causing excessive MRI-related heating of the implanted components.

If all of the instructions stated in this document are followed, MRI examinations of the head using a transmit/receive RF head coil of a 1.5-Tesla horizontal bore MRI system may be safely performed.

Boston Scientific has conducted testing consistent with ASTM standards F2052² for magnetically induced displacement force, F2213³ for magnetically induced torque, and F2182⁴ for RF heating as well as functional testing of the Precision™ SCS System to

² ASTM F2052-02, “Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment”

³ ASTM F2213-04, “Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment”

⁴ ASTM F2182-02a, “Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants during Magnetic Resonance Imaging”

define the recommendations and guidelines contained within this document.

The instructions in this document describe how to conduct a *head-only* MRI examination of a patient with the Precision SCS System, using a transmit/receive RF head coil of a 1.5-Tesla, horizontal bore MRI system.

WARNINGS

MRI Transmit or Transmit/Receive RF Body Coil. Patients implanted with the Precision System should not be subjected to a MRI transmit or transmit/receive RF body coil. MRI exposure may result in dislodgement of implanted components, heating of the neurostimulator, damage to the device electronics and/or voltage induction through the leads and Stimulator causing an uncomfortable or “jolting” sensation.

MRI Transmit/Receive RF Head Coil. An MRI examination of the head *only* (no other part of the body) can be conducted safely using *only* a transmit/receive RF head coil of a 1.5-Tesla horizontal bore MRI system when all instructions in this document are followed. The use of a transmit/receive RF head coil constrains the majority of the RF fields to the head region and away from the implanted components of the Precision System. MRI examinations performed on patients with the Precision System using open-sided MRI systems, or other types of MRI systems operating at other static magnetic field strengths (higher or lower) have not been assessed and, as such, should not be performed.

Limitations

- MRI should not be considered for patients with the Precision SCS System if other potentially safer diagnostic methods such as computed tomography, X-ray, or other methods will provide adequate diagnostic information for patient management.
- The instructions in this document apply to all implanted components of the Precision SCS System, consisting of Boston Scientific IPG, leads, surgical accessories, and extensions. Contact Boston Scientific at the telephone number or address listed at the end of this document for information about newer models or any updates.
- If the patient has any other *active* or *passive* medical implant from a manufacturer that prohibits or contraindicates an MRI examination, follow the instructions from the manufacturer. The instructions in this document apply only to Boston Scientific products described herein.
- Physicians should not prescribe MRI for patients undergoing trial neurostimulation and/or having systems that are not fully implanted.
- Do not remove the IPG and leave the lead system implanted as this can result in higher than expected lead heating. Excessive heating can cause tissue damage.
- Allow at least six weeks time from the date of Precision SCS System implantation to the time of the MRI examination.

Risks Associated with MRI Examination

Exposing a patient with an implanted neurostimulation system or component to MRI may potentially injure the patient or damage the Precision™ SCS System. The known potential risks are as follows:

- Induced currents associated with the MRI examination in a patient with the Precision SCS System or any of its components may cause heating, especially at the lead electrode site, resulting in tissue damage. Induced electrical currents may also stimulate or shock the patient.
- Factors that increase the risks of patient injury include, but are not limited to, the following:
 - Performing MRI scans outside Normal Operating Mode
 - The use of a transmit/receive RF head coil that extends over any component of the Precision SCS System (i.e. the implanted leads, IPG, etc).
- Exposure to gradients exceeding a dB/dt limit of 20 Tesla per second may result in overstimulation, shocking, or damage to the IPG.
- MRI may permanently damage the neurostimulator, requiring device explantation or replacement.
- MRI may affect the operation of the Precision SCS System. The MRI may also reset the parameters to power-on-reset settings, requiring reprogramming with the Precision SCS System Clinician's Programmer.

- The Precision IPG may move within the implant pocket which may cause patient discomfort or a recent IPG pocket incision to open.
- SCS implant geometries or locations outside the IFU have not been evaluated.
- Longer Lead lengths and use of extensions may increase the risk of overstimulation, shocking, or damage to the IPG.

MRI PROCEDURE USING A TRANSMIT/RECEIVE RF HEAD COIL

If all of the instructions stated in this document are followed, MRI examinations of the head using a transmit/receive RF head coil in a 1.5-Tesla horizontal bore MRI system may be safely performed.

Supervision

Patient interaction during MRI – The patient should be in a psychological condition and mental state where the patient is able to provide immediate feedback of any problems during the examination.

Monitor the patient during the MRI examination. Verify that the patient is feeling normal and is responsive during and between each individual scan sequence of the MRI examination. Discontinue the MRI immediately if the patient becomes unresponsive to questions or experiences any heating, pain, shocking sensations/ uncomfortable stimulation, or unusual sensations.

Prior to the MRI examination, an individual with the proper knowledge of MRI equipment such as an MRI-trained radiologist or MRI physicist must ensure the MRI examination will be conducted according to the information outlined in this document.

Note: Due to the additional requirements in these instructions, MRI examination time may be significantly extended.

MRI exposure requirements

Prior to an MRI examination, determine whether the patient has another active or passive medical implant. If the patient has any other *active* or *passive* medical implant from a manufacturer that prohibits or contraindicates an MRI examination then you should NOT conduct an MRI examination. The instructions in this document apply only to Boston Scientific products described herein.

If the following requirements cannot be met, **do NOT** proceed with the MRI examination.

- Use only a transmit/receive RF head coil.
 - **Important Note:** If you are unsure if your MRI system has a transmit/receive RF head coil capability or if you are unsure if your MRI system displays “head SAR”, check with your MRI system manufacturer.
- No component of the Precision SCS System (i.e. the implanted leads, IPG, etc) extends into the head coil.
- Use only a 1.5 Tesla, horizontal bore MRI system (i.e., do not use open-sided or MRI systems operating at other static magnetic field strengths).

- Ensure MRI system is in the Normal Operating Mode.
- The MRI pulse sequence parameters must meet the following requirements. If they do not, the parameters must be adjusted so that they comply with these requirements.
 - Ensure dB/dt is limited to no more than 20 Tesla per second.
 - Use MRI examination pulse sequence parameters that limit the MRI system-reported head SAR level to 1.5 W/kg or less for all sequences.
- Enter the correct patient weight into the MRI console to ensure that the head SAR level is estimated correctly by the MRI system.

Note: The requirements provided in this document are based on *in-vitro* testing and should result in a safe MRI examination of a patient with the Precision SCS System when all instructions in this document are followed. However, due to the many variables that affect safety, the safety of patients or continued functionality of the Precision SCS System exposed to MRI cannot be absolutely ensured.

Preparation for the MRI examination

The following procedures must be followed prior to performing an MRI examination on a patient with the Precision SCS System:

1. Determine if the patient has an implant or condition that would prohibit or contraindicate the MRI examination. Do NOT conduct an MRI examination if any implant or condition is present that would prohibit or contraindicate the procedure.
2. Allow at least six weeks from the date of Precision™ SCS System implantation to the time of the MRI examination.
3. Inform the patient of the risks of undergoing an MRI examination as stated in this document including implanted component movement and/or excessive heating, and that the metallic components of the neurostimulation systems can affect the MRI image, potentially impacting the diagnostic use of this modality.
4. The patient should be in a psychological condition and mental state where the patient is able to provide immediate feedback of any problems during the examination.
5. Instruct the patient to immediately inform the MRI system operator if any discomfort, stimulation, shocking, or heating occurs during the examination.
6. Carefully monitor the patient throughout the MRI examination using audio and visual techniques.
7. The radiologist or an MRI physicist must verify that all proposed MRI examination parameters comply with the “MRI exposure requirements” listed above. If not, the parameters must be modified to meet these requirements. If parameters cannot be modified, do NOT perform the MRI examination.
8. If the patient has implanted leads but does not have an IPG, do NOT perform an MRI examination.
9. If the patient has an implanted neurostimulator, perform the following steps:
 - a. Review the neurostimulator stimulation parameters with a clinician programmer and print out a copy of the programmed stimulation parameters for reference.
 - b. Ensure that no Precision SCS external components (i.e. Remote Control, or Charger) enters within the designated boundary of the MRI system.
10. The Precision SCS System must be programmed to stimulation OFF before allowing the patient within the designated safety boundary of the MRI system.

During the MRI examination

- Continuously monitor the patient both visually and audibly. Check the patient during and between each imaging sequence. Discontinue the MRI examination immediately if the patient is unable to respond to questions or reports any problems.
- Conduct the MRI examination using only the pulse sequences and conditions that the MRI-trained radiologist or MRI physicist has confirmed meets the “MRI exposure requirements” outlined in this document.

Post-MRI examination review

- Verify that the patient feels normal.
- Verify that the neurostimulator is functional.
- Reprogram the neurostimulator to pre-MRI settings.

AUS Australian Sponsor Address
Boston Scientific (Australia) Pty Ltd

PO Box 332
BOTANY
NSW 1455
Australia
Free Phone 1800 676 133
Free Fax 1800 836 666

EC REP EU Authorised Representative
Boston Scientific International S.A.

55 Avenue des Champs Pierreux
TSA 51101
92729 NANTERRE CEDEX
FRANCE

European Headquarters – Paris

T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

Argentina

T: +5411 4896 8556 F: +5411 4896 8550

Australia / New Zealand

T: 1800 676 133 F: 1800 836 666

Austria

T: +43 1 60 810 F: +43 1 60 810 60

Balkans

T: 0030 210 95 37 890 F: 0030 210 95 79 836

Belgium

T: 080094 494 F: 080093 343

Brazil

T: +55 11 5853 2244 F: +55 11 5853 2663

Bulgaria

T: +359 2 986 50 48 F: +359 2 986 57 09

Canada

T: +1 888 359 9691 F: +1 888 575 7396

Chile

T: +562 445 4904 F: +562 445 4915

China – Beijing

T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou

T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai

T: +86 21 6391 5600 F: +86 21 6391 5100

Colombia

T: +57 1 629 5045 F: +57 1 629 5082

Czech Republic

T: +420 2 3536 2911 F: +420 2 3536 4334

Denmark

T: 80 30 80 02 F: 80 30 80 05

Finland

T: 020 762 88 82 F: 020 762 88 83

France

T: +33 1 39 30 49 00 F: +33 1 39 30 49 01

Germany

T: 0800 072 3301 F: 0800 072 3319

Greece

T: +30 210 95 42401 F: +30 210 95 42420

Hong Kong

T: +852 2960 7100 F: +852 2563 5276

Hungary

T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore

T: +91 80 5112 1104/5 F: +91 80 5112 1106

India – Chennai

T: +91 44 2648 0318 F: +91 44 2641 4695

India – Delhi

T: +91 11 2618 0445/6 F: +91 11 2618 1024

India – Mumbai

T: +91 22 5677 8844 F: +91 22 2617 2783

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea

T: +82 2 3476 2121 F: +82 2 3476 1776

Malaysia

T: +60 3 7957 4266 F: +60 3 7957 4866

Mexico

T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East / Gulf / North Africa

T: +961 1 805 282 F: +961 1 805 445

The Netherlands

T: +31 30 602 5555 F: +31 30 602 5560

Norway

T: 800 104 04 F: 800 101 90

Philippines

T: +63 2 687 3239 F: +63 2 687 3047

Poland

T: +48 22 435 1414 F: +48 22 435 1410

Portugal

T: +351 21 3801243 F: +351 21 3801240

Singapore

T: +65 6418 8888 F: +65 6418 8899

South Africa

T: +27 11 840 8600 F: +27 11 463 6077

Spain

T: +34 901 111 215 F: +34 902 267 888

Sweden

T: 020 65 25 30 F: 020 55 25 35

Switzerland

T: 0800 826 786 F: 0800 826 787

Taiwan

T: +886 2 2747 7278 F: +886 2 2747 7270

Thailand

T: +66 2 2654 3810 F: +66 2 2654 3818

Turkey – Istanbul

T: +90 216 464 3666 F: +90 216 464 3677

Uruguay

T: +59 82 900 6212 F: +59 82 900 6212

UK & Eire

T: +44 844 800 4512 F: +44 844 800 4513

Venezuela

T: +58 212 959 8106 F: +58 212 959 5328

NOTE: Phone numbers and fax numbers may change. For the most current contact information, please refer to our website at <http://www.bostonscientific-international.com/> or write to the following address:

Boston Scientific Neuromodulation

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² ASTM F2052-02, „Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment“

³ ASTM F2213-04, „Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment“

⁴ ASTM F2182-02a, „Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants during Magnetic Resonance Imaging“

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² ASTM F2052-02, "Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment"

³ ASTM F2213-04, "Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment"

⁴ ASTM F2182-02a, "Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants during Magnetic Resonance Imaging"

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or its affiliates. All rights reserved.
© 2014 Boston Scientific Corporation