Vercise™ Neural Navigator Programming Manual

Directions for Use

Directions for Use	en
Instrucciones de uso	es
Mode d'emploi	fr
Gebrauchsanweisung	de
Istruzioni per l'uso	it
Gebruiksaanwijzing	nl
Bruksanvisning	sv
Käyttöohje	fi
Bruksanvisning	no
Brugsanvisning	da
Instruções de Utilização	pt
Instruções de uso	ptBR
Kullanım Talimatları	tr
Указания по использованию	ru
إرشادات الاستخدام	ar
Οδηγίες χρήσης	el
Návod k použití	cz
Návod na použitie	sk
Instrukcja obsługi	pl
Használati utasítások	hu
Указания за употреба	bg

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Trademarks

All trademarks are the property of their respective holders.

Additional Information

For indications, contraindications, warnings, precautions, adverse events summary, sterilization, component disposal, storage and handling, and warranty information, refer to the Information for Prescribers DFU for your Boston Scientific DBS System. For other device specific information not included in this manual or labeling symbols, refer to the appropriate DFU for your Boston Scientific DBS System as listed on your DBS Reference Guide.

Product Model Numbers

Model Number	Description
*DB-7161	Vercise™ DBS Clinician Programmer
*DB-7161-R	Vercise DBS Clinician Programmer (refurbished)
*DB-7161-203	Clinician Programmer
*DB-7161-203-R	Clinician Programmer (refurbished)
*NM-7161	Clinician Programmer
*NM-7161-R	Clinician Programmer (refurbished)
DB-7105-N2	Vercise Neural Navigator Software 2 Installer
NM-7190 and DB-7190	Programming Wand
NM-6316	International Power Converter

*Applicable after installation of Vercise Neural Navigator 2 (Software version 9028429-200).

Table of Contents

Introduction	1
Intended Use.....	1
Connecting the Programming Wand to the Clinician Programmer.....	1
Starting a Programming Session	2
Starting Vercise Neural Navigator.....	2
Connecting to the Stimulator.....	2
Programming the Stimulator	4
Configuring the Leads.....	4
Measuring Impedances.....	4
The Programming Screen.....	5
Creating or Modifying a Program.....	7
Selecting Stimulation Areas.....	7
Selecting Contacts.....	7
Steering Mode.....	7
Custom Mode.....	8
Turning Stimulation OFF for Individual Areas.....	8
Turning All Stimulation OFF.....	8
Increasing and Decreasing the Amplitude.....	8
Increasing and Decreasing the Pulse Width.....	9
Increasing and Decreasing the Rate.....	9
Programming Multiple Areas with Different Rates.....	9
Selecting Patient Amplitude Range.....	9
Viewing the Stimulation Field Model.....	9
Mapping the Patient's Clinical Effects of Stimulation.....	10
Reports	11
Tools and Additional Information	12
Backup Database.....	12
Update Leads.....	13
Deleting Clinical Effects Data.....	14
Stimulator Programmable Characteristics.....	14
Charge Density.....	15
Energy Use Index.....	15
Estimated Charge Time.....	16
Elective Replacement Indicator (ERI) Message.....	16
End of Service (EOS) Message.....	16
Ending a Programming Session.....	16
Adjusting CP Time and Date.....	17
Troubleshooting	18

This page intentionally left blank.

Introduction

Vercise™ Neural Navigator is a software program that is used to set and adjust stimulation parameters for the Vercise PC and Vercise Gevia™ DBS Systems.

An initial programming session may consist of the following steps:

1. Starting Vercise Neural Navigator
2. Connecting to the Stimulator
3. Configuring the Stimulator and Leads
4. Testing different Stimulation Settings

This manual will provide instruction on how to accomplish these steps as well as perform additional functions, such as exporting reports and backing up data.

If you have any issues, please contact Boston Scientific Customer Service.

Note: Screens depicted in this manual may differ slightly from the screens on your Vercise Neural Navigator.

Intended Use

Vercise Neural Navigator is a software program that is used to set and adjust stimulation parameters for the Vercise PC and Vercise Gevia DBS Systems.

Connecting the Programming Wand to the Clinician Programmer

The Clinician Programmer (CP) communicates with the Stimulator via a Programming Wand (See Figure 1). The Programming Wand uses a radio frequency (RF) link to communicate with the Stimulator.

Figure 1. Clinician Programmer (CP) and Programming Wand

CAUTION: Use only Vercise PC or Vercise Gevia DBS System components with Vercise Neural Navigator. Failure to do so may result in the inability to program the Stimulator.


CAUTION: The CP is not equipment for the patient environment as defined by IEC 60601-1. The CP and the person using the CP should not be in contact with the patient while programming.

Connecting the Programming Wand to the CP:

1. Power ON the CP.
2. Log in as ClinicUser. You will be prompted to setup a password the first time you log into the CP.
3. Connect the Programming Wand to the CP using the USB cable provided with the Programming Wand. Plug the Mini USB end of the cable into the USB port on the side of the Programming Wand and the Standard USB end of the cable into the USB port on the CP.
4. The Wand will perform a short self-test. At the end of the self-test, the Wand will beep.
 - (a). If the Power Light on the Wand turns green, the Wand is ready for use.
 - (b). If the Power Light on the Wand remains red, contact a Boston Scientific representative, or call Customer Service. For more information on the Programming Wand refer to the *Vercise DBS Programming Wand Manual*.
5. Place the Wand over the Stimulator.

Starting a Programming Session

Starting Vercise Neural Navigator

1. Power ON the CP and log in.
2. Click the Vercise Launcher icon  on the desktop.
3. Select  to start Vercise Neural Navigator.

Note: Multiple softwares should not be run simultaneously on the same CP.

Note: The Vercise Neural Navigator can also be launched in Demo Mode using the Vercise Launcher. Demo Mode is used for demonstration purposes only.

Connecting to the Stimulator

Upon starting Vercise Neural Navigator, the screen will show the **Connect Tab** and the software will automatically attempt to connect to a Stimulator (See Figure 2). If no Stimulator is found, move the Wand closer to the Stimulator that you are trying to connect to and click the **Rescan** button.


Figure 2. Connect Tab

#	Feature	Description
1	Connect Tab	Displays the connection status between the CP, Wand, and Stimulator.
2	Configure Tab	Configure leads and edit patient profile.
3	Program Tab	Adjust the Stimulator program settings.
4	Report Tab	Generate, print and export reports.
5	Tools Tab	Enable new leads, backup or delete patient data.
6	End Session Tab	Disconnect from the Stimulator or exit the application.
7	Battery Indicator	Displays the battery status of the CP.

Note: The CP cannot connect to the Vercise Gevia™ Stimulator in MRI Mode. Exit MRI Mode using the Remote Control and rescan to connect. For instructions on exiting MRI Mode, refer to the Remote Control Manual as listed in your DBS Reference Guide.

Once a connection has been established between the CP and the Stimulator, the following screen will appear (See Figure 3).


Figure 3. Connection established between CP and Stimulator

#	Feature	Description												
1	Patient ID	Displays the Patient ID number												
2	Connection Status	Displays the connection status between the CP, Wand and Stimulator along with the Model and Serial Number of each of the devices <table border="1" data-bbox="280 905 1484 1140"> <thead> <tr> <th>#</th> <th>Connection Status</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td>1</td> <td></td> <td>Connection between CP, Wand and Stimulator was established.</td> </tr> <tr> <td>2</td> <td></td> <td>Connection between Wand and Stimulator not established. Move the Wand closer to the Stimulator and Rescan.</td> </tr> <tr> <td>3</td> <td></td> <td>Connection between CP, Wand and Stimulator not established. Check the USB cable connections to the Wand and refer to the Vercise DBS Programming Wand Manual for more information.</td> </tr> </tbody> </table>	#	Connection Status	Description	1		Connection between CP, Wand and Stimulator was established.	2		Connection between Wand and Stimulator not established. Move the Wand closer to the Stimulator and Rescan.	3		Connection between CP, Wand and Stimulator not established. Check the USB cable connections to the Wand and refer to the Vercise DBS Programming Wand Manual for more information.
#	Connection Status	Description												
1		Connection between CP, Wand and Stimulator was established.												
2		Connection between Wand and Stimulator not established. Move the Wand closer to the Stimulator and Rescan.												
3		Connection between CP, Wand and Stimulator not established. Check the USB cable connections to the Wand and refer to the Vercise DBS Programming Wand Manual for more information.												
3	Rescan Button	Scan for available Stimulators. Disabled if CP is already connected to the Stimulator.												
4	Connection Status	Displays a green filled circle if the Stimulator is connected to the CP.												
5	Connect or Disconnect Button	Connect or disconnect from a Stimulator. When a Stimulator is not connected, this button reads 'Connect'. When a Stimulator is connected, this button reads 'Disconnect'.												

Programming the Stimulator

Configuring the Leads

Once a connection has been established between the CP, Wand and Stimulator, switch to the **Configure Tab** to configure the leads that are connected to the Stimulator (See Figure 4).

Note: You only need to perform this step the first time you connect to a new Stimulator. Once a Stimulator has been configured, you can directly switch to the **Program Tab** after connection has been established.


Figure 4. Configuring Leads

#	Feature	Description
1	Leads Configuration	Select the leads, brain hemisphere and target area. Use the arrow keys to access the drop down menus and identify the implanted lead(s) and brain hemisphere (Right Side or Left Side) Type in the Target Area using an on-screen or attached keyboard. Note: If available, additional lead options may be added via the Tools Tab . See the “Tools Tab” section for more information.
2	Measure Button	Measure Impedances. See the “Measuring Impedances” section for more information.
3	Patient Amplitude Control	Turn ON/OFF the patient’s ability to change stimulation amplitude. The range of the amplitude is adjusted in the Program Tab .
4	Stimulator Information	Displays the Stimulator Information including Serial number, Model number, Firmware Version and type of Stimulator.
5	Implant Date	The Implant Date is defined by the date on which a CP first connects to a new Stimulator. You may adjust it by selecting the Implant Date button.
6	Patient ID	The Patient ID is defined as the Stimulator serial number by default. You may edit it by typing into the Patient ID field.

Measuring Impedances

Lead Impedances can be measured using the **Measure** button on the **Configure** or **Program Tabs**. The impedances of each of the 16 contacts may be checked to verify electrical integrity. When an impedance measurement is taken, impedances are assessed between a contact and the case (monopolar), and between pairs of contacts (bipolar). Impedances over 8000 Ω may be the result of open or unconnected wires and are displayed in yellow on the **Impedance Measurement** window. Impedances less than 200 Ω may be the result of short circuits and are displayed in orange. Contacts that have impedances outside the acceptable range are marked with a Ω symbol on the programming screen. The last set of impedance measurements are included in a report that can be printed or exported in the **Report Tab**.

The Programming Screen

Once the leads have been configured, select the **Program Tab** to begin programming. The Programming screen is divided into the following sections and features as seen in Figure 5. The Clinical Effects Map (CEM) shown in Figure 5 can be viewed only for a Standard Lead. Programming features exclusive to the Directional Lead can be seen in Figure 6. The Stimulation Field Model (SFM) shown in Figure 6 can be viewed for both a Standard Lead and Directional Lead.


Figure 5. Programming Screen


Figure 6. Directional Lead Programming Screen

#	Feature	Description
1	Program Button	Select the Program that you would like to set up or adjust.
2	Program Options Button	View battery longevity estimate for non-rechargeable Stimulators or recharge estimate for rechargeable Stimulators, delete and copy programs and change ramp and cycle times for Programs.
3	Area Panel	Select the Area within a Program that you would like to set up or adjust.
4	Delete Button	Delete an Area within a Program.
5	+	Add an area. Select from one of the ports defined in the Configure Tab .
6	Stimulation ON/OFF Button	Turn Stimulation OFF for the Area Selected. Note: When the Amplitude is at 0mA, increase the Amplitude to turn ON stimulation.

#	Feature	Description				
7	Pulse Width Buttons	Adjust the pulse width. <table border="1"> <thead> <tr> <th>Default</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>60 μS</td> <td>20 μS - 450 μS</td> </tr> </tbody> </table>	Default	Range	60 μ S	20 μ S - 450 μ S
Default	Range					
60 μ S	20 μ S - 450 μ S					
8	Rate Button	Adjust the rate <table border="1"> <thead> <tr> <th>Default</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>130 Hz</td> <td>2 Hz - 255 Hz</td> </tr> </tbody> </table>	Default	Range	130 Hz	2 Hz - 255 Hz
Default	Range					
130 Hz	2 Hz - 255 Hz					
9	Units Button	Select the Units in which Amplitude is displayed on contacts and case <table border="1"> <thead> <tr> <th>Default</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>%</td> <td>mA/%</td> </tr> </tbody> </table>	Default	Range	%	mA/%
Default	Range					
%	mA/%					
10	Patient Amplitude Buttons	Adjust Maximum and Minimum Patient Amplitude.				
11	Stop All Button	Turns off all Stimulation.				
12	Contact and Case Configuration	Displays percentage of anodic (+) or cathodic (-) energy assigned to the lead contacts and Stimulator Case for a given area. See the "Selecting Contacts" section for more information.				
13	Measure Button	Measures Impedances of the contacts.				
14	Programming Modes	Allows you to choose Steering or Custom Programming Modes				
15	Step Size	Select the step size for Amplitude adjustments: 0.1 mA or 0.5 mA				
16	Level Up and Down Buttons	Steer stimulation focus along the lead.				
17	Total Amplitude	Increase or decrease the total Amplitude delivered for a given Area.				
18	Clinical Effects Panel	Make note of the Therapeutic Benefits and/or Side Effects for the current Stimulation settings.				
19	Text Notes	Capture text notes for a given lead (up to 250 characters per lead port).				
20	Clinical Effects Map	Graphical summary of assigned Therapeutic Benefits and/or Side Effects at a given position along the DBS Lead array and a stimulation amplitude. Note: <i>Clinical effects data is captured and listed in reports but not plotted on the CEM for the DB2202-Directional Lead and for configurations not possible in Steering Mode.</i>				
21	Display Toggle	Toggle between viewing the CEM and the SFM. Note: <i>Only applies to DB2201-Standard Lead.</i>				
22	Stimulation Field Model	Visual representation of the estimated stimulation field for the currently programmed stimulation parameters. Note: <i>The SFM is not displayed for the DB2202-Directional Lead in Custom Mode.</i>				
23	Directional Presets	Select one of five preset directions for the stimulation field. The directional presets will steer the fully focused stimulation field in one of four orthogonal directions or put the stimulation field into "ring mode." Ring Mode generates, from a segmented contact level, stimulation fields equivalent to those generated by a standard "ring" or cylindrical contact. Note: <i>Only applies to DB2202-Directional Lead.</i>				
24	Rotate Buttons	Steer the stimulation focus circumferentially around the lead. Note: <i>Only applies to DB2202-Directional Lead.</i>				
25	Spread/Focus Buttons	Radially spread or focus the stimulation field. Note: <i>Only applies to DB2202-Directional Lead.</i>				
26	Directional Indicator	Visual indicator of the orientation of the radiopaque Directional Marker band on the Directional Lead. The orange line and dot correlate to the center of the radiopaque Directional Marker.				
27	SFM View Controls	Adjust the view of the SFM using Zoom (+ or -), Rotate, or Reset to the original view. Both the Lateral and Axial views of the SFM will adjust in unison using these controls.				

Creating or Modifying a Program

To create a new program or modify an existing program, select the Program Button and choose one of the four programs from the drop down menu. The system allows you to configure up to four Programs on a Stimulator.

For a given program, you can view and/or adjust several options via the Program Options button. The Program Options include:

#	Feature	Description						
1	Battery	<p>For a non-rechargeable (primary cell) Stimulator, the Energy Use Index for the current program is displayed. This value is used to provide an estimate of battery longevity for the current program on a new non-rechargeable Stimulator. See the “<i>Energy Use Index</i>” section for more information.</p> <p>For a rechargeable Stimulator, an Estimated Charge Time for the current program is displayed. This value provides an estimate of the duration and frequency of charging necessary to maintain stimulation.</p>						
2	Ramp	<p>Time to gradually increase the stimulation from zero to the programmed amplitude when Stimulation is turned ON.</p> <table border="1"> <thead> <tr> <th>Default</th> <th>Options</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>ON</td> <td>ON/OFF</td> <td>1 sec - 10 sec</td> </tr> </tbody> </table>	Default	Options	Range	ON	ON/OFF	1 sec - 10 sec
Default	Options	Range						
ON	ON/OFF	1 sec - 10 sec						
3	Cycle	<p>The cycled on and off duration of stimulation delivery.</p> <table border="1"> <thead> <tr> <th>Default</th> <th>Options</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>ON</td> <td>ON/OFF</td> <td>1 sec - 90 min</td> </tr> </tbody> </table>	Default	Options	Range	ON	ON/OFF	1 sec - 90 min
Default	Options	Range						
ON	ON/OFF	1 sec - 90 min						
4	Copy to	Copy the current program settings to another program.						
5	Delete	Delete the settings for the current program.						

Selecting Stimulation Areas

For a given Program, you may configure up to four areas. With a new Program, an Area will be automatically assigned to each Lead Port and named based on the defined Target and side of the brain selected in the **Configure Tab**. You may add an additional Area by selecting an empty Area (+) and choosing a Lead Port configuration (e.g. Left STN). You may reassign an Area by first selecting Delete to clear the setting for that Area and choosing a different Lead Port configuration.

Selecting Contacts

You can manually assign anodes and cathodes in **Custom Mode** or incrementally steer a stimulation field along the lead in the **Steering Mode**. Steering Mode is limited to a monopolar configuration of either a single cathode or adjacent cathodes. You may assign the case and all the contacts as anode or cathode individually in Custom Mode. The External Trial Stimulator (ETS) is limited to **Custom Mode** since the case cannot be assigned as an anode.

Note: Switching from Custom Mode to Steering Mode will clear the contact and case assignments.

Steering Mode

Steering Mode is an enhanced way of conducting a monopolar review where the contact(s) act as the cathode(s) and the stimulator case acts as the anode. This mode allows you to steer a monopolar cathode along the lead, eliminating the need to turn ON and OFF individual contacts. Steering Mode incrementally shifts a percentage of the cathodic current to the adjacent contact(s) using current steering technology to create smooth transitions between contacts.

The DB2201-Standard Lead has eight contacts per lead, labeled 1 through 8 for one lead and 9-16 for the second lead.

To Steer along the DB2201-Standard Lead:

1. Select **Steering Mode**.
2. Select a contact to assign it as a 100% cathode.
3. Use the  and  buttons to incrementally steer the stimulation focus along the length of the lead. The amount of cathodic current will shift in 10% increments.

Note: The amplitude for the selected Area will drop to 0 mA when another contact is selected.

The DB-2202 Directional Lead has a total of eight contacts per lead, labeled 1-8 for one lead and 9-16 for the second lead. Contacts 1 and 8 are the distal and proximal contacts while contacts 2-7 are the small directional contacts for one lead. Contacts 9 and 16 are the distal and proximal contacts while contacts 10-15 are the small directional contacts for the second lead.

To Steer along the DB2202-Directional lead:

1. Select **Steering Mode**.
2. Select a contact to assign it as a cathode.
3. Use the  and  buttons to incrementally steer the stimulation focus along the length of the lead.
4. Select one of five preset directions  for the stimulation field. The directional presets will steer the fully focused stimulation field in one of four orthogonal directions or put the stimulation field into “ring mode.”

The following steps can be used to refine the applied directional preset or as an alternative to the directional presets.

5. Use the  and  buttons to rotate and steer the stimulation focus circumferentially around the lead.
6. Use the  and  buttons to radially spread or shrink the focus of the stimulation field.
7. To choose another starting point or to steer on another contact, click another contact.

Note: *The total amplitude for the selected Area will drop to 0 mA when another contact is selected.*

Custom Mode

Custom Mode allows you to assign a percentage of anodic or cathodic current to individual contacts and the case.

To program the DB2201-Standard Lead and the DB2202-Directional Lead in Custom Mode:

1. Select **Custom Mode**.
2. Select the case or contact that you want to adjust. If it was blank, one tap will assign it as an anode (+). Another tap will reassign it as a cathode (-). Another tap will reassign it as OFF (blank). Tapping on a contact that already is assigned will first select it without changing the polarity.

Note: *Changing the contact polarities will reset the amplitude to zero.*

3. Select the + and – buttons for the Contact to adjust the percentage of anodic or cathodic current assigned to the selected Contact.

Note: *When using the External Trial Stimulator (ETS), monopolar configurations are not possible since the ETS “case” cannot be assigned as a cathode or anode.*

Note: *When using the ETS, Clinical Effects data is recorded but not plotted on the CEM.*

Turning Stimulation OFF for Individual Areas

To turn stimulation OFF for Individual Areas:

1. Make sure the Area you wish to turn OFF is selected by clicking on the appropriate Area on the Area Panel.
2. Press the Stimulation OFF button to turn Stimulation OFF.

Note: *When the Amplitude is at 0 mA, increase the Amplitude to turn ON stimulation.*

Turning All Stimulation OFF

Selecting the  button will stop stimulation for all active Areas. This function is only meant for turning all Stimulation OFF. To turn Stimulation ON, select each Area that you want to turn ON and select the Stimulation ON/OFF switch.

Increasing and Decreasing the Amplitude

Amplitude is measured in milliamperes (mA). The default setting for Amplitude is 0 mA and the range is 0 - 20 mA.

To increase or decrease the Amplitude:

1. Use the + and – buttons labeled Total Amplitude to increase or decrease Amplitude.
2. The default step size for Amplitude changes is 0.1 mA. You can change the step size to 0.5 mA using the Step Size buttons.

Note: *High stimulation levels can cause permanent tissue damage. A message will pop-up notifying you if you attempt to exceed a stimulation limit and settings that exceed this limit are not allowed.*

Increasing and Decreasing the Pulse Width

The Pulse Width of the stimulation is the length of time a burst of energy is applied per pulse. The Pulse Width is measured in microseconds (μS). The default setting for Pulse Width is 60 μS and the range is 20 - 450 μS .

Use the + and – buttons labeled Pulse Width to increase or decrease Pulse Width.

Note: High stimulation levels can cause permanent tissue damage. A message will pop-up notifying you if you attempt to exceed a stimulation limit and settings that exceed this limit are not allowed.

Increasing and Decreasing the Rate

The pulse rate of the stimulation, often called the Rate or the Frequency, dictates how many stimulation pulses are delivered in a second, measured in Hertz (Hz) or pulses per second (pps). The default setting for Rate is 130 Hz and the range is 2 - 255 Hz.

To increase or decrease the Rate:

1. Select the **Rate** button.
2. Select the desired Rate from the table of available rates. Incompatible rates are greyed out.

Note: Areas that are assigned to the same Lead Port cannot have rates that sum up to greater than 255 Hz.

Programming Multiple Areas with Different Rates

The Vercise PC and Vercise Gevia DBS Systems allow Areas to be programmed with different rates. By default, the Multiple Rate option is disabled. When you enable multiple rates, only the rates that are compatible with the rates and pulse widths from the other active Areas are available.

Note: Modifying the rate of an Area will alter the available rates for the other Areas.

Note: If you disable Multiple Rates, the rate for all Areas will be reset to the rate selected for the current Area.

Selecting Patient Amplitude Range

By default, patients do not have the ability to adjust the amplitude of their stimulation.

However, in some cases, you may want to give a patient the ability to adjust the Amplitude of their stimulation using the Remote Control. To give patients Amplitude Control, first turn ON the Patient Amplitude Control in the **Configure Tab**. Once the Patient Amplitude Control is turned ON, you can set the allowable Amplitude range in the **Program Tab** for each Area by setting a Minimum and Maximum.

Note: High stimulation levels can cause permanent tissue damage. A message will pop-up notifying you if you attempt to exceed a stimulation limit and settings that exceed this limit are not allowed.

Viewing the Stimulation Field Model

The Stimulation Field Model (SFM) is a visual representation of the estimated stimulation field for the currently programmed stimulation parameters. The SFM includes both a visual representation of the DBS lead as well as the approximated stimulation field shown in the color red (Figure 6). As programming parameters are adjusted and the stimulation is steered along the lead, the SFM will adjust accordingly.

Adjust the view of the SFM using  or  to Zoom,  or  to Rotate, or  to Reset the original view. Both the Lateral and Axial views of the SFM will adjust in unison using these controls. These controls will not affect or adjust any programming parameters.

Note: The SFM is displayed for the Standard Lead DB-2201 and the Directional Lead DB-2202. The SFM is not displayed for the the Directional Lead DB-2202 when in Custom Mode.

Mapping the Patient's Clinical Effects of Stimulation


Figure 7. Clinical Effects Map

For a given stimulation setting, you may make note of a 0-4 rating of potential therapeutic benefits and a 0-4 rating of possible side effects. Buttons labeled with potential symptoms and side effects are also available to select from. If not selected, one click will select a Therapeutic benefit and/or Side Effect and another click will remove the selection. Each button that is selected is captured as data associated with that stimulation setting for that patient.

You may also select the Text Notes button to enter and save up to 250 characters of text associated with each lead port.

When capturing Clinical Effects in Steering Mode, a dot is plotted on the CEM at the lead position and amplitude. The rating scale of the therapeutic benefit determines the color saturation of the center of dot. If a side effect is selected, an orange ring is displayed around the dot. Selecting a dot will display a pop-up window containing the date and time at which the dot was captured along with the stimulation setting and effects details (See Figure 7).

All of this data is saved on the Stimulator and available for export in the Reports Tab.

Note: Clinical effects data is captured and listed in reports but not plotted on the CEM for the DB2202-Directional Lead and configurations that are not possible in Steering Mode.

Reports

The Report Tab allows you to generate a report for the current programming session.

To generate a report for the current programming session, click on the **Report Tab** (See Figure 8). A Report may be printed and exported as a pdf or Excel file.


Figure 8. Report Tab

Select  on the **Report Tab** and select the desired information you want to include in the report by checking on any of the following check boxes:

- Programs
- Configuration
- Clinical Effects Maps
- Clinical Effects Details

You can also view reports for all Stimulators that were connected to the CP. Reports can be viewed when the CP is not connected to a Stimulator.

To view reports when the CP is not connected to a Stimulator (See Figure 9):

1. Click on the **Report Tab**
2. Select the patient whose report you would like to view and click on **View**.


Figure 9. Viewing Reports when CP is not connected to Stimulator

Tools and Additional Information

The **Tools Tab** allows you to Backup or Delete Patient Data and Enable New Leads.

Backup Database

The Backup/DataLog feature allows you to backup the entire patient database on the CP to a specified location. The backup location can be a folder on the CP or an external storage drive (for example, USB flash drive). This feature can be accessed from the **Tools Tab**.


Figure 10. Backup Database

To create a backup (See Figure 10):

1. Select the **Tools Tab**.
2. Select **Browse** to choose a backup location.
3. Click the **Backup** button to perform the backup function.

Note: After the backup is completed, a pop-up window will confirm the location of the file and indicate if the backup was successful.

Update Leads

This feature allows you to enable supported Leads within the Leads Configuration Panel. The Update Leads feature can be accessed from the **Tools Tab**. The DB2201-Standard Lead is enabled by default.

Note: You can only update Leads when disconnected from the Stimulator.


Figure 11. Update Leads

To enable the DB2202-Directional Lead (See Figure 11):

1. Disconnect from the patient's Stimulator by clicking on the **Disconnect** Button in the **Connect Tab**.
2. Select the **Tools Tab**.
3. Select the **Update Leads Tab**.
4. Select the row for the desired Lead.
5. Select the column labeled Key and enter the key 'KEVHSK'.
6. Select **Verify Key(s)**.

Deleting Clinical Effects Data

All the Clinical Effects Data for a patient can be deleted from the **Tools Tab**.

Note: This feature is available only when the CP is linked to a Patient's Stimulator.


Figure 12. Delete all Clinical Effects Data

To delete the Clinical Effects Data (See Figure 12):

1. Go to the **Tools Tab**.
2. Select the **Patient Data Tab** and click **Delete**.
3. Select **OK**.

Stimulator Programmable Characteristics

If two Leads are implanted, the stimulation parameters are independent such that stimulation of two different brain targets can have different amplitudes, pulse widths, stimulation rates and contact configurations. It is possible to configure one Lead as monopolar, and one as multipolar. It is also possible to configure a single Lead with both monopolar and multipolar areas.

The programmable parameter ranges for the Stimulator are shown below.

#	Parameter	Parameter Range
1	Amplitude ¹	0.1 mA - 20 mA
2	Rate ^{2,3}	2 Hz - 255 Hz
3	Pulse Width ⁴	20 μ s - 450 μ s
4	Cycle On/Off	1 sec - 90 minutes
5	Ramp On	1 - 10 seconds
6	Contact Connections	16
7	Independent Areas of Stim (4 programs with 4 areas per Program)	16

¹ The programmable coverage for each individual contact is limited to 12.7 mA. A programming interlock is enforced to limit the total output current to 20 mA or less per coverage area. For example, a maximum current output of 12.7 mA on one contact would limit the total summed current output on the remaining contacts to 7.3 mA within one coverage area.

² The rate is limited to 255 Hz for a given area.

³ The use of higher frequencies than the ones that are established (250 Hz) is the sole responsibility of the user.

⁴ The use of lower pulse widths than established (60-450 μ sec) is the sole responsibility of the user.

Charge Density


Figure 13. Charge Density Limits for Boston Scientific DBS Leads

Figure 13 displays the recommended maximum charge density for different combinations of Amplitude (mA) and Pulse Width (μs). The solid black line (Limit: 6mm^2) refers to all contacts on DB2201-Standard Lead and the proximal and distal contacts of DB2202-Directional Lead. The dashed black line (Limit: 1.5mm^2) refers to the small directional contacts of DB2202-Directional Lead. These estimates of charge density are only for the Boston Scientific DBS Leads.

WARNING: Patients may have the ability to change the Amplitude with the Remote Control. The physician should set and verify the maximum and minimum Amplitude levels allowed by the Remote Control to ensure that current levels remain safe.

Energy Use Index

The Energy Use Index applies to non-rechargeable, primary cell (PC) stimulators only. It does not apply to rechargeable stimulators. The Energy Use Index gives you an estimate of longevity of the battery life on the program selected. After the optimal settings have been identified for a program, from the **Program Tab**, click on **Program Options** and then click on **Battery** to obtain the Energy Use Index.

Use Figure 14 to identify the longevity that corresponds to this Energy Use Index. The figure takes into account nominal non-therapy power consumption, including shelf-life and patient remote control use. If the estimate for longevity obtained is below 12 months, consider evaluating a Boston Scientific rechargeable system.


Figure 14. Longevity Estimates Based On 24 Hour Per Day Usage

Estimated Charge Time

The Estimated Charge Time applies to rechargeable stimulators only. It does not apply to non-rechargeable, primary cell stimulators. The Estimated Charge Times provides an estimate of the duration and frequency of charging necessary to maintain stimulation for the selected program. After the optimal settings have been identified for a program, from the **Program Tab**, click on **Program Options** and then click on **Battery** to obtain the Estimated Charge Time.

Elective Replacement Indicator (ERI) Message

You will not be able to connect to a non-rechargeable (PC) Stimulator that is nearing the end of its battery life. The CP will display the Stimulator with an ERI message and the Stimulator battery voltage as seen in Figure 15 on the **Connect Tab**. No changes can be made to the Stimulator settings once it is in ERI.


Figure 15. ERI Message displayed on Connect Tab

End of Service (EOS) Message

When the Stimulator has reached its end of service, stimulation can no longer be provided. The CP will display the message as seen in Figure 16 on the **Connect Tab**.


Figure 16. EOS Message displayed on Connect Tab

Ending a Programming Session


The Patient's Remote Control automatically syncs with the Stimulator that it has been linked to. No additional steps need to be performed to copy programs from the Stimulator to the Remote Control. Refer to the appropriate Vercise Remote Control Manual, as listed in your *DBS Reference Guide*, for more information.

To end a Programming Session on the CP:

1. Select the **End Session Tab**.
2. Select **Exit Application** to end the programming session and close the application.
3. Select **Disconnect from Stimulator** to end the programming session and disconnect from the patient's stimulator. This will take you back to the **Connect Tab**.

Adjusting CP Time and Date

If system startup or hibernation is detected, the CP provides a notification to verify that the system time and date is correct.

If the Time and Date are correct, select **Verify** to dismiss the notification bar.

If the Time and Date are incorrect select **Adjust** to modify the time and date and click **OK** to confirm changes, as seen in Figure 17.


Figure 17. Date and Time Screens

Note: You can also click the Time and Date icon  on the Desktop to launch the time and date adjustment window.

Troubleshooting

There are no user serviceable parts. If you have a specific question or issue and need to contact Boston Scientific, choose your locality from the following list:

Argentina

T: +5411 4896 8556 F: +5411 4896 8550

Australia / New Zealand

T: 1800 676 133 F: 1800 836 666

Austria

T: +43 1 60 810 F: +43 1 60 810 60

Balkans

T: 0030 210 95 37 890 F: 0030 210 95 79 836

Belgium

T: 080094 494 F: 080093 343

Brazil

T: +55 11 5853 2244 F: +55 11 5853 2663

Bulgaria

T: +359 2 986 50 48 F: +359 2 986 57 09

Canada

T: +1 888 359 9691 F: +1 888 575 7396

Chile

T: +562 445 4904 F: +562 445 4915

China – Beijing

T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou

T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai

T: +86 21 6391 5600 F: +86 21 6391 5100

Colombia

T: +57 1 629 5045 F: +57 1 629 5082

Czech Republic

T: +420 2 3536 2911 F: +420 2 3536 4334

Denmark

T: 80 30 80 02 F: 80 30 80 05

Finland

T: 020 762 88 82 F: 020 762 88 83

France

T: +33(0)1 39 30 97 00 F: +33(0)1 39 30 97 99

Germany

T: 0800 072 3301 F: 0800 072 3319

Greece

T: +30 210 95 42401 F: +30 210 95 42420

Hong Kong

T: +852 2960 7100 F: +852 2563 5276

Hungary

T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore

T: +91 80 5112 1104/5 F: +91 80 5112 1106

India – Chennai

T: +91 44 2648 0318 F: +91 44 2641 4695

India – Delhi

T: +91 11 2618 0445/6 F: +91 11 2618 1024

India – Mumbai

T: +91 22 5677 8844 F: +91 22 2617 2783

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea

T: +82 2 3476 2121 F: +82 2 3476 1776

Malaysia

T: +60 3 7957 4266 F: +60 3 7957 4866

Mexico

T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East / Gulf / North Africa

T: +961 1 805 282 F: +961 1 805 445

The Netherlands

T: +31 30 602 5555 F: +31 30 602 5560

Norway

T: 800 104 04 F: 800 101 90

Philippines

T: +63 2 687 3239 F: +63 2 687 3047

Poland

T: +48 22 435 1414 F: +48 22 435 1410

Portugal

T: +351 21 3801243 F: +351 21 3801240

Singapore

T: +65 6418 8888 F: +65 6418 8899

South Africa

T: +27 11 840 8600 F: +27 11 463 6077

Spain

T: +34 901 11 12 15 F: +34 902 26 78 66

Sweden

T: 020 65 25 30 F: 020 55 25 35

Switzerland

T: 0800 826 786 F: 0800 826 787

Taiwan

T: +886 2 2747 7278 F: +886 2 2747 7270

Thailand

T: +66 2 2654 3810 F: +66 2 2654 3818

Turkey – Istanbul

T: +90 216 464 3666 F: +90 216 464 3677

Uruguay

T: +59 82 900 6212 F: +59 82 900 6212

UK & Eire

T: +44 844 800 4512 F: +44 844 800 4513

Venezuela

T: +58 212 959 8106 F: +58 212 959 5328

Note: Phone numbers and fax numbers may change. For the most current contact information, please refer to our website at <http://www.bostonscientific-international.com/> or write to the following address:

Boston Scientific Neuromodulation Corporation
25155 Rye Canyon Loop
Valencia, CA 91355 USA

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2018 Boston Scientific Corporation