**Boston
Scientific**

Advancing science for life™

Vercise Gevia™
Deep Brain Stimulation System
Information for Patients

92753242-20

Content: MP92753242-20 REV A

en

zh

zh-t

ko

How to Use this Manual

This manual provides information about the Vercise Gevia™ Deep Brain Stimulation (DBS) System. Read this manual carefully before using the Vercise Gevia DBS System. If you have any questions or concerns about your DBS System, contact your healthcare provider.

More Information

Refer to your *Remote Control Handbook* and/or *Charging Handbook* for:

- More instructions
- Compatibility among DBS System components
- Maintenance and monitoring using external devices
- More information

Refer to your *Labeling Symbols* document for an explanation of labeling symbols.

More information may be available online, visit www.bostonscientific.com/patientlabeling.

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

Trademarks

Vercise™, Vercise Gevia™, Cartesia™, SureTek™, ImageReady™, and Precision Spectra™ are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

Warranty

For device warranty information, visit www.bostonscientific.com/warranty.

Summary of Safety and Clinical Performance

For customers in the European Union, this information is available on the European database on medical devices (Eudamed) website: <https://ec.europa.eu/tools/eudamed>. Use the manufacturer and device name found in the labeling to search for the device's Summary of Safety and Clinical Performance.

Technical Support

There are no user serviceable parts. If you have a specific question or issue, contact your healthcare provider. To contact Boston Scientific for any other reason, use the contact information provided for your locality via www.bostonscientific.com.

Patient Implant Card

Your doctor should provide you with an Implant Card that identifies your particular implant. Carry your Implant Card with you at all times. Present it to all your healthcare providers (doctors, dentists, technicians) so they know you have an implanted device.

If you have not received your Implant Card, contact your healthcare provider.

Important Contact Information

For easy reference, note the contact information for your healthcare providers below:

Neurosurgeon _____

Neurologist _____

Caregiver _____

Information for Patients, Family Members, and Caregivers

- Boston Scientific recommends that you read this entire patient manual. It is unsafe to start using the device before reading the whole manual. If you have any questions about the information provided in this manual, contact your healthcare provider.
- **Always inform your healthcare providers that you have been implanted with a brain stimulation device.** If your healthcare providers have questions about your DBS System, they should contact Boston Scientific Technical Support. The contact list for your locality is provided via www.bostonscientific.com.
- If you have any questions or concerns about your DBS therapy, contact your healthcare provider. If there is an emergency, call 1-1-2 or local emergency services.

Table of Contents

Introduction.....	1
Description of the System	2
Rechargeable IPG Battery Information	4
Safety Information	6
Indications for Use	6
Patient Population	6
Intended Use	7
Clinical Benefit	7
Device Lifetime.....	7
Contraindications (When the Vercise Gevia DBS System Should Not Be Used).....	8
Warnings	9
Precautions	16
Adverse Events	25
Electromagnetic Compatibility.....	30
EN 60601-1-2 Classification Information.....	30
Notice to the User Per Industry Canada Radio Standards Specifications:	35

Essential Performance	35
Telemetry Information.....	35
Quality of Wireless Service	36
Troubleshooting Wireless Coexistence Issues.....	37
Wireless Security	38
System Components and Materials.....	39
Implantable Components	39
Materials.....	41
Glossary	43

Introduction

The Vercise Gevia DBS System provides therapy by stimulating structures in the brain with electrical pulses (Figure 1). The Vercise Gevia DBS System includes a rechargeable Implantable Pulse Generator (IPG).

Note: *Some patients may have their IPG implanted in the abdomen.*


Figure 1. Illustration of the Implanted DBS System (Pectoral Region)

Description of the System

The Vercise Gevia DBS System includes both implantable and non-implantable components. There are three main implantable components:

- **DBS Leads:** DBS Leads are thin, insulated wires that deliver electrical pulses to the brain. The body of the DBS Lead is implanted under the scalp. The other end of the Lead is connected to the DBS Lead Extension, typically behind the ear.
- **DBS Lead Extensions:** Lead Extensions are thin, insulated wires that connect the Leads to the IPG. One end of the Lead Extension is connected to the Lead, typically behind the ear. The other end of the Lead Extension is connected to the IPG. The Lead Extension transfers electrical stimulation from the IPG to the Lead.
- **Implantable Pulse Generator (IPG):** The IPG sends small electrical pulses to the end of the DBS Lead that is implanted in the brain. This produces stimulation in the brain. The IPG may be placed under the skin in the chest area or in the abdominal area. After your DBS surgery, your healthcare provider will adjust the stimulation settings for your IPG.

The IPG contains a rechargeable battery. The battery supplies power to your DBS System. You can charge the battery with the DBS Charging System.

The Vercise Gevia DBS System has two main non-implantable components:

- **Remote Control:** The Remote Control is a hand-held programmer that controls the IPG.
- **Charging System:** The Charging System charges the battery of rechargeable DBS IPGs.

To make the most of your Vercise Gevia DBS System, it is important to learn the following:

- How to live safely with your Vercise Gevia DBS System.
- How to use the Remote Control.
- How to use the Charging System to charge the IPG.

Note: *The Vercise Gevia DBS System components were not made with natural latex.*

Rechargeable IPG Battery Information

For information about charging your IPG, refer to the *Charging Handbook*. For information about checking the status of your IPG battery and battery messages, refer to the *Remote Control Handbook*.

IPG Battery

The rechargeable IPG battery should provide at least 5 years of service. In many cases, the IPG battery should provide at least 25 years of service. Battery life is dependent on the stimulation settings and conditions.

Charging Estimate

Boston Scientific recommends any charging schedule that fits your schedule and lifestyle while maintaining sufficient charge to maintain stimulation. Expect a daily charging time of 5 to 30 minutes, or a periodic charging time of 30 minutes to 4 hours every 1 to 2 weeks. Your charging routine may vary depending on your stimulation settings. High power users will require more frequent charging. Creating a charging routine involves finding the right balance among the following:

- How much power is needed for effective therapy
- How often you want to charge
- How long you want to charge
- How you would like to manage your personal schedule

The Vercise Gevia DBS System programming software gives your healthcare provider a conservative estimate for how often to charge. This estimate assumes stimulation is on the default stimulation level 24 hours per day, 7 days per week. While you may want to follow these estimates, you and your healthcare provider can create a charging routine that best fits your schedule.

If you do not charge your IPG before it enters Hibernation Mode, stimulation will stop until you charge the IPG. Create a charging routine to prevent loss of stimulation due to a low battery.

After years of service, the IPG may require shorter periods between charges. The IPG will need replacement when stimulation cannot be maintained with routine charging.

Safety Information

Indications for Use

The Vercise Gevia DBS System is indicated for use in the following:

- Unilateral or bilateral stimulation of the subthalamic nucleus (STN) or internal globus pallidus (GPi) for treatment of levodopa-responsive Parkinson's disease that is not adequately controlled with medication, for persons 18 years of age and older.
- Unilateral or bilateral stimulation of the subthalamic nucleus (STN) or internal globus pallidus (GPi) for treatment of intractable primary and secondary dystonia, for persons 7 years of age and older.
- Thalamic stimulation for the suppression of tremor that is not adequately controlled by medications in patients diagnosed with Essential Tremor or Parkinson's disease, for persons 18 years of age and older.

Patient Population

Target groups of patients for DBS therapy are those who have:

- Parkinson's disease (PD) – levodopa-responsive and not adequately controlled with medication in patients 18 years of age and older
- Dystonia – primary or secondary intractable dystonia in patients 7 years of age and older
- Tremor – in PD or essential tremor (ET) inadequately controlled by medication in patients 18 years of age and older

Intended Use

The intended use of the Boston Scientific DBS System is to treat movement disorders by stimulating a specific brain target.

Clinical Benefit

The clinical benefit of the Boston Scientific DBS System is reduction of motor symptoms of the disease state (Parkinson's disease, dystonia, tremor). Other benefits known to occur from motor symptom reductions are increased functionality and quality of life. Patients often have decreased use of medications and therefore fewer side effects from medication.

Device Lifetime

The minimum lifetime of the implantable components is 5 years.

Contraindications (When the Vercise Gevia DBS System Should Not Be Used)

The Vercise Gevia DBS System is not recommended for patients that have the following conditions or who will be exposed to the following procedures:

- **Diathermy:** You should not have any form of diathermy either as treatment for a medical condition or as part of a surgical procedure. The energy generated by diathermy can be transferred to the Vercise Gevia DBS System, causing tissue damage in the brain, which can result in severe injury or death. The use of shortwave, microwave, and/or therapeutic ultrasound diathermy with the implanted Vercise Gevia DBS System may result in severe injury or death.
- **Patient Incapability:** If you are unable to properly operate the Remote Control and Charging System, then you should not be implanted with the Vercise Gevia DBS System.
- **Poor Surgical Candidates:** The Vercise Gevia DBS System is not recommended for patients who are poor surgical candidates.

Warnings

Automobiles and Equipment

Operate the following with caution after receiving the Vercise Gevia DBS System:

- Automobiles
- Other motorized vehicles
- Potentially dangerous machinery or equipment

Avoid activities that would be dangerous if treated symptoms were to return. Avoid actions that cause stimulation changes to occur. Impaired driving and an increased accident risk have been previously reported for patients with Parkinson's disease.

You are at an increased risk of causing injury or death to yourself and others if your DBS System stops treatment for any reason during the following:

- Driving a car.
- Driving any other motorized vehicle.
- Using potentially dangerous machinery or equipment.

Electromagnetic Interference

Strong electromagnetic fields can potentially:

- Turn stimulation off.
- Cause temporary unpredictable changes in stimulation.
- Interfere with the communication of the Remote Control.

Strong electromagnetic fields may turn stimulation off for a short period of time. Once the electromagnetic field is removed, stimulation may return on its own.

You should avoid or exercise care around the following:

- Theft detectors, tag deactivators, and Radio Frequency Identification (RFID) devices, such as those used at department stores, libraries, and other public establishments.
 - If you must proceed through the detector, proceed with caution and move through the center of the detector as quickly as possible. Interference from these devices should not cause permanent damage to the implanted device.
- Security screeners, such as those used in Airport Security or at entrances to government buildings, including hand-held scanners.
 - Request assistance to bypass the security screener and advise the security staff that you have an implanted medical device. If you must proceed through the security screener, move through the security screener quickly and stay as far as allowed from the screener. Interference from these devices should not cause permanent damage to the implanted device.

- Power lines or power generators.
- Electric steel furnaces and arc welders.
- Large magnetized stereo speakers.
- Strong magnets.
- Automobiles or other motorized vehicles using a LoJack system or other anti-theft systems that can broadcast a radio frequency (RF) signal. The high energy fields produced by these systems may interfere with the function of the Remote Control and its ability to control stimulation.
- Other sources of electromagnetic disturbance, such as RF transmitters at:
 - Television broadcast stations
 - Radio broadcast stations
 - Amateur Radio radio transceivers
 - Citizens Band radio transceivers
 - Family Radio Service band transceivers

Note: *When in close proximity, equipment that generates strong electromagnetic fields might cause unintended stimulation or interfere with wireless communication between the Remote Control and IPG. This may occur even if the equipment complies with International Special Committee on Radio Interference (CISPR) requirements.*

Heat Due to Charging

Handle the Charger with care. The Charger may become warm while charging your IPG.

- Use the Charger with the Charging Collar, Charger Belt, or Adhesive Patch to avoid a burn.
- Do not charge your device while sleeping. This may result in a burn.
- Stop charging and contact your healthcare provider if you feel pain or discomfort.

High Stimulation Levels

High levels of stimulation may damage brain tissue. Your healthcare provider programs your Remote Control stimulation settings to ensure that your stimulation levels remain safe.

Intracranial Hemorrhage

Placing the DBS Leads in the brain may increase the risk of bleeding in or around the brain. Notify your healthcare provider if:

- You are more prone to bleeding.
- You have trouble forming blood clots.
- You take medication that thins the blood.

These factors may increase your risk of bleeding in or around the brain.

IPG Damage

Chemical burns may result if the IPG is ruptured or pierced and your tissue is exposed to battery chemicals.

Magnetic Resonance Imaging

As a patient implanted with the Vercise Gevia DBS System, you can have an MRI exam when specific conditions are met. Conditions for an MRI are specified in the physician manual *ImageReady™ MRI Guidelines for Boston Scientific DBS Systems*. This manual is available at www.IFU-BSCI.com.

Your healthcare provider should read this information in its entirety and determine all conditions are met before conducting and recommending an MRI exam.

External Devices:  The external/non-implantable components of the Vercise Gevia DBS System (Charging System, Remote Control, and accessories) are MR Unsafe. Do not take them into any MRI environment.

Other Implanted Stimulation Devices

Contact your healthcare provider if you must use other implantable stimulation devices at the same time as the DBS IPG. Using other implantable stimulation devices at the same time as the DBS IPG may interfere with the function of the devices. Careful programming of each system is needed. Examples of other implantable stimulation devices are:

- Pacemakers
- Cardioverter defibrillators
- Medication delivery pumps

Pregnancy

It is unknown whether the DBS device may cause complications with pregnancy.

Status Dystonicus (Dystonia Patients)

If dystonia symptoms return or become more severe, contact your healthcare provider right away.

DBS patients with dystonia may rarely experience a life-threatening worsening of dystonia symptoms. These symptoms may be referred to as “status dystonicus,” “dystonic storm,” or “dystonic crisis.” This condition may be triggered by several factors, such as the loss of stimulation from your DBS device. This condition is characterized by severe dystonia associated with muscle breakdown, organ failure, and possibly, death.

It is important that you monitor your symptoms and maintain DBS as prescribed. If you suspect that your DBS System is no longer providing adequate stimulation for any reason, contact your healthcare provider right away. If your dystonia symptoms worsen, contact your healthcare provider right away.

You should know how to use your Remote Control to ensure that your DBS System is on and providing stimulation. You should know how to charge your IPG.

Suicide

New onset or worsening depression which may be temporary or permanent is a risk that has been reported with DBS therapy. Suicidal thoughts, suicide attempts, and suicide are events that have also been reported. Patients and caregivers should consider the following:

- Be sure you talk to your treating healthcare provider(s) if you have a history of depression, suicidal thoughts, or have attempted suicide. Be sure you understand the possible risks of new onset or worsening depression (including suicidal thoughts) as well as the potential clinical benefits of DBS therapy.
- After the procedure, if you notice unusual changes in mood or behavior (such as increased anxiety, sleeping problems, loss of interest in activities, feeling of hopelessness, mood swings, weight loss or weight gain), or impulse control, contact your healthcare provider. If you are having thoughts of suicide, contact your healthcare provider or emergency services immediately.
- It is important to attend on going follow-up visits with your healthcare provider to manage your therapy.

Therapeutic Ultrasound

Do not expose the implanted components of the DBS System to therapeutic ultrasound.

Unauthorized Modification

Unauthorized changes to the medical devices are not allowed. If unauthorized changes are made to the devices, the integrity of the devices could be at risk and you may be harmed or injured.

Precautions

Activities Requiring Coordination

Loss of coordination is a possible side effect of DBS therapy. You should use caution when doing activities that require coordination, even if you were able to do them before receiving therapy (for example, swimming).

Bathing

Patients should exercise reasonable caution when bathing.

Cell Phones and Other Portable RF Communication Devices

Interference caused by cell phones is not expected. The full effects of interaction with cell phones are unknown at this time. Portable RF communications equipment (for example, cell phones) should be kept at least 6 in (15 cm) from the area of the implanted device. If interference occurs, move the cell phone away from the implanted IPG or turn off the phone. If there is a concern or problem, contact your healthcare provider.

Cleaning the Remote Control and Charging System

The Remote Control and Charging System components can be cleaned using alcohol or a mild detergent applied with a cloth or tissue. Remove residue from soapy detergents with a damp cloth and then wipe dry. Do not use abrasive cleansers for cleaning. Do not clean any of the devices while they are directly or indirectly connected to a power outlet.

Remove the Charger, Counterweight, and, if applicable, Charger Spacer from the Charging Collar before washing. Remove the Charger and, if applicable, Charger Spacer from the Charger Belt before washing. Hand wash the Charging Collar or Charger Belt with mild soap and cool water. Do not machine wash the Charging Collar or Charger Belt. Let the Charging Collar or Charger Belt air dry. Do not use the Charging Collar or Charger Belt when it is damp or wet.

As an operator of these external devices, only perform the following service and maintenance tasks on the external devices:

- Charging the battery
- Cleaning

Ensure that the devices are not in use while performing service and maintenance tasks.

Component Disposal

DBS components removed from the body should be returned to Boston Scientific. The IPG should be removed from the body prior to cremation and returned to Boston Scientific. Cremation may cause the IPG battery to explode.

The Remote Control and Charging System contain batteries which may explode when exposed to fire and cause injury. Do not dispose of the Remote Control or Charging System by fire. Dispose of used batteries in accordance with local laws and regulations.

Delayed Effectiveness

It may take some time before DBS therapy becomes effective. To achieve the best symptom control, you may need frequent adjustment of stimulation settings over a period of several weeks to months.

Device Connections

Do not connect the Remote Control USB port to any other device except to the Remote Control charger provided with it. Do not connect the Charger to any other device except to the Base Station provided with it.

Device Failure

Implants can fail at any time due to:

- Random component failure
- Loss of battery function
- Lead breakage
- Lead migration

Stopping stimulation suddenly can cause serious reactions to develop. If the IPG stops working, even after complete charging, turn off stimulation with your Remote Control and contact your healthcare provider right away. Your healthcare provider will assess the DBS System and provide proper medical care to manage the return of symptoms.

IPG Orientation

Never attempt to change the orientation of the IPG or turn over the IPG. Avoid touching the incisions or IPG site. If the IPG flips over in your body, it may be unable to communicate with the Remote Control or Clinician Programmer. If the IPG flips over in your body, it cannot be charged. If stimulation cannot be turned on after charging, contact your healthcare provider to arrange an exam of the System.

If you notice a change in appearance of the skin at the IPG location, such as the skin becoming thin over time, contact your healthcare provider.

Massage Therapy

Avoid massage therapy near the implanted components of the DBS System. If you receive massage therapy, inform the massage therapist that you have an implanted device. Show the massage therapist where the IPG, DBS Lead Extension, and DBS Leads are located. Have the massage therapist avoid these areas and proceed with caution.

Medical Devices/Therapies

The following medical therapies or procedures may:

- Turn stimulation off.
- Cause permanent damage to the IPG.
- Cause you injury, especially if they are used near the IPG.
- Diagnostic ultrasonic scanning is unlikely to damage the IPG if stimulation is turned off. Testing has been completed to pertinent standards.
- Electrocautery – The use of a heated electric probe to stop bleeding during surgery. This can potentially damage the IPG or may cause injury. If you will have this procedure, your healthcare provider should review and follow the guidance provided in the *Information for Prescribers* manual.
- External Defibrillation – The use of electrically charged paddles to restart the heart in an emergency. Safe use has not been established in DBS patients. Testing has been completed to pertinent standards. Defibrillation is unlikely to permanently damage the implanted device if:
 - Stimulation is turned off.
 - The defibrillator electrode does not contact any component of the implanted device.

- Lithotripsy – High-output sound or shock waves often used to treat gall stones and kidney stones. High frequency signals directed near the IPG may damage circuitry.
- Radiation Therapy – Ionizing energy commonly used to treat cancer. Lead shielding should be used over the IPG to prevent damage from high radiation. Any damage to the DBS IPG by radiation may not be known right away.
- Transcranial Stimulation – The use of electric or magnetic fields to stimulate nerve cells in the brain. Safe use of these therapies in DBS patients has not been established.
- X-ray and CT scans may damage the IPG if stimulation is on. X-ray and CT scans are unlikely to damage the DBS IPG if stimulation is turned off.

If any of the above is needed, the procedure(s) should be performed as far from the implanted components as possible. The DBS System may need to be removed as a result of severe injury or damage to the device.

Before having these procedures, medical therapies, or diagnostics, have your healthcare provider call the Boston Scientific Technical Support department for proper instructions. Or have them refer to the guidance provided in the *Information for Prescribers* manual.

Operating Temperature

The operating temperature of the Remote Control is 5 °C to 40 °C (41 °F to 104 °F). Do not use the Charging System if the ambient temperature is above 35 °C (95 °F).

Other Models of External Devices

Only use the Remote Control and Charger provided with the Vercise Gevia DBS System. Other similar models of the Remote Control and Charger do not work with the Vercise Gevia DBS System.

Pediatric Patients

Caregivers should know how to use the Remote Control to ensure the system is on and providing stimulation. Caregivers must know how to charge the IPG. If symptoms return, the caregiver should contact the healthcare provider right away.

Post-Operative

After your surgery, the medical staff provides standard medical care:

- A CT scan may be taken to record the position of the DBS Leads and IPG.
- You and a family member are educated on the DBS System, including:
 - How to turn stimulation on and off.
 - How to charge the IPG battery.
 - Expectations of stimulation in the treatment of your disease.
- Antibiotics may be prescribed to prevent infection.
- Post-operative pain concerns are addressed by your healthcare provider before you are discharged from the hospital.
- A responsible adult that can fully understand the post-operative care instructions must drive you home after the surgery.

After surgery, use care so that healing secures the implanted components and closes the surgical incision:

- Restrict head movements, as instructed by your healthcare provider, until healed. This may include turning or moving your head up, down, and side to side.
- Do not try to move heavy objects.
- Do not shower until cleared by your healthcare provider. Surgical sutures and staples need to be removed by your healthcare provider in a follow-up visit.
- Follow the instructions provided by your healthcare provider on how to care for the dressing covering the area where the IPG is implanted.

Stimulation therapy may be delayed up to 2 months until swelling resolves. Therapy timing depends on the judgment of your healthcare provider. Expect post-surgical swelling to subside. If swelling is present at the IPG site once stimulation therapy has begun, swelling may lead to longer charging times or prevent charging. If swelling lasts, contact your healthcare provider.

There may be some pain in the area of the IPG until healing is complete. If pain lasts beyond 2 weeks, contact your healthcare provider.

If you notice excess redness or drainage around the wound areas, contact your healthcare provider. In rare cases, adverse tissue reaction to implanted materials can occur.

Storage, Handling, and Transport

Do not expose the Remote Control or Charging System to very hot or cold conditions. Do not leave these devices in your car or outdoors for long periods of time. The sensitive electronics can be damaged by temperature extremes, especially high heat. If the Remote Control or the Charging System are stored for a long period of time, ensure the storage temperature does not exceed -20 °C to 60 °C (-4 °F to 140 °F).

Handle the System components and accessories with care. Do not drop them or submerge them in water. Keep accessories, including the Remote Control and Charger, dry and protected from moisture. Keep the Remote Control and Charger away from pets, pests, and children to avoid damage to these devices. Reliability testing has been performed to ensure quality manufacturing and performance. However, the devices can be permanently damaged by:

- Dropping on hard surfaces
- Dropping in water
- Rough handling

Adverse Events

Not all risks with the use of deep brain stimulation are known. The following is a list of known risks. If any of these events occur, inform your healthcare provider as soon as possible. Report any serious incident that occurs in relation to this device to Boston Scientific and to the relevant local regulatory authority for medical devices in your country. For customers in Australia, report any serious incident that occurs in relation to this device to Boston Scientific and to the Therapeutic Goods Administration (<https://www.tga.gov.au>).

Risks Associated with Surgical Procedure and Post-Operative Period

- Blood clot forming, for example in the veins of the legs
- Blood clot or air in the blood stream, which can block blood flow to parts of the lungs or other tissue
- Cerebral spinal fluid (CSF) that normally surrounds the brain leaking outside the skull
- Confusion or problems with attention, thinking, or memory
- Death
- Infection
- Injury to areas next to the implant, such as blood vessels, nerves, the chest wall, and the brain
- Neurosurgery and anesthesia risks, including unsuccessful implant and pneumonia
- Pain, headache, or discomfort
- Seizures
- Status dystonicus
- Stroke, resulting in temporary or permanent problems, such as paralysis or difficulty speaking
- Swelling or bruising

Possible Side-Effects

- Confusion or problems with attention, thinking, or memory
- Pain, headache, or discomfort
- Psychiatric disturbances, such as:
 - Anxiety
 - Depression
 - Lessened interest or emotion
 - Mania
 - Emotional sensitivity
 - Sleep problems
 - Suicide
 - Suicidal thoughts or attempts
- Seizures
- Sensory changes
- Speech or language problems, difficulty swallowing, as well as complications such as pneumonia from inhaling fluid
- Systemic effects, such as:
 - Rapid heart beat
 - Sweating
 - Fever
 - Dizziness

- Changes in kidney function
- Difficulty passing urine
- Sexual effects
- Nausea
- Difficulty having bowel movements
- Bloating
- Weakness, muscle spasms, shaking, restlessness, or problems with movement, walking, coordination, or body positioning, and falls or injuries resulting from these problems
- Undesirable sensations, such as tingling
- Visual problems, eyelid or eye movement difficulties, or other eye-related symptoms
- Weight changes

Device-Related Risks

- Allergic or immune system response to implanted materials
- Failure or malfunction of any part of the device. This may or may not require device removal and/or replacement. This includes, but is not limited to:
 - Battery leakage
 - Battery failure
 - Lead or Lead Extension breakage
 - Hardware malfunctions
 - Loose connections
 - Electrical shorts or open circuits
 - Lead insulation breaches
- Implant site complications, such as:
 - Pain
 - Poor healing
 - Redness
 - Warmth
 - Swelling
 - Wound reopening
- Implanted device components (IPG, Lead, or Lead Extension) may move from the original implanted location or wear through the skin, which may require additional surgery

- Infection
- Interference from external electromagnetic sources
- Loss of adequate stimulation
- During an MRI exam, there are potential interactions with the implanted DBS Lead, Lead Extension, and IPG, and risk of patient harm. It is important that your healthcare provider follow the information provided in the *ImageReady™ MRI Guidelines for Boston Scientific DBS Systems* before giving an MRI exam. It is available at www.IFU-BSCI.com.
- Pain, headache, or discomfort
- Skin irritation or burns at the IPG site
- Status dystonicus
- Stiffness in muscles or with joint movement
- Sudden return of symptoms, if the stimulation or medication is stopped abruptly. In patients with Parkinson's disease, there have been rare cases of rapid symptom return progressing to inability to move.
- Swelling, including fluid collecting around the device

Electromagnetic Compatibility

EN 60601-1-2 Classification Information

- Internally Powered Equipment
- Continuous Operation
- Ordinary Equipment
- Class II

Table 1: Guidance and Manufacturer's Declaration Electromagnetic Emissions		
The Vercise Gevia DBS System is intended for use in the electromagnetic environment described below. The customer or user of the Vercise Gevia DBS System should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The Vercise Gevia DBS System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Vercise Gevia DBS System is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Table 2: Guidance and Manufacturer's Declaration
Electromagnetic Immunity**

The Vercise Gevia DBS System is intended for use in the electromagnetic environment described below. The customer or user of the Vercise Gevia DBS System should ensure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Air: $\pm 2 \text{ kV}, \pm 4 \text{ kV}, \pm 8 \text{ kV}, \pm 15 \text{ kV}$ Contact: $\pm 8 \text{ kV}$	Air: Remote Control and Charger: $\pm 2 \text{ kV}, \pm 4 \text{ kV}, \pm 8 \text{ kV}, \pm 15 \text{ kV}$ Contact: Remote Control and Charger: $\pm 8 \text{ kV}$	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Note: <i>Applies to external electrical components.</i>
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. Magnetic fields from common appliances are not expected to affect the device.

**Table 3: Guidance and Manufacturer's Declaration
Electromagnetic Immunity**

The Vercise Gevia DBS System is intended for use in the electromagnetic environment described below. The customer or user of the Vercise Gevia DBS System should ensure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m 80 MHz to 2,7 GHz	Professional healthcare facility environment and home healthcare environment. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the symbol shown below: 

Note: These guidelines may not apply in all situations. Absorption and reflection from structures, objects, and people affects electromagnetic propagation.

a Field strengths from fixed transmitters cannot be predicted theoretically with accuracy. This includes base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Vercise Gevia DBS System is used exceeds the applicable RF compliance level above, the Vercise Gevia DBS System should be observed to verify normal operation. If abnormal performance is observed, more measures may be needed, such as re-orienting or relocating the Vercise Gevia DBS System.

**Table 4: Immunity Testing
RFID Readers**

The external electrical components of the Vercise Gevia DBS System have been tested for immunity to interference from RFID readers per the following specifications.

RFID Spec Per AIM 7351731	Frequency	Test Level (RMS)
ISO 14223	134.2 kHz	65 A/m
ISO/IEC 14443-3 (Type A)	13.56 MHz	7.5 A/m
ISO/IEC 14443-4 (Type B)	13.56 MHz	7.5 A/m
ISOAEC 15693 (ISO 18000-3 Mode 1)	13.56 MHz	5 A/m
ISO 18000-3 Mode 3	13.56 MHz	12 A/m
ISO/IEC 18000-7	433 MHz	3 V/m
ISO/IEC 18000-63 Type C	860-960 MHz	54 V/m
ISO/IEC 18000-4 Mode 1	2.45 GHz	54 V/m

**Table 5: Manufacturer's Declaration
Proximity Fields**

The Vercise Gevia DBS System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The users of the Vercise Gevia DBS System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Vercise Gevia DBS System as recommended below, according to the maximum output power of the communications equipment.

Proximity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment Guide
IEC 61000-4-3	385 MHz: 27 V/m @ 18 Hz pulse modulation	27 V/m	Recommended separation distance d = 30 cm
	450 MHz: 28 V/m @ FM modulation	28 V/m	
	710 MHz, 745 MHz, 780 MHz: 9 V/m @ 217 Hz pulse modulation	9 V/m	
	1720 MHz, 1845 MHz, 1970 MHz: 28 V/m @ 217 Hz pulse modulation	28 V/m	
IEC 61000-4-3	1720 MHz, 1845 MHz, 1970 MHz: 28 V/m @ 217 Hz pulse modulation	28 V/m	Recommended separation distance d = 30 cm
	2450 MHz: 28 V/m @ 217 Hz pulse modulation	28 V/m	
	5240 MHz, 5500 MHz, 5785 MHz: 9 V/m @ 217 Hz pulse modulation	9 V/m	

Note: For the frequency bands in this table, use the recommended separation distance. The recommended minimum separation distance of 30 cm between portable and mobile RF communications equipment (transmitters) and the Vercise Gevia DBS System apply to all other frequencies within the stated ranges.

Note: These guidelines may not apply in all situations. Absorption and reflection from structures, objects, and people affects electromagnetic propagation.

Notice to the User Per Industry Canada Radio Standards Specifications:

The Boston Scientific DBS devices comply with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

1. The device may not cause interference.
2. The device must accept any interference, including interference that may cause undesired operation of the device.

Essential Performance

Failure of the external electrical components will not result in an unacceptable risk to the user.

Telemetry Information

The following values describe the wireless communication link between the IPG and Remote Control:

- Frequency Band: 119 kHz to 131 kHz
- Modulation Type: FSK
- Effective Radiated Power: 0.05 mW (-13 dBm) maximum
- Magnetic Field Strength (at 3 m distance): 46 μ A/m

Quality of Wireless Service

The Vercise Gevia DBS System uses a Half-Duplex, direct point-to-point, primary-secondary communication system with the characteristics listed in Table 6:

Table 6: Quality of Wireless Service of the Configure Tab
Typical Range Between the Remote Control and Implanted IPG
36 in (91.4 cm)

Timing

Once a command is initiated by the user, the system responds in less than 1.5 seconds.

Troubleshooting Wireless Coexistence Issues

Other wireless and RF technology based equipment running close to a similar frequency band may degrade the range and responsiveness of the Vercise Gevia DBS System. If you experience issues with the wireless communication behavior between the Remote Control and IPG, try the following steps to correct the behavior:

- Decrease the distance between the two devices if possible.
- Change the orientation of the communicating device.
- Move the communicating devices away from other equipment or devices that may cause interference, such as:
 - Television monitors
 - Computer monitors
 - Short range RFID electronic tracking systems, such as:
 - Badge scanners
 - Parking lot scanners
 - Power adapters
 - Wireless chargers

Wireless Security

The Vercise Gevia DBS System has a short range inductively coupled telemetry system. A Remote Control has to be linked with an IPG to allow communication. The IPG will not respond to any unlinked device. There are additional mechanisms that ensure the integrity of the communicated data.

System Components and Materials

Implantable Components

Your Boston Scientific DBS System may include the following implantable components:

Table 7: Implantable Components		
Device Type	Model Number	Description
<ul style="list-style-type: none"> • IPG • Port Plug 	DB-1200	Vercise Gevia™ Implantable Pulse Generator Kit, 16 Contact
<ul style="list-style-type: none"> • Lead • Lead Boot • Suture Sleeve 	DB-2201-30DC	Lead Kit, 8 Contact, 30cm
	DB-2201-45DC	Lead Kit, 8 Contact, 45cm
	DB-2202-30	Vercise™ Cartesia™ Directional Lead Kit, 8 Contact, 30cm
	DB-2202-45	Vercise™ Cartesia™ Directional Lead Kit, 8 Contact, 45cm
	DB-2500-C	Vercise™ Physician's Spare Kit, 8 Contact Lead
<ul style="list-style-type: none"> • Lead Extension • Port Plug 	NM-3138-55	Lead Extension Kit, 8 Contact, 55cm
	DB-3128-55	Lead Extension Kit, 2x8 Contact, 55cm
	DB-3128-95	Lead Extension Kit, 2x8 Contact, 95cm
<ul style="list-style-type: none"> • Burr Hole Cover 	DB-4600-C	SureTek™ Burr Hole Cover Kit
	DB-4605-C	SureTek™ Burr Hole Cover Spares Kit

Table 7: Implantable Components

Device Type	Model Number	Description
• Adapter	DB-9208-15	Vercise™ Adapter S8, 8 Contact, 15cm
	DB-9208-55	Vercise™ Adapter S8, 8 Contact, 55cm
	DB-9218-15	Vercise™ Adapter M8, 8 Contact, 15cm
	DB-9218-55	Vercise™ Adapter M8, 8 Contact, 55cm
• Medical Adhesive	SC-4320	Medical Adhesive
• Port Plug	SC-4401	Precision Spectra™ IPG Port Plug

Materials

The following materials and substances may contact the body during the DBS implant procedure.

Table 8: Materials		
Device Type	Materials	Percent Weight
IPG	Titanium	18% to 57%
	Epoxy	42% to 82%
	Silicone	Up to 3%
IPG Template	Stainless Steel	100%
Port Plug	Thermoplastic	100%
Lead	Thermoplastic	69% to 76%
	Platinum - Iridium	7% to 23%
	Epoxy	0% to 22%
8 Contact Extension	Silicone	60% to 66%
	Thermoplastic	34% to 40%
	Epoxy	< 1%
2x8 Contact Extension	Silicone	14% to 39%
	Thermoplastic	61% to 85%
	Gold	< 1%
Suture Sleeve	Silicone	100%

Lead Boot	Silicone	100%
Burr Hole Base, Cap, and Lead Clip	Thermoplastic	100%
Burr Hole Cover Screw	Titanium	100%
Tunneling Tool Shaft	Stainless Steel	100%
Tunneling Tool Straw	PTFE	100%
Medical Adhesive	Medical Adhesive Silicone	100%
Adapter	Silicone	37% to 74%
	Thermoplastic	26% to 63%
	Epoxy	< 1%
<p>The IPG Template and Tunneling Tool may contain cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.</p> <p>Note: <i>The IPG Template and Tunneling Tool are made with stainless steel which may contain cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.</i></p>		

Glossary

Adhesive Patch. A non-reactive skin patch that temporarily attaches the Charger to the skin over the IPG site.

Adverse Event. An undesirable effect.

Amplitude. The measure of strength of delivered stimulation.

Base Station. A power supply holder that charges the Charger.

Battery. The power source for your IPG.

Cardiac Pacemaker. A small implantable device that controls the rhythm of the heart.

Charger. A portable device that charges the battery of implanted rechargeable IPGs.

Charger Belt. An accessory that holds the Charger over rechargeable IPGs for proper charging.

Charger Spacer. A piece of material that may be placed, if directed by your healthcare provider, behind the Charger in the pocket of the Charging Collar or Charger Belt.

Charging Collar. An accessory that holds the Charger over rechargeable IPGs for proper charging.

Charging System. The Charging System charges the battery of rechargeable DBS IPGs.

Computerized Axial Tomography (CT or CAT) Scans. A procedure that creates a 3-D image of your brain or other parts of your body.

Contacts. Metal electrodes on the DBS Lead that deliver electrical stimulation pulses to the brain.

Contraindication. A condition under which the Vercise Gevia DBS System should not be used because the risks outweigh any possible benefit.

Control Buttons. Buttons located on the Remote Control that adjust stimulation settings.

Counterweight. A device placed in the Charging Collar on the opposite side of the Charger to balance the Charging Collar.

DBS Lead. An insulated wire that delivers electrical stimulation pulses from the IPG to the brain.

Deep Brain Stimulation (DBS). A method of applying electrical pulses to the brain to deliver therapy for various disorders.

Diathermy. A therapeutic procedure that heats body tissue with high-frequency electromagnetic currents or ultrasound.

Electrical Stimulation. Electrical pulses created by the DBS IPG.

Electromagnetic Disturbance. Any electromagnetic event that can degrade the performance of a DBS device or system.

Electromagnetic Interference. Reduced performance of a transmission channel or system caused by an electromagnetic disturbance.

Fluoroscopy. An X-ray procedure used during surgery.

Hibernation Mode. A state that the rechargeable IPG reaches when the battery level is too low to apply stimulation.

Idle Mode (Sleep Mode). A time-out period when the Remote Control is not being used.

Implantable Cardioverter Defibrillator (ICD). A small implantable device that treats sudden cardiac arrest and restores a normal heartbeat.

Implantable Pulse Generator (IPG). A device that sends electrical pulses to the brain. The IPG is also referred to as the “Battery” or “Stimulator.”

Incision. A small surgical cut in the skin.

Indicator Light. A signal light on the Charger that shows the status of the Charger.

Level. A term used on the Remote Control screen to identify the amplitude or strength of stimulation.

Magnetic Resonance Imaging (MRI). A technique that uses magnetic fields and radio waves to create pictures of areas inside the body.

Paresthesia. A tingling sensation.

Patient Implant Card. A wallet-size card that lists the patient name, physician names, and the model number and serial number of the implanted medical devices.

Precaution. Generally, situations that you should be aware of to avoid potentially undesirable stimulation effects and/or damage to your Vercise Gevia DBS System.

Program. Settings that define the pattern of your stimulation.

Rechargeable IPG. An IPG with a battery that can be charged. When the battery is depleted, it must be charged for stimulation to continue.

Remote Control. A battery-powered hand-held programmer that adjusts electrical stimulation.

Stimulation. Low level electrical pulses applied to the brain.

Stimulator. A device that sends electrical pulses to the brain. The Stimulator is also referred to as the “Battery” or “Implantable Pulse Generator” (IPG).

Ultrasound. A procedure that uses high frequency sound waves to visualize structures inside the body.

Warning. A potential hazard that you must be aware of to avoid serious situations that may cause injury or death.

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2024 Boston Scientific Corporation