RF 3000™ **Radio Frequency** **Generator**

OPERATION AND **SERVICE MANUAL**

Instructions for Use	3
Mode d'emploi.....	20
Gebrauchsanweisung	39
Istruzioni per l'uso	58
Instrucciones de uso.....	77
Gebruiksaanwijzing	96
Instruções de Utilização	115

ENGLISH

Français

DEUTSCH

ITALIANO

ESPAÑOL

NEDERLANDS

Português

Boston
Scientific

RF 3000™ RADIO FREQUENCY GENERATOR

R_x ONLY **CAUTION!** Federal (USA) law restricts this device to be used by or under the direction of a physician.

INTRODUCTION

The RF 3000 Radio Frequency (RF) Generator is intended for thermal coagulation of soft tissue with electrodes separately cleared by the United States FDA. The generator delivers an isolated RF output of up to 200 watts to the electrode. Full power is available in the impedance range of 25Ω to 100Ω at a constant RF voltage; lower powers are available outside this range. Power can be adjusted manually; however, changes in load impedance will change the actual power delivered (a safety feature described more fully in "Passive Power Limitation"). When power is on, actual delivered RF power and impedance are displayed.

The Operation and Service Manual is intended as an instructional guide only.

GENERAL PROCEDURE SAFETY

WARNING! Safe and effective electrosurgery is dependent not only on equipment design, but also, to a large extent, on factors under the control of the operator. It is important that the instructions supplied, both with this equipment and the accessories, be read, understood, and followed in order to ensure safety and effectiveness.

WARNING! No modification of this equipment is allowed

WARNING! The use of RF energy can produce unintended neuromuscular stimulation. Appropriate precautions, including the use of lower power settings and continuous monitoring of the patient during treatment, should be taken to minimize the risk of patient injury.

CAUTION! The safety of electrosurgery will be greatly enhanced by a thorough knowledge of the medical literature on the subject. Study of specific information on the hazards and complications of the procedure in question is especially recommended.

CAUTION! Whenever practical, the patient should not be allowed to come into contact with electrically conductive parts that are earth grounded or have an appreciable capacitance to earth (for example supports, etc.). The use of antistatic sheeting is recommended for this purpose.

CAUTION! The risk of igniting flammable gases or other materials is inherent in electrosurgery and cannot be eliminated by device design. Precautions must be taken to avoid contact of flammable materials and substances with electrosurgical electrodes. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of HF surgery.

CAUTION FOR MAGNETIC RESONANCE (MR) IMAGING! The RF 3000 Generator is not MR compatible and should be kept outside the shielded MR scan room. In addition, the RF Generator emits signals that may interfere with MR imaging. Therefore, the Generator should be "OFF" when the MR scanner is acquiring image data.

CAUTION! Apparent low output or failure of the RF 3000 Generator to function correctly at the normal operating setting may indicate faulty application of the electrode or poor contact in its connections. In this case, the application of the electrode and its connections should be checked before selecting a higher output power.

CAUTION! Skin-to-skin contact (for example between the arms and body of the patient) should be avoided. Insertion of dry gauze between the skin contact areas is recommended.

EQUIPMENT SETUP

INSTALLATION

Inspect the RF Generator for any signs of physical damage to the front panel, chassis, or cover. If physical damage is found, DO NOT USE THE GENERATOR; contact Boston Scientific. All returned units/parts must be approved by Boston Scientific prior to return.

PREPARATION FOR USE

The RF Generator may be placed on a mounting cart, a sturdy table or platform. Refer to hospital procedures or local codes for detailed information.

Provide at least four to six inches (10-15 cm) of space around the rear, sides and top of the generator for convection cooling. DO NOT obstruct the fan vents on the rear panel. If under continuous use for extended periods of time, it is normal for the top and rear panel to become warm.

CONNECTION TO MAINS POWER

A hospital-grade mains power cord is shipped with the RF Generator. Do not use extension cords or three-prong to two-prong adapters. Connect equipment to input power at 100-120V~ or 220-240V~ at 50-60Hz with medically approved ground. Ensure that the fuse block and voltage selector switch (see "Rear Panel" section) are in the proper positions with respect to mains voltage. Operating the generator with these items in the wrong position could damage the instrument and void the warranty.

WARNING! To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

GENERAL CLEANING INSTRUCTIONS

CAUTION! Do not spray or pour cleaners or any other liquids on the RF generator.

Use a mild detergent and damp cloth to clean the generator cover, front panel, controls and power cable. The RF generator cannot be sterilized. Do not allow fluids to enter the chassis. The generator may be disinfected using an alcohol solution applied with a cloth. Non flammable agents should be used for cleaning and disinfection wherever possible. Flammable agents used for cleaning and disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of RF 3000™ Generator surgery. Also, attention should be called to the danger of ignition of endogenous gases.

ELECTROMAGNETIC INTERFERENCE

WARNING! The use of electrosurgical devices can cause electromagnetic interference (EMI) in other devices, particularly cardiac pacemakers. Precautions should be taken to ensure that the patient's well-being is maintained in the event of such interference.

If EMI appears to be a problem, changing the exact routing of the active or dispersive electrode wires or re-applying the diagnostics sensors (such as ECG leads or oximetry sensors) at a better location (usually farther from the electrodes) can eliminate interference.

TECHNICAL SPECIFICATIONS

OUTPUT POWER:

Adjustable to 200W max. (subject to 2.83A_{RMS} and 141V_{RMS} limits); 2W min.; no output below 15Ω, no output if initial impedance is above 400Ω.

POWER ACCURACY:

Accuracy: ±5% ±2W
Resolution: 1W

OUTPUT FREQUENCY:

461 kHz ± 5 kHz

OUTPUT WAVEFORM:

Quasi Sinusoidal

ELECTRICAL ISOLATION:

4kVRMS @ 60Hz from patient-connected parts to ground.

IMPEDANCE:

Range: 15Ω to 999Ω

Rated Load: 50Ω

Accuracy: ±10% ±5Ω

(displays "LO" below 15Ω and "HI" above 999Ω)

TIME DISPLAY:

Range: 0 to 20 minutes

Accuracy: ±1%

Resolution: 1s

RETURN PAD CURRENT:

Current Limit to any pad: 0.79ARMS ±5% ±0.02ARMS max.

(P-X error, where X is the pad number)

MAINS INPUT VOLTAGE:

100-120V~ or 220-240V~ configured through the power entry module fuse block and a line voltage switch.

MAINS INPUT FREQUENCY:

50/60 Hz

INPUT CURRENT:

4.8 ARMS max @ 100-120V~

2.5 ARMS max @ 220-240V~

FUSES:

100-120V~: T6.3AL 250V, 5 x 20 mm

220-240V~: T3.15AL 250V, 5 x 20 mm

THERMAL SHUTOFF:

Excessive temperature (80 °C) inside the chassis inhibits RF output.

RS-232 PORT:

Female DB9 connector, with the following pinout:

Pin2: Rx/D

Pin3: Tx/D

Pin5: GND (serial communication reference)

Pins 1, 4, 6-9: n/c

RS-232 Communication protocol: 9600 Baud, No Parity, 8 bits, 1 stop bit (9600, N, 8, 1)

SIZE:

35.5 cm W x 41.9 cm L x 10.6 cm H (14" W x 16.5" L x 4.2" H) (including protrusions)

WEIGHT:

6.350 kg (14 lb)

OPERATING CONDITIONS:

Temperature: +10 °C to +40 °C

Humidity: 30% to 75% (non-condensing)

Atmospheric Pressure: ≥ 800 millibars (up to 2,000 meters)

CAUTION! For altitudes above 1,500 meters, the maximum operating temperature is 35 °C.

STORAGE/TRANSPORT CONDITIONS:

Temperature: 0 °C to +70 °C

Humidity: 0% to 95% (non-condensing)

Atmospheric Pressure: ≥ 800 millibars (up to 2,000 meters)


Figure A: Typical Maximum Output Power vs. Load Impedance


Figure 1B: Actual Output Power vs. Setpoint into Nominal, 50W Load


Figure 1C: Maximum Peak Output Voltage vs. Power Setpoint

OPERATION


Figure 2: RF 3000™ Radio Frequency Generator Front Panel

DISPLAYS

Item	Function
Time (min:sec)	<i>Time Display</i> – In Ready mode, displays the time duration for which RF energy will be applied, as entered by the operator. In Operating mode, displays the total elapsed time. Default time is 15:00.
Power (W)	<i>Power Display</i> – During Operating mode, displays the power output. The unit operates as a constant voltage source during Operating mode, so output power will vary as load impedance changes. During Ready mode, this display is blank.
Impedance (Ω)	<i>Impedance Display</i> – In Operating mode, this display shows both the currently measured impedance and a 10-element display of the value relative to the initially measured impedance. During Ready mode, this display is blank.
 Active	<i>Active Indicator</i> – Illuminates during Operating mode indicating active RF output (2W or more of RF power). When this indicator is on, measured power and elapsed time are also displayed.

CONTROLS

Symbol	Function
	<i>Increment</i> – Increases the value displayed above. A single depression of the key increases the value by one unit. Continuous depression increases to the maximum value, pausing at specified intervals.
	<i>Decrement</i> – Similar function to the Increment button, except that decreases in value are produced.
	<i>Reset</i> – Transfers the system from Standby mode or Fault mode into Ready mode.
	<i>Start</i> – Transfers the system from Ready mode into Operating mode.

Stop – Transfers the system from Operating mode into Standby mode.

RECEPTACLES

Symbol

Function

Electrode Receptacle – Accepts the plug on the electrode cord for energy delivery through the LeVeve™ electrode or Soloist™ Electrodes.

Return Receptacles – Accept the plugs on the Megadyne® 0850C Patient Return Electrodes. Also, uniquely identifies the four return pads for troubleshooting.

Safety Symbol – Indicates that the receptacles have CF (Cardiac, Floating) type electrical isolation from mains and ground and that the RF 3000™ Generator is protected from damage during a defibrillator discharge directly into the receptacles.

Electrode isolated from earth ground at high frequency.


REAR PANEL

Item


Function

Earth Ground Stud – Provides a means of securely linking the earth grounds of the RF generator to other grounded equipment.

RS-232 Port – For authorized service personnel only.

Power Entry Module – This module contains both the power ON/OFF switch and the fuses. The voltage is selected by the orientation of the fuse drawer, as marked.


Voltage Selector – This switch configures the generator to either 110V~ or 220V~ in conjunction with the Power Entry Module. The appropriate tool must be used to change configurations.

Volume Control Knob – This knob controls the tone volume from 40dBA @ 1 meter except pad-off tones which are fixed at 65dBA @ 1 meter.

LABELING

Symbol	Meaning
	Caution: Dangerous Voltage
	Caution / Warning
	Earth Ground
	Fuses
	Power ON
	Power OFF
	Radiofrequency energy source
	Apply 2 return pads to each leg (four pads total)
	Alternating Current

INITIAL POWER-ON MODE

CAUTION! Ensure that the Voltage Selector and the fuse block on the rear panel are set to the appropriate mains voltage before turning the mains power switch on.

The RF 3000™ Generator is switched on by setting the mains power switch, located on the rear of the unit on the power entry module, to the “I” position.

After mains power is applied, the RF 3000 Generator will briefly illuminate all displays and sound a short tone. A self test is performed with a brief application of RF power to an internal load. If any problems are identified, an error code will be displayed in the power window and the unit goes to Fault mode. If the self test is successful, the unit will proceed to the Ready mode.

READY MODE

In the Ready mode, the RF 3000 Generator is ready for ablation duration time to be adjusted. The TIME window will display the default time of 15:00 minutes, or a selected time in the acceptable range of values: 00:30 to 20:00 minutes. RF power is not applied in this mode of operation, therefore the ACTIVE light and the “RF On” tone are off. POWER and IMPEDANCE windows will be blank.

When START is pressed, the Operating mode is initiated.

OPERATING MODE

After START is pressed, the Operating mode is initiated. Before any significant RF power is applied between the output receptacles, the initial system impedance [Z(0)] is measured. If the initial impedance is between 15Ω and 400Ω, TIME starts incrementing from "00:00" and the POWER window will display the value of power selected. The maximum settable value for power is 200W. When power is being delivered, the ACTIVE light and "RF On" tone are on. The IMPEDANCE window gives an indication of both the currently measured system impedance and the impedance relative to Z(0). A greater number of elements illuminated indicates a higher value of impedance.

STANDBY MODE

In this mode, the final value of power and impedance, and the elapsed ablation time will be displayed. Settings cannot be changed in this mode and RF power, the ACTIVE light, and the "RF On" tone are all off.

Pressing the RESET button will clear the IMPEDANCE and POWER windows and initiate the Ready mode. Alternatively, pressing the START button will return the unit directly from the Standby mode to the Operating mode. In this case, ablation time will be the previously selected value, or 15:00 minutes if the default is used.

FAULT MODE

The Fault mode can be entered from any of the above modes. A series of rapid beeps indicates a fault, with the specific error code (tabulated below) displayed in the POWER window. "Recoverable" faults require only that the RESET button be pressed to clear the error code and return to Ready mode. For the "Unrecoverable" faults, sometimes cycling the mains power switch on the rear of the unit will bring about recovery. When this does not work, a more serious malfunction may be the cause. Boston Scientific should be contacted for assistance in further diagnosing the difficulty.

SHUT DOWN

The RF3000™ Generator is switched off by setting the mains power switch, located on the rear of the unit on the power entry module, to the "O" position

ERROR TABLE

Hardware Self-Test Errors:

Hxx: Unrecoverable errors due to problems with computer circuits or analog circuits

Measurement Errors:

Error Code	Description	Unrecoverable Error	Recoverable Error
E01	Low Impedance Measured (Displays "LO")		✓
E02	High Initial Impedance Measured		✓
E03	High Impedance Measured (Displays "HI")		✓
P-X	High current in Pad X		✓
E80	Pad Current Sum Error	✓	
E81	Hardware Shutdown:	✓	
	Over Voltage		
E82	Hardware Shutdown:	✓	
	Over Current		
E83	Hardware Shutdown:	✓	
	Over Power		
E84	Hardware Shutdown:	✓	
	Multiple		
E85	Reference Voltage Error	✓	
E86	Power Offset Error	✓	
E87	High Voltage Measured	✓	

E88	High Current Measured	☒
E89	High Power Measured	☒
E91	Software Shutdown	☒
E92	PAD Current Offset Error	☒
E93	V(rms) Offset Error	☒
E94	I(rms) Offset Error	☒

ABLATION

CAUTION! The output power selected should be as low as possible for the intended purpose.

CAUTION! Use of this device with laparoscopic instruments may result in a gas embolism. When this generator is used with a laparoscopic electrode, activation of the electrode, when it is not in contact with or in position to deliver energy to the target tissue, may cause capacitive coupling and result in patient burns.

CAUTION! Failure of the RF 3000™ Generator could result in an unintended increase of output power.

CAUTION! The RF3000 Generator is energized without continuous activation of a switch. Continuous monitoring of the procedure is recommended.

CAUTION! The cables of the surgical electrode should be positioned in such a way that contact with the patient or other leads is avoided. Temporarily unused active electrodes should be stored in a location that is isolated from the patient.

CAUTION! All Active accessories are required to have a minimum rated voltage of 200 Volts peak.

CAUTION! Electrodes and probes of monitoring, stimulating, and imaging devices can provide paths for high-frequency currents even if these devices are battery powered, insulated, or isolated at 60Hz (or 50Hz). The risk of burns can be reduced, but not eliminated, by placing the respective electrodes or probes as far away as possible from the electrosurgical site and from the dispersive electrode.

CAUTION! Apparent low output or failure of the RF 3000 Generator to function correctly at the normal operating settings may indicate faulty application of the electrode or poor contact in its connections. In this case, the application of the electrode and its connections should be checked before selecting a higher output power.

CAUTION! The user should regularly inspect the accessories and electrode cables for possible damage to the insulation. Discard if damage is detected.

CAUTION! When RF 3000 Generator and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended.

CONNECTION OF RETURN PADS

WARNING! The use and proper placement of dispersive electrodes (return pads) is a key element in the safe and effective use of monopolar electrosurgery, particularly in the prevention of burns. Follow the manufacturer's Instructions for Use.

Boston Scientific recommends the use of Megadyne® 0850C Disposable Electrosurgical Patient Return Electrodes (Pads). Use four pads and connect the return pad plugs to the RF 3000 Generator return receptacles, placing two pads on each leg, as shown in Figure 3. Be sure not to orient the pads improperly or to misalign the top of edge of the pads, as shown in Figure 4.

Figure 3: Proper Orientation and Alignment of Return Pads on Upper Mid-Thigh


Figure 4: Improper Orientation and Misalignment of Upper Edge

Surveillance of the dispersive electrodes is essential for safety, especially for lengthy ablations. This should occur at the beginning of the ablation, to prevent burns by early detection of a misapplied electrode. Surveillance is also important during ablation to detect a pad for which the interface quality is compromised, such as in peeling or tenting.

To aid detection of poorly adhered or improperly placed return pads, the RF 3000™ Generator has been equipped with the Pad-Guard™ Current Monitoring feature. Each of the four return pads is monitored for current flow. Current exceeding 0.8 amp results in an error message indicating the high current pad (P1, P2, P3, or P4) and halts the generator.

Causes for a P(x) error include improper location of pad placement (one pad placed closer to the active electrode, thus absorbing more current than the other pads) or poor pad contact with the skin (peeling or tenting).

In case of a P(x) error, all pads should be checked for location (see Figures 3 & 4) and proper skin contact.

Pad-Guard reduces the likelihood of a skin burn, but does not eliminate it completely. One should monitor both the pad contact and skin conditions periodically throughout the procedure.

CONNECTION OF ACTIVE ELECTRODE

WARNING! When not in use, the active electrodes should never touch the patient.

Boston Scientific recommends the use of the Boston Scientific LeVeen™ Electrode or Soloist™ Electrode with the RF 3000 Generator. Deploy active electrodes according to manufacturer's instructions. Ensure that the active electrode plug is connected securely to the RF 3000 Generator before proceeding.

PREPARING FOR ABLATION

- Follow all CAUTIONS, WARNINGS, and instructions in the following sections: Introduction, Operation, and Ablation.
- Ensure that the RF Generator is turned on by setting the mains power switch to "I". This switch is located on the rear panel on the power entry module (where the power cord attaches).

- Adjust the timer to the ablation time specified by the physician using the timer buttons.
- Ensure that the active and dispersive electrodes are properly deployed and attached both to the patient and the RF generator (see “Connection of Return Pads” and “Connection of Active Electrode”).

PERFORMING ABLATION

Ablation is initiated by pressing the START button, which causes the TIME window to display “00:00” and start incrementing. The POWER and IMPEDANCE windows are blank until the power is set to at least 2W. After pressing the START button, power should immediately be adjusted by the physician / user to the lowest value that delivers the power required for tissue effect within the desired ablation time.

If error code E01, E02, or E03 is displayed, check the connections between the generator and electrode and return pads.

If error code P-1, P-2, P-3, or P-4 is displayed, excessive RF current is flowing through the return pad indicated. Pads are identified next to the receptacle into which they are connected. Ensure that the pads are placed according to Figure 3 and not according to Figure 4.

CAUTION! The power applied to any electrode should be kept to the minimum necessary to achieve the desired surgical effect.

As RF ablation heats the tissue, impedance increases, resulting in more elements of the impedance display being lit. The “RF On” tone frequency increases as more elements are lit. When all ten elements are lit, impedance has roughly doubled. After this, impedance generally increases to a much higher value, and power decreases to less than 15W. When this happens the RF power should be turned off.

PASSIVE POWER LIMITATION

The RF 3000™ Generator is an adjustable, constant-voltage RF power supply. Delivering RF energy into target tissues at the desired rate produces heat, coagulation, and corresponding desiccation. This causes tissue impedance to rise. At a constant voltage, power is inversely proportional to the load or tissue impedance. For example, if impedance doubles, power decreases by 50%. RF power is therefore inherently limited by the response of the tissue which tends to protect the tissue from damage due to over-ablation.

TEST SPECIFICATIONS AND SERVICE

POWER-ON SELF-TEST

Every time the RF 3000 Generator mains power switch is turned on, a Power-On Self-Test is performed by the internal microprocessor that checks numerous systems and functions, including:

- the function of all displays and the “RF On” tone;
- the accuracy of power, voltage, and current measurements;
- the Reference Voltage accuracy;
- the Return Pad current-measurement accuracy;
- the A/D calibration accuracy;
- RAM functionality;
- firmware cyclical redundancy check on the software.

If any of these tests fail, then an error code is displayed and the radiofrequency generator (RFG) will not deliver energy. Sometimes, cycling mains power will clear the error, but if this does not work, then the RFG must be returned to Boston Scientific.

POWER ACCURACY TEST

CAUTION! To protect the test equipment, ensure that the electrosurgical analyzer can dissipate 200W at 100% duty cycle.

The only active socket in the electrode receptacle is the center one; the three peripheral sockets are for future use. Also, both pins of each of the four return receptacles (8 pins total) are electrically common with all the other pins.

For each return pad receptacle, the current is limited to 0.79 A_{RMS}. Therefore, measuring the maximum power output of the generator will require separate loads (100Ω to 400Ω each) connected in parallel to each return receptacle. If just one load is connected to all four return receptacles, then even a small difference in lead resistance between each return receptacle and the load can cause greatly different currents to flow which can create an over-current condition (P-x error) before maximum output is attained. One method for testing power accuracy at 200W is to use four non-inductive 100Ω loads connected between the active receptacle and each return receptacle. (The generator sees 25Ω in this configuration.)

Alternatively, a lower power level into a single 100Ω load can be used to test power accuracy. Refer to the Technical Specifications section for power accuracy specifications.

ERROR TESTING

The "Recoverable" errors (E01, 02, 03, P-x) may be tested, as outlined below. The "Unrecoverable" errors cannot be tested, since the conditions that lead to these errors are difficult to replicate.

"E 01" indicates that the initial impedance is below 15Ω. If a 10Ω load is connected between the active receptacle and one of the return receptacles, this error will be displayed in the Power window after the START button is pressed. The Impedance display will flash "LO" and the bar elements will flash until the RESET button is pressed.

"E 02" indicates that the initial impedance is above 400Ω. This error can be tested by pressing the START button with an open circuit between the active and return receptacles. The Impedance display will flash "HI" and the bar elements will flash until the RESET button is pressed.

"E 03" indicates that impedance is over 4000Ω. To test this error, connect a 50Ω load, increase power to 10W, then disconnect the load. The Impedance display will flash "HI" and the bar elements will flash until the RESET button is pressed.

"P-x" indicates that excessive current is flowing through return receptacle "x". To test this error, connect a short circuit from a return receptacle to the active receptacle and press the START button. The error will immediately be displayed since an RF current flows before the power increment button is pressed to measure the initial impedance. Repeat this test for each return receptacle.

LEAKAGE CURRENT TESTING

Every generator is tested for leakage currents at 132V~ for the 115V~ model or 264V~ for the 230V~ model.

CASE LEAKAGE

Case (or "enclosure") leakage current is measured between the metal enclosure and mains ground, with the safety ground open, at both normal and reverse polarity of mains voltage. Case leakage must be below 100μA.

SOURCE LEAKAGE, APPLIED PART

Source leakage current (or "patient leakage source current") is measured between mains ground and all output receptacles, with the safety ground open, at both normal and reverse polarity of mains voltage. Source leakage must be below 10μA.

SINK CURRENT, APPLIED PART

Sink current (or “isolated leakage sink current”) is measured between mains input and all output receptacles, at both normal and reverse polarity of mains voltage. Sink current must be below 50 μ A.

FIELD SERVICE PROCEDURES

The only field-replaceable components of the RFG are the fuses, which are identified in “Fuses”. Be sure to use the correct fuses and replace the fuse drawer in the proper position for the mains voltage of the RFG, as marked on the fuse drawer.

For all other service requirements, please contact Boston Scientific.

DISPOSAL

The user should follow local and national regulations for disposal of electronics when disposing of this unit.

RF 3000™ Radio Frequency Generator Safety Specifications

DEVICE DESCRIPTION
Class I, Defibrillation proof Type CF Equipment, IPX0, not AP/APG
Mode of Operation: Continuous
EMC Emissions and Susceptibility: The RF 3000 Radio Frequency Generator has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This system generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instruction given below, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

ELECTRICAL ISOLATION
Leakage current conforms to IEC 60601-1
Dielectric withstand conforms to IEC 60601-1

GUIDANCE AND MANUFACTURER’S DECLARATION – ELECTROMAGNETIC EMISSIONS		
The RF 3000 Radio Frequency Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the RF 3000 Radio Frequency Generator should assure that it is used in such an environment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
RF emissions CISPR 11	Group 2	The RF 3000 Radio Frequency Generator must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected. The RF 3000 Radio Frequency Generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The RF 3000™ Radio Frequency Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the RF 3000 Radio Frequency Generator should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	±2 kV for power supply lines ±1 kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_t (>95% dip in U_t) for 0.5 cycles 40% U_t (60% dip in U_t) for 5 cycles 70% U_t (30% dip in U_t) for 25 cycles <5% U_t (>95% dip in U_t) for 5 sec	<5% U_t (>95% dip in U_t) for 0.5 cycles 40% U_t (60% dip in U_t) for 5 cycles 70% U_t (30% dip in U_t) for 25 cycles <5% U_t (>95% dip in U_t) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the RF 3000 Radio Frequency Generator requires continued operation during power mains interruptions, it is recommended that the RF 3000 Radio Frequency Generator be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of typical location in a typical commercial or hospital environment.
NOTE - The RF 3000 does not possess a DC power port. The RF 3000 interface cables do not exceed 3 meters in length.			

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The RF 3000™ Radio Frequency Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the RF 3000 Radio Frequency Generator should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the RF 3000 Radio Frequency Generator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.33\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RF 3000 Radio Frequency Generator or any of its components are used exceeds the applicable RF compliance level above, the RF 3000 Radio Frequency Generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire RF 3000 Radio Frequency Generator.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The RF 3000™ Radio Frequency Generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RF 3000 Radio Frequency Generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RF 3000 Radio Frequency Generator as recommended below, according to the maximum output power of the communications equipment.			
RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER M		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.18	0.18	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.7	11.7	23.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

WARRANTY

Boston Scientific Corporation (BSC) warrants for one year from the date of purchase that reasonable care has been used in the design and manufacture of this product. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, cleaning and storage of the product as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Boston Scientific Corporation's control may directly affect the product and results obtained from it. Boston Scientific Corporation shall repair or replace, at its option, any part of the product that Boston Scientific Corporation determines was defective at time of shipment if notice thereof is received within one year of shipment. Boston Scientific Corporation shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of the product. Boston Scientific Corporation neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with the product.

Buyer shall be responsible for the ongoing support and maintenance of the product not covered by this one-year warranty and after the one year warranty period has expired. Buyer may, at its sole cost and expense, purchase an extended warranty from Boston Scientific Corporation (BSC) to extend the term of this warranty.

Megadyne is a trademark of Megadyne Medical Products, Incorporated.

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Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited