AngioJet™ Ultra Power Pulse™ Kit

Y-sets

Instructions for Use	2
Instrucciones de uso	5
Mode d'emploi	9
Gebrauchsanweisung	13
Istruzioni per l'uso	17
Instructies voor gebruik	21
Instruções de Utilização	25

TABLE OF CONTENTS

REUSE WARNING 3

DEVICE DESCRIPTION 3

 Figure 1. Power Pulse Kit (shown with protective caps removed)..... 3

 Contents 3

 Non-Pyrogenic..... 3

INTENDED USE/INDICATIONS FOR USE 3

CLINICAL BENEFIT STATEMENT..... 3

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE..... 3

CONTRAINDICATIONS..... 3

WARNINGS..... 3

PRECAUTIONS..... 3

POTENTIAL ADVERSE EVENTS 3

HOW SUPPLIED 3

 Device Details 4

 Handling and Storage..... 4

OPERATIONAL INSTRUCTIONS 4

 Required Items (not included in the package) 4

 Power Pulse Kit Preparation 4

 Enter Power Pulse Mode 4

 Power Pulse Spray Infusion Procedure..... 4

 Exit Power Pulse Mode 4

 Thrombectomy Set Operational Procedure 4

 Disposal 4

 Post-Procedure 4

INFORMATION TO BRIEF THE PATIENT..... 4

WARRANTY..... 4

SYMBOL DEFINITIONS..... 4

AngioJet™ Ultra Power Pulse™ Kit Y-sets

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

DEVICE DESCRIPTION

The AngioJet Ultra Power Pulse Kit Y-sets (Power Pulse Kit or Y-set) (see Figure 1) enable the AngioJet Thrombectomy Set to deliver a pulsed infusion of physician-specified fluid to a local treatment area. The Power Pulse Kit is a plastic Y-set with vented bag spikes that are bonded to tubing. Each upper arm of the Y-set contains a tubing clamp, one white and one red.

The vented bag spikes are used to access the standard intravenous saline solution bag used with the AngioJet System and a second intravenous bag containing physician-specified fluid. The tubing clamps are used to control the flow of the two fluids. The physician-specified fluid is delivered using the AngioJet System, which includes the AngioJet Console and Thrombectomy Set.

Figure 1. Power Pulse Kit (shown with protective caps removed)

Contents

The Power Pulse Kit is supplied as follows:

Quantity	Material
Five (5) each	Y-set

Non-Pyrogenic

This device meets pyrogen limit specifications.

INTENDED USE/INDICATIONS FOR USE

The Power Pulse Kit is intended for use only with AngioJet Thrombectomy Sets indicated for the control and selective infusion of physician-specified fluids, including thrombolytic agents, into the peripheral vascular system using the AngioJet System.

CLINICAL BENEFIT STATEMENT

The clinical benefit of the AngioJet Thrombectomy System is to break apart and remove thrombus from target vessels with or without infusion of physician-specified fluids resulting in improved patency with low rates of major complications.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website: <https://ec.europa.eu/tools/eudamed>

CONTRAINDICATIONS

Refer to the individual Thrombectomy Set Instructions for Use for specific contraindications.

WARNINGS

- Inspect sealed sterile package before opening. If packaging integrity has been compromised, do not use, reprocess or resterilize the contents, and contact a Boston Scientific representative.
- Refer to the individual AngioJet Thrombectomy Set Instructions for Use for specific warnings and precautions.
- Refer to the product insert supplied with the physician-specified fluid used during Power Pulse technique for contraindications, side effects, warnings and precautions.
- The Power Pulse Kit has NOT been tested for use in administering contrast medium. Do not administer contrast medium using this Y-set. Use of contrast medium may result in patient injury.
- Verify the Power Pulse Kit tubing clamps are in the required open or closed position prior to running in Thrombectomy Mode. Improper fluid delivery may result in patient injury.

PRECAUTIONS

- The Power Pulse Kit assembly should only be used with those AngioJet Thrombectomy Sets which are indicated for use with Power Pulse technique.
- Verify Y-set "Use By" date printed on the product label has not expired. Do not use product if "Use By" date has expired, as sterility or function may be compromised.
- Visually inspect the Power Pulse Kit and other AngioJet components prior to use to ensure that no damage has occurred during shipment. To avoid potential injury do not use a damaged device. If any damage is noted, contact a Boston Scientific representative.

POTENTIAL ADVERSE EVENTS

Potential adverse events, which may be associated with use of the AngioJet Thrombectomy Sets, including Power Pulse technique in peripheral vessels, are similar to those associated with other interventional procedures and include but are not limited to:

- Allergic reactions (to drug, contrast, device or other)
- Aneurysm
- Arrhythmia
- Bleeding/hemorrhage
- Cerebrovascular accident (CVA), transient ischemic attack (TIA)
- Death
- Electric shock
- Embolism (air, plaque, thrombus, device or other)
- Hematoma
- Hemolysis
- Hypotension/hypertension
- Infection/sepsis
- Myocardial ischemia/infarction
- Need for additional intervention or surgery
- Pain
- Pancreatitis
- Renal insufficiency/failure
- Thrombus/thrombosis
- Vasospasm
- Vessel trauma (dissection, injury, perforation, pseudoaneurysm, rupture)
- Vessel occlusion

HOW SUPPLIED

The Power Pulse Kit is supplied sterile inside a pouch. The device is sterilized using an ethylene oxide (EO) process.

Device Details

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

The AngioJet System should only be used by physicians who have a thorough understanding of thrombectomy treatments and percutaneous interventional procedures.

A thorough understanding of each component in the AngioJet System is required for proper operation. Read the Instructions for Use, the Instructions for Use packaged with the AngioJet Thrombectomy Set, as well as the Operator's Manual supplied with the AngioJet Console before attempting to use any part of the System.

Required Items (not included in the package)

- AngioJet Console
- AngioJet Thrombectomy Set
- Heparinized normal saline (HepNS) for injection
- Physician-specified fluid

Power Pulse Kit Preparation

Note: For easy visual identification, use the short line/white clamp for heparinized saline, and use the long line/red clamp for physician-specified fluid.

1. Prepare heparinized saline per standard procedure.
2. Prepare physician-specified fluid according to the medication package insert.

Note: Bag hooks are provided on the Console for hanging fluid sources.

3. Remove Y-set assembly from sterile package.
4. Select the saline line spike.
5. Close the tubing clamp, remove the spike cover and insert the spike into the saline bag.
6. Open the clamp and allow saline to fill the tubing; close clamp.
7. Select the physician-specified fluid line spike.
8. Close the tubing clamp, remove the spike cover and insert the spike into the physician-specified fluid bag.
9. Open the clamp and allow physician-specified fluid to fill the remainder of the tubing; close clamp.
10. Remove the vented cap from the end of the Y-set and connect Thrombectomy Set to the Y-set.
11. Open the clamp on the saline line to allow the fluid to flow.
12. Prepare and prime per the AngioJet Thrombectomy Set Instructions for Use.

Enter Power Pulse Mode

The following set-up instructions are specific to the Thrombectomy Set use with the AngioJet Ultra Console. If you are not using the AngioJet Ultra Console, do not use the following set-up instructions and instead refer to the Operator's Manual for the AngioJet Console you are using.

1. Press the CATHETER button on the Console.
2. Press the CATHETER button a second time, within 3 seconds, to advance to the next screen.
3. Press the arrow button to select YES on the decision (YES/NO) screen to confirm access to Power Pulse Mode.
4. Press the CATHETER button again to enter Power Pulse Mode.

Note: The Power Pulse volume and time are displayed on the main screen and "PP" is shown in the Console time display. The title on the screen should read INFUSION IN PROGRESS.

Power Pulse Spray Infusion Procedure

Note: The following infusion procedure is a guideline. The exact procedure and infused fluid are determined by the physician.

1. Open the clamp on the physician-specified fluid line. Close the clamp on the saline line.
2. Advance the catheter to the treatment site.
3. Press the foot switch to infuse the desired volume of physician-specified fluid.

Note: Each pump stroke delivers approximately 0.6 cm³ (0.6 cc) of fluid.

4. The catheter should be moved through the treatment site at the physician's discretion.
5. When the desired amount of physician-specified fluid has been delivered, release the foot switch.
6. Close the clamp on the physician-specified fluid line.

7. It is recommended to allow a dwell time after infusion, typically 30 minutes, before starting thrombectomy. Dwell times may vary.

Exit Power Pulse Mode

Method 1:

1. Select the ALARM RESET button to exit Power Pulse Mode immediately and return to Thrombectomy Mode.
2. Open the clamp to the heparinized saline.

Method 2:

1. Access the decision (YES/NO) screen by pressing the CATHETER button twice.
2. Select NO and press the CATHETER button again to advance to the READY screen.
3. Open the clamp to the heparinized saline.

Note: After exiting Power Pulse Mode, the total Power Pulse time and volume will be displayed on the screen. If Power Pulse Mode is re-entered, these values will not reset and will continue to increase.

Thrombectomy Set Operational Procedure

Note: When changing fluids, purge the lines for 45 seconds to ensure lines are cleared.

1. OPEN clamp to heparinized saline.
2. CLOSE clamp to physician-specified fluid.
3. Proceed with thrombectomy per AngioJet Thrombectomy Set Instructions for Use.

Disposal

Follow proper precautions for the handling of infectious waste. Reuse of Power Pulse Kit is prohibited due to risk of contamination by blood waste products. After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

- Assess patient for hematoma and/or other signs of bleeding at the puncture site.
- Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling the patient on the use of the AngioJet System in association with the interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document, both for AngioJet System and for other interventional treatments likely to be employed.
- Inform the patient that cardiac arrhythmias during catheter operation have been reported in a small number of patients and may require additional measures during the procedure such as temporary pacing.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and homecare, or rehabilitation guidelines.

No separate patient information is provided because the Power Pulse Kit is used in association with other interventional devices. Risks and benefits to the patient are similar and, from the patient's perspective, a part of those interventional procedures.

WARRANTY

For device warranty information, visit www.bostonscientific.com/warranty.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

The following are trademarks of Boston Scientific Corporation or its affiliates: AngioJet, Power Pulse. All other trademarks are the property of their respective owners.

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2023 Boston Scientific Corporation or its affiliates.