IceRod™

1.5 MRI STRAIGHT
1.5 MRI 90°

IceSeed™

1.5 MRI STRAIGHT
1.5 MRI 90°

Cryoablation Needle

en	Instructions for Use	2
es	Instrucciones de uso	9
fr	Mode d'emploi	16
de	Gebrauchsanweisung	23
it	Istruzioni per l'uso	30
nl	Instructies voor gebruik	38
pt-EU	Instruções de Utilização	45

TABLE OF CONTENTS

REUSE WARNING 3

DEVICE DESCRIPTION 3

 Contents 3

 Scan Card 3

 Figure 1. 1.5 MRI Cryoablation Needle – Components 3

 Figure 2. 1.5 MRI Cryoablation Needle – Needle Shaft Marks 3

 Table 1. 1.5 MRI Cryoablation Needle – Handle Colors 4

 Materials 4

 Operating Principle 4

 Needle Performance - Cryoablation 4

 Figure 3. 1.5 MRI Cryoablation Needle Isotherm Data 4

 Needle Performance - Active Thaw 4

 User Information 4

INTENDED USE 4

 Patient Groups 4

INDICATIONS FOR USE 5

 Clinical Benefit Statement 5

CONTRAINDICATIONS 5

WARNINGS 5

PRECAUTIONS 6

 General 6

 Handling 6

 Procedural 6

ADVERSE EVENTS 6

CONFORMANCE TO STANDARDS 7

HOW SUPPLIED 7

 Device Details 7

 Handling and Storage 7

MRI SAFETY INFORMATION 7

OPERATIONAL INSTRUCTIONS 7

 Needle Preparation 7

 Needle Handling and Insertion 7

 Notes for Conducting Freezing 8

 Notes for Conducting Thawing 8

 Needle Removal 8

 Disposal 8

 Post-Procedure 8

INFORMATION TO BRIEF THE PATIENT 8

WARRANTY 8

SYMBOL DEFINITIONS 8

IceRod™

1.5 MRI STRAIGHT
1.5 MRI 90°

IceSeed™

1.5 MRI STRAIGHT
1.5 MRI 90°

Cryoablation Needle

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

The Boston Scientific IceRod Cryoablation Needle and IceSeed Cryoablation Needle, compatible with Magnetic Resonance Imaging (MRI Needles), are disposable, single-use products. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

DEVICE DESCRIPTION

Contents


One (1) of the following MRI Needles:

- IceRod Cryoablation Needle, 1.5 MRI Straight
- IceRod Cryoablation Needle, 1.5 MRI 90°
- IceSeed Cryoablation Needle, 1.5 MRI Straight
- IceSeed Cryoablation Needle, 1.5 MRI 90°

The Boston Scientific MRI Needles (straight or angled 90°) have a sharp cutting tip, a shaft with circumferential marks, a color-coded handle, gas tubing, and a gas connector. All components are illustrated in Figure 1. Needle shaft marks are shown in Figure 2. Needle handle colors are listed in Table 1.


Scan Card

For use with the “Other Procedure Barcode Card” refer to the SeedNet Gold Cryoablation User Manual.

#	Component
1	Cutting Tip
2	Needle Shaft (175 mm; diameter 1.5 mm)
3	Needle Handle (color coded)
4	Gas Tubing(length ~2.5 m)
5	Gas Connector

Figure 1. 1.5 MRI Cryoablation Needle – Components

#	Shaft Marks	Description
1	Thick Marks	Start 20 mm from tip, in 10 mm intervals
2	Thin Marks	Start 25 mm from tip, in 5 mm intervals
3	Single-thick Mark	Located 50 mm from tip
4	Double-thick Mark	Located 100 mm from tip
5	Triple-thick Mark	Located 150 mm from tip

Figure 2. 1.5 MRI Cryoablation Needle – Needle Shaft Marks

Table 1. 1.5 MRI Cryoablation Needle – Handle Colors

Needle Brand Name	Configuration	Handle Color
IceSeed Cryoablation Needle 1.5 MRI	Straight	Black
IceSeed Cryoablation Needle 1.5 MRI	Angled 90°	Black
IceRod Cryoablation Needle 1.5 MRI	Straight	Red
IceRod Cryoablation Needle 1.5 MRI	Angled 90°	Red

Materials

Materials	MRI Needle Shaft	Nickel-Chromium superalloy
	Needle Handle	Brass (coated with heat shrink tubing)
	Gas Tubing	Polyurethane
	Gas Connector	Polyoxymethylene
Sterilization Method	Ethylene Oxide (EO)	

Operating Principle

The Boston Scientific Visual-ICE MRI Cryoablation System is a portable system intended for cryoablative tissue destruction using a minimally invasive procedure. The system is computer-controlled with a touch screen user interface that allows the user to control and monitor the procedure. The therapy delivered by the system is based on the Joule-Thomson effect displayed by compressed gases. The Joule-Thomson effect is a change in the temperature of a compressed gas as it flows through a narrow orifice and expands to a lower pressure. Certain gases, such as argon, decrease in temperature due to the Joule-Thomson effect, while other gases, such as helium, increase in temperature. The cryoablation system uses high-pressure argon gas that circulates through closed-tip cryoablation needles to induce tissue freezing. Freezing is due to the Joule-Thomson effect augmented by a recuperative heat exchanger within the cryoablation needle.

Details on how tissue thawing is achieved can be found in the **Needle Performance - Active Thaw** section.

Cryoablative treatments are performed using a small-diameter cryoablation needle(s) inserted into the target lesion, usually under image guidance. When cryoablation needle(s) are placed into or near the target tissue and freezing is initiated, an iceball grows around the distal end of the needle shafts. In time, the iceball completely engulfs and destroys the target tissue. The shape of the iceball (ellipsoidal or spherical) depends on the type of needle used.

Needle Performance - Cryoablation

In vivo iceball dimensions and resulting ablation zone are determined by the selected cryoablation needle, number of needles placed, tissue and tumor characteristics, thermal heat sink from surrounding vasculature, and treatment duration. Monitoring iceball formation provides direct control throughout the procedure and is key to cryoablation success.

NOTE: Using MRI visualization, monitor iceball formation throughout the cryoablation procedure.

Iceball shapes and dimensions are specific to the needle brand. The following laboratory model of iceball dimensions is provided to assist users in selecting the cryoablation needle(s) and needle placement to appropriately ablate the target area. Typically, in vivo dimensions are smaller than the dimensions generated in the laboratory. Laboratory testing was performed in room temperature (21 °C) gel; measurements were made after two (10 minute) Freeze cycles separated by a 5 minute passive Thaw cycle. Accuracy is ±3 mm width, ±4 mm length.


Figure 3. 1.5 MRI Cryoablation Needle Isotherm Data

Needle Performance - Active Thaw

The MRI Needles support active thawing using helium.

NOTE: In vivo thawing performance will be determined by numerous factors, including the brand (type) and number of cryoablation needles used, needle locations, tissue characteristics, and duration of thaw activation.

User Information

The intended users are medical professionals who have a thorough understanding of the technical principles, clinical applications, and risk associated with cryoablation procedures. Refer to the Visual-ICE MRI Cryoablation System User Manual for optional education.

INTENDED USE

The Boston Scientific disposable IceRod Cryoablation Needle and IceSeed Cryoablation Needle (MRI Needles) are meant to be connected to a Boston Scientific Visual-ICE MRI Cryoablation System when performing cryoablative tissue destruction through application of extremely cold temperatures. The needles are intended to convert high-pressure gas to either a very cold Freezing application or to a warm Thawing application.

Patient Groups

The intended population includes patients intended for cryoablative destruction of tissue during surgical procedures.

INDICATIONS FOR USE

The Boston Scientific disposable IceRod Cryoablation Needle and IceSeed Cryoablation Needle (MRI Needles), when used with a Boston Scientific Visual-ICE MRI Cryoablation System, are designed to destroy tissue by the application of extremely cold temperatures. A full list of specific indications can be found in the Boston Scientific Visual-ICE MRI Cryoablation System User Manual.

Clinical Benefit Statement

The Visual-ICE MRI Cryoablation System (used in conjunction with Boston Scientific accessory products) is intended to destroy tissue (including prostate and kidney tissue, liver metastases, tumors, and skin lesions) by the application of extremely cold temperatures during minimally invasive procedures.

The clinical benefit is measured by overall clinical outcomes, including acceptable safety and performance outcomes specific to the target anatomy.

CONTRAINDICATIONS

There are no known contraindications specific to use of MRI Needles.

WARNINGS

- A thorough understanding of the technical principles, clinical applications, and risk associated with cryoablation is necessary before using this product. Refer to the Visual-ICE MRI Cryoablation System User Manual for optional education.
- Do not use this device for any purpose other than the stated intended use.
- Observe the expiration date of this product. Do not use past the listed expiration date.
- Select MRI Needles appropriate for the application and tumor size. The iceball shape and size for MRI Needles are described in the **Device Description** Section.
- Do not use cryoablation needles that are not MRI compatible (labeled MR Unsafe) within an MRI suite or environment.
- Use special care to ensure that a cryoablation needle does not come into contact with an implanted device.
- Cryoablation needles have sharp tips. Use care to ensure safe handling of a cryoablation needle to eliminate the risk of injury or possible exposure to blood-borne pathogens.
- The sterile field and sterility of the cryoablation needle should be maintained at all times. Do not contaminate the distal end of the sterile cryoablation needle.
- MRI Needles are designed and indicated for Freezing and Thawing applications. These needles are not designed, indicated or tested for thermal protection. Serious injury to patient tissue may result if used for thermal protection.
- Inspect the packaging for damage. Do not use the needle if packaging appears opened or damaged; in the event of such occurrence, contact Boston Scientific Technical Assistance Center to arrange return of the complete package with the product.
- Do not use the needle if it is bent or damaged while attempting to unpack or use it. Never use a defective needle for a cryoablation procedure. A defective cryoablation needle that has a gas leak can cause a gas embolism or pneumatosis in the patient.
- Avoid bending the needle shaft. Do not grasp needles with auxiliary instruments as this may cause damage to the needle shaft.
- During use, avoid damage to the needle from other surgical instruments.
- Do not kink, pinch, cut or pull excessively on the needle tubing. Damage to needle handle or tubing may cause the needle to become unusable.

- Ensure appropriate stability of tubing to avoid inadvertent movement of tubing or needle shifting during a procedure.
- Each needle must be locked into a needle channel before initiating a cryoablation procedure to avoid the risk of forceful ejection of the needles while under gas pressure.
- Prior to starting a cryoablation procedure, set up the cryoablation system (reference the SYSTEM SETUP Section within the cryoablation system User Manual) and then perform the Needle Integrity and Functionality Test. This test must be successfully completed in order to begin the procedure.
- Do not use the needle if there is no ice formation during the Freeze phase while performing the Needle Integrity and Functionality Test. Refer to cryoablation system user manual for troubleshooting. If the issue does not resolve, obtain a new needle and repeat the testing procedure.
- Do not use the needle if bubbles are seen escaping from the needle during Needle Integrity and Functionality Test.
- Ensure adequate measures are taken to protect organs and structures adjacent to the targeted tissue.
- Continuously monitor the cryoablation procedure needle insertion, needle positioning, iceball formation and removal using direct visualization or MRI image guidance to ensure adequate tissue coverage and to avoid damage to adjacent structures.
- The MRI Needles may heat up when the scanner Radio Frequency (RF) field is active. Hazards associated with heating can be minimized by using scan sequences with low Specific Absorption Rate (SAR) and limiting the duration of the scan. It is recommended that the scanner be operated in Normal Operating Mode with a whole body averaged SAR of ≤ 1.5 watts per kilogram (W/kg). It is also recommended that the scan duration be limited to less than one minute per scan if the needle is not operating in freezing mode during the scan.
- In the rare event that a needle breaks while inserted in the tissue, act immediately to remove needle parts from the patient's body and report such event to Boston Scientific.
- If a needle unintentionally strikes bone, do not start or continue the Freezing process as needle integrity may be compromised. Replace the needle prior to continuing the procedure.
- Needle handles and the needle shaft may frost during freezing. Avoid prolonged contact with frosted portions of the needle handle to avoid unintended thermal tissue damage to the patient or clinician. The patient's skin surface should be protected by warm saline irrigation or other means as determined by the physician.
- Needle tubing may become extremely cold when conducting freeze cycles during a cryoablation procedure. It is important that a patient's skin is protected from direct contact with needle tubing to avoid the potential for thermal injury to the patient. Ensure an appropriate insulating barrier is placed as needed (such as towels) or other method is employed to prevent needle tubing from touching a patient's skin.
- The needle handle may become warm during active thawing. Prolonged contact with warm portions of the needle handle could cause unintended thermal tissue damage/burn to the patient or clinician.
- Active thawing produces heat along the distal needle shaft. Use care to avoid thermal injury/burn to nontargeted tissues.
- Ensure adequate thawing or cooling before attempting to remove needles from the patient.

- Discontinue all needle operation prior to needle removal to minimize risk of thermal injury and/or tissue injury.
- Remove needles from the patient prior to disconnecting needles from the Boston Scientific cryoablation system.
- Dispose of needle(s) in accordance with the **Disposal** Section.
- No data regarding cryoablation in combination with other therapies is available from Boston Scientific.

PRECAUTIONS

General

- Boston Scientific MRI Needles operate only with a Boston Visual-ICE MRI Cryoablation System.
- Boston Scientific recommends that only needles of the same type be placed together in a single channel. Using needles of differing types in a channel may affect the accuracy of the Gas Indicator.
- Confirm availability of sufficient gas to conduct the planned procedure: the number of needles, needle operations activated, gas cylinder size, pressure and gas flow affect the required gas volume.
- When performing MRI guided cryoablation procedures, ensure Boston Scientific MRI Needles are selected. (MRI text on the needle handle; packaging is labeled MR conditional).

Handling

- Use of multiple needles is recommended to fully cover a target site and provide a suitable margin.
- Multiple needles placed in an adjacent configuration will typically create a large, coalesced iceball. Iceball formation must be monitored using image guidance to optimize a successful ablation procedure.
- Availability of a back-up needle is recommended should a replacement or additional needle be required during a procedure.
- Ensure all connections between the Boston Scientific MRI mobile connection panel and the MRI Needles are tight.
- During a cryoablation procedure, do not immerse the proximal handle or gas tubing in fluids.
- Do not expose a cryoablation needle to organic solvents such as alcohol which may damage the needle.

Procedural

- In a safe and controlled fashion, twist the protective sleeve in a rotational motion while simultaneously pulling it away from the device handle. Take extra care when removing the protective sleeve from the device to avoid contact with the distal end of the needle.
- To minimize image artifacts, position needle tubing such that the tubing runs in the center of the MRI system. Avoid routing the tubing along the side of the MRI system and keep away from the RF coils. When routing the needle tubing, avoid looping or crossing of the tubing.

ADVERSE EVENTS

The potential adverse events associated with the device and/or cryoablation procedure include, but are not limited to:

- Allergic reaction (contrast, device, other)
- Angina
- Arrhythmia
- Atelectasis
- Bladder spasms
- Bleeding/hemorrhage
- Burn/frostbite
- Cerebrovascular accident (CVA)/stroke
- Cryoshock phenomenon (e.g. multi organ failure, severe coagulopathy, disseminated intravascular coagulation (DIC))
- Death
- Distension
- Edema/swelling
- Ejaculatory dysfunction
- Embolism (air, device, thrombus)
- Erectile dysfunction
- Fever
- Fistula
- Fracture
- Gastrointestinal symptoms (e.g. nausea, vomiting, diarrhea, constipation)
- Healing, impaired
- Hematoma
- Hematuria
- Hemothorax
- Hepatic dysfunction/failure
- Hernia
- Hypertension
- Hypotension
- Hypothermia
- Ileus
- Impotence
- Infection/abscess/sepsis
- Inflammation
- Muscle spasm
- Myocardial infarction
- Necrosis
- Need for additional intervention or surgery
- Nerve injury
- Neuropathy
- Obstruction
- Pain/discomfort

- Perforation (including organ and adjacent structures)
- Pericardial effusion
- Perirenal fluid collection
- Pleural effusion
- Pneumatosi s (air or gas in an abnormal quantity and/or place in the body)
- Pneumothorax
- Post ablation syndrome (e.g. fever, pain, nausea, vomiting, malaise, myalgia)
- Renal insufficiency/failure
- Renal parenchymal or capsule fracture
- Respiratory distress/insufficiency/failure
- Scrotal edema
- Stenosis/stricture
- Subcutaneous emphysema
- Thrombosis/thrombus
- Tissue damage
- Transient ischemic attack (TIA)
- Tumor cell seeding
- Urethral sloughing
- Urinary frequency/urgency
- Urinary incontinence
- Urinary retention
- Urinary tract infection
- Vasovagal response
- Vessel trauma (e.g. dissection, injury, perforation, pseudoaneurysm, rupture, or other)
- Wound infection

CONFORMANCE TO STANDARDS

EMF output conforms with IEC 60601-1-2 “General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests.”

HOW SUPPLIED

Device Details

Boston Scientific disposable MRI Needles are individually packaged in a sealed film-Tyvek Pouch. Each pouch is labeled Sterile, for SINGLE USE only.

NOTE: Boston Scientific cryoablation systems are supplied separately.

Needle connectors are protected by a rubber cap. This cap should be removed before connecting a needle to a cryoablation system.

Do not use if package is damaged or unintentionally opened before use.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the MRI Needles are MR Conditional. A patient undergoing a cryoablation procedure may safely undergo MRI scanning under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name(s)	• IceRod Cryoablation Needle, 1.5 MRI Straight • IceRod Cryoablation Needle, 1.5 MRI 90° • IceSeed Cryoablation Needle, 1.5 MRI Straight • IceSeed Cryoablation Needle, 1.5 MRI 90°
Static Magnetic Field Strength (B ₀)	1.5 T or 3.0 T
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	1.5 W/kg (Normal Operating Mode)
Scan Duration	1.5 W/kg whole-body average SAR for 1 minute of continuous RF (a sequence or back to back series/scan without breaks) if not operating the needle in freezing mode.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

OPERATIONAL INSTRUCTIONS

NOTE: For information regarding software control of needle performance, refer to the Boston Scientific Visual-ICE MRI Cryoablation System User Manual.

Needle Preparation

1. Using aseptic technique, carefully remove the cryoablation needle from the package and place in a sterile work area. In a safe and controlled fashion, remove the protective sleeve by twisting the protective sleeve in a rotational motion while simultaneously pulling it away from the device handle.
2. To conduct the Needle Integrity and Functionality Test, prepare a large basin containing sterile water or saline that is sufficiently able to accommodate the full needle shaft length.
3. Remove the connector cap, then connect the needle to the selected port.
4. Place the needles, individually or in groups, in the basin such that the full length of the needle is submerged in the sterile water or saline.

The needle is now ready for the Needle Integrity and Functionality Test.

NOTE: For detailed instructions on connecting needles to the system’s Needle Connection Panel and performing the Needle Integrity and Functionality Test, refer to the Boston Scientific Visual-ICE MRI Cryoablation System User Manual.

Needle Use

Needle Handling and Insertion

- Correct insertion of cryoablation needles into the target tissue is the responsibility of the physician.

NOTE: Although the needle has a sharp tip, a small skin incision may be made at the point of needle insertion.

- Always use two hands and support the needle mid-shaft with two fingers to eliminate the risk of bending. Do not insert the needle into tissue while holding the handle with one hand only.
- Insertion depth may be estimated using the needle marks on the shaft. Use image guidance as necessary to guide needle insertion and placement.
- Use image guidance to verify that the cryoablation needle is placed at the desired location prior to activating the needle.
- After needle insertion, route the needle tubing in a manner that ensures it is not in contact with the patient's skin.
- Place an appropriate insulating barrier (such as towels) or employ other methods to prevent needle tubing from touching a patient's skin.
- The needle handle becomes cold during cryoablation. If the handle is in contact with skin, the skin surface should be protected by warm saline irrigation or other means as determined by the physician.

Notes for Conducting Freezing

- Select Freeze intensity and initiate Freezing.
- Continue Freeze cycle for duration and intensity necessary to optimize cryoablation of the target site.

NOTE: See the Boston Scientific Visual-ICE MRI Cryoablation System User Manual for instructions on system controls available to manage each Freeze cycle.

WARNING: Continuously monitor iceball formation using direct visualization or MRI guidance to ensure adequate tissue coverage and to avoid damage to adjacent structures.

- The needle shaft becomes cold during cryoablation. Protect the patient's skin by using warm saline irrigation or other means as determined by the physician.
- Pay attention to needle handle position. Prolonged contact with frosted portions of the needle handle could cause unintended thermal tissue damage to the patient or clinician.
- Avoid allowing needle tubing to be in direct contact with a patient's skin. Direct contact with needle tubing could cause unintended thermal damage to the patient's skin. Place an appropriate insulating barrier (such as towels) or employ other methods to prevent needle tubing from touching a patient's skin.

Notes for Conducting Thawing

NOTE: See the Boston Scientific Visual-ICE MRI Cryoablation System User Manual for instructions on system controls available to manage each Thaw cycle.

Needle Removal

- Thaw thoroughly and stop all needle operation prior to removing needles to minimize the risk of tissue injury.

NOTE: Boston Scientific's MRI Needles are specially designed with a sharp, conical tip to minimize bleeding. However, some bleeding may occur. In the event of bleeding, apply treatment in accordance with good clinical practice and the hospital's treatment protocol. For example, following needle removal, hold compression until hemostasis is achieved; if necessary, place an appropriate dressing on the needle insertion site.

Disposal

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

For customers in Australia, report any serious incident that occurs in relation to this device to Boston Scientific and to the Therapeutic Goods Administration (<https://www.tga.gov.au>).

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling the patient on the use of the cryoablation system in association with the interventional procedure:

- Discuss the risks and benefits, including review of potential adverse events listed in this IFU and the cryoablation system User Manuals, of other interventional treatments likely to be employed.
- Discuss post-procedure instructions, including any lifestyle changes, medications, and home care and rehabilitation guidelines.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

IceRod, IceSeed, and Visual-ICE MRI are trademarks of Boston Scientific Corporation or its affiliates. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU in Adobe® Portable Document Format (PDF), and be sure to have the product label available for reference.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.