Innova™

OVER-THE-WIRE

Self-Expanding Stent System

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REUSE WARNING 3

DEVICE DESCRIPTION 3

 Table 1. Innova Self-Expanding Stent System 3

 Contents 3

 Operating Principle 3

 Device Component Description 3

 Product Matrix 3

 Materials 3

 Table 2. Innova Self-Expanding Stent System Models and Sizes 3

 Non-Pyrogenic 4

 User Information 4

INTENDED USE/INDICATIONS FOR USE 4

CLINICAL BENEFIT STATEMENT 4

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE 4

CONTRAINDICATIONS 4

WARNINGS 4

 Pre- and Post-Procedure Antiplatelet Therapy 4

PRECAUTIONS 4

MAGNETIC RESONANCE IMAGING (MRI) 5

 3.0 Tesla Temperature Information 5

 1.5 Tesla Temperature Information 5

 Image Artifact 5

 Recommendations 5

ADVERSE EVENTS 5

HOW SUPPLIED 5

 Device Details 5

 Handling and Storage 5

OPERATIONAL INSTRUCTIONS 6

 Additional Items for Safe Use 6

 Patient Preparation 6

 Inject Contrast Media 6

 Evaluate and Mark the Stenosis 6

 Select Proper Stent System 6

 Preparation of Stent Delivery System 6

 Figure 1. Stent Delivery System 6

 Delivery Procedures 7

 Stent Deployment Procedure (Reference Figure 1) 7

 Figure 2. Eliminate slack 7

Recommended Method of Deployment 7

 Figure 3. Long stents (150 mm, 180 mm, and 200 mm) require the pull grip to be retracted only after the white activation arrow is visible to complete deployment. 8

Disposal 8

Post-Procedure 8

Implantable Device Patient Information 8

INFORMATION TO BRIEF THE PATIENT 8

 Expected Lifetime 8

WARRANTY 8

Self-Expanding Stent System

ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Innova Self-Expanding Stent System is comprised of two components: the implantable endoprosthesis and the stent delivery system. The characteristics of the Innova Stent System are described in Table 1.

Table 1. Innova Self-Expanding Stent System

Innova Self-Expanding Stent System	
Available Stent Lengths (mm)	20, 40, 60, 80, 100, 120, 150, 180, 200
Available Stent Diameters (mm)	5, 6, 7, 8
Stent Material	Nickel Titanium Alloy (NiTi)
Delivery System Effective Length	75 cm, 130 cm

Contents

- One (1) Innova Stent with Delivery System

Operating Principle

Device Component Description

The stent system is comprised of: the implantable endoprosthesis and the stent delivery system. The stent is a laser cut self-expanding stent composed of a nickel titanium alloy (nitinol). On both the proximal and distal ends of the stent, radiopaque markers made of tantalum increase visibility of the stent to aid in placement. The stent is constrained within a 6F (2.1 mm maximum OD) delivery system. The stent is loaded on the delivery system. The delivery system is a triaxial design with an outer shaft to stabilize the stent delivery system, a middle shaft to protect and constrain the stent, and an inner shaft to provide a guidewire lumen. The delivery system is compatible with 0.035 in (0.89 mm) guidewires.

The Innova Stent System is provided sterile. When ready to be implanted, the stent is deployed by retracting the middle sheath of the delivery system via actuation of the thumbwheel. As the stent is exposed to body temperature, it expands to appose the vessel wall. The 150 mm, 180 mm, and 200 mm length stents require a two step deployment method whereas the shorter lengths can be fully deployed with the thumbwheel only. Refer to Figure 1. for product illustration.

Product Matrix

The Boston Scientific Innova Self-Expanding Stent System (Innova Stent System) is available in a variety of stent diameters and stent lengths, offered in 75 cm and 130 cm working lengths (the length of the finished delivery system that can be inserted into the body), and is compatible with stiff 0.035 in (0.89 mm) guide wires and 6F (2.0 mm) introducer sheaths.

Materials

The Innova Stent System is a self-expanding nitinol (nickel-titanium alloy) stent structure with tantalum radiopaque markers. Below are the materials used in the Innova Stent by percent weight.

Nitinol 99.99

Tantalum <0.1

Table 2. Innova Self-Expanding Stent System Models and Sizes

Stent Nominal Diameter (mm)	Unconstrained Length (mm)	Reference Vessel Diameter (mm)	Stent Length at Vessel Diameter (mm)
5	20	4.0	20
	40		41
	60		58
	80		78
	100		99
	120		118
	150		148
	180		174
6	200	4.0 – 5.0	194
	20		20
	40		40
	60		60
	80		79
	100		101
	120		119
	150		150
	180		178
	200		198

Stent Nominal Diameter (mm)	Unconstrained Length (mm)	Reference Vessel Diameter (mm)	Stent Length at Vessel Diameter (mm)
7	20	5.0 – 6.0	20
	40		39
	60		60
	80		79
	100		100
	120		118
	150		149
	180		176
	200		198
8	20	6.0 – 7.0	20
	40		39
	60		59
	80		80
	100		99
	120		120
	150		149
	180		179
	200		197

Non-Pyrogenic

Innova Self-Expanding Stent System meets pyrogen limit specifications.

User Information

Only physicians, technicians, and nurses experienced in preparing for and performing peripheral vascular procedures should use this device.

INTENDED USE/INDICATIONS FOR USE

The Innova Vascular Self-Expanding Stent System is indicated to improve luminal diameter in the treatment of adults with symptomatic de-novo or restenotic lesions in the native superficial femoral artery (SFA) and/or proximal popliteal artery (PPA) with reference vessel diameters from 4.0 mm to 7.0 mm and lesion lengths up to 190 mm.

CLINICAL BENEFIT STATEMENT

The Innova Self-Expanding Stent System is designed to improve luminal diameter in the treatment of symptomatic de-novo or restenotic lesions in the native superficial femoral artery (SFA) and/or proximal popliteal artery (PPA). The clinical benefit can be measured by overall clinical outcomes as demonstrated by primary patency rates, freedom from amputation, freedom from TLR, and overall survival with respect to other existing therapies.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (EUDAMED) website: (<https://ec.europa.eu/tools/eudamed>).

CONTRAINDICATIONS

The Innova Self-Expanding Stent System is contraindicated for use in any situation in which percutaneous transluminal angioplasty (PTA) is contraindicated.

WARNINGS

- The delivery system is not designed for use with power injection systems.
- Do not use a kinked delivery system.
- Only advance the stent delivery system over a recommended guidewire. Use of other guidewires may lead to deployment difficulties resulting in adverse event or need for urgent intervention/surgery.
- When catheters are in the body, they should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed.
- The stent delivery system is not intended for arterial blood monitoring.
- Prior to completion of the procedure, utilize fluoroscopy to ensure proper positioning of the stent. If the lesion is not fully covered, use additional stents as necessary to adequately treat the lesion.
- The minimally acceptable introducer or guide sheath size is printed on the package label. Do not attempt to pass the stent delivery system through a smaller size introducer or guide sheath than indicated on the label.
- In the event of thrombosis of the expanded stent, thrombolysis and/or PTA should be considered.
- In the event of complications such as infection or vessel trauma, surgical removal of the stent may be required.
- Recrossing a partially or fully expanded deployed stent with adjunct devices must be performed with extreme caution to ensure that the adjunct device does not get caught within previously placed stent struts.
- Do not remove the thumbwheel lock prior to deployment. Premature removal of the thumbwheel lock may result in an unintended deployment of the stent.
- In the event of stent deployment difficulties (e.g. partial deployment), manipulation of the device, removal/replacement of the delivery system, or urgent medical or surgical intervention may be required.
- This product should not be used in patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Persons with a known allergy to nickel or titanium may suffer an allergic response to this implant.

Pre- and Post-Procedure Antiplatelet Therapy

The device carries an associated risk of thrombosis. It is strongly advised that the treating physician follow the Inter-Society Consensus (TASC II) Guidelines recommendations (or other applicable country guidelines) for antiplatelet therapy pre- and post-procedure to reduce the risk of thrombosis.

PRECAUTIONS

- Do not use after the "Use By" date specified on the package. See How Supplied section prior to use.
- Stenting across a bifurcation or side branch could compromise future diagnostic or therapeutic procedures.
- The stent is not designed for repositioning.
- Once the stent is partially deployed, it cannot be "recaptured" or "reconstrained" using the stent delivery system.

- The stent may cause an embolism from the site of the implant down the arterial lumen.
- Do not use if the temperature exposure indicator dot on the pouch label is red, indicating that stent expansion may have been compromised.
- Do not use if the temperature exposure indicator dot on the pouch label is missing.
- Do not expose to organic solvents (e.g. alcohol).
- Should more than one stent be required, allow for at least 5 mm of stent overlap.
- When multiple stents are required, if placement results in metal to metal contact, stent materials should be of similar composition.

MAGNETIC RESONANCE IMAGING (MRI)

Non-clinical testing has demonstrated the Innova Stent is MRI Conditional. It can be scanned safely up to a total length of 200 mm and overlapping stents up to 200 mm under the following conditions:

- Static magnetic field of 3 Tesla and 1.5 Tesla.
- Static magnetic gradient field of 100 Tesla/meter (extrapolated).
- Normal operating mode only with a maximum whole body (WB) averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning for patient landmarks above the umbilicus (patient navel).
- Maximum whole body (WB) SAR of 0.41 W/kg for 15 minutes of scanning for patient landmarks below the umbilicus.
- Use whole body transmit coil only. Do not use local transmit coils. Local receive coils can be used. MRI at 3 T or 1.5 T may be performed immediately following the implantation of the Innova Stent. The Innova Stent should not migrate in this MRI environment. This stent has not been evaluated to determine if it is MRI Conditional beyond these conditions.

3.0 Tesla Temperature Information

In non-clinical testing, the Innova Stent at single lengths of 200 mm and overlapped lengths of 200 mm produced a maximum temperature rise of less than 6.58 °C at a maximum whole body averaged of 0.49 W/kg, that was determined by validated calculation for 15 minutes of MR scanning in a 3 Tesla Siemens Magnetom Trio, software version Numaris/4, Syngo MR A30, N4 VA30A Latest 20070315 PT, MR scanner. In this model, the reported temperatures are conservative as they do not take into account the cooling effects of perfusion and blood flow.

- For landmarks above the umbilicus the calculated temperature rise was 2.9 °C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.
- For landmarks below the umbilicus the calculated temperature rise was 6.58 °C for a whole body average SAR value of 0.49 W/kg and a continuous scan time of 15 minutes.

1.5 Tesla Temperature Information

In non-clinical testing, the Innova Stent at single lengths of 152 mm (as measured) and overlapped lengths of 200 mm produced a maximum temperature rise of less than 4.5 °C at a maximum whole body averaged specific absorption rate (SAR) of 0.41 W/kg, that was determined by validated calculation for 15 minutes of MR scanning in a 1.5 Tesla Philips Intera, software version Release 10.6.2.5, 2006-03-10 MR scanner. In this model, the reported temperatures are conservative as they do not take into account the cooling effects of perfusion and blood flow.

- For landmarks above the umbilicus the calculated temperature rise was 1.6 °C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

- For landmarks below the umbilicus the calculated temperature rise was 4.48 °C for a whole body average SAR value of 0.41 W/kg and a continuous scan time of 15 minutes.

Image Artifact

The image artifact extends approximately 5 mm from the perimeter of the device diameter and 1.5 mm beyond each end of the length of the stent when scanned in non-clinical testing using the sequence, Spin Echo. With a Gradient Echo sequence the image artifact extends 12 mm from the perimeter of the device diameter and 1.8 mm beyond each end of the length of the stent with both sequences partially shielding the lumen in a 3.0 Tesla Achieva (Achieva Upgrade), Philips Medical Solutions, software version Release 2.5.3.0 2007-09-28 MR system with a Quadrature transmit/receive head coil. Image artifacts in a body birdcage coil are similar to the image artifacts in the transmit/receive CP head coil.

Recommendations

It is recommended that patients register the conditions under which the implant can be scanned safely with the MedAlert Foundation (www.medicalert.org) or equivalent organization.

ADVERSE EVENTS

Potential adverse events which may be associated with the use of a peripheral stent include but are not limited to:

- Allergic reaction (to drug, contrast, device or other)
- Bleeding/Hemorrhage
- Death
- Embolism (air, plaque, thrombus, device, tissue, or other)
- Hematoma
- Hypotension/hypertension
- Ischemia/necrosis
- Need for additional intervention or surgery
- Renal insufficiency or failure
- Restenosis of stented artery
- Sepsis/infection
- Thrombosis/thrombus
- Vasospasm
- Vessel occlusion
- Vessel trauma (perforation, pseudoaneurysm, injury, rupture, and dissection)

HOW SUPPLIED

The Innova Self-Expanding Stent System is supplied sterile inside a pouch. The device is sterilized via Ethylene Oxide.

Device Details

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.

Handling and Storage

The Innova Stent is a nitinol stent that has an upper temperature limit of 51 °C (124 °F).

Precaution: Do not use if the temperature indicator dot on the carton or pouch is red indicating that stent expansion may have been compromised.

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

Recommended Materials (not included in stent system package)

- 0.035 in (0.89 mm) guidewire of appropriate length (300 cm length recommended for 130 cm length stent delivery systems)
- Introducer or guide sheath of appropriate size and length and equipped with a hemostatic valve
- Luer lock syringe 10 ml (10 cc) for prepping the stent delivery system

Patient Preparation

The percutaneous placement of a self-expanding stent in a stenotic or obstructed artery should be done in an angiography procedure room with the appropriate imaging equipment. Patient preparation and sterile precautions should be the same as for any angioplasty procedure. Appropriate antiplatelet and anticoagulation therapy must be administered pre- and post-procedure in accordance with standard practices. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. Access vessels must be sufficiently patent to proceed with further intervention. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice.

Inject Contrast Media

Perform angiogram using standard technique.

Evaluate and Mark the Stenosis

Observe fluoroscopically the most distal view of the stenotic or obstructed artery.

Obtain a road map image of the lesion area.

Select Proper Stent System

1. Measure the diameter of the reference vessel (proximal and distal to the lesion or obstruction). Select a stent based on **Table 2**, to ensure a secure placement.
2. Measure the entire length of the actual lesion and select the proper length of the stent(s) to be deployed. To help ensure adequate apposition, it is recommended

that the length of the stent be chosen so that the ends of the stent extend at least 5 mm beyond both ends of the lesion into healthy tissue.

Precaution: Should more than one stent be required to cover the lesion, allow for at least 5 mm of stent overlap. It is generally recommended that the distal stent be placed first. When multiple stents are required, if placement results in metal to metal contact, stent materials should be of similar composition.

3. Estimate the distance between the lesion and the entry site to select the proper stent delivery system length.

Preparation of Stent Delivery System

1. Open the outer box to reveal the pouch containing the stent delivery system.
2. Check the temperature exposure indicator on the pouch label to confirm that the product has not been compromised. See Precautions section.
3. After careful inspection of the pouch looking for damage to the sterile barrier carefully peel open the pouch and extract the stent delivery system tray.
4. Carefully withdraw the stent delivery system from the tray by grasping the handle of the delivery system.
5. Examine the stent delivery system for any damage. If it is suspected that the sterility or integrity of the device has been compromised (i.e. kinking or missing component), the device should not be used. The device should not be used if the device is kinked, or if the thumbwheel lock is not attached.
6. Do not remove the thumbwheel lock [2] prior to deployment. Premature removal of the thumbwheel lock may result in an unintended deployment of the stent.
7. Attach a 10 ml (10 cc) syringe filled with saline to the flushing luer [6] on the handle. Apply positive pressure. Continue to flush until saline appears at the distal end of the guidewire lumen.
8. Remove the flushing luer [6] (by pulling the syringe or by pulling flushing luer [6]) (Reference **Figure 1**).


Figure 1. Stent Delivery System

Delivery Procedures

1. Gain arterial access utilizing a 6F (2.1 mm) or larger sheath with a hemostatic valve.

Precaution: Always use an introducer or guide sheath for the implant procedure to protect the access site and prevent system damage.

Precaution: Do not use a kinked delivery system. Kinking of the introducer/guide sheath at the access site can restrict the movement of the delivery system during deployment.

2. Pass a 0.035 in (0.89 mm) guidewire of appropriate length (300 cm length recommended for 130 cm length stent delivery systems) across the target lesion or obstruction.

Note: A **stiff** 0.035 in guidewire is strongly recommended for deployment of the stent, especially for tortuous anatomy and contralateral approaches. Use of undersized guidewires may lead to insufficient support of the device which can compromise stent delivery.

Note: If using a hydrophilic guidewire, ensure that it is hydrated at all times.

3. Pre-dilate the lesion with a balloon dilatation catheter using conventional technique. After the lesion has been properly dilated, remove the dilatation catheter, leaving the guidewire with the tip distal to the lesion for stent system advancement.

Precaution: Physicians should use judgement based on experience in dilating arterial lesions and/or obstructions. Never force a balloon catheter to inflate to the point of risking dissection of the arterial wall.

4. Place the Innova Stent Delivery System over the guidewire. Advance the delivery system as a unit through the hemostatic valve of the introducer or guide sheath.

Note: Do not tighten trolley-burst such that it restricts the movement of the delivery system.

Stent Deployment Procedure (Reference Figure 1)

1. Remove slack from the system by advancing the system just beyond the target lesion, then, pulling the system back until stent radiopaque markers [1] are centered over the target lesion.

Note: Prior to deployment, ensure adequate distance between the proximal end of the stent and the introducer/guide sheath to prevent deployment within the introducer/guide sheath.

2. Remove the thumbwheel lock [2] by compressing the tabs and pulling. Confirm that the radiopaque markers are still properly positioned across the target lesion.

Precaution: If strong resistance is met with the introduction of the delivery system or if unable to initiate release of the stent, remove the entire system from the patient and introduce a new system.

Note: For optimal performance, keep the entire length of the delivery system that is outside the body as straight and stable as possible. To do so, remove slack from the system, maintain slight backward tension on the delivery system, and anchor the handle on the patient or operating table during deployment. Alternatively, the operator may straighten and stabilize the distal end of the blue outer shaft during deployment.

Note: Failure to eliminate slack (Reference **Figure 2**) and/or curvature of the delivery system catheter between the introducer/guide sheath and the delivery system handle during deployment may adversely affect deployment accuracy, especially in ipsilateral cases.

Note: If repositioning of the stent delivery system is required, reinserting the thumbwheel lock will prevent inadvertent deployment.


Figure 2. Eliminate slack

Note: Stabilize the deployment handle by keeping the handle firmly placed on the patient or operating table to prevent inadvertent pulling of the system or introduction of slack during deployment.

Recommended Method of Deployment

1. While using fluoroscopy maintain position of the distal and proximal stent radiopaque markers [1] relative to the targeted site. Roll the thumbwheel [4] of the deployment handle in the direction of the arrow indicated on the handle. Continue to roll thumbwheel until the middle shaft radiopaque marker band [5] passes the distal stent radiopaque markers. Watch for the distal stent radiopaque markers to begin separating: separation of the distal stent radiopaque markers signals that the stent is deploying.
2. Continue to roll thumbwheel until the middle shaft radiopaque marker band [5] passes the proximal radiopaque markers of the stent resulting in full deployment, or until the white activation arrow is visible on pull grip extension rod (for 150 mm, 180 mm, and 200 mm length stents), which indicates that pull grip activation is required to complete stent deployment (Reference **Figure 3**).

Note: When activating the pull grip, avoid rapid deployment.

Note: Do not restrict movement of the thumbwheel [4] otherwise deployment difficulties could be encountered. Do not attempt to pull a partially expanded stent back into introducer/guide sheath as dislodgement may occur.

Note: Do not push or pull the delivery system during deployment as this may compromise stent length.


Figure 3. Long stents (150 mm, 180 mm, and 200 mm) require the pull grip to be retracted only after the white activation arrow is visible to complete deployment.

3. Long stents (150 mm, 180 mm, and 200 mm) require pull grip deployment after the white activation arrow becomes visible on the pull grip extension rod. Grasp the manual pull grip [3] and gently pull away from the handle in the direction of the arrow. Slowly pull back until the middle shaft radiopaque marker band [5] passes the proximal radiopaque markers of the stent resulting in full deployment.
4. View the delivery system under fluoroscopy, ensuring that the middle shaft radiopaque marker band [5] has crossed the proximal stent markers. The delivery system can now be withdrawn.
5. Grasp the guidewire a short distance from handle and repeatedly retract the system over the wire until fully removed. Use caution when withdrawing the stent delivery system and always manipulate under fluoroscopy. If unusual resistance is felt, carefully readvance and rotate the delivery system in an attempt to center the delivery system within the vessel, then carefully attempt to repeat withdrawal.

Note: Avoid bending the guidewire excessively near handle when retracting device to aid removal and prevent guidewire kinking.

6. If incomplete expansion exists within the stent at any point along the lesion, balloon dilatation can be performed utilizing standard PTA technique.

Precaution: Never post-dilate the stent using a balloon that is larger in diameter than the nominal (labeled) diameter of the stent.

7. Withdraw guidewire and sheath from patient and establish hemostasis per conventional technique.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose of device and packaging as follows:

After use, device may contain biohazardous substances. Devices that contain biohazardous substances are considered biohazardous waste and should be stored in a biohazard container that is labeled with the biological hazard symbol. Untreated biohazardous waste should not be disposed of in the municipal waste system. Biohazardous substances should be removed via chemical disinfection prior to disposal.

Alternatively, biohazardous waste may be disposed of utilizing a certified facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Assess patient for hematoma and/or other signs of bleeding at the puncture site.

Determine appropriate antiplatelet therapy based on the Inter-Society Consensus (TASC II) Guidelines recommendations (or other applicable country guidelines) for antiplatelet therapy post-procedure to reduce the risk of thrombosis.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

Implantable Device Patient Information

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

INFORMATION TO BRIEF THE PATIENT

The physician should consider the following points while counseling patients on the use of the Innova Stent in association with the interventional procedure:

- Discuss the risks and benefits including review of potential adverse events listed in this document, both the Innova Stent and for other interventional treatments likely to be employed.
- Discuss patient allergies in particular the risk for patients who may be allergic to nickel and/or titanium.
- Discuss the risks and benefits of antiplatelet therapy including risk of thromboembolism should the patient be allergic or discontinue use.
- Discuss post procedure instructions, including any follow-up appointments, lifestyle changes, medications, and home-care or rehabilitation guidelines.
- Provide the patient with the completed implant card to carry and advise the patient that additional information, including MRI conditions, may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).
- Instruct the patient to present the implant card to their Healthcare professionals (doctors, dentist, technicians) so they can take the necessary precautions.

Expected Lifetime

Inform the patient that the stent is a permanent implantable and has been tested for structural integrity (fracture resistance) for a minimum of 10 years; however, the materials of the device are nonbiodegradable and are intended to last for the lifetime of the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

Innova is a registered trademark of Boston Scientific Corporation or its affiliates.

Magnetom Trio is a trademark of Siemens Aktiengesellschaft Corporation.

Syngo is a trademark of Siemens Aktiengesellschaft Corporation.

Intera is a trademark of Koninklijke Philips Electronics N.V. Corporation.

All other trademarks are the property of their respective owners.

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INSTRUCTIONS FOR USE

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU in Adobe® Portable Document Format (PDF), and be sure to have the product label available for reference.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
© 2022 Boston Scientific Corporation or its affiliates.