AdVance™ XP

Male Sling System

REF 720163-03

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2022-11

TABLE OF CONTENTS

REUSE WARNING	3
DEVICE DESCRIPTION	3
Contents	3
Operating Principle	4
Materials	4
Animal Tissues or Derivatives	4
MRI Safety Information	4
User Information	4
INTENDED USE/INDICATIONS FOR USE	5
Clinical Benefit Statement	5
Summary of Safety and Clinical Performance	5
CONTRAINDICATIONS	5
WARNINGS	5
Procedural	5
Device Related	5
PRECAUTIONS	5
Procedural	6
Post Procedure	6
Device Related	6
ADVERSE EVENTS	6
HOW SUPPLIED	7
Handling and Storage	7
OPERATIONAL INSTRUCTIONS	7
Procedure	7
Patient Preparation	7
Dissection	7
Placing the Sling	8
Attaching the Sling	9
Tensioning the Sling	9
Disposal	10
Post Procedure	10
IMPLANTABLE DEVICE PATIENT INFORMATION	10
IMPLANT CARD INSTRUCTIONS	10
PATIENT COUNSELING INFORMATION	10
Device Lifetime	11
Duration of Use	11
WARRANTY	11
SYMBOL DEFINITIONS	11

AdVance™ XP

Male Sling System

Rx ONLY

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The AdVance XP Male Sling System consists of two sterile, single-use surgical instruments called "needle passers" and a mesh sling with attached connectors. One end of each needle passer is keyed to allow a secure attachment of the sling connector. The sling mesh is constructed of polypropylene monofilament mesh that is pre-cut. The sling arms contain polypropylene tissue anchors molded directly to the sling mesh. Polypropylene tensioning fibers are integrated into the mesh weave. Plastic sheaths are placed over each arm of the mesh to aide in ease of placement. The sheaths contain liners to protect the plastic sheath from tearing upon removal. The connectors are attached to the ends of the needle passers during the procedure. The mesh and tissue anchors of the sling are intended to remain in the body as a permanent implant. The materials used in the sling are non-absorbable.

The AdVance XP Male Sling System is intended to provide a minimally-invasive surgical treatment for male stress urinary incontinence by providing additional support to the urethral bulb via an implanted mesh sling. The AdVance XP Male Sling System is centered under the corpus spongiosum, and secured by sutures at the corpus spongiosum with the arms and tissue anchors of the mesh sling embedded in the soft tissues of the pelvis.

Contents

- (1) Male Sling
- (2) Needle Passers


Figure 1. Needle Passers


Figure 2. AdVance XP Male Sling

Operating Principle

The AdVance XP Male Sling is a suburethral mesh implant designed to treat male incontinence by repositioning the urethra. AdVance XP needle passers assist in the placement of the male sling. AdVance XP is an implanted, surgical device that penetrates inside the body. Its intended purpose is for the treatment of male stress urinary incontinence (SUI) by the placement of a suburethral sling.

Materials

The permanent implantable portion of the device consists of a polypropylene monofilament mesh and polypropylene tissue anchors that are directly molded into the sling arms of the mesh. The material quantity of the mesh and tissue anchors calculated based on total weight percent is listed below:

Mesh and Tissue Anchor Permanent Implantable Material	Weight Percent
Polypropylene	98 %
Stabilizer	2 %

The AdVance XP needle passers are constructed of stainless steel and therefore may contain the following substance defined as Carcinogenic, Mutagenic and toxic to Reproduction (CMR 1B) in a concentration above 0.1 % by weight:

Contains Cobalt; CAS No. 7440-48-4; EC No.231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects. A device specific evaluation has determined that the presence of cobalt does not present a risk within the clinical use of this device.

This device meets pyrogen limit specifications.

Animal Tissues or Derivatives

The device and packaging contain tallow derivatives for which the transmissible spongiform encephalopathies and other pathogenic risk is deemed acceptable based on the standards that apply.

MRI Safety Information

The AdVance XP Male Sling is composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic and poses no known hazards resulting from exposure to any MR environment. Therefore, the AdVance XP Male Sling is considered MR Safe.

User Information

Precaution: Read the entire Instructions for Use in detail before using the AdVance XP Male Sling System. The system should only be used by or under the supervision of physicians with an advanced understanding of sub-urethral sling devices and disease states that result in incontinence. A thorough understanding of the technical principles, clinical applications, and risk associated with this procedure is necessary before using this device.

INTENDED USE/INDICATIONS FOR USE

The AdVance XP Male Sling System is intended for the treatment of male stress urinary incontinence (SUI) by the placement of a suburethral sling.

Clinical Benefit Statement

Improved Quality of Life for male patients with stress urinary incontinence (SUI), assessed using validated questionnaires.

Summary of Safety and Clinical Performance

For customers in the European Union, visit (<https://ec.europa.eu/tools/eudamed>) to view a summary of this device's safety and clinical performance information on the EUDAMED website.

CONTRAINDICATIONS

Contraindications associated with the use of the AdVance XP Male Sling System include patients with:

- Urinary tract infections.
- Blood coagulation disorders.
- A compromised immune system or any other condition that would compromise healing.
- Renal insufficiency and upper urinary tract relative obstruction.

WARNINGS

- Only physicians who have been trained as to its use should use the AdVance XP Male Sling System.
- The AdVance XP Male Sling System must not be re-sterilized or re-used.
- If AdVance XP Male Sling System sterile barrier package is found unintentionally opened or damaged before use, sterility may be compromised, rendering the device nonsterile.
- Do not use product beyond the indicated expiration date.
- Physicians should be familiar with surgical procedures and techniques involving non-absorbable meshes before using the AdVance XP Male Sling System device.
- Good surgical practice should be followed for management of contaminated or infected wounds.
- Urinary tract infection should be treated prior to implantation.

Procedural

- Do not insert the AdVance XP Male Sling System needles with an approach that deploys the needle from the perineum to the obturator.
- Take care to avoid vessel perforation. Observe patient for any signs of retropubic or periurethral bleeding.

Device Related

- If non-functional instruments, or potentially compromised sterile barrier is found; do not use the device, return it to Boston Scientific and call your Boston Scientific representative.

PRECAUTIONS

- The implantation of this device should only be considered for adult patients whom the physician determines are appropriate surgical candidates.
- The risks and benefits of using the AdVance XP Male Sling System procedure should be carefully considered on patients with:
 - Blood coagulation.
 - Compromised immune systems or any other conditions that would compromise healing.
 - Renal insufficiency and upper urinary tract obstruction.
- The possibility of urgency incontinence should be carefully considered before a sling implant is conducted. It is recommended that:
 - Good bladder function (bladder capacity > 250 ml and post void residual urine < 50 ml) be demonstrated by candidates for a male sling.

- The presence of bladder neck or urethral strictures be ruled out for male sling candidates.
- A condition involving cystitis, urethritis or prostatitis be ruled out for male sling candidates.
- Detrusor instability of a neurological origin be ruled out for male sling candidates.
- A 6 month period of non-invasive treatment (e.g., behavior modification, bladder exercises, biofeedback, extra corporeal magnetic stimulation of the pelvic floor, or drug therapy) is recommended before a sling implant is considered for males with stress urinary incontinence.

Procedural

- Check for bleeding and pelvic organ dysfunction due to intraoperative vessel or nerve damage associated with anomalous location, during needle passage through the medial area of the obturator foramen membrane.
- It is important to verify the tension and placement of the sling prior to closure.
- Do not remove the plastic sheath until the sling is in its desired position.
- If during the procedure it becomes necessary to remove the sling, cut the sling near its midline and draw the sling arms out of the patient by pulling on the distal (connector) ends.
- Do not contact the sling with any staples, clips, or other instruments which may damage the mesh.

Post Procedure

- Acute inflammatory tissue reaction and transitory local irritation has been reported with the use of the absorbable suture.
- After several (2 weeks - 3 weeks) weeks post-implant, the mesh becomes fixated by in-growth of surrounding tissues and cannot be entirely removed. However, the urethral repositioning may be reversed by cutting through the center of the mesh to release the tension.
- Patients should be counseled on abstaining from heavy lifting, strenuous exercise, and intercourse for a minimum of 6 weeks. Patients may return to other normal daily activities at the physician's discretion, often one to two weeks post procedure.
- Advise the patient that additional information may be available to them on the Boston Scientific website: www.bostonscientific.com/patientlabeling
- If dysuria, bleeding, or other problems occur, the patient should be instructed to call the surgeon immediately.

Device Related

- To maintain sterility, only the contents of the tray package should be introduced into the sterile field.

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

ADVERSE EVENTS

Applicable to males with stress urinary incontinence

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the implant could include extrusion, erosion through the urethra or other surrounding tissue, migration of the device from the desired location, fistula formation and inflammation. The occurrence of these responses may require removal of the entire sling.
- Like all foreign bodies, the sling may allow an existing infection to propagate.
- Incorrect sling tension may cause temporary or permanent lower urinary tract obstruction and retention of urine which may require surgical intervention.
- Unintended damage to muscle, tissue, nerves, or vessels.

Known risks of surgical procedures for the treatment of incontinence include: Allergic reaction, prolonged procedure, scar tissue, toxicity, hematoma, perforation, exposure to biohazardous material, unretrieved device fragment, pain, infection, erosion, device migration, bleeding, pelvic organ dysfunction, injury (urethra, corpus spongiosum, nerve) and malposition, and complete failure of the procedure, resulting in continued incontinence and mild to moderate urinary incontinence due to incomplete support or overactive bladder.

HOW SUPPLIED

The Male Sling System containing two primary devices: (1) sling and (2) needle passers. The Male System is provided in the AdVance XP Kit, which is a procedure pack.

Contents supplied STERILE using an ethylene oxide (EO) process.

Do not use if package is unintentionally opened or damaged before use.

Do not use if labeling is incomplete or illegible.

Do not use product beyond the indicated expiration date.

Handling and Storage

Recommended storage at $\leq 25^{\circ}\text{C}$ (77°F); excursions permitted to 55°C (131°F).

OPERATIONAL INSTRUCTIONS

The AdVance XP Male Sling System procedure can be carried out under spinal or general anesthesia. A small midline perineal incision is made followed by periurethral dissection. Two small stab incisions are also made above the obturator foramen for needle passer entry.

Procedure

Open the box. Remove all tray packages from the box. Verify that the package integrity has not been compromised during shipping. Inspect the device for damage.

Patient Preparation

1. Place the patient in a dorsal lithotomy position, with legs bent at approximately 90° and slightly spread apart.
2. Shave the genital area.
3. After shaving, scrub the area following approved hospital pre-operative scrub procedure.
4. Ensure the bladder is empty. A catheter is not required but may aid in identifying the urethra during the procedure.
5. An adhesive antimicrobial film may be used to protect the surgical site from post-operative stool contamination.

Dissection

1. Elevate the scrotum and make a low, midline perineal incision.
2. Carry the incision through Colles' fascia.
3. Open the bulbospongiosus muscle in the midline and reflect it laterally. Mobilize the corpus spongiosum proximally to the central tendon. Identify the central tendon by lifting the proximal bulb anteriorly. Mark or place a biodegradable tack suture at the origin of the central tendon on the proximal bulb. Sharply release the firm, fibrous portion of the central tendon.
4. The needle passer will be inserted into the obturator foramen at a point bordering the inferior pubic ramus defining the foramen which lies approximately one-third of the distance below the foraminal apex (about one finger breadth below the insertion of the adductor longus tendon). Palpate the inferior pubic ramus and feel for the bony landmarks to locate the proper position. A needle through the skin can be used to probe the bone to help confirm that the correct location for the needle passer entry point is found, but it is not required.
5. Make a small stab incision at the correct location over each obturator foramen. Confirm that both incisions lie in a straight line at the level shown below.


Figure 3. Needle Insertion Locations

6. Tissue retraction may need to be removed to avoid interfering with the needle passers. The patient is now ready for sling placement.

Placing the Sling

1. Open the AdVance XP Male Sling System tray package.
2. Remove the AdVance XP Male Sling System needle passer that is designated for the patient's left side and the male sling from the package.
3. Hold the needle passer handle at a 45° angle to the midline incision as shown below.


Figure 4. Needle Approach

4. Place a finger of the opposite hand into the perineal incision to receive the needle and protect the urethra. The finger should be at the apex of the triangle formed by the bulbar urethra medially and the ischiaopubic ramus laterally.
5. Insert and advance the needle tip through the obturator foramen along the lateral edge of the pubic ramus. The thumb of the opposite hand may be placed on the outside curve of the needle to help push it through.
6. Once two "pops" are felt, stop advancing the needle and drop the needle handle toward the midline to direct the needle path towards the fingertip placed in the anterior aspect of the triangle between the corpus spongiosum and ischiopubic ramus.
7. Rotate the needle until the tip meets your finger (placed at the anterior apex of the triangle). If the needle tip cannot be located, withdraw the needle passer and carefully advance again. Using the fingertip to guide the needle and protect the urethra, complete the rotation to bring the needle tip into the perineal incision.
8. Connect the sling to the needle by pressing the white connector onto the needle tip until it clicks into place.
9. Make sure that the blue printing on the sling sleeves is facing away from the corpus spongiosum and that the sling is not twisted. Rotate the needle back along the insertion pathway to pull the sling through the obturator foramen and out the stab incision. Pull only enough to bring the center of the sling to the patient's midline.
10. Cut the insertion sheath and mesh just above the tab on the white sheath liner.

Note: Leave enough sling material above the level of the skin so that the sling can be tensioned and the sheath can later be removed.

11. Using your opposite hand and the other needle passer, repeat the needle passage and sling insertion procedure on the patient's contralateral side to loosely position the sling with the blue marks on the sling sleeves facing outward, away from the urethra.

Attaching the Sling

1. Position the center of the sling flat against the corpus spongiosum with the edge of the proximal flap at the previously marked point where the central tendon was attached.
2. Use tack sutures (resorbable sutures recommended) to securely attach the center portion of the sling to the corpus spongiosum proximally and distally. The sutures should take an adequate bite of spongiosum and pass through the mesh at least two pores from the edge of the sling to ensure sufficient fixation.


Figure 5. Suture

Tensioning the Sling

If a catheter was used, it may be removed prior to tensioning but will not interfere if left in-situ. The sling is then tensioned by pulling equally on the distal end of each sling arm. It is important to exert sufficient force to reposition the bulbar urethra 2 cm - 4 cm.

Note: Do not over-tighten the sling. The goal is to reposition the urethra without causing obstruction or compression of the urethral lumen. Apply only enough tensioning force to the sling to achieve a 2 cm - 4 cm displacement of the bulbar urethra. Stop tensioning the sling upon achieving this level of displacement. Tensioning force will be dependent on patient anatomy.

1. Clamp a device such as a hemostat across each sling-arm, at the lateral stab incisions. Be sure that the complete width of the sling is captured within the clamp. The sling may be rolled around the clamp to improve the grip.
2. Pull firmly on both ends of the sling to observe approximately 2 cm - 4 cm proximal movement of the urethra. The patient's legs may be partially lowered to achieve complete repositioning.
3. Cystoscopy may be used to confirm coaptation after tensioning. A flexible cystoscope should be used if cystoscopy is performed. If coaptation is not seen, reposition the sling (usually further distal) by loosening it with an instrument placed between the sling and the urethra, removing the tack sutures, and re-tacking and re-tensioning the sling in the correct position. Loosening the sling may also be accomplished by gripping the sling arms with a clamp at a point deep to the wide center section of the mesh and pulling forward.

4. Remove the plastic sheaths from the sling by gripping both the white sheath liner and the sheath with a clamp. Be careful not to grip the mesh. Confirm the correct tension of the sling after the sheath has been removed. Give the sling ends a final pull if needed.
5. If desired, the ends of the sling may be tunneled subcutaneously back to the perineal incision.
6. Trim the sling at the level of the subcutaneous tissue. In the event that a tissue anchor becomes exposed at the sling exit, the sling arms may still be tunneled if desired or the sling may be trimmed beneath the exposed anchor.
7. Close the bulbospongiosus muscle then the rest of the perineal incision and the stab incisions in several layers after irrigation and hemostasis.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging that came into contact with bodily fluid should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

Post Procedure

- A catheter may be used at the discretion of the surgeon.
- Antibiotic prophylaxis should be given.
- The ability of the patient to empty his bladder should be confirmed.

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

IMPLANTABLE DEVICE PATIENT INFORMATION

Advise the patient that additional information may be available to them on the Boston Scientific website (www.bostonscientific.com/patientlabeling).

IMPLANT CARD INSTRUCTIONS

To complete the implant card that came with the AdVance XP product, peel off label from the product labeling and apply peel-off label onto the supplied patient implant card. Fill in the implant date, patient name, healthcare institution and/or physician information. Give the completed card to the patient.

PATIENT COUNSELING INFORMATION

Prior to the use of this device, brief the patient on any contraindications, warnings, precautions, adverse events, and post-procedure instructions outlined in this document that are of relevance to the patient. Explain that AdVance XP is a permanent implant and is MR Safe.

Patients should be counseled on abstaining from heavy lifting, strenuous exercise, and intercourse for a minimum of 6 weeks. Patients may return to other normal daily activities at their physician's discretion. A stool softener may be prescribed to prevent straining with bowel movement. A pain medication may be prescribed for short term pain management. If discomfort while urinating, unable to urinate, bleeding, or other problems occur, the patient should contact their physician.

Device Lifetime

The mesh and tissue anchors are a permanent implant and are intended to treat Stress Urinary Incontinence (SUI) over the patient's lifetime.

Although the outcome of the surgery is intended to be permanent, the functional life expectancy of the mesh may be impacted by unpredictable future factors such as: changes in pelvic anatomy, changes in body structure, new pelvic floor conditions such as urge incontinence, and/or surgeries that may occur after implantation to address potential mesh related complications.

Note: Complete removal of the mesh may not be possible and additional surgeries may not always fully correct the complications as the mesh becomes fixated by in-growth of surrounding tissues after several weeks post-implant.

The materials used in the sling are non-absorbable. Biocompatibility testing was completed to support material biocompatibility for AdVance XP's usage as a long term implant device.

Duration of Use

AdVance XP is a long term implant (>30 days); the needle passers and sheath assemblies are considered limited contact devices (≤24 hours).

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: AdVance.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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or its affiliates. All rights reserved.
©2022 Boston Scientific Corporation