**Boston
Scientific**

AdVance™ XP

Male Sling System

Directions for Use	2
Instrucciones de uso	11
Mode d'emploi	20
Gebrauchsanweisung	30
Istruzioni per l'uso	40
Gebruiksaanwijzing	50
Instruções de Utilização	60

51117738-01

2020-04

TABLE OF CONTENTS

WARNING	3
DEVICE DESCRIPTION.....	3
User Information	3
MR Safe	3
Contents.....	3
Figure 1. Needle Passers	4
Figure 2. AdVance™ XP Male Sling.....	4
INTENDED USE/INDICATIONS FOR USE	4
CONTRAINDICATIONS.....	4
WARNINGS.....	4
Procedural.....	4
Device Related	5
PRECAUTIONS.....	5
Procedural.....	5
Post Procedure.....	5
Device Related	6
ADVERSE EVENTS	6
HOW SUPPLIED.....	6
Handling and Storage.....	6
OPERATIONAL INSTRUCTIONS.....	6
PROCEDURE.....	6
Patient Preparation.....	7
Dissection.....	7
Figure 3. Needle Insertion Locations	7
Placing the Sling	7
Figure 4. Needle Approach	8
Attaching the Sling	8
Figure 5. Sutures.....	9
Tensioning the Sling.....	9
Immediate Post-Operative Care	10
WARRANTY	10

AdVance™ XP

Male Sling System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The AdVance XP Male Sling System consists of two sterile, single-use surgical instruments called “needle passers” and a mesh sling with attached connectors. One end of each needle passer is keyed to allow a secure attachment of the sling connector. The sling mesh is constructed of polypropylene monofilament mesh that is precut. The sling arms contain polypropylene tissue anchors molded directly to the sling mesh. Polypropylene tensioning fibers are integrated into the mesh weave. Plastic sheaths are placed over each arm of the mesh to aide in ease of placement. The sheaths contain liners to protect the plastic sheath from tearing upon removal. The connectors are attached to the ends of the needle passers during the procedure. The mesh and tissue anchors of the sling are intended to remain in the body as a permanent implant. The materials used in the sling are not absorbed or degraded by the action of tissue in-growth or tissue enzymes.

User Information

Precaution: Read the entire Directions for Use in detail before using the AdVance XP Male Sling System. The system should only be used by or under the supervision of physicians thoroughly trained in sling placement technique. A thorough understanding of the technical principles, clinical applications, and risk associated with this procedure is necessary before using this device.

MR Safe

The AdVance XP Male Sling is composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic and poses no known hazards resulting from exposure to any MR environment. Therefore, the AdVance XP Male Sling is considered MR Safe.

Contents

- (1) Male Sling
- (1) Needle Passers
- (1) Retractor Ring
- (1) Blunt Stay Hooks


Figure 1. Needle Passers


Figure 2. AdVance™ XP Male Sling

INTENDED USE/INDICATIONS FOR USE

The AdVance XP Male Sling System is intended for the treatment of male stress urinary incontinence (SUI) by the placement of a suburethral sling.

CONTRAINDICATIONS

Contraindications associated with the use of the AdVance XP Male Sling System include patients with:

- urinary tract infections;
- blood coagulation disorders;
- a compromised immune system or any other condition that would compromise healing;
- renal insufficiency and upper urinary tract relative obstruction.

WARNINGS

- Only physicians who have been trained as to its use should use the AdVance XP Male Sling System.
- The AdVance XP Male Sling System must not be re-sterilized or re-used.
- Inspect the packaging for visible damage prior to use. Any damage to the sterile barriers renders the device non sterile.
- Do not use product beyond the indicated expiration date.
- Only the innermost pouch may be introduced into the sterile field.
- Physicians should be familiar with surgical procedures and techniques involving non-absorbable meshes before using the AdVance XP Male Sling System device.
- Good surgical practice should be followed for management of contaminated or infected wounds.
- Urinary tract infection should be treated prior to implantation.

Procedural

- Do not insert the AdVance XP Male Sling System needles with an approach that deploys the needle from the perineum to the obturator.
- Take care to avoid vessel perforation. Observe patient for any signs of retropubic or periurethral bleeding.

Device Related

- Non-functional instruments should not be used and should be returned to Boston Scientific.
 - Do not use the AdVance™ XP Male Sling System if the package is opened or damaged, as sterility may be compromised.
-

PRECAUTIONS

- The implantation of this device should only be considered for adult patients whom the physician determines are appropriate surgical candidates.
- The risks and benefits of using the AdVance XP Male Sling System procedure should be carefully considered on patients with:
 - blood coagulation;
 - compromised immune systems or any other conditions that would compromise healing;
 - renal insufficiency and upper urinary tract obstruction.
- The possibility of urgency incontinence should be carefully considered before a sling implant is conducted. It is recommended that:
 - good bladder function (bladder capacity > 250 ml and post void residual urine < 50 ml) be demonstrated by candidates for a male sling;
 - the presence of bladder neck or urethral strictures be ruled out for male sling candidates;
 - a condition involving cystitis, urethritis or prostatitis be ruled out for male sling candidates;
 - detrusor instability of a neurological origin be ruled out for male sling candidates.
- A 6 month period of non-invasive treatment (e.g., behavior modification, bladder exercises, biofeedback, extra corporeal magnetic stimulation of the pelvic floor, or drug therapy) is recommended before a sling implant is considered for males with stress urinary incontinence.

Procedural

- Check for bleeding and pelvic organ dysfunction due to intraoperative vessel or nerve damage associated with anomalous location, during needle passage through the medial area of the obturator foramen membrane.
- It is important to verify the tension and placement of the sling prior to closure.
- Do not remove the plastic sheath until the sling is in its desired position.
- If during the procedure it becomes necessary to remove the sling, cut the sling near its midline and draw the sling arms out of the patient by pulling on the distal (connector) ends.
- Do not contact the sling with any staples, clips, or other instruments which may damage the mesh.

Post Procedure

- Acute inflammatory tissue reaction and transitory local irritation has been reported with the use of the absorbable suture.
- After several (2-3) weeks post-implant, the mesh becomes fixated by in-growth of surrounding tissues and cannot be entirely removed. However, the urethral repositioning may be reversed by cutting through the center of the mesh to release the tension.
- Patients should be counseled on abstaining from heavy lifting, strenuous exercise, and intercourse for a minimum of 6 weeks. Patients may return to other normal daily activities at the physician's discretion, often one to two weeks post procedure.
- If dysuria, bleeding, or other problems occur, the patient should be instructed to call the surgeon immediately.

Device Related

- Do not use any part of the AdVance™ XP Male Sling System beyond the indicated expiration date.
 - To maintain sterility, only the innermost package of the sling should be introduced into the sterile field.
 - Store the AdVance XP Male Sling System in a clean, dry, area away from direct sunlight and at room temperature.
-

ADVERSE EVENTS

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the implant could include extrusion, erosion through the urethra or other surrounding tissue, migration of the device from the desired location, fistula formation and inflammation. The occurrence of these responses may require removal of the entire sling.
- Like all foreign bodies, the sling may allow an existing infection to propagate.
- Incorrect sling tension may cause temporary or permanent lower urinary tract obstruction and retention of urine which may require surgical intervention.
- Unintended damage to muscle, tissue, nerves, or vessels.

Known risks of surgical procedures for the treatment of incontinence include:

- Allergic reaction
- Prolonged procedure
- Scar tissue
- Toxicity
- Hematoma
- Perforation
- Exposure to biohazardous material
- Unretrieved device fragment
- Pain
- Infection
- Erosion
- Device migration
- Bleeding
- Pelvic organ dysfunction
- Injury (urethra, corpus spongiosum, nerve) and malposition
- Complete failure of the procedure, resulting in incontinence and mild to moderate urinary incontinence due to incomplete support or overactive bladder.

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Do not use product beyond the indicated expiration date.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

The AdVance XP Male Sling System procedure can be carried out under spinal or general anesthesia. A small midline perineal incision is made followed by periurethral dissection. Two small stab incisions are also made above the obturator foramen for needle passer entry.

PROCEDURE

Verify that the package integrity has not been compromised during shipping. Inspect the device for damage.

Patient Preparation

1. Place the patient in a dorsal lithotomy position, with legs bent at approximately 90° and slightly spread apart.
2. Shave the genital area.
3. After shaving, scrub the area following approved hospital pre-operative scrub procedure.
4. Ensure the bladder is empty. A catheter is not required but may aid in identifying the urethra during the procedure.
5. An adhesive antimicrobial film may be used to protect the surgical site from post-operative stool contamination.

Dissection

1. Elevate the scrotum and make a low, midline perineal incision.
2. Carry the incision through Colles' fascia.
3. Open the bulbospongiosus muscle in the midline and reflect it laterally. Mobilize the corpus spongiosum proximally to the central tendon. Identify the central tendon by lifting the proximal bulb anteriorly.

Mark or place a biodegradable tack suture at the origin of the central tendon on the proximal bulb. Sharply release the firm, fibrous portion of the central tendon.

4. The needle passer will be inserted into the obturator foramen at a point bordering the inferior pubic ramus defining the foramen which lies approximately one-third of the distance below the foraminal apex (about one finger breadth below the insertion of the adductor longus tendon). Palpate the inferior pubic ramus and feel for the bony landmarks to locate the proper position. A spinal or percutaneous needle through the skin can be used to probe the bone to help confirm that the correct location for the needle passer entry point is found, but it is not required.
5. Make a small stab incision at the correct location over each obturator foramen. Confirm that both incisions lie in a straight line at the level shown below.


Figure 3. Needle Insertion Locations

6. Tissue retraction may need to be removed to avoid interfering with the needle passers. The patient is now ready for sling placement.

Placing the Sling

1. Remove the AdVance™ XP Male Sling System needle passer that is designated for the patient's left side from the package.
2. Hold the needle handle at a 45° angle to the midline incision as shown below.


Figure 4. Needle Approach

3. Place a finger of the opposite hand into the perineal incision to receive the needle and protect the urethra. The finger should be at the apex of the triangle formed by the bulbar urethra medially and the ischiaopubic ramus laterally.
4. Insert and advance the needle tip through the obturator foramen along the lateral edge of the pubic ramus. The thumb of the opposite hand may be placed on the outside curve of the needle to help push it through.
5. Once two “pops” are felt, stop advancing the needle and drop the needle handle toward the midline to direct the needle path towards the fingertip placed in the anterior aspect of the triangle between the corpus spongiosum and ischiopubic ramus.
6. Rotate the needle until the tip meets your finger (placed at the anterior apex of the triangle). If the needle tip cannot be located, withdraw the needle passer and carefully advance again. Using the fingertip to guide the needle and protect the urethra, complete the rotation to bring the needle tip into the perineal incision.
7. Connect the sling to the needle by pressing the white connector onto the needle tip until it clicks into place.
8. Make sure that the blue printing on the sling sleeves is facing away from the corpus spongiosum and that the sling is not twisted. Rotate the needle back along the insertion pathway to pull the sling through the obturator foramen and out the stab incision. Pull only enough to bring the center of the sling to the patient’s midline.
9. Cut the insertion sheath and mesh just above the tab on the white sheath liner.

Note: Leave enough sling material above the level of the skin so that the sling can be tensioned and the sheath can later be removed.

10. Using your opposite hand and the other needle, repeat the needle passage and sling insertion procedure on the patient’s contralateral side to loosely position the sling with the blue marks on the sling sleeves facing outward, away from the urethra.

Attaching the Sling

1. Position the center of the sling flat against the corpus spongiosum with the edge of the proximal flap at the previously marked point where the central tendon was attached.
2. Use tack sutures (resorbable sutures recommended) to securely attach the center portion of the sling to the corpus spongiosum proximally and distally. The sutures should take an adequate bite of spongiosum and pass through the mesh at least two pores from the edge of the sling to ensure sufficient fixation.


Figure 5. Sutures

Tensioning the Sling

If a catheter was used, it may be removed prior to tensioning but will not interfere if left in-situ.

The sling is then tensioned by pulling equally on the distal ends of each sling arm. It is important to exert sufficient force to reposition the bulbar urethra 2-4 cm.

Note: Do not over-tighten the sling. The goal is to reposition the urethra without causing obstruction or compression of the urethral lumen. Apply only enough tensioning force to the sling to achieve a 2-4 cm displacement of the bulbar urethra. Stop tensioning the sling upon achieving this level of displacement.

1. Clamp a device such as a hemostat across each sling-arm, at the lateral stab incisions. Be sure that the complete width of the sling is captured within the clamp. The sling may be rolled around the clamp to improve the grip.
2. Pull firmly on both ends of the sling to observe approximately 2–4 centimeters proximal movement of the urethra. The patient's legs may be partially lowered to achieve complete repositioning.
3. Cystoscopy may be used to confirm coaptation after tensioning. A flexible cystoscope should be used if cystoscopy is performed. If coaptation is not seen, reposition the sling (usually further distal) by loosening it with an instrument placed between the sling and the urethra, removing the tack sutures, and re-tacking and re-tensioning the sling in the correct position. Loosening the sling may also be accomplished by gripping the sling arms with a clamp at a point deep to the wide center section of the mesh and pulling forward.
4. Remove the plastic sheaths from the sling by gripping both the white sheath liner and the sheath with a clamp. Be careful not to grip the mesh. Confirm the correct tension of the sling after the sheath has been removed. Give the sling ends a final pull if needed.
5. If desired, the ends of the sling may be tunneled subcutaneously back to the perineal incision.
6. Trim the sling at the level of the subcutaneous tissue. In the event that a tissue anchor becomes exposed at the sling exit, the sling arms may still be tunneled if desired or the sling may be trimmed beneath the exposed anchor.
7. Close the bulbospongiosus muscle then the rest of the perineal incision and the stab incisions in several layers after irrigation and hemostasis.

Immediate Post-Operative Care

- A catheter may be used at the discretion of the surgeon.
- Antibiotic prophylaxis should be given.
- The ability of the patient to empty his bladder should be confirmed.
- **Patients should be counseled on abstaining from heavy lifting, strenuous exercise, and intercourse for a minimum of 6 weeks. Patients may return to other normal daily activities at their physician's discretion.**
 - **A stool softener may be prescribed to prevent straining with bowel movement**
 - **A pain medication may be prescribed for short term pain management**

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
**Boston Scientific Limited
**Boston Scientific (Australia) Pty Ltd
© 2020 Boston Scientific Corporation or its affiliates.