


CXR Rear Tip Extender (RTE) Kit

Instructions for Use

English	CXR Rear Tip Extender (RTE) Kit Instructions for Use.....	3
Deutsch	Proximale, nicht verkettungsfähige CXR Verlängerungskappen Gebrauchsanleitung.....	5
Español	Kit de extensores traseros CXR Instrucciones de uso.....	7
Français	Kit de prolongateurs arrière CXR Mode d'emploi.....	9
Italiano	Kit di estensori posteriori (RTE) CXR Istruzioni per l'uso.....	11
Português	Conjunto de Extensão de Ponta Traseira (RTE) CXR Instruções de Utilização.....	13

R_X ONLY

AMSTM

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CXR Rear Tip Extender (RTE) Kit Instructions for Use

Contents

This Rear Tip Extender (RTE) kit contains 2 sets of RTEs of the following sizes:

4 cm, 5 cm, 6 cm

Indication for Use

These RTEs are intended only for use with the AMS 700™ CXR inflatable penile prosthesis cylinders.

Instructions

1. Follow the instructions in the Operating Room Manual for implanting AMS 700 CXR Cylinders.
2. Select the proper length RTE for each cylinder that will best fit the patients anatomy. Sizing to the 1/2cm can be achieved, if desired, by attaching the 0.5cm, the 1.5cm, or for lengths longer than 1.5cm, by stacking one additional whole centimeter length RTE to the 1.5cm RTE.
3. Firmly press the RTE down on to the proximal end of the AMS 700 CXR cylinder. The RTE is correctly attached when the groove on the cylinder engages the corresponding locking ring on the RTE, resulting in a smooth outer profile. The same principle applies when attaching an additional RTE to the 1.5cm stacking RTE to achieve 1/2cm sizing lengths.
4. Slightly twist each RTE to verify that the locking ring is engaged.

Caution

1. Do not stack CXR Rear Tip Extenders other than the 1,5cm RTE. If other sizes of RTEs are stacked, then the locking mechanism will not engage and rear tip extenders may not stay connected to each other.
2. CXR RTEs are not compatible with CX or LGX cylinders.

How Supplied and Storage

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Store product in a clean, dry, dark area at room temperature.

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All rights reserved.
Boston Scientific (Australia) Pty Ltd
American Medical Systems, Inc.
U.S.A.
American Medical Systems Europe B.V.
©2016 Boston Scientific Corporation or its affiliates.