# Sterilization Instructions for AMS Tools

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## Sterilization Instructions for AMS Tools

### How to Sterilize AMS Tools prior to use

American Medical Systems does not sterilize AMS Tools.

Before using, refer to the product packaging label to determine if the tool should be removed from the packaging prior to sterilization.

Sterilize these stainless-steel instruments by steam autoclave using a gravity cycle at 250-254°F (121-123°C) for 20 minutes or a Prevac cycle at 273-279°F (134-137°C) for 3 1/2 minutes.

### Cleaning and Resterilization Instructions

Before next implant procedure, all components of AMS Tools must be thoroughly cleaned prior to initiation of the sterilization process. This cleaning process must include the following steps:

1. Rinse the device/components including lumens (if applicable) with warm running tap water for approximately two minutes or until all gross debris is removed.
2. If applicable, disassemble the device/components using the instructions provided. Use a cotton gauze pad moistened with warm tap water and wipe each piece of the device/components.
3. Rinse the device/components, including lumens (if applicable). Using a cleaning brush (CMS catalog #039-941 or equivalent), remove any visible debris from all areas of the device/components (including lumens).  
NOTE: It is the responsibility of the health care facility to ensure adequate cleaning of the brushes between each use.
4. Place the device/components into an ultrasonic cleaning bath containing ENZOL™ enzymatic detergent or equivalent cleaner diluted according to the manufacturer's instructions. Allow the device/components to remain in the ultrasonic bath for ten minutes.

5. Remove device/components and rinse for two minutes using warm running tap water.
6. Visually inspect the device/components to verify that all soil has been removed. If it has not, repeat the cleaning process.
7. Allow the device/components to dry. Double wrap the device/components in disposable sterilization pouches.

**WARNING:**

**Do not use an automated cleaning method, as it may be ineffective. The thorough manual cleaning process outlined above is required.**

Sterilize the device/components using a gravity steam autoclave cycle with the following parameters:

Time: 20 minutes

Temperature: 250-254°F (121-123°C)

Or a Prevac steam autoclave cycle with the following parameters:

Time: 3 1/2 minutes

Temperature: 273-279°F (134-137°C)

**Notices**

This document is written for professional medical audiences.

American Medical Systems, Inc. periodically updates product literature. If you have any questions regarding the currency of this information, please contact American Medical Systems.

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©2016 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are the property of the respective owners.
Pty Ltd
U.S.A.
**Europe B.V.**