# AMS Ambicor™

## Inflatable Penile Prosthesis

---

### Directions for Use

---

2

---

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---

11

---

### Mode d'emploi

---

20

---

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---

29

Oстаряла версия. Да не се използва.  
Zastaralá verzia. Nepoužívať.  
Forældet version. Ikke anvendes.  
Version überholte. Nicht verwenden.  
Aegunud versioon. Muy antigua. No utilizar.  
Палідійова версія. Не використовуйте.  
Outdated version. No longer in use.  
Versão obsoleta. Non utilizzate.  
Úrejt útgáta. Ne használja!  
Versione obsoleta. Neizmantoš.  
Novecojusi versija. Neizmantot.  
Pasenusi versió. Skal ikke brukes.  
Elavult verzió. Ne használja!  
Dit is een verouderde versie. Niet gebruiken.  
Utdatert versjon. Skal ikke brukes.  
Wersja przeterminowana. Nie używaj.  
Versão obsoleta. Não utilize.  
Versiune expirată. A nu se utiliza.  
Zastaraná verzia. Nepoužívať.  
Zastarela različica. Ne uporabite.  
Vanhentunut versio. Använd ej.  
Föråldrad version. Använd ej.  
Güncel olmayan sürüm. Kullanmayın.

**TABLE OF CONTENTS**

<b>WARNING .....</b>	<b>3</b>
<b>DEVICE DESCRIPTION .....</b>	<b>3</b>
User Information .....	3
Contents .....	3
Material of Animal Origin .....	3
<b>INTENDED USE.....</b>	<b>3</b>
<b>INDICATIONS FOR USE.....</b>	<b>3</b>
<b>CONTRAINdications.....</b>	<b>4</b>
<b>WARNINGS.....</b>	<b>4</b>
<b>PRECAUTIONS.....</b>	<b>4</b>
<b>PATIENT COUNSELING INFORMATION.....</b>	<b>4</b>
<b>SILICONE INFORMATION.....</b>	<b>4</b>
<b>MAGNETIC RESONANCE IMAGING (MRI).....</b>	<b>4</b>
MRI-Related Heating.....	4
Artifact Information .....	4
<b>ADVERSE EVENTS .....</b>	<b>5</b>
<b>CLINICAL STUDIES .....</b>	<b>6</b>
<b>HOW SUPPLIED.....</b>	<b>6</b>
Handling and Storage .....	6
<b>OPERATIONAL INSTRUCTIONS.....</b>	<b>6</b>
Pre-procedural Instructions .....	6
Procedural Instructions .....	6
Post-procedure Information .....	9
<b>PRODUCT REPLACEMENT INFORMATION.....</b>	<b>10</b>
<b>WARRANTY .....</b>	<b>10</b>

# AMS Ambicor<sup>TM</sup>

## Inflatable Penile Prosthesis

### Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

### WARNING

Contents supplied STERILE using steam or an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

### DEVICE DESCRIPTION

The AMS Ambicor Inflatable Penile Prosthesis is a closed, fluid-filled system consisting of a pair of cylinders, which are implanted in the corpora cavernosa, and a pump which is implanted in the scrotum (Figure 1). Rear Tip Extenders (RTEs) are provided for additional length, if needed. All components are preconnected by kink-resistant tubing and prefilled with normal saline. The cylinders are inflated as fluid is pumped from the reservoirs, which are located in the proximal ends of the cylinders, into the main cylinder body, creating an erection. They are deflated as fluid is transferred back to the reservoirs, making the penis flaccid once again.


Figure 1. AMS Ambicor Inflatable Penile Prosthesis

The AMS Ambicor Penile Prosthesis cylinder specifications are shown in Table 1.

Table 1. AMS Ambicor Cylinder Specifications

Cylinder Diameter	Cylinder Length	Kink-Resistant Tubing Length
12.5 mm	14 cm	9 cm
	16 cm	
	18 cm	
14 mm	16 cm	9 cm
	18 cm	
	20 cm	
15.5 mm	18 cm	9 cm
	20 cm	
	22 cm	

### Device Function

#### Inflation

The inflation pump is placed in the scrotum. To inflate the cylinders, the patient stabilizes the pump with one hand and uses the other hand to squeeze and release the pump bulb several times to make the cylinders stiff (Figure 2). When the cylinders are fully inflated, the pump bulb will be hard and can no longer be compressed.


Figure 2. Inflating the AMS Ambicor Cylinders

#### Deflation

To deflate the cylinders, the patient places his thumb under the shaft of the penis to act as a fulcrum and places his fingers on top of the shaft of the penis. Using the fingers of the same or opposite hand, the patient bends his penis down over his thumb (fulcrum) towards his scrotum at an approximate 65 degree angle, making sure that both cylinders are bent (Figure 3). The cylinders are held in this position for approximately 10-12 seconds and then released. This opens the valves allowing the fluid to flow back into the reservoirs and pump. The patient may also deflate the cylinders by bending them upwards; this can be done by placing his thumb on top of the shaft of his penis and following the same instructions as above.


Figure 3. Deflating the AMS Ambicor Cylinders

### User Information

The implantation of this device is to be performed only by physicians who are knowledgeable and have received appropriate training regarding the use of inflatable penile prostheses. This manual is not intended to be a complete reference.

### Contents

- (2) Cylinders
- (1) Pump
- (2) 0.5 cm RTE
- (2) 1 cm Stackable RTE
- (2) 2 cm Stackable RTE
- (2) 3 cm Stackable RTE
- (2) Keith Needles

### Material of Animal Origin

An internal component of the Ambicor Penile Prosthesis cylinders contains a non-viable tallow derivative. The component does not contact blood or tissue.

### INTENDED USE

The AMS Ambicor Penile Prosthesis is designed to provide the patient with control over the erect and flaccid states of his penis.

### INDICATIONS FOR USE

The AMS Ambicor Penile Prosthesis is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

## CONTRAINDICATIONS

The implantation of this device is contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery.

## WARNINGS

- The implantation of this device will make latent natural or spontaneous erections as well as other interventional treatment options, impossible.
- Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.
- Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.
- Implantation of a penile prosthesis may result in penile curvature or scarring.
- This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
- Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.

## PRECAUTIONS

### Surgery-Related

- Spontaneous unintended inflation or deflation of the cylinders may occur and may result in unintended partial or full erections.
- Migration of the device components can occur if the cylinders are improperly sized, or if the pump is not positioned properly.
- Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or may make it impossible.
- Improper measurement technique, positioning, or sizing may reduce cylinder life.
- Unsuccessful outcomes have been reported due to improper surgical technique, anatomical misplacement of components, improper sizing of components, or tubing kinks.

### Device-Related

- Some of the materials used in the construction of this device have been shown to cause minor irritation when implanted in animals. Therefore, implantation of this device may cause minor irritation or discomfort in some patients.
- Do not use product that has damaged or open packaging, as sterility may be compromised.

### Patient-Related

- Adequate patient manual dexterity and strength are required for proper device inflation and deflation.
- Mental or psychological conditions, such as senile dementia, may inhibit the patient's successful operation of the prosthesis.
- Trauma to the pelvic area, such as impact injuries associated with sports (e.g., bicycle riding), can result in damage of the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction, including replacement of the device.
- The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.
- Use of injection therapy concurrently with the penile prosthesis can damage the prosthesis. Patients should not use injection therapy after they receive their implant.

## PATIENT COUNSELING INFORMATION

Patients should be counseled in order to have a realistic expectation of the physical (e.g., length and girth), psychological, and functional outcome of the implantation. The risks, benefits, and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option. A discussion concerning the patient's pump placement preference may be appropriate.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process. Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal.

Implantation of a penile prosthesis may result in penile curvature or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations. Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring, scrotal deformity, pump bulge in the scrotum, lack of concealability, and other possible adverse events. Patients should also be aware that penile prostheses are not considered to be lifetime implants.

Improper implantation of a penile prosthesis may not provide rigidity to the glans, which may result in a floppy glans and in a lack of rigidity of the corpus spongiosum. Penile flaccidity may be less than it was before implantation.

Patients who undergo revision surgery may notice a change in the character of their erection compared to their previous implant, which may include differences in sensation, length, girth, rigidity, and/or flaccidity. It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (See Silicone Information).

## SILICONE INFORMATION

This device is composed of a number of materials, including solid silicone elastomers. Silicone gel is not a component in the materials of this device. Silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate "allergic-like" symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer.

There are reports of malignant tumor formation in laboratory animals only, not humans, associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS Ambicor™ Device. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

## MAGNETIC RESONANCE IMAGING (MRI)

### Important Safety Information

Nonclinical testing demonstrated that the AMS Ambicor Penile Prosthesis is MR Conditional. The device can be safely scanned under the conditions shown in Table 2.

**Table 2. Scanning Conditions**

Static Magnetic Field	1.5 Tesla <sup>a</sup>	3.0 Tesla <sup>b</sup>
Spatial Gradient Field	450 Gauss/cm or less	720 Gauss/cm or less
Maximum whole body averaged Specific Absorption Rate (SAR)	1.5 W/kg for 15 minutes of scanning as assessed by calorimetry	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry

### MRI-Related Heating

Nonclinical testing demonstrated that the AMS Ambicor Penile Prosthesis produced the temperature rises shown in Table 3 during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the patient.

**Table 3. Temperature Rises**

Static Magnetic Field	1.5 Tesla <sup>a</sup>	3.0 Tesla <sup>b</sup>
Highest Temperature Change	≤ + 0.4° C	≤ + 1.9° C

(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI)

(b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5

### Artifact Information

Nonclinical testing demonstrated that the AMS Ambicor Penile Prosthesis may compromise the MR image quality if the area of interest is relatively close to the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a "moderate" localized signal void in the size and shape of the implant (Table 4). Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

**Table 4. Signal Voids**

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	6,244 mm <sup>2</sup>	1,589 mm <sup>2</sup>	10,295 mm <sup>2</sup>	2,779 mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

## ADVERSE EVENTS

A clinical trial was conducted to determine the safety and effectiveness of the AMS Ambicor™ Inflatable Penile Prosthesis. A total of 179 patients were implanted with follow-up out to 18 months or longer for 141 patients. There were four patient deaths during the course of the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

The Adverse Device Events (ADE), detailed in Tables 5 and 6, were noted during the duration of this clinical trial for all enrolled patients.

**Table 5: AMS Ambicor Penile Prosthesis Clinical Trial Device Related Adverse Events**

ADE	# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Urogenital Pain (Typically Associated with Healing Process)	38 (21.2%)	41	83 (0 – 548)
Patient Dissatisfaction (With Length, Ability to Use and Nonspecific Reasons)	22 (12.3%)	25	269 (0 – 1105)
Mechanical Malfunction	11 (6.1%)	11	394 (0 – 1508)
Penile Sensation Change	11 (6.1%)	11	266 (10 – 573)
Urogenital Hematoma	9 (5.0%)	9	9 (2 – 34)
Urogenital Edema	8 (4.5%)	9	9 (0 – 29)
Infection/Device Infection	6 (3.4%)	8	207 (32 – 606)
Auto-inflation	4 (2.2%)	4	110 (27 – 224)
Penile Curvature	4 (2.2%)	4	323 (0 – 999)
Abnormal Ejaculation	3 (1.7%)	3	256 (40 – 487)
Application Site Reaction (Wound Separation, Delay in Cutaneous Closure)	3 (1.7%)	3	26 (0 – 46)
Phimosis/Paraphimosis	3 (1.7%)	3	21 (3 – 38)
Post-Operative Bleeding	3 (1.7%)	3	4 (2 – 6)
Seroma	3 (1.7%)	3	242 (34 – 596)
Urinary retention	3 (1.7%)	3	121 (0 – 361)
Urogenital Ecchymosis	3 (1.7%)	3	123 (6 – 351)
Adhesion	2 (1.1%)	2	42 (34 – 50)
Auto-deflation	2 (1.1%)	2	486 (212 – 760)
Nocturia	2 (1.1%)	2	245 (105 – 385)
Penile Rash	2 (1.1%)	2	272 (43 – 501)
Renal Failure	2 (1.1%)	2	476 (396 – 555)
Urethral Laceration	2 (1.1%)	2	1 (0 – 2)
Urinary Tract Infection	2 (1.1%)	2	283 (4 – 561)
Urination Impaired	2 (1.1%)	2	225 (64 – 385)
Urogenital Erythema	2 (1.1%)	2	193 (13 – 373)
Other Urogenital	23 (12.8%)	24	97 (0 – 727)
Other Non-urogenital	10 (5.6%)	14	178 (0 – 556)

\* Percentages based on total number of patients implanted (n=179).

The following "other" urogenital adverse events each occurred in one patient: abnormal sexual function, allergic reaction, BPH, burning urination, crossover from left to right corpora, device migration, dribbling, dysuria, epididymitis, extrusion, incontinence, induration in penile area, lower urinary tract symptoms, malposition, mental excoriation, micro-hematuria, penile scar tissue, prerenal azotemia, scab on corpora, soft glans, tunica perforation, urosepsis, wound dehiscence and wound infection.

**Table 6: Resolution of Device Related Adverse Events**

ADE	# Patients (%)	# Events	Surgical	Medication	Other**	No Medical Intervention	Continuing
Urogenital Pain	38 (21.2 %)	41	3	9	1	29	8
Patient Dissatisfaction	22 (12.3 %)	25	3	3	3	16	14
Mechanical Malfunction	11 (6.1 %)	11	7	0	1	3	1
Penile Sensation Change	11 (6.1 %)	11	0	2	0	10	4
Urogenital Hematoma	9 (5.0 %)	9	0	1	3	6	0
Urogenital Edema	8 (4.5 %)	9	0	2	2	5	0
Infection/Device Infection	8 (4.5 %)	8	6	1	0	1	1
Auto-inflation	4 (2.2 %)	4	0	0	0	4	2
Penile Curvature	4 (2.2 %)	4	2	0	0	2	2
Abnormal Ejaculation	3 (1.7 %)	3	0	1	0	2	1
Application Site Reaction	3 (1.7 %)	3	0	1	0	2	0
Phimosis/Paraphimosis	3 (1.7 %)	3	1	2	0	0	0
Post-Operative Bleeding	3 (1.7 %)	3	0	0	0	3	0
Seroma	3 (1.7 %)	3	0	0	1	2	0
Urinary retention	3 (1.7 %)	3	1	1	2	0	1
Urogenital Ecchymosis	3 (1.7 %)	3	0	0	0	3	0
Adhesion	2 (1.1 %)	2	1	0	0	1	0
Auto-deflation	2 (1.1 %)	2	0	0	1	1	0
Nocturia	2 (1.1 %)	2	0	1	0	1	2
Penile Rash	2 (1.1 %)	2	0	1	1	0	1
Renal Failure	2 (1.1 %)	2	0	0	1	1	1
Urethral Laceration	2 (1.1 %)	2	0	0	1	1	0
Urinary Tract Infection	2 (1.1 %)	2	0	2	0	0	0
Urination Impaired	2 (1.1 %)	2	0	1	1	0	2
Urogenital Erythema	2 (1.1 %)	2	0	2	0	0	0

\*\* Other methods of resolution consist of standard medical treatments.

A total of 14 patients underwent revision surgeries in the study. Information on device revisions is described in the Clinical Studies section.

The following events occurred during this clinical study but were classified by the Investigators as not related to the device or the implantation surgery: urethral stricture (9 cases), prostate cancer (4 cases), BPH (3 cases), incontinence (3 cases), urinary retention (3 cases), urinary tract infection (3 cases), abnormal prostate exam/PSA (2 cases), abnormal sexual function (2 cases), balanitis (2 cases), hematuria (2 cases), pain (2 cases), renal calculus (2 cases), abnormal ejaculation (1 case), bladder instability (1 case), bladder spasms (1 case), blood in urine (1 case), dysuria (1 case), edema (1 case), epididymal cyst (1 case), eroded artificial sphincter (1 case), erythema (1 case), infection (1 case), inguinal hernia (1 case), nocturia (1 case), prostate nodule (1 case), prostatitis (1 case), pyuria (1 case), renal colic (1 case), renal insufficiency (1 case), testicular cyst (1 case), testicular hydrocele (1 case), trichomonas (1 case), urethral mucosal atrophy (1 case), urination impaired (1 case), urgency (1 case), and urosepsis (1 case).

The following adverse events may be associated with this product: abscess, bleeding, erosion, excessive fibrous capsular growth, exposure to biohazardous material, genital changes, granulomas, hemorrhage, immune-related connective tissue disorders, improper size, inguinal hernia, ischemia, necrosis, pain (which may be prolonged or severe), perforation or injury of bladder, corpus cavernosum, nerve, tunica, or urethra, prolonged procedure, Super Sonic Transporter (SST) (as a result of improper sizing), ulceration, untrieved device fragment, vascular compromise, and ventral chordee.

## CLINICAL STUDIES

A clinical trial was undertaken to demonstrate that the AMS Ambicor™ Inflatable Penile Prosthesis provides an erection that is suitable for intercourse and has acceptable rates of surgical revision and of significant clinical events associated with the implantation and use of these devices. This trial was also designed to demonstrate the implantation of these devices does not negatively impact the sexual satisfaction, psychological well-being, self-esteem, or quality of life of patients who receive these devices. It was a prospective, multi-center cohort trial in which the patients served as their own control.

One hundred and seventy-nine (179) male patients, over 21 years of age (range 31 – 81 years), were enrolled and implanted in this study. The majority of study participants were White (66.5%), with African Americans making up 27.4% of the study population followed by Hispanics (3.4%), Indians (1.1%), Asians (0.6%), Native Americans (0.6%), and participants from the Middle East (0.6%). All patients with diagnosed organic erectile dysfunction were eligible for enrollment, if they did not present with a history of allergy/sensitivity to silicone, preexisting autoimmune or connective tissue diseases or active urogenital infection.

All safety-related data, diagnoses and health status evaluations were captured on detailed case report forms. The Investigators' professional evaluation of the erections provided by the device after implantation and its suitability for intercourse was the primary efficacy endpoint. The number of surgical revisions performed and reported by the Investigators was the primary safety endpoint. Patient self-evaluations on four validated outcome instruments were the secondary efficacy endpoints (concerning quality of life, self-esteem, psychological well-being, and sexual satisfaction and functioning).

This clinical trial provided the following results through the 18-month evaluation for the first 141 patients to reach this post-surgical follow-up. Some patients were studied for longer than 18 months.

### Physician Assessment of Device Function

The function of 141 implanted devices was evaluated on or after the 18-month follow-up exam. Of these 141 devices, 140 could be inflated (99.3%). The device that could not be inflated was surgically replaced. All of the 140 devices that were successfully inflated were determined to provide an erection suitable for intercourse. This assessment of device function does not include the following information regarding device malfunctions: Of the 140 devices found to be functioning properly on or after 18 months, 2 required revision surgery prior to this evaluation to correct a mechanical malfunction. These and other cases of revision surgery are discussed in the next section.

### Surgical Revisions

The incidence of revisions was evaluated in the 141 patients with follow-up out to at least 18 months, and an additional 6 patients who were explanted prior to 18 months and subsequently discontinued from the study. A revision is considered any urogenital surgical intervention that is related to the function, placement, or site reaction to the implanted device. Of these 147 patients, 14 (9.5%) experienced a total of 15 revision surgeries, and 133 (90.5%) have not been revised.

The average time to the first revision surgery was 12 months (ranging from 0 to 50 months). Of the fifteen (15) revision surgeries, there were six (6) revisions due to "Mechanical Malfunction", five (5) due to "Infection", one (1) due to "Migration/Extrusion" (pump migration and cylinder aneurysm/bulging/extrusion), one (1) due to "Malposition", one (1) due to "Phimosis", and one (1) due to "Patient Dissatisfaction".

### Patient Evaluation of Quality of Life, Self-esteem, Psychological and Sexual Well-Being

No decline in overall patient quality of life, self-esteem, psychological well-being, sexual functioning, or sexual satisfaction was demonstrated during this clinical trial. The psychological well-being of these patients, as judged by the Brief Symptom Inventory, was determined to be equivalent to the pre-implant state. The self-esteem of the patients, as evaluated by the Rosenberg Self-Esteem Scale, was also equivalent to the pre-implant state. The sexual functioning and sexual satisfaction did improve from the pre-implant state, as measured by the Sexual History Form.

### HOW SUPPLIED

The AMS Ambicor Penile Prosthesis and RTEs are provided sterile in a fluid-filled pouch. The prosthesis is preconnected and prefilled. The RTEs are packaged with the prosthesis.

- Do not use if package is open or damaged.
- Do not use if labeling is incomplete or illegible.
- Do not use the device if past the "Use By" date.

### Handling and Storage

Store in a cool, dry, dark place at room temperature.

## Supplemental Surgical Tools

The following Boston Scientific surgical tools used to help facilitate implantation of the AMS Ambicor Penile Prosthesis may be ordered separately.

The following tools are provided non-sterile and must be sterilized before use (for reprocessing information, refer to the instructions manual provided with the tools):

- Sizer

### Contents:

- (1) 9.5 mm Sizer

- (1) 12 mm Sizer

- Furlow Insertion Tool

- Closing Tool

The following are provided sterile and are intended for single use only:

- Proximal Tool (dual use for cylinder insertion and closing)
- Disposable Dilators
- SKW Deep Scrotal Retraction System

### Required Materials

Instruments normally required for a urological surgical procedure are recommended. In addition to the AMS Ambicor Penile Prosthesis, RTEs, and Keith Needles, the following sterile setup is needed:

- Corporal length measuring tool
- Furlow Insertion Tool for passing traction sutures through the glans
- Disposable Dilators or a set of corporal dilators
- Sterile surgical stand (Mayo stand) or a stainless steel tray
- Basin of sterile normal saline
- Antibiotic solution for irrigation

## OPERATIONAL INSTRUCTIONS

### Pre-procedural Instructions

#### Prepare the Patient

Before surgery, take adequate steps to limit the risk of postoperative infection.

Once the patient is in the operating room, shave the abdominal and genital area. Following the shave, scrub the area using the approved hospital preoperative scrub procedure.

Establish the sterile field, drape, position, and prepare the patient so that a penoscrotal incision can be made.

### Procedural Instructions

The following description is an overview of the penoscrotal surgical approach. This approach leaves the incision well-hidden and provides convenient access to the corpora cavernosa. Throughout the procedure, the surgical site may be flushed with copious amounts of broad-spectrum antibiotics.

#### Dilate and Size

1. Place a Foley catheter to facilitate identification of the urethra.
2. Some physicians use the SKW Deep Scrotal Retraction System to place the penis on "stretch" and to provide exposure of the corpora.
3. Make a 2 cm to 3 cm incision through the median raphe of the scrotum at the penoscrotal junction. Some physicians prefer a high scrotal incision for better proximal corporal access. Laterally retract the corpus spongiosum to avoid damaging the urethra.
4. Dissect through to Buck's fascia to expose the tunica albuginea. Place stay sutures to use as a reference point when measuring the corpora. Make an incision into one of the corpora cavernosa.
5. Dilate the proximal corpus (crus) and the distal corpus to create a space for inserting a penile cylinder. Dilate the corpus cavernosum to 13 mm-16 mm.

---

**Note:** The available AMS Ambicor Penile Prosthesis cylinder diameters are 12.5 mm, 14 mm, and 15.5 mm.

---

After dilating one corpus cavernosum, incise and dilate the adjacent corpus cavernosum following the same procedure. If the corpora cavernosa are unevenly sized or the measuring tool cannot reach mid glans, consider dissecting through the fibrotic tissue as needed.

6. To select the cylinders and rear tip extenders (if needed) to fit the patient's anatomy, measure each corpus length proximally and distally using a corporal length measuring device (See the Sizing section for choosing the cylinder size). As a general rule, the corporotomy is best placed when two-thirds of the total corporal measurement is distal to the incision and one-third is proximal. This facilitates the placement of the cylinders, and may avoid the need to extend the corporotomy during the procedure.

## Sizing Method

Select the appropriate size cylinders and if needed, RTEs.

This method allows the tubing to exit directly from the corporotomy. Follow the formula described below to select the appropriate cylinder length and number of rear tip extenders. If necessary, extend the length of the corporotomy.

### 1. Calculate the Total Corporal Length (Distal + Proximal)

Example

Distal Corporal Length	11 cm
Proximal Corporal Length	+8 cm
Total Corporal Length	19 cm

### 2. Subtract 2 cm from the Total Corporal Length to obtain an Adjusted Measurement.

Example

Total Corporal Length	19 cm
	-2 cm
Adjusted Measurement	17 cm

### 3. Select the closest cylinder size that is shorter than or equal to the Adjusted Measurement.

Example

Adjusted Measurement	17 cm
Selected Cylinder Length	16 cm

### 4. Subtract the Selected Cylinder Length from the Total Corporal Length to determine the length of rear tip extenders required to fit the patient.

Example

Total Corporal Length	19 cm
Selected Cylinder Length	-16 cm
Rear Tip Extender Length	3 cm

## Unpack the Device

**NOTE:** Do not open the prosthesis package until the surgeon dilates and sizes both corpora cavernosa intraoperatively. The adhesive label at one end of the box and the small, peelable labels on the foil pouches show the part and serial/lot numbers and the size of the components. Record this information on and affix this label to the Patient Information Form.

Keep the AMS Ambicor™ Device in the box until it is in the operating room.

## Remove the AMS Ambicor Penile Prosthesis from the Packaging

1. Remove the pouch from the box. Remove the inner fluid-filled foil pouch from the outer pouch.

**CAUTION:** To avoid the possibility of cutting the AMS Ambicor Penile Prosthesis or RTEs, do not open the package with scissors.

2. Place the sterile, fluid-filled foil pouch onto a sterile, lint free Mayo stand.

**CAUTION:** Cloth towels placed on the Mayo stand may transfer lint to the AMS Ambicor Components.

3. Hold the foil pouch over a basin containing sterile normal saline. Carefully tear the foil pouch at the side notch (Figure 4). Empty the AMS Ambicor Device, its RTEs, and the fluid into the basin. Add sterile normal saline as needed to cover the components in the basin (Figure 5).


Figure 4. Tear the pouch at the notch

Figure 5. Empty the AMS Ambicor Device and RTEs into the basin

**NOTE:** Do NOT leave the prosthesis exposed to air. Submerge the prosthesis in the fluid from the foil pouch or in sterile normal saline immediately after removing it from the foil pouch. Keep the prosthesis submerged until implantation to prevent air from entering the device.

4. Remove the tips of the AMS Ambicor Penile Prosthesis cylinders from the protective end cushion.

5. Discard the end cushion.

6. Leave the device and RTEs in the basin until implantation.

## Apply the RTEs if Needed

The RTEs are supplied in two stacks of four each in the following lengths: (2) 0.5 cm, (2) 1 cm, (2) 2 cm, and (2) 3 cm RTEs (Figure 6).


Figure 6. RTEs: 0.5 cm, 1 cm, 2 cm, and 3 cm

1. Gently pull the RTEs apart (Figure 7).

The 0.5 cm RTEs are not stackable and if used, must be the last RTE placed on the cylinders. The 1 cm, 2 cm, and 3 cm are stackable and can be combined in lengths appropriate for the patient's anatomy. Stack the RTEs from longest to shortest on the cylinder.


Figure 7. Gently pull the RTEs apart

2. Attach the RTEs to the proximal end of the cylinders using the same size RTEs on each cylinder (Figure 8).


Figure 8. Attach the RTEs

## Implant the Cylinders

Use the Furlow Insertion Tool to introduce the cylinders into the corpora cavernosa, as follows:

1. When inserting the cylinder distally, use the Furlow Insertion Tool to put the penis on a mild stretch. The tool should be palpable in the glans.
2. To help prevent cross over through the intercarvernosus septum when placing the first cylinder, a dilator or other surgical tool may be placed in the contralateral corpus cavernosum.
3. Withdraw the obturator of the Furlow Insertion Tool into the retracted, or locked position.
4. Thread both ends of the traction suture (3-0 braided, nonabsorbable, PTFE coated polyester) at the distal tip of the cylinder through a Keith Needle (Figure 9), and place the blunt end of the needle in the insertion tool (Figure 10). Place the suture in the tool slot, completely retract the suture into the slot and fully draw the needle into the barrel of the tool.


Figure 9. Thread the traction suture through the Keith Needle


Figure 10. Place the Keith Needle into the Furlow Insertion Tool

5. Hold the four strands of suture against the tool and insert the tool into the distal portion of the corporal body until the front tip is at mid glans.
  - a. During insertion, be sure that the Furlow Insertion Tool is in the ipsilateral corpora at the distal penis. Because the intercavernous septum may be inconsistent distally, it is easy to cross over to the contralateral side.
  - b. If you think you may have crossed over, remove the Furlow Insertion Tool, place a dilator into the contralateral side and reposition the Furlow Insertion Tool on the ipsilateral side. No repair is necessary.

6. Pass the needle through the glans by fully inserting the obturator into the barrel (Figure 11). Grasp the needle with a needle holder and pull it completely through the glans.


Figure 11. Pass the Keith Needle through the glans using the Furlow Insertion Tool

7. Detach the needle from the suture, and remove the needle from the area to prevent any accidental cylinder punctures. Attach a hemostat to the remaining suture to prevent retraction through the glans.
8. Repeat the procedure using the remaining cylinder.
9. Inflate the cylinders fully, then insert the proximal end of the first cylinder into the crus (Figure 12). The Proximal Tool may be used to place the proximal cylinder end into the corpora.


Figure 12. Insert the proximal end of the cylinder into the crus.

10. Deflate the cylinder by bending it to approximately 65 degrees and holding it for 10-12 seconds (Figure 13).



Figure 13. Deflate the cylinder before placing the distal portion in the corpus cavernosum

11. Position the distal end of the cylinder by pulling on the traction sutures (Figure 14) to guide the cylinder until the front tip is placed at mid glans.


**Figure 14. Use the traction suture to position the distal cylinder end at mid glans**

12. Repeat Steps 9 through 11 to place the other cylinder.

13. Assess each cylinder length for satisfactory fit within the corpora cavernosa by ensuring that the distal tip is snugly at mid glans, the cylinder lies within the corporotomy, and the proximal end is firmly against the crus.

**NOTE:** The corporal incision may need to be extended to ensure the input tubing exits directly from the corporotomy.

14. Once both cylinders are implanted, close the tunica albuginea. Place the Closing Tool or the closing end of the Proximal Tool (or other suitable instruments) over the cylinder to protect it from inadvertent needle injury, so that the corporal body can be closed (Figure 15). This should be done with meticulous attention to hemostasis.

**Figure 15. Place an instrument over the cylinder to protect it while closing the corporotomy**


#### Implant the Pump

1. Use blunt dissection to form a pocket in the most lateral and dependent portion of the scrotum (Figure 16). The pump may also be placed in a dependent position between the two testicles.


**Figure 16. Use blunt dissection to form a pocket**

2. Insert the pump into the scrotal pocket (Figures 17 and 18). The tubing between the pump and cylinders should not be palpable to the patient. A suture may be placed to secure the pump in the scrotum, if desired.

**Figure 17. Insert the pump in the scrotal pocket**

**Figure 18. Final pump placement**


3. After both the cylinders and pump are implanted, inflate the device to check the quality of the erection and deflate to evaluate flaccidity. The penis should lie close to the body when deflated. There may be some swelling that precludes a good flaccid result.

4. Close the incision.

**Note:** Some physicians close the dartos in two layers with a running 2-0 absorbable suture.

5. Remove the traction sutures from the cylinders by cutting one end of the traction suture approximately 2 cm from the glans. Pull it out slowly to minimize trauma to the glans and the distal tip of the cylinder. Repeat with the other traction suture.

6. Apply a wound dressing (Figure 19) and leave the cylinders partially inflated. Support the penis on the lower abdomen. A drain may be placed for 12 to 24 hours, if desired.

**Figure 19. Apply a wound dressing and support the penis on the lower abdomen (example)**

#### Post-procedure Information

##### Immediately Postoperative

After the surgery, some physicians partially inflate the cylinders for the first 24 hours. This will aid hemostasis. A closed system drain may be placed in the abdomen to drain excess fluid from the incision site.

After 24 hours, remove the dressing and completely deflate the cylinders. Support the penis on the abdomen for four to six weeks to obtain a straight erection.

##### After the Patient is Released from the Hospital

After the patient has returned home and the swelling from the surgery has subsided, the physician may ask the patient to pull down on the pump located in the scrotum to properly position it. Positioning the pump makes it easier for the patient to locate the pump.

The frequency of positioning the pump is the physician's decision. Some physicians have their patients position the pump several times daily.

To position the pump in the scrotum, the patient should be told to:

1. Locate the pump in the scrotum.
2. Grasp the pump firmly and carefully pull the pump down in the scrotum. The patient should gently pull the pump into a position close to the outer scrotal wall.

After three to six weeks, the physician may instruct the patient to begin cycling the device for the first time. To cycle the device, the patient inflates and deflates the prosthesis several times. The cylinders must be fully inflated before they can be deflated.

It may be painful for the patient the first few times that he inflates and deflates the device. However, after the postoperative healing period, the pain should subside. Instruct the patient to inflate and deflate the prosthesis several times daily.

Four to six weeks postoperatively determine if the patient is ready to use the device:

1. Check the incision site to be sure that it has healed properly. There should be no redness, swelling, or drainage. Any of these findings may indicate that an infection is present and it should be treated promptly with antibiotics.
2. Ask the patient about pain when cycling the device and observe the patient inflating and deflating the device.

After determining that the patient knows how to operate the device and that the device is functioning correctly, inform the patient that it is possible to have gentle intercourse. Advise the patient not to inflate the device for long periods of time when not in use.

If the patient is familiar with injection therapies for erectile dysfunction, remind the patient that such therapies can cause damage to the penile prosthesis and should not be used.

#### **Evaluating Long-Term Function and Placement**

After the postoperative healing period, continue to have contact with the patient at least on an annual basis to evaluate the function of the device and to check for signs of infection or erosion. Ask the patient how the device is functioning and if he has noticed any changes in the function, for example, the cylinders losing rigidity. If the patient is having mechanical difficulty with the device, or if there is infection or erosion present, revision surgery may be necessary.

#### **PRODUCT REPLACEMENT INFORMATION**

For more information about the replacement of a device, component or accessory, please contact a Boston Scientific representative or local distributor. Replacement may require the return of product to Boston Scientific in accordance with the applicable Product Replacement Policy for Implanted Devices.

#### **WARRANTY**

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling and storage of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control may directly affect the device and the results obtained from its use. BSC's obligation under this warranty is limited to the replacement of this device and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. **BSC assumes no liability with respect to devices that are reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

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© 2019 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd