GreenLight HPS™

Laser Fiber

REF 0010-2090

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebrauchsanweisung	14
it	Istruzioni per l'uso	18
nl	Instructies voor gebruik	22
pt-EU	Instruções de Utilização	26

GreenLight HPS™

Laser Fiber

ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

These instructions detail the care and use of the GreenLight HPS Laser Fiber. These instructions are not recommendations for the medical or surgical application of the devices. The laser fibers are designed for use with compatible American Medical Systems/Boston Scientific Corporation (AMS/ BSC) GreenLight XPS™ laser systems. See the Laser System User's Manual, Safety and Professional Information Sections, for specific instructions concerning warnings, cautions, precautions, contraindications, and adverse events.

DEVICE DESCRIPTION

The GreenLight HPS Laser Fiber is sterile, single-use, disposable delivery system with a 600-micron fiber. Each device is designed to deliver laser energy at the 532 nm wavelength. The proprietary connector permits the laser to recognize the delivery device and activates the software safety features. The GreenLight HPS Laser Fiber is designed to deliver laser energy from the tip laterally at 70 degrees. Each device is intended to be used in surgical applications where lateral delivery of the laser energy is needed (i.e. procedures where lateral delivery of laser energy is normally accomplished with fibers alone or in conjunction with mirrors or fiber diverters).

The GreenLight HPS Laser Fiber is designed to accommodate a variety of surgical applications and endoscopic instrumentation. Each device can be used in non-contact and rotated 360 degrees, which allows tissue access in multiple planes and is compatible with standard endoscopes, cystoscopes and Albarran bridges. The device will cut, vaporize, or coagulate tissue in liquid endoscopic surgical procedures.


Figure 1. GreenLight HPS Laser Fiber

1. Connector (Laser firing input)
2. Control Knob (with alignment indicator)
3. Fiber Cap (Laser firing output including alignment indicators)

For correct alignment, the Blue Triangle (Alignment Arrow) must be facing up. If the Red Octagon is facing up, alignment is incorrect

Operating Principle

The GreenLight HPS Laser Fiber is a single-use disposable device intended for delivering 532 nm laser energy to the tissue. This laser energy transmitted within the 600-micron fiber core is delivered from the tip laterally at 70 degrees and upon impact cuts, vaporizes, ablates, or coagulates the tissue. The HPS laser fiber delivers 120 Watts maximum power.

It can be used for any application requiring the emission of side firing laser energy from the forward pointing direction of the axis of the laser fiber.

This device shall be compatible with the GreenLight HPS Laser System and the GreenLight XPS Laser System.

The laser fibers are compatible with most cystoscopes.

Table 1. Product Specifications

Product Description	GreenLight HPS Laser Fiber
Environment	Liquid
Wavelength	532 nm
Maximum Power (watts)*	120
Maximum Energy / Time	400 kJ
Limit Alert	<50,000 J/ <5 min
Time Out (Fiber card expires)	150 min. after 1st exposure
Pause Limit (Fiber card expires)	30 min. pause limit after 1st exposure
Time Out Warning	5 min. remaining
Overall length	3048 mm
Operating length (knob to distal tip)	357 mm
Core diameter	600 µm

*The recommended maximum powers shown in "Table 1. Product Specifications" are based on results obtained during clinical evaluations. They are not necessarily indicative of clinical performance.

USER INFORMATION

Physicians using the device should be appropriately trained in the use of medical lasers. Other users (prep nurse, operating tech) should have reasonable experience with common medical equipment and the procedure outlined in the Instructions for use (IFU). Only persons trained in the use of the laser fibers and who have sufficient medical knowledge (skilled urologist along with their team of technicians and/or nurses of varying experience and technical skill levels) should use the GreenLight HPS Laser Fiber. The user should have a demonstrated understanding of the laser fibers' operation and the risk associated with improper use.

CONTENTS

- Laser fiber
- Fiber card

MATERIALS

This device does not contain any regulated materials.

INTENDED USE

The GreenLight HPS Laser Fiber is indicated for endoscopic (cystoscopic) 532 nm laser procedures involving surgical incision, excision, vaporization, ablation, hemostasis and coagulation of urological soft tissue.

INDICATIONS FOR USE

- For the treatment of benign prostatic hypertrophy/hyperplasia (BPH).
- Other endoscopic procedures involving surgical incision, excision, vaporization, ablation, hemostasis, and coagulation of urological soft tissue in adults.

CLINICAL BENEFIT STATEMENT

The clinical benefits of HPS Fibers include improved urodynamics, improved quality of life (QoL) and improved International Prostate Symptoms Score (IPSS).

CONTRAINDICATIONS

The GreenLight HPS Laser Fiber should only be used by a qualified and trained surgeon. The use of the laser fiber is contraindicated for patients:

- Whose general medical condition contraindicates surgical intervention
- Where appropriate anesthesia is contraindicated by patient history
- Where tissue (especially tumors) has calcified
- For hemostasis of vessels over approximately two millimeters in diameter
- Where laser therapy is not considered the treatment of choice
- Who have bleeding disorders and coagulopathy
- Who have prostate cancer
- Who have acute urinary tract infection (UTI)
- Who have severe urethral stricture: a severe stricture shows visible narrowing on urethrography or ultrasonography, with near total obstruction that makes passage of instruments difficult or dangerous.

WARNINGS

BSC has no clinical information or experience concerning the use of the laser fiber in pediatric patients, pregnant women or nursing mothers.

Procedural Warnings

To avoid GreenLight HPS Laser Fiber damage or breakage:

- If using a deflecting mechanism, such as an Albarran bridge or a fiber diverter, do not attempt to remove the GreenLight HPS Laser Fiber while the bridge is deflected.

- In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of fiber or other material. If the aim beam or working beam exit the end of the fiber, discontinue use immediately and discard the fiber.
- Check the fiber package for damage. If damage is noted or sterility is in question, DO NOT USE.
- Reuse or misuse of a GreenLight HPS Laser Fiber will result in an increased potential for damage to the GreenLight Laser Fiber and ancillary equipment.

Usage Warnings

- Do not press the footswitch without having the appropriate protective eyewear. Pressing the footswitch may activate the laser beam. Activation of the laser beam when not wearing protective eyewear may result in eye injury to the user, personnel in the room, or the patient. Refer to the laser system's User Manual for requirements concerning eye protection.
- Reuse of a GreenLight HPS Laser Fiber can pose health and safety risks to patients and may compromise treatment outcomes.
- Do not fire the laser when the red octagon is visible as this may result in damage to the cystoscope or unintended tissue.
- Do not use if the fiber appears damaged. If damaged, replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.
- Never activate the laser energy if the fiber tip does not extend beyond the end of the endoscope.
- Care should be taken to prevent bladder perforation by avoiding excessive application of laser energy to the bladder wall. Low power and brief exposures should be employed until the effect of the laser on tissue is determined.
- Do not press the footswitch while checking the aim beam. Pressing the footswitch will activate the laser beam. Inappropriate activation of the laser beam may result in injury to the user, personnel in the room, or the patient.

Fire Hazard Warnings

- The GreenLight HPS Laser Fiber is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy.

To prevent drape fires or burns:

- Do not wrap any portion of the fiber in drapes or cloth.
- Do not attach the fiber directly to surgical drapes with a crushing clamp such as a hemostat.
- Do not place or drop instrumentation onto the fiber.
- Do not step on the fiber.
- Do not bend the fiber where it connects to the laser console.
- Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.

PRECAUTIONS

- When using accessories, tools, disposables, or materials that were in contact with the patient, take protective measures to prevent cross-contamination.
- The fiber and GreenLight XPS Laser System can be used in the treatment of urethral strictures with proper cautions. Use caution to avoid injury to urethral tissue. Refer to contraindications regarding severe urethral strictures.
- The fiber card must remain inserted until the procedure is complete. Premature removal of the card could result in the discontinued use of the fiber during the procedure.
- Do not activate the laser if the blue triangle is not visible. Damage may occur to the fiber or endoscope and require replacement of the damaged devices.

To avoid GreenLight HPS Laser Fiber damage or breakage:

- Do not bury the tip in tissue.
- Do not retract or probe tissue with the device.
- Do not bend at sharp angles.
- A sudden increase in bubbles at the probe tip indicates insufficient irrigation, tissue contamination or device failure. If increased irrigation does not eliminate the bubbling then stop, clean and inspect the device. Discontinue use and discard the fiber if damage is evident or abnormal bubbling occurs upon continuation of lasing.

ADVERSE EVENTS

Potential adverse events include, but are not limited to, the following:

- Abdominal bloating (intestinal gas)
- Acute renal failure
- Allergic reaction
- Aspiration
- Bladder neck contracture
- Bladder spasm
- Bleeding
- Burn
- Chills

- Clot retention
- Contamination of the device may lead to injury, illness, or death of the patient
- Deep venous thrombosis
- Delay in healing
- Dysuria
- Edema
- Embolism
- Epididymitis
- Erectile dysfunction (ED)
- Fatigue or weakness
- Fever
- Fluid overload/hyponatremia
- Gas over-distension
- Hematospermia
- Hematuria
- Infection
- Inflammation
- Leukocytosis
- Malfunction of laser fiber or console resulting in an injury or prolonged procedure
- Nocturia
- Overactive bladder
- Pain:
 - Abdominal pain unresponsive to Nonsteroidal Anti-inflammatory Drugs (NSAIDs)
 - Arm or leg pain
 - Headache
 - Back/low back pain
 - Body aches
 - Pelvic
 - Penile
 - Pelvic hematoma
 - Penile urethral injury
 - Perforation
 - Pneumothorax
 - Profuse perspiration (not fever-related)
 - Prostatitis
 - Pulmonary embolus
 - Retrograde ejaculation
 - Sepsis
 - Stricture
 - Tissue damage
 - Tissue sloughing
 - Ulceration
 - Unretrieved device fragment
 - Ureteral orifice injury
 - Urethral stricture
 - Urgency
 - Urinary frequency
 - Urinary incontinence
 - Urinary retention
 - Urinary tract infection

As with conventional endoscopic treatment, in extreme cases, death may occur due to procedural complications, concurrent illness, or the application of the laser. Use caution when treating patients who had difficulty with previous endoscopic procedures.

See the GreenLight XPS Laser System User's Manual for additional complications.

HOW SUPPLIED

GreenLight HPS Laser Fiber is supplied sterile as a single-use product not intended or validated for reprocessing or reuse.

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.
- Do not use past the use by date displayed on product label. See Disposal to properly dispose of this device.

HANDLING AND STORAGE

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Laser Console Set-up

1. Before use of the laser console, see the laser system User's Manual for information and instructions for use. Ensure the user, the patient, and all personnel in the room are wearing appropriate protective eyewear before using the laser console.
2. Turn the laser console on. Wait for self-test to be completed and the prompt "Insert fiber card" to appear.
3. Remove the sterile pouch containing the fiber card from the box. Do not open the pouch.
4. Remove fiber card from the pocket on the outside of the pouch. Insert card into the reader with chip side facing operator.
5. After the fiber card is recognized by the laser console, there will be a prompt to "Attach a device."

Before the Procedure

WARNING: BSC has no clinical information or experience concerning the use of the laser fiber in pediatric patients, pregnant women or nursing mothers.

PRECAUTION: The fiber and GreenLight XPS Laser System can be used in the treatment of urethral strictures with proper cautions. Use caution to avoid injury to urethral tissue. Refer to contraindications regarding severe urethral strictures.

1. Before use, inspect the fiber package for visible damage.

WARNING: Check the fiber package for damage. If damage is noted or sterility is in question, DO NOT USE.

WARNING: The GreenLight HPS Laser Fiber is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy.

WARNING: To prevent drape fires or burns:

- Do not wrap any portion of the fiber in drapes or cloth.
- Do not attach the fiber directly to surgical drapes with a crushing clamp such as a hemostat.
- Do not place or drop instrumentation onto the fiber.
- Do not step on the fiber.
- Do not bend the fiber where it connects to the laser console.
- Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.

WARNING: Reuse or misuse of a GreenLight HPS Laser Fiber will result in an increased potential for damage to the GreenLight HPS Laser Fiber and ancillary equipment.

2. Using aseptic technique, open the fiber pouch.
3. Remove the fiber from the pouch and the protective tube. Place the fiber onto the sterile field.

PRECAUTION: When using accessories, tools, disposables, or materials that were in contact with the patient, take protective measures to prevent cross-contamination.

4. Inspect the fiber for visible damage (e.g. broken, crushed or kinked).

WARNING: Do not use if the fiber appears damaged. If damaged, replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.

PRECAUTION: To avoid GreenLight HPS Laser Fiber damage or breakage do not bend at sharp angles.

5. Insert the connector of the fiber into the fiber port on the laser console and turn clockwise until it locks (1/4 turn).

WARNING: Do not press the footswitch while checking the aim beam. Pressing the footswitch will activate the laser beam. Inappropriate activation of the laser beam may result in injury to the user, personnel in the room, or the patient.

Performing the Procedure

PRECAUTION: The fiber card must remain inserted until the procedure is complete. Premature removal of the card could result in the discontinued use of the fiber during the procedure.

1. Ensure user, patient, and personnel in the room are wearing appropriate protective eyewear before using the laser.

WARNING: Do not press the footswitch without having the appropriate protective eyewear. Pressing the footswitch may activate the laser beam. Activation of the laser beam when not wearing protective eyewear may result in eye injury to the user,

personnel in the room, or the patient. Refer to the laser system's User Manual for requirements concerning eye protection.

2. Connect the tubing from the sterile saline for irrigation solution to the endoscope luer lock connector. Elevate the saline solution to ensure proper flow.

PRECAUTION: A sudden increase in bubbles at the probe tip indicates insufficient irrigation, tissue contamination or device failure. If increased irrigation does not eliminate the bubbling then stop, clean and inspect the device. Discontinue use and discard the fiber if damage is evident or abnormal bubbling occurs upon continuation of lasing.

3. Open the flow valve on the endoscope.

WARNING: If using a deflecting mechanism, such as an Albarran bridge or a fiber diverter, do not attempt to remove the GreenLight HPS Laser Fiber while the bridge is deflected.

4. Insert and advance the endoscope to the treatment area using standard medical procedure.
5. Insert the fiber into the endoscope until the distal end of the fiber is visualized. A blue triangle and a red octagon are located on fiber to provide an extension reference.
6. Laser light is emitted opposite the blue triangle. The blue triangle must be visible when firing the laser.

WARNING: Do not fire the laser when the red octagon is visible as this may result in damage to the endoscope or unintended tissue.

WARNING: Care should be taken to prevent bladder perforation by avoiding excessive application of laser energy to the bladder wall. Low power and brief exposures should be employed until the effect of the laser on tissue is determined.

7. Place the laser in Ready mode. This will enable the footswitch control.

WARNING: Never activate the laser energy if the fiber tip does not extend beyond the end of the endoscope.

8. Determine the direction of the laser beam by using the red colored aiming beam. The laser beam is inline with the raised indicator on the Control knob.

PRECAUTION: Do not activate the laser if the blue triangle is not visible. Damage may occur to the fiber or endoscope and require replacement of the damaged devices.

9. With the aim beam on the targeted tissue, depress the footswitch and begin lasing. The treatment times will vary based on distance to tissue, the power settings, duration of beam application, and other factors. Use settings that are as low as possible for the desired treatment effect.

PRECAUTION: To avoid GreenLight HPS Laser Fiber damage or breakage do not bury the tip in tissue.

PRECAUTION: To avoid GreenLight HPS Laser Fiber damage or breakage do not retract or probe tissue with the device.

PRECAUTION: To avoid GreenLight HPS Laser Fiber damage or breakage do not bend at sharp angles.

10. The laser fiber can be used up to an energy or time limit. The user is alerted by an audible tone and a screen message per the Product Specifications. Once the energy or time limit is reached, to continue the procedure, a new fiber card and fiber must be inserted.

If, during the procedure, the tip accumulates debris, turn the laser to the STANDBY mode and carefully wipe the tip clean with a wet sterile gauze or towel. Wipe the laser fiber, starting at the tip and moving down the laser fiber. Frequent cleaning of the laser fiber will help maintain vaporization efficiency and increase the laser fiber life.

If excessive tissue or debris is allowed to accumulate on the laser fiber cap, over heating of the laser fiber cap will occur leading to premature laser fiber degradation and/or laser fiber failure.

WARNING: In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of fiber or other material. If the aim beam or working beam exit the end of the fiber, discontinue use immediately and discard the fiber.

After the Procedure

1. With the endoscope still in place, pull the fiber into the scope approximately 5 cm.
2. Confirm the amount of desired tissue is removed.
3. Observe for bleeding and restart irrigation and coagulate, as necessary.
4. After use, dispose of product and packaging in accordance with hospital, administrative, and local government policy.

WARNING: Reuse of a GreenLight HPS Laser Fiber can pose health and safety risks to patients and may compromise treatment outcomes.

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

DISPOSAL

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

PATIENT COUNSELING INFORMATION

Prior to use of this device, discuss the risks and benefits of Greenlight Laser Therapy and associated fibers, alternative therapies, and obtain the patient's informed consent. Brief the patient on any contraindications, warnings, precautions, potential adverse events, undesirable side effects and relevant post-procedure instructions. Advise the patient regarding the expected post-operative course and when return to regular physical activities, as well as sexual intercourse, may be resumed. Inform the patient to contact their physician with any problems or concerns.

WARRANTY

For device warranty information, visit www.bostonscientific.com/Warranty.

The following are trademarks of Boston Scientific Corporation or its affiliates: GreenLight HPS, GreenLight XPS.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary. Additional symbols are defined at the end of this document.

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet. Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
© 2022 Boston Scientific Corporation or its affiliates.