MoXy™

LIQUID COOLED

Laser Fiber

REF 0010-2400

en	Instructions for Use	2
es	Instrucciones de uso	6
fr	Mode d'emploi	10
de	Gebrauchsanweisung	14
it	Istruzioni per l'uso	19
nl	Instructies voor gebruik	23
pt-EU	Instruções de Utilização	27

MoXy™

LIQUID COOLED

Laser Fiber

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

These instructions detail the care and use of the MoXy Laser Fiber. These instructions are not recommendations for the medical or surgical application of the device. The MoXy Laser Fiber is designed for use with the GreenLight XPS™ Laser System. See the GreenLight XPS Laser System User's Manual for specific instructions concerning warnings, precautions, contraindications, and adverse events.

DEVICE DESCRIPTION

The MoXy Laser Fiber features a side firing mechanism delivering up to 180 W of 532 nm light to tissue. The fiber is designed/intended to be used with the Greenlight XPS Laser System and delivers laser energy to the desired treatment site via an endoscope/cystoscope. The fiber can be used for the surgical incision/excision, vaporization, ablation, hemostasis, and coagulation of urological soft tissue.

The fiber is supplied as a single-use, sterile product and is compatible with the GreenLight XPS Laser System. It can access the tissue in multiple planes. The fiber is a liquid cooled delivery device, which enables the delivery of up to 180 W of power and aids in maintaining a clear environment at the fiber cap.


Figure 1. MoXy Laser Fiber

- 1. Metal cap
- 2. Laser firing point
- 3. Blue triangle; red octagon located on opposite side
- 4. Outer flow tubing
- 5. Control knob
- 6. Beam alignment indicator
- 7. Fiber handle
- 8. Inlet flow tubing
- 9. Luer lock connector
- 10. Connector

The GreenLight™ MoXy Laser Fiber was designed to deliver 50% more power (180 W) and maintain the safety profile of the predicate GreenLight HPS™ fibers 0010-2090 (120 W) and 0010-2092 (40 W). By increasing the core diameter of the GreenLight MoXy Laser Fiber, the beam area is increased targeting more tissue than the GreenLight HPS Laser Fibers. The average power density of 180 W of laser power delivered using this fiber is comparable to the 120 W power density of laser power delivered using the GreenLight HPS Laser Fiber. With comparable power densities, the tissue is removed to a similar depth and the heat affect zone (HAZ) for both devices falls within the acceptable safety range of 0.5 mm to 5 mm. Therefore, tissue continues to be vaporized safely and effectively.

The active Cooling Cap™ Feature of the GreenLight MoXy Laser Fiber prevents overheating and beam devitrification at the laser firing point. The internal fluid channel is created via a metal cap that is bonded to fiber's glass cap, and nylon outer flow tubing that surrounds and travels the length of the laser fiber. This fluid channel provides saline flow down the length of the laser fiber in addition to the saline irrigation flow between the cystoscope and the fiber. The saline flowing down both fiber and the cystoscope helps expel tissue debris from the firing point during the procedure. In addition to removing the tissue debris, the internal and external irrigation fluid reduces heat generated on the laser firing point through turbulent heat convection flow.

OPERATING PRINCIPLE

The MoXy Laser Fiber is a single-use, disposable device intended for delivering up to 180 W of 532 nm laser energy to the tissue. This laser energy transmitted within the 750-micron fiber core is delivered from the tip laterally at 70 to 80 degrees and upon impact cuts, vaporizes, ablates or coagulates the tissue. It can be used for any application requiring the side fire emission of laser energy from the forward pointing direction of the axis of the fiber.

This device will be compatible with the GreenLight XPS Laser System. The laser fiber is compatible with most cystoscopes. The laser fiber is cooled during the procedure with a saline solution increasing the life of device

The laser fiber and cap assembly have an outer diameter of ≤2.3 mm. The laser fiber is designed to be compatible with cystoscopes which have a minimum working channel of 7.5 Fr.

Table 1. Product Specifications

Product Description	MoXy Laser Fiber
Environment	Liquid
Wavelength	532 nm
Maximum Power (watts)*	180
Maximum Energy / Time	650 kJ
Limit Alert	<50,000 J/ <5 min
Time Out (Fiber card expires)	150 min. after 1st exposure
Pause Limit (Fiber card expires)	30 min. pause limit after 1st exposure
Time Out Warning	5 min. remaining
Overall length	3048 mm
Operating length (knob to distal tip)	360 mm
Core diameter	750 µm

*The recommended maximum powers shown in "Table 1. Product Specifications" are based on results obtained during clinical evaluations. They are not necessarily indicative of clinical performance.

USER INFORMATION

Only persons trained in the use of the laser fibers and who have sufficient medical knowledge (skilled urologist along with their team of technicians and/or nurses of varying experience and technical skill levels) should use the MoXy Laser Fiber. The user should have a demonstrated understanding of the laser fibers' operation and the risk associated with improper use.

CONTENTS

- Laser fiber
- Fiber card

MATERIALS

Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

Note: Metal Cap is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

INTENDED USE

The GreenLight Moxy Laser Fiber is indicated for endoscopic (cystoscopic) 532 nm laser procedures involving surgical incision, excision, vaporization, ablation, hemostasis and coagulation of urological soft tissue.

INDICATIONS FOR USE

- For the treatment of benign prostatic hypertrophy/hyperplasia (BPH).
- Other endoscopic procedures involving surgical incision, excision, vaporization, ablation, hemostasis, and coagulation of urological soft tissue in adults.

CLINICAL BENEFIT STATEMENT

The clinical benefits of Moxy Laser Fiber include improved urodynamics, improved quality of life (QoL) and improved International Prostate Symptoms Score (IPSS).

CONTRAINDICATIONS

The use of the fiber and GreenLight XPS Laser System are contraindicated for patients with the following conditions:

- Whose general medical condition contraindicates surgical intervention
- Where appropriate anesthesia is contraindicated by patient history
- Who have calcified tissue (especially tumors)
- For hemostasis of vessels over approximately two millimeters in diameter
- Where laser therapy is not considered the treatment of choice
- Who have uncontrolled bleeding disorders and coagulopathy
- Who have prostate cancer
- Who have acute Urinary Tract Infection (UTI)
- Who have severe urethral stricture: a severe stricture shows visible narrowing on urethrography or ultrasonography, with near total obstruction that makes passage of instruments difficult or dangerous.

WARNINGS

BSC has no clinical information or experience concerning the use of the laser fiber in pediatric patients, pregnant women or nursing mothers.

Procedural Warnings

To avoid fiber damage or breakage:

- Do not use if the package is damaged. Damage to the package may result in damage to the fiber or compromised sterility.
- Do not depress the footswitch to activate the working beam, unless the aim beam is visible on the intended tissue.
- Reuse or misuse of the fiber will result in an increased potential for damage to the fiber and ancillary equipment.
- In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of fiber or other material. If the aim beam or working beam exit the end of the fiber, discontinue use immediately and discard the fiber.

Usage Warnings

- Do not fire the laser when the red octagon is visible as this may result in damage to the cystoscope or unintended tissue.
- The laser beam is emitted at an angle. Avoid forward deflection of the beam at the bladder neck as damage to the ureteral orifices may occur.
- Do not press the footswitch while checking the aim beam. Pressing the footswitch will activate the laser beam. Inappropriate activation of the laser beam may result in injury to the user, personnel in the room, or the patient.
- Do not use if the red colored aim beam is not visible. Replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.
- Do not press the footswitch without having the appropriate protective eyewear. Pressing the footswitch may activate the laser beam. Activation of the laser beam when not wearing protective eyewear may result in eye injury to the user, personnel in the room, or the patient. Refer to the laser system's User Manual for requirements concerning eye protection.
- Do not use if the fiber appears damaged. If damaged, replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.
- Reuse of a fiber can pose health and safety risks to patients and may compromise treatment outcomes.
- Never activate the laser energy if the fiber tip does not extend beyond the end of the endoscope.

Fire Hazard Warnings

- The fiber is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy. Do not excessively manipulate or bend the fiber.

To prevent drape fires or burns:

- Do not wrap portion of the fiber in drapes or cloth.
- Do not attach the fiber directly to surgical drapes with a crushing clamp such as a hemostat.
- Do not place or drop instrumentation onto the fiber.
- Do not step on the fiber.
- Do not bend the fiber where it connects to the laser console.
- Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.

PRECAUTIONS

- When using accessories, tools, disposables, or materials that were in contact with the patient, take protective measures to prevent cross-contamination.
- The fiber and GreenLight XPS Laser System can be used in the treatment of urethral strictures with proper cautions. Use caution to avoid injury to urethral tissue. Refer to contraindications regarding severe urethral strictures.

- The fiber card must remain inserted until the procedure is complete. Premature removal of the card will result in the discontinued use of the fiber during the procedure.
- Failure to position the saline solution at least 106 cm or not opening the flow valve will decrease the amount of fluid across the fiber end. Decreased fluid flow may limit the lasing time of the fiber.
- Do not depress the footswitch to activate the working beam, if the aim beam exits straight from the end of the fiber. Replace the fiber.
- Do not activate the laser if the blue triangle is not visible. Damage may occur to the fiber or endoscope and require replacement of the damaged devices.
- Do not bury the tip of the fiber into the tissue. Reduced irrigation fluid flow may result in damage to the fiber.
- Do not retract, dissect, or probe tissue with the tip of the fiber. Damage may occur to the fiber tip.
- Do not press the fiber end into the tissue, which will create stress between the fiber and the opening (edge) of the endoscope. Damage to the fiber may result.
- Do not bend the fiber at sharp angles. Damage may occur to the fiber

ADVERSE EVENTS

Potential adverse events include, but are not limited to, the following:

- Abdominal bloating (intestinal gas)
- Acute renal failure
- Allergic reaction
- Aspiration
- Bladder neck contracture
- Bladder spasm
- Bleeding
- Burn
- Chills
- Clot retention
- Contamination of the device may lead to injury, illness, or death of the patient
- Deep venous thrombosis
- Delay in healing
- Dysuria
- Edema
- Embolism
- Epididymitis
- Erectile dysfunction (ED)
- Extravasation
- Fatigue or weakness
- Fever
- Fluid overload/hyponatremia
- Gas over-distension
- Hematospermia
- Hematuria
- Infection
- Inflammation
- Leukocytosis
- Malfunction of laser fiber or console resulting in an injury or prolonged procedure
- Nocturia
- Overactive bladder
- Pain:
 - Abdominal pain unresponsive to Nonsteroidal Anti-inflammatory Drugs (NSAIDs)
 - Arm or leg pain
 - Headache
 - Back/low back pain
 - Body aches
 - Pelvic
 - Penile
- Pelvic hematoma
- Penile urethral injury
- Perforation
- Pneumothorax
- Profuse perspiration (not fever-related)
- Prostatitis

- Pulmonary embolus
- Retrograde ejaculation
- Sepsis
- Stricture
- Tissue damage
- Tissue sloughing
- Ulceration
- Unretrieved device fragment
- Ureteral orifice injury
- Urethral stricture
- Urgency
- Urinary frequency
- Urinary incontinence
- Urinary retention
- Urinary tract infection

As with conventional endoscopic treatment, in extreme cases, death may occur due to procedural complications, concurrent illness, or the application of the laser. Use caution when treating patients who had difficulty with previous endoscopic procedures. See the GreenLight XPS Laser System User's Manual for additional complications.

HOW SUPPLIED

MoXy Laser Fiber is supplied sterile as a single-use product not intended or validated for reprocessing or reuse.

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.
- Do not use past the use by date displayed on product label. See Disposal to properly dispose of this device.

HANDLING AND STORAGE

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Laser Console Set-up

1. Before use of the laser console, see the GreenLight XPS Laser System User's Manual for information and instructions for use. Ensure the user, the patient, and all personnel in the room are wearing appropriate protective eyewear before using the laser console.
2. Turn the laser console on. Wait for self-test to be completed and the prompt 'Insert Fiber Card' to appear.
3. Remove the sterile pouch containing the fiber card from the box. Do not open the pouch.
4. Remove fiber card from the pocket on the outside of the pouch. Insert card into the reader with chip side facing operator.
5. After the fiber card is recognized by the laser console, there will be a prompt to 'Attach a Fiber'.

PRECAUTION: The fiber card must remain inserted until the procedure is complete. Premature removal of the card will result in the discontinued use of the fiber during the procedure.

Before the Procedure

WARNING: BSC has no clinical information or experience concerning the use of the laser fiber in pediatric patients, pregnant women or nursing mothers.

PRECAUTION: The fiber and GreenLight XPS Laser System can be used in the treatment of urethral strictures with proper cautions. Use caution to avoid injury to urethral tissue. Refer to contraindications regarding severe urethral strictures.

1. Before use, inspect the fiber package for visible damage.

WARNING: Do not use if the package is damaged. Damage to the package may result in damage to the fiber or compromised sterility.

WARNING: The fiber is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy. Do not excessively manipulate or bend the fiber.

WARNING: To prevent drape fires or burns:

- Do not wrap portion of the fiber in drapes or cloth.
- Do not attach the fiber directly to surgical drapes with a crushing clamp such as a hemostat.
- Do not place or drop instrumentation onto the fiber.
- Do not step on the fiber.
- Do not bend the fiber where it connects to the laser console.
- Do not use this laser fiber in the presence of flammable anesthetics, combustible materials, or in an oxygen enriched environment due to risk of fire or explosion.

WARNING: Reuse or misuse of the fiber will result in an increased potential for damage to the fiber and ancillary equipment.

2. Using aseptic technique, open the fiber pouch.
3. Remove the fiber from the pouch and the protective tube. Place the fiber onto the sterile field.

PRECAUTION: When using accessories, tools, disposables, or materials that were in contact with the patient, take protective measures to prevent cross-contamination.

4. Inspect the fiber for visible damage (e.g. broken, crushed, or kinked).

WARNING: Do not use if the fiber appears damaged. If damaged, replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.

PRECAUTION: Do not bend the fiber at sharp angles. Damage may occur to the fiber.

5. Insert the connector of the fiber into the fiber port on the laser console and turn clockwise until it locks (1/4 turn).

WARNING: Do not press the footswitch while checking the aim beam. Pressing the footswitch will activate the laser beam. Inappropriate activation of the laser beam may result in injury to the user, personnel in the room, or the patient.

WARNING: Do not use if the red colored aim beam is not visible. Replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.

Performing the Procedure

1. Ensure user, patient, and personnel in the room are wearing appropriate protective eyewear before using the laser.

WARNING: Do not press the footswitch without having the appropriate protective eyewear. Pressing the footswitch may activate the laser beam. Activation of the laser beam when not wearing protective eyewear may result in eye injury to the user, personnel in the room, or the patient. Refer to the laser system's User Manual for requirements concerning eye protection.

2. Connect the tubing from the sterile saline for irrigation solution to the fiber luer lock connector. Position the saline solution at least 106 cm higher than the level of the endoscope. Increased irrigation flow by means of a saline pressure bag (set to 250 mmHg – 300 mmHg) will further increase liquid cooling effect and may reduce fiber tip damage.
3. Open the flow valve. Fluid should be observed flowing out of the metal cap window.

PRECAUTION: Failure to position the saline solution at least 106 cm or not opening the flow valve will decrease the amount of fluid across the fiber end. Decreased fluid flow may limit the lasing time of the fiber.

4. Insert and advance the endoscope to the treatment area using standard medical procedure.
5. Insert the fiber into the endoscope until the distal end of the fiber is visualized. A blue triangle and a red octagon are located on fiber to provide an extension reference.
6. Laser light is emitted opposite the blue triangle. The blue triangle must be visible when firing the laser.

WARNING: Do not fire the laser when the red octagon is visible as this may result in damage to the endoscope or unintended tissue.

WARNING: The laser beam is emitted at an angle. Avoid forward deflection of the beam at the bladder neck as damage to the ureteral orifices may occur.

7. Place the console in Ready mode to enable the footswitch.

WARNING: Never activate the laser energy if the fiber tip does not extend beyond the end of the endoscope.

8. Determine the direction of the laser beam by using the red colored aiming beam. The laser beam is inline with the raised indicator on the Control knob.

WARNING: Do not depress the footswitch to activate the working beam unless the aim beam is visible on the intended tissue.

PRECAUTION: Do not depress the footswitch to activate the working beam if the aim beam exits straight from the end of the fiber. Replace the fiber.

9. With the aim beam on the targeted tissue, depress the footswitch and begin lasing. The treatment times will vary based on distance to tissue, the power settings, duration of beam application, and other factors. Use settings that are as low as possible for the desired treatment effect.

PRECAUTION: Do not activate the laser if the blue triangle is not visible. Damage may occur to the fiber or endoscope and require replacement of the damaged devices.

PRECAUTION: Do not bury the tip of the fiber into the tissue. Reduced irrigation fluid flow may result in damage to the fiber.

PRECAUTION: Do not retract, dissect, or probe tissue with the tip of the fiber. Damage may occur to the fiber tip.

PRECAUTION: Do not press the fiber end into the tissue, which will create stress between the fiber and the opening (edge) of the endoscope. Damage to the fiber may result.

PRECAUTION: Do not bend the fiber at sharp angles. Damage may occur to the fiber.

10. The fiber can be used up to an energy or time limit. The user is alerted by an audible tone and a screen message when less than 50,000 Joules or 5 minutes remaining. Once the energy or time limit is reached, to continue the procedure, a new fiber card and fiber must be inserted.

WARNING: In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of fiber or other material. If the aim beam or working beam exit the end of the fiber, discontinue use immediately and discard the fiber.

After the Procedure

1. With the endoscope still in place, pull the fiber into the scope approximately 5 cm.
2. Confirm the amount of desired tissue is removed.
3. Observe for bleeding and restart irrigation and coagulate, as necessary.
4. After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

WARNING: Reuse of a fiber can pose health and safety risks to patients and may compromise treatment outcomes.

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

Table 2. Troubleshooting and Procedural Considerations

Condition	Steps
A sudden increase in bubbles at the fiber tip	1. Increase the amount of irrigation through the endoscope. 2. If increased irrigation does not eliminate the bubbling, repeat the step above. 3. If the problem persists, discontinue use and replace with a new fiber.
Lack of bubbles inefficient vaporization due to the sweep speed and distance from tissue	1. Adjust sweep rate of the laser beam for even tissue removal. 2. Ensure the distance from the fiber to the tissue is approximately 2 mm (1 to 3 mm). 3. Adjust the output of the laser.
Accumulation of debris in the metal cap window or on the fiber	NOTE: The liquid cooling design of the MoXy Laser Fiber eliminates the need to clean the fiber during recommended use. CAUTION: Clean the fiber only when tissue or other foreign material is forced into the window. Unnecessary cleaning may damage the fiber. 1. Turn the laser console to the Standby mode. 2. Remove the fiber from the endoscope. 3. Gently clean the metal window opening with wet sterile gauze or pointed swab. 4. Insert the fiber into the endoscope and resume the procedure. WARNING: Failure to place the laser console in the Standby mode may result in user, personnel, or patient exposure to laser energy and potential injury. CAUTION: Excessive or rough cleaning of the fiber tip may result in damage to the fiber. Accumulation of debris may result in over heating of the fiber cap leading to premature fiber degradation or fiber failure.

DISPOSAL

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

PATIENT COUNSELING INFORMATION

Prior to use of this device, discuss the risks and benefits of Greenlight Laser Therapy and associated fibers, alternative therapies, and obtain the patient's informed consent. Brief the patient on any contraindications, warnings, precautions, potential adverse events, undesirable side effects and relevant post-procedure instructions. Advise the patient regarding the expected post-operative course and when return to regular physical activities, as well as sexual intercourse, may be resumed. Inform the patient to contact their physician with any problems or concerns.

WARRANTY

For device warranty information, visit www.bostonscientific.com/Warranty

The following are trademarks of Boston Scientific Corporation or its affiliates: GreenLight, GreenLight XPS, MoXy, Cooling Cap, GreenLight HPS. All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary

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INSTRUCTIONS FOR USE

en

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet. Visit www.IFU-BSCI.com to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, please contact Boston Scientific Customer Service or your local country contact. A copy will be sent to you at no charge and should arrive within seven days.

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All rights reserved.
Boston Scientific Limited
© 2023 Boston Scientific Corporation or its affiliates.