






IVUS Set-Up and Workflow

Imaging Catheters Preparation Instructions

- 1** **Wet assemble** the accessories (ensure air is removed from all the components).



- 2** **Fully open the telescope.**



- 3** **Flush the first 3 cc syringe slowly (first 1 cc is the most important to flush slowly).** Refill the 3 cc syringe and flush a second time with the telescope open.



- 4** If not already done, connect the MDU Fast Pullback Motor Drive to the sled. Position the MDU in the back end of the Pullback Sled.

- 5** Cover the MDU Fast Pullback Motor Drive and the pre-connected Permanent Sled with the **Sterile Bag**.



- 6** **Connect the catheter** to the MDU Fast Pullback.










IVUS Set-Up and Workflow

Imaging Catheters Preparation Instructions

7

a. **Turn on « Imaging »** via the console , or on the MDU Fast Pullback Motor Drive by pressing on the Imaging button. 




b. **Close telescope.** Press the **Brake Release Trigger**  and **push forward the MDU at the front** of the Permanent Sled. **Lock the red anchor**  of the catheter to the Permanent Sled U-shaped Metallic Connector and its appropriate position on the Permanent Sled Bag, identified with a „2“.




8

Flush again and confirm image quality.



9

Turn off the « Imaging » MDU Fast Pullback Motor Drive by re-clicking on the Imaging on the console or the MDU button. 

To prevent catheter damaging and/or complications, do not advance the IVUS catheter whilst imaging is activated.

10

Following standard best practices: introduce the imaging catheter on the guidewire and into the patient. It's recommended to position the IVUS probe past the zone of interest if anatomy allows safely.





IVUS Set-Up and Workflow

IVUS Software Workflow




IVUS 123 Essential Workflow


- 1 **Turn Imaging On**
 - + Click IMAGING to start imaging and to confirm the sensor is at the right position.
- 2 **Record Image**
 - + Click REC to start and stop recording.
 - + The pullback method is auto-selected, slide to manual to continue with manual pullback.
 - + Bookmarks can be added while recording. To stop recording, click the stop button.
- 3 **Record New Run**
 - + Reposition the sensor to the desired area of interest.
 - + Click on LIVE at the top of the screen to be able to turn on imaging and/or recording.



PRE-PCI WORKFLOW



POST-STENT WORKFLOW



We are Modern PCI

SEE

PREP

TREAT

All trademarks are property of their respective owner. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.