# AMS 700™ with MS Pump™

## *Penile Prosthesis Product Line*

### Instructions for Use

English	AMS 700™ with MS Pump™ <i>Penile Prosthesis Product Line</i> Instructions for Use .....	3
Deutsch	AMS 700™ mit MS Pump™ <i>Produktreihe der Schwellkörperimplantate</i> Gebrauchsanweisung .....	11
Español	AMS 700™ con bomba MS Pump™ <i>Línea de prótesis de pene</i> Instrucciones de uso .....	21
Français	AMS 700™ avec MS Pump™ <i>Gamme de prothèses péniennes</i> Mode d'emploi.....	33
Italiano	AMS 700™ con MS Pump™ <i>Linea di prodotti di protesi peniene</i> Istruzioni per l'uso .....	45
Nederlands	AMS 700™ met MS Pump™ <i>Productlijn penisprothesen</i> Gebruiksaanwijzing .....	55
Português	AMS 700™ com MS Pump™ <i>Linha de produtos de próteses penianas</i> Instruções de utilização .....	65
Svenska	AMS 700™ med MS Pump™ <i>Penisprotesproduktssortiment</i> Bruksanvisning .....	75

**Rx ONLY**

**AMS™**

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## *Penile Prosthesis Product Line*

### Instructions For Use

**NOTE:** Refer to the Operating Room Manual for further information on the AMS 700™ product line and their implantation.

#### Device Description

The AMS 700 Series Inflatable Penile Prosthesis product line includes the AMS 700 LGX™ Preconnect, AMS 700 CX Preconnect, AMS 700 LGX, AMS 700 CX, and the AMS 700 CXR Preconnect Penile Prosthesis. These configurations are available with and without InhibiZone™, an antibiotic impregnation of rifampin (rifampicin) and minocycline hydrochloride (HCl).<sup>\*</sup> Use of an InhibiZone treated penile prosthesis results in a significant reduction in the rate of revisions due to infection in patients receiving a first time or a revision AMS 700 device compared to an untreated AMS 700 device (this benefit is also seen among high risk diabetic patients receiving a first time AMS 700 implant). This clinical data demonstrates that the use of InhibiZone reduces device-related infections. These prostheses are closed fluid-filled systems consisting of a pair of cylinders, optional rear tip extenders (RTEs), a Momentary Squeeze (MS) pump, and a fluid reservoir. All components are connected by kink-resistant tubing (KRT). The cylinders are inflated as fluid is pumped from the reservoir, creating an erection. They are deflated as fluid is transferred back to the reservoir, making the penis flaccid once again. This device contains solid silicone elastomer. These devices are for men who, after appropriate patient history and diagnostic evaluations as well as discussions with the urologist about other alternative treatment methods, are determined to be suitable candidates for implantation surgery.

<sup>\*</sup> not available in all markets.

#### Indications for Use

The AMS 700 Series Inflatable Penile Prosthesis product line is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

#### Contraindications

The implantation of this device is contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery.

The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin (rifampicin), or to minocycline HCl, or other tetracyclines.

The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to aggravate this condition.

#### Warnings

1. Implantation of the device will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible.
2. Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.
3. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.
4. Implantation of a penile prosthesis may result in penile shortening, curvature or scarring.
5. This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
6. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.
7. If a hypersensitivity reaction develops to a device coated with InhibiZone, the penile prosthesis should be removed and the patient treated appropriately.

#### Precautions

##### **Surgery Related**

1. Improper reservoir placement or filling technique can result in spontaneous unintended inflation or deflation of the cylinders that may result in unintended partial or full erections.
2. Migration of the device components can occur if the cylinders are improperly sized, if the pump or the reservoir is not positioned properly, or if the tubing lengths are incorrect.
3. Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or may make it impossible.
4. Improper measurement technique, positioning or sizing may reduce cylinder life.
5. Unsuccessful outcomes have been reported due to improper surgical technique, anatomical misplacement of components, improper sizing and filling of components, or tubing kinks.
6. Implantation of AMS 700 LGX cylinders in patients with Peyronie's disease may not provide a satisfactory result.

##### **Device Related**

1. AMS Quick Connect Sutureless Window Connectors, provided in the AMS 700 Accessory Kit, should not be used in revision procedures involving previously implanted component tubing. In this situation the Quick Connect Sutureless Window Connectors may be less effective.
2. Some of the materials used in the construction of this device have been shown to cause minor irritation when implanted in animals. Therefore, implantation of this device may cause minor irritation or discomfort in some patients.
3. Devices in the AMS 700 product line should be filled with sterile, normal saline. Some patients may have a hypersensitivity to contrast media.
4. Do not use product that has damaged or open packaging, as sterility may be compromised.
5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped in a sterile solution immediately prior to implant, if desired.
7. CXR RTEs are not compatible with CX or LGX cylinders.
8. Verify proper attachment of RTEs by spinning them once seated. Properly attached RTEs should spin freely without accidental disengagement or material bulging.

- Do not stack the CX, LGX or CXR Snapcone RTEs with the exception of the 1.5 cm. The locking ring will not engage with the smooth outer surface of the RTE, which may result in the RTE disconnecting.

### Patient Related

- A thorough preoperative consultation should include a discussion between patient and physician of all available treatment options and their risks and benefits.
- Adequate patient manual dexterity and strength are required for proper device inflation and deflation.
- Mental or psychological conditions, such as senile dementia, may inhibit the patient's successful operation of the prosthesis.
- Trauma to the pelvic or abdominal areas, such as impact injuries associated with sports (eg, bicycle riding), can result in damage of the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction, including replacement of the device.
- The contour, elasticity, and dimension of the tunica albuginea may limit the length and/or diameter expansion of the AMS 700 cylinders.
- The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.
- Use of injection therapy concurrently with the penile prosthesis can damage the prosthesis. Patients should not use injection therapy after they receive their implant.

### InhibiZone Related

- InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used when implanting an inflatable penile prosthesis.
- Because the products with InhibiZone are impregnated with a combination of rifampin (rifampicin) (a derivative of rifamycin B) and minocycline HCl (a derivative of tetracycline), the contraindications, warnings, and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline HCl and rifampin (rifampicin) in patients receiving this device are unlikely to be detected.
- Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin (rifampicin) and minocycline HCl can cause additional stress on the hepatic and renal systems.
- Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
- Patients who receive a device with InhibiZone and are also taking warfarin should have their prothrombin time monitored, because tetracyclines have been reported to slow coagulation.
- Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid, and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin (rifampicin).

### Adverse Events

A clinical trial was conducted to determine the safety and effectiveness of the AMS 700 Series of inflatable penile prostheses. This trial involved only devices without InhibiZone. A total of 300 patients were enrolled with follow-up out to 5 years for 126 patients. The Adverse Device Effects (ADE), detailed in the following table, were noted during the duration of this clinical trial for all enrolled patients.

**Table 1: AMS Clinical Trial Adverse Device Effects**

ADE	# Patient (%**)	Mean Onset Time in Days (Range in Days)
Urogenital Pain (Typically Associated with Healing Process)	160 (53.3%)	21 (0 – 876)
Urogenital Edema	106 (35.3%)	8 (0 – 722)
Urogenital Ecchymosis	30 (10.0%)	4 (0 – 150)
Reservoir Encapsulation (persistent in 11/19 cases)	19 (6.3%)	275 (38 – 1731)
Patient Dissatisfaction (With Length, Ability to Use and Nonspecific Reasons)	18 (6.0%)	384 (0 – 1830)
Auto-Inflation	17 (5.7%)	141 (0 – 608)
Mechanical Malfunction (Leaks, Incomplete Inflation/Deflation, Kinking)	13 (4.3%)	905 (0 – 1915)
Urination Impaired (Slow Stream, Split Stream, Voiding Difficulties or Obstructive Symptoms)	11 (3.7%)	239 (0 – 930)
Urogenital Erythema	10 (3.3%)	36 (0 – 320)
Joint Pain, Swelling, or Stiffness	9 (3.0%)	609 (1 – 1592)
Decreased Penile Sensation	7 (2.3%)	124 (0 – 214)
Urogenital Hematoma	7 (2.3%)	4 (0 – 25)
Abnormal Ejaculation (Delayed, Burning, or General Nonspecific Problems)	6 (2.0%)	409 (40 – 1797)
Infection	6 (2.0%)	216 (9 – 716)
Dysuria	5 (1.7%)	231 (2 – 684)
Penile Curvature	5 (1.7%)	144 (0 – 257)
Application Site Reaction (Wound Separation, Delay in Cutaneous Closure)	4 (1.3%)	14 (0 – 30)
Erosion/Extrusion (Pump/Cylinder)	4 (1.3%)	425 (72 – 1066)
Paresthesia	4 (1.3%)	490 (0 – 1897)
Urogenital Inflammation	4 (1.3%)	12 (0 – 27)
Adhesion of the Pump/Scrotum	3 (1.0%)	13 (10 – 19)
Device Malposition	3 (1.0%)	278 (43 – 574)
Device Migration (Pump/Cylinder)	3 (1.0%)	210 (40 – 548)
Transient Urinary Retention	3 (1.0%)	85 (3 – 248)
Urinary Frequency	3 (1.0%)	277 (99 – 409)
Weakness	3 (1.0%)	1072 (519 – 1592)
Abnormal Sexual Function	2 (0.7%)	239 (128 – 349)
Device Cylinder Aneurysm/Bulge	2 (0.7%)	945 (110 – 1780)
Dizziness	2 (0.7%)	929 (7 – 1850)
Dry Mouth	2 (0.7%)	1721 (1592 – 1850)
Hematuria	2 (0.7%)	902 (13 – 1791)

ADE	# Patient (%**)	Mean Onset Time in Days (Range in Days)
Low Grade Fever	2 (0.7%)	13 (7 – 18)
Memory Difficulties	2 (0.7%)	1318 (1107 – 1592)
Pelvic Pain	2 (0.7%)	270 (42 – 498)
Rheumatoid Arthritis	2 (0.7%)	281 (189 – 372)
Other	22 (7.0%)	N/A

\*\* Percentages based on total number of patients implanted (300).

**Table 2: Resolution of Adverse Device Effects**

ADE	Method of Resolution					
	# Patient (%)	Surgical	Medication	Other <sup>1</sup>	No Medical Intervention	Continuing
Urogenital Pain	160 (53.3%)	1%	31%	0%	68%	11%
Urogenital Edema	106 (35.3%)	0%	3%	11%	86%	2%
Urogenital Ecchymosis	30 (10.0%)	0%	0%	0%	100%	0%
Reservoir Encapsulation	19 (6.3%)	5%	0%	0%	95%	58%
Patient Dissatisfaction	18 (6.0%)	0%	0%	0%	100%	56%
Auto-Inflation	17 (5.7%)	0%	0%	0%	100%	35%
Mechanical Malfunction	13 (4.3%)	46%	0%	8%	46%	62%
Urination Impaired	11 (3.7%)	0%	64%	9%	27%	0%
Urogenital Erythema	10 (3.3%)	10%	30%	0%	60%	0%
Joint Pain, Swelling, or Stiffness	9 (3.0%)	0%	11%	11%	78%	67%
Decreased Penile Sensation	7 (2.3%)	0%	0%	0%	100%	72%
Urogenital Hematoma	7 (2.3%)	0%	0%	0%	100%	0%
Abnormal Ejaculation	6 (2.0%)	0%	17%	0%	83%	17%
Infection	6 (2.0%)	67%	33%	0%	0%	17%
Dysuria	5 (1.7%)	0%	60%	0%	40%	0%
Penile Curvature	5 (1.7%)	0%	0%	0%	100%	60%
Application Site Reaction	4 (1.3%)	0%	25%	0%	75%	25%
Erosion/Extrusion	4 (1.3%)	100%	0%	0%	0%	0%
Paresthesia	4 (1.3%)	0%	0%	0%	100%	50%
Urogenital Inflammation	4 (1.3%)	0%	50%	0%	50%	0%
Adhesion of Pump/Scrotum	3 (1.0%)	0%	0%	0%	100%	33%
Device Malposition	3 (1.0%)	67%	0%	0%	33%	0%
Device Migration	3 (1.0%)	100%	0%	0%	0%	0%
Transient Urinary Retention	3 (1.0%)	0%	0%	100%	0%	0%
Urinary Frequency	3 (1.0%)	0%	33%	0%	67%	67%
Weakness	3 (1.0%)	0%	0%	67%	33%	67%
Abnormal Sexual Function	2 (0.7%)	0%	0%	0%	100%	100%
Device Cylinder Aneurysm/Bulge	2 (0.7%)	50%	0%	0%	50%	50%
Dizziness	2 (0.7%)	0%	0%	0%	100%	50%
Dry Mouth	2 (0.7%)	0%	0%	0%	100%	100%
Hematuria	2 (0.7%)	0%	50%	0%	50%	50%
Low Grade Fever	2 (0.7%)	50%	50%	0%	0%	0%
Memory Difficulties	2 (0.7%)	0%	0%	0%	100%	0%
Pelvic Pain	2 (0.7%)	0%	0%	0%	100%	50%
Rheumatoid Arthritis	2 (0.7%)	0%	0%	0%	100%	100%
Other	22 (7.0%)	N/A	N/A	N/A	N/A	N/A

1 Other treatments included back brace, physical therapy, urine culture, ice packs, elevation, hot soaks, hot sitz bath, manual manipulation, patient education, fillforms and followers, foley catheter, ultrasound/CT scan, and cystoscopy.

The following "Other" adverse device effects (in alphabetical order) each occurred in less than 0.5% of the patients: Alopecia, Back Pain, Cellulitis, Depression, Diabetes Mellitus, Epigastric Pain, Eye Disorder, Eye Pain, Fecal Incontinence, Fibrosis, Glans Hypermobile Dorsally, Kidney Calculus, Libido Decreased, Migraine, Necrosis, Phimosis, Photosensitivity Reaction, Pump Fixation, Thickening of the Skin, Urinary Tract Infection, Urinary Urgency, and Vertigo.

The following adverse events (in alphabetical order) have been associated with the use of this product: Bleeding, Cavemosal Crossover, Exposure to Biohazardous Material, Granuloma Formation, Ischemia, Non-Rheumatoid Arthritis Immune-Related Tissue Disorders, Pain (which may be prolonged or severe), Perforation or Injury of Bladder, Corpus Cavernosum, Nerve, Tunica, or Urethra, Prolonged Procedure, Seroma, Thrombosis, Ulceration, Unretrieved Device Fragment, Vascular Compromise, Ventral Chordee, and Vessel Trauma.

There were 18 patient deaths during the course of the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

A total of 22 patients underwent revision surgeries in the 5 year study period. Information on device revisions is described in the Clinical Studies section.

## Clinical Studies

A clinical trial was undertaken to demonstrate that the AMS 700 product line provides an erection that is suitable for intercourse and has acceptable rates of surgical revision and of significant clinical events associated with the implantation and use of these devices. This trial included only devices **without** InhibiZone. This trial was also designed to demonstrate the implantation of these devices does not negatively impact the sexual satisfaction, psychological well-being, self-esteem or quality of life of patients who receive these devices. It was a prospective, multicenter cohort trial in which the patients served as their own control. The choice of device model implanted (ie, 700 CX, 700 CX Preconnect, 700 CXM, 700 Ultrex™, 700 Ultrex Plus) was at the discretion of the patient and implanting physician.

**NOTE:** The AMS 700 MS Pump was not available at the time the clinical study was conducted. However, based on the similarities between the AMS 700 MS Pump and the AMS 700 Inflate/Deflate Pump, the clinical results also apply to this new model.

**NOTE:** The AMS 700 LGX Preconnect was not available at the time the clinical study was conducted. However, since providing the AMS 700 LGX in preconnected form is not expected to affect the safety and effectiveness of the prosthesis, these clinical results also apply to the new model.

**NOTE:** The AMS 700 CXR was not available at the time the clinical study was conducted. However, based on the similarities between the AMS 700 CXR and the AMS 700 CXM models, the clinical results also apply to this new model.

**NOTE:** The AMS Conceal™ Low Profile Reservoir\* was not available at the time the clinical study was conducted. However, based on similarities between the AMS Conceal Low Profile Reservoir and the spherical reservoir, the clinical results also apply to this new model.

\* not available in all markets.

Three hundred male patients, over 21 years of age, were enrolled in this study. All patients with diagnosed organic erectile dysfunction were eligible for enrollment, if they did not present with a history of allergy/sensitivity to silicone, preexisting autoimmune or connective tissue diseases, or active urogenital infection.

All safety-related data, diagnoses, and health status evaluations were captured on detailed case report forms. The investigators' professional evaluation of the erections provided by the IPPs (Inflatable Penile Prostheses) after implantation and their suitability for intercourse was the primary efficacy endpoint. The number of surgical revisions performed and reported by the investigators was the primary safety endpoint. Patient self-evaluations on 4 validated outcome instruments were the secondary efficacy endpoints (concerning quality of life, self-esteem, and sexual satisfaction and functioning).

This clinical trial provided the following results through the 5 year evaluation for the first 126 patients to reach this postsurgical follow-up.

### Physician Assessment of Device Function

One hundred twenty-six devices were evaluated at the 5 year follow-up, of which 123 (97.6%) could be inflated. Of these 123 devices, all (100%) were determined to provide an erection suitable for intercourse. However, it should be noted that this analysis does not include the following information regarding device malfunctions: (i) 3 of the 123 devices found to be functioning properly at the 5 year exam were surgically revised prior to this exam to correct a mechanical malfunction and (ii) 3 additional devices not evaluated at the 5 year follow-up exam were also surgically revised due to mechanical malfunction. These cases of device revision are discussed further in the next section.

### Surgical Revisions

The incidence of revisions was evaluated in the 126 patients with follow-up out to 5 years, as well as 16 additional patients who experienced one or more revision surgeries and did not reach the 5-year follow-up exam. (A revision is considered any urogenital surgical intervention that is related to the function, placement, or site reaction to the implanted device.) Of these 142 patients, 22 (15.5%; 95% confidence interval = 21.5%) experienced a total of 26 revision surgeries, and 120 (84.5%) were not revised.

The average time to the first revision surgery was 15 months (ranging from 0.9 months to 60.1 months). Of the 26 revision surgeries, there were 5 revisions due to "infection," 2 for "infection/ Erosion," 2 for "Migration/Malposition," 2 for "Erosion," 2 for "Malposition," 7 for "Mechanical Malfunction," 2 for "Fibrous Capsular Complication," 2 for "Reimplantation Following Previous Revision," and 2 due to reasons listed as "Other." The "Other" reasons for revision surgeries included 1 for Cylinder kink/auto-inflation and 1 for Corporal body aneurysm. In 5 of these revision surgeries, no device components were explanted or replaced. The components were manipulated/repositioned but were not removed.

### Patient Evaluation of Quality of Life, Self-esteem, Psychological and Sexual Well-Being

In accordance with the study protocol, overall health-related quality of life (using the Medical Outcomes Study Health Survey, MOS-20), self-esteem (using the Rosenberg Self-esteem Scale), psychological well-being (using the Brief Symptom Inventory), and sexual functioning and satisfaction (using the Sexual History Form) were evaluated in patients through 2 years post-implantation. Throughout the 2-year follow-up period, patient quality of life, self-esteem, and psychological well-being were determined to be equivalent to the preimplant state. Sexual functioning and sexual satisfaction, on the other hand, was significantly improved over the preimplant state.

### Supplementary Clinical Information

Although it is not feasible to predict exactly how long an implanted penile prosthesis will function in a particular patient, American Medical Systems (AMS), Inc., has assembled a set of data on device removals and revisions to help gain insight into the product performance over time.

The following 2 tables provide an estimate of the long-term rates of device removals and revisions for Ultrex and CX models. The first data set comes from Patient Information Forms (PIFs) submitted to AMS by physicians for surgical procedures requiring parts replacement under AMS warranty (Table 3). All forms reporting devices implanted between January 1993 and December 2000 were included in a life table analysis that was used to calculate the revision rates for each category. Revision surgeries may not be reported to AMS. Therefore, the incidence of surgeries after original implant would likely be underreported, if one were to rely solely on the PIF data. AMS also assembled a second set of data directly from a retrospective review of physicians' medical records (Table 4). These medical records capture each surgery performed by that physician after original implant for any reason.

**NOTE:** These analyses did not include the AMS 700 with MS Pump, snapcone cylinders, or parylene coated components.

**NOTE:** These analyses did not include the AMS Conceal Low Profile Reservoir.

### PIF Study

Table 3: Revision rates based on PIF data:\*

Reason for Removal or Replacement Surgery	700 CX Revision Rate (5 YRS) n=12,080	700 Ultrex Revision Rate (5 YRS) n=20,438
Mechanical Revision	5.6%	4.8%
Removal for Infection	2.2%	2.0%
Removal for Erosion	1.4%	1.0%
Migration/Malposition of Component	0.5%	0.5%
Cylinder Aneurysm	0.5%	0.5%
Other Reasons	3.6%	2.9%
<b>ALL REASONS</b>	<b>11.6%</b>	<b>9.8%</b>

\* Interpretations of the PIF data may be limited by a number of factors:

- Statistics are based only on surgery data voluntarily reported to AMS by hospitals and physicians in the U.S. as part of the AMS Product Replacement. Because surgeries may not be reported to AMS, the number of patients implanted and the incidence of removal/replacement surgery may actually be higher.
- These statistics pertain only to the incidence of removal/replacement surgery and not to the current functioning of devices that have not been removed.

**NOTE:** The amount for "ALL REASONS" is lower than the total of the individual percentages due to the incidence of multiple reasons for removal/replacement surgery.

## Medical Records Study

**Table 4: Revision rates based on medical records data\***

Reason for Removal or Replacement Surgery	700 CX Revision Rate (5 YRS) n=512	700 Ultrex Revision Rate (3 YRS) n=155
Mechanical Revision	9.4%	3.2%
Removal for Infection	2.9%	2.0%
Removal for Erosion	2.4%	1.6%
Migration/Malposition of Component	2.8%	0.7%
Cylinder Aneurysm	0.6%	0.0%
Other Reasons	2.1%	4.0%
<b>ALL REASONS</b>	<b>18.1%</b>	<b>10.4%</b>

\* Interpretations of the medical record data may be limited by several factors:

- These percentages reflect the known revisions performed by the original implant physicians.
- These statistics pertain only to the incidence of removal/replacement surgery and not to the current functioning of devices which have not been removed.

**NOTE:** The amount for "ALL REASONS" is lower than the total of the individual percentages due to the incidence of multiple reasons for removal/replacement surgery.

## Antibiotic Information

The antibiotics present in InhibiZone, minocycline HCl and rifampin (rifampicin), are well characterized and have been in use for years. The dosage present on the penile prostheses is intended to act on organisms that attempt to colonize the device. The AMS 700 components are treated with very low levels of antibiotics. AMS provides numerous completed configurations of the AMS 700 to individualize treatment; however, a complete device (reservoir, pump, and 2 cylinders), regardless of configuration, contains  $\leq 33$  mg rifampin and  $\leq 12$  mg minocycline HCl which represents less than 2% of oral dose exposure for a complete course of rifampin or minocycline HCl with the maximum dose calculated for the most common device configuration's average concentration plus (1) standard deviation.

The following in vitro data describe the Zones of Inhibition provided by samples from InhibiZone treated devices.

**Table 5: In vitro zones of inhibition for device samples\* with InhibiZone treatment**

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
<i>Staphylococcus epidermidis</i>	22.6	2.9	21
<i>Staphylococcus aureus</i>	17.5	5.0	25
<i>Escherichia coli</i> **	6.5	2.6	24
<i>Enterococcus faecalis</i> **	4.8	6.7	21
<i>Candida albicans</i> **	0.1	0.4	21
<i>Proteus mirabilis</i> **	0.6	1.0	17

\* Obtained using standardized KRT test samples containing 12  $\mu$ g minocycline and 26  $\mu$ g rifampin (rifampicin)

\*\* The isolates tested were not susceptible to rifampin (rifampicin) and/or minocycline control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and 5 rabbits were implanted subcutaneously with 6 control samples each. One rabbit received 3 test samples and 3 control samples. The test samples were portions of an InhibiZone treated AMS 700 Pump and the control samples were portions of a standard AMS 700 Pump without InhibiZone. All samples were soaked in a 103-104 colony-forming units solution of *Staphylococcus aureus*, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

## Antibiotic Human Clinical Experience

A post marketing study was conducted to evaluate the effect of the antibiotic surface treatment on reducing the incidence of penile implant infections, which consisted of a comprehensive review of more than 43,000 patients in the AMS PIF database. This study shows a significant improvement in the rate of revisions due to infection for patients with original AMS 700 InhibiZone (IZ) implants or revision AMS 700 IZ implants compared to those who received AMS 700 devices without the IZ treatment. Further, patients with diabetes who received original AMS 700 IZ implants also demonstrate significant improvements in the rate of revisions due to infections, compared to those who received AMS 700 devices without the IZ treatment. The data set included PIF data for patients who underwent an AMS 700 implant surgery in the United States from May 1, 2001, through September 1, 2007. The most common reason for implantation of an AMS 700 device was recorded as "organic" on the PIFs (Table 6). Post prostatectomy, vascular disease and diabetes were also commonly reported as the etiology leading to the need for the AMS 700 implants.

**Table 6: Etiology for AMS 700 Implant Patients**

Description	IZ		Non-IZ		Mixed	
	#	%	#	%	#	%
Original AMS 700 Implant						
Diabetes Mellitus*	4,977	17.7%	627	17.7%	177	18.1%
Peyronie's Disease	2,553	9.1%	328	9.2%	102	10.4%
Post Prostatectomy	7,215	25.7%	819	23.1%	242	24.7%
Vascular Disease	4,037	14.4%	576	16.2%	120	12.3%
Other	1,107	3.9%	151	4.3%	50	5.1%
Organic	8,185	29.2%	1,046	29.5%	288	29.4%
Total	28,074	100.0%	3,547	100.0%	979	100.0%
Revision or Replacement AMS 700 Implant						
Diabetes Mellitus	1,979	16.2%	349	15.4%	63	15.2%
Peyronie's Disease	1,243	10.2%	208	9.2%	40	9.6%
Post Prostatectomy	2,016	16.5%	330	14.5%	68	16.4%
Vascular Disease	2,566	21.0%	599	26.4%	101	24.3%
Other	1,034	8.5%	184	8.1%	34	8.2%
Organic	3,378	27.7%	602	26.5%	109	26.3%
Total	12,216	100.0%	2,272	100.0%	415	100.0%
TOTAL	40,290		5,819		1,394	

**NOTE:** None or more than one etiology could be selected for each patient; therefore, the numbers of patients within each subgroup by etiology above may not be identical to the number of patients in each life table analysis.

\* The diabetes subgroup with an original implant was identified and analyzed as a specific high-risk subgroup.

Three patient populations were analyzed in this PIF study: the first group included all patients with an AMS 700 original implant between May 1, 2001, and September 1, 2007. A subgroup of these AMS 700 original implant patients included all patients with diabetes. The third group included all patients with an AMS 700 revision implant (removing or replacing all or part of an earlier implanted device) between May 1, 2001, and September 1, 2007. An event was recorded for any revision surgery with a surgery reason recorded as infection, mechanical malfunction, fluid loss, erosion, patient dissatisfaction, or other. Revision surgeries were tracked through December 1, 2007, when the reporting period was closed. When multiple reasons were listed for 1 revision surgery (eg, infection, device malfunction, erosion, etc), the individual revision surgery was counted and included in the data for each of the indicated revision reasons. Thus, all revisions that report some presence of infection will be captured in the infection category presented in this study (even if multiple other reasons for the revision surgery are listed).

Life tables were constructed and the logrank test was used to compare the device survival distributions between the Treatment (InhibiZone) and Control (non-InhibiZone) groups. All significance tests were performed at an alpha level of 0.05. For original implants, device survival distributions were defined as the time between the implantation of the original device in the study period and the first replacement, revision, or removal of any component of the device due to infection, malfunction, fluid loss, erosion, patient dissatisfaction, or other reason (the Original Implant Group). For revision surgery patients, device survival distributions were created for the time between the implantation of the revision device and the next revision device due to infection, malfunction, fluid loss, erosion, patient dissatisfaction, or other reason (the Revision Surgery Group). The study also included a subgroup of patients who received a penile implant after spinal cord injury. However, this aspect of the study continues to be unknown, since there were too few patients in this group to be able to determine valid conclusions.

### Original Implant Patient Group

Among the original implant patients, a total of 31,328 (98.7%) patient devices did not have any events reported due to infection and were censored during the timeframe of this report. A total of 338 (1.2%) device events were due to infection in the IZ treated group while 84 (2.5%) were due to infection in the non-IZ treated group. Life table analyses compared the rates of revision surgery due to infection for original IZ or non-IZ treated implants. In this device survival analysis with infection as the reason for the revision surgery, InhibiZone treated AMS 700 devices had a significantly better survival than non-InhibiZone treated devices among original implant patients ( $p < 0.0001$ ) (Figure 1).


Figure 1: Device survival—freedom from revision due to infection among original implant patients

### Diabetic Original Implant Patient Group

Among the subgroup of original implant patients who had diabetes, a total of 5604 diabetic patients included 4977 (88.8%) with IZ treated devices and 627 (11.2%) with non-IZ treated devices. A total of 5508 (98.3%) diabetic patient devices did not have any events due to infection during the timeframe of this report. Only 71 (1.4%) infection device events were in the IZ treated group while 25 (4.0%) events were in the non-IZ treated group. Life table analyses compared rates of revision surgery due to infection for diabetic patients. InhibiZone treated devices had significantly better survivals than non-InhibiZone treated devices among diabetic patients ( $p < 0.0001$ ) (Figure 2).



Figure 2: Device survival—freedom from revision due to infection among diabetic patients

### Revision Surgery Implant Group

Among the revision surgery implant patients, a total of 10,769 (97.3%) patient devices did not have any events due to infection and were censored during the timeframe of this report. Only 229 (2.5%) infection device events occurred in the IZ treated group while 66 (3.7%) infection events occurred in the non-IZ treated group. Life table analysis compared the rates of revision surgery due to infection for IZ or non-IZ treated implants among revision surgery implant patients. In this device survival analysis with infection as the reason for the revision surgery, InhibiZone treated AMS 700 devices had a significantly better survival than non-InhibiZone treated devices among revision surgery implant patients ( $p=0.0252$ ) (Figure 3).

**Figure 3: Device survival—freedom from revision due to infection among revision surgery implant patients**

In conclusion, this human clinical evidence from the InhibiZone PIF study indicates the revisions due to infection were significantly lower among patients who received an AMS 700 InhibiZone treated device when compared to a nontreated AMS 700 device. This reduced rate of revisions subsequent to infection for the AMS 700 InhibiZone treated device when compared to a nontreated AMS 700 device was reported for original as well as revision implants and included the high-risk subgroup of original implant patients who had diabetes. The study supports the following conclusions:

- No significant differences were found between InhibiZone treated and untreated AMS 700 penile prostheses for overall rates of revision due to mechanical malfunction, fluid loss, erosion, patient dissatisfaction, or all other reasons.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection among penile prosthesis implants.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection among penile prosthesis implants in diabetic patients.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection in original and revision penile prosthesis implants.

### Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation. The risks, benefits, and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders and a diagnostic work-up, should be a part of the patient decision-making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal.

Implantation of a penile prosthesis may result in penile shortening, curvature, or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations. Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring, scrotal deformity, pump bulge in the scrotum, lack of concealability, and other possible adverse events. Patients should also be aware that penile prostheses are not considered to be lifetime implants.

Improper implantation of a penile prosthesis may not provide rigidity to the glans, which may result in a floppy glans and may result in a lack of rigidity of the corpus spongiosum. Penile flaccidity may be less than it was before implantation.

Patients who undergo revision surgery may notice a change in the character of their erection compared to their previous implant, which may include differences in sensation, length, girth, rigidity, and/or flaccidity.

It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (see Silicone Information section).

### Silicone Information

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate "allergic-like" symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

There are reports of malignant tumor formation in laboratory animals only, not humans, associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the prostheses in the AMS 700. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

## Magnetic Resonance Imaging (MRI) Important Safety Information

### MR Conditional

The AMS 700 Series Three-Piece Inflatable Penile Implant was determined to be MR Conditional.

Nonclinical testing demonstrated that the AMS 700 Series Three-Piece Inflatable Penile Implant is MR Conditional. A patient with this device can be safely scanned immediately after placement under the following conditions:

#### Static Magnetic Field

-Static magnetic field of 3 T or less

-Maximum spatial gradient magnetic field of 720 Gauss/cm or less

#### MRI-Related Heating

In nonclinical testing, the AMS 700 Series Three-Piece Inflatable Penile Implant produced the following temperature rise during MRI performed for 15 minutes of scanning (ie, per pulse sequence) in the 3-Tesla (3 T/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) magnetic resonance (MR) system:

**Highest temperature change** +1.9°C

Therefore, the MRI-related heating experiments for the AMS 700 Series Three-Piece Inflatable Penile Implant at 3 T using a transmit/receive radio frequency (RF) body coil at an MR system reported whole body averaged specific absorption rate (SAR) of 2.9 W/kg (ie, associated with a calorimetry measured whole body averaged value of 2.7 W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.9°C.

#### Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the AMS 700 Series Three-Piece Inflatable Penile Implant. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	6,244 mm <sup>2</sup>	1,589 mm <sup>2</sup>	10,295 mm <sup>2</sup>	2,779 mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

## Inventory Returns and Product Replacement Information

### In the United States

Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the PIF.

Follow **all** of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Product Warranty Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

### Outside the United States

Customers outside of the United States should contact their local AMS Representative prior to returning any product.

This document is written for professional medical audiences. Contact American Medical Systems for lay publications.

American Medical Systems periodically updates product literature. If you have questions about the currency of this information, contact American Medical Systems.

## How Supplied and Storage

**WARNING:** Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, or if labeling information appears missing, illegible, or compromised, call your AMS representative for clarification, and do not use the AMS 700 device until integrity of the product has been verified.

**For single patient use only.** Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the AMS 700 device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the AMS 700 device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the AMS 700 device may lead to injury, illness, or death of the patient. After use, dispose of the AMS 700 device and packaging in accordance with hospital, administrative, and/or local government policy.

Store product in a clean, dry, dark area at room temperature.

Do not use device beyond the indicated expiration date.

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©2017 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are the property of the respective owners.
Boston Scientific (Australia) Pty Ltd
American Medical Systems, Inc.
U.S.A.
Europe B.V.