AMS 800TM

Artificial Urinary Sphincter for Male and
Female Patients

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WARRANTY18

AMS 800™

Artificial Urinary Sphincter for Male and Female Patients

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using steam or an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.


Figure 1. AMS 800 AUS implanted in a male patient


Figure 2. AMS 800 AUS implanted in a female patient

DEVICE DESCRIPTION

System Description

The AMS 800 Artificial Urinary Sphincter (AUS) is an implantable, fluid filled device (Figures 1, 2). The device consists of three components: an occlusive cuff (cuff), a control pump, and a pressure regulating balloon (PRB). Kink-resistant tubing (KRT) connects these components to form a closed system. The AMS 800 AUS simulates sphincter function by applying circumferential pressure to aid in opening and closing the urethra. When the cuff is inflated, the urethra is closed and urine stays in the bladder (Figure 3). When the patient wishes to void, he or she cycles the device by squeezing and releasing the pump several times, moving the fluid from the cuff to the PRB (Figure 4). The cuff deflates and urine passes through the open urethra. Pressure from the PRB pushes fluid back into the cuff, occluding the urethra and restoring continence (Figure 5).

The AMS 800 AUS cuffs and pumps are available with InhibiZone™, an antibiotic treatment of rifampin (rifampicin) and minocycline hydrochloride (minocycline HCl)¹.

¹ InhibiZone Antibiotic Treatment is not available in all markets.


Figure 3. Urethra closed by cuff


Figure 4. Opening of cuff and urination by manipulation of pump


Figure 5. Automatic cuff inflation and urethral coaptation

Occlusive Cuff

The cuff (Figures 6, 7) occludes the urethra by applying pressure circumferentially. It is made of silicone elastomer and is available in thirteen sizes in length: 3.5 cm, 4.0 cm, 4.5 cm, 5.0 cm, 5.5 cm, 6.0 cm, 6.5 cm, 7.0 cm, 7.5 cm, 8.0 cm, 9.0 cm, 10.0 cm, and 11.0 cm. All cuffs are approximately 1.8 cm wide when deflated. The proper size to be used in the patient is determined by the circumference of the tissue around the urethra intraoperatively. The cuff tubing is joined to the control pump tubing with a connector.


Figure 6. Occlusive cuff


Figure 7. Occlusive cuff, side view

Control Pump

The control pump (Figure 8) is implanted in the soft tissue of the scrotum or labium. It is made of silicone elastomer and is approximately 1.3 cm wide and 3.5 cm long. The upper part of the pump (the valve block) contains the resistor and valves needed to transfer fluid between the components. The bottom half of the pump is a bulb that the patient squeezes and releases to transfer fluid out of the cuff to void. The deactivation button is located on the valve block. Pressing the deactivation button prevents fluid from being transferred between the components. This feature allows the physician to leave the cuff deflated for transurethral procedures, and during the postoperative period.


Figure 8. Control pump

Pressure Regulating Balloon (PRB)

The pressure regulating balloon (Figure 9), implanted in the prevesical space, controls the amount of pressure exerted by the cuff. It is made of silicone elastomer and is provided in the following three pressure ranges:

- 51-60 cm H₂O (5.0-5.9 kPa)
- 61-70 cm H₂O (6.0-6.9 kPa)
- 71-80 cm H₂O (7.0-7.8 kPa)

The PRB tubing is joined to the pump tubing by a connector.


Figure 9. Pressure regulating balloon

Kink-Resistant Tubing (KRT)

Kink-resistant tubing (Figure 10) is reinforced with suture material that is either clear or black. The tubing is color-coded to help surgeons make the correct connections between components.


Figure 10. Pump with clear and black KRT

- Clear tubing connects the pump to the cuff.
- Black tubing connects the pump to the PRB.

KRT is connected with connectors, of which there are two types: Quick Connect Window Connectors and Suture-Tie Connectors.

AMS 800™ Accessory Kit

The AMS 800 Accessory Kit contains the materials needed for one implant procedure. The kit contains the following disposable items:

Cuff Sizer

- One disposable Cuff Sizer, which aids in selection of a cuff

Special Purpose Needles

- Two 15 gauge disposable, Blunt Tip Needles
- Two 22 gauge disposable, Blunt Tip Needles

Hemostat Shods

- Two 30 cm lengths of blue silicone tubing

Tubing Connection Accessories

- Three Straight Quick Connect Window Connectors
- Two Right Angle Quick Connect Window Connectors
- One Y Quick Connect Window Connector (a 3-way connector)
- One Collet Holder with eight Collets
- Three Straight Suture-Tie Connectors
- Two Right Angle Suture-Tie Connectors
- One Y Suture-Tie Connector (3-way connector)

InhibiZone™ Antibiotic Treatment

The InhibiZone antibiotic treatment process is intended to elute antibiotics from AMS 800 AUS component surfaces when they are exposed to a warm, moist environment. The dosage present on the AMS 800 AUS is intended to act on organisms that attempt to colonize the surface of the device. InhibiZone is a formulation of minocycline hydrochloride and rifampin (rifampicin). BSC provides various component configurations of the AMS 800 AUS to individualize treatment. However, while the AMS 800 AUS pressure regulating balloon (PRB) is not InhibiZone treated, a complete device (PRB, pump, and one or two cuffs), regardless of configuration, contains ≤ 6.5 mg rifampin and ≤ 8 mg minocycline HCl. This represents less than 2% of oral dose exposure for a complete course of rifampin or minocycline HCl with the maximum dose calculated from the means and 95% tolerance intervals.

In vitro data presented in Table 1 are available, but their clinical significance is unknown. No clinical studies have been performed to evaluate the effect of the antibiotic treatment on reducing the incidence of artificial urinary sphincter implant infections.

Table 1. In Vitro Zone of Inhibition for Device Samples* with InhibiZone Treatment

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
Staphylococcus epidermidis	22.6	2.9	21
Staphylococcus aureus	17.5	5.0	25
Escherichia coli**	6.5	2.6	24
Enterococcus faecalis**	4.8	6.7	21
Candida albicans**	0.1	0.4	21
Proteus mirabilis**	0.6	1.0	17

*obtained using standardized kink-resistant tubing test samples containing approximately 12 µg minocycline HCl and 26 µg rifampin
**the isolates tested were not susceptible to rifampin and/or minocycline HCl control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and five rabbits were implanted subcutaneously with 6 control samples each. One rabbit received three test samples and three control samples. The test samples were portions of an InhibiZone treated AMS 700™ Inflatable Penile Prosthesis pump and the control samples were portions of a standard AMS 700 pump without InhibiZone. The AMS 700 pumps used in the AMS 700 device are similar to AMS 800 AUS pumps used in the AMS 800 Artificial Urinary Sphincter device in regards to material composition, adhesive, and InhibiZone application process. All samples were soaked in a 103-104 CFU solution of Staphylococcus aureus, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

Existing prophylactic antibiotic protocols should be maintained as determined by the physician and/or institution.

User Information

Only physicians with advanced training in urology, and management of patients in need of an artificial urinary sphincter, shall implant the AMS 800 AUS. This manual is not intended to be a complete reference.

Contents

The following is a list of the AMS 800™ AUS components and accessories.

- (1) Occlusive Cuff
- (1) Control Pump
- (1) Pressure Regulating Balloon
- (1) AMS 800 Artificial Urinary Sphincter Accessory Kit:
 - (3) Straight Quick Connect Window Connectors
 - (2) Right Angle Quick Connect Window Connectors
 - (1) Y Quick Connect Window Connector
 - (3) Straight Suture-Tie Connectors
 - (2) Right Angle Suture-Tie Connectors
 - (1) Y Suture-Tie Connector
 - (8) Collets
 - (1) Collet Holder
 - (1) Cuff Sizer
 - (2) 22 gauge Blunt Tip Needles
 - (2) 15 gauge Blunt Tip Needles
 - (2) 30 cm lengths of tubing

Materials of Animal Origin

The following components of the AMS 800 Artificial Urinary Sphincter contain materials of animal origin:

- InhibiZone™ treated components may contain a non-viable derivative of porcine origin or a bovine milk derivative.
- The Quick Connect Window Connectors and the collets contain a non-viable tallow derivative.

INTENDED USE

The AMS 800 AUS inflated cuff applies circumferential pressure on the urethra to increase urethral resistance to urinary flow.

INDICATIONS FOR USE

The AMS 800 AUS is used to treat urinary incontinence due to reduced urethral/bladder outlet resistance (intrinsic sphincter deficiency) in males and females.

CONTRAINDICATIONS

- This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
- This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
- This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
- The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin or to minocycline HCl or other tetracyclines.
- The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to aggravate this condition.

WARNINGS

Surgery-Related

- Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component malposition. Failure to evaluate and promptly treat erosion may result in substantial worsening of the condition, leading to infection and/or loss of tissue.
- The implanter should check that there is an adequate amount of bulbospongiosus muscle to surround and support a bulbous urethral cuff implant. Thinner spongiosum typically occurs toward the distal end of the bulbous urethra, and implantation of the cuff where the spongiosum is thin increases the chance of erosion and other complications. This warning is especially important for double cuff implants, where the second cuff is placed distal to the first implanted cuff.

Patient-Related

- Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection. Infection that fails to respond to antibiotic therapy may result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.

- Poor bladder compliance or a small, fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
- Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia, or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
- This device contains solid silicone elastomers. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
- Previous patient history of an adverse reaction to any radiopaque solution precludes their use as a filling medium for the prosthesis. Instead, sterile normal saline should be used to fill the device.

Post-Procedural

- Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device, to prevent potential damage to the urethra or the AMS 800 AUS.
- Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device. Surgical complications, or physical or psychological complications (such as any progressively degenerative disease), if they occur, may affect the use of the system. Continued improper use may lead to revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
- Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in outflow obstruction. Mechanical events should be evaluated carefully by the treating physician, and the patient should consider risks and benefits of treatment options, including removal or revision surgery.
- If a hypersensitivity reaction develops to a device treated with InhibiZone, the cuff and pump should be removed and the patient treated appropriately.
- Female patients with persistent incontinence should be evaluated to rule out vesicovaginal fistula, which may have resulted from an unrecognized iatrogenic injury.

PRECAUTIONS

Surgery-Related

- Improper cuff sizing, improper balloon selection, or other causes may result in tissue erosion, migration of components, or continued incontinence.
- Component migration can occur if the cuff is sized improperly, if the pump or balloon is not positioned correctly, or if the tubing lengths are incorrect. Migration can result in pain, complications, device malfunction, and surgical revision.
- Unsuccessful outcomes have been reported due to improper surgical technique, such as anatomical misplacement of components, improper sizing and/or filling of components, as well as tubing kinks. Product wear, component disconnection, or other mechanical problems may lead to malfunctioning of the components and leakage of fluid.
- Tubing kinks may result if the connecting tubing is cut to an improper length during the implant procedure.

Device-Related

- Quick Connect Window Connectors should not be used in revision procedures involving previously implanted component tubing. Suture-Tie Connectors should be used in revision procedures. The Quick Connect system may be used in revision surgeries when all previously implanted components are removed and replaced with new components.
- The stainless steel tubing plugs in the Deactivation Package contain nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may occur.

In the event of large pressures within the bladder, automatic pressure relief that normally occurs with the device would be prevented. Cycling the device can relieve the outflow obstruction.

Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the device, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.

Release of the deactivation valve may require greater pressure than that used to cycle the device.

- System pressure changes may occur over time if the balloon is filled with radiopaque solution of incorrect concentration. Follow the instructions in the Filling Solutions section to prepare the radiopaque solution in the correct concentration.

Patient-Related

- The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.
- Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff.
- Trauma or injury to the pelvic, perineal, or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the device.
- Female patients of child-bearing age must be forewarned that pregnancy is acceptable, but caesarean section may be indicated to minimize the risk of damage to the bladder neck and its surrounding cuff. For those patients who become pregnant, deactivation of the device during the third trimester is recommended to reduce the risk of erosion. Patients contemplating pregnancy should consider delaying implantation.
- Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5 cm) generally exceeds 28F. Additional clearance is required to accommodate the patient’s urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician’s assessment to determine its impact on sizing.

InhibiZone™-Related

- InhibiZone does not replace normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
- Components with InhibiZone are treated with a combination of rifampin and minocycline HCl. The contraindications, warnings, and precautions regarding these antimicrobial agents apply, and should be followed for the use of this device, although systemic levels of minocycline HCl and rifampin in patients receiving this device are unlikely to be detected.
- Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin and minocycline HCl can cause additional stress on the hepatic and renal systems.
- Patients who receive a device with InhibiZone and are also taking methoxyflurane should be carefully monitored for signs of renal toxicity.
- Patients who receive a device with InhibiZone and are also taking warfarin should have their prothrombin time monitored because tetracyclines have been reported to slow coagulation.
- Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid, and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin.
- Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, or acetone, or other non-polar solvents. These solvents may remove the antibiotics from the device.
- InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.

PATIENT COUNSELING INFORMATION

Although the prosthesis is designed to close the urethra, some patients continue to have a degree of incontinence after this procedure.

Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcomes of the implantation. The risks, benefits, and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process. A discussion concerning the patient’s preference for which side of their body the pump should be placed on may be appropriate.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal. Patients should also be aware that the AMS 800™ is not considered to be a lifetime implant.

It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device.

Patients should be counseled on expected postoperative pain including severity and duration. Patients may experience pain when the device is activated in the postoperative period and during periods of initial use. Cases of chronic pain associated with the presence of the AMS 800 AUS device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention.

Silicone Information

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Solid silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

There are reports of malignant tumor formation in laboratory animals only associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS 800 AUS. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

MAGNETIC RESONANCE IMAGING (MRI)

Important Safety Information

Non-clinical testing has demonstrated that the AMS 800 AUS is MR Conditional. The device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 Tesla (T) or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the AMS 800 AUS produced the following temperature rise during MRI performed for 15 minutes of scanning (i.e., per pulse sequence) in the 3T magnetic resonance (MR) system (3T/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI):

Highest temperature change +2.0°C.

Therefore, the MRI-related heating experiments for the AMS 800 AUS at 3T using a transmit/receive radio frequency (RF) body coil reported a whole body averaged specific absorption rate (SAR) of 2.9 W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7 W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +2.0°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the AMS 800 AUS. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Table 2. Artifact Information

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	5,800 mm ²	1,956 mm ²	6,096 mm ²	2,650 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

ADVERSE EVENTS

Adverse Events – Males

A prospective clinical study was conducted to demonstrate the safety and efficacy of the AMS 800 AUS. A total of 87 patients were enrolled in the study and 85 patients were implanted with the device. The total device exposure was 143.1 patient years with a mean of 1.7 patient years and a range of 0.03 to 3.4 patient years.

During the study, 26 patients experienced 43 device related adverse events. There were 3 patient deaths during the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

Table 3 lists the device related adverse events reported during the study. This trial involved only devices without InhibiZone™.

Table 3. AMS 800™ Prospective Clinical Study – Device Related Adverse Events

Adverse Event Category				Interventions*		
	Total Events	Patients with Adverse Event	Events Resolved	None Reported	Medical**	Surgical
Impaired Device Function	7	6	4	2	2	4
Pain/Discomfort	6	5	4	3	3	1
Delayed Wound Healing	5	5	5	2	3	0
Bladder Spasms	2	2	0	0	2	0
Difficult Activation	2	2	2	1	1	0
Migration	3	3	1	2	0	1
Tissue Erosion	2	2	2	0	0	2
Difficult Deactivation	1	1	1	0	1	0
Infection	2	2	2	0	0	2
Recurrent Incontinence	3	3	3	1	0	2
Fistula Formation	1	1	1	0	0	1
Hematoma	1	1	1	0	1	0
Swelling	2	2	2	0	2	1
Hydrocele	1	1	1	0	1	1
Tissue Erosion/Infection	1	1	1	0	0	1
Patient Dissatisfaction	1	1	1	0	0	1
Positional Incontinence	1	1	0	1	0	0
Wound Infection	1	1	1	0	1	0
Urinary Retention	1	1	1	0	1	0

*Events may have been addressed with more than one type of intervention.

** Medical interventions included: medication, education, frequent device deactivation, dressing changes and catheterization.

For the 85 patients implanted with the device followed under the prospective clinical study, 14 patients (16.5%) experienced a total of 15 revisions up to 24 months following implant. Information on device revisions is described in the Clinical Studies – Males section.

Adverse Events – Males and Females

Adverse events potentially related to the implantation and use of the AMS 800 AUS include: adhesion, allergic reaction, contracture, deep vein thrombosis, dehiscence, dysuria, edema, erosion, exposure to biohazardous material, fibrosis, fistula, foreign body/unretrieved device fragment, hematoma, hematuria, hemorrhage/bleeding, herniation of the device through tissue, infection, limited urethral coaptation (may be due to device leak, sizing, malfunction, placement, or other causes), malfunction, malposition of components, mechanical malfunction or mechanical difficulty, migration of component, overactive bladder, pain/discomfort, prolonged procedure, unintended surgical damage (perforation or injury to the bladder, urethra, nerves, vessels, or other local structures), urethral atrophy, urethral stricture, urinary retention, and urinary urgency.

Adverse Events – Males

Additional adverse events reported specifically for males include: bladder spasms, delayed wound healing, hydrocele, and patient dissatisfaction.

Adverse Events – Females

Additional adverse events reported specifically for females with additional etiologies include: device damage – unrelated procedure, early post-operative morbidity, incisional hernia, phlebitis, pulmonary embolus, upper urinary tract deterioration, and wound seroma.

CLINICAL STUDIES

Clinical Studies – Males

A prospective, multi-center, non-randomized clinical study was undertaken to demonstrate that the AMS 800 AUS can be surgically implanted without serious adverse sequelae, provides an acceptable level of continence, and enhances quality of life. Each patient served as their own control. Efficacy data and safety data related to adverse events, revision surgery, diagnoses, and health status evaluations were captured on case report forms. Patient self-evaluations related to health status and non-illness specific quality of life were measured on two validated outcome instruments. Patient and physician assessments of continence were measured on a recognized, standardized, non-validated instrument.

Eighty-seven (87) male patients were enrolled in the study of which 85 patients were implanted with the device during the study. The patients ranged in age from 48.5 to 96.1 years at time of enrollment. The majority of the participants were White (97.70%), 1.15% were African American, and 1.15% were Asian. Patients available at the follow-up intervals were 6-months (n=67), 12-months (n=60), 18-months (n=55), and 24-months (n=41). Patients diagnosed with intrinsic sphincter deficiency (ISD) resulting from prostate surgery were eligible for enrollment. Patients with a history of allergy/sensitivity to silicone, pre-existing autoimmune or connective tissue disease, or active urogenital infection were excluded from the study.

This trial involved only devices without InhibiZone.

Endpoints

The primary effectiveness endpoint evaluated the effect of the prosthesis on patient quality of life using the *Incontinence Impact Questionnaire*, an incontinence-specific quality of life questionnaire. The primary safety endpoint evaluated the five-year revision-free rate using a Bayesian hierarchical model. The safety endpoint was a five-year revision-free rate equivalent to 75% using a 10% delta with a two-sided 95% lower bound greater than 65%.

Incontinence Impact Scores

The primary effectiveness endpoint was a reduction in the Incontinence Impact Score from pre- to post-implant status. Incontinence impact was measured pre- and post-implant at 6, 12, 18, and 24 months. Thirty-nine (39) patients answered the *Incontinence Impact Questionnaire* (IIQ) at 24-month follow-up. The IIQ is a 30-item, self-administered questionnaire designed to assess the impact of urinary incontinence on several subscales including physical, emotional, and social. The IIQ used in the study was developed from a validated instrument¹. The mean pre-implant score was significantly higher ($p<0.0001$) than the mean scores at all follow-up visits. Therefore, the impact of incontinence was reduced for patients following AMS 800 AUS implantation and the primary objective was met.

Physician and Patient Assessment of Continence

Physician assessed continence was 63.6% dry and 34.1% required some additional protection at one-year follow-up (n=43). At two-year follow-up (n=30), 73.3% were dry and 23.3% required some additional protection. Patient assessed continence was 61.7% dry and 36.7% required some additional protection at one-year follow-up (n=60). At two-year follow-up (n=41), 65.9% were dry and 31.7% required some additional protection. No significant difference existed between the physicians' assessment and the patients' assessment of their incontinence.

Patient Evaluation of Health Status and Self-esteem

General Quality of Life as measured by the *Health Status Questionnaire* and the *Rosenberg Self-esteem Questionnaire* was evaluated at pre- and post-implant at 6, 12, 18, and 24 months. Thirty-eight (38) patients answered the *Health Status* and *Rosenberg Self-esteem Questionnaires* at 24-month follow-up. The self-administered *Health Status Questionnaire*² was used to assess non-illness specific parameters such as physical functioning, social functioning, energy/fatigue, pain, health perception, and emotional problems. A high score indicates that overall health was perceived to be high. The mean score was 596 at pre-implant and 612 at two-year follow-up. No significant difference in health status scores was observed during the study. The self-administered *Rosenberg Self-esteem Questionnaire*³ was used to assess changes in patient self-esteem. The range of possible scores is 0-6, with a score of 6 indicating high self-esteem. The mean score at implant was 3.5 and at two-year follow-up was 4.1. The increase in mean score indicates a more positive self-esteem following implantation of the AMS 800 AUS. The device did not have an adverse effect on sexual function. Some patients with improved continence following implantation also reported increased sexual activity. The positive impact of the device on patient's lives measured in the clinical study is consistent with results obtained by other authors.^{4,5,6}

Surgical Revisions

A revision is a surgical intervention related to the function, placement, or site reaction to the implanted device. For the 85 patients implanted with the device followed under the prospective clinical study, 14 patients (16.5%) experienced a total of 15 revisions up to 24 months following implant. One patient had two revisions. Three (3) revisions were due to mechanical malfunction. Two (2) revisions were due to recurrent incontinence. Two (2) revisions were due to erosion. Two (2) revisions were due to infection. One (1) revision each (total = 6) was due to migration, pain, erosion/infection, persistent incontinence/patient dissatisfaction, recurring incontinence/malfunction, and infection/pain/urethrocutaneous fistula. Multiple reasons were provided for some revisions. Four of the 14 patients who experienced revisions elected to have the device removed and 10 elected to have the device replaced. The probability of remaining revision-free 24 months following implant was 79.5% (95% CI with 95% lower confidence bound 69.8%) based on the prospective clinical study.

Additional data on the number of surgical revisions and their reasons were collected under two retrospective studies. Each of these studies are briefly described below and comparisons of revision data of these two retrospective studies and the prospective study are presented in Tables 4 and 5.

These studies involved only devices without InhibiZone™.

Patient Information Form (PIF) Analysis—A retrospective analysis of patients implanted (n=12,713) in the U.S. during the five-year period 1995-1999. The analysis examined PIF data voluntarily sent to the manufacturer by the implanting physician for original implants and revisions. PIF data is required to be on file with the manufacturer in order to be eligible for product replacement. Revision data presented in Table 4 and Table 5 below are based on a total of 2,116 revisions reported for 2,014 patients that required one or more revisions during the five-year period of the analysis.

Retrospective Clinical Study—The retrospective clinical study was an analysis of patients implanted (n=390) in the U.S. between 1987-1990. The study examined pre- and post-implant medical records and follow-up data collected by questionnaires and physician examinations. Post-implant data was available for 356 patients. The study followed patients for up to ten years (mean: 4.1 years; range: 0.03-10.3 years). The revision data presented in Table 4 and Table 5 below are based on a total of 317 revisions reported on 164 patients that required one or more revisions during the ten-year period of the study.

The data in Table 4 presents the percentage of patients revised during the specified follow-up period, the average number of revisions conducted on patients requiring a revision, and the number of revisions expected per 100 patients for these studies in comparison with the data of the prospective clinical study.

Table 4. Comparison of Revision Data from Three Different Clinical Studies

	Prospective Study (85 pts. over 24 months)	PIF Analysis (12,713 pts. over 5 years)	Retrospective Study (356 pts. over 9 years)
% patients (pts.) revised	16.5% (14/85)	15.8% (2,014/12,713)	46.1% (164/356)
avg. # of revisions per patients revised	1.07 (15/14)	1.05 (2,116/2,014)	1.93 (317/164)
# of revisions per 100 patients	18 (15/85)	17 (2,116/12,713)	89 (317/356)

Table 5 shows revision data stratified by each reported reason for revision from three different studies of male patients implanted with the AMS 800™ AUS. Under the PIF Analysis and Retrospective Study, multiple reasons were sometimes provided for a single revision. Therefore, in order to stratify this revision data by reason, all occurrences were included and presented as “% reason.” The total number of reasons therefore exceeds the total number of revisions reported for these studies.

Table 5. Reasons for Revision in Three Different Studies

Revision Reason ^a	Prospective Study (n=85)		PIF Analysis (n=12,713)		Retrospective Study (n=356)	
	revisions %	#	reasons ^b %	#	reasons %	#
Infection	2.4%	(2)	2.3%	(297)	8.1%	(29)
Infection/erosion	1.2%	(1)	-	-	-	-
Erosion	2.4%	(2)	3.6%	(451)	22.5%	(80)
Recurring Incontinence	2.4%	(2)	5.7%	(724)	42.4%	(151)
Fluid Loss	-	-	2.3%	(298)	9.3%	(33)
Fluid Transfer Impaired	-	-	0.3%	(38)	-	-
Pressure Too Low	-	-	1.1%	(140)	-	-
Mechanical Malfunction	3.5%	(3)	0.7%	(89)	13.8%	(49)
Migration/Malposition	3.5%	(3)	0.4%	(46)	4.8%	(17)
Iatrogenic Complications	-	-	0.4%	(51)	0.6%	(2)
Reimplantation/ Replacement	-	-	-	-	3.1%	(11)
Pain	1.2%	(1)	0.2%	(22)	1.4%	(5)
Patient Dissatisfaction	1.2%	(1)	0.2%	(27)	1.7%	(6)
Other ^c	-	-	2.4%	(305)	-	-
Not Indicated	-	-	1.9%	(242)	-	-

a Note that some adverse events in the table such as fluid loss, pressure too low, fluid transfer impaired and malposition could fall into the category of mechanical malfunction or iatrogenic error. Since information is not available to place them in either category, they are listed separately.

b Numbers of reasons can vary for the same percentage due to rounding.

c Other includes: double cuff, pressure too high, unable to activate, unable to deactivate, atrophy, difficult to operate, urinary retention, air in the system, hematoma.

Device Survival

Although it is not possible to predict exactly how long an implanted prosthesis will function in a particular patient, data was gathered from two sources on device removals and revisions to help gain insight into product performance over time. Although these sources included information gathered from male patients, it is expected that similar data would result from analyzing data exclusively from female patients. Figure 11 presents device survival results from the prospective clinical study and a Bayesian analysis that uses data from the prospective clinical study and the PIF Analysis to estimate device survival at five years.

Prospective Clinical Study—A device survival curve was calculated from data collected during a prospective clinical study (n=85) with two-year follow-up. Using Kaplan-Meier analysis, the two-year revision-free rate for the AMS 800 AUS was 79.5% (95% CI with 95% lower confidence bound 69.8%).

Bayesian Analysis—A Bayesian hierarchical model was used to evaluate device safety in the prospective clinical study. The Bayesian model estimated device survival using historical data (PIF Analysis n=12,713) and prospective clinical study data (n=85) on the AMS 800 AUS. A log-normal distribution fit the AMS 800 AUS historical data. Based on the log-normal hierarchical model, it was estimated that the five-year revision-free rate for the AMS 800 AUS is approximately 73.8% with 95% CI ranging from 67.3% to 79.6%. The results met the primary safety endpoint for the clinical study of a five-year revision free rate at 75% using a 10% delta with two-sided 95% lower bound greater than 65% (Figure 11).


Figure 11. Estimated device survival rates for AMS 800 AUS

Device Use

Retrospective Clinical Study—The study included male AMS 800 AUS patients (n=390) implanted between 1987-1990. Data from this study was used to estimate the device use for patients (n=356) with available data through 1997 (range: 0.03-10.3 years). For the retrospective study, “device use” means the span of time from implant to removal, including revisions. Life table analysis indicated that the probability for a 9-year span of device use was 83.9%. Thirty-three (33) of the 356 patients had their device removed. The remaining 323 patients had a functioning device at last contact.

Clinical Information – Females

The AMS 800 AUS is indicated for the treatment of urinary incontinence in females. A literature review was performed to evaluate the safety and efficacy of the AMS 800 AUS in female patients. A total of seventeen documents, specifically reporting clinical information for female patients only, were analyzed. Clinical data, including patient demographics, etiology, complications, and success criteria and rates are reported below.

Patients

Clinical efficacy has been reported for 938 female patients, ages ranging from 9-85, implanted with the AMS 800 AUS⁷⁻²³. The mean reported follow-up on 785 AMS 800 patients was 3.6 years⁷⁻²³. Various patient etiologies were reported, including neuropathic bladder dysfunction^{7,8,9,10,19}, stress incontinence^{7,10,12,13,14,15,18,21,22,23}, pelvic/perineal trauma^{7,10,18,19}, congenital defects (exstrophy)⁷, hysterectomy⁷, radiotherapy⁷, cystectomy⁷, nonneurogenic bladder^{8,19}, intrinsic sphincter deficiency^{8,12}, myelomeningocele^{19,20}, spinal injury¹⁹, and cerebral palsy²⁰. Prior surgery or radiotherapy were reported in 354 patients^{8,9,11,12,14,16,18,21,22}.

These studies involved only devices without Inhibizone.

Adverse Events–Females

Table 6 summarizes the reported adverse events for 938 patients receiving an AMS 800 AUS. The most common reported complication was bladder injury, which occurred in 5.8% of patients (n=54). There may have been more than one reported complication for each patient.

Table 6. Adverse Events in Patients Studied

Complication	Percentage	Reference
Bladder injury	5.8%	8,14, 17, 23
Vaginal injury	4.2%	8,10, 14, 17, 19
Cuff erosion	3.4%	9, 12, 13, 16, 19, 23
Infection/erosion	3.3%	7
Device fluid leak	2.8%	8, 9, 15, 18, 19, 21, 22, 23
Incisional hernia	1.5%	8
Urinary retention	1.3%	8, 22
Urethral injury	1.2%	8, 10, 17
Urinary urgency	1.1%	8, 12
Labia hematoma	1.0%	8, 20
Infection	0.6%	7, 10, 14, 16, 23
Pelvic abscess	0.4%	22
Early post-operative morbidity	0.4%	8
Pressure atrophy at cuff site	0.3%	23
Cutaneous erosion	0.3%	10
Bothered by device	0.3%	8
Pulmonary embolus	0.1%	7
Enterocutaneous fistula	0.1%	7
Phlebitis	0.1%	8
Upper urinary tract deterioration	0.1%	19
Superficial wound dehiscence	0.1%	22
Device malfunction	0.1%	7, 8, 19, 22, 23
Device displacement	0.1%	9
Device damage - unrelated procedure	0.1%	12
Wound seroma	0.1%	15
Pump erosion	0.1%	19
Sphincter related problem	0.1%	16
Loose cuff	0.1%	21, 22

Success Rates

Patient success rates were reported and analyzed. For those defining clinical success^{8,10,11,12,13,18,19,22,23}, rates of success were categorized as follows:

- Fully and socially continent = 0-1 pad used per day
- Occasional incontinence = Consistent use of 1-2 pads per day
- Fair/improved continence = Use of 3 pads per day
- Not improved = Use of more than 3 pads per day

Success rates are summarized in Table 7.

Table 7. Success Rates

Success Criteria	% patients
Fully and socially continent	86.9
Occasional incontinence	4.4
Fair/improved	4.3
Not improved	7.1

Patient Information Form (PIF) Analysis

A retrospective analysis was performed for female patients implanted (n=637) with the AMS 800™ AUS during the five-year period 2000-2005. The analysis examined Patient Information Form (PIF) data sent to the manufacturer by the implanting physician for original implants and revision surgeries. Of the 637 female patients implanted during the period studied, 76 patients had undergone a total of 88 revision surgeries (1.16 revisions per patient). There were a total of 121 reported reasons for revisions. Note that more than one reason may be listed for the same revision. The data in Table 8 presents the reported reasons for revision and the percentage of each reason for revision.

This analysis did not include AMS 800 AUS devices treated with InhibiZone™.

Table 8. PIF Analysis Revision Data

Reason for revision	% revisions
Incontinence	19.8
Not specified	19.8
Infection	13.2
Fluid loss	11.6
Erosion	10.7
Other medical condition	8.3
Device out of specification/malfunction	5.8
Improper size of device	3.3
Other*	3.3
Urinary retention	1.7
Pain	1.7

*Other includes: hematoma, device migration/malposition, patient dissatisfaction.

HOW SUPPLIED

The components of the AMS 800 AUS and Accessory Kit are supplied sterile.

- Do not use if package is open or damaged.
- Do not use if labeling is incomplete or illegible.
- Do not use if the device is past the "Use By" date.

Handling and Storage

InhibiZone treated components: Keep dry and protect from light. Recommended storage at 25°C (77°F); excursions permitted to 15°C-40°C (59°F-104°F).

Non-InhibiZone treated components: Store in a cool, dry, dark place.

Supplemental Surgical Tools

The following BSC surgical tools which are used to help facilitate implantation of the AMS 800 AUS may be ordered separately.

Non-Sterile Tools

The following tools are provided non-sterile and must be sterilized before use. For reprocessing information, refer to the instructions manual provided with the tools.

- Quick Connect Assembly Tool

In order to use the Quick Connect Window Connectors in the AMS 800 Accessory Kit, the Quick Connect Assembly Tool must be ordered. This is a reusable stainless steel instrument used to close the connectors (Figure 12).


Figure 12. Quick Connect Assembly Tool

- Insertion Package (2 stainless steel Tubing Passers)

Two curved tubing passers are used to route the tubing of the components from incision site to incision site. The machined end of the passer provides a snug fit in the tubing lumen (Figure 13).


Figure 13. Tubing Passers

Sterile Tools

- SKW Deep Scrotal Retraction System
- Deactivation Package

Contents:

(3) Tubing plugs

(1) Straight Suture-Tie Connector

The Deactivation Package is not required for an initial implant. The Package contains 3 stainless steel tubing plugs and 1 straight Suture-Tie Connector. The tubing plugs may be used to prevent fluid from entering or leaving the prosthesis during revision surgeries.

Required Materials

Supplies and Instrument Requirements

Instruments normally required for a urological procedure are recommended. In addition to the AMS 800™ AUS components and the AMS 800 Accessory Kit, the following sterile setup should be available.

- Sterile stainless steel tray
- 1000 ml graduate container
- 500 ml graduate container
- Sponge bowl
- Medicine cup
- Emesis basin
- Two 30 ml disposable syringes
- One 10 ml disposable syringe
- Eight shod hemostats (for clamping tubing)
- One straight pair of clean, sharp scissors for trimming tubing
- Hegar dilators
- Babcock clamps
- Asepto™ syringe
- Antibiotic solution (for irrigation)
- Catheter
- Penrose drain
- Umbilical tape
- Vaginal pack
- Sterile normal saline (filling and flushing solution)
- Retractor
- Extended nasal speculum (optional)
- Rectal tube (optional)
- Centimeter ruler (optional)

Use a plastic-draped surgical stand or a stainless steel tray as a station for handling and filling the components of the prosthesis.

Position splash basins so that the surgeons can conveniently clean their gloves during the surgical implantation procedure, especially before they make the tubing connections.

Filling Solutions

The fluid used to fill the prosthesis must be sterile and free of blood and particulate matter. The presence of any foreign materials in the fluid can affect the operation of the prosthesis. The solution must also be isotonic to minimize the transfer of fluid across the silicone membrane, which is semi-permeable. **Normal saline is the recommended isotonic solution to use when filling the prosthesis.**

However, if contrast media is preferred, one of the tested solutions in Table 9 may be used for filling. If you do not use the contrast media in the mixture proportions, you may alter the isotonicity of the mixtures and promote the formation of particulate matter.

Note: The products listed below are some of the radiographic solutions tested by Boston Scientific for use in the AMS 800 AUS; only sterile water should be used for dilution. For a complete list, contact Boston Scientific.

CAUTION: Do not use sterile saline or lactated Ringer's solution to dilute the contrast solutions.

WARNING: Contrast media are contraindicated if the patient has an iodine allergy.

Table 9. Filling Solution Dilution

Contrast Media	Dilution	Validated for InhibiZone™ Use
Conray™ 43	30 ml Conray 43 + 60 ml sterile H ₂ O	Yes
Cysto-Conray™ II	60 ml Cysto-Conray II + 15 ml sterile H ₂ O	Yes
Isovue™ 200	60 ml Isovue 200 + 23 ml sterile H ₂ O	No
Isovue™ 300	57 ml Isovue 300 + 60 ml sterile H ₂ O	No
Isovue™ 370	38 ml Isovue 370 + 60 ml sterile H ₂ O	No
Omnipaque™ 180	60 ml Omnipaque 180 + 14 ml sterile H ₂ O	No
Omnipaque™ 240	60 ml Omnipaque 240 + 38 ml sterile H ₂ O	No
Omnipaque™ 300	57 ml Omnipaque 300 + 60 ml sterile H ₂ O	Yes
Omnipaque™ 350	48 ml Omnipaque 350 + 60 ml sterile H ₂ O	No
Telebrix 12	53 ml Telebrix 12 + 47 ml sterile H ₂ O	Yes
Use an equivalent ratio of contrast media and sterile water for a larger total volume.		

OPERATIONAL INSTRUCTIONS

Surgical Approaches

There are several approaches to implant the AMS 800 AUS components. It is important for the surgical staff to know which approach the surgeon intends to use because the approach will influence the position of the patient, instrumentation, and surgical procedure.

The following are brief descriptions of surgical approaches that can be used to implant the AMS 800 AUS.

Male Perineal Approach

The patient is placed in the lithotomy position. The cuff is placed around the bulbous urethra (Figure 14) by making a midline perineal incision followed by dissection of the bulbocavernosus muscle from around the urethra. The PRB is placed in the prevesical space, by making a transverse suprapubic incision through the rectus fascia, followed by dissection of the linea alba and prevesical space. The pump is placed within a subdartos pouch in the scrotum through the suprapubic incision, using blunt dissection. A second cuff can be placed using this surgical approach; for additional information refer to the Troubleshooting section.

Transverse Scrotal Approach

The patient is placed in the supine position with legs gently abducted with spreaders. The cuff is placed around the proximal bulbar urethra (Figure 14) by making a transverse scrotal incision followed by dissection of the subcutaneous tissue and scrotal septum. The PRB is placed in the prevesical space by displacing the transverse scrotal incision and dissection of the transversalis fascia. Alternatively, to avoid the necessity of piercing the fascia in patients with scarred retroperitoneum from radiation or radical surgery, create a pouch beneath the rectus abdominus and anterior to the transversalis fascia. The pump is placed within a subdartos pouch in the scrotum using blunt dissection. A second cuff can be placed using this surgical approach; for additional information refer to the Troubleshooting section.


Figure 14. Bulbous urethra placement

Male Bladder Neck Approach

The patient is placed in the supine position with legs gently abducted with spreaders. The cuff is placed around the bladder neck (Figure 15) by making a suprapubic incision followed by dissection around the bladder neck. The PRB is placed in the prevesical space by blunt dissection. The pump is placed within a subdartos pouch in the scrotum using blunt dissection.


Figure 15. Bladder neck cuff placement, male

Female Single Incision Bladder Neck Approach

The patient is placed in the supine position with legs gently abducted with spreaders. The cuff is placed around the bladder neck (Figure 16) by making a suprapubic incision followed by dissection around the bladder neck. The PRB is placed in the prevesical space by blunt dissection. The pump is placed within a pouch in the labium using blunt dissection.

Female Double Incision Bladder Neck Approach (Transvaginal)

The patient is placed in the lithotomy position. The cuff is placed around the bladder neck by making a transvaginal incision in the shape of an inverted "U" on the anterior wall of the vagina, followed by dissection from around the bladder neck. The PRB is placed in the prevesical space by making a suprapubic incision followed by blunt dissection. The pump is placed within a subdartos pouch in the labium through the suprapubic incision, using blunt dissection.


Figure 16. Bladder neck cuff placement, female

PRE-PROCEDURAL INSTRUCTIONS

Use appropriate surgical and sterile techniques per the hospital's standard procedures.

CAUTION: Be certain that the components do not come into contact with paper or cloth drapes. Fragments of these can cause possible obstruction of fluid flow if they enter the device.

Open the AMS 800™ Accessory Kit

1. Remove the tray from the box in the operating room.
2. Remove inner tray from outer tray and place the inner tray on a sterile, lint free or plastic draped surgical stand.
3. Open inner tray.

Note: Record the part and serial/lot numbers of the AMS 800 Accessory Kit on the Patient Information Form (PIF). The adhesive label at one end of the outer box and the small, removable labels on the outer trays contain the name of the component, the part and serial/lot numbers, and the size of the components. This information is also listed on the Tyvek™ lid of the outer tray.

Prepare the Hemostats

Use the following procedure to cover the hemostat jaws with the blue tubing provided in the AMS 800 Accessory Kit:

1. Place blue tubing on both jaws of hemostats to completely cover the serrated surfaces.
2. Close the jaws until the first notch to prevent excessive pressure on the blue tubing.
3. Trim blue tubing at jaw tip with sharp, clean scissors.
4. Reserve scissors for use as tubing scissors throughout procedure.

CAUTION: When clamping KRT, do not advance a hemostat's ratchet more than one notch. Excessive pressure will permanently damage the tubing.

Prepare the Patient

Before surgery, take adequate steps to limit the risk of postoperative infection.

CAUTION: Using a device treated with InhibiZone™ does not change the need to follow normal hospital protocols for prophylactic antibiotic administration.

Once the patient is in the operating room, shave the abdominal and genital area. Following the shave, scrub the area using the approved hospital preoperative scrub procedure.

Establish the sterile field, drape, and prepare the patient. Throughout the procedure, flush the surgical site with copious amounts of broad-spectrum antibiotics. Position the patient according to the preferred surgical approach.

PROCEDURAL INSTRUCTIONS

Make the Incision and Dissect

1. Place a Foley catheter or a 20 Fr sound to facilitate urethra identification. The catheter will help decompress the bladder and help avoid bladder injury during PRB placement.
2. Make the appropriate incision for the surgical approach chosen.
3. Bluntly dissect to the urethra (Figures 17, 18, 19).

CAUTION: Ensure that there is an adequate amount of muscle tissue at the implant site.


Figure 17. Dissect bulbocavernosus muscle

Figure 18. Posterior dissection of the urethra


Figure 19. Dissect around bladder neck

4. Assess the size of the urethra and select the appropriate cuff size (Figures 20, 21).
 - a. Remove the catheter or sound from the urethra. Place the cuff sizer (or Penrose drain) around the urethra where the cuff is to be implanted. It should fit snugly without constricting the urethra.

NOTE:

- A cuff length of 4.0 cm to 4.5 cm is typically required for bulbous urethra placement.
- Do not stretch the cuff sizer before use.


Figure 20. Assess the size of the urethra


Figure 21. Assess the size of the urethra

- b. The cuff size only provides approximate size of urethral circumference. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on cuff sizing. Cuff length size describes the outside circumference of cuff when it encircles the urethra.

Cuff Preparation, Placement, and Filling

1. Prepare the cuff so that the cuff is aspirated of fluid and air.
 - a. Select appropriate cuff.
 - b. Attach a 15 gauge, blunt tip needle to a 30 ml syringe.
 - c. Fill the syringe with approximately 10 ml of filling solution.
 - d. Hold the cuff with one hand and squeeze it until the air is removed.
 - e. Insert the needle into the end of the cuff tubing.
 - f. Fill with solution and aspirate air from the cuff until air bubbles are removed. Manipulate the cuff as needed, taking care not to overfill or stretch the cuff. (Figures 22, 23).

Note: The amount of fluid needed to fill a cuff will vary depending on the cuff size (approximately 1-5 ml). Larger cuffs require more fluid than smaller cuffs.


Figure 22. Manipulate cuff


Figure 23. Do not overfill the cuff

- g. Inspect for bubbles. If air remains in either the cuff or the tubing, repeat the above steps.
- h. Pull back on the plunger and squeeze the cuff to aspirate the fluid.

CAUTION: Do not over-aspirate the cuff because air can be drawn into the system through the semi-permeable silicone membrane.

- i. With the air and fluid removed from the cuff, and keeping the needle within the tubing, clamp one shod hemostat (one notch only) 3 cm below the needle. Clamp a second shod hemostat (one notch only) 3 cm below the first (Figure 24).


Figure 24. Aspirate air from cuff, then clamp

- j. Remove the 15 gauge blunt tip needle from the cuff tubing.
- k. For an InhibiZone™ treated cuff, place the prepared cuff onto an empty sterile tray or kidney basin and cover with a sterile drape.

CAUTION: Soaking antibiotic treated devices in saline will cause the antibiotics to diffuse off the device into the solution. This will cause the solution to turn orange and will reduce the concentration of antibiotics on the device.

- l. For a non-InhibiZone treated cuff, submerge the cuff in a kidney basin of filling solution until the surgeon is ready to implant the cuff.
2. Implant the cuff.
 - a. Prior to implanting, the cuff should be inspected for entrapped air.
 - b. Position the cuff at the implant site with the mesh backing toward the outside and the inflatable side toward the urethra (Figures 25, 26).

Note: The tubing should be on the patient's preferred pump side.


Figure 25. Implant cuff


Figure 26. Implant cuff

- c. Pass the prepared cuff, tab first, under the urethra and grasp the cuff tab with a shod hemostat to wrap the cuff around the urethra.

CAUTION: To avoid damage to the cuff, grasp the cuff tab with a shod hemostat.

- d. Pass the end of the cuff tubing through the hole until the hemostat meets the hole. Clamp a hemostat onto the cuff tubing on the opposite side of the hole, and then release the first hemostat, so air does not enter the cuff.

Note: Flushing the tubing before clamping the second hemostat prevents air from entering the cuff.

- e. Pull the remainder of the tubing through the hole and close the cuff by pulling the tab over the tubing adaptor (button). Ensure that the edges of the hole fit into the slot of the adaptor.
- f. Rotate the cuff so the adaptor tubing avoids contact with the cuff.

PRB Preparation, Placement, and Filling

1. Select the appropriate PRB.

Note: The surgeon usually selects the lowest PRB pressure needed to maintain closure of the bladder neck or bulbous urethra.

2. Prepare the PRB.
 - a. Attach a 15 gauge, blunt tip needle to a 30 ml syringe.
 - b. Fill the syringe with approximately 25 ml of filling solution.

- c. Hold the PRB with one hand and squeeze it until the air is removed.
- d. Insert the needle into the end of the PRB tubing.
- e. Fill PRB with 20 ml of solution, then aspirate the air from the PRB until air bubbles are removed. Manipulate the PRB as needed, taking care not to stretch the PRB (Figures 27, 28, 29).


Figure 27. Rotate the balloon


Figure 28. Collect air bubbles into one bubble


Figure 29. Aspirate air from PRB, then clamp

CAUTION: Do not over-aspirate the PRB because air can be drawn into the system through the semi-permeable silicone membrane.

- f. Inspect for bubbles, if air remains in either the PRB or tubing, repeat the steps above.
- g. Pull back on the plunger and squeeze the PRB to aspirate the fluid.
- h. With the air and fluid removed from the PRB, and keeping the needle within the tubing, clamp one shod hemostat (one notch only) 3 cm below the needle (Figure 29). Clamp a second shod hemostat (one notch only) 3 cm below the first.
- i. Remove the 15 gauge blunt tipped needle from the PRB tubing.
- j. Submerge the PRB into a kidney basin of filling solution until the surgeon is ready to implant the PRB.

CAUTION: Do not place any hemostats on top of the PRB. Any instruments on the PRB can damage it.

3. Dissect. Then, position the PRB within the prevesical space or anterior to the transversalis fascia on the patient's preferred pump side.
 - a. Make the incision and/or bluntly dissect to create a space for the PRB.
 - b. Prior to implanting, the PRB should be inspected for trapped air.
 - c. Place the PRB in the appropriate space per surgical approach.
 - d. Route the cuff and PRB tubing for later connections as needed (Figure 30).


Figure 30. Route the balloon KRT and the cuff KRT

4. Flush and fill the PRB.
 - a. Flush the PRB tubing end using a 22 gauge needle attached to a 10 ml syringe filled with filling solution.
 - b. Attach a 15 gauge needle to a 30 ml syringe filled with filling solution to PRB and unclamp the tubing.
 - c. Fill the PRB with approximately 22 ml of the chosen filling solution.

Note: Larger cuff sizes may require more filling solution. See Cuff Pressurization Option.

- d. Clamp (one notch only) the tubing approximately 3 cm below the needle with a shod hemostat.

Cuff Pressurization Option:

1. Flush the PRB and the cuff tubing.
2. Make a temporary connection with a Suture-Tie Connector between the PRB and the cuff.
3. Unclamp; wait one minute; then re-clamp.
4. Disconnect PRB and cuff.
5. Flush PRB tubing. Unclamp and aspirate filling solution.
6. Refill with 20 ml of filling solution and clamp.
7. Secure PRB within tissue by narrowing the opening with absorbable suture.

Pump Preparation and Positioning

Prepare Control Pump

1. Submerge both ends of the tubing in a kidney basin with the appropriate filling solution.
2. Hold the pump at a 45 degree angle with the black tubing on top (Figure 31).


Figure 31. Place ends of tubing in filling solution

3. Squeeze and release pump bulb repeatedly until all air in the pump and tubing has been displaced with fluid.
4. Inspect for bubbles, if air remains, repeat above steps.
5. While keeping tubing submerged, use shod hemostats to clamp (one notch only) each tube 4-5 cm from the end (Figure 32).


Figure 32. Clamp tubing with ends in filling solution

- For an InhibiZone™ treated pump, place the pump onto empty sterile tray or kidney basin and cover with a sterile drape.

CAUTION: Soaking antibiotic treated devices in saline will cause the antibiotics to diffuse off the device into the solution. This will cause the solution to turn orange and will reduce the concentration of antibiotics on the device.

- For a non-InhibiZone treated pump, submerge the pump into a kidney basin of filling solution until the surgeon is ready to implant the pump.

Pump Placement

- Use blunt dissection to create a dependent pouch in the scrotum or labium on the patient's preferred side (Figures 33, 34).
- Prior to implanting, the pump and tubing should be inspected for trapped air.
- Place pump into pouch with deactivation button facing outward so that it is palpable (Figures 35, 36).
- Route the tubing to appropriate surgical incision.
- Loosely tie a suture around the opening of the pouch to secure the pump position.


Figure 33. Create pouch


Figure 34. Create pouch


Figure 35. Place pump into pouch


Figure 36. Place pump into pouch

CONNECT TUBING

Quick Connect Window Connectors or Suture-Tie Connectors are used to connect the tubing between components.

Note: Clear tubing connects the control pump to the occlusive cuff. Black tubing connects the control pump to the pressure regulating balloon.

Using Quick Connect Window Connectors

Note: The Quick Connect System components cannot be resterilized. Conventional hospital sterilization will damage the connector components. However, the Quick Connect Assembly Tool may be resterilized. For reprocessing information, refer to the instructions manual provided with the tool.

- Cut the tubing length to fit the patient's anatomy with clean, straight scissors, making sure the cut end is square (Figure 37).


Figure 37. Cut the tubing squarely

- Insert the small diameter portion of the collet holder into the tubing.
- Slide a collet onto the tubing, ensuring that the teeth of the collet face toward the open tubing end (Figure 38).


Figure 38. Slide a collet onto each end of the tubing

- Use a 22 gauge needle attached to a 10 ml syringe and filled with filling solution to flush the connector and tubing to remove particulate matter and air.
- Firmly push one end of the tubing to the middle wall of the connector and check the tubing placement through the connector window (Figure 39).


Figure 39. Firmly push the first end of the tubing onto the connector

Note: Three types of Quick Connect Window Connectors are available depending on the location of the tubing: Straight Connector, Right Angle Connector, and Y Connector.

CAUTION:

- Right Angle Connectors should always be used when the tubing makes a sharp curve at the point of connection.
- The long branch of the Y Connector must be connected to the control pump.

- Flush tubing and connector before firmly pushing the other end of the tubing to the middle wall (Figure 40). Check the connector window to make certain both tubing ends are still touching the middle walls of the connector.


Figure 40. Firmly push the remaining tubing onto the connector

- Slide the collets towards the connector until the collet teeth touch the connector (Figure 41).


Figure 41. Slide the collets towards the connector until the teeth touch

8. Place the tubing connection in the Quick Connect Assembly tool jaw and squeeze the tool handles until the closure stop touches the opposite handle. Do not hold the tubing connection while squeezing the tool jaw as this could cause misalignment. Check that all the collet teeth are entering the connector before completing the connection.

- If using a Straight Connector, use the tool once (Figure 42).
- If using a Right-Angle Connector, use tool twice, once on each end of the connector (Figure 43).
- If using a Y Connector, the assembly tool must be used 3 times - once on each end of the connector (Figure 44).


Figure 42. Place the connection in the tool jaw with the tubing exiting straight from the connector ends


Figure 43. Use the tool twice when using a right angle connector – once on each connector end


Figure 44. Quick Connect Assembly Tool (approach from the side of the Y Connector)

CAUTION:

- Check the tubing before closing the assembly tool. Do not trap tubing between the assembly tool jaw and the connector. The tubing must exit straight from the ends of the connector, through the slots in the assembly tool.
- After using assembly tool, tubing should bulge through connector window. This indicates that tubing is firmly against middle of connector wall.

9. After all connections have been made, cycle the device to confirm function and urethral coaptation, and deactivate the device (see Deactivate Cuff section).

Suture-Tie Connectors

All connections using Suture-Tie Connectors are tied with 3-0 non-absorbable polypropylene.

1. Cut the tubing length to fit the patient's anatomy with clean, straight scissors, making sure the cut end is square (Figure 37).
2. Use a 22 gauge needle attached to a 10 ml syringe filled with filling solution to flush connector and tubing to remove particulate matter and air.

3. Push the tubing over the ends of the connector so that they meet at the center hub of the connector.

Note: Make sure the tubing is on the connector straight.

4. Use a double-throw overhand surgeon's knot followed by a minimum of two single throws to attach the tubing to the connector (Figure 45).


Figure 45. Tie suture

Note: The suture should crimp, but not cut the tubing.

5. Pass the suture 180 degrees and use the same tying technique on the opposite side of the connector. Use another suture and repeat on the opposite end of the connector.

6. After all connections have been made, cycle the device to confirm function and urethral coaptation, and deactivate the device (see Deactivate Cuff section).

Deactivate Cuff

The device system must be left in the deactivated mode for four to six weeks following implantation.

1. Squeeze and release the pump bulb several times to remove all fluid from cuff (Figure 46).


Figure 46. Squeeze and release pump bulb

Note: The cuff will be empty when the pump remains flat.

2. Allow pump bulb to partially refill (approximately 30 to 60 seconds).

Note: It is recommended to record the time needed to refill the pump and the number of pump squeezes needed to empty the cuff and pump as this information is helpful postoperatively.

3. When slight indentation in pump bulb is felt, press the deactivation button (Figure 47). After the deactivation button is pressed, the pump bulb may feel firmer than usual.


Figure 47. Press deactivation button when slight indentation is felt

Note:

- It is important to leave a slight indentation in the pump bulb to ensure there is enough fluid in the pump to activate the device later.
- You will always feel the deactivation button, even when the device is deactivated.
- The pump bulb will remain partially filled until activated. When the device is deactivated, the cuff will not be inflated (urethra not closed) and your patient will be incontinent. No fluid will move into the cuff or pump when it is deactivated.

POST-PROCEDURE INFORMATION

Some surgeons use a prophylactic antibiotic prior to surgery and intravenous antibiotics immediately postoperatively. Most choose to send their implanted patients home with a five to ten day course of antibiotics. The following paragraphs provide additional details on postoperative care.

Immediately Postoperative

After the surgery, deactivate the cuff and introduce a catheter into the urethra prior to closing. The length of time the catheter is left in place is at the discretion of the physician.

After 24 hours, the nursing staff may place ice packs in the region of the pump to reduce postoperative edema. The patient should be advised about using absorbent pads or condom catheters until the device is activated four to six weeks after the surgery. The patient should be advised to avoid undue compression of the cuff area.

After Release from Hospital

The patient is usually discharged within one to four days after surgery.

The patient must return to the physician's office to activate the device prior to use of the AMS 800™ AUS. The device is normally activated four to six weeks postoperatively. To determine that the patient is ready to use the device, check the incision site to be sure that it has healed properly.

At that time, instruct the patient that it is possible to begin using the prosthesis to urinate. The patient may experience some mild discomfort the first few times the prosthesis is being used.

Note: During voiding, patients with a slow or weak stream may find that voiding has not been completed before the cuff automatically closes and restores continence. In this case, the patient should be advised to deflate the cuff a second time to finish voiding.

Activating the device may be difficult if deactivation occurred when the pump bulb was deflated. If unable to cycle the prosthesis, squeezing the sides of the control pump adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally. Refer to the Activate (Reactivate) Cuff: Normal Method section.

There should be no redness, swelling, or drainage. Any of these things may indicate that an infection is present and the infection should be treated appropriately. Ask the patient about tenderness and/or discomfort when cycling the device.


The physician may want to observe the patient for up to an hour in the office to determine if sufficient continence is achieved with the device activated.

Provide patient education about the operation. Patients should carry a permanent patient ID card to inform others of the device in case of emergency and to prevent catheterization without device deactivation (which could damage the urethra or the device).

Activate (Reactivate) Cuff: Normal Method

To activate (reactivate) the device, complete the following instructions. See the Troubleshooting section for optional methods.

Push the deactivation button a few times to loosen the poppet inside the valve block. Then give the pump bulb a quick, forceful squeeze (Figure 48).

**Figure 48. Activate (reactivate) the cuff:
Normal method**


Figure 49. Cuff and pump refilled

Note:

- This will move the deactivation poppet back to activated position. After the device is activated, the pump will fill first and then cuff will refill (Figure 49). It will take a few minutes for the device to refill and for the cuff to close off the urethra or bladder neck. When the system is activated, the pump may become less firm.
- If it is difficult to activate the device, there may not be enough fluid remaining in the pump to push the deactivation button to its activated position. Use an optional method described in the Troubleshooting section if this happens.

Evaluating Long-Term Function and Placement

After the postoperative healing period, the surgeon should continue to have contact with the patient at least on an annual basis to evaluate the function of the device. During the annual evaluation, the surgeon should ask the patient about how the device is functioning and if he or she has noticed any changes in the function.

If the patient is having mechanical difficulty with the device, or there is infection or erosion present, revision or removal surgery may be necessary. In the event of a revision surgery, follow the same preparation and implantation techniques outlined herein. Suture-Tie Connectors must be used to connect any original tubing to new tubing.

CAUTION: In a revision case, when radiopaque solution was used for the primary case, and one or more primary components are preserved, use the same radiopaque solution at the same density. Device performance may deteriorate if different types or densities of filling solutions are mixed and the isotonicity may be lost or particles may form.

Deactivation Package

The Deactivation Package is useful during revision surgical procedures. The following example describes explanting only the cuff and preserving the other components.

When the cuff needs to be explanted due to erosion, first clamp the tubing (clear KRT) with shod hemostats at both sides of the connector. Then, cut the tubing (clear KRT) and explant the cuff. Unless infected, the pressure regulating balloon and control pump may be left in the body, with the tubing plugs inserted into the tubing ends while the tissue around the urethra heals after the cuff is explanted.

By inserting a tubing plug into the end of the control pump tubing (clear KRT), the surgeon can protect the filling solution from contamination by blood or other materials during deactivation.

1. Flush the inside of the tubing end with a 22 gauge blunt needle attached to a syringe, then insert a tubing plug from the Deactivation Package.
2. Suture-tie the tubing plug with 3-0 non-absorbable suture. When tying, use one double-throw, overhand surgical knot followed by at least two single throw knots to attach the connector to the tubing. The suture should clamp the tubing, but not so tight that it cuts the tubing.
3. The tubing plug connected to the control pump tubing is placed in the same superficial location as the original connectors, and the incision is closed.
4. After healing takes place, when the new cuff is inserted, clamp the tubing beneath the tubing plug with a shod hemostat. Remove the plug. Make sure the ends of the tubing have been cut squarely. Using a 22 gauge needle attached to a 10 ml syringe, flush the end of the tubing with system fluid. Reconnect to the component using a Suture-Tie Connector. No additional filling or priming is necessary as all the original fluid is maintained within the system.

TROUBLESHOOTING**Attaching a Second Cuff to the System**

Clinical literature reports that some patients with severe stress incontinence may continue to report some degree of incontinence even after the placement of an artificial urinary sphincter. In such cases, two cuffs may be implanted around the bulbous urethra in males. The second cuff is connected to the device system using a Y Connector included in the AMS 800 Accessory Kit. If the patient is experiencing incontinence secondary to urethral atrophy at the site of the first cuff, the first cuff may also require replacement with a smaller cuff or cuff may need to be repositioned. Usable cuff combinations are shown in Table 10.

Table 10. Double Cuff Combinations

Cuff Sizes	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	9.0	10.0	11.0
3.5	●	●	●	●	●	●	●	●	●	●	●	●	
4.0	●	●	●	●	●	●	●	●	●	●	●	●	
4.5	●	●	●	●	●	●	●	●	●	●	●	●	
5.0	●	●	●	●	●	●	●	●					
5.5	●	●	●	●	●	●	●	●					
6.0	●	●	●	●	●	●	●	●					
6.5	●	●	●	●	●	●	●	●					
7.0	●	●	●	●	●	●	●	●					
7.5	●	●	●										
8.0	●	●	●										
9.0	●	●	●										
10.0	●	●	●										
11.0													

● Indicates cuff combinations that may be used.

Place a Second Cuff

- Deactivate the AMS 800™ AUS (see Deactivate Cuff section) and place a Foley catheter or a 20 Fr sound to facilitate urethra identification. The catheter will help decompress the bladder and help avoid bladder injury during PRB placement.
- Follow the chosen surgical approach to make an incision and dissect, to expose the urethra, cuff, and existing clear and black tubing.
- Clamp the clear and black tubing on each side of the current connectors.
- Cut the tubing on both sides of the connector(s), between the clamp and the connector. Remove the cuff if replacing.
- Determine appropriate cuff size combination.
- Prepare the cuff(s) as described in the Cuff Preparation, Placement, and Filling section.
- Implant the cuff(s) with a 1 to 2 cm gap between the cuffs, to prevent them from rubbing against each other, and to maintain vascularization.
- Fill the cuff using one of the following methods:
 - Flush and fill the new cuff with 1-2 ml of filling solution and re-clamp the cuff, or
 - Aspirate all fluid and refill the PRB with 24-26 ml filling solution, then re-clamp the PRB, or
 - Use the cuff pressurization option (see the PRB Preparation, Placement, and Filling, Cuff Pressurization Option section), and use the Y connector for the temporary connection.
- Connect the system using an appropriate connection method: clear cuff tubing to clear pump tubing, black PRB tubing to black pump tubing.

Note: The short ends of the Y Connector are attached to the cuffs.

Problem	What to Do
Total Device	
The device fails to cycle.	Check connections between components. If they are correct, change the entire device.
Leak in any of the components.	By pump manipulation or by ultrasound, determine whether there is a leak. If there is a leak, replace all the components (because body fluid has entered into the system).
Occlusive Cuff	
Too tight or too loose around the urethra.	Remove improperly-sized cuff. Re-size the urethra with the cuff sizer and implant the proper size.
Punctured or damaged.	Remove and replace with new cuff.
Pressure Regulating Balloon	

Punctured during filling.	Remove and replace with new pressure regulating balloon.
Control Pump	
Difficulty activating (or reactivating) the device.	Squeeze and release the sides of the control pump adjacent to the deactivation button to allow fluid to fill the pump bulb. When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze.

Activate (Reactivate) Cuff: Optional Methods

If the normal activation method does not work, use one of the following optional methods.

Side Squeeze Method

- Squeeze sides of control pump adjacent to deactivation button to allow fluid to fill the pump bulb (Figure 50).

Note: It may take several minutes for the pump to refill.

- When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze to reactivate the system.


Figure 50. Side squeeze method

Cotton Swab Method

- Feel control pump to locate deactivation button.
- Take a cotton tip swab and apply pressure to area directly behind deactivation button (Figure 51).

Note: This should unseat the poppet and allow fluid to fill pump and then cuff.

- When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze to reactivate the system.


Figure 51. Cotton swab method

Fulcrum Method

- Feel the control pump, locate the deactivation button, and place your index finger above it (at tubing side) (Figure 52).
- Place tip of your thumb below deactivation button on the opposite side.
- Place index finger of your other hand on firm portion of pump (valve block portion) in front of deactivation button (toward the pump bulb).
- Firmly bend pump end down to activate, by using thumbs as a fulcrum.
- Release after bending.
- When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze to reactivate the system.


Figure 52. Fulcrum method

REFERENCES

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*. 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. *Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up*. *Journal of Urology* 1997; 158:435-439.
- Litwiler SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S., *The Artificial Urinary Sphincter for Post-radical*

Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. Journal of Urology 1996; 155:1260.

7 Thomas K, Venn SN, Mundy AR, *Outcome of the Artificial Urinary Sphincter in Female Patients, Journal of Urology*, Apr 2002, Vol. 167, pp. 1720-1722

8 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L, *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test, Journal of Urology*, Apr 2001, Vol. 165:4, pp. 1172-1176

9 Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H, *Explantations and Revisions with Artificial Urinary Sphincter AMS 800™ in Females, Journal of Urology*, May 1996, Vol. 155:5, Supplement

10 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C, *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up, Journal of Urology*, 1996, Vol. 155:5, supplement

11 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H, *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse, Acta Urol Bel*, May 1995, Vol.63(2), pp. 45-46

12 Hadley R, Lioses P, Dickinson M, *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon, Journal of Urology*, Apr 1995, Vol. 153, supplement

13 Stone KT, Diokn AC, Mitchell BA, *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females, Journal of Urology*, Apr 1995, Vol. 153, supplement

14 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H, *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence, Journal of Urology*, May 1994, Vol. 151:5, supplement

15 Karram MM, Rosenzweig, BA, Bhatia, NN, *Artificial Urinary Sphincter for Recurrent/ Severe Incontinence in Women: Urogynecologic Perspective, Journal of Reproductive Medicine*, Oct 1993, Vol. 38, pp. 791-794

16 Duncan HJ, Nurse DE, Mundy AR, *Role of the Artificial Urinary Sphincter in Women, British Journal of Urology*, Feb 1992, Vol. 69, pp. 141-143

17 Salisz JA, Diokno AC, *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient, Journal of Urology*, Nov 1992, Vol. 148, pp. 1528-1530

18 Webster GD, Perez LM, Khoury JM, Timmons SL, *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter, Urology*, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503

19 Parulkar BG, Barrett DM, *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females, Gynecology and Obstetrics*, Aug 1990, Vol. 171, pp. 131-138

20 Abbassian A, *A New Operation for the Insertion of the Artificial Urinary Sphincter, Journal of Urology*, Sept 1988, Vol. 140, pp. 512-513

21 Appell RA, *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach, Neurourology and Urodynamics*, 1998, Vol. 7, pp. 613-619

22 Diokno A., Hollander JB, Alderson TP, *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results, Journal of Urology*, Oct 1987, Vol. 138, pp. 778-780

23 Light JK, Scott FB, *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter, Journal of Urology*, Sept 1985, Vol. 134, pp. 476-478

24 Shellock F, *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature, AJR*, October 1988

25 Shellock F, *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review, Radiology*, 1991, Vol. 180, pp. 541-550

26 Shellock F, *MR Procedures and Biomedical Implants, Materials and Devices: 1993 Update, Radiology*, 1993, Vol. 189, pp. 587-599

27 Shellock F, *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

PRODUCT REPLACEMENT INFORMATION

For more information about the replacement of a device, component or accessory, please contact a Boston Scientific representative or local distributor. Replacement may require the return of product to Boston Scientific in accordance with the applicable Product Replacement Policy for Implanted Devices.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling and storage of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control may directly affect the device

and the results obtained from its use. BSC's obligation under this warranty is limited to the replacement of this device and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. **BSC assumes no liability with respect to devices that are reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

Tyvek is a trademark of E.I. DuPont De Nemours and Company Corporation.

Asepto is a trademark of Becton, Dickson & Co. Corporation.

Conray and Cysto-Conray are trademarks of Mallinckrodt Medical, Inc.

Isovue is a trademark of E.R. Squibb & Sons, Inc.

Omnipaque is a trademark of Sterling Drug Inc.

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Tamaños de manguito	3,5	4,0	4,5	5,0	5,5	6,0	6,5	7,0	7,5	8,0	9,0	10,0	11,0
3,5	●	●	●	●	●	●	●	●	●	●	●	●	
4,0	●	●	●	●	●	●	●	●	●	●	●	●	
4,5	●	●	●	●	●	●	●	●	●	●	●	●	
5,0	●	●	●	●	●	●	●	●					
5,5	●	●	●	●	●	●	●	●					
6,0	●	●	●	●	●	●	●	●					
6,5	●	●	●	●	●	●	●	●					
7,0	●	●	●	●	●	●	●	●					
7,5	●	●	●										
8,0	●	●	●										
9,0	●	●	●										
10,0	●	●	●										
11,0													
● Indica las combinaciones de manguitos que se pueden utilizar.													

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

REFERENCIAS

1 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.

2 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.

3 Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.

4 Haab F, Trockman BA, Zimmern PE, Leach GE. *Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up*. *Journal of Urology* 1997; 158:435-439.

5 Litwiler SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.

6 Fleshner N, Herschorn S., *The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life*. *Journal of Urology* 1996; 155:1260.

7 Thomas K, Venn SN, Mundy AR. *Outcome of the Artificial Urinary Sphincter in Female Patients*, *Journal of Urology*, Apr 2002, Vol. 167, pp. 1720-1722

8 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L, *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, *Journal of Urology*, Apr 2001, Vol. 165:4, pp. 1172-1176

9 Costa P, Mottet N, Lanfey P, Sahwi A, Dagues F, Louis JF, Navratil H, *Explantations and Revisions with Artificial Urinary Sphincter AMS 800™ in Females*, *Journal of Urology*, May 1996, Vol. 155:5, Supplement

10 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C, *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, *Journal of Urology*, 1996, Vol. 155:5, supplement

11 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H, *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, *Acta Urol Bel*, May 1995, Vol.63(2), pp. 45-46

12 Hadley R, Lioses P, Dickinson M, *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, *Journal of Urology*, Apr 1995, Vol. 153, supplement

13 Stone KT, Diokno AC, Mitchell BA, *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, *Journal of Urology*, Apr 1995, Vol. 153, supplement

14 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H, *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, *Journal of Urology*, May 1994, Vol. 151:5, supplement

15 Karram MM, Rosenzweig, BA, Bhatia, NN, *Artificial Urinary Sphincter for Recurrent/ Severe Incontinence in Women: Urogynecologic Perspective*, *Journal of Reproductive Medicine*, Oct 1993, Vol. 38, pp. 791-794

16 Duncan HJ, Nurse DE, Mundy AR, *Role of the Artificial Urinary Sphincter in Women*, *British Journal of Urology*, Feb 1992, Vol. 69, pp. 141-143

17 Salisz JA, Diokno AC, *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, *Journal of Urology*, Nov 1992, Vol. 148, pp. 1528-1530

18 Webster GD, Perez LM, Khoury JM, Timmons SL, *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, *Urology*, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503

- 19 Parulkar BG, Barrett DM, *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, *Gynecology and Obstetrics*, Aug 1990, Vol. 171, pp. 131-138
- 20 Abbassian A, *A New Operation for the Insertion of the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1988, Vol. 140, pp. 512-513
- 21 Appell RA, *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, *Neurourology and Urodynamics*, 1998, Vol. 7, pp. 613-619
- 22 Diokno A., Hollander JB, Alderson TP, *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, *Journal of Urology*, Oct 1987, Vol. 138, pp. 778-780
- 23 Light JK, Scott FB, *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1985, Vol. 134, pp. 476-478
- 24 Shellock F, *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature*, *AJR*, October 1988
- 25 Shellock F, *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review*, *Radiology*, 1991, Vol. 180, pp. 541-550
- 26 Shellock F, *MR Procedures and Biomedical Implants, Materials and Devices: 1993 Update*, *Radiology*, 1993, Vol. 189, pp. 587-599
- 27 Shellock F, *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

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RÉFÉRENCES

- 1 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- 2 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 3 Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- 4 Haab F, Trockman BA, Zimmern PE, Leach GE. *Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up*. *Journal of Urology* 1997; 158:435-439.
- 5 Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.
- 6 Fleshner N, Herschorn S., *The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life*. *Journal of Urology* 1996; 155:1260.
- 7 Thomas K, Venn SN, Mundy AR. *Outcome of the Artificial Urinary Sphincter in Female Patients*, *Journal of Urology*, Apr 2002, Vol. 167, pp. 1720-1722
- 8 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L. *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, *Journal of Urology*, Apr 2001, Vol. 165:4, pp. 1172-1176
- 9 Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H. *Explantations and Revisions with Artificial Urinary Sphincter AMS 800™ in Females*, *Journal of Urology*, May 1996, Vol. 155:5, Supplement
- 10 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C. *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, *Journal of Urology*, 1996, Vol. 155:5, supplement
- 11 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H. *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, *Acta Urol Bel*, May 1995, Vol.63(2), pp. 45-46
- 12 Hadley R, Lioses P, Dickinson M. *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, *Journal of Urology*, Apr 1995, Vol. 153, supplement
- 13 Stone KT, Diokno AC, Mitchell BA. *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, *Journal of Urology*, Apr 1995, Vol. 153, supplement
- 14 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H. *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, *Journal of Urology*, May 1994, Vol. 151:5, supplement
- 15 Karram MM, Rosenzweig BA, Bhatia, NN. *Artificial Urinary Sphincter for Recurrent/Severe Incontinence in Women: Urogynecologic Perspective*, *Journal of Reproductive Medicine*, Oct 1993, Vol. 38, pp. 791-794
- 16 Duncan HJ, Nurse DE, Mundy AR. *Role of the Artificial Urinary Sphincter in Women*, *British Journal of Urology*, Feb 1992, Vol. 69, pp. 141-143
- 17 Salisz JA, Diokno AC. *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, *Journal of Urology*, Nov 1992, Vol. 148, pp. 1528-1530
- 18 Webster GD, Perez LM, Khoury JM, Timmons SL. *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, *Urology*, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503
- 19 Parulkar BG, Barrett DM. *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, *Gynecology and Obstetrics*, Aug 1990, Vol. 171, pp. 131-138
- 20 Abbassian A. *A New Operation for the Insertion of the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1988, Vol. 140, pp. 512-513
- 21 Appell RA. *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, *Neurourology and Urodynamics*, 1998, Vol. 7, pp. 613-619
- 22 Diokno A., Hollander JB, Alderson TP. *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, *Journal of Urology*, Oct 1987, Vol. 138, pp. 778-780
- 23 Light JK, Scott FB. *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1985, Vol. 134, pp. 476-478
- 24 Shellock F. *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature*, *AJR*, October 1988
- 25 Shellock F. *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review*, *Radiology*, 1991, Vol. 180, pp. 541-550

26 Shellock F. *MR Procedures and Biomedical Implants, Materials and Devices: 1993 Update*, *Radiology*, 1993, Vol. 189, pp. 587-599

27 Shellock F. *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

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REFERÊNCIAS

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. *Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up*. *Journal of Urology* 1997; 158:435-439.
- Litwiler SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S., *The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life*. *Journal of Urology* 1996; 155:1260.
- Thomas K, Venn SN, Mundy AR. *Outcome of the Artificial Urinary Sphincter in Female Patients*, *Journal of Urology*, Apr 2002, Vol. 167, pp. 1720-1722
- Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L. *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, *Journal of Urology*, Apr 2001, Vol. 165:4, pp. 1172-1176
- Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H. *Explantations and Revisions with Artificial Urinary Sphincter AMS 800™ in Females*, *Journal of Urology*, May 1996, Vol. 155:5, Supplement
- Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C. *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, *Journal of Urology*, 1996, Vol. 155:5, supplement
- Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H. *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, *Acta Urol Bel*, May 1995, Vol. 63(2), pp. 45-46
- Hadley R, Lioses P, Dickinson M. *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, *Journal of Urology*, Apr 1995, Vol. 153, supplement
- Stone KT, Diokno AC, Mitchell BA. *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, *Journal of Urology*, Apr 1995, Vol. 153, supplement
- Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H. *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, *Journal of Urology*, May 1994, Vol. 151:5, supplement
- Karram MM, Rosenzweig BA, Bhatia, NN. *Artificial Urinary Sphincter for Recurrent/ Severe Incontinence in Women: Urogynecologic Perspective*, *Journal of Reproductive Medicine*, Oct 1993, Vol. 38, pp. 791-794
- Duncan HJ, Nurse DE, Mundy AR. *Role of the Artificial Urinary Sphincter in Women*, *British Journal of Urology*, Feb 1992, Vol. 69, pp. 141-143
- Salisz JA, Diokno AC. *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, *Journal of Urology*, Nov 1992, Vol. 148, pp. 1528-1530
- Webster GD, Perez LM, Khoury JM, Timmons SL. *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, *Urology*, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503
- Parulkar BG, Barrett DM. *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, *Gynecology and Obstetrics*, Aug 1990, Vol. 171, pp. 131-138
- Abbassian A. *A New Operation for the Insertion of the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1988, Vol. 140, pp. 512-513
- Appell RA. *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, *Neurourology and Urodynamics*, 1998, Vol. 7, pp. 613-619
- Diokno A., Hollander JB, Alderson TP. *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, *Journal of Urology*, Oct 1987, Vol. 138, pp. 778-780
- Light JK, Scott FB. *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, *Journal of Urology*, Sept 1985, Vol. 134, pp. 476-478
- Shellock F. *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature*, *AJR*, October 1988
- Shellock F. *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review*, *Radiology*, 1991, Vol. 180, pp. 541-550

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28 Shellock F, *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

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© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd