AMS 800™

Urinary Control System For Male Patients

Instructions For Use

English	AMS 800™ <i>Urinary Control System For Male Patients</i> Instructions For Use.....	1
Français	AMS 800™ <i>Système de contrôle urinaire pour les hommes</i> Mode d'emploi.....	13
Deutsch	AMS 800™ <i>System zur Blasenkontrolle für den Mann</i> Gebrauchsanweisung.....	25
Italiano	AMS 800™ <i>Sistema di controllo urinario per pazienti maschi</i> Istruzioni per l'uso.....	37
Español	AMS 800™ <i>Sistema de control urinario para varones</i> Instrucciones de uso.....	49
Português	AMS 800™ <i>Sistema de controlo urinário para doentes do sexo masculino</i> Instruções de utilização.....	61

Rx ONLY

AMS™

	<p>en Catalog number de Katalognummer es Número de catálogo fr Référence catalogue it Numero di catalogo pt Número de catálogo</p>
	<p>en Lot number de Chargennummer es N.º de lote fr Numéro de lot it Numero di lotto pt Número de lote</p>
	<p>en Date of manufacture de Herstellungsdatum es Fecha de fabricación fr Date de fabrication it Data di fabbricazione pt Data de fabrico</p>
	<p>en Use by YYYY-MM-DD de Verwendbar bis JJJJ-MM-TT es Usar antes de AAAA-MM-DD fr Utiliser avant le AAAA-MM-JJ it Usare entro AAAA-MM-GG pt Utilizar até AAAA-MM-DD</p>
	<p>en Do not reuse de Nicht wiederverwenden es No reutilizar fr Ne pas réutiliser it Non riutilizzare pt Não reutilizar</p>
	<p>en Do not reesterilize de Nicht reesterilisieren es No volver a esterilizar fr Ne pas résteriliser it Non risterilizzare pt Não reesterilizar</p>
	<p>en Sterilized Using Steam de Mit Dampf sterilisiert es Esterilizar con vapor fr Stérilisé à la vapeur it Sterilizzato a vapore pt Esterilizado por vapor</p>
	<p>en Sterilized using ethylene oxide de Mit Ethylenoxid sterilisiert es Esterilizado con óxido de etileno fr Stérilisé à l'oxyde d'éthylène it Sterilizzato con ossido di etilene pt Esterilizado por óxido de etileno</p>
	<p>en Upper Limit of Temperature de Obere Temperaturgrenze es Límite superior de temperatura fr Limite supérieure de température it Limite superiore di temperatura pt Limite superior de temperatura</p>
	<p>en Non-sterile de Nicht steril es No estéril fr Non stérile it Non sterile pt Não esterilizado</p>
	<p>en Consult Instructions for Use de Gebrauchsanweisung beachten es Consultar las instrucciones de uso fr Consulter le mode d'emploi it Consultare le istruzioni per l'uso pt Consultar as instruções de utilização</p>
	<p>en Manufacturer de Hersteller es Fabricante fr Fabricant it Produttore pt Fabricante</p>
	<p>en Authorized Representative in the European Community de Autorisierte Vertretung in der EU es Representante autorizado en la Comunidad Europea fr Représentant autorisé dans la Communauté Européenne it Rappresentante autorizzato della comunità europea pt Representante autorizado na comunidade europeia</p>
	<p>en Recyclable packaging de Wiederverwertbare Verpackung es Embalaje reciclable fr Emballage recyclable it Confezione riciclabile pt Embalagem Reciclável</p>
	<p>en CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician. de VORSICHT: Nach US-amerikanischem Recht darf dieses Gerät nur von einem Arzt oder auf ärztliche Anordnung verkauft werden. es PRECAUCIÓN: las leyes federales de los EE. UU. limitan la venta de este dispositivo a médicos o bajo prescripción facultativa. fr MISE EN GARDE : Selon la législation fédérale américaine, ce dispositif ne peut être vendu que sur prescription d'un médecin. it ATTENZIONE: le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo ai medici o dietro prescrizione medica. pt ATENÇÃO: a lei federal dos Estados Unidos só permite a venda deste dispositivo mediante receita médica.</p>

AMS 800™

Urinary Control System

For Male Patients

Instructions For Use

NOTE: For implant procedure information, consult the AMS 800 Urinary Control System Operating Room Manual.

Brief Device Description

The AMS 800 Urinary Control System is an implantable, fluid filled, solid silicone elastomer device used to treat stress urinary incontinence. It is designed to restore the natural process of urinary control. The device simulates normal sphincter function by opening and closing the urethra at the control of the patient. The AMS 800 consists of three interconnected components: a cuff, a pump, and a pressure-regulating balloon. The three components are connected with kink-resistant tubing. The AMS 800 can be implanted at either the bulbous urethra or the bladder neck.

The AMS 800 Urinary Control System Cuff and Pump are available with InhibiZone™, an antibiotic coating of rifampin (rifampicin) and minocycline.*

A System implant contains less than 2% of the oral dose exposure during a complete course of rifampin and minocycline, even if two of the largest size Cuffs are implanted. The average amounts of rifampin and minocycline contained on a prosthesis implant are represented by the means and 95% tolerance intervals of the following implant configurations:

- 1.9 mg rifampin (0.7 – 3.1 mg) and 2.8 mg minocycline (2.1 – 3.5 mg) for the implant configuration with the lowest drug levels, i.e., single 4.0 cm cuff + control pump.
- 3.7 mg rifampin (0.9 – 6.5 mg) and 6.3 mg minocycline (4.7 – 8.0 mg) for the implant configuration with the highest drug levels, i.e., double 11.0 cm cuff + control pump.

Indications for Use

The AMS 800 is used to treat urinary incontinence due to reduced outlet resistance (intrinsic sphincter deficiency).

Contraindications

1. This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
2. This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
3. This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
4. The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin or to minocycline or other tetracyclines.

* not available in all markets

5. The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline has been reported to aggravate this condition.

Warnings

1. Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection.
Infection that fails to respond to antibiotic therapy may result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.
2. Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component misplacement. The cuff may erode around the urethra or bladder neck. The control pump may erode through the scrotum. The pressure-regulating balloon may erode into the bladder. Acute urinary tract infection can interfere with proper functioning of the device and may lead to erosion of the urethra in the cuff area. Failure to evaluate and promptly treat the erosion may result in a substantial worsening of the condition leading to infection and/or loss of tissue.
3. Poor bladder compliance or a small fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
4. Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
5. Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device to prevent potential damage to the urethra or the AMS 800.
6. This device contains solid silicone elastomers. This device does not contain silicone gel. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
7. Surgical, physical, psychological, or mechanical complications, if they occur, may necessitate revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
8. Product wear, component disconnection or other mechanical problems may lead to surgical intervention. Mechanical complications may include malfunctioning of the components and leakage of fluid. Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in overflow obstruction. Mechanical events should be evaluated carefully by the treating physician and the patient should consider risks and benefits of treatment options, including revision surgery.
9. Previous patient history of adverse reaction(s) to radiopaque solution precludes its use as a filling medium for the prosthesis. Instead, saline should be used to fill the device.

10. The implanter should check that there is an adequate amount of bulbospongiosus muscle to surround and support a bulbous urethral cuff implant. Thinner spongiosum typically occurs toward the distal end of the bulbous urethra, and implantation of the cuff where the spongiosum is thin increases the chance of erosion and other complications. This warning is especially important for double cuff implants, where the second cuff is placed distal to the first implanted cuff.
11. If a hypersensitivity reaction develops to a device coated with InhibiZone, the Cuff and Pump should be removed and the patient treated appropriately.

Precautions

Patient Related

1. Patient selection requires thorough preoperative consultation and evaluation by the physician.
2. Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation of an AMS 800. Although the prosthesis is designed to restore urinary control, some patients continue to have a degree of incontinence after this procedure.
3. Patients may experience pain when the device is activated in the postoperative period and during the period of initial use. Cases of chronic pain associated with device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention. Patients should be counseled on expected postoperative pain including severity and duration.
4. Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff at the bulbous urethra or bladder neck.
5. Any progressively degenerative disease, e.g. multiple sclerosis, may limit the future usefulness of the implanted prosthesis as a treatment for the patient's urinary incontinence.
6. Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device.
7. Trauma or injury to the pelvic, perineal or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the device. The physician should advise patients of these possibilities and warn them to avoid trauma to these areas.
8. Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5cm) generally exceeds 28F. Additional clearance is required to accommodate the patient's urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on sizing.

InhibiZone™ Related

1. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin (rifampicin) and minocycline can cause additional stress on the hepatic and renal systems.

2. Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
3. Patients who receive a device with InhibiZone and are also taking Warfarin should have their prothrombin time monitored because tetracyclines have been reported to slow coagulation.
4. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses or rifampin (rifampicin).
5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.
7. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
8. Because products with InhibiZone are impregnated with a combination of rifampin (rifampicin) and minocycline, the contraindications, warnings and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline and rifampin (rifampicin) in patients receiving this device are unlikely to be detected.

Surgery Related

1. Improper cuff sizing, improper balloon selection, or other causes may result in tissue erosion, migration of components, or continued incontinence.
2. Component migration can occur if the cuff is sized improperly, if the pump or balloon is not positioned correctly, or if the tubing lengths are incorrect. Migration can result in pain, complications, device malfunction and surgical revision.
3. Unsuccessful outcomes may result from improper surgical technique, improper sterile technique, anatomical misplacement of components, improper sizing and/or filling of components.
4. Although reinforced tubing has been designed to be more resistant to tubing kinks, tubing kinks may still result from tailoring the connecting tubing to an improper length during the implant procedure.

Device Related

1. If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may arise as a result:
 - a. In the event of large pressures within the bladder, automatic pressure relief that normally occurs with the device would be prevented. Cycling the device can relieve the outflow obstruction.
 - b. Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the device, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.

- c. Release of the deactivation valve may require greater pressure than that used to cycle the device.
2. System pressure changes may occur over time if you fill the balloon with radiopaque solution of incorrect concentration. Follow the instructions in the Operating Room Manual to prepare the radiopaque solution with the correct concentration.

Adverse Events

A prospective clinical study was conducted to demonstrate the safety and efficacy of the AMS 800 Urinary Control System. A total of 87 patients were enrolled in the study and 85 patients were implanted with the device. During the study, 26 patients experienced 43 device related adverse events. Table 1 lists the device related adverse events reported during the study.

This trial involved only devices without InhibiZone.

Table 1: AMS 800 Prospective Clinical Study Device Related Adverse Events

Adverse Event Category	Total Events	Patients with AE	Events Resolved	Interventions*		
				None Reported	Medical**	Surgical
Impaired Device Function	7	6	4	2	2	4
Pain/Discomfort	6	5	4	3	3	1
Delayed Wound Healing	5	5	5	2	3	0
Bladder Spasms	2	2	0	0	2	0
Difficult Activation	2	2	2	1	1	0
Migration	3	3	1	2	0	1
Tissue Erosion	2	2	2	0	0	2
Difficult Deactivation	1	1	1	0	1	0
Infection	2	2	2	0	0	2
Recurrent Incontinence	3	3	3	1	0	2
Fistula Formation	1	1	1	0	0	1
Hematoma	1	1	1	0	1	0
Swelling	2	2	2	0	2	1
Hydrocele	1	1	1	0	1	1
Tissue Erosion/Infection	1	1	1	0	0	1
Patient Dissatisfaction	1	1	1	0	0	1
Positional Incontinence	1	1	0	1	0	0
Wound Infection	1	1	1	0	1	0
Urinary Retention	1	1	1	0	1	0

*Events may have been addressed with more than one type of intervention.

**Medical interventions included: medication, education, frequent device deactivation, dressing changes and catheterization.

How Supplied and Storage

Store device in a clean, dry, dark area at room temperature.

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Store device in a clean, dry, dark area at room temperature.

Clinical Studies

A prospective, multi-center, non-randomized clinical study was undertaken to demonstrate that the AMS 800 can be surgically implanted without serious adverse sequelae, provides

an acceptable level of continence and enhances quality of life. Each patient served as their own control. Efficacy data and safety data related to adverse events, revision surgery, diagnoses and health status evaluations were captured on case report forms. Patient self-evaluations related to health status and non-illness specific quality of life were measured on two validated outcome instruments. Patient and physician assessments of continence were measured on a recognized, standardized non-validated instrument.

Eighty-seven (87) male patients were enrolled in the study of which 85 patients were implanted with the device during the study. Patients available at the follow-up intervals were 6-months (n=67), 12-months (n=60), 18-months (n=55), and 24-months (n=41). Patients diagnosed with intrinsic sphincter deficiency (ISD) resulting from prostate surgery were eligible for enrollment. Patients with a history of allergy/sensitivity to silicone, pre-existing autoimmune or connective tissue disease or active urogenital infection were excluded from the study.

This trial involved only devices without InhibiZone.

Endpoints

The primary effectiveness endpoint evaluated the effect of the prosthesis on patient quality of life using the *Incontinence Impact Questionnaire*, an incontinence-specific quality of life questionnaire. The primary safety endpoint evaluated the five-year revision-free rate using a Bayesian hierarchical model. The safety endpoint was a five-year revision-free rate equivalent to 75% using a 10% delta with a two-sided 95% lower bound greater than 65%.

Incontinence Impact Scores

The primary effectiveness endpoint was a reduction in Incontinence Impact Score from pre- to post-implant status. Incontinence impact was measured pre- and post-implant at 6, 12, 18, and 24 months. Thirty-nine (39) patient answered the *Incontinence Impact Questionnaire* (IIQ) at 24-month follow-up. The IIQ is a 30-item, self-administered questionnaire designed to assess the impact of urinary incontinence on several subscales including physical, emotional, and social. The IIQ used in the study was developed from a validated instrument.⁴ The mean pre-implant score was significantly higher ($p<0.0001$) than mean scores at all follow-up visits. Therefore, the impact of incontinence was reduced for patients following AMS 800 implantation and the primary objective was met.

Physician and Patient Assessment of Continence

Physician assessed continence was 63.6% dry and 34.1% required some additional protection at one-year follow-up (n=43). At two-year follow-up (n=30), 73.3% were dry and 23.3% required some additional protection. Patient assessed continence was 61.7% dry and 36.7% required some additional protection at one year follow-up (n=60). At two-year follow-up (n=41), 65.9% were dry and 31.7% required some additional protection. No significant difference existed between physicians' assessment and the patients' assessment of their incontinence.

Patient Evaluation of Health Status, and Self-esteem

General Quality of Life as measured by the *Health Status Questionnaire* and the *Rosenberg Self-esteem Questionnaire* was evaluated at pre- and post-implant at 6, 12, 18, and 24 months. Thirty-eight (38) patients answered the *Health Status* and *Rosenberg Self-esteem Questionnaires* at 24-month follow-up. The self-administered *Health Status Questionnaire*² was used to assess

non-illness specific parameters such as physical functioning, social functioning, energy/fatigue, pain, health perception, and emotional problems. A high score indicates that overall health was perceived to be high. The mean score was 596 at pre-implant and 612 at two-year follow-up. No significant difference in health status scores was observed during the study. The self-administered *Rosenberg Self-esteem Questionnaire*³ was used to assess changes in patient self-esteem. The range of possible scores is 0-6, with a score of 6 indicating high self-esteem. The mean score at implant was 3.5 and at two-year follow-up was 4.1. The increase in mean score indicates a more positive self-esteem following AMS 800 implant. The device did not have an adverse effect on sexual function. Some patients with improved continence following implant also reported increased sexual activity. The positive impact of the device on patient's lives measured in the clinical study is consistent with results obtained by other authors.^{4,5,6}

Surgical Revisions

A revision is a surgical intervention related to the function, placement, or site reaction to the implanted device. For the 85 patients implanted with the device followed under the prospective clinical study, 14 patients (16.5%) experienced a total of 15 revisions up to 24 months following implant. One patient had two revisions. Three (3) revisions were due to mechanical malfunction. Two (2) revisions were due to recurrent incontinence. Two (2) revisions were due to erosion. Two (2) revisions were due to infection. One (1) revision each (total = 6) was due to migration, pain, erosion/infection, persistent incontinence/patient dissatisfaction, recurring incontinence/malfunction, infection/pain/urethrocuteaneous fistula. Multiple reasons were provided for some revisions. Four of the 14 patients who experienced revisions elected to have the device removed and 10 elected to have the device replaced. The probability of remaining revision-free 24 months following implant was 79.5% (95% CI with 95% lower confidence bound 69.8%) based on the prospective clinical study.

Additional data on the number of surgical revisions and their reasons were collected under two retrospective studies. Each of these studies are briefly described below and comparisons of revision data of these two retrospective studies and the prospective study are presented in Tables 2 and 3.

Patient Information Form (PIF) Study – The PIF study was a retrospective analysis of patients implanted (n=12,713) in the U.S during the five-year period 1995-1999. The study examined PIF data voluntarily sent to the manufacturer by the implanting physician for original implants and revisions. PIF data is required to be on file with the manufacturer in order to be eligible for product replacement. Revision data presented in Table 2 and Table 3 below are based on a total of 2,116 revisions reported for 2,014 patients that required one or more revisions during the five-year period of the study.

Retrospective Clinical Study – The retrospective clinical study was an analysis of patients implanted (n=390) in the U.S. between 1987-1990. The study examined pre- and post-implant medical records and follow-up data collected by questionnaires and physician examinations. Post-implant data was available for 356 patients. The study followed patients for up to ten years (mean: 4.1 years; range: 0.03-10.3 years). The revision data presented in Table 2 and Table 3 below are based on a total of 317 revisions reported on 164 patients that required one or more revisions during the ten-year period of the study.

The data in Table 2 presents the percentage of patients revised

during the specified follow-up period, the average number of revisions conducted on patients requiring a revision and the number of revisions expected per 100 patients for these studies in comparison with the data of the prospective clinical study.

Table 2: Comparison of Revision Data from Three Different Clinical Studies

	Prospective Study (85 pts. over 24 months)	PIF Study (12,713 pts. over 5 years)	Retrospective Study (356 pts. over 9 years)
% pts. revised	16.5% (14/85)	15.8% (2014/12713)	46.1% (164/356)
avg. # of revisions per pts. revised	1.07 (15/14)	1.05 (2116/2014)	1.93 (317/164)
# of revisions per 100 pts.	18 (15/85)	17 (2116/12713)	89 (317/356)

Table 3 shows revision data stratified by each reported reason for revision from three different studies of male patients implanted with the AMS 800. Under the PIF Study and Retrospective Study multiple reasons were sometimes provided for a single revision. Therefore, in order to stratify this revision data by reason, all occurrences were included and presented as “% reason.” The total number of reasons therefore exceeds the total number of revisions reported for these studies.

Table 3: Reasons for Revision in Three Different Studies

Revision Reason ^a	Prospective Study (n=85)		PIF Study (n=12713)		Retrospective Study (n=356)	
	revisions		reasons ^b		reasons	
	%	#	%	#	%	#
Infection	2.4%	(2)	2.3%	(297)	8.1%	(29)
Infection/erosion	1.2%	(1)	----	----	----	----
Erosion	2.4%	(2)	3.6%	(451)	22.5%	(80)
Recurring Incontinence	2.4%	(2)	5.7%	(724)	42.4%	(151)
Fluid Loss	----	----	2.3%	(298)	9.3%	(33)
Fluid Transfer Impaired	----	----	0.3%	(38)	----	----
Pressure too low	----	----	1.1%	(140)	----	----
Mechanical Malfunction	3.5%	(3)	0.7%	(89)	13.8%	(49)
Migration/ Malposition	3.5%	(3)	0.4%	(46)	4.8%	(17)
Iatrogenic Complications	----	----	0.4%	(51)	0.6%	(2)
Reimplantation/ Replacement	----	----	----	----	3.1%	(11)
Pain	1.2%	(1)	0.2%	(22)	1.4%	(5)
Patient Dissatisfaction	1.2%	(1)	0.2%	(27)	1.7%	(6)
Other ^c	----	----	2.4%	(305)	----	----
Not indicated	----	----	1.9%	(242)	----	----

a Note that some adverse events in the table such as fluid loss, pressure too low, fluid transfer impaired and malposition could fall into the category of mechanical malfunction or iatrogenic error. Since information is not available to place them in either category, they are listed separately.

b Numbers of reasons can vary for the same percentage due to rounding.

c Other includes: double cuff, pressure too high, unable to activate, unable to deactivate, atrophy, difficult to operate, urinary retention, air in the system, hematoma.

Device Survival

Although it is not possible to predict exactly how long an implanted prosthesis will function in a particular patient, American Medical Systems, Inc. gathered data from two sources on device removals and revisions to help gain insight into product performance over time. Figure 1 presents device survival results from the prospective clinical study and a Bayesian analysis that uses data from the prospective clinical study and the PIF Study to estimate device survival at five years.

Prospective Clinical Study – A device survival curve was calculated from data collected during a prospective clinical study (n=85) with two-year follow-up. Using Kaplan-Meier analysis, the two-year revision-free rate for the AMS 800 was 79.5% (95% CI with 95% lower confidence bound 69.8%).

Bayesian Analysis – A Bayesian hierarchical model was used to evaluate device safety in the prospective clinical study. The Bayesian model estimated device survival using historical data (PIF Study n=12,713) on the AMS 800 and prospective clinical study data (n=85) on the AMS 800. A log-normal distribution fit the AMS 800 historical data. Based on the log-normal hierarchical model, it was estimated that the five-year revision-free rate for the AMS 800 is approximately 73.8% with 95% CI ranging from 67.3% to 79.6%. The results met the primary safety endpoint for the clinical study of a five-year revision free rate at 75% using a 10% delta with two-sided 95% lower bound greater than 65%.


Figure 1: Estimated Survival Rates for the AMS 800

Device Use

Retrospective Clinical Study – The study included male AMS 800 patients (n=390) implanted between 1987-1990. Data from this study was used to estimate the device use for patients (n=356) with available data through 1997 (range: 0.03-10.3 years). For the retrospective study, “device use” means the span of time from implant to removal, including revisions. Life table analysis indicated that the probability for a 9-year span of device use was 83.9%. Thirty-three (33) of the 356 patients had their device removed. The remaining 323 patients had a functioning device at last contact.

Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological and functional outcome of the implantation. The risks, benefits and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal. Patients should also be aware that the AMS 800 is not considered to be a lifetime implant.

It is also important that the physician discusses with the patient

the possibility of an allergic reaction to the materials in the device (See Silicone Information).

Antibiotic Information

The antibiotics present in InhibiZone, minocycline and rifampin (rifampicin), are well characterized and have been in use for years. The dosage present on the artificial urinary sphincter is intended to act on organisms that attempt to colonize on the device.

The average amounts of rifampin and minocycline contained on a prosthesis implant are represented by the means and 95% tolerance intervals of the following implant configurations:

- 1.9 mg rifampin (0.7 – 3.1 mg) and 2.8 mg minocycline (2.1 – 3.5 mg) for the implant configuration with the lowest drug levels, i.e., single 4.0 cm cuff + control pump.
- 3.7 mg rifampin (0.9 – 6.5 mg) and 6.3 mg minocycline (4.7 – 8.0 mg) for the implant configuration with the highest drug levels, i.e., double 11.0 cm cuff + control pump.

The following in vitro data are available, but their clinical significance is unknown. No clinical studies have been performed to evaluate the effect of the antibiotic surface treatment on reducing the incidence of artificial urinary sphincter implant infections.

Table 4: In vitro Zone of Inhibition for Device Samples* with InhibiZone Treatment

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
Staphylococcus epidermidis	22.6	2.9	21
Staphylococcus aureus	17.5	5.0	25
Escherichia coli**	6.5	2.6	24
Enterococcus faecalis**	4.8	6.7	21
Candida albicans**	0.1	0.4	21
Proteus mirabilis**	0.6	1.0	17

*obtained using standardized KRT test samples containing approximately 12µg minocycline and 26 µg rifampin (rifampicin)

**the isolates tested were not susceptible to rifampin (rifampicin) and/or minocycline control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and five rabbits were implanted subcutaneously with 6 control samples each. One rabbit received three test samples and three control samples. The test samples were portions of an InhibiZone treated AMS 700 pump and the control samples were portions of a standard AMS 700 pump without InhibiZone. The AMS 700 Pumps used in the 700 Series of Inflatable Penile Prosthesis is similar to AMS 800 Pumps used in the 800 Series of Artificial Urinary Sphincters in regards to material composition, adhesive, and InhibiZone application process. All samples were soaked in a 10^3 - 10^4 CFU solution of staphylococcus aureus, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

Silicone Information

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Solid silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluoro-silicone lubricant.

There are reports of malignant tumor formation in laboratory animals only associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS 800. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

Magnetic Resonance Imaging (MRI) Important Safety Information

Non-clinical testing has demonstrated the AMS 800 product line is MR Conditional. The device can be scanned safely under the following conditions:

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Spatial Gradient Field	450 Gauss/cm or less	720 Gauss/cm or less
Maximum whole body averaged Specific Absorption Rate (SAR)	1.5 W/kg for 15 minutes of scanning as assessed by calorimetry	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry

MRI-Related Heating

Non-clinical testing has demonstrated the AMS 800 product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Highest Temperature Change	≤ +0.4°C	≤ +2.0°C
(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5		

Artifact Information

Non-clinical testing has demonstrated that the AMS 800 product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	5,800 mm ²	1,956 mm ²	6,096 mm ²	2,650 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Inventory Returns and Product Replacement Information

A Patient Information Form (PIF) must be filled out and filed with AMS at the time of implant to activate the product warranty. Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form.

Follow all of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS. Request an AMS Product Return Kit from the AMS Customer Service Department to return any explanted components to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Limited Warranty Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

Outside the United States

Customers outside of the United States should contact their local AMS Representative prior to returning any product.

This document is written for professional medical audiences. Contact American Medical Systems for lay publications.

American Medical Systems periodically updates product literature. If you have questions about the currency of this information, contact American Medical Systems.

References

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiler SE, Kim KB, Fone PD, Devere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Références

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVeré White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Literatur

- 01 Shumaker SA, Wyman JF, Übersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306, 1994.
- 02 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 03 Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- 04 Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- 05 Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- 06 Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Bibliografia

- 01 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- 02 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 03 Rosenberg M. *Society and the Adolescent Self-Image*, Princeton, New Jersey: Princeton University Press, 1965.
- 04 Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- 05 Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- 06 Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Referencias

- 01 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- 02 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 03 Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- 04 Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- 05 Litwiller SE, Kim KB, Fone PD, De Vere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- 06 Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Referências

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. *Quality of Life Research*, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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**American Medical Systems
Europe B.V.**

Haarlerbergweg 23 G
1101 CH Amsterdam Zuid-Oost
The Netherlands

©2017 Boston Scientific Corporation or its affiliates.
All rights reserved.
1003273 (P/N)
1003181 (A/W Rev A) (2017-02)

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Boston Scientific (Australia) Pty Ltd
U.S.A.