# AMS 800™

## *Urinary Control System*

### *For Female and Pediatric Patients*

## Instructions for Use

English	AMS 800™ <i>Urinary Control System</i> <i>For Female and Pediatric Patients</i> Instructions For Use.....	1
Français	AMS 800™ <i>Système de contrôle urinaire</i> <i>pour les femmes et les enfants</i> Mode d'emploi.....	11
Deutsch	AMS 800™ <i>System zur Blasenkontrolle</i> <i>für Frauen und Kinder</i> Gebrauchsanweisung.....	23
Italiano	AMS 800™ <i>Sistema di controllo urinario</i> <i>per pazienti femmine e pediatrici</i> Istruzioni per l'uso.....	35
Español	AMS 800™ <i>Sistema de control urinario</i> <i>para mujeres y niños</i> Instrucciones de uso.....	47
Português	AMS 800™ <i>Sistema de controlo urinário</i> <i>Para doentes do sexo feminino e pediátricos</i> Instruções de utilização.....	59
	References / Références / Literatur / Bibliografia / Referencias / Referências.....	71

**Rx ONLY**

**AMS™**

	<p>en Catalog number de Katalognummer es Número de catálogo fr Référence catalogue it Numero di catalogo pt Número de catálogo</p>
	<p>en Lot number de Chargennummer es N.º de lote fr Numéro de lot it Numero di lotto pt Número de lote</p>
	<p>en Date of manufacture de Herstellungsdatum es Fecha de fabricación fr Date de fabrication it Data di fabbricazione pt Data de fabrico</p>
	<p>en Use by YYYY-MM-DD de Verwendbar bis JJJJ-MM-TT es Usar antes de AAAA-MM-DD fr Utiliser avant le AAAA-MM-JJ it Usare entro AAAA-MM-GG pt Utilizar até AAAA-MM-DD</p>
	<p>en Do not reuse de Nicht wiederverwenden es No reutilizar fr Ne pas réutiliser it Non riutilizzare pt Não reutilizar</p>
	<p>en Do not resterilize de Nicht resterilisieren es No volver a esterilizar fr Ne pas restériliser it Non risterilizzare pt Não reesterilizar</p>
	<p>en Sterilized Using Steam de Mit Dampf sterilisiert es Esterilizar con vapor fr Stérilisé à la vapeur it Sterilizzato a vapore pt Esterilizado por vapor</p>
	<p>en Sterilized using ethylene oxide de Mit Ethylenoxid sterilisiert es Esterilizado con óxido de etileno fr Stérilisé à l'oxyde d'éthylène it Sterilizzato con ossido di etilene pt Esterilizado por óxido de etileno</p>
	<p>en Upper Limit of Temperature de Obere Temperaturgrenze es Límite superior de temperatura fr Limite supérieure de température it Limite superiore di temperatura pt Limite superior de temperatura</p>
	<p>en Non-sterile de Nicht steril es No estéril fr Non stérile it Non sterile pt Não esterilizado</p>
	<p>en Consult Instructions for Use de Gebrauchsanweisung beachten es Consultar las instrucciones de uso fr Consulter le mode d'emploi it Consultare le istruzioni per l'uso pt Consultar as instruções de utilização</p>
	<p>en Manufacturer de Hersteller es Fabricante fr Fabricant it Produttore pt Fabricante</p>
	<p>en Authorized Representative in the European Community de Autorisierte Vertretung in der EU es Representante autorizado en la Comunidad Europea fr Représentant autorisé dans la Communauté Européenne it Rappresentante autorizzato della comunità europea pt Representante autorizado na comunidade europeia</p>
	<p>en Recyclable packaging de Wiederverwertbare Verpackung es Embalaje reciclable fr Emballage recyclable it Confezione riciclabile pt Embalagem Reciclável</p>
	<p>en CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician. de VORSICHT: Nach US-amerikanischem Recht darf dieses Gerät nur von einem Arzt oder auf ärztliche Anordnung verkauft werden. es PRECAUCIÓN: las leyes federales de los EE. UU. limitan la venta de este dispositivo a médicos o bajo prescripción facultativa. fr MISE EN GARDE: Selon la législation fédérale américaine, ce dispositif ne peut être vendu que sur prescription d'un médecin. it ATTENZIONE: le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo ai medici o dietro prescrizione medica. pt ATENÇÃO: a lei federal dos Estados Unidos só permite a venda deste dispositivo mediante receita médica.</p>

# AMS 800™

## Urinary Control System For Female and Pediatric Patients

### Instructions For Use

**NOTE:** For implant procedure information, consult the AMS 800 Urinary Control System Operating Room Manual.

### Brief Device Description

The AMS 800 Urinary Control System is an implantable, fluid filled, solid silicone elastomer device used to treat urinary incontinence. It is designed to restore the natural process of urinary control. The device simulates normal sphincter function by opening and closing the urethra at the control of the patient. The AMS 800 consists of three interconnected components: a cuff, a pump, and a pressure-regulating balloon. The three components are connected with kink-resistant tubing. The AMS 800 can be implanted at the bladder neck.

The AMS 800 Urinary Control System Cuff and Pump are available with InhibiZone™, an antibiotic coating of rifampin (rifampicin) and minocycline.\*

A System implant contains less than 2% of the oral dose exposure during a complete course of rifampin and minocycline, even if two of the largest size Cuffs are implanted. The average amounts of rifampin and minocycline contained on a prosthesis implant are represented by the means and 95% tolerance intervals of the following implant configurations:

- 1.9 mg rifampin (0.7 – 3.1 mg) and 2.8 mg minocycline (2.1 – 3.5 mg) for the implant configuration with the lowest drug levels, i.e., single 4.0 cm cuff + control pump.
- 3.7 mg rifampin (0.9 – 6.5 mg) and 6.3 mg minocycline (4.7 – 8.0 mg) for the implant configuration with the highest drug levels, i.e., double 11.0 cm cuff + control pump.

### Indications for Use

The AMS 800 is used to treat urinary incontinence due to reduced urethral/bladder outlet resistance (intrinsic sphincter deficiency) in males, females and children. The information presented in this document focuses on implantation of the AMS 800 in females and children. For information relating to the implantation of the device in males, see the AMS 800 Urinary Control System Instructions for Use.

### Contraindications

1. This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
2. This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
3. This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
4. The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin or to minocycline or other tetracyclines.

\* not available in all markets

5. The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline has been reported to aggravate this condition.

## Warnings

1. Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection. Infection that fails to respond to antibiotic therapy may result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.
2. Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component misplacement. The cuff may erode around the urethra. The control pump may erode through the labia or scrotum. The pressure-regulating balloon may erode into the bladder. Acute urinary tract infection can interfere with proper functioning of the device and may lead to erosion of the urethra in the cuff area. Failure to evaluate and promptly treat the erosion may result in a substantial worsening of the condition leading to infection and/or loss of tissue.
3. Poor bladder compliance or a small fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
4. Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
5. Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device to prevent potential damage to the urethra or the AMS 800.
6. This device contains solid silicone elastomers. This device does not contain silicone gel. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
7. Surgical, physical, psychological, or mechanical complications, if they occur, may necessitate revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
8. Product wear, component disconnection or other mechanical problems may lead to surgical intervention. Mechanical complications may include malfunctioning of the components and leakage of fluid. Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in overflow obstruction. Mechanical events should be evaluated carefully by the treating physician and the patient should consider risks and benefits of treatment options, including revision surgery.
9. Previous patient history of adverse reaction(s) to radiopaque solution precludes its use as a filling medium for the prosthesis. Instead, saline should be used to fill the device.
10. Female patients with persistent incontinence should be evaluated to rule out vesicovaginal fistula, which may have resulted from an unrecognized iatrogenic injury.

11. If a hypersensitivity reaction develops to a device coated with InhibiZone, the Cuff and Pump should be removed and the patient treated appropriately.

## Precautions

### *Patient Related*

1. Patient selection requires thorough preoperative consultation and evaluation by the physician.
2. Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation of an AMS 800. Although the prosthesis is designed to restore urinary control, some patients continue to have a degree of incontinence after this procedure.
3. Patients may experience pain when the device is activated in the postoperative period and during the period of initial use. Cases of chronic pain associated with device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention. Patients should be counseled on expected postoperative pain including severity and duration.
4. Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff at the urethra.
5. Any progressively degenerative disease, e.g. multiple sclerosis, may limit the future usefulness of the implanted prosthesis as a treatment for the patient's urinary incontinence.
6. Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device.
7. Trauma or injury to the pelvic, perineal or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the device. The physician should advise patients of these possibilities and warn them to avoid trauma to these areas.
8. Female patients of child-bearing age must be forewarned that pregnancy is acceptable, but cesarean section may be indicated to minimize the risk of damage to the bladder neck and its surrounding cuff. For those patients who become pregnant, deactivation of the device during the third trimester is recommended to reduce the risk of erosion. Those patients contemplating pregnancy should consider delaying implantation.
9. Children receiving device should be evaluated at regular intervals. Lifelong radiological and urodynamic surveillance of the urinary tract is crucial. Before implantation, the patient and their family should be informed of the complication rate and the need for long-term follow-up.
10. Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5cm) generally exceeds 28F. Additional clearance is required to accommodate the patient's urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on sizing.

### *InhibiZone™ Related*

1. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin (rifampicin) and minocycline can cause additional stress on the hepatic and renal systems.
2. Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
3. Patients who receive a device with InhibiZone and are also taking Warfarin should have their prothrombin time monitored because tetracyclines have been reported to slow coagulation.
4. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses or rifampin (rifampicin).
5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.
7. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
8. Because products with InhibiZone are impregnated with a combination of rifampin (rifampicin) and minocycline, the contraindications, warnings and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline and rifampin (rifampicin) in patients receiving this device are unlikely to be detected.

### *Surgery Related*

1. Improper cuff sizing, improper balloon selection, or other causes may result in tissue erosion, migration of components, or continued incontinence.
2. Component migration can occur if the cuff is sized improperly, if the pump or balloon is not positioned correctly, or if the tubing lengths are incorrect. Migration can result in pain, complications, device malfunction and surgical revision.
3. Unsuccessful outcomes may result from improper surgical technique, improper sterile technique, anatomical misplacement of components, improper sizing and/or filling of components.
4. Although reinforced tubing has been designed to be more resistant to tubing kinks, tubing kinks may still result from tailoring the connecting tubing to an improper length during the implant procedure.

### *Device Related*

1. If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may arise as a result:
  - a. In the event of large pressures within the bladder, automatic pressure relief that normally occurs with

the device would be prevented. Cycling the device can relieve the outflow obstruction.

- b. Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the device, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.
  - c. Release of the deactivation valve may require greater pressure than that used to cycle the device.
2. System pressure changes may occur over time if you fill the balloon with radiopaque solution of incorrect concentration. Follow the instructions in the Operating Room Manual to prepare the radiopaque solution with the correct concentration.

## How Supplied and Storage

Store device in a clean, dry, dark area at room temperature.

**WARNING:** Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## Clinical Information – Females

The AMS 800 is indicated for the treatment of urinary incontinence in females and in children. A literature review was performed to evaluate the safety and efficacy of the AMS 800 in female patients. A total of seventeen documents, specifically reporting clinical information for female patients only, were analyzed. Clinical data, including patient demographics, etiology, complications, and success criteria and rates are reported below.

### Patients

Clinical efficacy has been reported for 938 female patients, ages ranging from 9-85, implanted with the AMS 800<sup>1-17</sup>. The mean reported follow-up on 785 AMS 800 patients was 3.6 years<sup>1-17</sup>. Various patient etiologies were reported, including neuropathic bladder dysfunction<sup>1,2,3, 4,13</sup>, stress incontinence<sup>1,4,6,7,8,9,12,15,16,17</sup>, pelvic/perineal trauma<sup>1,4,12,13</sup>, congenital defects (exstrophy)<sup>1</sup>, hysterectomy<sup>1</sup>, radiotherapy<sup>1</sup>, cystectomy<sup>1</sup>, nonneurogenic bladder<sup>2,13</sup>, intrinsic sphincter deficiency<sup>2,6</sup>, myelomeningocele<sup>13,14</sup>, spinal injury<sup>13</sup>, and cerebral palsy<sup>14</sup>. Prior surgery or radiotherapy were reported in 354 patients<sup>2,3,5,6,8,10,12,15,16</sup>.

### Complications

Table 1 summarizes the reported complications for 938 patients receiving an AMS 800 implant. The most common reported complication was bladder injury, which occurred in 5.8% of patients (n=54). There may have been more than one reported complication for each patient.

This trial did not include AMS 800 devices impregnated with InhibiZone.

**Table 1. Complications in Patients Studied**

Complication	Percentage	Reference
Bladder injury	5.8%	2,8,11,17
Vaginal injury	4.2%	2,4,8,11,13
Cuff Erosion	3.4%	3,6,7,10,13,17
Infection/erosion	3.3%	1
Device Fluid leak	2.8%	2,3,9,12,13,15,16,17
Incisional hernia	1.5%	2
Urinary retention	1.3%	2,16
Urethral injury	1.2%	2,4,11
Urinary urgency	1.1%	2,6
Labia hematoma	1.0%	2,14
Infection	0.6%	1,4,8,10,17
Pelvic abscess	0.4%	16
Early post-operative morbidity	0.4%	1
Pressure atrophy at cuff site	0.3%	17
Cutaneous erosion	0.3%	4
Bothered by device	0.3%	2
Pulmonary embolus	0.1%	1
Enterocutaneous fistula	0.1%	1
Phlebitis	0.1%	2
Upper urinary tract deterioration	0.1%	13
Superficial wound dehiscence	0.1%	16
Device malfunction	0.1%	1,2,13,16,17
Device displacement	0.1%	3
Device damage – unrelated procedure	0.1%	6
Wound seroma	0.1%	9
Pump erosion	0.1%	13
Sphincter-related problem	0.1%	10
Loose cuff	0.1%	15,16

### Success Rates

Patient success rates were reported and analyzed. For those defining clinical success<sup>2,4,5,6,7,12,13,16,17</sup>, rates of success were categorized into the following criteria:

- Fully and socially continent = 0-1 pad use per day
- Occasional incontinence = Consistent use of 1-2 pads per day
- Fair/improved continence = Use of 3 pads per day
- Not improved = more than 3 pads per day

Success rates are summarized in Table 2.

**Table 2. Success Rates**

Success Criteria	% patients
Fully and socially continent	86.9%
Occasional incontinence	4.4%
Fair/improved	4.3%
Not improved	7.1

### Patient Information Form (PIF) Study

A retrospective analysis was performed for female patients implanted (n=637) with the AMS 800 during the five-year period 2000- 2005. The study examined Patient Information Form (PIF) data sent to the manufacturer by the implanting physician for original implants and revision surgeries. A revision is a surgical intervention related to the function, placement, or site reaction to the implanted device. PIF data is required to be on file by the manufacturer in order to be eligible for product replacement. Of the 637 female patients implanted during the period studied, 76 patients had undergone a total of 88 revision surgeries (1.16 revisions per patient). There were a total of 121 reported reasons for revisions. Note that more than one reason may be listed for the same revision. The data in Table 3 presents the reported reasons for revision and the percentage of each reason for revision.

This trial did not include AMS 800 devices impregnated with InhibiZone.

**Table 3. PIF Study Revision Data**

Reason for revision	% revisions
Incontinence	19.8%
Not specified	19.8%
Infection	13.2%
Fluid loss	11.6%
Erosion	10.7%
Other medical condition	8.3%
Device out of specification/malfunction	5.8%
Improper size of device	3.3%
Other*	3.3%
Urinary retention	1.7%
Pain	1.7%

\* Other includes: hematoma, device migration/malposition, patient dissatisfaction.

### Clinical Information – Children

The AMS 800 is indicated for the treatment of urinary incontinence in children of both sexes. The safety and efficacy of the AMS 800 in children was evaluated by conducting a literature review. Four documents, reporting clinical information for patients under the age of 18, were analyzed and a summary of the clinical evidence is reported here.

Clinical efficacy has been reported for 322 children who have undergone implantation of the AMS 800. The mean age at the time of implant ranged from 9.9 to 14 years, and mean follow-up time reported was 5 to 15.4 years. Etiologies were similar to those reported for female patients, with the addition of sacral agenesis<sup>18,19,21</sup>, lipoma<sup>18,21</sup>, spinal cord abnormalities<sup>18,20,21</sup>, sacralcoccygeal teratoma<sup>18</sup> or cyst<sup>21</sup>, accidental injection of bicarbonate into the umbilical artery<sup>18</sup>, urethral injury<sup>20</sup>, idiopathic<sup>19</sup>, lipomeningocele<sup>21</sup>, posterior urethral valves<sup>21</sup>, interstitial cystitis<sup>21</sup>, congenital absence of bladder neck<sup>21</sup>, and ureter malfunction<sup>21</sup>.

Reported complications were similar to those reported for female patients. Device revisions and bladder augmentation were the most commonly reported reason for additional surgeries<sup>18,19,21</sup>. Reported success rates were comparable to those reported for female patients, with a reported range of 67% - 86%<sup>18,20,21</sup>.

This trial did not include AMS 800 devices impregnated with InhibiZone.

### Device Survival

Although it is not possible to predict exactly how long an implanted prosthesis will function in a particular patient, American Medical Systems, Inc. gathered data from two sources on device removals and revisions to help gain insight into product performance over time. Although these sources included information gathered from male patients, it is expected that similar data would result from analyzing data exclusively from female patients and children. Figure 1 presents device survival results from a prospective clinical study and a Bayesian analysis that uses data from a prospective clinical study and a PIF Study to estimate device survival at five years.

*Prospective Clinical Study* – A device survival curve was calculated from data collected during a prospective clinical study (n=85) with two-year follow-up. Using Kaplan-Meier analysis, the two-year revision-free rate for the AMS 800 was 79.5% (95% CI with 95% lower confidence bound 69.8%).

*Bayesian Analysis* – A Bayesian hierarchical model was used to evaluate device safety in the prospective clinical study. The Bayesian model estimated device survival using historical data (PIF Study n=12,713) on the AMS 800 and prospective clinical study data (n=85) on the AMS 800. A log-normal distribution fit the AMS historical data. Based on the log-normal model, it

was estimated that the five-year revision-free rate for the AMS 800 is approximately 73.8% with 95% CI ranging from 67.3% to 79.6%. The results met the primary safety endpoint for the clinical study of a five-year revision free rate at 75% using a 10% delta with two-sided 95% lower bound greater than 65%.


Figure 1: Estimated Survival Rates for the AMS 800

## Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological and functional outcome of the implantation. The risks, benefits and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal. Patients should also be aware that the AMS 800 is not considered to be a lifetime implant.

It is important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (See Silicone Information).

## Antibiotic Information

The antibiotics present in InhibiZone, minocycline and rifampin (rifampicin), are well characterized and have been in use for years. The dosage present on the artificial urinary sphincter is intended to act on organisms that attempt to colonize on the device. The average amounts of rifampin and minocycline contained on a prosthesis implant are represented by the means and 95% tolerance intervals of the following implant configurations:

- 1.9 mg rifampin (0.7 – 3.1 mg) and 2.8 mg minocycline (2.1 – 3.5 mg) for the implant configuration with the lowest drug levels, i.e., single 4.0 cm cuff + control pump.
- 3.7 mg rifampin (0.9 – 6.5 mg) and 6.3 mg minocycline (4.7 – 8.0 mg) for the implant configuration with the highest drug levels, i.e., double 11.0 cm cuff + control pump.

The following in vitro data are available, but their clinical significance is unknown. No clinical studies have been performed to evaluate the effect of the antibiotic surface treatment on reducing the incidence of artificial urinary sphincter implant infections.

**Table 4. In Vitro Zone of Inhibition for Device Samples\* with InhibiZone Treatment**

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
Staphylococcus epidermidis	22.6	2.9	21
Staphylococcus aureus	17.5	5.0	25
Escherichia coli**	6.5	2.6	24
Enterococcus faecalis**	4.8	6.7	21
Candida albicans**	0.1	0.4	21
Proteus mirabilis**	0.6	1.0	17

\* obtained using standardized KRT test samples containing approximately 12µg minocycline and 26 µg rifampin(rifampicin)

\*\* the isolates tested were not susceptible to rifampin (rifampicin) and/or minocycline control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and five rabbits were implanted subcutaneously with 6 control samples each. One rabbit received three test samples and three control samples. The test samples were portions of an InhibiZone treated AMS 700 pump and the control samples were portions of a standard AMS 700 pump without InhibiZone. The AMS 700 Pumps used in the 700 Series of Inflatable Penile Prosthesis is similar to AMS 800 Pumps used in the 800 Series of Artificial Urinary Sphincters in regards to material composition, adhesive, and InhibiZone application process. All samples were soaked in a  $10^3$ - $10^4$  CFU solution of staphylococcus aureus, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

### **Silicone Information**

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Solid silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

There are reports of malignant tumor formation in laboratory animals only associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS 800. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

## Magnetic Resonance Imaging (MRI) Important Safety Information

Non-clinical testing has demonstrated the AMS 800 product line is MR Conditional. The device can be scanned safely under the following conditions:

Static Magnetic Field	1.5 Tesla <sup>a</sup>	3.0 Tesla <sup>b</sup>
Spatial Gradient Field	450 Gauss/cm or less	720 Gauss/cm or less
Maximum whole body averaged Specific Absorption Rate (SAR)	1.5 W/kg for 15 minutes of scanning as assessed by calorimetry	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry

### MRI-Related Heating

Non-clinical testing has demonstrated the AMS 800 product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	1.5 Tesla <sup>a</sup>	3.0 Tesla <sup>b</sup>
Highest Temperature Change	≤ +0.4 °C	≤ +2.0 °C
(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5		

### Artifact Information

Non-clinical testing has demonstrated that the AMS 800 product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	5,800 mm <sup>2</sup>	1,956 mm <sup>2</sup>	6,096 mm <sup>2</sup>	2,650 mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

### Inventory Returns and Product Replacement Information

A Patient Information Form (PIF) must be filled out and filed with AMS at the time of implant to activate the product warranty. Contact your local AMS Representative prior to returning any product. Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form. Follow **all** of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS. Request an AMS Product Return Kit from the AMS Customer Service Department to return any explanted components to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Limited Warranty Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

This document is written for professional medical audiences. Contact American Medical Systems for lay publications.

American Medical Systems periodically updates product literature. If you have questions about the currency of this information, contact American Medical Systems.

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- 1 Thomas K, Venn SN, Mundy AR, *Outcome of the Artificial Urinary Sphincter in Female Patients*, Journal of Urology, Apr 2002, Vol. 167, pp. 1720-1722
- 2 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L, *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, Journal of Urology, Apr 2001, Vol. 165:4, pp. 1172-1176
- 3 Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H, *Explantations and Revisions with Artificial Urinary Sphincter AMS 800 in Females*, Journal of Urology, May 1996, Vol. 155:5, Supplement
- 4 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C, *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, Journal of Urology, 1996, Vol. 155:5, supplement
- 5 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H, *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, Acta Urol Bel, May 1995, Vol.63(2), pp. 45-46
- 6 Hadley R, Lioses P, Dickinson M, *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 7 Stone KT, Diokn AC, Mitchell BA, *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 8 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H, *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, Journal of Urology, May 1994, Vol. 151:5, supplement
- 9 Karram MM, Rosenzweig, BA, Bhatia, NN, *Artificial Urinary Sphincter for Recurrent/Severe Incontinence in Women: Urogynecologic Perspective*, Journal of Reproductive Medicine, Oct 1993, Vol. 38, pp. 791-794
- 10 Duncan HJ, Nurse DE, Mundy AR, *Role of the Artificial Urinary Sphincter in Women*, British Journal of Urology, Feb 1992, Vol. 69, pp. 141-143
- 11 Salisz JA, Diokno AC, *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, Journal of Urology, Nov 1992, Vol. 148, pp. 1528-1530
- 12 Webster GD, Perez LM, Khoury JM, Timmons SL, *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, Urology, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503
- 13 Parulkar BG, Barrett DM, *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, Gynecology and Obstetrics, Aug 1990, Vol. 171, pp. 131-138
- 14 Abbassian A, *A New Operation for the Insertion of the Artificial Urinary Sphincter*, Journal of Urology, Sept 1988, Vol. 140, pp. 512-513
- 15 Appell RA, *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, Neurourology and Urodynamics, 1998, Vol. 7, pp. 613-619
- 16 Diokno A., Hollander JB, Alderson TP, *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, Journal of Urology, Oct 1987, Vol. 138, pp. 778-780
- 17 Light JK, Scott FB, *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, Journal of Urology, Sept 1985, Vol. 134, pp. 476-478
- 18 Simeoni J, Guys JM, Mollard P, Buzelin JM, Moscovici J, Bondonny JM, Melin Y, Lortat-Jacob S, Aubert D, Costa F, Galifer B, Debeugny P, *Artificial Urinary Sphincter Implantation for Neurogenic Bladder: A Multi-institutional Study in 107 Children*, British Journal of Urology, 1996, Vol. 78, pp. 287-293
- 19 Kryger JV, Barthold JS, Fleming P, Gonzalez R, *The Outcome of Artificial Sphincter Placement After a 15- Year Follow-up in a Paediatric Population*, British Journal of Urology International, 1999, Vol. 83, pp. 1026-1031
- 20 Castera R, Podesta ML, Ruarte A, Herrera M, Medel R, *10-Year Experience with Artificial Urinary Sphincter in Children and Adolescents*, Journal of Urology, June 2001, Vol. 165, pp. 2373-2376

- 21 Herndon, CDA, Rink RC, Shaw MBK, Simmons GR, Cain MP, Keafer M, Casale AJ, *The Indiana Experience with Artificial Sphincters in Children and Young Adults*, Journal of Urology, Feb 2003, Vol. 169, pp. 650-654
- 22 Shellock F, *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature*, AJR, October 1988
- 23 Shellock F, *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review*, Radiology, 1991, Vol. 180, pp. 541-550
- 24 Shellock F, *MR Procedures and Biomedical Implants, Materials and Devices: 1993 Update*, Radiology, 1993, Vol. 189, pp. 587-599
- 25 Shellock F, *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

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All trademarks are the property of the respective owners.

# AMS<sup>TM</sup>

**American Medical Systems, Inc.**

10700 Bren Road West  
Minnetonka, MN 55343  
U.S.A.  
U.S. Toll Free: 800 328 3881  
Tel: +1 952 930 6000  
Tel: +31 20 593 8800

STERILE EO

STERILE

**American Medical Systems  
Europe B.V.**

Haarlerbergweg 23 G  
1101 CH Amsterdam Zuid-Oost  
The Netherlands

©2017 Boston Scientific Corporation or its affiliates.  
All rights reserved.  
1003276 (P/N)  
1003182 (A/W Rev A) (2017-02)

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Boston Scientific (Australia) Pty Ltd