
SYNERGY MEGATRON™

Everolimus-Eluting Platinum Chromium Coronary Stent System

**Purpose-built
for Large Vessels**

SYNERGY MEGATRON™ & SYNERGY™ XD LV Model


SYNERGY
MEGATRON™
Everolimus-Eluting Platinum Chromium Coronary Stent System

STENT Design

Purpose built stent architecture SYNERGY MEGATRON™

12 Crown Design with Shorter Strut Length


Platinum Chromium(PtCr) Alloy

4 Connectors on Proximal Two Segments

3 Connectors Throughout the Body

Strut Thickness and Width

SYNERGY MEGATRON™
BP Stent


89 μm

92 μm

SYNERGY™ XD : LV Model

10 Crown Design with Shorter Strut Length¹


Platinum Chromium(PtCr) Alloy

5 Connectors on Proximal Two Segments¹

2 Connectors Throughout the Body

Strut Thickness and Width

SYNERGY™ XD BP Stent
Large Vessel Sizes

81 μm

94 μm

※SYNERGY™ XD 4.0mm 径

REFERENCE 1. The Small Vessel and Workhorse sizes have 4 connectors on the two most proximal segments. The Large Vessel sizes have 5 connectors on the two most proximal segments.

Purpose Built for Large Proximal Vessels

SYNERGY
MEGATRON™
Percutaneous-Deployable Platinum-Chromium Coronary Stent System

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

Uniform Lesion

Scaffolding To Reduce tissue prolapse

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
PtCr for

Visibility to aid accurate stent placement

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


Radial Compression Resistance Test

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
Stent

REFERENCE 1. Bench tests performed by Boston Scientific Corporation. Data on file. 1. All stents 4.0mm. N=3 minimum. Bench test and computational model results may not necessarily be indicative of clinical performance.

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Axial Compression Resistance Test

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Uniform Lesion Scaffolding with a 12-Crown Design

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
5.0mm



3.5mm


Minimal tissue prolapse (smoother edges)

0.20 Max Displacement

Visible tissue prolapse ("pillowing")

0.30 Max Displacement


9-Crown Prototype DES

Purpose Built for Large Proximal Vessels

Maximum Expanded Cell Diameter¹

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Side Cell Expansion


5.0mm

Maximum expanded cell diameter (MECD)¹

Additional connector on the Proximal 2 Segments for Strength

3.8mm MECD in the first two proximal segments¹

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
REFERENCE 1. Based on 5.0mm diameter. Testing Completed by Boston Scientific data on file. Computational model results may not necessarily be indicative of clinical performance.

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





REFERENCE 1. Based on 3.5mm diameter. Testing Completed by Boston Scientific data on file. Bench test and computational model results may not necessarily be indicative of clinical performance.

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


Bench test and computational model results may not necessarily be indicative of clinical performance.

Purpose Built for Large Proximal Vessels

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


Cardiac History & Lesion Characteristics

Primary Endpoint : 1 year TLF

10.8% 95% CI [5.3%, 18.9%]

TLF: composite of ischemia-driven TLR, MI related to the target vessel, or cardiac death;
Denominator includes patients with sufficient follow-up or having a TLF event

Cardiac Death and MI at 1 year

Cardiac Death

6.5% 95% CI [2.4%, 13.5%]


Six cardiac-related deaths unrelated to the study stent;
None of these were related to stent thrombosis

Myocardial Infarction

6.5% 95% CI [2.4%, 13.5%]

Six subjects experienced target vessel related MI,
3 of which were periprocedural MI

Additional Outcomes at 1 year


SYNERGY Product Family Specifications

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SYNERGY MEGATRON™ Stent System Key SDS Specifications

Lesion entry profile	0.017 in (0.43 mm)
Average stent profile	3.50mm 0.049" (1.24mm), 4.00mm 0.049" (1.24mm), 4.50mm 0.053" (1.34mm), 5.00mm 0.058" (1.47mm)
Total catheter length	144 cm
Distal shaft coating	Bioslide™ hydrophilic coating
Hypotube coating	PTFE
SDS shaft outer diameter (OD)	≤0.70 mm (2.1 F) proximally and 2.6F or 2.7F** (≤ 0.92 mm) distally
Guide catheter compatibility	5F (≥0.056" / 1.42mm) for 3.50 - 4.00 mm sizes
	6F (≥0.070" / 1.78 mm) for 4.50 - 5.00 mm sizes
Delivery system ports	Single access port to inflation lumen.
	Guidewire exit port is located approximately 25 cm from tip.
	Designed for guidewire ≤ 0.014 in (0.36 mm)
Balloon material	PEBAX®
Balloon fold	5 wings
Marker band material	Platinum Iridium
Marker band length	1 mm
Marker band placement	Nominally spaced 0.4 mm away from the end of the stent

*3.50 mm device. **Dependent on diameter and length. Boston Scientific data on file.

PEBAX is registered trademarks of Arkema.

Purpose Built for Large Proximal Vessels

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Nominal pressure
Rated Burst pressure*

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製品の詳細に関しては添付文書等でご確認いただくか、弊社営業担当へご確認ください。
©2023 Boston Scientific Corporation or its affiliates. All rights reserved.
All trademarks are the property of their respective owners.

**Boston
Scientific**
Advancing science for life™

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