Advancing science for life™

GUIDE™ DBS Software Programming Manual

Directions for Use

Directions for Use

en

Mode d'emploi

fr

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es

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ptBR

91062299-07

Content: 91072429 REV AB

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

Trademarks

All trademarks are the property of their respective holders.

Additional Information

For indications for use, contraindications, warnings, precautions, adverse events, storage and handling, sterilization, labeling symbols, limited warranty, or other device-specific information not included in this manual, refer to the appropriate DFU for your DBS System as listed on your DBS *Reference Guide*.

Product Model Numbers

Model Numbers	Description
DB-7154	GUIDE™ DBS Station
DB-7154R	GUIDE DBS Station (R stands for Refurbished Products)

Glossary

Anterior Commissure

A bundle of nerve fibers (white matter) connecting the two cerebral hemispheres across the midline.

Atlas

See Morel Atlas.

Axial Plane

A horizontal plane dividing the brain into top and bottom sections.

Coronal Plane

A vertical plane dividing the brain into anterior and posterior sections.

DICOM

Digital Imaging and Communications in Medicine (DICOM) is a standard for handling, storing, printing, and transmitting information in medical imaging. It includes a file format definition and a network communications protocol. DICOM is known as NEMA Standard PS3 and as ISO Standard 12052.

Encapsulation Layer

A thin layer of glial scar tissue which forms around a foreign body, in this case, the implanted lead.

Fusing/Fusion

The process of aligning one or more DICOM images or a model (for example, atlas) to one another so that they can both be referred to in the same coordinate system.

GUIDE™ DBS

Refers to GUIDE™ DBS Software.

Hemisphere

A cerebral hemisphere is one of the two regions of the brain that are delineated by the median plane (medial longitudinal fissure). The brain can thus be described as being divided into left and right cerebral hemispheres.

IPG

Implantable Pulse Generator.

Level Control

The medical imaging term for contrast control. See Leveling.

Leveling

The process of defining the window/level (brightness/contrast) of a DICOM image.

Mid-Commissural Point (MCP)

The midpoint along the line between the Anterior Commissure (AC) and Posterior Commissure (PC) points.

Mid-Sagittal Plane

A vertical plane through the midline of the brain, dividing the brain into right and left hemispheres.

Morel Atlas

A collection of three-dimensional neuroanatomical volumes rendered from histological slicing of six cadaver brains. Work was done by Anne Morel et al.¹

Posterior Commissure (PC)

A bundle of nerve fibers (white matter) connecting the two cerebral hemispheres across the midline.

Reference Volume

The Reference Volume is designed to visually indicate a volume of tissue which, if stimulated, has been shown in peer-reviewed literature to give favorable clinical outcomes with minimal clinical side effects.²

Registration Process

The steps performed to register a Patient's images with the implanted Lead and Morel Atlas.

Sagittal Plane

The vertical plane that passes from front to rear dividing the brain into right and left sections.

Stimulation Field Model (SFM)

Displays the volume of tissue surrounding the Lead that is likely to be activated for the given stimulation parameters. This volume is based on neuronal and finite element modeling.

Transformation

A mapping from the coordinate system of one image into the coordinate system of another image.

Window Control

The medical imaging term for brightness control. See Leveling.

¹ Krauth A, Blanc R, Poveda A, Jeanmonod D, Morel A, Székely G. A mean three-dimensional atlas of the human thalamus: generation from multiple histological data. *Neuroimage*. 2010 Feb 1;49(3):2053-62.

² The Reference Volume (RV) was derived by Dr. Cameron McIntyre and Dr. Chris Butson specifically for DBS treatment of Parkinson's disease when the Lead is implanted in the neighborhood of the subthalamic nucleus. The RV is designed to visually indicate a volume of tissue which, if stimulated, has been shown in peer-reviewed literature to give favorable clinical outcomes with minimal clinical side effects. Dr. McIntyre and Dr. Butson integrated the RV into their custom Cicerone software package (Miocinovic, 2007), and its use for programming human subjects has been documented in the peer reviewed literature (Frankemolle et al., 2010). Subsequently, the RV has been registered by Dr. Butson to the anatomical atlas that will be used in the GUIDE™ DBS software.

Miocinovic, S., Noecker, A.M., Maks, C.B., Butson, C.R., McIntyre, C.C., 2007. Cicerone: stereotactic neurophysiological recording and deep brain stimulation electrode placement software system. *Acta Neurochir Suppl.* 97, 561-7.

en

This page intentionally left blank

Table of Contents

Glossary iii

Introduction..... 1

Creating a Visual Model 2

 Creating a New Patient Record 2

 Registering the Patient’s MR and CT..... 3

 Configuring a Stimulation Field Model 5

Indications for Use 6

Device Description and Intended Use 6

GUIDE™ DBS Kit Contents* 6

Software 6

 ⚠ **Warnings** 7

Setting Up and Shutting Down..... 8

Setting Up GUIDE™ DBS..... 8

Creating a New User..... 9

Setting Up the Mouse..... 9

 Troubleshooting Mouse Pairing 10

Connecting the Optical Drive 11

Setting Up the Flash Drive..... 11

Launching GUIDE™ DBS..... 12

Exiting GUIDE™ DBS 14

Setting the Time Zone 15

Monitoring the Battery Level 16

Restoring the Software Display 16

Shutting Down GUIDE™ DBS..... 16

Patient Browser 17

Adding a Patient 19

Ending a Session..... 21

Editing a Patient 22

Deleting a Patient 23

Viewing Patient Information 24

Searching for a Patient 25

Sorting the Patient List 25

Exporting Patient Files..... 26

Importing Patient Files..... 28

en

Image Registration 30

- Optimizing Registration 30**
 - Minimum Scan Requirements..... 30
 - Recommended Scan Requirements..... 30
 - Rescanning..... 31
- Image Registration Process Flow 33**
- Loading MR Series 34**
- Loading Post-Op CT Series 37**
- Verifying Fusion of MR and CT Scans..... 40**
 - Using the Image Registration Tools 43
 - Selecting a Registration View Mode 47
 - Using the Viewports..... 52
- Selecting the Mid-Sagittal Plane 54**
- Selecting the Anterior Commissure (AC) and Posterior Commissure (PC)..... 57**
 - Using the AC/PC Markers..... 58
- Verifying/Locating the Lead Tips 60**
 - Verifying the Lead Tips 60
 - Locating the Lead Tips Manually 64
- Verifying/Locating the Lead Shafts..... 65**
 - Verifying the Lead Shafts..... 65
 - Locating the Lead Shafts Manually..... 67
- Previewing Left Lead Registration..... 68**
- Previewing Right Lead Registration 70**
- Ending a Registration Session..... 72**

Programming Simulation 73

- Creating Programming Simulations 73**
- Using the Simulation Panel 78**
 - Adding a Simulation Slot..... 78
 - Deleting a Simulation Slot..... 78
 - Duplicating a Simulation Slot..... 79
 - Configuring the Programming of a Simulation Slot..... 80
- Using the 2D Parameters Panel 81**
 - Programming the Contact Configuration 81
 - Programming the Pulse Width..... 87
 - Programming the Rate..... 87
 -  Charge Density Warning..... 88
 -  Display Limit Reached Warning..... 89
- Using the 3D Visualization Panel 90**
- Creating a Simulation Report..... 92**
- Exporting a Simulation Report..... 93**


en

Introduction

GUIDE™ DBS Software (GUIDE DBS) is a simulation system that allows you to plan the programming of a Patient with a Boston Scientific Deep Brain Stimulation (DBS) System. GUIDE DBS provides:

- A view of the Lead(s) relative to the deep brain structures associated with treatment.
- The tools for simulating neuron response to stimulation.
- A visual model of the stimulation field in the context of an anatomical atlas.

Using GUIDE DBS, you can determine a start point for the monopolar review that may reduce programming time.

CAUTION: Screens depicted in this manual may differ slightly from the screens on your GUIDEDBS.

Creating a Visual Model

To create a visual model of the stimulation field in the context of a Patient's anatomical atlas, you will use GUIDE™ DBS to do the following:

- Create a New Patient record.
- Register the Patient's Pre-Op Magnetic Resonance (MR) and Post-Op Computerized Tomography (CT) images.
- Configure one or more stimulation settings for each DBS Lead.

Each of these steps is described in detail in this manual. The following is a brief overview to introduce you to the process.

Creating a New Patient Record

To create a new Patient Record, from the **Patient Browser** screen you will select  **New** and enter the following Patient information and hardware information on the form shown below:

- **Patient Name (Last, First), Patient ID, and Date of Birth**
- **IPG (Implantable Pulse Generator) Date of Implant (optional)**
- For each Lead implanted, the **Lead Model, Lead Location (Left or Right), Electrode #'s (1-8 or 9-16), and Date of Implant (optional)**

Fields marked with an "*" are required.

* Patient Name

Last

First

* Patient ID

* Date of Birth

DD

MM

YYYY

+ * IPG Model

Vercise™

Date of Implant

DD

MM

YYYY

Lead 1

* Lead Model

1x8

* Lead Location

Select one...

* Electrode #'s

E1-E8

Date of Implant

DD

MM

YYYY

+ Lead 2

* Lead Model

Select one...

* Lead Location

Select one...

* Electrode #'s

Select one...

Date of Implant

DD

MM


YYYY

GUIDE™ DBS Software Programming Manual
91062299-07 2 of 407


Registering the Patient's MR and CT

Once the Patient record is created, you will align the Patient's Pre-Op MR and Post-Op CT, align the MR with the atlas, and locate the Lead(s) in the CT. The Image Registration process is summarized below:


1. Select and load the appropriate Pre-Op MR.
2. Select and load the appropriate Post-Op CT.
3. Verify the automatic fusion of the MR and CT or manually adjust the image alignment.


4. Align the atlas to MR by identifying the Mid-Sagittal Plane (orange line) in the MR.

5. Finish the alignment by selecting the  Anterior and  Posterior Commissures in the MR.


6. Verify the automatic location of the Lead Tip(s) in the CT  or manually locate the Lead Tip(s).


7. Verify the automatic location of the Lead Shaft(s) in the CT or manually locate the Lead Shaft(s).

8. Preview the Lead(s) in the atlas, in 3D, to verify the location(s).

Configuring a Stimulation Field Model

Once the DBS Lead(s) has been located, you will use GUIDE™ DBS to visualize the Stimulation Field Model (SFM) for the Patient's therapy. Up to four SFMs can be stored for each Lead.

CAUTION: GUIDE DBS is a simulation tool and should not be a substitute for standard clinical practice.

Indications for Use

GUIDE DBS is intended to display medical images and simulate stimulation output. It includes functions for image manipulation and 3D visualization (reconstructions and volume rendering). Features include display of a simulated DBS Lead from a Patient's CT scan compared with an anatomical atlas.


Device Description and Intended Use

GUIDE DBS consists of a GUIDE DBS Station with a touch pad, Mouse, keyboard, USB Optical Drive capable of reading DICOM files from removable optical media, and Flash Drive. The system is intended to be used in an office environment as a planning tool for the programming of the Boston Scientific DBS system.

GUIDE™ DBS Kit Contents*

- (1) GUIDE DBS Station with Power Supply
- (1) Optical Drive
- (1) Mouse
- (1) Flash Drive
- (1) Manual

* Contents packaged and labeled separately.

Software

GUIDE DBS uses a Patient's post-operative imaging to display the simulated DBS Lead(s) relative to an anatomical atlas. GUIDE DBS reads Magnetic Resonance (MR) and Computerized Tomography (CT) data from any removable media and loads selected images into a Patient profile. Using the Patient's imaging, GUIDE DBS allows functionality for the following:

- 3D visualization of the relative position of the DBS lead in an anatomical atlas.
- 3D visualization of the theoretical neural activation field.
- 3D visualization of the reference target stimulation volume.
- 3D visualization of the theoretical neural activation field relative to the reference target stimulation volume.

Note: *GUIDE DBS displays a warning when the current density of the output goes above the allowed warning threshold.*

GUIDE™ DBS is only compatible with the Boston Scientific Deep Brain Stimulation (DBS) System components and parameters. Compatible leads are as follows:

Model Number	Product/Description
DB-2201-xx-C	Vercise Lead Kit, 30cm
DB-2201-xxAC	Vercise Lead Kit, 30cm
DB-2201-xx-C	Vercise Lead Kit, 45cm
DB-2201-xxBC	Vercise Lead Kit, 45cm


For more information on programming the Boston Scientific DBS System, refer to the appropriate Directions for Use as listed on your Reference Guide.

 Warnings

- GUIDE DBS is to be used with Boston Scientific products only.
- Do not install any other software on GUIDE DBS.
- Do not change the GUIDE DBS computer hardware or software configuration.
- If you have a visual impairment, we advise against your using GUIDE DBS as a tool with your standard clinical care.


Setting Up and Shutting Down

GUIDE DBS Station

Note: The above illustration is a sample set-up. The USB port on the right is the 3.0 port and may be faster for data import/export.

Setting Up GUIDE™ DBS

1. Attach the power cable to the port on the right side of GUIDE DBS.
2. Slide the prongs on the power cable to the open position.
3. Plug the prongs into a power outlet or appropriate adapter.

Creating a New User

1. Power on GUIDE DBS.
2. Select "ClinicUser".
3. The following message will display:

The User's password must be changed before logging on the first time.

4. Select  to change the User's password.
5. Create a password. The minimum length is eight characters.
6. Leave the existing password section blank.

Flash Drive

Setting Up the Mouse

1. Turn off any other Bluetooth-enabled mice. During setup, the receiver may detect such devices.
2. Insert two AA batteries into the Mouse. The battery compartment is on the bottom of the Mouse.
3. Turn on the Mouse. The on/off button is on the bottom of the Mouse. Slide the button to the Green (On) position.
4. In the Windows taskbar at the bottom of the screen, select Start.
5. Select **Devices and Printers**.
6. In the Devices and Printers window that appears, select **Add a device**.
7. Press and *hold* the Blue Bluetooth button on the bottom of the Mouse for 2 seconds. The Blue LED on the top of the Mouse should start blinking.
8. In the **Add a device window**, wait for a minimum of 20 seconds until the Mouse appears in the list. Select the Mouse, and select **Next**.
9. GUIDE DBS should connect with the Mouse.
10. If the Mouse has been successfully added to GUIDE DBS, a verification message will appear. Select **Close**.

11. If the Mouse was not successfully added and the Mouse is listed in the devices window, delete the Mouse entry and repeat the procedure.

Note: For more information on the Mouse, refer to the Mouse instruction booklet.


Troubleshooting Mouse Pairing

Description	Action
Message: “Adding this device to this computer failed. It took too long to add this device.”	Select “Try again.”
Message: “Select a pairing option.”	Choose: “Pair without using a code.”
Message: “Adding this device failed because of an authentication error.”	Message appears if you choose the wrong pairing option or wait too long to select an option. Select “Try again.”
Mouse is connected but does not work.	Uninstall the Mouse and reinstall it. 1. From the Devices and Printers menu, select the Mouse. 2. Select “Remove device.” 3. Select “Next.” 4. When prompted to enter an Admin password, press the blue button on the bottom of the Mouse.

Connecting the Optical Drive

The Optical Drive is used to import Patient MR and CT data. The Optical Drive is read-only.

Note: You will need the Flash Drive to export Patient data and Patient reports.

To set up the Optical Drive, connect the Optical Drive cable to the Optical Drive and to a USB port on the GUIDE™ DBS Station.


To open the Optical Drive:

1. Slide the Lever on the front to the right. The lid will open automatically.
2. Insert the removable media into the Optical Drive and close the lid.

Note: Refer to the manufacturer's instructions for more information on the Optical Drive.

Setting Up the Flash Drive

The Flash Drive is used to export Patient data and Patient reports.

1. Remove the cap from the end of the Flash Drive.
2. Slide the Flash Drive into the USB port on the right side of GUIDE DBS.

Note: The USB port on the right is the 3.0 port and will be faster for data import/export.

Launching GUIDE™ DBS

1. Double select the GUIDE DBS icon on the desktop to launch GUIDE DBS.


2. The GUIDE DBS Splash screen displays the following information:
 - Software (**SW**) Version
 - **Model Number**
 - **Software Product Number**

Note: The GUIDE DBS Splash screen can be accessed from any screen by selecting

3. After the Splash screen, the **Patient Browser** screen displays. See “Patient Browser” on page 17.


The screenshot shows the 'Patient Browser' interface. At the top, there is a 'Menu' button and a title bar 'Patient Browser'. Below the title bar, there are buttons for '+ New', 'Edit', 'Delete', 'Register', and 'Program'. A search bar is located on the right side of the top bar. The main area contains a table with the following columns: 'PATIENT NAME', 'PATIENT ID', 'DATE OF BIRTH', 'LEAD 1 IMPLANT DATE', 'LEAD 2 IMPLANT DATE', 'REGISTRATION', and 'LAST MODIFIED'. The table lists ten patients, with the first row, 'Chung, Kristina H.', highlighted in purple. The 'Boston Scientific' logo and a battery status indicator (50%) are visible in the top right corner.

PATIENT NAME	PATIENT ID	DATE OF BIRTH	LEAD 1 IMPLANT DATE	LEAD 2 IMPLANT DATE	REGISTRATION	LAST MODIFIED
Chung, Kristina H.	12345678	04 JAN 1957	09APR 09	09 APR 09	4. Verify Fusion	09 MAR 09
Feng, Jane	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	5. Locate Mid-Sagittal Plane	08 MAR 09
Kenton, Paige H.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	7. Locate Lead Tip	04 MAR 09
Melton, Sherri E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	8. Locate Lead Shaft	04 MAR 09
Plunkett, Paul	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	2. Load MR	06 MAR 09
Weissenburgh, Sullivan	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	4. Verify Fusion	08 MAR 09
Wikowsky, Karen U.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09
Williamson, Patrick O.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09
Willow, Elsie A.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	3. Load CT	08 MAR 09
Zachson, Hazel E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	1. Skipped Registration	08 MAR 09
Zuruchy, Malcolm A.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09

Exiting GUIDE™ DBS


To exit Guide DBS at any time:

- 1. From the **Patient Browser** screen select  > **Exit**. See “Patient Browser” on page 17.

- 2. The following message displays:

Are you sure you want to exit GUIDE DBS?

- 3. Select  to exit the GUIDE DBS software.

Or

Select  to return to the **Patient Browser** screen.

Note: To return to the **Patient Browser** Screen from a Patient Session, see “Ending a Session” on page 21. To return to the **Patient Browser** screen during Image Registration, see “Ending a Registration Session” on page 72.


Setting the Time Zone

At the top right of all screens, GUIDE™ DBS displays the Station's battery level, date, and time.


The time zone can be changed to the local time zone if needed.

1. Select  **Menu** > **Exit**, to exit GUIDE DBS. See “Exiting GUIDE™ DBS” on page 14.
2. Select the time and date in the lower right corner of the screen to open the **Date and Time** window.

3. Select “**Change date and time settings...**”.

4. Select “**Change time zone...**”.
5. Customize the time zone and press **OK**.
6. Close the **Date and Time** window.
7. Double select the GUIDE DBS icon on the desktop to launch GUIDE DBS.

Monitoring the Battery Level

A sampling of the battery level and charging status icons is shown below:

-  Critical
-  25%
-  100%
-  Critical, Charging
-  25%, Charging
-  100%, Charging

Note: Low battery level may affect screen brightness. Ensure GUIDE DBS has adequate power during use.

Restoring the Software Display

Occasionally, when GUIDE DBS is launched and left unattended, the software display may change. To restore the software display:

1. Select  **Menu** > **Exit**, to exit GUIDE DBS. See “Exiting GUIDE™ DBS” on page 14.
2. Double select the GUIDE DBS icon on the desktop to re-launch GUIDE DBS.

Shutting Down GUIDE™ DBS

To shut down GUIDE DBS:

1. Exit GUIDE DBS. See “Exiting GUIDE™ DBS” on page 14.
2. Power down GUIDE DBS.
3. Turn off the Mouse. Slide the on/off button to the Red (Off) position.
4. Unplug and stow GUIDE DBS along with the Mouse, Optical Drive, Flash Drive, and power supply.

Patient Browser

The **Patient Browser** is for managing Patient data for Image Registration and Program Simulation. From the **Patient Browser** you can:

- Add new Patients.
- Enter Patient information and hardware information (Implantable Pulse Generator (IPG) and Lead(s)).
- Modify existing Patient data.
- Delete Patients.
- Import/Export Patient Files.
- Select Image Data Registration or Skip Registration.
- Start Image Registration.
- Start Program Simulation.

The **Patient Browser** screen initially displays existing Patients sorted by **Patient Name**. Patients can be sorted by any column by selecting the column header. A scrollbar appears on the right when more than 11 Patients are in the Patient List.

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

Menu

New

Edit

Delete

Register

Program



Search...

Patient Selection Header Row

Patient Selection List


PATIENT NAME	PATIENT ID	DATE OF BIRTH	LEAD 1 IMPLANT DATE	LEAD 2 IMPLANT DATE	REGISTRATION	LAST MODIFIED
Chung, Kristina H.	12345678	04 JAN 1957	09 APR 09	09 APR 09	4. Verify Fusion	09 MAR 09
Feng, Jane	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	5. Locate Mid-Sagittal Plane	08 MAR 09
Kenton, Paige H.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	7. Locate Lead Tip	04 MAR 09
Melton, Sherri E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	8. Locate Lead Shaft	04 MAR 09
Plunkett, Paul	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	2. Load MR	06 MAR 09
Weissenburgh, Sullivan	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	4. Verify Fusion	08 MAR 09
Wikowsky, Karen U.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09
Williamson, Patrick O.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09
Willow, Elsie A.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	3. Load CT	08 MAR 09
Zachson, Hazel E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	1. Skipped Registration	08 MAR 09
Zuruchy, Malcolm A.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	9. Complete	08 MAR 09

The components on the **Patient Browser** screen are described in the table below.

Patient Browser Screen		
	Component	Description
1		Allows you access to the following Menu of options: ● Patient Information (if a Patient is selected; otherwise disabled) ● Export Patient Files ● Import Patient Files ● About ● Exit
2		Allows you to enter New Patient information.
3		Allows you to Edit an existing Patient's information.
4		Allows you to Delete a Patient.
5		Allows you to Register a Patient.
6		Allows you to create Program simulations for a selected Patient.
7		Allows you to Search for a specific Patient.
8	Patient Selection List	Display of existing Patients.
9	Patient Selection Header Row	Headings for Patient data as entered by the user. Select a column header to sort the Patient List by that column.
10	Patient Name	Last Name, First Name of Patient.
11	Patient ID	ID assigned by the user to the Patient.
12	Date of Birth	Patient's Date of Birth.
13	Lead 1 Implant Date	Date of implant for Lead 1.
14	Lead 2 Implant Date	Date of implant for Lead 2.
15	Registration	Status of Current Image Registration step (for example, Complete , or Verify Fusion).
16	Last Modified	Date on which the Patient data (information, registration, or simulations) were last modified.
17	 Scrollbar	Allows you to scroll through the Patient List. It appears on the right of the screen when more than 11 Patients are in the Patient List.

Adding a Patient

1. Select . The **Enter Patient Information** screen displays.

Note: An Asterisk (*) denotes a required field. Press **Tab** on the keyboard to navigate to the next field. Press **Shift + Tab** to navigate in reverse.

2. Enter the **Patient Name (Last, First)**.
3. Enter the **Patient ID** (alphanumeric).
4. Enter the **Date of Birth** in DD MMM YYYY format, using the drop-down list for month.

Note: To unselect a month, select the “blank” option at the top of the drop-down list.

5. In the **IPG Panel**, enter the **Date of Implant** in DD MMM YYYY format. (Optional)
6. In the **Lead 1** panel:
- Select **Lead Location (Left or Right)**.
 - Select **Electrode #'s (E1-E8 or E9-E16)**.
 - Enter the **Date of Implant** in DD MMM YYYY format. (Optional)

7. If there is a second Lead, select  to enable the **Lead 2** panel.

Note: Selecting  will close the **Lead 2** Panel.

Or

If there is only one Lead, continue to the **Select or Skip Registration** panel.

en

- en
8. In the **Lead 2** panel:

• Select **Lead Location**. (Selection of **Right** or **Left** is linked to selection for **Lead 1**. When **Left** is selected for Lead 1, **Right** is selected by default for Lead 2, and vice versa.)

• Select **Electrode #'s**. (Selection of **E1-E8** or **E9-E16** is linked to selection for Lead 1. When **E1-E8** is selected for Lead 1, **E9-E16** is selected by default for Lead 2, and vice versa.)


• Enter the **Date of Implant** in DD MMM YYYY format. (Optional)

9. In the **Select or Skip Registration** panel (right side of screen), select one of the following buttons:
- | | |
|--|--|
| <div><div>Image Data Registration
Requires:
Pre-op MR
Post-op CT</div></div> | Image Data Registration button
Requires:
Pre-op MR and Post-op CT |
| <div><div>Skip Registration
Warning:
Atlas will not be displayed during
Program Simulation</div></div> | Skip Registration button
Warning:
Atlas will not be displayed during Program Simulation |
10. If you selected the **Image Data Registration** button, select **Next** at the top right of the screen. The **Image Registration** screen displays. Refer to the Image Registration section of this manual for further instructions.
- Note:** If **Next** is disabled, verify that all required fields are properly completed.
- Or
- If you selected the **Skip Registration** button, because MR or CT images are not available, select **Finish** at the top right of the screen. The **Patient Browser** screen displays.
- Note:** If **Finish** is disabled, verify that all required fields are properly completed.
- GUIDE™ DBS Software Programming Manual
91062299-07 20 of 407

Ending a Session

1. To return to the **Patient Browser** screen, select  > **End Session**.
2. If all required data have *not* been entered, the following message will appear:

Do you want to exit without adding the patient? Not all required fields have data so the patient cannot be added.


Select  to exit without adding the Patient.

Or

Select  to stay in the session.

3. If all required data have been entered, the following message will appear:

Would you like to save changes before ending the session?

Select  to save changes.

Or

Select  to end the session without saving changes.


Or

Select  to stay in the session.

Editing a Patient

1. Select a Patient on the **Patient Browser** screen.
2. Select .
3. The following message will display:


Changes made to the lead configuration may affect registration and simulations.


Press .

4. The **Enter Patient Information** screen displays. Edit the information as desired.
5. Select  to return to the **Patient Browser** screen.

Or

Select  > **End Session**. The following message will appear:

Would you like to save changes before ending the session?

Select  to save changes.

Or

Select  to end the session without saving changes.


Or

Select  to stay in the session.

Deleting a Patient

1. Select a Patient on the **Patient Browser** screen.
2. Select  **Delete**. The following message appears:

Are you sure you want to delete this patient?

3. Select  to delete the Patient.
- Or
- Select  to cancel the action and not delete the Patient.

Viewing Patient Information

1. On the **Patient Browser** screen, select a Patient, then press  **Menu** .
2. Select **Patient Information** from the drop-down list. The **Patient Information** window displays.

Patient Information

Patient information for John Doe.

Patient Information			
Patient Name	Patient ID	Date of Birth	
John Doe	A005	17 NOV 1956	
Loaded MR Series	Date	Images	Comments
		0	
Loaded CT Series	Date	Images	Comments
		0	

Hardware Information			
	Model	Location	Date of Implant
 IPG 1	Vercise™		22 JUN 2012
 E1-E8	1x8	Right	22 JUN 2012

Close

The **Patient Information** window displays the following **Patient Information** and **Hardware Information**:

- Patient Name
- Patient ID
- Date of Birth
- Loaded MR Series (Date, Images, Comments)
- Loaded CT Series (Date, Images, Comments)
- IPG
 - Model
 - Date of Implant
- Lead 1
 - Model
 - Location
 - Date of Implant

- **Lead 2**
 - **Model**
 - **Location**
 - **Date of Implant**

3. Select  to exit the **Patient Information** window.

Searching for a Patient

You can search for a Patient by **Patient Name**, **Patient ID**, **Date of Birth**, **Lead 1 Implant Date**, **Lead 2 Implant Date**, or **Last Modified** date. When one or more characters is entered in the  **Search** text box, the Patient List will display only Patients whose data *begin* with the characters entered.

Note: *Search does not find characters within words.*

Clear the  **Search** text box to view all Patients in the Patient List.

Sorting the Patient List

To sort the Patient List by a particular column, select the column header. Select the header again to toggle between ascending and descending order:

-  = The files sort in ascending order.
-  = The files sort in descending order.

Exporting Patient Files

The **Export Patient Files** function allows you to select Patient files to export to the removable media (i.e., Flash Drive). The **Export Patient Files** screen provides the following:

- Capability to sort by column.
- Indication of whether or not the removable media is available to write to (see bottom right of screen).
- Indication of available storage room on the provided removable media.

Note: *The Optical Drive is Read Only.*

To export Patient data:

 **Warning:** *Exporting will delete any existing patients on the removable media.*

1. Insert the removable media (Flash Drive) to which you would like to export files.
2. From the Patient Browser, select  **Menu** > **Export Patient Files**. The **Export Patient Files** screen displays.

Export List Header

Export List

Menu

Patient Browser

New

Edit

Delete

Register

Program

Search

Export Patient Files

Select patients you would like to export to the removable media. *

☒	PATIENT NAME	PATIENT ID	DATE OF BIRTH	LEAD 1 IMPLANT DATE	LEAD 2 IMPLANT DATE	LAST MODIFIED
☒	Chung, Kristina H.	12345678	04 JAN 1957	09 APR 09	09 APR 09	12 MAR 09
☒	Kenton, Paige H.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	12 MAR 09
☒	Melton, Sherri E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	11 MAR 09
☒	Needren, Tom J.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	10 MAR 09
☒	Orion, Bill T.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09
☒	Plunket, Roger D.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09
☒	Shun, Phi	12345678	04 JAN 1957	12 MAR 09	12 MAR 09	12 MAR 09
☒	Sho, Pho	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09

* Warning: Exporting will delete any existing patients on the removable media.

OK

Cancel

Export Disk Space Indicator

Note: *All Patient files are selected by default on the **Export Patient Files** screen, and the files are initially sorted by **Patient Name**. To sort by a different column, select the column header.*

3. Select Patients you would like to export to the removable media.

Note: To de-select a Patient, uncheck the box next to the Patient's name. To de-select all Patients, uncheck the box next to **Patient Name**. To select all Patients, check the box next to **Patient Name**.

4. Select  to export the files.

Or

Select  to cancel the export function.

Note: If you first search in the **Patient Browser** for a subset of Patients, only the subset appears in the Export Patient Files List. Clear **Search** to export all Patients.

Importing Patient Files

The **Import Patient Files** function allows you to select Patient files to import from the removable media. The **Import Patient Files** screen provides the following:

- Capability to sort Patient files by column.
- Indication of whether or not the removable media is available to import from (see bottom right of screen).
- Indication of available storage room on the GUIDE™ DBS hard drive following importing of all selected Patients.

To import Patient Files:

1. Insert the removable media (Flash Drive) containing the Patient files.
2. Select  **Menu** > **Import Patient Files** from the drop-down list. The **Import Patient Files** screen displays.

Note: All Patient files are selected by default on the **Import Patient Files** screen, and the files are initially sorted by **Patient Name**. To sort by a different column, select the column header.

3. Select Patients you would like to import from the removable media.

Import List Header

Import List

Menu

New

Edit

Delete

Register

Program

Patient Browser

Search

Import Patient Files

Select patients you would like to import from the removable media.

	PATIENTNAME	PATIENT ID	DATE OF BIRTH	LEAD 1 IMPLANT DATE	LEAD 2 IMPLANT DATE	LAST MODIFIED
☒	Chung, Kristina H.	12345678	04 JAN 1957	09 APR 09	09 APR 09	12 MAR 09
☒	Kenton, Paige H.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	12 MAR 09
☒	Melton, Sherri E.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	11 MAR 09
☒	Needren, Tom J.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	10 MAR 09
☒	Orion, Bill T.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09
☒	Plunket, Roger D.	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09
☒	Shun, Phi	12345678	04 JAN 1957	12 MAR 09	12 MAR 09	12 MAR 09
☒	Sho, Pho	12345678	04 JAN 1957	09 MAR 09	09 MAR 09	09 MAR 09

Local Disk (C:)

OK

Cancel

Import Disk Space Indicator

Note: To de-select a Patient, uncheck the box next to the Patient's name. To de-select all Patients, uncheck the box next to **Patient Name**. To select all Patients, check the box next to **Patient Name**.

4. Select  to import the files.

Or

Select  to cancel the import function.

Image Registration

Optimizing Registration

GUIDE™ DBS Image Registration registers an atlas to the Patient's anatomy and locates the Lead relative to this.

The Image Registration:

- Fuses the Patient's MR and CT Volumes.
- Allows the user to define the Patient's AC, PC, and Mid-Sagittal plane.
- Transforms and fuses the atlas with the Patient's MR.
- Allows the user to define the Lead(s) within the Patient's post-op CT.
- Displays a preview of the Lead(s) in an atlas of the Patient's brain.

Minimum Scan Requirements

Minimum scan requirements are the least that are required to allow the GUIDE DBS software to be used.

- MR data sets must have at least 20 slices. They can be either T1 or T2 scans.
- CT data sets must have at least 20 slices. They must be taken at least two weeks and preferably at least four weeks following the last Lead implant date.

Recommended Scan Requirements

- MR data sets can be either T1 or T2 scans. They should be obtained at 1.5T with a field of view of at least 200 mm. The slices should be axial with a slice thickness of at most 2 mm. A submillimeter slice thickness is preferred.
- CT data sets should be obtained on a Helical scanner using a whole head protocol. The field of view should be at least 200 mm. The slices should be axial and the slice thickness should be at most 2 mm. A submillimeter slice thickness is preferred. The data set should be comprised of slices that include the Lead as well as slices that are distal to the Lead.
- Image scans should have consistent slice settings. When selecting a series, make certain the series has consistent slice settings.


Rescanning

Select **Rescan** any time a new removable media containing Patient imaging has been connected to the GUIDE DBS Station. Selecting the Rescan button will examine all removable media and generate a list of all imaging data found.

To rescan the removable media:

1. From the **Load MR Series** screen or the **Load Post-Op CT Series** screen, select **Rescan**.
2. The following message displays:

Are you sure you want to scan the removable media? Scanning may take a few minutes to complete.


3. Select **Yes** to scan the removable media.
Or
Select **No** to return to the **Load MR Series** or **Load Post-Op CT Series** screen without rescanning the removable media.
4. Verify that the MR and CT Patient names match the name on the Patient record.
5. If the names do not match, replace the removable media and select **Rescan**.

Cancel Rescanning

To cancel the rescan process:

1. Select .
2. The following message displays:

Are you sure you want to stop scanning the removable media? DICOM files listed may not be complete.

3. Select  to cancel the rescan process and return to the **Load MR Series** or **Load Post-Op CT Series** screen.

Or

Select  to continue the rescan process.


Image Registration Process Flow

The **Image Registration** process flows in the following order:

1. Load MR Series
2. Load Post-Op CT Series
3. Verify Fusion of MR and CT Scans
4. Select Mid-Sagittal Plane
5. Select Anterior Commissure (AC) and Posterior Commissure (PC)
6. Verify/Locate the Lead Tip
7. Verify/Locate the Lead Shaft
8. Left Lead Registration Preview (if applicable)
9. Right Lead Registration Preview (if applicable)

Nine separate Image Registration screens perform the various functions listed above. Use the  and  buttons to navigate between the **Image Registration** screens.

The  and  buttons are displayed only when going forward or backward is possible. The  button displays in a disabled state if all requirements for the current screen have not been met.

Note: You can end an Image Registration session from any screen in the Image Registration process by selecting  > **End Session**.

Note: You can view the Patient Information window from any screen by selecting  > **Patient Information**.

Loading MR Series


Image Data

- 1. From the **Enter Patient Information** screen, select the **Registration** button and select **Next** .


Or

From the **Patient Browser** screen, select a Patient and select **Register** . If the Patient record's **Registration** status is "**Load MR**", the **Load MR Series** screen will display. The components of the **Load MR Series** screen are described in the table below.

- 2. Select an item in the Load MR Series image list. The image displays in the Image Display Panel.

- 3. Select the Image Scrub Bar arrows to view the next and previous images in the series. Drag the Scrub Thumb Handle to view the image associated with the position of the Scrub Thumb Handle in the Image Scrub Area. Or, hold the Mouse cursor over the scrub bar to use the scroll wheel to scrub through the slices.

Thumb Handle

Note: For **Next** to be enabled, the series selected must contain at least 20 images. A series with fewer than 20 images will be grayed.

Load MR Series Screen		
	Component	Description
1		Allows you access to the following Menu options: ● Patient Information ● About ● End Session
2	Patient Name	Displayed at top of screen, Patient name as entered by the user.
3	 Registration Step Indicator	Indicates the current step in the registration process. Eight step icons display if one Lead is defined. Nine step icons display if two Leads are defined.
4	Load MR Series List Panel	Lists Patient name, date, images and comments.
5	Name	Patient name as it appears in the DICOM header.
6	Date	Date for the image series as it appears in the DICOM header.
7	Images	Number of images in the series.
8	Comments	Series descriptions as they appear in the DICOM header.
9		Allows you to Rescan the removable media.
10	Load MR Series Image Display Panel	2D display of the images in the selected series.
11		Allows you to choose from the following Alignment options: ● Flip Left and Right – flips the image across the sagittal plane. ● Flip Anterior and Posterior – flips the image across the coronal plane. ● Flip Superior and Inferior – flips the image across the axial plane.
12	 Image Scrub Bar	The Image Scrub Bar arrows allow you to view the next and previous images in the series. The Scrub Thumb Handle allows you to view the image associated with the position of the Scrub Thumb Handle in the Image Scrub Area.
13		Selects the highlighted image set and proceeds to the Next step in the Image Registration process. For  to be enabled, the series selected must contain at least 20 images.


en

4. If an image is not aligned correctly, select **Alignment** to open the **Alignment** menu.


⚠ Warning: Only use the Alignment button when an image is not oriented correctly.

5. Select from the following alignment options:
- **Flip Left and Right** – flips the image across the sagittal plane.
 - **Flip Anterior and Posterior** – flips the image across the coronal plane.
 - **Flip Superior and Inferior** – flips the image across the axial plane.

Before


After

Loading Post-Op CT Series


1. From the **Load MR Series** screen, select **Next** to select the highlighted MR image series and move to the **Load Post-Op CT Series** screen. The components of the **Load Post-Op CT Series** screen are described in the table below.

Note: For **Next** to be enabled, the series selected must contain at least 20 images. A series with fewer than 20 images will be grayed.

2. Select an item in the Load Post-Op CT Series image list to display it in the Image Display Panel.
3. Select the Image Scrub Bar arrows to view the next and previous images in the series. Drag the Scrub Thumb Handle to view the image associated with the position of the Scrub Thumb Handle in the Image Scrub Area. Or, hold the Mouse cursor over the scrub bar to use the scroll wheel to scrub through the slices.

Load Post-Op CT Series Screen		
	Component	Description
1		Allows you to select from the following Menu options: ● Patient Information ● About ● End Session
2	Patient Name	Displayed at top of screen, Patient name as entered by the user.
3	 Registration Step Indicator	Indicates the current step in the registration process. Eight step icons display if one Lead is defined. Nine step icons display if two Leads are defined.
4	Load Post-Op CT Series List Panel	Lists Patient name, date, images, and comments.
5	Name	Patient name as it appears in the DICOM header.
6	Date	Date for the image series as it appears in the DICOM header.
7	Images	Number of images in the series.
8	Comments	Series descriptions as they appear in the DICOM header.
9		Allows you to Rescan the removable media.
10	Load Post-Op CT Series Image Display Panel	2D display of the images in the selected series.
11		Allows you to choose from the following Alignment options: ● Flip Left and Right – flips the image across the sagittal plane. ● Flip Anterior and Posterior – flips the image across the coronal plane. ● Flip Superior and Inferior – flips the image across the axial plane.
12	 Image Scrub Bar	The Image Scrub Bar arrows allow you to view the next and previous images in the series. The Scrub Thumb Handle allows you to view the image associated with the position of the Scrub Thumb Handle in the Image Scrub Area.
13		The Back button returns you to the previous step in the Image Registration process.
14		The Next button advances you to the next step in the Image Registration process.

4. If an image is not aligned correctly, select **Alignment** to open the **Alignment** menu.

CAUTION: Only use the Alignment button when an image is not oriented correctly.

5. Select from the following alignment options:
- **Flip Left and Right** – flips the image across the sagittal plane.
 - **Flip Anterior and Posterior** – flips the image across the coronal plane.
 - **Flip Superior and Inferior** – flips the image across the axial plane.
6. Select **Next** to proceed to the **Verify Fusion of MR and CT Scans** screen.
7. The **Verify Patient Information** Dialog will display.

The dialog box is titled "Verify Patient Information" and contains two main sections: "Patient Information" and "Hardware Information".

Patient Information

Patient Name	Patient ID	Date of Birth
Sullivan Weissenburgh	201	05 FEB 1975

Loaded MR Series	Date	Images	Comments
Sullivan Weissenburgh	03 SEP 2009	166	AX T2 GRE

Loaded CT Series	Date	Images	Comments
Sullivan Weissenburgh	05 OCT 2009	98	1 mm Axial for DBS

Hardware Information

Model	Location	Date of Implant
IPG Vercise™		21 DEC 2011
E1-8 In-Line 1x8	Left	21 DEC 2011
E9-E16 In-Line 1x8	Left	21 DEC 2011

At the bottom are "Continue" and "Cancel" buttons.

8. Verify that the MR and CT Patient names match the name on the Patient record.
9. Verify that the correct number of Leads is listed under “**Hardware Information**”.

Note: If only one Lead is selected while two Leads are implanted, GUIDE™ DBS will not be able to automatically identify the Lead Tips. You may have to start over again from this step.

10. Select **Cancel** to cancel automated fusion and return to the **Load Post-Op CT Series** screen.
- Or
- Select **Continue** to proceed with automated fusion.


11. The **Verify Fusion of MR and CT Scans** screen will display.

en

Verifying Fusion of MR and CT Scans

CAUTION: Verifying the fusion of MR and CT scans requires MR and CT image interpretation skills and anatomical knowledge of brain structures.

1. From the **Load Post-Op CT Series** screen, select **Next** to begin automated fusion. The process will load the selected MR and CT images from the external drive.
2. The **Verify Fusion of MR and CT Scans** screen will display. The components of the **Verify Fusion of MR and CT Scans** screen are described in the table below.
3. Use either the Merged or Overlay Modes (described in the following pages) to verify the alignment of the MR and CT images. You can make any necessary adjustments to the alignment using the XForm tool.

Verify Fusion of MR and CT Series Screen		
	Component	Description
1		Allows you to access the following Menu options: ● Patient Information ● About ● End Session
2	 Registration Step Indicator	Indicates the current step in the registration process. Eight step icons display if one Lead is defined. Nine step icons display if two Leads are defined.
3	Toolbar	GUIDE™ DBS Image Registration tools.
4	 DRAG	Select to Drag the MR/CT images in the Viewports.
5	 ZOOM	Select to Zoom in/out on the MR/CT images in the Viewports.
6	 LEVELS	Select to change the window/level of the MR/CT images.
7	 XFORM (Transform)	Select to translate and rotate the MR relative to the CT.
8	 AUTO FUSE	Select to re-perform automatic fusion.
9	 RESET	Select to return the Markers to the Marker Drawer.
10	 UNDO	Select to Undo the last action performed that has not been undone.
11	 REDO	Select to Redo an action. Only actions that have been undone with the Undo button can be redone.
12	 View Mode	When the MR Only button is selected, the Axial, Coronal, and Sagittal Viewports display only the respective MR Viewport Images.
13	 View Mode	When the CT Only button is selected, the Axial, Coronal, and Sagittal Viewports display only the respective CT Viewport Images.
14	 View Mode	When the Merged button is selected, the Axial, Coronal, and Sagittal Viewports display both the respective MR Viewport Images and the respective CT Viewport Images. (The images are displayed overlapping.)

Verify Fusion of MR and CT Series Screen		
15	 View Mode	When the Overlay button is selected, the Axial, Coronal, and Sagittal Viewports display both the MR Viewport Images and CT Viewport Images. (The CT images are displayed in a vertical strip that is overlaid on the MR images.)
16	Axial Viewport	Displays a 2D image representing a plane from the Axial anatomical perspective. Each Viewport contains: • Orientation labels. • Volume Location Indicator (red line) when one of the other Viewports is being scrubbed. • Volume Scrub Bar Control.
17	Coronal Viewport	Displays a 2D image representing a plane from the Coronal anatomical perspective.
18	Sagittal Viewport	Displays a 2D image representing a plane from the Sagittal anatomical perspective.
19	 Viewport Control	Changes the viewable size of the Viewport Images.
20	 Blend Control	Allows you to change the transparency of the MR or CT in the fused image. (Only visible in Merged View.)
21	 Volume Scrub Bar	Allows you to modify which 2D image of the respective Volume is displayed.
22		Select the Back button to return to the previous Registration step.
23		Select to proceed to the Next step in the Image Registration process.

Using the Image Registration Tools

Drag Tool

Select  **DRAG** to move an image. Place the cursor on the image. Left click the Mouse and drag the cursor to move the image. Alternatively, use the Mouse scroll wheel to move the image up/down, and with Shift + scroll wheel, move the image right/left.


Zoom Tool

Select  **ZOOM** to zoom in and out on an image. Place the cursor on the image. Left click the Mouse and drag the cursor downward to zoom out. Drag the cursor upward to zoom in. Alternatively, use the Mouse scroll wheel to zoom in/out.

Levels Tool

The terms “**window**” and “**level**” are related to contrast and brightness adjustment. Making the **window** wider or narrower decreases or increases the display contrast. Moving the **level** left or right changes the display brightness.

Select  **LEVELS** to open the CT and MR Window Level Controls window.


Note: If the MR Viewport Images are not displayed, the MR Window Level Controls will not display. If the CT Viewport Images are not displayed, the CT Window Level Controls will not display.

CT Level Elevator

Drag the CT Level Elevator to the left to decrease the level of the CT Viewport Image, and to the right to increase the level.

CT Window Elevator

Drag the CT Window Elevator to the left to decrease the window of the CT Viewport Image, and to the right to increase the window.

CT Leveling Technique

The dynamic range of some CT images when imported makes leveling more difficult. To optimize leveling of these CT images:


1. From the **CT ONLY** screen, select **LEVELS** to open the CT Levels Control.
2. Slide the CT Level Control and the CT Window Control to the far left, so that the CT image is fully white.
3. Slowly slide the CT Level Control to the right until the image inside the skull is approximately 50% black and 50% white.
4. Slowly slide the CT Window Control to the right until the brain structures (ventricles, sulci, gyri, cerebellum, brainstem, etc.) are visible.
5. To fine-tune leveling, click in the black bars, to the left and right of the slider, or place the Mouse above one of the sliders and use the scroll wheel to finely move the bar.

MR Level Elevator

Drag the MR Level Elevator to the left to decrease the level of the MR Viewport Image, and to the right to increase the level.

MR Window Elevator


Drag the MR Window Elevator to the left to decrease the window for the MR Viewport Image, and to the right to increase the window.

en

Transform Tool

Select  **XFORM** to open the **Transform** tool and rotate the MR image relative to the CT image. If the view mode is **MR Only** or **CT Only**, the mode will change to **Merged**.

- Select and drag the Rotation Handle  counterclockwise to rotate the MR Volume counterclockwise on an axis perpendicular to the plane of the Viewport and through the Rotation Anchor.
- Select and drag the Rotation Handle  clockwise to rotate the MR Volume clockwise on an axis perpendicular to the plane of the Viewport and through the Rotation Anchor.
- Select and drag the Translation Handle  to translate the Volume in the plane of the Viewport. The Rotation Handle and Rotation Anchor will track the Translation Handle.

Note: The Rotation Anchor is used as the center of rotation and translation. All Viewports will be updated during the rotation and translation to reflect the changes in the Volume.

- Select  **UNDO** to reverse a rotation or translation.

Auto-Fuse Tool

Select  **FUSE** to automatically fuse the MR and CT scans.

Undo Tool

Select  **UNDO** to reverse the last action performed.

Redo Tool

Select  **REDO** to redo the last action undone.

Selecting a Registration View Mode


GUIDE™ DBS has four Registration View Modes:

- **MR ONLY** MR Only
- **CT ONLY** CT Only
- **MERGED** Merged
- **OVERLAY** Overlay

The Registration View Mode buttons are located in the top right of the **Verify Fusion of MR and CT Scans** screen. Select a button to change the View Mode.


MR Only

When the **MR ONLY** View Mode is selected, the Axial, Coronal, and Sagittal Viewports display only the respective MR Viewport Images.

CT Only

When the **CT ONLY** View Mode is selected, the Axial, Coronal, and Sagittal Viewports display only the respective CT Viewport Images.


Merged and Overlay

The **Verify Fusion of MR and CT Scans** screen provides two methods of comparing the MR to CT fusion: **Merged** and **Overlay**.

Note: If the quality of the images is sufficient but the fusion is not acceptable, you can use the tools on the screen to modify the alignment of the MR and CT images. If the fusion is not successful because the MR or CT image is not of sufficient quality, press **Back** to return to the **Load MR Series** and **Load Post-Op CT Series** screens to select different MR or CT images. Then repeat the **Load MR Series** and **Load Post-Op CT Series** registration steps.

Merged

When **MERGED** is selected, the Axial, Coronal, and Sagittal Viewports display both the respective MR Viewport Images and the respective CT Viewport Images. Merged View shows transparent images with the CT on top of the MR. The transparency of the CT and MR can be controlled to compare the two modalities.

Blend Control


The Blend Control allows you to change the transparency of the MR or CT in the fused image. Sliding the Blend Control to the left makes the CT more transparent, showing more of the MR. Sliding the control to the right makes the MR more transparent, showing more of the CT. When the Blend Control is at the half-way point, the CT image and MR image are both at 50% transparency.

Note: The Blend Control is only visible in Merged View.

Overlay


Overlay View, seen below, shows a slice of the CT displayed fully opaque over the MR. The width of the slice and its location can be varied. Setting the Level and Window for the MR and CT images allows you to compare the images at the brightness and contrast stretch level that is most effective.

Overlay Handle


The viewable area of the Overlaying Image is constrained by the size of the Overlay View Area. The viewable area can be changed using the Overlay Handle Control.

- Drag the Overlay Left Handle left/right to view more/less of the left side of the CT Image.
- Drag the Overlay Right Handle right/left to view more/less of the right side of the CT Image.
- Drag the Overlay Thumb Handle to move the Overlay View Area.

Using the Viewports

The Viewport Image is a 2D image representing a plane from the Volume. The Viewport Image can be from an Axial, Coronal, or Sagittal perspective. In the following Viewport Images, the Volume Location Indicator (red line) shows the relative plane in the two orthogonal views.


The components of the Viewports are shown in the table below.

Viewports		
	Component	Description
1	Viewport Image	2D image representing a plane from the Volume.
2	Viewport Handle Control	Allows you to change the viewable size of the Viewports.
3	Volume Location Indicator (red line)	Relative plane in the two orthogonal views.
4	Orientation Label	Left. Right. Anterior. Posterior. Superior.
5	 Volume Scrub Bar Control	Allows you to modify which 2D image of the respective Volume is displayed. The Volume Location is displayed above the handle.
6	View Icon	Axial. Coronal. Sagittal.

Viewport Handle Control

The Viewports are separated by a Viewport Handle Control. The Viewport Handle is used to change the viewable size of the Viewport Images. When one of the Viewport Images is increased in size (zooming in), the adjacent Viewport Image(s) will decrease in size (zoom out) to maintain the Image within the Viewport. Drag the Viewport Handle to the left or right to adjust the Viewport Image size. Double select any Viewport Handle to return all of the Viewports to their original size.


Volume Scrub Bar Control

The Volume Scrub Bar Control Right and Left Arrows are used to navigate through the Viewport Images. Select the Right Arrow to display the next image in the series. Select the Left Arrow to display the previous image in the series.

The Volume Location represents the depth within the Volume of the Viewport Image displayed.

The Volume Scrub Bar Thumb Handle is used to display the plane of the Volume that is parallel to the Viewport and located at the depth displayed by the Volume Location Indicator. Drag the Thumb Handle left or right to view adjacent slices.

When the Volume Scrub Bar Thumb Handle is selected, the Volume Location Indicator will display as a red line in the other Viewports. The red line represents the Volume Location from the Volume Scrub Thumb Handle that is selected.

Selecting the Mid-Sagittal Plane


CAUTION: *Selecting the Mid-Sagittal Plane of the Patient requires anatomical knowledge of brain structures. The instructions in this manual describe how to use GUIDE™ DBS to specify the location of the Mid-Sagittal Plane in the Patient’s MR. They are not meant to describe how to locate the Mid-Sagittal Plane in a Patient’s MR.*

From the **Verify Fusion of MR and CT Scans** screen, select **Next** to open the **Select Mid-Sagittal Plane** screen.

The **Select Mid-Sagittal Plane** screen allows you to identify and specify the Mid-Sagittal Plane (MSP) in the Patient’s MR. The MSP is the imaginary plane that separates the left and right hemispheres. In the Axial and Coronal Viewports, the brain is approximately mirrored around it.

At any time, the Sagittal Viewport displays the MR slice through the selected MSP.


Using the tools on the screen you can align the MSP with the Mid-Sagittal Indicator (dashed orange lines running vertically in the Axial and Coronal Viewports).

The tools available for moving and rotating the Viewport Images are the same as those in the **Verify Fusion of MR and CT Scans** screen, with some modifications:

- The  **FUSE** button is disabled.
- The  **DRAG** button in the Axial and Coronal Viewports will not move the images to the left or right, only up or down. Use the center point on the  **XFORM** (Transform) control to move images left and right.
- In the Sagittal Viewport, the Volume Scrub Bar is disabled.

All actions that change the selection of the MSP must be done using the  **XFORM** (Transform) control.

To identify and align the MSP:

1. Begin in either the Axial or Coronal Viewport.
2. Zoom in so that only the deep structures are visible within the Viewport.
3. Select  **XFORM** (Transform).
4. Use the Rotation Arm  to rotate the MR scan so that the MSP lines up parallel to the dotted line.
5. Use the center point  to translate the MR left or right to align the dotted line with the MSP.

6. Repeat steps 4 and 5 in the other Viewport.
7. In both Viewports, scrub through the slices to verify that the MSP lines up with the dotted line within the deep structures of the brain.

This alignment is an iterative process, where you move and rotate the MR to bring it closer in line with the Mid-Sagittal indicator, until you are satisfied with the alignment.

Selecting the Anterior Commissure (AC) and Posterior Commissure (PC)

CAUTION: Selecting the Anterior and Posterior Commissures requires anatomical knowledge of brain structures on the part of the user. The instructions in this manual describe how to use the GUIDE™ DBS tools to specify where the AC and PC are located in a Patient's MR and are not meant to describe how to identify them.


From the **Select Mid-Sagittal Plane** screen, select **Next** to open the **Select Anterior Commissure (AC) and Posterior Commissure (PC)** screen.

The **Select Anterior Commissure (AC) and Posterior Commissure (PC)** screen allows you to specify the location of the Anterior and Posterior Commissures in the Patient's MR.

The tools available on the **Select Anterior Commissure (AC) and Posterior Commissure (PC)** screen are the same as those in the **Verify Fusion of MR and CT Scans** screen, with the following modifications:

- The AC/PC Markers and Marker Drawer are added.
- The  **RESET** button is enabled. When selected, the button returns the Markers to the Marker Drawer.

When the **Select Anterior Commissure (AC) and Posterior Commissure (PC)** screen opens, the AC/PC Markers are in the Marker Drawer.

Using the AC/PC Markers


The Markers in the three Viewports represent the same point in space from different perspectives. When a Marker has been placed in one Viewport, the Marker will automatically be placed in the other Viewports. When a Marker is moved in one Viewport, the Markers in the other Viewports move as well.

When a Marker is placed and the Viewport Image is scrubbed using the Volume Scrub Bar, the Marker will become more transparent as the Volume Location Indicator changes. This represents the Marker moving through the slices. When the Volume Location Indicator is 5 mm or more from the value that the Marker was dropped on, the Marker will become hidden. Select the Marker to snap the view in all three Viewports to the slice containing the Marker.

To place the AC/PC Markers:

1. From the AC/PC Marker Drawer, select and drag the  AC Marker or  PC Marker and place it on the appropriate location of any Viewport Image.
2. Select  **RESET** to return the Markers to the Marker Drawer.
3. When satisfied with the location of the AC and PC Markers, select  to complete this registration step and move on to the next step: Verify/Locate the Lead Tip.

Note: *The AC Marker should not be placed more posterior than the PC Marker. If it is,  is disabled and the AC-PC distance indicator will be hidden.*

Note: AC/PC Distance displays the straight line distance between the placed AC and PC Markers.

Verifying/Locating the Lead Tips

CAUTION: *Identifying and verifying the locations of the Lead Tips requires the capability to interpret CT imagery. The instructions in this manual describe how to use the GUIDE™ DBS tools to verify and adjust the location of the Lead Tips and are not meant to describe how to identify them.*

The tools available for moving and rotating the Viewport Images are the same as those in the **Verify Fusion of MR and CT Scans** screen. In addition:

- The  **LOCATE** button is enabled. Select  **LOCATE** to repeat the auto locate process.

Verifying the Lead Tips

When the **Locate the Lead Tip** screen is displayed and the Lead Tip Markers have not been placed, GUIDE DBS will attempt to automatically locate the Lead Tips and Shafts. Leads appear as bright spots in CT images.

Note: *The number of Leads defined on the **Patient Information** screen and the Lead location must match the number and location of the Leads actually in the CT image or the auto locate process will fail.*


1. The **Locating Lead Tips and Shafts** message will display automatically.
2. GUIDE DBS will locate the Lead Tips and Shafts within one minute.
3. If GUIDE DBS locates the Lead Tips and Shafts:
 - The **Verify the Lead Tip** screen Viewport Images will be updated with the Lead Tip 1 Marker Control placed.
 - If two Leads are defined, the Lead Tip 2 Marker Control will also be placed in the Viewport Images.
4. Select each Marker individually to verify placement. Make any necessary adjustments to ensure the proper positioning using the manual method detailed on page 64.
 - To move a Marker, select and drag a Marker in one of the Viewports. The Markers in the other two Viewports will update automatically.
5. Select  to accept the location of the Lead Tip Marker(s) and complete this registration step.


Failure: Number of Leads not Found

If GUIDE™ DBS fails to locate the Lead(s) the following message displays:

Failed to locate lead. The defined number of leads was not found.

This message may display because the number of Leads defined on the Patient Information screen does not match the number of Leads actually in the CT image.

1. Select  .
2. Select  > **End Session** to end the Registration Session. See “Ending a Registration Session” on page 72.
3. Edit the Patient Information. See “Editing a Patient” on page 22.
4. Repeat the Image Registration Process.

This message may also display if the Mid-Sagittal Plane is poorly aligned.

1. Select  .
2. Select  to return to the **Select Mid-Sagittal Plane** screen and re-align the MSP.


This message may also display if you selected a pre-op CT that does not have a Lead in it.

1. Select  .
2. Select  to return to the **Load Post-Op CT Series** screen and select a CT with a Lead in it.

Failure: No Leads Found

If the location of the Lead defined on the Patient Information screen does not match the location of the Lead actually in the CT image (for instance, if the Lead location is set to “Left” but the Lead is actually in the right hemisphere), the following message displays:

Failed to locate lead. No leads were found.

1. Select  .
2. Select  > **End Session** to end the Registration Session. See “Ending a Registration Session” on page 72.
3. Edit the Patient Information so that the Lead’s specified location matches the CT. See “Editing a Patient” on page 22.

If the location of the Lead on the Patient Information screen and the CT image matches:

1. Select  until the **Load Pre-op CT Series** screen displays.
2. Adjust image alignment, if needed.
3. Repeat the Image Registration Process.


Locating the Lead Tips Manually

If GUIDE™ DBS fails to locate the Lead Tips, and the number and location of Leads on the Patient Information screen match the number and location of Leads actually in the CT Image, use the Lead Tip Markers to manually locate the Lead Tips. Leads appear as bright spots in CT images.

Note: The number of Lead Tip Markers in the Marker Drawer will match the number of Leads defined.

To place the Lead Tip Marker for a given Lead:

1. Determine the distal-most slice in which the Lead is visible (in the Sagittal View).
2. Zoom in on the bright spot corresponding to the Lead in the slice.
3. Drag the Lead Tip Marker from the Lead Tip Marker Drawer into place on the appropriate Viewport Image.
4. Place the Marker at the center of the Lead Tip. Adjust as needed in the other two views to center the Marker on the center of the Lead Tip in all perspectives.

5. Select **Next** to accept the location of the Lead Tip Marker(s) and complete this registration step.

Note: Select a Marker to move the scrub level of all Viewports to align to the location of that Marker.

Verifying/Locating the Lead Shafts

CAUTION: Identifying and verifying the locations of the Lead Shafts requires the capability to interpret CT imagery. The instructions in this manual describe how to use the GUIDE™ DBS tools to verify and adjust the location of the Lead Shafts. They are not meant to describe how to identify the Lead Shafts.


The tools available for moving and rotating the Viewport Images are the same as those in the **Verify Fusion of MR and CT Scans** screen.

Verifying the Lead Shafts

If GUIDE DBS automatically locates the Lead Tips and Lead Shafts:

1. From the **Verify the Lead Tip** screen select  to display the **Verify the Lead Shaft** screen.
 - The **Verify the Lead Shaft** screen Viewport Images will be updated with the Lead Shaft 1 Marker Control placed.
 - If two Leads are defined, the Lead Shaft 2 Marker Control will also be placed in the Viewport Images.
2. Select each Marker individually to verify placement. Make any necessary adjustments to ensure the proper positioning using the manual method detailed on page 67.
 - To move a Marker, select and drag a Marker in one of the Viewports. The Markers in the other two Viewports will update automatically.
3. Select  to accept the location of the Lead Shaft Marker(s) and complete this registration step.

en


Locating the Lead Shafts Manually

If GUIDE™ DBS fails to locate the Leads, drag the Lead Shaft Markers from the Lead Shaft Marker Drawer into place on the appropriate Viewport Image. Lead Shaft Markers should be placed at least 10 mm superior to the Lead Tip Markers.

Note: If the Lead Shaft Markers are not placed superior to the Lead Tip Markers, the **Next** button will be disabled.

To place the Lead Shaft Markers:

1. Identify the correct slice (as with the Lead Tip Marker placement).
2. Identify the bright spot corresponding to the Lead within this slice.
3. Zoom in on this bright spot and place the Marker at its center.
4. Repeat for second Lead Shaft, if appropriate.
5. Select **Next** to accept the location of the Lead Shaft Marker(s) and complete this registration step.

Note: Select a Marker to move the scrub level of all Viewports to align to the location of that Marker.


Previewing Left Lead Registration

CAUTION: *Verifying the overall Lead to atlas registration requires anatomical knowledge of brain structures and subject matter expertise on DBS Leads and their placement. The instructions in this manual describe how to use the GUIDE™ DBS tools to verify the registration and are not meant to describe a valid registration.*

If there is a Left Lead, registration preview will start in the **Left Lead Registration Preview** screen. Otherwise, it will start in the **Right Lead Registration Preview** screen.

The **Left Lead Registration Preview** screen allows you to see where the Left DBS Lead is located with respect to the various brain structures and reference volumes in the GUIDE DBS atlas. From the **Left Lead Registration Preview** screen you can:

- Change the view using the  Vertical and  Horizontal Viewport Slider controls, or left select the image and drag, to view from any angle looking in towards the center of the contact areas on the Lead.
- Zoom in using the  button or by hovering the Mouse on the image and scrolling up.
- Zoom out using the  button or by hovering the Mouse on the image and scrolling down.
- Open the **View**  menu to enable or disable the display of objects in the atlas:
 - **Anterior Commissure (AC)**
 - **Encapsulation Layer**
 - **Posterior Commissure (PC)**
 - **Red Nucleus (RN)**
 - **Reference Volume**
 - **Stimulation Field Model (SFM)**
 - **Subthalamic Nucleus (STN)**
 - **Thalamus**

To assist with orientation, a reference head displays in the lower right corner. The head rotates with the displayed perspective.

Select **Next** to display the **Right Lead Registration Preview** screen, or select **Finish** if there is no Right Lead. Pressing **Finish** will complete the registration process and display the Patient List.


Previewing Right Lead Registration

CAUTION: *Verifying the overall Lead to atlas registration requires anatomical knowledge of brain structures on the part of the user and subject matter expertise on DBS Leads and their placement. The instructions in this manual describe how to use the GUIDE™ DBS tools to verify the registration and are not meant to describe a valid registration.*

The **Right Lead Registration Preview** screen opens only if a Right Lead is defined for the Patient. If both a Left and a Right Lead are defined for the Patient, the registration preview will start in the **Left Lead Registration Preview** screen. If only a Right Lead is defined for the Patient, registration preview will start in the **Right Lead Registration Preview** screen.

The **Right Lead Registration Preview** screen allows you to see where the Right DBS Lead is located with respect to the various brain structures and reference volumes in the GUIDE DBS atlas. From the **Right Lead Registration Preview** screen you can:

- Change the view using the  vertical and  horizontal Viewport Slider Controls, or left select the image and drag, to view from any angle looking in towards the center of the contact areas on the Lead.
- Zoom in using the  button or by hovering the Mouse on the image and scrolling up.
- Zoom out using the  button or by hovering the Mouse on the image and scrolling down.
- Open the  menu to enable or disable the display of objects in the atlas:
 - **Anterior Commissure (AC)**
 - **Encapsulation Layer**
 - **Posterior Commissure (PC)**
 - **Red Nucleus (RN)**
 - **Reference Volume**
 - **Stimulation Field Model (SFM)**
 - **Subthalamic Nucleus (STN)**
 - **Thalamus**

To assist with orientation, a reference head displays in the lower right corner. The head rotates with the displayed perspective.


When you are done previewing the Right Lead, select **Finish** to complete the registration process and display the Patient List.

Ending a Registration Session

1. Select  **Menu** > **End Session**. A dialog box displays the following question:

Would you like to save changes to the patient record before ending the session?

2. Select one of the following buttons:
-  to end the session and save changes to the Patient record.
 -  to end the session without saving changes to the Patient record.
 -  to cancel ending the session and return to the registration process.


Programming Simulation



On the **Simulation** screen you can create and save four simulations for each Lead implanted. For each simulation, the Lead and Stimulation Field Model (SFM) created are displayed relative to the atlas that has been registered to the Patient's anatomy.


CAUTION: SFMs generated without Image Registration (i.e., when “Skip Registration” is selected) are not suitable for use in a Patient.



Creating Programming Simulations



1. On the **Patient Browser** screen, select a Patient and select **Program** . (Or, if the registration step is listed as “Complete”, double select the Patient name.) The **Simulation** screen displays. The components on the **Simulation** screen are described in the table below.

Programming Simulation Screen		
	Component	Function
1		Allows you to select from the following Menu options: ● Patient Information ● About ● End Session
2		Displays a simulation of the Left Lead Program in the atlas; may be unavailable if there is no Left Lead.
3		Displays a simulation of the Right Lead Program in the atlas; may be unavailable if there is no Right Lead.
4		Displays a Simulation Report . See “Creating a Simulation Report” on page 92 for details. May be unavailable if “Skip Registration” was selected.
5	Simulation Panel	GUIDE™ DBS allows you to access, add, and save up to four simulated programs for each Lead. The simulated programs are displayed in the Simulation Panel . See “Using the Simulation Panel” on page 78.
6	2D Parameters Panel	Defines the requirements for the programming simulation. Contains the settings for the current programming simulation. See “Creating Programming Simulations” on page 73.
7	Lead Contact Amplitude  % Mode  mA Mode	An indication of the current amplitude being simulated to a given contact. Use the arrows to increase/decrease the setting. For % Mode When the system is in % Mode, the overall amount of current being applied in mA to all active contacts is displayed in the Amplitude box. The fraction of this overall current which is shuttled to individual contacts is displayed next to each contact as a percentage of 100%. In a bipolar configuration, the anodic contacts and the cathodic contacts each add up to 100%. For mA Mode When the system is in mA Mode, the absolute amount of current shuttled to each contact is displayed next to that contact. The sum of the current being sent to all cathodic contacts is displayed, in a lighter gray, in the Amplitude box.

Programming Simulation Screen																																																																																																																																				
8	Lead Contacts	Each Lead Contact displays a number corresponding to the physical Lead contact it represents and an image representing the Lead Contact State: Selecting the Lead Contact toggles the contact between Open State, Cathode State, and Anode State, in that order. 4 = Open State (the initial value) 4 — = Cathode State 4 + = Anode State																																																																																																																																		
9		The Amplitude button displays the Current in milliamps when the system is in % Mode. When the Current Adjusters are selected, the Current changes to the next/previous valid current value. Alternatively, select the Amplitude button to display an array of values to choose from. 0.6 mA <table><tr><td>0.0 mA</td><td>1.0 mA</td><td>2.0 mA</td><td>3.0 mA</td><td>4.0 mA</td><td>5.0 mA</td><td>6.0 mA</td><td>7.0 mA</td><td>8.0 mA</td><td>9.0 mA</td><td>10.0 mA</td><td>11.0 mA</td><td>12.0 mA</td></tr><tr><td>0.1 mA</td><td>1.1 mA</td><td>2.1 mA</td><td>3.1 mA</td><td>4.1 mA</td><td>5.1 mA</td><td>6.1 mA</td><td>7.1 mA</td><td>8.1 mA</td><td>9.1 mA</td><td>10.1 mA</td><td>11.1 mA</td><td>12.1 mA</td></tr><tr><td>0.2 mA</td><td>1.2 mA</td><td>2.2 mA</td><td>3.2 mA</td><td>4.2 mA</td><td>5.2 mA</td><td>6.2 mA</td><td>7.2 mA</td><td>8.2 mA</td><td>9.2 mA</td><td>10.2 mA</td><td>11.2 mA</td><td>12.2 mA</td></tr><tr><td>0.3 mA</td><td>1.3 mA</td><td>2.3 mA</td><td>3.3 mA</td><td>4.3 mA</td><td>5.3 mA</td><td>6.3 mA</td><td>7.3 mA</td><td>8.3 mA</td><td>9.3 mA</td><td>10.3 mA</td><td>11.3 mA</td><td>12.3 mA</td></tr><tr><td>0.4 mA</td><td>1.4 mA</td><td>2.4 mA</td><td>3.4 mA</td><td>4.4 mA</td><td>5.4 mA</td><td>6.4 mA</td><td>7.4 mA</td><td>8.4 mA</td><td>9.4 mA</td><td>10.4 mA</td><td>11.4 mA</td><td>12.4 mA</td></tr><tr><td>0.5 mA</td><td>1.5 mA</td><td>2.5 mA</td><td>3.5 mA</td><td>4.5 mA</td><td>5.5 mA</td><td>6.5 mA</td><td>7.5 mA</td><td>8.5 mA</td><td>9.5 mA</td><td>10.5 mA</td><td>11.5 mA</td><td>12.5 mA</td></tr><tr><td>0.6 mA</td><td>1.6 mA</td><td>2.6 mA</td><td>3.6 mA</td><td>4.6 mA</td><td>5.6 mA</td><td>6.6 mA</td><td>7.6 mA</td><td>8.6 mA</td><td>9.6 mA</td><td>10.6 mA</td><td>11.6 mA</td><td>12.6 mA</td></tr><tr><td>0.7 mA</td><td>1.7 mA</td><td>2.7 mA</td><td>3.7 mA</td><td>4.7 mA</td><td>5.7 mA</td><td>6.7 mA</td><td>7.7 mA</td><td>8.7 mA</td><td>9.7 mA</td><td>10.7 mA</td><td>11.7 mA</td><td>12.7 mA</td></tr><tr><td>0.8 mA</td><td>1.8 mA</td><td>2.8 mA</td><td>3.8 mA</td><td>4.8 mA</td><td>5.8 mA</td><td>6.8 mA</td><td>7.8 mA</td><td>8.8 mA</td><td>9.8 mA</td><td>10.8 mA</td><td>11.8 mA</td><td></td></tr><tr><td>0.9 mA</td><td>1.9 mA</td><td>2.9 mA</td><td>3.9 mA</td><td>4.9 mA</td><td>5.9 mA</td><td>6.9 mA</td><td>7.9 mA</td><td>8.9 mA</td><td>9.9 mA</td><td>10.9 mA</td><td>11.9 mA</td><td></td></tr></table> Select  or  to navigate through the array of values. Select the desired value.	0.0 mA	1.0 mA	2.0 mA	3.0 mA	4.0 mA	5.0 mA	6.0 mA	7.0 mA	8.0 mA	9.0 mA	10.0 mA	11.0 mA	12.0 mA	0.1 mA	1.1 mA	2.1 mA	3.1 mA	4.1 mA	5.1 mA	6.1 mA	7.1 mA	8.1 mA	9.1 mA	10.1 mA	11.1 mA	12.1 mA	0.2 mA	1.2 mA	2.2 mA	3.2 mA	4.2 mA	5.2 mA	6.2 mA	7.2 mA	8.2 mA	9.2 mA	10.2 mA	11.2 mA	12.2 mA	0.3 mA	1.3 mA	2.3 mA	3.3 mA	4.3 mA	5.3 mA	6.3 mA	7.3 mA	8.3 mA	9.3 mA	10.3 mA	11.3 mA	12.3 mA	0.4 mA	1.4 mA	2.4 mA	3.4 mA	4.4 mA	5.4 mA	6.4 mA	7.4 mA	8.4 mA	9.4 mA	10.4 mA	11.4 mA	12.4 mA	0.5 mA	1.5 mA	2.5 mA	3.5 mA	4.5 mA	5.5 mA	6.5 mA	7.5 mA	8.5 mA	9.5 mA	10.5 mA	11.5 mA	12.5 mA	0.6 mA	1.6 mA	2.6 mA	3.6 mA	4.6 mA	5.6 mA	6.6 mA	7.6 mA	8.6 mA	9.6 mA	10.6 mA	11.6 mA	12.6 mA	0.7 mA	1.7 mA	2.7 mA	3.7 mA	4.7 mA	5.7 mA	6.7 mA	7.7 mA	8.7 mA	9.7 mA	10.7 mA	11.7 mA	12.7 mA	0.8 mA	1.8 mA	2.8 mA	3.8 mA	4.8 mA	5.8 mA	6.8 mA	7.8 mA	8.8 mA	9.8 mA	10.8 mA	11.8 mA		0.9 mA	1.9 mA	2.9 mA	3.9 mA	4.9 mA	5.9 mA	6.9 mA	7.9 mA	8.9 mA	9.9 mA	10.9 mA	11.9 mA	
0.0 mA	1.0 mA	2.0 mA	3.0 mA	4.0 mA	5.0 mA	6.0 mA	7.0 mA	8.0 mA	9.0 mA	10.0 mA	11.0 mA	12.0 mA																																																																																																																								
0.1 mA	1.1 mA	2.1 mA	3.1 mA	4.1 mA	5.1 mA	6.1 mA	7.1 mA	8.1 mA	9.1 mA	10.1 mA	11.1 mA	12.1 mA																																																																																																																								
0.2 mA	1.2 mA	2.2 mA	3.2 mA	4.2 mA	5.2 mA	6.2 mA	7.2 mA	8.2 mA	9.2 mA	10.2 mA	11.2 mA	12.2 mA																																																																																																																								
0.3 mA	1.3 mA	2.3 mA	3.3 mA	4.3 mA	5.3 mA	6.3 mA	7.3 mA	8.3 mA	9.3 mA	10.3 mA	11.3 mA	12.3 mA																																																																																																																								
0.4 mA	1.4 mA	2.4 mA	3.4 mA	4.4 mA	5.4 mA	6.4 mA	7.4 mA	8.4 mA	9.4 mA	10.4 mA	11.4 mA	12.4 mA																																																																																																																								
0.5 mA	1.5 mA	2.5 mA	3.5 mA	4.5 mA	5.5 mA	6.5 mA	7.5 mA	8.5 mA	9.5 mA	10.5 mA	11.5 mA	12.5 mA																																																																																																																								
0.6 mA	1.6 mA	2.6 mA	3.6 mA	4.6 mA	5.6 mA	6.6 mA	7.6 mA	8.6 mA	9.6 mA	10.6 mA	11.6 mA	12.6 mA																																																																																																																								
0.7 mA	1.7 mA	2.7 mA	3.7 mA	4.7 mA	5.7 mA	6.7 mA	7.7 mA	8.7 mA	9.7 mA	10.7 mA	11.7 mA	12.7 mA																																																																																																																								
0.8 mA	1.8 mA	2.8 mA	3.8 mA	4.8 mA	5.8 mA	6.8 mA	7.8 mA	8.8 mA	9.8 mA	10.8 mA	11.8 mA																																																																																																																									
0.9 mA	1.9 mA	2.9 mA	3.9 mA	4.9 mA	5.9 mA	6.9 mA	7.9 mA	8.9 mA	9.9 mA	10.9 mA	11.9 mA																																																																																																																									

Programming Simulation Screen																																																														
10		The Pulse Width Button displays the Pulse Width in microseconds. When the Pulse Width Adjusters are selected, the Pulse Width changes to the next/previous valid value. Alternatively, select the Pulse Width button to display the valid values (shown below), then select the desired value. Valid Pulse Width values range from 10 to 450 μs in steps of 10 μs. <table><tr><td></td><td>100 μs</td><td>200 μs</td><td>300 μs</td><td>400 μs</td></tr><tr><td>10 μs</td><td>110 μs</td><td>210 μs</td><td>310 μs</td><td>410 μs</td></tr><tr><td>20 μs</td><td>120 μs</td><td>220 μs</td><td>320 μs</td><td>420 μs</td></tr><tr><td>30 μs</td><td>130 μs</td><td>230 μs</td><td>330 μs</td><td>430 μs</td></tr><tr><td>40 μs</td><td>140 μs</td><td>240 μs</td><td>340 μs</td><td>440 μs</td></tr><tr><td>50 μs</td><td>150 μs</td><td>250 μs</td><td>350 μs</td><td>450 μs</td></tr><tr><td>60 μs</td><td>160 μs</td><td>260 μs</td><td>360 μs</td><td></td></tr><tr><td>70 μs</td><td>170 μs</td><td>270 μs</td><td>370 μs</td><td></td></tr><tr><td>80 μs</td><td>180 μs</td><td>280 μs</td><td>380 μs</td><td></td></tr><tr><td>90 μs</td><td>190 μs</td><td>290 μs</td><td>390 μs</td><td></td></tr></table>		100 μ s	200 μ s	300 μ s	400 μ s	10 μ s	110 μ s	210 μ s	310 μ s	410 μ s	20 μ s	120 μ s	220 μ s	320 μ s	420 μ s	30 μ s	130 μ s	230 μ s	330 μ s	430 μ s	40 μ s	140 μ s	240 μ s	340 μ s	440 μ s	50 μ s	150 μ s	250 μ s	350 μ s	450 μ s	60 μ s	160 μ s	260 μ s	360 μ s		70 μ s	170 μ s	270 μ s	370 μ s		80 μ s	180 μ s	280 μ s	380 μ s		90 μ s	190 μ s	290 μ s	390 μ s											
	100 μ s	200 μ s	300 μ s	400 μ s																																																										
10 μ s	110 μ s	210 μ s	310 μ s	410 μ s																																																										
20 μ s	120 μ s	220 μ s	320 μ s	420 μ s																																																										
30 μ s	130 μ s	230 μ s	330 μ s	430 μ s																																																										
40 μ s	140 μ s	240 μ s	340 μ s	440 μ s																																																										
50 μ s	150 μ s	250 μ s	350 μ s	450 μ s																																																										
60 μ s	160 μ s	260 μ s	360 μ s																																																											
70 μ s	170 μ s	270 μ s	370 μ s																																																											
80 μ s	180 μ s	280 μ s	380 μ s																																																											
90 μ s	190 μ s	290 μ s	390 μ s																																																											
11		The Rate Button displays the Rate in hertz. When the Rate Adjusters are selected, the Total Rate changes to the next/previous value of the defined Valid Rate values. The Valid Rate values range from 2 to 255 Hz in varying increments. <table><tr><td>2 Hz</td><td>37 Hz</td><td>74 Hz</td><td>116 Hz</td><td>154 Hz</td><td>208 Hz</td></tr><tr><td>6 Hz</td><td>40 Hz</td><td>76 Hz</td><td>119 Hz</td><td>159 Hz</td><td>223 Hz</td></tr><tr><td>10 Hz</td><td>45 Hz</td><td>79 Hz</td><td>128 Hz</td><td>170 Hz</td><td>227 Hz</td></tr><tr><td>15 Hz</td><td>50 Hz</td><td>85 Hz</td><td>130 Hz</td><td>174 Hz</td><td>231 Hz</td></tr><tr><td>17 Hz</td><td>53 Hz</td><td>89 Hz</td><td>132 Hz</td><td>176 Hz</td><td>238 Hz</td></tr><tr><td>20 Hz</td><td>57 Hz</td><td>95 Hz</td><td>136 Hz</td><td>179 Hz</td><td>255 Hz</td></tr><tr><td>22 Hz</td><td>60 Hz</td><td>99 Hz</td><td>139 Hz</td><td>185 Hz</td><td></td></tr><tr><td>26 Hz</td><td>63 Hz</td><td>104 Hz</td><td>143 Hz</td><td>190 Hz</td><td></td></tr><tr><td>30 Hz</td><td>66 Hz</td><td>112 Hz</td><td>146 Hz</td><td>198 Hz</td><td></td></tr><tr><td>34 Hz</td><td>71 Hz</td><td>113 Hz</td><td>149 Hz</td><td>204 Hz</td><td></td></tr></table>	2 Hz	37 Hz	74 Hz	116 Hz	154 Hz	208 Hz	6 Hz	40 Hz	76 Hz	119 Hz	159 Hz	223 Hz	10 Hz	45 Hz	79 Hz	128 Hz	170 Hz	227 Hz	15 Hz	50 Hz	85 Hz	130 Hz	174 Hz	231 Hz	17 Hz	53 Hz	89 Hz	132 Hz	176 Hz	238 Hz	20 Hz	57 Hz	95 Hz	136 Hz	179 Hz	255 Hz	22 Hz	60 Hz	99 Hz	139 Hz	185 Hz		26 Hz	63 Hz	104 Hz	143 Hz	190 Hz		30 Hz	66 Hz	112 Hz	146 Hz	198 Hz		34 Hz	71 Hz	113 Hz	149 Hz	204 Hz	
2 Hz	37 Hz	74 Hz	116 Hz	154 Hz	208 Hz																																																									
6 Hz	40 Hz	76 Hz	119 Hz	159 Hz	223 Hz																																																									
10 Hz	45 Hz	79 Hz	128 Hz	170 Hz	227 Hz																																																									
15 Hz	50 Hz	85 Hz	130 Hz	174 Hz	231 Hz																																																									
17 Hz	53 Hz	89 Hz	132 Hz	176 Hz	238 Hz																																																									
20 Hz	57 Hz	95 Hz	136 Hz	179 Hz	255 Hz																																																									
22 Hz	60 Hz	99 Hz	139 Hz	185 Hz																																																										
26 Hz	63 Hz	104 Hz	143 Hz	190 Hz																																																										
30 Hz	66 Hz	112 Hz	146 Hz	198 Hz																																																										
34 Hz	71 Hz	113 Hz	149 Hz	204 Hz																																																										
12	 Open State  Anode State	IPG (Internal Pulse Generator) icon. Open State: Bipolar, Inactive. Anode State: Monopolar.																																																												


Programming Simulation Screen		
13		The Current Equalize button is enabled when the Total Current is not distributed equally across the Lead Contacts in either the Anode or Cathode State. When selected, the Total Current will equalize across all active Lead Contacts.
14		The Current Mode button allows you to toggle between % Mode and mA Mode . The Current Mode button displays the Current Mode that will be in effect if selected, not the Current Mode that is displayed.
15	 E-Troll Control  Initiate E-Troll Button	The E-Troll Control up/down arrows are used to methodically move stimulation along the main axis of the Lead. The E-Troll Control is accessible after selecting the Initiate E-Troll Button. The Initiate E-Troll button is visible in mA mode, and in % Mode when either bipolar stimulation is ongoing or when multiple monopoles are not adjacent. Refer to “Using the 3D Visualization Panel” on page 90.
16	3D Visualization Panel	If registration is complete, the 3D Visualization Panel will display the Lead relative to the atlas, otherwise a vertical 3D Lead will be used with no atlas shown.
17		The View menu opens to allow you to toggle viewing of the 3D objects.
18	Stimulation Field Model (SFM)	Displays the volume of tissue surrounding the Lead which is likely to be activated for the given stimulation parameters. This volume is based on neuronal and finite element modeling.

en

Using the Simulation Panel

The Simulation Panel can store up to four Programmed Simulations. Each Programmed Simulation is stored within a Slot, and each Slot displays:

- A plus sign when the IPG is set to an anode.
- The Current value for each contact.
- The total Current value.
- The Pulse Width value.
- The Rate value.
- The **Simulation** number (1 through 4).

Simulation Slots can be added, deleted, and duplicated.

Adding a Simulation Slot


1. Select  **Add Simulation Slot**.

Deleting a Simulation Slot

1. Select the desired  Simulation Slot if not selected already.
2. Select the  Simulation Slot again to open the Simulations menu. The menu will include the following options:
 - **Duplicate** – when there is an Empty Simulation Slot displayed.
 - **Delete** – when there are two or more Programmed Simulation Slots displayed.
3. Select **Delete** from the **Simulations** menu.

Duplicating a Simulation Slot

1. Select the desired  Simulation Slot if not selected already.
2. Select the  Simulation Slot again to open the Simulations menu.

3. Select **Duplicate** from the **Simulations** menu.
4. The duplicated simulation will show up in the next open simulation slot.

Configuring the Programming of a Simulation Slot


1. Select

Left Simulations

 to view the **Left Lead Programing Simulation** screen or select

Right Simulations

 to view the **Right Lead Programing Simulation** screen.

2. Select a Simulation Slot.
3. Configure the programming for the Simulation Slot using the programming tools on the **2D Parameters Panel**. See “Using the 2D Parameters Panel” on page 81.
4. View the simulation in the **3D Visualization Panel**. See “Using the 3D Visualization Panel” on page 90.

Using the 2D Parameters Panel

For each Simulation you will need to program the contact configuration, pulse width, and rate.

Programming the Contact Configuration

Using % Mode

From % Mode, you can adjust individual contact percentages and overall Amplitude.


1. While in  **% Mode**, select each Active Contact until it is set to the desired state:

-  Off
-  Cathode
-  Anode.


Each successive selection will cycle the state from Off to Cathodic to Anodic and back to Off.


2. For each Active Contact, select the  **Amplitude** button to open the Amplitude Values Table. Select the desired total amplitude from the table of values.

Or

Select the **Amplitude** up/down adjustment arrows  to increase/decrease the Amplitude.

Or

Hover the Mouse cursor over the  **Amplitude** button and use the scroll wheel to increase/decrease the Amplitude.

3. For each Active Contact, select the arrows next to the contact  to increase/ decrease the percentage of the current Amplitude being simulated to the contact.

Or

Hover the Mouse cursor over the contact and use the scroll wheel to increase/decrease the percentage of the current amplitude being simulated to the contact.

Using mA Mode

From **mA Mode**, you can adjust individual contact Amplitudes.

1. Select the  **mA Mode** button to enter mA Mode.

Note: *The Amplitude button is disabled in mA Mode.*

2. Select each Active Contact until it is set to the desired state:

-  Off
-  Cathode
-  Anode.

Each successive selection will cycle the state from Off to Cathodic to Anodic and back to Off.

3. For each Active Contact, select the arrows next to the contact  to increase/decrease the Amplitude.

Or

Hover the Mouse cursor over the desired contact and use the scroll wheel to increase/decrease the Amplitude.

Using E-Troll™ Mode


E-Troll Mode is used to move the current up/down the Lead in small increments.

To use E-Troll Mode:

- Initiate

E-Troll


1. Select the the **Initiate E-Troll** button to enter E-Troll Mode.

- Initiate

Initiate E-Troll

2. When the **Initiate E-Troll** button is selected:
- The current Mode is set to %.
 - Lead Contact 1 is set to the Cathode State.
 - All other Lead contacts are set to Open State.
 - The IPG is set to the Anode State.
 - If the Total Current is greater than 12.7 mA, it will be set to 12.7 mA.
 - If the Total Current is equal to 0.0 mA, it will be set to 0.1 mA.
 - The E-Troll Control is displayed.

3. Adjust the Amplitude. The Amplitude is initially set at 0.1 mA. Select the

Amplitude button to open the Amplitude Values Table. Select the desired total amplitude from the table of values.

Or

Select the **Amplitude** up/down adjustment arrows

to increase/decrease the Amplitude.

Or

Hover the Mouse cursor over the

Amplitude button and use the scroll wheel to increase/decrease the Amplitude.


Note: The E-Troll Control also displays while in % Mode with a single monopole.

4. Select the  E-Troll™ Control up/down arrows to move the percentages of the current up/down the contact in small increments.

Or

Hover the Mouse cursor over the 

E-Troll Control and use the scroll wheel to move the percentages up/down the contact.

Notes:

- For any programming method, the total Amplitude across all Active Contacts cannot exceed 20 mA, and the total on any individual Active Contact cannot exceed 12.7 mA.
- For any programming mode, the total of all Active Anodes must equal the total of all Active Cathodes.

- At any time, the **Equalize** button can be selected to equally distribute the total percentage of Amplitude across all Active Anodes and Cathodes, depending on the mode.

- The **IPG** icon automatically will be set as the Anode in a monopolar configuration.

Using the E-Troll Control

Up Arrow:

When the **E-Troll™** Control is enabled and the up arrow selected, the upper active E-Troll Lead Contact increases and the lower active Lead Contact decreases according to the E-Troll Lead Contact Current Rules. When the upper active E-Troll Lead Contact equals the Total Current the lower active contact changes to the Open state. If the Lead Contact equals the Total Current, with the next up arrow press the Lead Contact above it is set to a Cathode State with a current value set according to the E-Troll Lead Contact Current Rules.




Down Arrow:

When the  **E-Troll™** Control is enabled and the down arrow selected, the lower active E-Troll Lead Contact increases and the upper active Lead Contact decreases according to the E-Troll Lead Contact Current Rules. When the lower active E-Troll Lead Contact equals the Total Current the upper active contact changes to the Open State. If the Lead Contact equals the Total Current, with the next down arrow press the Lead Contact below it is set to a Cathode State with a current value set according to the E-Troll Lead Contact Current Rules.




Lead Contact Current Rules

Total Current	inc/dec mA Value
≤ 2.0 mA	0.1 mA
> 2.0 mA and < 4.1 mA	0.2 mA
> 4.0 mA and < 6.1 mA	0.3 mA
> 6.0 mA and < 8.1 mA	0.4 mA
> 8.0 mA and < 10.1 mA	0.5 mA
> 10.0 mA and < 12.8 mA	0.6 mA

Programming the Pulse Width

- 
1. Select the **Pulse Width** button to open the Pulse Width Values Table. Select the desired Pulse Width from the table of values.
- Or
- 
- Select the **Pulse Width** up/down adjustment arrows to increase/decrease the Pulse Width.
- Or
- 
- Hover the Mouse cursor over the **Pulse Width** button and use the scroll wheel to increase/decrease the Pulse Width.

Programming the Rate

- 
1. Select the **Rate** button to open the Rate Values Table. Select the desired Rate from the table of values.
- Or
- 
- Select the **Rate** up/down adjustment arrows to increase/decrease the Rate.
- Or
- 
- Hover the Mouse cursor over the **Rate** button and use the scroll wheel to increase/decrease the Rate.


en

⚠ Charge Density Warning

WARNING: High levels of stimulation may damage brain tissue. The GUIDE™ DBS Software includes a warning when the stimulation results in a charge density of 30 $\mu\text{C}/\text{cm}^2/\text{phase}$ or greater to minimize the likelihood of injury. However, the software does allow the physician to exceed this limit at his/her discretion.


The following warning message displays when the charge density is exceeded:

Charge Density Exceeded. Parameter settings for the simulation may cause tissue damage and should be lowered.

Lower current on contacts or decrease pulse width.

The graph below displays the recommended maximum charge density of 30 $\mu\text{C}/\text{cm}^2/\text{phase}$ for different combinations of current amplitude (mA) and pulse width (μs) for Boston Scientific DBS Leads.

Display Limit Reached Warning

The Stimulation Field Models (SFMs) have only been verified within certain spatial extents. When the stimulation parameters create an SFM that goes beyond these extents, the display will be hidden and the following warning message will display:

Display Limit Reached. The selected parameters result in SFMs that exceed the verified SFM boundaries. Please reduce amplitude and/or pulse width to return to display.


Display Limit Reached

The selected parameters result in SFMs that exceed the verified SFM boundaries. Please reduce amplitude and/or pulse width to return to display.

Using the 3D Visualization Panel

- en
1. Select **Left Simulations** to view the **Left Lead Programing Simulation** screen or select **Right Simulations** to view the **Right Lead Programing Simulation** screen.

2. Select the **View** menu to toggle viewing of the following structures:
 - Anterior Commissure (AC)
 - Encapsulation Layer
 - Posterior Commissure (PC)
 - Red Nucleus (RN)
 - Reference Volume
 - Stimulation Field Model (SFM)
 - Subthalamic Nucleus (STN)
 - Thalamus

3. Change the view using the  vertical and  horizontal Viewport Slider Controls, or by grabbing and dragging the image, to view from any angle looking in towards the center of the contact areas on the Lead.

To assist with orientation, a reference head displays in the lower right corner. The head rotates with the displayed perspective.

4. Zoom in using the  button . Zoom out using the  button. Or, zoom in/out by hovering the Mouse over the view and scrolling the wheel up/down.

Creating a Simulation Report

Select **Report** on the **Simulations** screen. The Simulation Report displays.

Report Patient Information

Report Hardware Information

Report Simulation Settings

Menu

Roy Jenkins

Left Simulations

Right Simulations

Report

Patient Information

Hardware Information

IPG Simulation Settings

Report Export Button

Patient Information		Date of Birth	Patient ID
Roy Jenkins		01 MAY 1963	201
Loaded MR Series		MR Images	MR Comments
Roy Jenkins		01 JUL 2009	124
Loaded CT Series		CT Images	Post SPGR AX
Roy Jenkins		20 JUL 2009	136
		CT Comments	
Hardware Information			
Model	Date of Implant	Location	Electrode Pz
IPG	Vertex™	-	
Lead 1	1x8	Left	E1-E8
Lead 2	1x8	Right	E9-E16
IPG Simulation Settings			
Left Lead-Simulation 1		Left Lead-Simulation 2	
Contacts	Amplitude	Contacts	Amplitude
8	3.0 mA	8	2.3 mA
7	Pulse Width	7	8 + 54 %
6	60 µs	6	2.8 mA
5	Rate	5	16 + 58 %
4	130 Hz	4	1.2 mA
3	3 + 48 %	3	1.9 mA
2	2	2	1.2 mA
1	1	1	1.2 mA
Left Lead-Simulation 3		Left Lead-Simulation 4	
Contacts	Amplitude	Contacts	Amplitude
8	3.0 mA	8	2.3 mA
7	Pulse Width	7	8 + 54 %
6	60 µs	6	2.8 mA
5	Rate	5	16 + 58 %
4	130 Hz	4	1.2 mA
3	3 + 48 %	3	1.9 mA
2	2	2	1.2 mA
1	1	1	1.2 mA
Right Lead-Simulation 1		Right Lead-Simulation 2	
Contacts	Amplitude	Contacts	Amplitude
16	1.3 mA	16	1.3 mA
15	15 - 62 %	15	15 - 62 %
14	60 µs	14	60 µs
13	Rate	13	Rate
12	130 Hz	12	130 Hz
11	11	11	11
10	10	10	10
9	9	9	9
Right Lead-Simulation 3		Right Lead-Simulation 4	
Contacts	Amplitude	Contacts	Amplitude
16	1.3 mA	16	1.3 mA
15	15 - 62 %	15	15 - 62 %
14	60 µs	14	60 µs
13	Rate	13	Rate
12	130 Hz	12	130 Hz
11	11	11	11
10	10	10	10
9	9	9	9

Note: GUIDE™ DBS is a simulation tool and should not be a substitute for standard clinical practice.

Note: **Report** may not be visible on the **Simulations** screen if Image Registration was skipped. To register a Patient's MR and CT, see "Image Registration" on page 30.

Exporting a Simulation Report


1. Ensure the removable media (Flash Drive) is connected (to the right USB port).
2. From the Simulation **Report** screen, select **Export**.
3. The Report is saved to the connected removable media.

The format of the report name is LASTNAME_FIRSTNAME_YYYY_MM_DD_HH_MM_SS.PDF.

Or

If the removable media is not connected to GUIDE™ DBS, the following message displays:

The removable media could not be found and the report was not exported.


Select **Close**.

Or

If the removable media disk is full, the following message displays:

The removable media is full and the report was not exported.


Select **Close**.

Or

If export fails for another reason than the two listed above, the following message displays:

There was an export error and the report was not exported.

Select .

4. The Simulation Report may display the following information:

Patient Information

- Name
- Date of Birth
- Patient ID
- Loaded MR Series (Name)
- MR Date
- MR Images (Count)
- MR Comments
- Loaded CT Series (Name)
- CT Date
- CT Images (Count)
- CT Comments

Hardware Information

- IPG Model, Date of Implant
- Lead 1 Model, Date of Implant, Location, Electrode #s
- Lead 2 Model, Date of Implant, Location, Electrode #s

IPG Simulation Settings

- Left Lead-Simulation # (where # equals the number of the simulation 1-4)
Contacts, Amplitude, Pulse Width, Rate
- Right Lead-Simulation # (where # equals the number of the simulation 1-4)
Contacts, Amplitude, Pulse Width, Rate

CAUTION: GUIDE™ DBS is a simulation tool and should not be a substitute for standard clinical practice.

Roy Jenkins

Patient Information

Name		Date of Birth	Patient ID
Roy Jenkins		01 MAY 1963	201
Loaded MR Series	MR Date	MR Images	MR Comments
Roy Jenkins	01 JUL 2009	124	Post SPGR AX
Loaded CT Series	CT Date	CT Images	CT Comments
Roy Jenkins	20 JUL 2009	136	

Hardware Information

	Model	Date of Implant	Location	Electrode #s
IPG	Vercise™	-		
Lead 1	1x8	-	Left	E1-E8
Lead 2	1x8	-	Right	E9-E16

IPG Simulation Settings

Left Lead-Simulation 1		Left Lead-Simulation 2		Left Lead-Simulation 3		Left Lead-Simulation 4	
Contacts	Amplitude	Contacts	Amplitude	Contacts	Amplitude	Contacts	Amplitude
8	3.0 mA	8	2.3 mA	8 + 54 %	2.8 mA	8	1.9 mA
7	Pulse Width	7	Pulse Width	7 + 46 %	Pulse Width	7	Pulse Width
6	60 µs	6	60 µs	6	60 µs	6	60 µs
5 - 100 %	Rate	5 + 52 %	Rate	5	Rate	5 - 37 %	Rate
4	130 Hz	4 - 100 %	130 Hz	4	130 Hz	4 - 21 %	130 Hz
3		3 + 48 %		3		3 - 42 %	
2		2		2 - 36 %		2	
1		1		1 - 64 %		1	

Right Lead-Simulation 1		Right Lead-Simulation 2		Right Lead-Simulation 3		Right Lead-Simulation 4	
Contacts	Amplitude	Contacts	Amplitude	Contacts	Amplitude	Contacts	Amplitude
16	1.3 mA	16	1.5 mA	16	2.1 mA	16 + 58 %	1.2 mA
15 - 62 %	Pulse Width	15	Pulse Width	15	Pulse Width	15 - 42 %	Pulse Width
14 - 38 %	60 µs	14	60 µs	14	60 µs	14 - 58 %	60 µs
13	Rate	13 - 100 %	Rate	13 + 48 %	Rate	13 + 42 %	Rate
12	130 Hz	12 + 100 %	130 Hz	12 - 33 %	130 Hz	12	130 Hz
11		11		11 + 52 %		11	
10		10		10 - 67 %		10	
9		9		9		9	

Note: GUIDE™ DBS is a simulation tool and should not be a substitute for standard clinical practice.

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The screenshot shows a dark-themed form titled "Fields marked with an '*' are required." The form is organized into sections. The first section contains fields for "Patient Name" (split into "Last" and "First"), "Patient ID", and "Date of Birth" (with a date picker). The second section, marked with a "+" icon, contains "IPG Model" (with a dropdown menu showing "Vercise™") and "Date of Implant" (with a date picker). The third section, marked with a "1" icon, is for "Lead 1" and includes "Lead Model" (dropdown with "1x8"), "Lead Location" (dropdown with "Select one..."), "Electrode #'s" (dropdown with "E1-E8"), and "Date of Implant" (date picker). The fourth section, marked with a "+" icon, is for "Lead 2" and includes "Lead Model" (dropdown with "Select one..."), "Lead Location" (dropdown with "Select one..."), "Electrode #'s" (dropdown with "Select one..."), and "Date of Implant" (date picker).

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or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
© 2018 Boston Scientific Corporation