# **Intrasept™**

Intraosseous Nerve Ablation System

INSTRUCTIONS FOR USE

|    |                      |    |
|----|----------------------|----|
| en | Instructions for Use | 2  |
| es | Instrucciones de uso | 10 |
| de | Gebrauchsanweisung   | 18 |
| it | Istruzioni per l'uso | 26 |
| fr | Mode d'emploi        | 34 |

TABLE OF CONTENTS

INDICATIONS FOR USE..... 3

INTENDED PURPOSE..... 3

PATIENT POPULATION ..... 3

INTENDED USER..... 3

INTENDED USE ENVIRONMENT..... 3

DEVICE DESCRIPTION ..... 3

REQUIRED EQUIPMENT ..... 3

HOW SUPPLIED..... 3

    INSTRUCTIONS FOR USE..... 4

MATERIALS ..... 4

STORAGE ..... 4

CONTRAINDICATIONS..... 4

WARNINGS..... 4

PRECAUTIONS ..... 4

CLINICAL BENEFITS ..... 4

POTENTIAL ADVERSE EVENTS ..... 5

DIRECTIONS FOR USE ..... 5

DISPOSAL..... 6

POST PROCEDURE..... 6

LIMITATION OF LIABILITY ..... 7

DISCLAIMER OF LIABILITY FOR DAMAGES..... 7

FIGURES ..... 7

    FIGURE 1..... 7

    FIGURE 2 ..... 7

    FIGURE 3 ..... 7

    FIGURE 4..... 8

    FIGURE 5 ..... 8

    FIGURE 6..... 8

SYMBOLS USED ON LABELING..... 9

    Symbol..... 9

    Definition ..... 9

# Intrasept™

## Intraosseous Nerve Ablation System

### INSTRUCTIONS FOR USE

**Caution:** This device is restricted to sale, distribution, and use by or on the order of a physician.

#### INDICATIONS FOR USE

The Intrasept™ Intraosseous Nerve Ablation System is intended to be used in conjunction with radiofrequency (RF) generators for the ablation of basivertebral nerves of the L3 through S1 vertebrae for the relief of chronic low back pain of at least six months duration that has not responded to at least six months of conservative care, and is also accompanied by features consistent with Type 1 or Type 2 Modic changes on an MRI such as inflammation, edema, vertebral endplate changes, disruption and fissuring of the endplate, vascularized fibrous tissues within the adjacent marrow, hypointensive signals (Type 1 Modic change), and changes to the vertebral body marrow including replacement of normal bone marrow by fat, and hyperintensive signals (Type 2 Modic change).

#### INTENDED PURPOSE

The intended purpose of Intrasept Intraosseous Nerve Ablation System (RF Probes and Access Instruments) is to treat pain by performing RF ablation procedures of the basivertebral nerve within the vertebral body.

#### PATIENT POPULATION

The Intrasept Intraosseous Nerve Ablation System is intended to be used on patients with chronic low back pain, specifically with Type 1 or Type 2 Modic changes.

#### INTENDED USER

The intended users of the Intrasept System are physicians. Before using the Intrasept System, users should read all product labeling. If additional support is needed, call your Relievant Medsystems representative.

#### INTENDED USE ENVIRONMENT

The Intrasept System shall be used in healthcare facilities equipped to provide a sterile environment for standard minimally invasive spine treatments or surgeries such as operating rooms, outpatient or inpatient surgical suites, procedure rooms or equivalent.

#### DEVICE DESCRIPTION

The Intrasept Intraosseous Nerve Ablation System ("Intrasept System") is designed to ablate the basivertebral nerve ("BVN") within a vertebral body. The Intrasept Access Instruments are used to create a path to a treatment location within a vertebral body to allow for targeted deployment of the Intrasept RF Probe for the ablation of the basivertebral nerve. The Intrasept System consists of the following components:

- Intrasept Access Instruments ("Instruments" - see Figure 1):
  - Introducer made up of the Introducer Cannula and the Stylet (either Diamond or Bevel)
  - Curved Cannula Assembly made up of the Curved Cannula and J-Stylet
  - Straight Stylet
  - Drill
- Intrasept RF Probe ("Probe" - see Figure 2)

The Stylets (Diamond and Bevel) are labeled with "1", the J-Stylet is labeled with "2", and the Straight Stylet is labeled with "3" all on their respective blue handle surfaces; the numbers indicate the sequence of use. The Introducer Cannula is packaged pre-assembled with a Stylet (Diamond or Bevel) to form the Diamond or Bevel Introducer, and one Curved Cannula is packaged pre-assembled with a J-Stylet to form the Curved Cannula Assembly. The Drill is individually packaged. The Bevel Stylet, Introducer Cannula, and J-Stylet all exhibit arrows that indicate the directionality of the Instrument. The Stylets (Diamond and Bevel), Introducer Cannula, J-Stylet, Straight Stylet, and Drill are designed to be visible under fluoroscopic or CT imaging.

#### REQUIRED EQUIPMENT

In addition to the Instruments and Probe, the Intrasept System requires the use of the Intrasept Radiofrequency Generator (RFG) to provide the energy needed to ablate the tissue. The RFG provides the following treatment parameters when used with the Intrasept RF Probe:

| Ablation Parameter Set #1                | Ablation Parameter Set #2                 |
|------------------------------------------|-------------------------------------------|
| • Temperature: 75°C                      | • Temperature: 85°C                       |
| • Temperature ramp: 0.5°C/second         | • Temperature ramp: 1°C/second            |
| • Ablation time: 7 minutes (420 seconds) | • Ablation time: 15 minutes (900 seconds) |
| • Ablation zone: ~5 mm radially          | • Ablation zone: ~6 mm radially           |

#### HOW SUPPLIED

The Access Instruments and Probe are packaged in a dual barrier tray system within a carton.

The outer tray has a single sterile barrier. The inner tray is sterile and contains the Access Instruments and Probe.

The Access Instruments, Probe, and inner tray are sterilized with irradiation and intended for single patient use only.

Do not use if package is opened or damaged, as product integrity and/or sterility may be compromised. If damage is found or if the outer tray is opened, call your Relievant Medsystems representative.

The Intrasept Access Instruments are available in the following kitting configurations:

- |               |                                       |               |                                       |
|---------------|---------------------------------------|---------------|---------------------------------------|
| • RLV 0045-EU | Intrasept Access Instrument kit, 2 VB | • RLV 0055-EU | Intrasept Access Instrument kit, 3 VB |
|               | One (1) Bevel Introducer              |               | One (1) Bevel Introducer              |
|               | One (1) Diamond Introducer            |               | One (1) Diamond Introducer            |
|               | One (1) Curved Cannula Assembly       |               | One (1) Curved Cannula Assembly       |
|               | One (1) Curved Cannula                |               | Two (2) Curved Cannula                |
|               | One (1) Straight Stylet               |               | One (1) Straight Stylet               |
| • RLV 0075-EU | Additional Vertebral Body (VB) kit    | • RLV 0035-EU | Intrasept Access Instruments – Drill  |
|               | One (1) Introducer Cannula            |               | One (1) Drill                         |
|               | One (1) Curved Cannula                |               |                                       |
| • RLV 0015-EU | Intrasept RF Probe                    |               |                                       |
|               | One (1) Probe                         |               |                                       |

The Intrasept Access Instruments kit may be used to access up to two or three vertebral bodies (depending on kit configuration) in a single patient. The Introducer Cannula, Diamond Stylet, Bevel Stylet, J-Stylet, Straight Stylet, and Drill may be used to treat up to a maximum of four (4) vertebrae (L3, L4, L5, S1) in a single patient, while the Curved Cannula may only be used on one (1) vertebra in a single patient.

The Probe may be used to treat up to a maximum of four (4) vertebrae (L3, L4, L5, S1) on a single patient.

### INSTRUCTIONS FOR USE

The Instructions for Use (IFU) for this product are supplied in electronic form over the Internet.

Visit [www.IFU-BSCI.com](http://www.IFU-BSCI.com) to access the IFU using your internet browser. Information from the product labeling such as Product Code, UDI-DI or Product Name may be used to search for the specific IFU.

If you have difficulty accessing the IFU online, or would prefer to receive a paper copy, call your Relevant Medsystems representative. A copy will be sent to you at no charge and should arrive within seven days.

### MATERIALS

The Probe and Access Instruments may contain cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1% weight by weight.

**Note:** These devices are made with stainless steel which may contain cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

### STORAGE

The Intrasept Access Instruments and Intrasept RF Probe should be stored in a cool, dry place.

### CONTRAINDICATIONS

Use of the Intrasept Intraosseous Nerve Ablation System is contraindicated in:

- Patients with severe cardiac or pulmonary compromise
- Patients with active implantable pulse generators (e.g. pacemakers, defibrillators)
- Patients where the targeted ablation zone is < 10 mm away from a sensitive structure not intended to be ablated, including the vertebral foramen (spinal canal)
- Patients with active systemic infection or local infection in the area to be treated
- Patients who are pregnant
- Skeletally immature patients (generally ≤ 18 years of age)

---

### WARNINGS

- Do not use if package is opened or damaged as product integrity and/or sterility may be compromised. If damage is found or if the outer tray is opened, call your Relevant Medsystems representative.
- Do not use after the expiration date has passed as product integrity and/or sterility may be compromised.
- For SINGLE PATIENT USE ONLY. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Reconditioning, refurbishing, repair or modification of the device to enable further use is prohibited.
- Do not use this product if you have not been properly trained in its use. Physicians using the device should be familiar with the physiology and pathology of the selected anatomy to be treated and be trained in the performance of the chosen surgical technique. Improper surgical use and technique may lead to suboptimal clinical outcomes.
- As with any surgical instrument, careful attention must be exercised to ensure that excessive force is not placed on the Instruments or Probe. Excessive force may lead to product failure, causing device breakage, fragmentation, or loss of functionality.
- Read and understand the Instructions For Use (IFU) completely prior to use.
- The Intrasept System must be used with the Intrasept Radiofrequency Generator in order to provide the required treatment parameters.
- Failure of the RFG could result in an unintended increase of output power to the Probe.
- The Intrasept RF Probe may interfere and adversely influence the operation of other electronic equipment.

---

### PRECAUTIONS

- The Intrasept Access Instruments and the Intrasept RF Probe are single patient use only in a single procedure. The Introducer Cannula, Diamond Stylet, Bevel Stylet, J-Stylet, Straight Stylet, Drill, and Probe may be used to treat up to a maximum of four (4) vertebrae (L3, L4, L5, S1) on a single patient, while the Curved Cannula may only be used on one (1) vertebra on a single patient.
- Prior to the procedure, CT or MRI imaging must be utilized to help define the desired treatment site and define working anatomical landmarks that can be used for mapping access to the treatment site.
- The device should be manipulated only while under fluoroscopic or CT observation.
- DO NOT use this device in the presence of flammable anesthetics, other flammable gases or objects, near flammable fluids such as skin prepping agents and tinctures, or oxidizing agents. Observe appropriate fire precautions at all times. There is a risk of pooling of flammable solutions under the patient or in body depressions and cavities.
- Fluids pooled in the body depressions and cavities of the patient should be mopped up before RFG is used.
- There is a danger of ignition of endogenous gases (e.g., cotton and gauze saturated with oxygen may be ignited by sparks produced during normal use of the RFG).
- Operator may choose to use smoke-plume extraction apparatus, though the indicated procedure is unlikely to produce noticeable smoke.
- The Probe is for use WITHOUT a neutral electrode (i.e., a grounding pad).
- Discontinue use if inaccurate, erratic or sluggish temperature readings are observed. Use of damaged equipment may cause patient injury.
- Prior to operation, visually inspect the Probe, Instruments, and RFG for physical damage, obvious cracks in insulation or loose parts.
- During power delivery, the cable should not come in direct contact with the patient's skin or other patient leads.
- Safety and effectiveness in patients with conditions that are associated with poor bone quality such as osteoporosis have not been established.

---

### CLINICAL BENEFITS

The clinical benefit of the Intrasept System includes improvement in pain with a potential reduction of medications. These benefits can relate to improved functional capacity and an improved quality of life.

## POTENTIAL ADVERSE EVENTS

The following adverse events or complications associated with this procedure include, but are not limited to:

- Nerve injury including puncture of the spinal cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis
- Breakage, slippage, misuse or mishandling of Instruments or Probe, such as on sharp edges, may cause burns or injury to the patient and surgical staff, as well as produce device fragments that may result in adverse effects such as nerve injury or allergic reactions if not retrieved
- Unintended puncture wounds including vascular puncture and dural tear
- Embolism of fat, thrombus or other materials resulting in symptomatic pulmonary embolism or other clinical sequelae
- Deep or superficial wound infection
- Hemorrhage
- Hematoma
- Pedicle or vertebral body fracture
- Risks related to having a surgical procedure/anesthesia, e.g. respiratory or cardiovascular problems, exposure to biohazardous material, death, radiation exposure
- Allergic reaction or inflammation
- Pain or discomfort

Refer to the Incept Radiofrequency Generator User Manual for additional adverse events that have been reported for specific radiofrequency ablation applications.

## DIRECTIONS FOR USE

1. Using preoperative MRI or CT imaging, identify a target treatment site in tissue for forming an ablation zone of up to 12 mm in diameter (6 mm radially from the Probe tip) at the terminus of the basivertebral foramen. Intraoperatively, the target will be identified by the measurements made in this planning process, noting the distances from the vertebra bony landmarks in all three axes. These measurements should be converted to percentages that can be used to determine location when using fluoroscopic imaging. In order to determine the safety boundary for the vertebral body, 10 mm is divided by the vertebral body length in the sagittal view.

---

**Warning:** The Incept System is contraindicated for targeted treatment sites < 10 mm from boundary of the vertebral body, or other sensitive structures not intended to be ablated, including the vertebral foramen (spinal canal).

---

2. Open the carton and the outer tray.

3. Transfer the inner tray or devices to the sterile field.

**Note:** If the surgery is interrupted, call your Relevant Medsystems representative, as the inner tray cannot be stored.

4. Patient must be positioned in prone recumbency. Prepare and drape the skin at the intended access site using standard sterile technique.

**Note:** The Incept procedure should be conducted under at least moderate conscious patient sedation.

5. Make a 3 - 5 mm skin incision with a scalpel blade over the area selected for initial transpedicular or extrapedicular access.

6. Using manual control, advance the Diamond or Bevel Introducer through the incision and dock onto the selected bone surface. With appropriate imaging observation, advance the Introducer using controlled manual force or a small mallet to the desired depth. For a transpedicular approach, the Introducer is advanced through the pedicle until the scalloped section in the Bevel or Diamond Stylet reaches the posterior cortical wall of the vertebral body as viewed in a lateral fluoroscopic image. The Drill may be used to aid pedicle access.

---

**Precaution:** Fluoroscopic or CT imaging must be utilized frequently throughout the procedure to ensure proper entry and placement of the Instruments and Probe. Improper entry and placement of the devices may result in patient injury and/or sub-optimal clinical treatment effects.

---

**Note:** Care should be taken when handling the sharp tip of the Bevel and Diamond Stylets to avoid injury.

**Note:** Graduated markings on the Introducer Cannula may be used as reference only, they are not intended to replace the use of fluoroscopic or CT observation.

7. Holding the Introducer, remove the Stylet by simultaneously pressing the release button and rotating the Stylet 180 degrees counter-clockwise (see Figure 3). Set the Stylet aside within the sterile field in case it is needed again.
8. Ensure the Gear Wheel on the Curved Cannula Assembly is held by the friction lock at the lowest thread position (rotate clockwise if necessary). Depress (open) the Bail tab (red indicator is visible) ensuring the J-Stylet tip is retracted into the Curved Cannula (see Figure 4). While maintaining Introducer Cannula position, grip the handle of the Curved Cannula Assembly and carefully insert the distal tip into the Introducer Cannula handle while aligning the tip of the Curved Cannula Assembly with the inner lumen of the Introducer Cannula. In one motion, feed the Curved Cannula Assembly into the Introducer Cannula while raising it upright as it slides into the Introducer Cannula (see Figure 5). Continue insertion until the Curved Cannula Assembly is fully seated and the arrow on the J-Stylet aligns with the arrows on the Introducer Cannula. Ensure the Bail tab is in the closed position (see Figure 5) prior to Curved Cannula Assembly deployment.

**Note:** Prior to inserting the Curved Cannula Assembly into the Introducer Cannula, ensure the Bail tab is in the open position (visible red indicator).

**Note:** Prior to deployment of the Curved Cannula Assembly into the vertebra, ensure the Bail tab is in the closed position.

9. Under fluoroscopic or CT observation determine whether the direction of the distal tip of the J-Stylet is directed toward the treatment location. The direction of distal tip deployment is indicated by directional arrows on the Introducer Cannula handle and Curved Cannula Assembly handle. If necessary, reposition the instruments to obtain proper trajectory. To deploy the Curved Cannula Assembly, rotate the Gear Wheel counterclockwise until it is frictionally locked in its highest thread position and in contact with the Curved Cannula handle.

**Note:** Ensure that the Introducer Cannula does not advance once the Curved Cannula Assembly has begun deployment. If inadvertent Introducer Cannula advancement is observed, it must be retracted prior to repositioning.

10. Begin deployment of the Curved Cannula Assembly using a small mallet to advance the assembly.

---

**Precaution:** In the event of inadvertent advancement of the Introducer Cannula during deployment of the Curved Cannula Assembly, care must be taken to prevent potential further device damage. This can occur if mallet strikes continue to be delivered when the Gear Wheel is in contact with the Introducer Cannula. This scenario can create a condition in which the Introducer Cannula advances along its straight trajectory, while the Curved Cannula is held in the bone along its curved trajectory. This can lead to the creation of a kink in the Curved Cannula material. If it is suspected that this occurred, or if retraction of the Curved Cannula is difficult, it is recommended that Introducer Cannula and Curved Cannula be removed as a single unit, rather than retracting the Curved Cannula by itself.

---

11. Continue advancement of the Curved Cannula Assembly using controlled mallet strikes and small rotational adjustments to adjust the trajectory as needed to guide the tip towards the target site; use both Anterior-Posterior (AP) and Lateral imaging to confirm trajectory. Continue advancement until the proximal section of the J-Stylet tip crosses the mid-sagittal plane of the vertebral body or the Curved Cannula Assembly is fully deployed, denoted by the Gear Wheel seated within the Introducer Cannula.

**Note:** If the Curved Cannula Assembly needs to be removed during deployment (in cases without inadvertent advancement of the Introducer Cannula as noted above), retract the Curved Cannula Assembly into the Introducer Cannula Assembly by turning the gearwheel in a clockwise direction until the Gear Wheel has traversed the threaded section of the Curved Cannula. Once fully retracted, remove the Curved Cannula Assembly and Introducer Cannula from the patient.

12. Rotate the Gear Wheel clockwise until it contacts the Introducer Cannula handle (if not already in contact). This will serve as a stop and will not allow further advancement of the Curved Cannula. Once the Gear Wheel is properly seated, with one hand hold the handle of the Introducer Cannula and the Gear Wheel together. Partially retract the J-Stylet, rotate it 180 degrees, and then fully remove from the Curved Cannula. Set aside the J-Stylet within the sterile field in case it's needed again.

**Note:** The Bail tab may be placed in the open position to help initiate removal of the J-Stylet.

---

**Precaution:** Care should be taken when removing the J-Stylet to prevent blood/tissue from ejecting from the tip of the J-Stylet as it comes free of the Curved Cannula. In addition to rotating the J-Stylet 180 degrees as it is removed, gauze may be placed over the slot in the Curved Cannula handle to prevent this action.

---

13. If additional deployment depth is required after fully deploying the Curved Cannula Assembly, the Straight Stylet may be used to extend the channel. Under imaging observation, insert the Straight Stylet into the Curved Cannula and use a small mallet to slowly advance the Straight Stylet to the target treatment site. Extend the Straight Stylet about 5 mm beyond the mid-sagittal plane of the vertebral body. This will create a channel for the Probe to traverse.
14. Withdraw the Straight Stylet from the Curved Cannula while maintaining the Curved Cannula and Introducer Cannula at the access site. Set aside the Straight Stylet within the sterile field in case it's needed again.
15. Connect the Probe to the RFG, matching the blue connector to the RFG's Instrument Connection Port (also blue) and aligning the arrow vertically (12 o'clock position).
16. Introduce the Probe into the Curved Cannula and advance the Probe distally beyond the Curved Cannula's tip to the targeted treatment site. The Probe should be advanced by hand, requiring minimal force to reach the desired position. If resistance is met during insertion, the Probe should be removed and the Straight Stylet introduced to clear the channel. Verify the Probe position using fluoroscopic or CT imaging in conjunction with the measurements taken during preoperative imaging, making sure that the probe remains  $\geq 10$  mm from boundary of the vertebral body in all directions.

---

**Precaution:** Both Anterior-Posterior (A-P) and lateral imaging should be reviewed to ensure that the Probe is positioned so that its tip transects the basivertebral foramen at the target treatment site, while remaining  $\geq 10$  mm from the boundary of the vertebral body, or other sensitive structures not intended to be ablated. Unintended damage to tissue can occur, with potential harm to the patient, if advanced beyond this target site.

---

17. Once target location has been confirmed, turn the Gear Wheel clockwise to retract the Curved Cannula until the proximal marker band on the Probe shaft has retracted below the surface of the Curved Cannula handle. This ensures that both Probe electrodes are properly exposed during ablation.

**Note:** If Probe electrodes are not properly exposed from the Curved Cannula, the size of the ablation zone may be reduced.

18. Re-confirm the placement of the Probe using fluoroscopic or CT imaging, as the Probe may move during the retraction of the Curved Cannula.
19. The nerve ablation process requires use of the Incept Radiofrequency Generator. Refer to the Incept Radiofrequency Generator User Manual for detailed instructions, warnings, precautions, and contraindications.
20. If the Probe is  $\leq 5$  mm from the intended target location, Ablation Parameter Set #1 may be selected; if  $\leq 6$  mm, Ablation Parameter Set #2 may be selected. (See examples in Figure 6.) The targeted vertebral body will be treated with the physician-selected ablation parameters:

**Ablation Parameter Set #1**

- Temperature: 75°C
- Temperature ramp: 0.5°C/second
- Ablation time: 7 minutes (420 seconds)
- Ablation zone:  $\sim 5$  mm radially

**Ablation Parameter Set #2**

- Temperature: 85°C
- Temperature ramp: 1°C/second
- Ablation time: 15 minutes (900 seconds)
- Ablation zone:  $\sim 6$  mm radially

**Note:** The RFG includes feedback controls to automatically control the ablation therapy via the Probe. If specific operating parameters fall outside of pre-determined ranges, the RFG may automatically suspend therapy. If this occurs, the event should be evaluated. In some cases, the Probe may be repositioned, or removed from the patient and wiped down with gauze and sterile saline before being placed back at the treatment site. If the cause of the event, or solution consistent with these instructions, cannot be determined, therapy should be aborted, and the Probe and Instruments should be removed from the patient, followed by appropriate wound care and post-operative patient management according to generally accepted standards of care. If the RFG does not reach the set temperature, then therapy should be manually aborted.

---

**Precaution:** Do not remove or withdraw the Probe while energy is being delivered, due to risk of burn, injury, or incomplete tissue ablation.

---

21. Upon completion of treatment, carefully remove the Probe from the Curved Cannula. Retract the Curved Cannula into the Introducer Cannula by rotating the Gear Wheel clockwise. When the Gear Wheel is at its lowest threaded position and frictionally locked preventing further rotation, the curved portion of the Curved Cannula will be inside the Introducer Cannula lumen. Remove the Introducer Cannula together with the Curved Cannula from the patient using a gauze pad at the skin incision site to prevent blood/tissue displacement. Separate the Introducer Cannula from the Curved Cannula; set aside the Introducer Cannula within the sterile field in case it's needed again. Dispose of the Curved Cannula per standard practice; DO NOT REUSE the Curved Cannula for treatment of additional vertebral bodies.

**Note:** Care should be taken when retracting the Curved Cannula to prevent the Curved Cannula from kinking against the Introducer Cannula. If resistance is felt with Gear Wheel rotation, it is recommended to halt further retraction and remove the Introducer Cannula and Curved Cannula as a single unit.

22. If treating an additional vertebral body, wipe down the Probe with gauze and sterile saline (to prevent blood from drying on the electrodes) and set within the sterile field before proceeding to step 23 below. If procedure has been completed, proceed to step 26.

**Note:** The Probe should remain plugged into the RFG until all vertebral bodies have been treated.

To treat an additional vertebral body:

23. Assemble an Introducer by using an Introducer Cannula from the Incept Access Instruments kit, with either the Bevel Stylet or Diamond Stylet. Orient the alignment tab of the Stylet to the slot of the Introducer Cannula and insert until it securely clicks into position. The Stylet is now locked to the Introducer Cannula.
24. Assemble a new Curved Cannula Assembly by using an unused Curved Cannula with the J-Stylet from step 12. Grip the shaft of the J-Stylet just proximal to the bend and carefully insert the distal tip into the slotted channel of the Curved Cannula handle until it is seated in the slot. Simultaneously raise the J-Stylet upright and feed the instrument into the Curved Cannula. Visually locate the alignment tab underneath the J-Stylet handle and align it to the slot on the top of the Curved Cannula handle. Complete insertion until the J-Stylet handle touches the Curved Cannula handle.
25. Repeat steps 5 to 21 to create a path to the treatment location within the next targeted vertebral body to allow for basivertebral nerve ablation.
26. Upon completion of the procedure, disconnect the Probe from the RFG by pulling on the connector body to release the captive latch. Do not pull directly on the wire, as this can pull the generator from the platform.
27. Dispose of product in accordance with accepted medical practice and applicable local and national regulations.
28. After skin closure using standard operating procedures, patient should receive appropriate wound care and post-operative patient management according to generally accepted standards of care.

**Note:** The Incept Access Instruments and Incept RF Probe are not implantable. Patients treated shall not be contraindicated from future imaging using traditional diagnostic imaging modalities such as CT, ultrasound, and MRI.

## DISPOSAL

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

This device contains a sharps hazard. Take precautions to ensure that sharps are handled properly. Dispose of all sharps directly into a sharps disposal container labeled with a biological hazard symbol. Sharps waste should be safely disposed of using available sharps waste channels in accordance with hospital, administrative, and/or local government policy.

## POST PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.


**LIMITATION OF LIABILITY**

Relievent Medsystems, Inc. ("Relievent") warrants that reasonable care has been used in the manufacture of this device. There is no express or implied warranty, including fitness for a particular purpose, on the Relievent Intercept Intraosseous Nerve Ablation System. The description or specifications are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Relievent Medsystems, Inc. is not responsible for any direct, incidental, special, or consequential loss, damage, or expense based on any defect, failure, or malfunction of this product, or arising from or related to its use, other than as expressly provided by mandatory provisions of applicable law. No person has the authority to bind Relievent Medsystems, Inc. to any representation or warranty except as indicated in this Limited Warranty.


**DISCLAIMER OF LIABILITY FOR DAMAGES**

In any claim or lawsuit for damages arising from alleged breach of warranty, breach of contract, negligence, product liability or any other legal or equitable theory, the buyer specifically agrees that neither Distributor nor Manufacturer, nor any agent thereof, shall be liable for damages for loss of profits or claims of buyer's customers for any such damages. The sole liability of Distributor or Manufacturer, or any agent thereof, for damages shall be limited to the cost to buyer of the specified goods sold by Distributor to buyer which give rise to the claim for liability.


**FIGURES**

**FIGURE 1**

**FIGURE 2**

**FIGURE 3**


FIGURE 4


FIGURE 5


FIGURE 6

SYMBOLS USED ON LABELING

| Symbol                                                                             | Definition                                                      |
|------------------------------------------------------------------------------------|-----------------------------------------------------------------|
|    | Sterilized using irradiation                                    |
|    | Single sterile barrier system                                   |
|    | Single sterile barrier system with protective packaging outside |
|    | Do not re-use                                                   |
|    | Do not re-sterilize                                             |
|    | Catalogue number                                                |
|    | Batch code                                                      |
|    | Serial number                                                   |
|     | Use by date                                                     |
|    | Consult instructions for use                                    |
|    | Unique device identifier                                        |
|    | Medical Device under EU Legislation                             |
|    | Contents                                                        |
|    | Manufacturer                                                    |
|   | Date of Manufacture                                             |
|  | Global Trade Item Number                                        |
|  | Contains hazardous substances                                   |
|  | Do not use if package is damaged                                |
|  | Keep dry                                                        |
|   | Upper limit of temperature                                      |
|  | Authorized representative in the European Community             |

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Intracept is a trademark of Relievable Medsystems Inc.
© 2025 Relievable Medsystems, Inc. All rights reserved.
Relievable Medsystems, Inc.
Boston Scientific Limited