Vercise™ Adapter S8 Instructions for Use

Instructions for Use

Mode d'emploi

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Kasutusjuhised

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How to Use this Manual

This manual describes the usage of the Vercise™ Adapter S8. This product is a component of the Boston Scientific Deep Brain Stimulation (DBS) System. Read all instructions carefully before using the DBS System.

For indications for use, contraindications, warnings, precautions, adverse events, storage and handling, sterilization, device disposal, and post-procedure information, refer to the *DBS Information for Prescribers IFU*. For other device-specific information not included in this manual, refer to the appropriate Instructions for Use (IFU) for your Boston Scientific DBS System as listed in your *DBS Reference Guide*. Advise the patient that additional information may be available to them on the Boston Scientific website www.bostonscientific.com/patientlabeling.

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

Trademarks

Vercise™ is a trademark of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

Warranty

For device warranty information, visit www.bostonscientific.com/warranty.

Summary of Safety and Clinical Performance

For customers in the European Union: visit <https://ec.europa.eu/tools/eudamed> to view a summary of this device's safety and clinical performance information on the EUDAMED website.

Technical Support

There are no user serviceable parts. If you have a specific question or issue, contact your sales representative. To contact Boston Scientific for any other reason, use the contact information provided for your locality via www.bostonscientific.com.

Product Model Numbers

Model Number	Description
DB-9208-15	Vercise™ Adapter S8, 8 Contact, 15cm
DB-9208-55	Vercise™ Adapter S8, 8 Contact, 55cm

Table of Contents

Device Description	1
Package Contents	1
Specifications and Technical Data	1
Materials List	2
Device Compatibility	2
Safety Information	3
Intended Use	3
Warnings	3
Cautions	3
Connecting the Adapter to the IPG and Lead Extensions	4
Explanting or Replacing the Adapter	6

This page intentionally left blank.

Device Description

The Adapters are thin, insulated wires that connect Abbott (St. Jude Medical) Lead Extensions to the Stimulator. The Adapters transfer the electrical stimulation from the Stimulator to the Lead.

Package Contents

The package contents of the Vercise Adapter S8 are provided in Table 1.

Table 1: DB-9208-15 and DB-9208-55 Adapter Kits	
Description	Quantity
Vercise Adapter S8	1
Hex Wrench	1
Note: All contents of the inner package or tray are sterile.	

Specifications and Technical Data

The physical characteristics of the Vercise Adapter S8 are provided in Table 2.

Table 2: Technical Specifications Adapter S8 (DB-9208-15 and DB-9208-55)	
Feature	Specification
Adapter Length	15 cm or 55 cm
Adapter Diameter	1.3 mm
Number of Contacts	8
Contact Material	Platinum/Iridium
Insulation Material	Polyurethane, Silicone
Conductor Material	MP35N-DFT-28% Ag

Materials List

Materials and substances that may contact the body during the DBS implant procedure are provided in Table 3.

Table 3: Materials	
Materials	Percent Weight
Silicone	37% to 74%
Thermoplastic	26% to 63%
Epoxy	< 1%

Device Compatibility

The Vercise Adapter S8 is compatible with the following Abbott (St. Jude Medical) 8-Channel Lead Extensions:

- Model 6371
- Model 6372
- Model 6373
- Model 6377
- Model 6378
- Model 6379

These Abbott (St. Jude Medical) Lead Extensions should be connected to the following Abbott (St. Jude Medical) Leads for programming:

- Model 6170
- Model 6171
- Model 6172
- Model 6173
- Model 6178
- Model 6179
- Model 6180
- Model 6181

Safety Information

Intended Use

The Vercise™ Adapter S8 is intended to connect compatible Abbott (St. Jude Medical) Lead Extensions to the compatible Boston Scientific Implantable Pulse Generator (IPG).

Warnings

Use of the Vercise Adapter S8 as part of the DBS System has not been tested for compatibility with MRI. Patients with the Vercise Adapter S8 should not be subjected to an MRI.

Cautions

Do not bend, kink, or stretch the Lead Extension or Vercise Adapter S8; this may break its wires or connections. Broken wires or connections may create an open circuit, resulting in component failure, loss of stimulation, and potential surgical intervention.

Do not handle the Vercise Adapter S8 with instruments that were not supplied with the Adapter. Other instruments may cause damage to the device, including broken wires or connections.

The Setscrew in the Vercise Adapter S8 must be retracted to insert the Lead Extension. There is no stop feature to prevent excessive retraction of the Setscrew and disassembly of the Adapter. To avoid damage to the sealing grommet, do not retract the Setscrew assembly more than is needed to insert the Lead Extension.

Minimize the number of Hex Wrench insertions into the Vercise Adapter S8 to prevent damage to the septum seal plug.

Connecting the Adapter to the IPG and Lead Extensions

1. Refer to the appropriate Abbott (St. Jude Medical) product labeling for instructions on explanting the Abbott (St. Jude Medical) IPG. Explant the Abbott (St. Jude Medical) IPG.
2. Carefully open the Vercise Adapter S8 package and transfer the contents of the inner tray to the sterile field.
3. Ensure that the Abbott (St. Jude Medical) Lead Extension is clean and not damaged.
4. Retract the Setscrews in all Adapters. To do so, insert the Hex Wrench into the septum containing the Setscrew and turn the Hex Wrench counterclockwise (Figure 1).


Figure 1. The Setscrew Assembly in the Vercise Adapter S8

5. Ensure that stimulation is OFF before connecting to the Vercise Adapter S8.
6. Ensure that the Lead Extension is clean and dry before inserting into the Vercise Adapter S8.
7. Gently guide the proximal end of the Abbott (St. Jude Medical) Lead Extension into the Vercise Adapter S8 until the Lead Extension reaches the end stop of the Adapter. Use one Vercise Adapter S8 for each Abbott (St. Jude Medical) Lead Extension (Figure 2).

Caution: Take care not to bend or kink the Abbott (St. Jude Medical) Lead Extension when inserting into the Vercise Adapter S8. Take care not to over insert the Abbott (St. Jude Medical) Lead Extension into the Vercise Adapter S8.


Figure 2. Vercise Adapters S8 Connected to Abbott (St. Jude Medical) Lead Extensions and a Boston Scientific DBS IPG

8. Ensure that the proximal end of the Vercise Adapter S8 is clean and insert the proximal end of the Adapter into the Boston Scientific DBS Stimulator.
9. Perform the Stimulator implantation procedure. For instructions on implanting the Stimulator, refer to the *Surgical Implant Manual* for your Boston Scientific DBS System as listed in your *DBS Reference Guide*. Complete Step a prior to closing the incisions.
 - a. Tighten the Setscrew in the Vercise Adapter S8 until the Hex Wrench clicks, indicating the Setscrew is fully secured.

Note: Do not place the distal segment of the Vercise Adapter S8 between its proximal end and the Boston Scientific DBS Stimulator. Placement of wires between these two devices could result in damage to the wire insulation resulting in loss of stimulation.

Note: Do not wrap or coil the excess Lead Extension wires or the distal Adapter segment around the perimeter or in front of the Stimulator. Wrapping around the Stimulator may cause wires to slip in front of the Stimulator. Placing wires in front of the Stimulator may cause kinking of the Lead Extension and Adapter wires, potential interference with telemetry and/or recharge operation, and damage during future Stimulator replacement surgery.

Note: If the Vercise Adapter S8 is placed or moves superficial to the Stimulator, it could cause communication or charging difficulties.

10. To program Leads, refer to the *Programming Manual* for your Boston Scientific DBS System, as listed in your *DBS Reference Guide*.
11. Place the device peel-off stickers, where necessary.

Explanting or Replacing the Adapter

To disconnect and remove the Vercise Adapter S8 from the IPG or Lead Extension, follow these steps:

1. Turn off the IPG.
2. Palpate the implanted area to locate the IPG.
3. Surgically open the pocket where the IPG is located. The incision should be large enough to remove the IPG from the pocket. Be careful not to damage the implanted components to allow for proper analysis following explant.
4. Remove the IPG from the IPG pocket.
5. Use the Hex Wrench to retract the IPG Setscrew connecting the Vercise Adapter S8 to the IPG header without fully removing the Setscrew.
6. Gently pull the Vercise Adapter S8 out of the IPG header.
7. Use the Hex Wrench to retract the Setscrew on the Vercise Adapter S8 without fully removing the Setscrew.
8. Gently pull the Lead Extension out of the Vercise Adapter S8.
9. To replace the Vercise Adapter S8, see the “Connecting the Adapter to the IPG and Lead Extensions” section of this manual and repeat all steps.
10. Ship the explanted Vercise Adapter S8 to Boston Scientific.

This page intentionally left blank.

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or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
Neuromodulation Corporation
© 2024 Boston Scientific Corporation