



Refer to the device directions for use for complete instructions on device use.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Intended Use

The Pneumatic Inflator is intended for inflation and deflation of the Rigiflex II Single Use Achalasia Balloon Dilator while simultaneously monitoring the balloon pressure.

Indications For Use

The Pneumatic Inflator is only to be used with the Rigiflex II Single Use Achalasia Balloon Dilator for dilatation of the cardia (LES) in patients with achalasia.

Contraindications

None known.

Warnings

This product contains natural rubber latex which may cause allergic reactions.