

Extractor™ Pro RX-S

Retrieval Balloon Catheter

Boston
Scientific

Advancing science for life™




RX Biliary System™
We're in!™

Extractor™ Pro RX-S

Retrieval Balloon Catheter

The Extractor Pro RX-S Retrieval Balloon Catheter is indicated for use endoscopically to remove stones from the biliary system or to facilitate the injection of contrast medium while occluding the duct with the balloon.






Enhanced Design

- Increased catheter stiffness designed for improved pushability and kink-resistance.
- Each of the three balloon sizes can be inflated to two distinct diameters: 9–12mm; 12–15mm; 15–18mm to better accommodate anatomical variations.
- **Squared-shoulder balloon** is designed to facilitate stone removal.

Procedural Efficiency and Control

- The RX Biliary System locking device allows physicians to maintain guidewire position and helps facilitate multiple sweeps of the duct during stone retrieval.
- Option to front load or back load the guidewire.
- Injection-above or injection-below options to support physician preference and facilitate procedural needs.
- Designed for excellent injection flow rate.



Ordering Information: Extractor Pro RX-S Retrieval Balloon Catheter

Order Number	Balloon O.D. (mm)	Injection
M00547300	9-12mm	Above
M00547330	9-12mm	Below
M00547310	12-15mm	Above
M00547340	12-15mm	Below
M00547320	15-18mm	Above
M00547350	15-18mm	Below

All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: The law restricts this device to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

RX Biliary System™

We're in!™

The RX Biliary System helps put the physician in control of the guidewire. Using the RX Locking Device, the physician can secure the guidewire to maintain access while exchanging devices. The RX Biliary System and its compatible devices offer a broad portfolio of products for pancreatico-biliary stone and stricture management.

Boston Scientific
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com/gastro
www.EndoSuite.com

©2016 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-396406-AA June 2016