



# Enteral Stents

## WALLSTENT™ Colonic and Duodenal Stents (Continued)



 Magnetic Resonance Conditional<sup>1</sup>

### WALLSTENT Single-Use Colonic & Duodenal Endoprosthesis with UniStep™ Plus Delivery System\*



Order Number	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565560	20	60	230	255	10	Each
M00565570	20	90	230	255	10	Each
M00565580	22	60	230	255	10	Each
M00565590	22	90	230	255	10	Each
M00565620	20	60	135	160	10	Each
M00565630	20	90	135	160	10	Each
M00565640	22	60	135	160	10	Each
M00565650	22	90	135	160	10	Each


\*Recommended Guidewire .035" WALLSTENT Super Stiff Guidewire, Order #18001, (see below).

<sup>1</sup>Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

### WALLSTENT Super Stiff Guidewire



Order Number	Description	Distal Tip O.D. (in)	Length (cm)	Tip Style	Packaging
H965180011	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5



 Magnetic Resonance Conditional<sup>1</sup>

### Ultraflex™ Precision Single-Use Colonic Stent System\*



Order Number	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery Tip Diameter (Fr)	Packaging
M00557330	25	30	57	100	105	22	Each
M00557360	25	30	87	100	105	22	Each
M00557380	25	30	117	100	105	22	Each

\*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, (see page 4).

<sup>1</sup>Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.