



**INTERLOCK™** Coils

# MORE FIBER. MORE POWER.



Interlock Coils  
**PERFORM  
NEARLY  
3X FASTER\***



Interlock Coils feature dense networks of fibers to facilitate shutting down vessels faster, with fewer coils.

\* Bench testing performed by independent laboratory. Data on file. n=7. Bench test results may not necessarily be indicative of clinical performance.

# Interlock™ -18 Fibered IDC™ Occlusion System

## Standard Length 2D Configurations



UPN	Description	Diameter (mm)	Length (cm)	Shape
M001361480	Interlock - 18 Coil	2	4	2D
M001361490	Interlock - 18 Coil	2	6	2D
M001361500	Interlock - 18 Coil	3	6	2D
M001361510	Interlock - 18 Coil	3	12	2D
M001361520	Interlock - 18 Coil	4	8	2D
M001361530	Interlock - 18 Coil	4	15	2D
M001361540	Interlock - 18 Coil	5	8	2D
M001361550	Interlock - 18 Coil	5	15	2D
M001361560	Interlock - 18 Coil	6	10	2D
M001361570	Interlock - 18 Coil	6	20	2D
M001361580	Interlock - 18 Coil	8	20	2D
M001361590	Interlock - 18 Coil	10	20	2D
M001361600	Interlock - 18 Coil	10	30	2D
M001361610	Interlock - 18 Coil	12	20	2D
M001361620	Interlock - 18 Coil	12	30	2D
M001361630	Interlock - 18 Coil	14	20	2D
M001361640	Interlock - 18 Coil	14	30	2D

## Long Length 2D Configurations



UPN	Description	Diameter (mm)	Length (cm)	Shape
M001361920	Interlock - 18 Coil	10	50	2D
M001361930	Interlock - 18 Coil	14	50	2D
M001361940	Interlock - 18 Coil	18	50	2D
M001361950	Interlock - 18 Coil	20	50	2D
M001361960	Interlock - 18 Coil	22	60	2D

## Diamond Configurations



UPN	Description	Diameter (mm)	Length (cm)	Shape
M001361740	Interlock - 18 Coil	2/3	2.3	Diamond
M001361750	Interlock - 18 Coil	2/4	4.1	Diamond
M001361760	Interlock - 18 Coil	2/5	5.8	Diamond
M001361770	Interlock - 18 Coil	2/6	8	Diamond

**Boston Scientific**  
Advancing science for life™

### INTERLOCK™ FIBERED IDC™ OCCLUSION SYSTEM

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions. **INTENDED USE/INDICATIONS FOR USE:** The Interlock Fibered IDC Occlusion System is a modified interlocking detachable coil indicated to obstruct or reduce rate of blood flow in the peripheral vasculature. This device is not intended for neurovascular use. **CONTRAINDICATIONS:** None known. **GENERAL PRECAUTIONS:** Do not advance the Interlock Fibered IDC Occlusion System if it becomes lodged within the microcatheter. Determine the cause of the resistance and replace the microcatheter and coil if necessary.

**ADVERSE EVENTS:** The complications that may result from a peripheral embolization procedure include, but are not limited to: • Complications related to catheterization (e.g., hematoma at the site of entry, vessel injury, etc.) • Death • Emboli • Foreign body reactions necessitating medical intervention • Hemorrhage • Infection necessitating medical intervention • Ischemia • Pain • Recanalization • Temporary neurological deficit • Tissue necrosis • Undesirable clot formation of the vasculature • Vasospasm (REV AA)

Interlock and Fibered IDC are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

### Peripheral Interventions

300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2014 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-276902-AA DEC2014