

Endoscopy Product Catalogue

**Boston
Scientific**
Advancing science for life™

Produktkatalog – Catalogue Produits
Catálogo de Productos – Catalogo dei Prodotti

2017



Updated 01/2017



Contents



Alphabetical
Index

Product Code
Index

Contents

RX Biliary System™

Cannulas	Tandem™ RX Cannula	2
	RX Cannulas	2
Sphincterotomes	Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	3
	Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	3
	Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	4
	Autotome™ RX Cannulating Sphincterotomes	5
	RX Needle Knife XL Triple Lumen Needle Knife	5
Stone Retrieval	Extractor™ Pro RX Triple-lumen Retrieval Balloons – Distal Injection	6
	Extractor™ Pro RX Triple-lumen Retrieval Balloons – Proximal Injection	6
	Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection	7
	Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection	7
	Trapezoid™ RX Wire-guided Retrieval Baskets	7
Dilation	Hurricane™ RX Biliary Dilation Balloons	8
Cytology	RX Cytology Brush	8
Accessories	RX Accessories	9
	Encore™ 26 Inflation Device	9

Biliary

Cannulas	Tandem™ XL Triple-lumen ERCP Cannula	11
	Contour™ ERCP Cannulas	11
	Fluoro Tip™ ERCP Cannulas	12
Sphincterotomes	Stonetome™ Stone Removal Device	12
	TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	13
	TRUEtome™ Dreamwire Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	13
	TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	13
	TRUEtome™ Cannulating Sphincterotomes	14
	Ultratome XL™ Triple-lumen Sphincterotomes	15
	Ultratome™ Double-lumen Sphincterotomes	15
	Microknife™ XL Triple-lumen Needleknife	15
Stone Retrieval	Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection	16
	Extractor™ Pro XL Retrieval Balloon Catheter – Proximal Injection	16
	Extractor™ Pro DL Retrieval Balloon Catheter	16
Dilation	Max Force™ Biliary Dilation Balloons	17
Accessories	Biliary Accessories	17

SpyGlass™ Direct Visualisation Systems

SpyGlass™ DS	SpyGlass™ DS Capital Equipment	19
	SpyGlass™ DS Disposable Equipment	19
SpyGlass™	SpyGlass™ Capital Equipment	20
	SpyGlass™ Reusable Equipment	20
	SpyGlass™ Disposable Equipment	20
	SpyGlass™ Optional Accessories	20



Contents (CONTINUED)

Endoscopic Ultrasound (EUS)

Fine Needle Aspiration (FNA)	Expect™ Standard Handle Needles	22
	Expect™ Slimline Handle Needles	22
Fine Needle Biopsy (FNB)	Acquire™ EUS-FNB Device	23
Stents and Delivery Systems	Hot AXIOS™ Stents and Electrocautery Enhanced Delivery Systems	23
	AXIOS™ Stents and Delivery Systems	23

Stents

RX Biliary Metal Stents	WallFlex™ Biliary RX Fully Covered Stent System RMV with Permalume™ Covering (including benign indication)	25
	WallFlex™ Biliary RX Partially Covered Stents with Permalume™ Covering	25
	WallFlex™ Biliary RX Non-covered Stents	26
	Wallstent™ RX Biliary Stents with Permalume™ Covering	27
	Wallstent™ RX Biliary Stents	27
Biliary Metal Stents	Wallstent™ Biliary Stents with Permalume™ Covering	28
	Wallstent™ Biliary Stents (Uncovered)	28
RX Biliary Plastic Stents	Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	29
	Advanix™ Stents – Singles	30–31
	NaviFlex™ RX Delivery System	31
	RX Biliary Plastic Stents	32
Biliary Plastic Stents	Flexima™ Biliary Stents	33
	Flexima™ Nasobiliary Stents	33
	Percuflex™ Amsterdam Biliary Stents	34
	Percuflex™ Amsterdam Introducer Kits	34
	C-flex™ Pigtail Biliary Stents	34
Pancreatic Plastic Stents	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail	35
	Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight	35
	Advanix™ Pancreatic Single Stents – Pigtail	36–37
	Advanix™ Pancreatic Single Stents – Straight	37–39
	NaviFlex™ RX Pancreatic Delivery System & NaviFlex Pushers	40
Esophageal Stents	WallFlex™ Fully Covered RMV Esophageal Stents with Permalume™ Covering (including refractory benign indication)	41
	WallFlex™ Fully Covered RMV Esophageal Stents Longer Loop with Permalume™ Covering (including refractory benign indication)	41
	WallFlex™ Partially Covered Esophageal Stents with Permalume™ Covering	41
	Ultraflex™ Covered Esophageal NG Stents – Proximal Release	42
	Ultraflex™ Covered Esophageal NG Stents – Distal Release	42
	Ultraflex™ Non-covered Esophageal NG Stents – Proximal Release	42
	Ultraflex™ Non-covered Esophageal NG Stents – Distal Release	42
	Polyflex Stent™ Self-expanding Silicone Esophageal Stents	43
Colonic and Duodenal Stents	WallFlex™ Colonic Stents	44
	WallFlex™ Duodenal Stents	44
	Wallstent™ Colonic and Duodenal Stents	45
	Ultraflex Precision™ Colonic Stents	45



Contents (CONTINUED)

Dilation

CRE™ Controlled Radial Expansion Balloon Dilators	CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
	CRE™ PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
	CRE™ Fixed Wire Balloon Dilators	48
Max Force™ TTS Balloon Dilators	Max Force™ TTS Esophageal Balloon Dilators	49
	Rigiflex™ II Achalasia Balloon Dilators	49
Inflation Systems	Alliance™ II Inflation System	50
Dilation Accessories	Inflation/Deflation Devices and Accessories	50

Guidewires

High Performance	Dreamwire™ High Performance Guidewires	52
	Hydra Jagwire™ High Performance Guidewires	52
	Jagwire™ High Performance Guidewires	53
	Jagwire™ High Performance Guidewires – Extendable Jagwire	53
	NovaGold™ High Performance Guidewires	54
	NaviPro™ 100% Hydrophilic Guidewires	54
Standard	Amplatz™ Super Stiff Guidewire	54
	Wallstent™ Super Stiff Guidewire	54

Hemostasis

Clipping Devices	Resolution 360™ Clip	56
	Resolution™ Clip	56
Multiple Band Ligators	Speedband Superview Super 7™ Multiple Band Ligators	56
Bipolar Electrohemostasis Catheters	Injection Gold Probe™ Electrohemostasis/Injection Catheters	57
	Gold Probe™ Electrohemostasis Catheters	57
Needles	Interject™ Sclerotherapy Needles	58
	Interject™ Clear Sclerotherapy Needles	58
	Optiflo™ Hemostasis Catheters	59
Electrosurgical Bipolar Cable Adaptors	Electrosurgical Bipolar Cable Adaptors	59

Tissue Resection

EMR Devices for Upper GI	Captivator™ EMR	61
Snares	Rotatable Snares	61
	Sensation™ Short Throw Snares	62
	Captiflex™ Snares	62
	Captivator™ Snares	63
	Captivator™ II Snares	63
	Profile™ Pediatric Snares	64
Hot Biopsy Forceps	Radial Jaw™ 4 Hot Biopsy Forceps	64
	Active Cords	64



Contents (CONTINUED)

Biopsy

Forceps	Radial Jaw™ 4 Standard Capacity Biopsy Forceps	66
	Radial Jaw™ 4 Large Capacity Biopsy Forceps	66
	Radial Jaw™ 4 Jumbo Biopsy Forceps	67
	Radial Jaw™ 4 Gastropediatric Biopsy Forceps	67
	Multibite™ Multiple Sample Biopsy Forceps	67

EndoVive™ Enteral Access Devices

Initial Placement	Standard PEGs	69
	Standard PEGs with ENFit	69
	PEG Accessories	69
	PEG Accessories with ENFit	69
	One Step Button Low Profile PEG Kits	70
	One Step Button and Low Profile Button Accessories	71
	Through-The-PEG (TTP) Jejunal Feeding Tube	71
Replacements	Standard Balloon Replacements	72
	Replacement Gastrostomy Feeding Tube with ENFit	73
	Low Profile Button Replacements	73

Pulmonary

Alair™ Bronchial Thermoplasty (BT) System	Alair™ Bronchial Thermoplasty (BT) System	75
	Alair™ BT Controller Accessory Kits ATS 201	75
Endoscopic Bronchial Ultrasound	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle – Olympus Scope Compatible	76
Airway Stents	Dynamic™ (Y) Stents	77
	Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents	77
	Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents	77
	Ultraflex™ Non-covered Tracheobronchial NG Stents – Proximal Release – Tracheal Stents	78
	Ultraflex™ Non-covered Tracheobronchial NG Stents – Proximal Release – Bronchial Stents	78
	Ultraflex™ Non-covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents	78
	Ultraflex™ Non-covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents	78
	Polyflex Stent™ Self-expanding Silicone Airway Stents	79
Biopsy Forceps	Radial Jaw™ 4 Pulmonary Biopsy Forceps	80
Transbronchial Aspiration Needles	eXcelon™ Transbronchial Aspiration Needle with Syringe	80
Cytology Brushes	Cellebri™ Endoscopic Cytology Brushes – Bronchoscope	80
Airway Retrieval Baskets	ZeroTip™ Airway Retrieval Baskets	81
Dilation Balloons	CRE™ Pulmonary Dilation Balloons	81
Guidewires	Jagwire™ Pulmonary Guidewire	81
Inflation System	Alliance™ II Inflation System	81

Alphabetical Index

Acquire™ EUS-FNB Device	23
Active Cords	17, 64
Advanix™ Pancreatic Single Stents – Pigtail	36–37
Advanix™ Pancreatic Single Stents – Straight	37–39
Advanix™ Pancreatic Stent Kits (Stent+ RX Pusher) – Pigtail	35
Advanix™ Pancreatic Stent Kits (Stent+ RX Pusher) – Straight	35
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System	29
Advanix™ Stents – Singles	30–31
Alair™ Bronchial Thermoplasty (BT) System	75
Alair™ BT Controller Accessory Kits ATS 201	75
Alliance™ II Inflation System	50, 81
Amplatz™ Super Stiff Guidewire	54
Autotome™ RX Cannulating Sphincterotomes	5
AXIOS™ Stents and Delivery Systems	23
Biliary Accessories	17
C-flex™ Pigtail Biliary Stents	34
Captiflex™ Snares	62
Captivator™ EMR	61
Captivator™ II Snares	63
Captivator™ Snares	63
Celebrity™ Endoscopic Cytology Brushes – Bronchoscope	80
Contour™ ERCP Cannulas	11
CRE™ Fixed Wire Balloon Dilators	48
CRE™ PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
CRE™ Pulmonary Dilation Balloons	81
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	3
Dreamwire™ High Performance Guidewires	52
Dynamic™ (Y) Stents	77
Electrosurgical Bipolar Cable Adaptors	59
Encore™ 26 Inflation Device	9, 17
eXcel™ Transbronchial Aspiration Needle with Syringe	80
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle – Olympus Scope Compatible	76
Expect™ Slimline Handle Needles	22
Expect™ Standard Handle Needles	22
Extractor™ Pro DL Retrieval Balloon Catheter	16
Extractor™ Pro RX Triple-lumen Retrieval Balloons – Distal Injection	6
Extractor™ Pro RX Triple-lumen Retrieval Balloons – Proximal Injection	6
Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection	7
Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection	7
Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection	16
Extractor™ Pro XL Retrieval Balloon Catheter – Proximal Injection	16
Flexima™ Biliary Stents	33
Flexima™ Nasobiliary Stents	33
Fluoro Tip™ ERCP Cannulas	12
Gold Probe™ Electrohemostasis Catheters	57
Hot AXIOS™ Stents and Electrocautery Enhanced Delivery Systems	23
Hurricane™ RX Biliary Dilation Balloons	8
Hydra Jagwire™ High Performance Guidewires	52
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	3



Alphabetical Index (CONTINUED)

Inflation/Deflation Devices and Accessories	50
Injection Gold Probe™ Electrohemostasis/Injection Catheters	57
Interject™ Clear Sclerotherapy Needles	58
Interject™ Sclerotherapy Needles	58
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	4
Jagwire™ High Performance Guidewires	53
Jagwire™ High Performance Guidewires – Extendable Jagwire	53
Jagwire™ Pulmonary Guidewire	81
Low Profile Button Replacements	73
Max Force™ Biliary Dilation Balloons	17
Max Force™ TTS Esophageal Balloon Dilators	49
Microknife™ XL Triple-lumen Needleknife	15
Multibite™ Multiple Sample Biopsy Forceps	67
NaviFlex™ RX Delivery System	31
NaviFlex™ RX Pancreatic Delivery System & NaviFlex Pushers	40
NaviPro™ 100% Hydrophilic Guidewires	54
NovaGold™ High Performance Guidewires	54
One Step Button and Low Profile Button Accessories	71
One Step Button Low Profile PEG Kits	70
Optiflo™ Hemostasis Catheters	59
PEG Accessories	69
PEG Accessories with ENFit	69
Percuflex™ Amsterdam Biliary Stents	34
Percuflex™ Amsterdam Introducer Kits	34
Polyflex Stent™ Self-expanding Silicone Airway Stents	79
Polyflex Stent™ Self-expanding Silicone Esophageal Stents	43
Profile™ Pediatric Snares	64
Radial Jaw™ 4 Gastropediatric Biopsy Forceps	67
Radial Jaw™ 4 Hot Biopsy Forceps	64
Radial Jaw™ 4 Jumbo Biopsy Forceps	67
Radial Jaw™ 4 Large Capacity Biopsy Forceps	66
Radial Jaw™ 4 Pulmonary Biopsy Forceps	80
Radial Jaw™ 4 Standard Capacity Biopsy Forceps	66
Replacement Gastrostomy Feeding Tube with ENFit	73
Resolution 360™ Clip	56
Resolution™ Clip	56
Rigiflex™ II Achalasia Balloon Dilators	49
Rotatable Snares	61
RX Accessories	9
RX Biliary Plastic Stents	32
RX Cannulas	2
RX Cytology Brush	8
RX Needle Knife XL Triple Lumen Needle Knife	5
Sensation™ Short Throw Snares	62
Speedband Superview Super 7™ Multiple Band Ligators	56
SpyGlass™ Capital Equipment	20
SpyGlass™ Disposable Equipment	20
SpyGlass™ DS Capital Equipment	19
SpyGlass™ DS Disposable Equipment	19
SpyGlass™ Optional Accessories	20



Alphabetical Index (CONTINUED)

SpyGlass™ Reusable Equipment	20
Standard Balloon Replacements	72
Standard PEGs	69
Standard PEGs with ENFit	69
Stonetome™ Stone Removal Device	12
Tandem™ RX Cannula	2
Tandem™ XL Triple-lumen ERCP Cannula	11
Through-The-PEG (TTP) Jejunal Feeding Tube	71
Trapezoid™ RX Wire-guided Retrieval Baskets	7
TRUEtome™ Cannulating Sphincterotomes	14
TRUEtome™ Dreamwire Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire	13
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire	13
TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire	13
Ultraflex Precision™ Colonic Stents	45
Ultraflex™ Covered Esophageal NG Stents – Distal Release	42
Ultraflex™ Covered Esophageal NG Stents – Proximal Release	42
Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents	77
Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents	77
Ultraflex™ Non-covered Esophageal NG Stents – Distal Release	42
Ultraflex™ Non-covered Esophageal NG Stents – Proximal Release	42
Ultraflex™ Non-covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents	78
Ultraflex™ Non-covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents	78
Ultraflex™ Non-covered Tracheobronchial NG Stents – Proximal Release – Bronchial Stents	78
Ultraflex™ Non-covered Tracheobronchial NG Stents – Proximal Release – Tracheal Stents	78
Ultratome XL™ Triple-lumen Sphincterotomes	15
Ultratome™ Double-lumen Sphincterotomes	15
WallFlex™ Biliary RX Fully Covered Stent System RMV with Permalume™ Covering (including benign indication)	25
WallFlex™ Biliary RX Non-covered Stents	26
WallFlex™ Biliary RX Partially Covered Stents with Permalume™ Covering	25
WallFlex™ Colonic Stents	44
WallFlex™ Duodenal Stents	44
WallFlex™ Fully Covered RMV Esophageal Stents Longer Loop with Permalume™ Covering (including refractory benign indication)	41
WallFlex™ Fully Covered RMV Esophageal Stents with Permalume™ Covering (including refractory benign indication)	41
WallFlex™ Partially Covered Esophageal Stents with Permalume™ Covering	41
Wallstent™ Biliary Stents (Uncovered)	28
Wallstent™ Biliary Stents with Permalume™ Covering	28
Wallstent™ Colonic and Duodenal Stents	45
Wallstent™ RX Biliary Stents	27
Wallstent™ RX Biliary Stents with Permalume™ Covering	27
Wallstent™ Super Stiff Guidewire	54
ZeroTip™ Airway Retrieval Baskets	81



Product Code Index

Order Number	Description	Page Number
Wallstent™ Super Stiff Guidewire		
H965180011	Wallstent Super Stiff Guidewire	54
Wallstent™ Biliary Stents (Uncovered)		
H965430100	Wallstent Biliary	28
H965430200	Wallstent Biliary	28
H965430300	Wallstent Biliary	28
H965430400	Wallstent Biliary	28
H965430500	Wallstent Biliary	28
H965430600	Wallstent Biliary	28
Wallstent™ Biliary Partially Covered Stents with Permalume™ Covering		
H965430700	Wallstent Biliary w/Permalume	28
H965430800	Wallstent Biliary w/Permalume	28
H965430900	Wallstent Biliary w/Permalume	28
H965431000	Wallstent Biliary w/Permalume	28
H965431100	Wallstent Biliary w/Permalume	28
H965431200	Wallstent Biliary w/Permalume	28
Alair™ Bronchial Thermoplasty (BT) System		
M005ATS20000	Alair RF Controller – monopolar (return electrode not included)	75
Alair™ BT Controller Accessory Kits ATS 201		
M005ATS20120	British Type G	75
M005ATS20130	EU Hybrid Type E/F	75
M005ATS20140	Australian Type I	75
M005ATS20150	Danish Type K	75
M005ATS20160	Brazil IEC 60906-1	75
M005ATS20180	Swiss Type J	75
Alair™ Bronchial Thermoplasty (BT) System		
M005ATS25020	Alair BT Catheter – Global	75
PEG Accessories with ENFit		
M00501261	EndoVive Y-Port w/ENFit 20Fr	69
M00501271	EndoVive Y-Port w/ENFit 24Fr	69
Standard PEGs with ENFit		
M00509081	EndoVive StdPEG Pull w/ENFit 20Fr	69
M00509101	EndoVive StdPEG Pull w/ENFit 24Fr	69
M00509091	EndoVive StdPEG Push w/ENFit 20Fr	69
M00509111	EndoVive StdPEG Push w/ENFit 24Fr	69
Replacement Gastrostomy Feeding Tube with ENFit		
M00509881	12Fr (4.00mm) Repl. G-Tube, Straight	73
M00509891	14Fr (4.67mm) Repl. G-Tube, Straight	73
M00509901	16Fr (5.33mm) Repl. G-Tube, Straight	73
M00509911	18Fr (6.00mm) Repl. G-Tube, Straight	73
M00509921	20Fr (6.67mm) Repl. G-Tube, Straight	73
M00509931	22Fr (7.33mm) Repl. G-Tube, Straight	73
M00509941	24Fr (8.00mm) Repl. G-Tube, Straight	73
M00509961	12Fr (4.00mm) Repl. G-Tube, Right Angle	73



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Replacement Gastrostomy Feeding Tube with ENFit (CONTINUED)		
M00509971	14Fr (4.67mm) Repl. G-Tube, Right Angle	73
M00509981	16Fr (5.33mm) Repl. G-Tube, Right Angle	73
M00509991	18Fr (6.00mm) Repl. G-Tube, Right Angle	73
M00510001	20Fr (6.67mm) Repl. G-Tube, Right Angle	73
M00510011	22Fr (7.33mm) Repl. G-Tube, Right Angle	73
M00510021	24Fr (8.00mm) Repl. G-Tube, Right Angle	73
Multibite™ Multiple Sample Biopsy Forceps		
M00510101	Multibite Multiple Sample Biopsy Forceps	67
M00510121	Multibite Multiple Sample Biopsy Forceps	67
Trapezoid™ RX Wire-guided Retrieval Baskets		
M00510860	Trapezoid RX Wire-guided Retrieval Basket	7
M00510870	Trapezoid RX Wire-guided Retrieval Basket	7
M00510880	Trapezoid RX Wire-guided Retrieval Basket	7
M00510890	Trapezoid RX Wire-guided Retrieval Basket	7
Optiflo™ Hemostasis Catheters		
M00511731	Optiflo Hemostasis Catheter	59
M00511751	Optiflo Hemostasis Catheter	59
ZeroTip™ Airway Retrieval Baskets		
M00513200	ZeroTip Airway Retrieval Basket	81
M00513210	ZeroTip Airway Retrieval Basket	81
Radial Jaw™ 4 Large Capacity Biopsy Forceps		
M00513301	Radial Jaw 4 Large Capacity	66
M00513302	Radial Jaw 4 Large Capacity	66
M00513303	Radial Jaw 4 Large Capacity	66
M00513311	Radial Jaw 4 Large Capacity w/Needle	66
M00513312	Radial Jaw 4 Large Capacity w/Needle	66
M00513313	Radial Jaw 4 Large Capacity w/Needle	66
M00513321	Radial Jaw 4 Large Capacity	66
M00513322	Radial Jaw 4 Large Capacity	66
M00513323	Radial Jaw 4 Large Capacity	66
M00513331	Radial Jaw 4 Large Capacity w/Needle	66
M00513332	Radial Jaw 4 Large Capacity w/Needle	66
M00513333	Radial Jaw 4 Large Capacity w/Needle	66
Radial Jaw™ 4 Jumbo Biopsy Forceps		
M00513343	Radial Jaw 4 Jumbo	67
M00513353	Radial Jaw 4 Jumbo w/Needle	67
M00513361	Radial Jaw 4 Jumbo	67
M00513362	Radial Jaw 4 Jumbo	67
M00513363	Radial Jaw 4 Jumbo	67
M00513371	Radial Jaw 4 Jumbo w/Needle	67
M00513372	Radial Jaw 4 Jumbo w/Needle	67
M00513373	Radial Jaw 4 Jumbo w/Needle	67



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Radial Jaw™ 4 Standard Capacity Biopsy Forceps		
M00513381	Radial Jaw 4 Standard Capacity	66
M00513382	Radial Jaw 4 Standard Capacity	66
M00513383	Radial Jaw 4 Standard Capacity	66
M00513391	Radial Jaw 4 Standard Capacity w/Needle	66
M00513392	Radial Jaw 4 Standard Capacity w/Needle	66
M00513393	Radial Jaw 4 Standard Capacity w/Needle	66
M00513401	Radial Jaw 4 Standard Capacity	66
M00513402	Radial Jaw 4 Standard Capacity	66
M00513403	Radial Jaw 4 Standard Capacity	66
M00513411	Radial Jaw 4 Standard Capacity w/Needle	66
M00513412	Radial Jaw 4 Standard Capacity w/Needle	66
M00513413	Radial Jaw 4 Standard Capacity w/Needle	66
Radial Jaw™ 4 Gastropediatric Biopsy Forceps		
M00513441	Radial Jaw 4 Gastropediatric	67
M00513443	Radial Jaw 4 Gastropediatric	67
M00513451	Radial Jaw 4 Gastropediatric w/Needle	67
M00513453	Radial Jaw 4 Gastropediatric w/Needle	67
Ultraflex™ Non-covered Esophageal NG Stents – Distal Release		
M00513700	Ultraflex Esophageal NG Stent	42
M00513710	Ultraflex Esophageal NG Stent	42
M00513720	Ultraflex Esophageal NG Stent	42
Ultraflex™ Covered Esophageal NG Stents – Distal Release		
M00513730	Ultraflex Esophageal NG Stent	42
M00513740	Ultraflex Esophageal NG Stent	42
M00513750	Ultraflex Esophageal NG Stent	42
Ultraflex™ Non-covered Esophageal NG Stents – Proximal Release		
M00513800	Ultraflex Esophageal NG Stent	42
M00513810	Ultraflex Esophageal NG Stent	42
M00513820	Ultraflex Esophageal NG Stent	42
M00513830	Ultraflex Esophageal NG Stent	42
Ultraflex™ Covered Esophageal NG Stents – Proximal Release		
M00513840	Ultraflex Esophageal NG Stent	42
M00513850	Ultraflex Esophageal NG Stent	42
M00513860	Ultraflex Esophageal NG Stent	42
Ultraflex™ Covered Esophageal NG Stents – Distal Release		
M00514200	Ultraflex Esophageal NG Stent	42
M00514210	Ultraflex Esophageal NG Stent	42
Ultraflex™ Non-covered Esophageal NG Stents – Proximal Release		
M00514220	Ultraflex Esophageal NG Stent	42
M00514230	Ultraflex Esophageal NG Stent	42
Ultraflex™ Covered Esophageal NG Stents – Proximal Release		
M00514240	Ultraflex Esophageal NG Stent	42
M00514250	Ultraflex Esophageal NG Stent	42



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Polyflex Stent™ Self-expanding Silicone Esophageal Stents		
M00514270	Polyflex Esophageal Stent	43
M00514280	Polyflex Esophageal Stent	43
M00514290	Polyflex Esophageal Stent	43
M00514300	Polyflex Esophageal Stent	43
M00514310	Polyflex Esophageal Stent	43
M00514320	Polyflex Esophageal Stent	43
M00514330	Polyflex Esophageal Stent	43
M00514340	Polyflex Esophageal Stent	43
M00514350	Polyflex Esophageal Stent	43
Radial Jaw™ 4 Hot Biopsy Forceps		
M00515031	Radial Jaw 4 Hot Biopsy Forceps	64
M00515032	Radial Jaw 4 Hot Biopsy Forceps	64
M00515033	Radial Jaw 4 Hot Biopsy Forceps	64
Jagwire™ Pulmonary Guidewire		
M00515171	Jagwire	81
Radial Jaw™ 4 Pulmonary Biopsy Forceps		
M00515181	Radial Jaw 4 Standard Capacity	80
M00515182	Radial Jaw 4 Standard Capacity	80
M00515191	Radial Jaw 4 Standard Capacity w/Needle	80
M00515192	Radial Jaw 4 Standard Capacity w/Needle	80
M00515201	Radial Jaw 4 Large Capacity	80
M00515202	Radial Jaw 4 Large Capacity	80
Celebrity™ Endoscopic Cytology Brushes – Bronchoscope		
M00516001	Celebrity Cytology Brush	80
M00516011	Celebrity Cytology Brush	80
M00516071	Celebrity Cytology Brush	80
M00516151	Celebrity Cytology Brush	80
WallFlex™ Fully Covered RMV Esophageal Stents with Permalume™ Covering (including refractory benign indication)		
M00516210	WallFlex Fully Covered RMV Esophageal Stent	41
M00516220	WallFlex Fully Covered RMV Esophageal Stent	41
M00516230	WallFlex Fully Covered RMV Esophageal Stent	41
M00516240	WallFlex Fully Covered RMV Esophageal Stent	41
M00516250	WallFlex Fully Covered RMV Esophageal Stent	41
M00516260	WallFlex Fully Covered RMV Esophageal Stent	41
WallFlex™ Fully Covered RMV Esophageal Stents Longer Loop with Permalume™ Covering (including refractory benign indication)		
M00516270	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
M00516280	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
M00516290	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
M00516300	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
M00516310	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
M00516320	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	41
WallFlex™ Partially Covered Esophageal Stents with Permalume™ Covering		
M00516900	WallFlex Partially Covered Esophageal Stent	41



Product Code Index (CONTINUED)

Order Number	Description	Page Number
WallFlex™ Partially Covered Esophageal Stents with Permalume™ Covering (CONTINUED)		
M00516910	WallFlex Partially Covered Esophageal Stent	41
M00516920	WallFlex Partially Covered Esophageal Stent	41
M00516930	WallFlex Partially Covered Esophageal Stent	41
M00516940	WallFlex Partially Covered Esophageal Stent	41
M00516950	WallFlex Partially Covered Esophageal Stent	41
Interject™ Clear Sclerotherapy Needles		
M00518111	Interject Sclerotherapy Needle	58
Interject™ Sclerotherapy Needles		
M00518151	Interject Sclerotherapy Needle	58
M00518161	Interject Sclerotherapy Needle	58
M00518251	Interject Sclerotherapy Needle	58
M00518261	Interject Sclerotherapy Needle	58
Interject™ Clear Sclerotherapy Needles		
M00518301	Interject Sclerotherapy Needle	58
M00518311	Interject Sclerotherapy Needle	58
Interject™ Sclerotherapy Needles		
M00518351	Interject Sclerotherapy Needle	58
M00518361	Interject Sclerotherapy Needle	58
Resolution 360™ Clip		
M00521230	Resolution 360 Clip	56
M00521231	Resolution 360 Clip	56
M00521232	Resolution 360 Clip	56
M00521233	Resolution 360 Clip	56
M00521240	Resolution 360 Clip	56
M00521241	Resolution 360 Clip	56
M00521242	Resolution 360 Clip	56
M00521243	Resolution 360 Clip	56
Resolution™ Clip		
M00522600	Resolution Clip	56
M00522601	Resolution Clip	56
M00522602	Resolution Clip	56
M00522610	Resolution Clip	56
M00522611	Resolution Clip	56
M00522612	Resolution Clip	56
Ultratome™ Double-lumen Sphincterotomes		
M00530800	Ultratome Short Nose	15
M00530810	Ultratome Long Nose	15
M00530820	Ultratome Short Nose	15
M00530830	Ultratome Long Nose	15
Contour™ ERCP Cannulas		
M00530850	Contour ERCP Cannula	11
M00530860	Contour ERCP Cannula	11
M00530870	Contour ERCP Cannula	11




Product Code Index (CONTINUED)

Order Number	Description	Page Number
Contour™ ERCP Cannulas (CONTINUED)		
M00530880	Contour ERCP Cannula	11
M00530890	Contour ERCP Cannula	11
Fuoro Tip™ ERCP Cannulas		
M00530960	Fuoro Tip ERCP Cannula, Standard Tip	12
M00530970	Fuoro Tip ERCP Cannula, Tapered Tip	12
M00530980	Fuoro Tip ERCP Cannula, Ultra Tapered Tip	12
C-flex™ Pigtail Biliary Stents		
M00532010	C-flex Pigtail Stent	34
M00532020	C-flex Pigtail Stent	34
M00532030	C-flex Pigtail Stent	34
M00532040	C-flex Pigtail Stent	34
M00532050	C-flex Pigtail Stent	34
M00532060	C-flex Pigtail Stent	34
Advanix™ Stents – Singles		
M00532160	Advanix Stent	31
M00532180	Advanix Stent	31
M00532190	Advanix Stent	31
M00532210	Advanix Stent	31
M00532220	Advanix Stent	31
M00532230	Advanix Stent	31
M00532240	Advanix Stent	31
M00532260	Advanix Stent	31
M00532270	Advanix Stent	31
M00532290	Advanix Stent	31
M00532300	Advanix Stent	31
M00532310	Advanix Stent	31
Microknife™ XL Triple-lumen Needleknife		
M00532810	Microknife XL Triple-lumen Needleknife	15
Advanix™ Stents – Singles		
M00532870	Advanix Stent	30
M00532880	Advanix Stent	30
M00532890	Advanix Stent	30
M00532900	Advanix Stent	30
M00532910	Advanix Stent	30
M00532920	Advanix Stent	30
M00532930	Advanix Stent	30
M00532940	Advanix Stent	30
M00532950	Advanix Stent	30
M00532960	Advanix Stent	30
M00532970	Advanix Stent	30
M00532980	Advanix Stent	30
M00532990	Advanix Stent	30
M00533000	Advanix Stent	30



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Stents – Singles (CONTINUED)		
M00533010	Advanix Stent	30
M00533020	Advanix Stent	30
M00533030	Advanix Stent	30
M00533040	Advanix Stent	30
M00533230	Advanix Stent	30
M00533240	Advanix Stent	30
M00533250	Advanix Stent	30
M00533260	Advanix Stent	30
M00533270	Advanix Stent	30
M00533280	Advanix Stent	30
M00533290	Advanix Stent	30
M00533300	Advanix Stent	30
M00533310	Advanix Stent	30
M00533320	Advanix Stent	30
M00533330	Advanix Stent	30
M00533340	Advanix Stent	30
M00533350	Advanix Stent	30
M00533360	Advanix Stent	30
M00533370	Advanix Stent	30
M00533380	Advanix Stent	30
M00533390	Advanix Stent	30
M00533400	Advanix Stent	30
NaviFlex™ RX Delivery System		
M00533560	RX Delivery System	31
M00533570	RX Delivery System	31
M00533580	RX Delivery System	31
Percuflex™ Amsterdam Biliary Stents		
M00533610	Percuflex Amsterdam Stent	34
M00533620	Percuflex Amsterdam Stent	34
M00533630	Percuflex Amsterdam Stent	34
M00533640	Percuflex Amsterdam Stent	34
M00533650	Percuflex Amsterdam Stent	34
M00533660	Percuflex Amsterdam Stent	34
M00533670	Percuflex Amsterdam Stent	34
M00533680	Percuflex Amsterdam Stent	34
M00533690	Percuflex Amsterdam Stent	34
M00533700	Percuflex Amsterdam Stent	34
Percuflex™ Amsterdam Introducer Kits		
M00533910	Percuflex Amsterdam Stent Kit w/push catheter for use with 7Fr (2.33mm) O.D. stent	34
M00533920	Percuflex Amsterdam Stent Kit w/push and guide catheter for use with 10Fr (3.33mm) O.D. stent	34
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System		
M00534200	Advanix Pre-loaded RX Stent	29
M00534210	Advanix Pre-loaded RX Stent	29



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System <small>(CONTINUED)</small>		
M00534220	Advanix Pre-loaded RX Stent	29
M00534230	Advanix Pre-loaded RX Stent	29
M00534240	Advanix Pre-loaded RX Stent	29
M00534250	Advanix Pre-loaded RX Stent	29
M00534260	Advanix Pre-loaded RX Stent	29
M00534270	Advanix Pre-loaded RX Stent	29
M00534280	Advanix Pre-loaded RX Stent	29
M00534290	Advanix Pre-loaded RX Stent	29
M00534300	Advanix Pre-loaded RX Stent	29
M00534310	Advanix Pre-loaded RX Stent	29
M00534320	Advanix Pre-loaded RX Stent	29
M00534330	Advanix Pre-loaded RX Stent	29
M00534340	Advanix Pre-loaded RX Stent	29
M00534350	Advanix Pre-loaded RX Stent	29
M00534360	Advanix Pre-loaded RX Stent	29
M00534370	Advanix Pre-loaded RX Stent	29
M00534560	Advanix Pre-loaded RX Stent	29
M00534570	Advanix Pre-loaded RX Stent	29
M00534580	Advanix Pre-loaded RX Stent	29
M00534590	Advanix Pre-loaded RX Stent	29
M00534600	Advanix Pre-loaded RX Stent	29
M00534610	Advanix Pre-loaded RX Stent	29
M00534620	Advanix Pre-loaded RX Stent	29
M00534630	Advanix Pre-loaded RX Stent	29
M00534640	Advanix Pre-loaded RX Stent	29
M00534650	Advanix Pre-loaded RX Stent	29
M00534660	Advanix Pre-loaded RX Stent	29
M00534670	Advanix Pre-loaded RX Stent	29
M00534680	Advanix Pre-loaded RX Stent	29
M00534690	Advanix Pre-loaded RX Stent	29
M00534700	Advanix Pre-loaded RX Stent	29
M00534710	Advanix Pre-loaded RX Stent	29
M00534720	Advanix Pre-loaded RX Stent	29
M00534730	Advanix Pre-loaded RX Stent	29
NaviFlex™ RX Pancreatic Delivery System & NaviFlex Pushers		
M00535000	Pancreatic Stent RX Pusher	40
M00535010	Pancreatic Stent RX Pusher	40
M00535020	Pancreatic Stent Pusher	40
M00535030	Pancreatic Stent Pusher	40
M00535060	NaviFlex RX Pancreatic Stent Delivery System	40
M00535080	NaviFlex RX Pancreatic Stent Delivery System	40
Stonetome™ Stone Removal Device		
M00535110	Stonetome (above cut wire)	12



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Stonetome™ Stone Removal Device (CONTINUED)		
M00535130	Stonetome (above cut wire)	12
M00535150	Stonetome (below cut wire)	12
M00535170	Stonetome (below cut wire)	12
M00535190	Stonetome (below cut wire)	12
M00535210	Stonetome (below cut wire)	12
Tandem™ XL Triple-lumen ERCP Cannula		
M00535700	Tandem XL Triple-lumen ERCP Cannula	11
Ultratome XL™ Triple-lumen Sphincterotomes		
M00535900	Ultratome XL Short Nose	15
M00535910	Ultratome XL Long Nose	15
M00535920	Ultratome XL Short Nose	15
M00535930	Ultratome XL Long Nose	15
Advanix™ Pancreatic Single Stents – Straight		
M00536000	Advanix Pancreatic Stent	37
M00536010	Advanix Pancreatic Stent	37
M00536020	Advanix Pancreatic Stent	37
M00536030	Advanix Pancreatic Stent	37
M00536040	Advanix Pancreatic Stent	37
M00536050	Advanix Pancreatic Stent	37
M00536060	Advanix Pancreatic Stent	37
M00536070	Advanix Pancreatic Stent	37
M00536080	Advanix Pancreatic Stent	38
M00536090	Advanix Pancreatic Stent	38
M00536100	Advanix Pancreatic Stent	38
M00536110	Advanix Pancreatic Stent	38
M00536120	Advanix Pancreatic Stent	38
M00536130	Advanix Pancreatic Stent	38
M00536140	Advanix Pancreatic Stent	38
M00536150	Advanix Pancreatic Stent	38
M00536160	Advanix Pancreatic Stent	38
M00536170	Advanix Pancreatic Stent	38
M00536180	Advanix Pancreatic Stent	38
M00536190	Advanix Pancreatic Stent	38
M00536200	Advanix Pancreatic Stent	38
M00536210	Advanix Pancreatic Stent	38
M00536220	Advanix Pancreatic Stent	38
M00536230	Advanix Pancreatic Stent	38
M00536240	Advanix Pancreatic Stent	38
M00536250	Advanix Pancreatic Stent	38
M00536260	Advanix Pancreatic Stent	38
M00536270	Advanix Pancreatic Stent	38
M00536280	Advanix Pancreatic Stent	38
M00536290	Advanix Pancreatic Stent	38



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Straight <small>(CONTINUED)</small>		
M00536300	Advanix Pancreatic Stent	38
M00536310	Advanix Pancreatic Stent	38
M00536320	Advanix Pancreatic Stent	38
M00536330	Advanix Pancreatic Stent	38
M00536340	Advanix Pancreatic Stent	38
M00536350	Advanix Pancreatic Stent	38
M00536360	Advanix Pancreatic Stent	38
M00536370	Advanix Pancreatic Stent	38
M00536380	Advanix Pancreatic Stent	38
M00536390	Advanix Pancreatic Stent	38
M00536400	Advanix Pancreatic Stent	38
M00536410	Advanix Pancreatic Stent	38
M00536420	Advanix Pancreatic Stent	38
M00536430	Advanix Pancreatic Stent	38
M00536440	Advanix Pancreatic Stent	38
M00536450	Advanix Pancreatic Stent	39
M00536460	Advanix Pancreatic Stent	39
M00536470	Advanix Pancreatic Stent	39
M00536480	Advanix Pancreatic Stent	39
M00536490	Advanix Pancreatic Stent	39
M00536500	Advanix Pancreatic Stent	39
M00536510	Advanix Pancreatic Stent	39
M00536520	Advanix Pancreatic Stent	39
M00536530	Advanix Pancreatic Stent	39
M00536540	Advanix Pancreatic Stent	39
M00536550	Advanix Pancreatic Stent	39
M00536560	Advanix Pancreatic Stent	39
M00536570	Advanix Pancreatic Stent	39
M00536580	Advanix Pancreatic Stent	39
M00536590	Advanix Pancreatic Stent	39
M00536600	Advanix Pancreatic Stent	39
M00536610	Advanix Pancreatic Stent	39
M00536620	Advanix Pancreatic Stent	39
M00536630	Advanix Pancreatic Stent	39
M00536640	Advanix Pancreatic Stent	39
M00536650	Advanix Pancreatic Stent	39
M00536660	Advanix Pancreatic Stent	39
M00536670	Advanix Pancreatic Stent	39
M00536680	Advanix Pancreatic Stent	39
M00536690	Advanix Pancreatic Stent	39
M00536700	Advanix Pancreatic Stent	39
M00536710	Advanix Pancreatic Stent	39
M00536720	Advanix Pancreatic Stent	39



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	Advanix™ Pancreatic Single Stents – Straight (CONTINUED)	
M00536730	Advanix Pancreatic Stent	39
	Advanix™ Pancreatic Single Stents – Pigtail	
M00536740	Advanix Pancreatic Stent	36
M00536750	Advanix Pancreatic Stent	36
M00536760	Advanix Pancreatic Stent	36
M00536770	Advanix Pancreatic Stent	36
M00536780	Advanix Pancreatic Stent	36
M00536790	Advanix Pancreatic Stent	36
M00536800	Advanix Pancreatic Stent	36
M00536810	Advanix Pancreatic Stent	36
M00536820	Advanix Pancreatic Stent	36
M00536830	Advanix Pancreatic Stent	36
M00536840	Advanix Pancreatic Stent	36
M00536850	Advanix Pancreatic Stent	36
M00536860	Advanix Pancreatic Stent	36
M00536870	Advanix Pancreatic Stent	36
M00536880	Advanix Pancreatic Stent	36
M00536890	Advanix Pancreatic Stent	36
M00536900	Advanix Pancreatic Stent	36
M00536910	Advanix Pancreatic Stent	36
M00536920	Advanix Pancreatic Stent	36
M00536930	Advanix Pancreatic Stent	36
M00536940	Advanix Pancreatic Stent	36
M00536950	Advanix Pancreatic Stent	36
M00536960	Advanix Pancreatic Stent	36
M00536970	Advanix Pancreatic Stent	36
M00536980	Advanix Pancreatic Stent	36
M00536990	Advanix Pancreatic Stent	36
M00537000	Advanix Pancreatic Stent	36
M00537010	Advanix Pancreatic Stent	36
M00537020	Advanix Pancreatic Stent	36
M00537030	Advanix Pancreatic Stent	36
M00537040	Advanix Pancreatic Stent	36
M00537050	Advanix Pancreatic Stent	36
M00537060	Advanix Pancreatic Stent	36
M00537070	Advanix Pancreatic Stent	36
M00537080	Advanix Pancreatic Stent	36
M00537090	Advanix Pancreatic Stent	36
M00537100	Advanix Pancreatic Stent	37
M00537110	Advanix Pancreatic Stent	37
M00537120	Advanix Pancreatic Stent	37
M00537130	Advanix Pancreatic Stent	37
M00537140	Advanix Pancreatic Stent	37



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Single Stents – Pigtail <small>(CONTINUED)</small>		
M00537150	Advanix Pancreatic Stent	37
M00537160	Advanix Pancreatic Stent	37
M00537170	Advanix Pancreatic Stent	37
M00537180	Advanix Pancreatic Stent	37
M00537190	Advanix Pancreatic Stent	37
M00537200	Advanix Pancreatic Stent	37
M00537210	Advanix Pancreatic Stent	37
M00537220	Advanix Pancreatic Stent	37
M00537230	Advanix Pancreatic Stent	37
M00537240	Advanix Pancreatic Stent	37
M00537250	Advanix Pancreatic Stent	37
M00537260	Advanix Pancreatic Stent	37
M00537270	Advanix Pancreatic Stent	37
M00537280	Advanix Pancreatic Stent	37
M00537290	Advanix Pancreatic Stent	37
M00537300	Advanix Pancreatic Stent	37
M00537310	Advanix Pancreatic Stent	37
M00537320	Advanix Pancreatic Stent	37
M00537330	Advanix Pancreatic Stent	37
M00537340	Advanix Pancreatic Stent	37
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail		
M00537350	Advanix Pancreatic Stent Kit	35
M00537360	Advanix Pancreatic Stent Kit	35
M00537370	Advanix Pancreatic Stent Kit	35
M00537380	Advanix Pancreatic Stent Kit	35
M00537390	Advanix Pancreatic Stent Kit	35
M00537400	Advanix Pancreatic Stent Kit	35
M00537410	Advanix Pancreatic Stent Kit	35
M00537420	Advanix Pancreatic Stent Kit	35
M00537430	Advanix Pancreatic Stent Kit	35
M00537440	Advanix Pancreatic Stent Kit	35
M00537450	Advanix Pancreatic Stent Kit	35
Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight		
M00537460	Advanix Pancreatic Stent Kit	35
M00537470	Advanix Pancreatic Stent Kit	35
M00537480	Advanix Pancreatic Stent Kit	35
M00537490	Advanix Pancreatic Stent Kit	35
M00537500	Advanix Pancreatic Stent Kit	35
M00537510	Advanix Pancreatic Stent Kit	35
M00537520	Advanix Pancreatic Stent Kit	35
M00537530	Advanix Pancreatic Stent Kit	35
M00537540	Advanix Pancreatic Stent Kit	35
M00537550	Advanix Pancreatic Stent Kit	35



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Advanix™ Pancreatic Stent Kits (Stent+ RX Pusher) – Straight (CONTINUED)		
M00537560	Advanix Pancreatic Stent Kit	35
Flexima™ Biliary Stents		
M00539200	Flexima Stent	33
M00539210	Flexima Stent	33
M00539220	Flexima Stent	33
M00539230	Flexima Stent	33
M00539240	Flexima Stent	33
M00539250	Flexima Stent	33
M00539260	Flexima Stent	33
M00539270	Flexima Stent	33
M00539280	Flexima Stent	33
M00539290	Flexima Stent	33
M00539300	Flexima Stent	33
M00539310	Flexima Stent	33
M00539320	Flexima Stent	33
M00539330	Flexima Stent	33
M00539340	Flexima Stent	33
M00539350	Flexima Stent	33
M00539360	Flexima Stent	33
M00539370	Flexima Stent	33
M00539380	Flexima Stent	33
M00539390	Flexima Stent	33
Flexima™ Nasobiliary Stents		
M00540120	Flexima Nasal Stent	33
M00540130	Flexima Nasal Stent	33
M00540140	Flexima Nasal Stent	33
M00540160	Flexima Nasal Stent	33
M00540170	Flexima Nasal Stent	33
M00540180	Flexima Nasal Stent	33
Speedband Superview Super 7™ Multiple Band Ligators		
M00542250	Super 7 Band Ligator	56
M00542251	Super 7 Band Ligator	56
M00542253	Super 7 Band Ligator	56
RX Cytology Brush		
M00545000	RX Cytology Brush	8
Autotome™ RX Cannulating Sphincterotomes		
M00545150	Autotome RX 49	5
M00545160	Autotome RX 49	5
M00545170	Autotome RX 44	5
M00545180	Autotome RX 44	5
M00545190	Autotome RX 39	5
M00545200	Autotome RX 39	5



Product Code Index (CONTINUED)

Order Number	Description	Page Number
RX Accessories		
M00545211	RX Locking Device – Olympus/Fujinon Brand Endoscope Compatible	9
M00545221	RX Locking Device – Pentax Brand Endoscope Compatible	9
M00545261	RX Locking Device and Biopsy Cap System – Olympus Compatible	9
M00545271	RX Locking Device and Biopsy Cap System – Pentax Compatible	9
M00545281	RX Locking Device and Biopsy Cap System – Fujinon Compatible	9
M00545471	RX Locking Device	9
RX Biliary Plastic Stents		
M00545550	RX Biliary Plastic Stent	32
M00545560	RX Biliary Plastic Stent	32
M00545570	RX Biliary Plastic Stent	32
M00545580	RX Biliary Plastic Stent	32
M00545590	RX Biliary Plastic Stent	32
M00545600	RX Biliary Plastic Stent	32
M00545610	RX Biliary Plastic Stent	32
M00545620	RX Biliary Plastic Stent	32
M00545630	RX Biliary Plastic Stent	32
M00545640	RX Biliary Plastic Stent	32
M00545650	RX Biliary Plastic Stent	32
M00545660	RX Biliary Plastic Stent	32
M00545670	RX Biliary Plastic Stent	32
M00545680	RX Biliary Plastic Stent	32
M00545690	RX Biliary Plastic Stent	32
M00545700	RX Biliary Plastic Stent	32
M00545710	RX Biliary Plastic Stent	32
M00545720	RX Biliary Plastic Stent	32
M00545730	RX Biliary Plastic Stent	32
M00545740	RX Biliary Plastic Stent	32
RX Cannulas		
M00545760	RX Standard ERCP Cannula	2
M00545770	RX Tapered ERCP Cannula	2
M00545780	RX Ball Tip ERCP Cannula	2
Tandem™ RX Cannula		
M00545820	Tandem RX Cannula	2
RX Needle Knife XL Triple Lumen Needle Knife		
M00545840	RX Needle Knife XL Triple Lumen Needle Knife	5
Hurricane™ RX Biliary Dilation Balloons		
M00545890	Hurricane RX Dilation Balloon	8
M00545900	Hurricane RX Dilation Balloon	8
M00545910	Hurricane RX Dilation Balloon	8
M00545920	Hurricane RX Dilation Balloon	8
M00545930	Hurricane RX Dilation Balloon	8
M00545940	Hurricane RX Dilation Balloon	8
M00545950	Hurricane RX Dilation Balloon	8



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Hurricane™ RX Biliary Dilation Balloons (CONTINUED)		
M00545960	Hurricane RX Dilation Balloon	8
SpyGlass™		
M00546030	SpyGlass Direct Visualisation Probe	20
M00546040	SpyGlass Ocular	20
M00546050	SpyGlass Storage Tray – Small	20
M00546060	SpyGlass Storage Tray – Large	20
M00546070	3-Joint Arm with Clamp	20
M00546110	SpyGlass Camera PAL	20
M00546120	SpyGlass Video Monitor	20
M00546140	Irrigation Pump with Foot Switch	20
M00546160	SpyGlass Travel Cart	20
M00546190	SpyGlass Lightsource	20
M00546200	SpyGlass Light Source Replacement Bulb	20
M00546210	SpyGlass Light Cable	20
M00546230	SpyScope™ Access and Delivery Catheter	20
M00546240	Isolation Transformer 230 V	20
M00546250	Power Cable Pack	20
M00546270	SpyBite™ Biopsy Forceps	19
SpyGlass™ DS		
M00546270	SpyBite™ Biopsy Forceps	19
M00546451	SpyGlass Irrigation Tube Set	19
SpyGlass™		
M00546451	SpyGlass Irrigation Tube Set	20
SpyGlass™ DS		
M00546600	SpyScope™ DS Access and Delivery Catheter	19
M00546650	SpyGlass DS Digital Controller	19
Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection		
M00547000	Extractor Pro RX Retrieval Balloon Catheter	6
M00547010	Extractor Pro RX Retrieval Balloon Catheter	6
M00547020	Extractor Pro RX Retrieval Balloon Catheter	6
Extractor™ Pro RX Retrieval Balloon Catheter – Proximal Injection		
M00547030	Extractor Pro RX Retrieval Balloon Catheter	6
M00547040	Extractor Pro RX Retrieval Balloon Catheter	6
M00547050	Extractor Pro RX Retrieval Balloon Catheter	6
Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection		
M00547100	Extractor Pro XL Retrieval Balloon Catheter	16
M00547110	Extractor Pro XL Retrieval Balloon Catheter	16
M00547120	Extractor Pro XL Retrieval Balloon Catheter	16
Extractor™ Pro XL Retrieval Balloon Catheter – Proximal Injection		
M00547130	Extractor Pro XL Retrieval Balloon Catheter	16
M00547140	Extractor Pro XL Retrieval Balloon Catheter	16
M00547150	Extractor Pro XL Retrieval Balloon Catheter	16



Product Code Index (CONTINUED)

Order Number	Description	Page Number
	Extractor™ Pro DL Retrieval Balloon Catheter	
M00547200	Extractor Pro DL Retrieval Balloon Catheter	16
M00547210	Extractor Pro DL Retrieval Balloon Catheter	16
M00547220	Extractor Pro DL Retrieval Balloon Catheter	16
	Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection	
M00547300	Extractor Pro RX-S Retrieval Balloon Catheter	7
M00547310	Extractor Pro RX-S Retrieval Balloon Catheter	7
M00547320	Extractor Pro RX-S Retrieval Balloon Catheter	7
	Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection	
M00547330	Extractor Pro RX-S Retrieval Balloon Catheter	7
M00547340	Extractor Pro RX-S Retrieval Balloon Catheter	7
M00547350	Extractor Pro RX-S Retrieval Balloon Catheter	7
	Expect™ Standard Handle Needle	
M00550000	Expect Standard Handle Needle	22
M00550001	Expect Standard Handle Needle	22
M00550010	Expect Standard Handle Needle	22
M00550011	Expect Standard Handle Needle	22
M00550020	Expect Standard Handle Needle	22
M00550021	Expect Standard Handle Needle	22
M00550040	Expect Standard Handle Needle	22
M00550041	Expect Standard Handle Needle	22
	Amplatz™ Super Stiff Guidewire	
M00550090	Amplatz Super Stiff Guidewire	54
	CRE™ Pulmonary Dilation Balloons	
M00550300	CRE Balloon	81
M00550310	CRE Balloon	81
M00550320	CRE Balloon	81
M00550330	CRE Balloon	81
M00550340	CRE Balloon	81
M00550350	CRE Balloon	81
	Alliance™ II Inflation System	
M00550601	60ml Syringe/Gauge Assembly	50, 81
M00550620	Alliance™ II Inflation Handle	7, 50, 81
	NovaGold™ High Performance Guidewires	
M00552001	NovaGold	54
M00552001	NovaGold	54
	Inflation/Deflation Devices and Accessories	
M00553200	Achalasia Pneumatic Hand Pump and Monitor	50
	AXIOS™ Stents and Delivery Systems	
M00553440	AXIOS Stent and Delivery System	23
M00553450	AXIOS Stent and Delivery System	23
	Hot AXIOS™ Stents and Electrocautery Enhanced Delivery Systems	
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	23
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	23



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Hot AXIOS™ Stents and Electrocautery Enhanced Delivery Systems <small>(CONTINUED)</small>		
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	23
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	23
Rigiflex™ II Achalasia Balloon Dilators		
M00554500	Rigiflex II Achalasia Balloon Dilator	49
M00554510	Rigiflex II Achalasia Balloon Dilator	49
M00554520	Rigiflex II Achalasia Balloon Dilator	49
Expect™ Slimline Handle Needle		
M00555500	Expect Slimline Handle Needle	22
M00555501	Expect Slimline Handle Needle	22
M00555510	Expect Slimline Handle Needle	22
M00555511	Expect Slimline Handle Needle	22
M00555520	Expect Slimline Handle Needle	22
M00555521	Expect Slimline Handle Needle	22
M00555530	Expect Slimline Handle Needle	22
M00555531	Expect Slimline Handle Needle	22
Acquire™ EUS-FNB Device		
M00555540	22ga Acquire EUS-FNB Needle	23
M00555541	22ga Acquire EUS-FNB Needle	23
M00555560	25ga Acquire EUS-FNB Needle	23
M00555561	25ga Acquire EUS-FNB Needle	23
Hydra Jagwire™ High Performance Guidewires		
M00556001	Hydra Jagwire Standard	52
M00556011	Hydra Jagwire Standard	52
M00556021	Hydra Jagwire Stiff Shaft	52
M00556031	Hydra Jagwire Stiff Shaft	52
M00556041	Hydra Jagwire Standard	52
M00556051	Hydra Jagwire Standard	52
M00556061	Hydra Jagwire Stiff Shaft	52
M00556071	Hydra Jagwire Stiff Shaft	52
Dreamwire™ High Performance Guidewires		
M00556101	Dreamwire Standard	52
M00556111	Dreamwire Standard	52
M00556121	Dreamwire Stiff Shaft	52
M00556131	Dreamwire Stiff Shaft	52
M00556141	Dreamwire Standard	52
M00556151	Dreamwire Standard	52
M00556161	Dreamwire Stiff Shaft	52
M00556171	Dreamwire Stiff Shaft	52
NaviPro™ 100% Hydrophilic Guidewires		
M00556181	NaviPro Standard	54
M00556191	NaviPro Standard	54
M00556201	NaviPro Standard	54
M00556211	NaviPro Standard	54



Product Code Index (CONTINUED)

Order Number	Description	Page Number
NaviPro™ 100% Hydrophilic Guidewires (CONTINUED)		
M00556221	NaviPro Standard	54
M00556231	NaviPro Standard	54
M00556241	NaviPro Stiff Shaft	54
M00556251	NaviPro Stiff Shaft	54
Jagwire™ High Performance Guidewires		
M00556461	Jagwire Standard	53
M00556481	Jagwire Stiff Shaft	53
M00556491	Jagwire Stiff Shaft	53
M00556511	Jagwire Standard	53
M00556521	Jagwire Stiff Shaft	53
M00556531	Jagwire Stiff Shaft	53
M00556561	Jagwire Standard	53
M00556581	Jagwire Standard	53
M00556591	Jagwire Standard	53
M00556601	Jagwire Stiff Shaft	53
M00556611	Jagwire Stiff Shaft	53
M00556621	Jagwire Stiff Shaft	53
Jagwire™ High Performance Guidewires – Extendable Jagwire		
M00556841	Jagwire Standard	53
M00556851	Jagwire Standard	53
M00556861	Jagwire Stiff Shaft	53
M00556871	Jagwire Stiff Shaft	53
M00556902	Jagwire Jag Tail	53
Ultraflex Precision™ Colonic Stents		
M00557330	Ultraflex Precision Colonic Stent	45
M00557360	Ultraflex Precision Colonic Stent	45
M00557380	Ultraflex Precision Colonic Stent	45
Jagwire™ High Performance Guidewires		
M00557641	Jagwire Standard	53
M00557651	Jagwire Standard	53
M00557661	Jagwire Stiff Shaft	53
M00557671	Jagwire Stiff Shaft	53
Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle – Olympus Scope Compatible		
M00558220	Expect Pulmonary 22ga – Olympus	76
M00558221	Expect Pulmonary 22ga – Olympus	76
M00558250	Expect Pulmonary 25ga – Olympus	76
M00558251	Expect Pulmonary 25ga – Olympus	76
M00558231	Expect Pulmonary Needle Adaptor – Olympus	76
CRE™ Fixed Wire Balloon Dilators		
M00558330	CRE Fixed Wire Balloon Dilator	48
M00558340	CRE Fixed Wire Balloon Dilator	48
M00558350	CRE Fixed Wire Balloon Dilator	48
M00558360	CRE Fixed Wire Balloon Dilator	48



Product Code Index (CONTINUED)

Order Number	Description	Page Number
CRE™ Fixed Wire Balloon Dilators (CONTINUED)		
M00558361	CRE Fixed Wire Balloon Dilator	48
M00558370	CRE Fixed Wire Balloon Dilator	48
M00558371	CRE Fixed Wire Balloon Dilator	48
M00558380	CRE Fixed Wire Balloon Dilator	48
M00558381	CRE Fixed Wire Balloon Dilator	48
CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)		
M00558600	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558610	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558620	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558630	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558640	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558650	CRE PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
CRE™ PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)		
M00558660	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558670	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558680	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558690	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558700	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
M00558710	CRE PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)	47
Gold Probe™ Electrohemostasis Catheters		
M00560070	Gold Probe Catheter	57
M00560071	Gold Probe Catheter	57
M00560100	Gold Probe Catheter	57
M00560101	Gold Probe Catheter	57
Injection Gold Probe™ Electrohemostasis/Injection Catheters		
M00560150	Injection Gold Probe Catheter	57
M00560160	Injection Gold Probe Catheter	57
Gold Probe™ Electrohemostasis Catheters		
M00560220	Gold Probe Catheter	57
Sensation™ Short Throw Snares		
M00560311	Medium Oval	62
M00560321	Small Oval	62
M00560331	Medium Crescent	62
Captivator™ II Snares		
M00561191	25 mm Rounded	63
M00561192	25 mm Rounded	63
M00561193	25 mm Rounded	63
M00561221	10 mm Rounded	63
M00561222	10 mm Rounded	63
M00561223	10 mm Rounded	63
M00561231	15 mm Rounded	63
M00561232	15 mm Rounded	63
M00561233	15 mm Rounded	63



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Captivator™ II Snares (CONTINUED)		
M00561241	20 mm Rounded	63
M00561242	20 mm Rounded	63
M00561243	20 mm Rounded	63
Accessories		
M00561250	Active Cord, Fits Olympus Brand Generator	17, 64
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, ValleyLab®, Aspen Labs and ERBE® Brand Generators	17, 64
Captivator™ II Snares		
M00561291	33 mm Rounded	63
Captivator™ Snares		
M00561311	Large Oval	63
Electrosurgical Bipolar Cable Adaptors		
M00561440	Bipolar Cable Adaptor	59
M00561450	Bipolar Cable Adaptor	59
M00561460	Bipolar Cable Adaptor	59
M00561470	Bipolar Cable Adaptor	59
M00561480	Bipolar Cable Adaptor	59
M00561490	Bipolar Cable Adaptor	59
M00561500	Bipolar Cable Adaptor	59
Captivator™ EMR		
M00561600	Captivator EMR for Standard Gastroscopes	61
M00561610	Captivator EMR for Large Gastroscopes	61
Rotatable Snares		
M00561821	Small Oval	61
M00561831	Medium Oval	61
Captivator™ Snares		
M00562301	Small Oval	63
M00562321	Medium Oval	63
M00562341	Medium Hexagonal	63
M00562371	Medium Crescent	63
M00562391	Large Oval (thin wire)	63
Captiflex™ Snares		
M00562401	Medium Oval	62
M00562402	Medium Oval	62
M00562421	Small Oval	62
M00562422	Small Oval	62
Captivator™ Snares		
M00562451	Small Hexagonal	63
Captiflex™ Snares		
M00562471	Extra Small Oval	62
Profile™ Pediatric Snares		
M00562531	Extra Small Oval	64
M00562551	Small Oval	64
M00562571	Medium Oval	64



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Sensation™ Short Throw Snares		
M00562651	Large Oval	56
M00562652	Large Oval	56
M00562671	Medium Oval	56
M00562672	Medium Oval	56
M00562673	Medium Oval	56
M00562691	Small Oval	56
M00562692	Small Oval	56
M00562693	Small Oval	56
eXcelon™ Transbronchial Aspiration Needle with Syringe		
M00564101	eXcelon Transbronchial Aspiration Needle w/Syringe	72
M00564111	eXcelon Transbronchial Aspiration Needle w/Syringe	72
M00564121	eXcelon Transbronchial Aspiration Needle w/Syringe	72
Ultraflex™ Non-Covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents		
M00564500	Ultraflex Bronchial NG Stent	70
M00564510	Ultraflex Bronchial NG Stent	70
M00564520	Ultraflex Bronchial NG Stent	70
Ultraflex™ Non-Covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents		
M00564530	Ultraflex Tracheal NG Stent	70
M00564540	Ultraflex Tracheal NG Stent	70
M00564560	Ultraflex Tracheal NG Stent	70
M00564570	Ultraflex Tracheal NG Stent	70
M00564590	Ultraflex Tracheal NG Stent	70
M00564600	Ultraflex Tracheal NG Stent	70
Ultraflex™ Non-Covered Tracheobronchial NG Stents – Proximal Release – Bronchial Stents		
M00564640	Ultraflex Bronchial NG Stent	70
M00564650	Ultraflex Bronchial NG Stent	70
M00564660	Ultraflex Bronchial NG Stent	70
M00564670	Ultraflex Bronchial NG Stent	70
M00564680	Ultraflex Bronchial NG Stent	70
M00564690	Ultraflex Bronchial NG Stent	70
M00564700	Ultraflex Bronchial NG Stent	70
Ultraflex™ Non-Covered Tracheobronchial NG Stents – Proximal Release – Tracheal Stents		
M00564710	Ultraflex Tracheal NG Stent	70
M00564720	Ultraflex Tracheal NG Stent	70
Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Bronchial Stents		
M00564740	Ultraflex Bronchial NG Stent	69
M00564750	Ultraflex Bronchial NG Stent	69
M00564760	Ultraflex Bronchial NG Stent	69
M00564770	Ultraflex Bronchial NG Stent	69
M00564780	Ultraflex Bronchial NG Stent	69
M00564790	Ultraflex Bronchial NG Stent	69
M00564800	Ultraflex Bronchial NG Stent	69
M00564810	Ultraflex Bronchial NG Stent	69



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release – Tracheal Stents		
M00564820	Ultraflex Tracheal NG Stent	69
M00564830	Ultraflex Tracheal NG Stent	69
M00564840	Ultraflex Tracheal NG Stent	69
M00564850	Ultraflex Tracheal NG Stent	69
M00564860	Ultraflex Tracheal NG Stent	69
M00564870	Ultraflex Tracheal NG Stent	69
M00564880	Ultraflex Tracheal NG Stent	69
M00564890	Ultraflex Tracheal NG Stent	69
M00564900	Ultraflex Tracheal NG Stent	69
M00564910	Ultraflex Tracheal NG Stent	69
WallFlex™ Duodenal Stents		
M00565010	WallFlex Duodenal Stent	40
M00565020	WallFlex Duodenal Stent	40
M00565030	WallFlex Duodenal Stent	40
WallFlex™ Colonic Stents		
M00565040	WallFlex Colonic Stent	40
M00565050	WallFlex Colonic Stent	40
M00565060	WallFlex Colonic Stent	40
M00565070	WallFlex Colonic Stent	40
M00565080	WallFlex Colonic Stent	40
M00565090	WallFlex Colonic Stent	40
M00565100	WallFlex Colonic Stent	40
M00565110	WallFlex Colonic Stent	40
M00565120	WallFlex Colonic Stent	40
M00565130	WallFlex Colonic Stent	40
M00565140	WallFlex Colonic Stent	40
M00565150	WallFlex Colonic Stent	40
Wallstent™ Colonic and Duodenal Stents		
M00565560	Wallstent Colonic and Duodenal Stent	41
M00565570	Wallstent Colonic and Duodenal Stent	41
M00565580	Wallstent Colonic and Duodenal Stent	41
M00565590	Wallstent Colonic and Duodenal Stent	41
M00565620	Wallstent Colonic and Duodenal Stent	41
M00565630	Wallstent Colonic and Duodenal Stent	41
M00565640	Wallstent Colonic and Duodenal Stent	41
M00565650	Wallstent Colonic and Duodenal Stent	41
Accessories		
M00566001	Rapid Refill Continuous Injection System	8, 15
Through-The-PEG (TTP) Jejunal Feeding Tube		
M00566320	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	63
M00566330	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	63
M00566340	12Fr (4.00mm) TTP Jejunal Feeding Tube	63
M00566350	12Fr (4.00mm) TTP Jejunal Feeding Tube	63



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Accessories		
M00566670	Encore™ 26 Inflation Device	8, 15
Max Force™ Biliary Dilation Balloons		
M00567330	Max Force Dilation Balloon	15
M00567340	Max Force Dilation Balloon	15
M00567370	Max Force Dilation Balloon	15
M00567380	Max Force Dilation Balloon	15
M00567410	Max Force Dilation Balloon	15
M00567420	Max Force Dilation Balloon	15
Low Profile Button Replacements		
M00568170	28Fr (9.33mm) Low Profile Button Replacement	65
M00568180	28Fr (9.33mm) Low Profile Button Replacement	65
M00568190	28Fr (9.33mm) Low Profile Button Replacement	65
M00568260	18Fr (6.00mm) Low Profile Button Replacement	65
M00568270	18Fr (6.00mm) Low Profile Button Replacement	65
M00568280	18Fr (6.00mm) Low Profile Button Replacement	65
M00568290	18Fr (6.00mm) Low Profile Button Replacement	65
M00568300	18Fr (6.00mm) Low Profile Button Replacement	65
M00568320	24Fr (8.00mm) Low Profile Button Replacement	65
M00568330	24Fr (8.00mm) Low Profile Button Replacement	65
M00568340	24Fr (8.00mm) Low Profile Button Replacement	65
M00568350	24Fr (8.00mm) Low Profile Button Replacement	65
Standard PEGs		
M00568381	20Fr Standard PEG Kit, Pull	61
M00568391	24Fr Standard PEG Kit, Pull	61
M00568401	20Fr Standard PEG Kit, Push	61
M00568411	24Fr Standard PEG Kit, Push	61
One Step Button Low Profile PEG Kits		
M00568450	18Fr (6.00mm) One Step Button Kit, Pull	62
M00568460	18Fr (6.00mm) One Step Button Kit, Pull	62
M00568470	18Fr (6.00mm) One Step Button Kit, Pull	62
M00568480	18Fr (6.00mm) One Step Button Kit, Pull	62
M00568500	24Fr (8.00mm) One Step Button Kit, Pull	62
M00568510	24Fr (8.00mm) One Step Button Kit, Pull	62
M00568520	24Fr (8.00mm) One Step Button Kit, Pull	62
M00568530	24Fr (8.00mm) One Step Button Kit, Pull	62
M00568560	18Fr (6.00mm) One Step Button Kit, Push	62
M00568570	18Fr (6.00mm) One Step Button Kit, Push	62
M00568580	18Fr (6.00mm) One Step Button Kit, Push	62
M00568600	24Fr (8.00mm) One Step Button Kit, Push	62
M00568610	24Fr (8.00mm) One Step Button Kit, Push	62
M00568620	24Fr (8.00mm) One Step Button Kit, Push	62
M00568630	24Fr (8.00mm) One Step Button Kit, Push	62



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Wallstent™ RX Biliary Stents		
M00569610	Wallstent RX Biliary	23
M00569620	Wallstent RX Biliary	23
M00569630	Wallstent RX Biliary	23
M00569640	Wallstent RX Biliary	23
M00569650	Wallstent RX Biliary	23
M00569660	Wallstent RX Biliary	23
Wallstent™ RX Biliary Stents with Permalume™ Covering		
M00569670	Wallstent RX w/Permalume	23
M00569680	Wallstent RX w/Permalume	23
M00569690	Wallstent RX w/Permalume	23
M00569700	Wallstent RX w/Permalume	23
M00569710	Wallstent RX w/Permalume	23
M00569720	Wallstent RX w/Permalume	23
Wallstent™ RX Biliary Stents		
M00569800	Wallstent RX Biliary	23
M00569810	Wallstent RX Biliary	23
Polyflex Stent™ Self-expanding Silicone Airway Stents		
M00570000	Polyflex Airway Stent	71
M00570010	Polyflex Airway Stent	71
M00570020	Polyflex Airway Stent	71
M00570030	Polyflex Airway Stent	71
M00570040	Polyflex Airway Stent	71
M00570050	Polyflex Airway Stent	71
M00570060	Polyflex Airway Stent	71
M00570070	Polyflex Airway Stent	71
M00570080	Polyflex Airway Stent	71
M00570090	Polyflex Airway Stent	71
M00570100	Polyflex Airway Stent	71
M00570110	Polyflex Airway Stent	71
M00570120	Polyflex Airway Stent	71
M00570130	Polyflex Airway Stent	71
M00570140	Polyflex Airway Stent	71
M00570150	Polyflex Airway Stent	71
M00570160	Polyflex Airway Stent	71
M00570170	Polyflex Airway Stent	71
M00570180	Polyflex Airway Stent	71
M00570190	Polyflex Airway Stent	71
M00570200	Polyflex Airway Stent	71
M00570210	Polyflex Airway Stent	71
M00570220	Polyflex Airway Stent	71
M00570230	Polyflex Airway Stent	71
M00570240	Polyflex Airway Stent	71
M00570250	Polyflex Airway Stent	71



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Polyflex Stent™ Self-expanding Silicone Airway Stents (CONTINUED)		
M00570260	Polyflex Airway Stent	71
M00570270	Polyflex Airway Stent	71
M00570280	Polyflex Airway Stent	71
M00570290	Polyflex Airway Stent	71
M00570300	Polyflex Airway Stent	71
M00570310	Polyflex Airway Stent	71
M00570320	Polyflex Airway Stent	71
M00570330	Polyflex Airway Stent	71
WallFlex™ Biliary RX Fully Covered Stent System RMV with Permalume™ Covering (including benign indication)		
M00570450	WallFlex RX Biliary w/Permalume	21
M00570460	WallFlex RX Biliary w/Permalume	21
M00570470	WallFlex RX Biliary w/Permalume	21
M00570480	WallFlex RX Biliary w/Permalume	21
M00570490	WallFlex RX Biliary w/Permalume	21
WallFlex™ Biliary RX Non-covered Stents		
M00570600	WallFlex RX Biliary	22
M00570610	WallFlex RX Biliary	22
M00570620	WallFlex RX Biliary	22
M00570630	WallFlex RX Biliary	22
M00570640	WallFlex RX Biliary	22
M00570650	WallFlex RX Biliary	22
M00570660	WallFlex RX Biliary	22
Dynamic™ (Y) Stents		
M00570670	Dynamic (Y) Stent	69
M00570680	Dynamic (Y) Stent	69
M00570690	Dynamic (Y) Stent	69
WallFlex™ Biliary RX Partially Covered Stents with Permalume™ Covering		
M00570700	WallFlex RX Biliary w/Permalume	21
M00570710	WallFlex RX Biliary w/Permalume	21
M00570720	WallFlex RX Biliary w/Permalume	21
M00570730	WallFlex RX Biliary w/Permalume	21
M00570740	WallFlex RX Biliary w/Permalume	21
WallFlex™ Biliary RX Non-covered Stents		
M00570890	WallFlex RX Biliary	22
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire		
M00573000	Jagtome RX 49	3
M00573010	Jagtome RX 49	3
M00573020	Jagtome RX 49	3
M00573030	Jagtome RX 49	3
M00573040	Jagtome RX 44	3
M00573050	Jagtome RX 44	3
M00573060	Jagtome RX 44	3
M00573070	Jagtome RX 44	3



Product Code Index (CONTINUED)

Order Number	Description	Page Number
Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire (CONTINUED)		
M00573080	Jagtome RX 39	3
M00573090	Jagtome RX 39	3
M00573100	Jagtome RX 39	3
M00573110	Jagtome RX 39	3
M00573500	Jagtome RX 49	3
M00573510	Jagtome RX 49	3
M00573520	Jagtome RX 49	3
M00573530	Jagtome RX 49	3
M00573540	Jagtome RX 44	3
M00573550	Jagtome RX 44	3
M00573560	Jagtome RX 44	3
M00573570	Jagtome RX 44	3
WallFlex™ Biliary RX Fully Covered Stent System RMV with Permalume™ Covering (including benign indication)		
M00574210	WallFlex RX Biliary w/Permalume	21
M00574230	WallFlex RX Biliary w/Permalume	21
M00574240	WallFlex RX Biliary w/Permalume	21
WallFlex™ Biliary RX Partially Covered Stents with Permalume™ Covering		
M00576740	WallFlex RX Biliary w/Permalume	21
M00576750	WallFlex RX Biliary w/Permalume	21
M00576760	WallFlex RX Biliary w/Permalume	21
M00576770	WallFlex RX Biliary w/Permalume	21
WallFlex™ Biliary RX Non-covered Stents		
M00576780	WallFlex RX Biliary	22
M00576790	WallFlex RX Biliary	22
WallFlex™ Biliary RX Fully Covered System RMV with Permalume™ Covering (including benign indication)		
M00576800	WallFlex RX Biliary w/Permalume	21
M00576810	WallFlex RX Biliary w/Permalume	21
M00576820	WallFlex RX Biliary w/Permalume	21
M00576830	WallFlex RX Biliary w/Permalume	21
One Step Button and Low Profile Button Accessories		
M00580171	18Fr (6.00mm) Bolus Straight Feeding Adapter	63
M00580181	18Fr (6.00mm) Continuous Right Angle Feeding Adapter	63
M00580231	24Fr (8.00mm) Bolus Straight Feeding Adapter	63
M00580241	24Fr (8.00mm) Continuous Right Angle Feeding Adapter	63
M00580271	28Fr (9.33mm) Bolus Straight Feeding Adapter	63
M00580281	28Fr (9.33mm) Continuous Right Angle Feeding Adapter	63
M00580401	18Fr (6.00mm) Decompression Tube	63
M00580411	18Fr (6.00mm) Decompression Tube	63
M00580421	18Fr (6.00mm) Decompression Tube	63
M00580431	18Fr (6.00mm) Decompression Tube	63
M00580441	18Fr (6.00mm) Decompression Tube	63
M00580461	24Fr (8.00mm) Decompression Tube	63
M00580471	24Fr (8.00mm) Decompression Tube	63



Product Code Index (CONTINUED)


Order Number	Description	Page Number
One Step Button and Low Profile Button Accessories (CONTINUED)		
M00580481	24Fr (8.00mm) Decompression Tube	63
M00580491	24Fr (8.00mm) Decompression Tube	63
PEG Accessories		
M00580641	Y-Port 20Fr	61
M00580651	Y-Port 24Fr	61
M00580731	20Fr (6.67mm) Round Bolsters	61
M00580781	20Fr (6.67mm) and 24F (8.00mm) C-Clamps	61
Standard Balloon Replacements		
M00582040	12Fr (4.00mm) Replacement G-Tube, Straight	64
M00582050	14Fr (4.67mm) Replacement G-Tube, Straight	64
M00582060	16Fr (5.33mm) Replacement G-Tube, Straight	64
M00582070	18Fr (6.00mm) Replacement G-Tube, Straight	64
M00582080	20Fr (6.67mm) Replacement G-Tube, Straight	64
M00582090	24Fr (8.00mm) Replacement G-Tube, Straight	64
M00582100	28Fr (9.33mm) Replacement G-Tube, Straight	64
M00582110	12Fr (4.00mm) Replacement G-Tube, Right Angle	64
M00582120	14Fr (4.67mm) Replacement G-Tube, Right Angle	64
M00582130	16Fr (5.33mm) Replacement G-Tube, Right Angle	64
M00582140	18Fr (6.00mm) Replacement G-Tube, Right Angle	64
M00582150	20Fr (6.67mm) Replacement G-Tube, Right Angle	64
M00582160	24Fr (8.00mm) Replacement G-Tube, Right Angle	64
M00582170	28Fr (9.33mm) Replacement G-Tube, Right Angle	64
M00582860	22Fr (7.33mm) Replacement G-Tube, Straight	64
M00582870	22Fr (7.33mm) Replacement G-Tube, Right Angle	64
Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire		
M00583000	Hydratome RX 49	2
M00583010	Hydratome RX 49	2
M00583020	Hydratome RX 49	2
M00583030	Hydratome RX 49	2
M00583040	Hydratome RX 44	2
M00583050	Hydratome RX 44	2
M00583060	Hydratome RX 44	2
M00583070	Hydratome RX 44	2
TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire		
M00583080	TRUEtome 39 Sphincterotome	11
M00583090	TRUEtome 44 Sphincterotome	11
M00583100	TRUEtome 44 Sphincterotome	11
M00583110	TRUEtome 44 Sphincterotome	11
M00583120	TRUEtome 39 Sphincterotome	11
M00583130	TRUEtome 44 Sphincterotome	11
M00583140	TRUEtome 44 Sphincterotome	11
M00583150	TRUEtome 44 Sphincterotome	11



Product Code Index (CONTINUED)


Order Number	Description	Page Number
TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire		
M00583160	TRUEtome 44 Sphincterotome	13
M00583170	TRUEtome 44 Sphincterotome	13
TRUEtome™ Dreamwire Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00583190	TRUEtome 44 Sphincterotome	13
M00583210	TRUEtome 44 Sphincterotome	13
Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire		
M00584000	Dreamtome RX 49	3
M00584010	Dreamtome RX 49	3
M00584020	Dreamtome RX 49	3
M00584030	Dreamtome RX 49	3
M00584040	Dreamtome RX 44	3
M00584050	Dreamtome RX 44	3
M00584060	Dreamtome RX 44	3
M00584070	Dreamtome RX 44	3
TRUEtome™ Cannulating Sphincterotomes		
M00584120	TRUEtome 39	14
M00584121	TRUEtome 39	14
M00584122	TRUEtome 39	14
M00584150	TRUEtome 49	14
M00584151	TRUEtome 49	14
M00584152	TRUEtome 49	14
M00584160	TRUEtome 49	14
M00584161	TRUEtome 49	14
M00584162	TRUEtome 49	14
M00584170	TRUEtome 44	14
M00584171	TRUEtome 44	14
M00584172	TRUEtome 44	14
M00584180	TRUEtome 44	14
M00584181	TRUEtome 44	14
M00584182	TRUEtome 44	14
M00584190	TRUEtome 39	14
M00584191	TRUEtome 39	14
M00584192	TRUEtome 39	14
Max Force™ TTS Esophageal Balloon Dilators		
M00589060	Max Force TTS Balloon Dilator	49
M00589080	Max Force TTS Balloon Dilator	49
M00589100	Max Force TTS Balloon Dilator	49
M00589120	Max Force TTS Balloon Dilator	49
M00589140	Max Force TTS Balloon Dilator	49
M00589150	Max Force TTS Balloon Dilator	49
M00589160	Max Force TTS Balloon Dilator	49
M00589180	Max Force TTS Balloon Dilator	49





RX Biliary System™


CANNULAS



Tandem™ RX Cannula

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545820	Tandem RX Cannula	7-5.5 / 2.33-1.84	210	.035 / 0.89

Packaged one per box





RX Cannulas

Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00545760	RX Standard ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545770	RX Tapered ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89
M00545780	RX Ball Tip ERCP Cannula	8.5-5 / 2.83-1.67	210	.035 / 0.89

Packaged one per box

RX Biliary System™

SPHINCTEROTOMES





Hydratome™ RX Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Pre-loaded Guidewire** Length (cm)	O.D. (in) / (mm)
M00583000	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583010	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	260	.035 / 0.89
M00583020	Hydratome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583030	Hydratome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	450	.035 / 0.89
M00583040*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583050*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	260	.035 / 0.89
M00583060*	Hydratome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89
M00583070*	Hydratome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	450	.035 / 0.89

Packaged one per box
 * Includes handle brake.
 ** Pre-loaded with .035in (0.89mm) Hydra Jagwire™ Guidewire.



Dreamtome™ RX Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Pre-loaded Guidewire** Length (cm)	O.D. (in) / (mm)
M00584000	Dreamtome RX 49	20	5	4.9 / 1.63	260	.035 / 0.89
M00584010	Dreamtome RX 49	30	5	4.9 / 1.63	260	.035 / 0.89
M00584020	Dreamtome RX 49	20	5	4.9 / 1.63	450	.035 / 0.89
M00584030	Dreamtome RX 49	30	5	4.9 / 1.63	450	.035 / 0.89
M00584040*	Dreamtome RX 44	20	5	4.4 / 1.47	260	.035 / 0.89
M00584050*	Dreamtome RX 44	30	5	4.4 / 1.47	260	.035 / 0.89
M00584060*	Dreamtome RX 44	20	5	4.4 / 1.47	450	.035 / 0.89
M00584070*	Dreamtome RX 44	30	5	4.4 / 1.47	450	.035 / 0.89

Packaged one per box
 * Includes handle brake.
 ** Pre-loaded with Straight-tip guidewires only.

RX Biliary System™

SPHINCTEROTOMES (CONTINUED)



Jagtome™ RX Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Pre-loaded Guidewire**		
							Shaft Stiffness	Length (cm)	O.D. (in) / (mm)
M00573000	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573020	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	260	.035 / 0.89
M00573010	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573030	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Standard	450	.035 / 0.89
M00573500	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573520	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	260	.035 / 0.89
M00573510	Jagtome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573530	Jagtome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	Stiff	450	.035 / 0.89
M00573040†	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573060†	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	260	.035 / 0.89
M00573050†	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573070†	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Standard	450	.035 / 0.89
M00573540†	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573560†	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	260	.035 / 0.89
M00573550†	Jagtome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573570†	Jagtome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	Stiff	450	.035 / 0.89
M00573080†	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573100†	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Stiff	260	.025 / 0.64
M00573090†	Jagtome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64
M00573110†	Jagtome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	Standard	450	.025 / 0.64



Packaged one per box

* Includes handle brake.

** Pre-loaded with the high performance Jagwire™ Guidewire (Jagtome RX 49 and RX 44 contain the Stiff Shaft Jagwire and Jagtome RX 39 the Standard Jagwire).

RX Biliary System™

SPHINCTEROTOMES (CONTINUED)



Autotome™ RX Cannulating Sphincterotomes

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00545150	Autotome RX 49	20	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545160	Autotome RX 49	30	7-5.5 / 2.33-1.83	5	4.9 / 1.63	200	.035 / 0.89
M00545170*	Autotome RX 44	20	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545180*	Autotome RX 44	30	7-5.5 / 2.33-1.83	5	4.4 / 1.47	200	.035 / 0.89
M00545190*	Autotome RX 39	20	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64
M00545200*	Autotome RX 39	30	7-5.5 / 2.33-1.83	5	3.9 / 1.30	200	.025 / 0.64

Packaged one per box
Recommended Guidewires: For Autotome RX 44 and 49, recommended guidewire is .035in (0.89mm) Dreamwire™ High Performance Guidewire
 For Autotome RX 39, recommended guidewire is .025in (0.69mm) Jagwire™ Guidewire
 * Includes handle brake.




RX Needle Knife XL Triple Lumen Needle Knife

Order Number	Description	Cut Wire (mm)	Shaft O.D. (F) / (mm)	Tip Length (mm)	Tip Diameter (F) / (mm)
M00545840	RX Needle Knife XL Triple Lumen Needle Knife	5	7-5.5 / 2.3-1.83	5	5.5 / 1.8

Packaged one per box
Recommended Guidewire: .035 in (0.89mm) Dreamwire™ High Performance Guidewire

RX Biliary System™

STONE RETRIEVAL




Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547000	Extractor Pro RX Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547010	Extractor Pro RX Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547020	Extractor Pro RX Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro RX Retrieval Balloon Catheter – Proximal Injection




Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547030	Extractor Pro RX Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547040	Extractor Pro RX Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547050	Extractor Pro RX Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

RX Biliary System™

STONE RETRIEVAL




NEW Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547300	Extractor Pro RX-S Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547310	Extractor Pro RX-S Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547320	Extractor Pro RX-S Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




NEW Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection



Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547330	Extractor Pro RX-S Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547340	Extractor Pro RX-S Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547350	Extractor Pro RX-S Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Trapezoid™ RX Wire-guided Retrieval Baskets



Order Number	Description	Open Basket Diameter (cm)	Required Working Channel (mm)
M00510860	Trapezoid RX Wire-guided Retrieval Basket	1.5	3.2
M00510870	Trapezoid RX Wire-guided Retrieval Basket	2.0	3.2
M00510880	Trapezoid RX Wire-guided Retrieval Basket	2.5	3.2
M00510890*	Trapezoid RX Wire-guided Retrieval Basket	3.0	3.2
M00550620	Alliance™ II Inflation Handle**		

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* The Trapezoid 3 x 6cm basket is designed for crushing calculus larger than 1.5cm in diameter.

** Alliance II recommended for use as Mechanical Lithotripter. Alliance II is reusable.



Contents




Alphabetical Index

Product Code Index

RX Biliary System™

DILATION



Hurricane™ RX Biliary Dilation Balloons

Order Number	Description	Balloon O.D. (mm)	Balloon Length (cm)	Catheter O.D. (F) / (mm)	Maximum Inflation Pressure (ATM) / (kPa)	Minimum Channel Size (mm)	Recommended Guidewire (in) / (mm)
M00545890	Hurricane RX Dilation Balloon	4	2	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545900	Hurricane RX Dilation Balloon	4	4	5.8 / 1.93	12 / 1216	3.2	.035 / 0.89
M00545910	Hurricane RX Dilation Balloon	6	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545920	Hurricane RX Dilation Balloon	6	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545930	Hurricane RX Dilation Balloon	8	2	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545940	Hurricane RX Dilation Balloon	8	4	5.8 / 1.93	11 / 1114	3.2	.035 / 0.89
M00545950	Hurricane RX Dilation Balloon	10	2	5.8 / 1.93	8 / 810	3.2	.035 / 0.89
M00545960	Hurricane RX Dilation Balloon	10	4	5.8 / 1.93	8 / 810	3.2	.035 / 0.89

Packaged one per box
Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670

CYTOLOGY



RX Cytology Brush

Order Number	Description	Brush Width (mm)	Catheter O.D. (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)
M00545000	RX Cytology Brush	2.1	8 / 2.67	200	.035 / 0.89

Packaged one per box

RX Biliary System™


ACCESSORIES

RX Accessories



Order Number	Description	Units
M00545471	RX Locking Device	box 10
M00545211†	RX Locking Device – Olympus/Fujinon Brand Endoscope Compatible	box 10
M00545221†	RX Locking Device – Pentax Brand Endoscope Compatible	box 10
M00545261	RX Locking Device and Biopsy Cap System – Olympus Compatible	box 10
M00545271	RX Locking Device and Biopsy Cap System – Pentax Compatible	box 10
M00545281	RX Locking Device and Biopsy Cap System – Fujinon Compatible	box 10
M00566001	Rapid Refill Continuous Injection System	box 10

† Product not available in the European Economic Area (EEA), except in Greece. Not for use or distribution in the EEA, except in Greece. Please check availability of the product with your local sales representative or your local customer service.





Encore™ 26 Inflation Device

Order Number	Description	Units
M00566670	Encore 26 Inflation Device	box 1


RX Biliary Metal and Plastic Stents are located in the Stents Section on pages 25 – 28 and 29 – 32

Biliary Guidewires are located in the Guidewires Section on page 52



Biliary

CANNULAS




Tandem™ XL Triple-lumen ERCP Cannula

Order Number	Description	Shaft O.D. (F) / (mm)	Distal O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
--------------	-------------	-----------------------	------------------------	-------------	-----	-----------------------------------

M00535700	Tandem XL Triple-lumen ERCP Cannula	7.0-5.5 / 2.33-1.83	5.5 / 1.83	210	Tapered	.035 / 0.89
-----------	-------------------------------------	---------------------	------------	-----	---------	-------------

Packaged one per box



Contour™ ERCP Cannulas


Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Tip	Recommended Guidewire (in) / (mm)
--------------	-------------	-----------------------	-------------	-----	-----------------------------------

M00530850	Contour ERCP Cannula	5 / 1.67	210	Standard	.035 / 0.89
M00530860	Contour ERCP Cannula	5 / 1.67	210	Tapered	.035 / 0.89
M00530870	Contour ERCP Cannula	5 / 1.67	210	Ultra Tapered	.025 / 0.64
M00530880	Contour ERCP Cannula	5 / 1.67	210	5-4-3	.018 / 0.46
M00530890	Contour ERCP Cannula	5 / 1.67	210	Ball Tip	.035 / 0.89

Packaged one per box

Biliary

CANNULAS (CONTINUED)




Fluoro Tip™ ERCP Cannulas



Order Number	Description	Shaft O.D. (F) / (mm)	Length (cm)	Recommended Guidewire (in) / (mm)
M00530960	Fluoro Tip ERCP Cannula, Standard Tip	5.5 / 1.83	210	.035 / 0.89
M00530970	Fluoro Tip ERCP Cannula, Tapered Tip	5.5 / 1.83	210	.035 / 0.89
M00530980	Fluoro Tip ERCP Cannula, Ultra Tapered Tip	5.5 / 1.83	210	.025 / 0.64

Packaged one per box

SPHINCTEROTOMES



Stonetome™ Stone Removal Device




Order Number	Description	Cut Wire (mm)	Balloon Diameter (mm)	Tip Length (mm)	Catheter Length (cm)
M00535150	Stonetome (below cut wire)	20	11.5	5	200
M00535170	Stonetome (below cut wire)	20	11.5	20	200
M00535190	Stonetome (below cut wire)	30	11.5	5	200
M00535210	Stonetome (below cut wire)	30	11.5	20	200
M00535110	Stonetome (above cut wire)	20	11.5	20	200
M00535130	Stonetome (above cut wire)	30	11.5	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

SPHINCTEROTOMES (CONTINUED)



TRUEtome™ Jag Cannulating Sphincterotomes with Pre-loaded Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in/cm)	Units
M00583080	TRUEtome 39 Sphincterotome	20	3.9 / 1.30	5	.025/450 Jagwire Guidewire	box 1
M00583090	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.025/450 Jagwire Guidewire	box 1
M00583100	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035/450 Jagwire Guidewire	box 1
M00583110	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035/450 Jagwire Guidewire	box 1
M00583120	TRUEtome 39 Sphincterotome	30	3.9 / 1.30	5	.025/450 Jagwire Guidewire	box 1
M00583130	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.025/450 Jagwire Guidewire	box 1
M00583140	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035/450 Jagwire Guidewire	box 1
M00583150	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035/450 Jagwire Guidewire	box 1

TRUEtome™ Dreamwire Cannulating Sphincterotomes with Pre-loaded Dreamwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in/cm)	Units
M00583190	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035/450 Dreamwire Guidewire	box 1
M00583210	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035/450 Dreamwire Guidewire	box 1


TRUEtome™ Hydra Cannulating Sphincterotomes with Pre-loaded Hydra Jagwire™ Guidewire



Order Number	Description	Cut Wire (mm)	Tip O.D. (F) / (mm)	Tip Length (mm)	Pre-loaded Guidewire (in/cm)	Units
M00583160	TRUEtome 44 Sphincterotome	20	4.4 / 1.47	5	.035/450 Hydra Jagwire Guidewire	box 1
M00583170	TRUEtome 44 Sphincterotome	30	4.4 / 1.47	5	.035/450 Hydra Jagwire Guidewire	box 1

Biliary

SPHINCTEROTOMES (CONTINUED)



TRUEtome™ Cannulating Sphincterotomes




Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)	Units
M00584190	TRUEtome 39	20	3.9 / 1.30	5	200	box 1
M00584120	TRUEtome 39	30	3.9 / 1.30	5	200	box 1
M00584191	TRUEtome 39	20	3.9 / 1.30	5	200	box 10
M00584121	TRUEtome 39	30	3.9 / 1.30	5	200	box 10
M00584192	TRUEtome 39	20	3.9 / 1.30	5	200	box 20
M00584122	TRUEtome 39	30	3.9 / 1.30	5	200	box 20
M00584170	TRUEtome 44	20	4.4 / 1.47	5	200	box 1
M00584180	TRUEtome 44	30	4.4 / 1.47	5	200	box 1
M00584171	TRUEtome 44	20	4.4 / 1.47	5	200	box 10
M00584181	TRUEtome 44	30	4.4 / 1.47	5	200	box 10
M00584172	TRUEtome 44	20	4.4 / 1.47	5	200	box 20
M00584182	TRUEtome 44	30	4.4 / 1.47	5	200	box 20
M00584150	TRUEtome 49	20	4.9 / 1.63	5	200	box 1
M00584160	TRUEtome 49	30	4.9 / 1.63	5	200	box 1
M00584151	TRUEtome 49	20	4.9 / 1.63	5	200	box 10
M00584161	TRUEtome 49	30	4.9 / 1.63	5	200	box 10
M00584152	TRUEtome 49	20	4.9 / 1.63	5	200	box 20
M00584162	TRUEtome 49	30	4.9 / 1.63	5	200	box 20

Recommended Guidewire: .025in (0.64mm) for TRUEtome 39; .035in (0.89mm) for TRUEtome 44 and 49

Biliary

SPHINCTEROTOMES (CONTINUED)




Ultratome XL™ Triple-lumen Sphincterotomes



Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00535900	Ultratome XL Short Nose	20	5.5 / 1.83	5	200
M00535910	Ultratome XL Long Nose	20	5.5 / 1.83	20	200
M00535920	Ultratome XL Short Nose	30	5.5 / 1.83	5	200
M00535930	Ultratome XL Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Ultratome™ Double-lumen Sphincterotomes



Order Number	Description	Cut Wire (mm)	Distal O.D. (F) / (mm)	Tip Length (mm)	Catheter Length (cm)
M00530800	Ultratome Short Nose	20	5.5 / 1.83	5	200
M00530810	Ultratome Long Nose	20	5.5 / 1.83	20	200
M00530820	Ultratome Short Nose	30	5.5 / 1.83	5	200
M00530830	Ultratome Long Nose	30	5.5 / 1.83	20	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Microknife™ XL Triple-lumen Needleknife




Order Number	Description	Shaft O.D. (F) / (mm)	Catheter Length (cm)
M00532810	Microknife XL Triple-lumen Needleknife	7-5.5 / 2.33-1.83	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

STONE RETRIEVAL




Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547100	Extractor Pro XL Retrieval Balloon Catheter	9.0-12.0	Above	7-6 / 2.33-2.00	200
M00547110	Extractor Pro XL Retrieval Balloon Catheter	12.0-15.0	Above	7-6 / 2.33-2.00	200
M00547120	Extractor Pro XL Retrieval Balloon Catheter	15.0-18.0	Above	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire




Extractor™ Pro XL Retrieval Balloon Catheter – Proximal Injection

Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547130	Extractor Pro XL Retrieval Balloon Catheter	9.0-12.0	Below	7-6 / 2.33-2.00	200
M00547140	Extractor Pro XL Retrieval Balloon Catheter	12.0-15.0	Below	7-6 / 2.33-2.00	200
M00547150	Extractor Pro XL Retrieval Balloon Catheter	15.0-18.0	Below	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Extractor™ Pro DL Retrieval Balloon Catheter


Order Number	Description	Balloon O.D. (mm)	Injection	Catheter O.D. (F) / (mm)	System Length (cm)
M00547200	Extractor Pro DL Retrieval Balloon Catheter	9.0-12.0	N/A	7-6 / 2.33-2.00	200
M00547210	Extractor Pro DL Retrieval Balloon Catheter	12.0-15.0	N/A	7-6 / 2.33-2.00	200
M00547220	Extractor Pro DL Retrieval Balloon Catheter	15.0-18.0	N/A	7-6 / 2.33-2.00	200

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Biliary

DILATION



Max Force™ Biliary Dilation Balloons

Order Number	Description	Inflated O.D. (mm)	Balloon Length (cm)	Catheter Size (F) / (mm)	Required Working Channel (mm)
M00567330	Max Force Dilation Balloon	4	2	5 / 1.67	2.8
M00567340	Max Force Dilation Balloon	4	4	5 / 1.67	2.8
M00567370	Max Force Dilation Balloon	6	2	5 / 1.67	2.8
M00567380	Max Force Dilation Balloon	6	4	5 / 1.67	2.8
M00567410	Max Force Dilation Balloon	8	2	5 / 1.67	2.8
M00567420	Max Force Dilation Balloon	8	3	5 / 1.67	2.8

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Recommended Inflation Device: Encore™ 26 Inflation Device, M00566670


ACCESSORIES

Accessories

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, ValleyLab®, Aspen Labs and ERBE® Brand Generators	box 1
M00566670	⊗ Encore 26 Inflation Device	box 1
M00566001	⊗ Rapid Refill Continuous Injection System	box 10



Biliary Metal and Plastic Stents are located in the Stents Section on pages 28 and 33 – 34

Biliary Guidewires are located in the Guidewires Section on page 52



SpyGlass™ Direct Visualisation Systems

SPYGLASS™ DS



SpyGlass™ DS Capital Equipment

Order Number	Description	Units
M00546650	SpyGlass DS Digital Controller	box 1


SpyGlass™ DS Disposable Equipment



Order Number	Description	Units
M00546600	SpyScope™ DS Access and Delivery Catheter	box 1
M00546270	SpyBite™ Biopsy Forceps	box 1
M00546451	SpyGlass Irrigation Tube Set	box 10

SpyGlass™ Direct Visualisation Systems

SPYGLASS™



SpyGlass™ Capital Equipment

Order Number	Description	Units
M00546110	SpyGlass Camera PAL*	box 1
M00546120	SpyGlass Video Monitor	box 1
M00546160	SpyGlass Travel Cart ¹	box 1
M00546190	SpyGlass Lightsource	box 1
M00546210	SpyGlass Light Cable	box 1
M00546040	SpyGlass Ocular	box 1
M00546240	Isolation Transformer 230 V**	box 1
M00546140	Irrigation Pump with Foot Switch	box 1
M00546060	SpyGlass Storage Tray – Large ²	box 1
M00546050	SpyGlass Storage Tray – Small ²	box 1
M00546250	Power Cable Pack	box 1

- 3-Joint Arm (without clamp) is included with the Travel Cart (M00546160). The 3-Joint Arm with Clamp (M00546070) is useful if the account does not want to purchase the Travel Cart as it allows the account to mount the arm to another cart or table.
 - Purchase a minimum of two storage trays. Large storage trays are recommended unless they are to be reprocessed in an Olympus® Medivator Machine, into which the large storage trays will not fit. In this case, two small storage trays (M00546050) should be ordered.
- * For countries providing mains voltage of 120V please order M00546100 Spyglass Camera NTSC.
 ** For countries providing mains voltage of 120V please order M00546260 Isolation Transformer 120V.

SpyGlass™ Reusable Equipment

Order Number	Description	Units
M00546030	SpyGlass Direct Visualisation Probe ¹	box 1

- SpyGlass Probes have been validated to 20 cycles of reprocessing without degradation of image. The true number of procedural uses, however, will depend on how the SpyGlass Probe is handled. Fiberoptics are very sensitive and need delicate handling.

SpyGlass™ Disposable Equipment

Order Number	Description	Units
M00546230	SpyScope™ Access and Delivery Catheter	box 1
M00546270	SpyBite™ Biopsy Forceps	box 1
M00546451	SpyGlass Irrigation Tube Set	box 10

SpyGlass™ Optional Accessories


Order Number	Description	Units
M00546070	3-Joint Arm with Clamp ¹	box 1
M00546200	SpyGlass Light Source Replacement Bulb	box 1

- 3-Joint Arm (without clamp) is included with the Travel Cart (M00546160). The 3-Joint Arm with Clamp (M00546070) is useful if the account does not want to purchase the Travel Cart as it allows the account to mount the arm to another cart or table.



Endoscopic Ultrasound (EUS)

FINE NEEDLE ASPIRATION (FNA)



Expect™ Standard Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00550000	Expect Standard Handle Needle	19 / 1.10	2.8	1.83	box 1
M00550001	Expect Standard Handle Needle	19 / 1.10	2.8	1.83	box 5
M00550040	Expect Standard Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00550041	Expect Standard Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00550010	Expect Standard Handle Needle	22 / 0.72	2.4	1.65	box 1
M00550011	Expect Standard Handle Needle	22 / 0.72	2.4	1.65	box 5
M00550020	Expect Standard Handle Needle	25 / 0.52	2.4	1.52	box 1
M00550021	Expect Standard Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock.
Working length: 137.5cm to 141.5cm, adjustable.
Needle length: 0cm to 8cm, adjustable.



Expect™ Slimline Handle Needles*



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555500	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 1
M00555501	Expect Slimline Handle Needle	19 / 1.10	2.8	1.83	box 5
M00555530	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 1
M00555531	Expect Slimline Handle Needle	19 Flex / 1.14	2.8	1.73	box 5
M00555510	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 1
M00555511	Expect Slimline Handle Needle	22 / 0.72	2.4	1.65	box 5
M00555520	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 1
M00555521	Expect Slimline Handle Needle	25 / 0.52	2.4	1.52	box 5

* Includes a 20cc syringe and one-way stopcock.
Working length: 137.5cm to 141.5cm, adjustable.
Needle length: 0cm to 8cm, adjustable.

Endoscopic Ultrasound (EUS)

FINE NEEDLE BIOPSY (FNB)



NEW Acquire™ EUS-FNB Device



Order Number	Description	Needle Size (gauge) / (mm)	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00555540	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 1
M00555541	22ga Acquire EUS-FNB Needle	22 / 0.72	2.4	1.65	box 5
M00555560	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 1
M00555561	25ga Acquire EUS-FNB Needle	25 / 0.52	2.4	1.52	box 5

STENTS AND DELIVERY SYSTEMS

Hot AXIOS™ Stents and Electrocautery Enhanced Delivery Systems



Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553520	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	14	6	8	9	138	146	box 1
M00553530	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	17	8	8	9	138	146	box 1
M00553540	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	21	10	10	10.8	138	146	box 1
M00553550	Hot AXIOS Stent and Electrocautery Enhanced Delivery System	24	15	10	10.8	138	146	box 1

AXIOS™ Stents and Delivery Systems



Order Number	Description	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter Outer Diameter (F)	Catheter Working Length (cm)	Catheter Total Length (cm)	Units
M00553440	AXIOS Stent and Delivery System	21	10	10	10.8	138	146	box 1
M00553450	AXIOS Stent and Delivery System	24	15	10	10.8	138	146	box 1




WALLFLEX™ ESOPHAGEAL
 FULLY COVERED STENT
 APPROVED FOR
 REFRACTORY
BENIGN
 INDICATION

WALLFLEX™ BILIARY RX
 FULLY COVERED STENT
 ALSO
 APPROVED FOR
BENIGN
 INDICATION

Stents

RX BILIARY METAL STENTS




WallFlex™ Biliary RX Fully Covered Stent System RMV with Permalume™ Covering (including benign indication)



Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00574210†	WallFlex RX Biliary w/Permalume	6	40	8.5 / 2.83	.035 / 0.89
M00574230†	WallFlex RX Biliary w/Permalume	6	60	8.5 / 2.83	.035 / 0.89
M00574240†	WallFlex RX Biliary w/Permalume	6	80	8.5 / 2.83	.035 / 0.89
M00570450†	WallFlex RX Biliary w/Permalume	8	60	8.5 / 2.83	.035 / 0.89
M00570460†	WallFlex RX Biliary w/Permalume	8	80	8.5 / 2.83	.035 / 0.89
M00576800†	WallFlex RX Biliary w/Permalume	8	100	9.0 / 3.00	.035 / 0.89
M00576810†	WallFlex RX Biliary w/Permalume	8	120	9.0 / 3.00	.035 / 0.89
M00570470†	WallFlex RX Biliary w/Permalume	10	40	8.5 / 2.83	.035 / 0.89
M00570480†	WallFlex RX Biliary w/Permalume	10	60	8.5 / 2.83	.035 / 0.89
M00570490†	WallFlex RX Biliary w/Permalume	10	80	8.5 / 2.83	.035 / 0.89
M00576820†	WallFlex RX Biliary w/Permalume	10	100	9.0 / 3.00	.035 / 0.89
M00576830†	WallFlex RX Biliary w/Permalume	10	120	9.0 / 3.00	.035 / 0.89

† Product commercialised only in countries where benign indication is approved.



WallFlex™ Biliary RX Partially Covered Stents with Permalume™ Covering





Order Number	Description	Diameter (mm)	Length (mm)	Covered Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570700	WallFlex RX Biliary w/Permalume	8	60	48	8.5 / 2.83	.035 / 0.89
M00570710	WallFlex RX Biliary w/Permalume	8	80	68	8.5 / 2.83	.035 / 0.89
M00576740	WallFlex RX Biliary w/Permalume	8	100	88	9.0 / 3.00	.035 / 0.89
M00576750	WallFlex RX Biliary w/Permalume	8	120	108	9.0 / 3.00	.035 / 0.89
M00570720	WallFlex RX Biliary w/Permalume	10	40	28	8.5 / 2.83	.035 / 0.89
M00570730	WallFlex RX Biliary w/Permalume	10	60	48	8.5 / 2.83	.035 / 0.89
M00570740	WallFlex RX Biliary w/Permalume	10	80	68	8.5 / 2.83	.035 / 0.89
M00576760	WallFlex RX Biliary w/Permalume	10	100	88	9.0 / 3.00	.035 / 0.89
M00576770	WallFlex RX Biliary w/Permalume	10	120	108	9.0 / 3.00	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

Stents

RX BILIARY METAL STENTS (CONTINUED)



WallFlex™ Biliary RX Non-covered Stents


Order Number	Description	Diameter (mm)	Length (mm)	Catheter Diameter (F) / (mm)	Guidewire Diameter (in) / (mm)
M00570600	WallFlex RX Biliary	8	40	8 / 2.67	.035 / 0.89
M00570610	WallFlex RX Biliary	8	60	8 / 2.67	.035 / 0.89
M00570620	WallFlex RX Biliary	8	80	8 / 2.67	.035 / 0.89
M00570630	WallFlex RX Biliary	8	100	8 / 2.67	.035 / 0.89
M00576780	WallFlex RX Biliary	8	120	8 / 2.67	.035 / 0.89
M00570890	WallFlex RX Biliary	10	40	8 / 2.67	.035 / 0.89
M00570640	WallFlex RX Biliary	10	60	8 / 2.67	.035 / 0.89
M00570650	WallFlex RX Biliary	10	80	8 / 2.67	.035 / 0.89
M00570660	WallFlex RX Biliary	10	100	8 / 2.67	.035 / 0.89
M00576790	WallFlex RX Biliary	10	120	8 / 2.67	.035 / 0.89

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

Stents

RX BILIARY METAL STENTS (CONTINUED)



Wallstent™ RX Biliary Stents with Permalume™ Covering



Order Number	Description	Covered Length (mm)	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
M00569670	Wallstent RX w/Permalume	30	8	40	8 / 2.67	194
M00569680	Wallstent RX w/Permalume	50	8	60	8 / 2.67	194
M00569690	Wallstent RX w/Permalume	70	8	80	8 / 2.67	194
M00569700	Wallstent RX w/Permalume	30	10	40	8 / 2.67	194
M00569710	Wallstent RX w/Permalume	50	10	60	8 / 2.67	194
M00569720	Wallstent RX w/Permalume	70	10	80	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

Wallstent™ RX Biliary Stents




Order Number	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
M00569610	Wallstent RX Biliary	8	40	8 / 2.67	194
M00569620	Wallstent RX Biliary	8	60	8 / 2.67	194
M00569630	Wallstent RX Biliary	8	80	8 / 2.67	194
M00569800	Wallstent RX Biliary	8	100	8 / 2.67	194
M00569640	Wallstent RX Biliary	10	40	8 / 2.67	194
M00569650	Wallstent RX Biliary	10	60	8 / 2.67	194
M00569660	Wallstent RX Biliary	10	80	8 / 2.67	194
M00569810	Wallstent RX Biliary	10	100	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm or .035in (0.89mm) Dreamwire™ High Performance Guidewire, 260cm

Stents

BILIARY METAL STENTS



Wallstent™ Biliary Partially Covered Stents with Permalume™ Covering



Order Number	Description	Covered Length (mm)	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
H965430700	Wallstent Biliary w/Permalume	30	8	40	8 / 2.67	194
H965430800	Wallstent Biliary w/Permalume	50	8	60	8 / 2.67	194
H965430900	Wallstent Biliary w/Permalume	70	8	80	8 / 2.67	194
H965431000	Wallstent Biliary w/Permalume	30	10	40	8 / 2.67	194
H965431100	Wallstent Biliary w/Permalume	50	10	60	8 / 2.67	194
H965431200	Wallstent Biliary w/Permalume	70	10	80	8 / 2.67	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft or .035in (0.89mm) Dreamwire™ High Performance Guidewire

Wallstent™ Biliary Stents (Uncovered)




Order Number	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (F) / (mm)	Working Length (cm)
H965430100	Wallstent Biliary	8	40	7.5 / 2.5	194
H965430200	Wallstent Biliary	8	60	7.5 / 2.5	194
H965430300	Wallstent Biliary	8	80	7.5 / 2.5	194
H965430400	Wallstent Biliary	10	40	7.5 / 2.5	194
H965430500	Wallstent Biliary	10	60	7.5 / 2.5	194
H965430600	Wallstent Biliary	10	80	7.5 / 2.5	194

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft or .035in (0.89mm) Dreamwire™ High Performance Guidewire

Stents

RX BILIARY PLASTIC STENTS



Advanix™ Pre-loaded RX Stents with NaviFlex™ RX Delivery System



Order Number	Description	Shape	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00534200	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	5
M00534210	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	7
M00534220	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	9
M00534230	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	12
M00534240	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	15
M00534250	Advanix Pre-loaded RX Stent	Duodenal Bend	7 / 2.3	18
M00534260	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	5
M00534270	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	7
M00534280	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	9
M00534290	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	12
M00534300	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	15
M00534310	Advanix Pre-loaded RX Stent	Duodenal Bend	8.5 / 2.83	18
M00534320	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	5
M00534330	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	7
M00534340	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	9
M00534350	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	12
M00534360	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	15
M00534370	Advanix Pre-loaded RX Stent	Duodenal Bend	10 / 3.3	18
M00534560	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	5
M00534570	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	7
M00534580	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	9
M00534590	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	12
M00534600	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	15
M00534610	Advanix Pre-loaded RX Stent	Center Bend	7 / 2.3	18
M00534620	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	5
M00534630	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	7
M00534640	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	9
M00534650	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	12
M00534660	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	15
M00534670	Advanix Pre-loaded RX Stent	Center Bend	8.5 / 2.83	18
M00534680	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	5
M00534690	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	7
M00534700	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	9
M00534710	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	12
M00534720	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	15
M00534730	Advanix Pre-loaded RX Stent	Center Bend	10 / 3.3	18

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



Contents

Alphabetical Index

Product Code Index

Stents

RX BILIARY PLASTIC STENTS (CONTINUED)



Duodenal Bend



Center Bend



Double Pigtail



Advanix™ Stents – Singles

Order Number	Description	Shape	O.D. (F)/(mm)	Length (cm) Barb to Barb
M00532870	Advanix Stent	Duodenal Bend	7 / 2.3	5
M00532880	Advanix Stent	Duodenal Bend	7 / 2.3	7
M00532890	Advanix Stent	Duodenal Bend	7 / 2.3	9
M00532900	Advanix Stent	Duodenal Bend	7 / 2.3	12
M00532910	Advanix Stent	Duodenal Bend	7 / 2.3	15
M00532920	Advanix Stent	Duodenal Bend	7 / 2.3	18
M00532930	Advanix Stent	Duodenal Bend	8.5 / 2.83	5
M00532940	Advanix Stent	Duodenal Bend	8.5 / 2.83	7
M00532950	Advanix Stent	Duodenal Bend	8.5 / 2.83	9
M00532960	Advanix Stent	Duodenal Bend	8.5 / 2.83	12
M00532970	Advanix Stent	Duodenal Bend	8.5 / 2.83	15
M00532980	Advanix Stent	Duodenal Bend	8.5 / 2.83	18
M00532990	Advanix Stent	Duodenal Bend	10 / 3.3	5
M00533000	Advanix Stent	Duodenal Bend	10 / 3.3	7
M00533010	Advanix Stent	Duodenal Bend	10 / 3.3	9
M00533020	Advanix Stent	Duodenal Bend	10 / 3.3	12
M00533030	Advanix Stent	Duodenal Bend	10 / 3.3	15
M00533040	Advanix Stent	Duodenal Bend	10 / 3.3	18
M00533230	Advanix Stent	Center Bend	7 / 2.3	5
M00533240	Advanix Stent	Center Bend	7 / 2.3	7
M00533250	Advanix Stent	Center Bend	7 / 2.3	9
M00533260	Advanix Stent	Center Bend	7 / 2.3	12
M00533270	Advanix Stent	Center Bend	7 / 2.3	15
M00533280	Advanix Stent	Center Bend	7 / 2.3	18
M00533290	Advanix Stent	Center Bend	8.5 / 2.83	5
M00533300	Advanix Stent	Center Bend	8.5 / 2.83	7
M00533310	Advanix Stent	Center Bend	8.5 / 2.83	9
M00533320	Advanix Stent	Center Bend	8.5 / 2.83	12
M00533330	Advanix Stent	Center Bend	8.5 / 2.83	15
M00533340	Advanix Stent	Center Bend	8.5 / 2.83	18
M00533350	Advanix Stent	Center Bend	10 / 3.3	5
M00533360	Advanix Stent	Center Bend	10 / 3.3	7
M00533370	Advanix Stent	Center Bend	10 / 3.3	9
M00533380	Advanix Stent	Center Bend	10 / 3.3	12
M00533390	Advanix Stent	Center Bend	10 / 3.3	15
M00533400	Advanix Stent	Center Bend	10 / 3.3	18

Stents

RX BILIARY PLASTIC STENTS (CONTINUED)

Advanix™ Stents – Singles (CONTINUED)



Order Number	Description	Shape	O.D. (F) / (mm)	Length (cm) Between Pigtails
M00532160	Advanix Stent	Double Pigtail	7 / 2.3	3
M00532180	Advanix Stent	Double Pigtail	7 / 2.3	5
M00532190	Advanix Stent	Double Pigtail	7 / 2.3	7
M00532210	Advanix Stent	Double Pigtail	7 / 2.3	10
M00532220	Advanix Stent	Double Pigtail	7 / 2.3	12
M00532230	Advanix Stent	Double Pigtail	7 / 2.3	15
M00532240	Advanix Stent	Double Pigtail	10 / 3.3	3
M00532260	Advanix Stent	Double Pigtail	10 / 3.3	5
M00532270	Advanix Stent	Double Pigtail	10 / 3.3	7
M00532290	Advanix Stent	Double Pigtail	10 / 3.3	10
M00532300	Advanix Stent	Double Pigtail	10 / 3.3	12
M00532310	Advanix Stent	Double Pigtail	10 / 3.3	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Scope Compatibility: 7F (2.3mm) and 8.5F (2.83mm) stents require a 3.2mm working channel; 10F (3.3mm) stents require a 4.2mm working channel



NaviFlex™ RX Delivery System



Order Number	Description	Catheter Diameter (mm)	Working Length (cm)
M00533560	RX Delivery System	7 / 2.3	202.5
M00533570	RX Delivery System	8.5 / 2.83	202.5
M00533580	RX Delivery System	10 / 3.3	202.5

Packaged one per box

Stents

RX BILIARY PLASTIC STENTS (CONTINUED)



RX Biliary Plastic Stents

Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00545550	RX Biliary Plastic Stent	7 / 2.33	5
M00545560	RX Biliary Plastic Stent	7 / 2.33	7
M00545570	RX Biliary Plastic Stent	7 / 2.33	10
M00545580	RX Biliary Plastic Stent	7 / 2.33	12
M00545590	RX Biliary Plastic Stent	7 / 2.33	15
M00545650	RX Biliary Plastic Stent	8.5 / 2.83	5
M00545660	RX Biliary Plastic Stent	8.5 / 2.83	7
M00545670	RX Biliary Plastic Stent	8.5 / 2.83	10
M00545680	RX Biliary Plastic Stent	8.5 / 2.83	12
M00545690	RX Biliary Plastic Stent	8.5 / 2.83	15
M00545600	RX Biliary Plastic Stent	10 / 3.33	5
M00545610	RX Biliary Plastic Stent	10 / 3.33	7
M00545620	RX Biliary Plastic Stent	10 / 3.33	10
M00545630	RX Biliary Plastic Stent	10 / 3.33	12
M00545640	RX Biliary Plastic Stent	10 / 3.33	15
M00545700	RX Biliary Plastic Stent	11.5 / 3.83	5
M00545710	RX Biliary Plastic Stent	11.5 / 3.83	7
M00545720	RX Biliary Plastic Stent	11.5 / 3.83	10
M00545730	RX Biliary Plastic Stent	11.5 / 3.83	12
M00545740	RX Biliary Plastic Stent	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire



Stents

BILIARY PLASTIC STENTS



Flexima™ Biliary Stents



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00539200	Flexima Stent	7 / 2.33	5
M00539210	Flexima Stent	7 / 2.33	7
M00539220	Flexima Stent	7 / 2.33	10
M00539230	Flexima Stent	7 / 2.33	12
M00539240	Flexima Stent	7 / 2.33	15
M00539250	Flexima Stent	8.5 / 2.83	5
M00539260	Flexima Stent	8.5 / 2.83	7
M00539270	Flexima Stent	8.5 / 2.83	10
M00539280	Flexima Stent	8.5 / 2.83	12
M00539290	Flexima Stent	8.5 / 2.83	15
M00539300	Flexima Stent	10 / 3.33	5
M00539310	Flexima Stent	10 / 3.33	7
M00539320	Flexima Stent	10 / 3.33	10
M00539330	Flexima Stent	10 / 3.33	12
M00539340	Flexima Stent	10 / 3.33	15
M00539350	Flexima Stent	11.5 / 3.83	5
M00539360	Flexima Stent	11.5 / 3.83	7
M00539370	Flexima Stent	11.5 / 3.83	10
M00539380	Flexima Stent	11.5 / 3.83	12
M00539390	Flexima Stent	11.5 / 3.83	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Flexima™ Nasobiliary Stents



Order Number	Description	O.D. (F) / (mm)	Catheter Length (cm)	Recommended Guidewire (in) / (mm)
M00540120†	Flexima Nasal Stent	6 / 2.00	250	.035 / 0.89
M00540130†	Flexima Nasal Stent	7.5 / 2.50	250	.035 / 0.89
M00540140†	Flexima Nasal Stent	8.5 / 2.83	250	.035 / 0.89
M00540160†	Flexima Nasal Stent	6 / 2.00	250	Steel Teflon® coated wire included*
M00540170†	Flexima Nasal Stent	7.5 / 2.50	250	Steel Teflon® coated wire included*
M00540180†	Flexima Nasal Stent	8.5 / 2.83	250	Steel Teflon® coated wire included*

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

† Product available in the European Economic Area (EEA) only. Please check availability of the product with your local sales representative or your local customer service.

* Teflon® is a registered trademark of E. I. DuPont Nemours & Co. Inc.



Stents

BILIARY PLASTIC STENTS (CONTINUED)



Percuflex™ Amsterdam Biliary Stents



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb
M00533610	Percuflex Amsterdam Stent	7 / 2.33	5
M00533620	Percuflex Amsterdam Stent	7 / 2.33	7
M00533630	Percuflex Amsterdam Stent	7 / 2.33	10
M00533640	Percuflex Amsterdam Stent	7 / 2.33	12
M00533650	Percuflex Amsterdam Stent	7 / 2.33	15
M00533660	Percuflex Amsterdam Stent	10 / 3.33	5
M00533670	Percuflex Amsterdam Stent	10 / 3.33	7
M00533680	Percuflex Amsterdam Stent	10 / 3.33	10
M00533690	Percuflex Amsterdam Stent	10 / 3.33	12
M00533700	Percuflex Amsterdam Stent	10 / 3.33	15

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

Percuflex™ Amsterdam Introducer Kits*



Order Number	Description
M00533910	Percuflex Amsterdam Stent Kit w/push catheter for use with 7F (2.33mm) O.D. stent
M00533920	Percuflex Amsterdam Stent Kit w/push and guide catheter for use with 10F (3.33mm) O.D. stent

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* Stents sold separately.



C-flex™ Pigtail Biliary Stents*



Order Number	Description	O.D. (F) / (mm)	Length (cm) Barb to Barb	Colour Code
M00532010	C-flex Pigtail Stent	7 / 2.33	5	Blue
M00532020	C-flex Pigtail Stent	7 / 2.33	10	Blue
M00532030	C-flex Pigtail Stent	7 / 2.33	15	Blue
M00532040	C-flex Pigtail Stent	10 / 3.33	5	Grey
M00532050	C-flex Pigtail Stent	10 / 3.33	10	Grey
M00532060	C-flex Pigtail Stent	10 / 3.33	15	Grey


Packaged one per box

Recommended Guidewire: .035in (0.89mm) Dreamwire™ High Performance Guidewire

* For use with Percuflex™ Amsterdam Single-use Introducer Kit.

Stents

PANCREATIC PLASTIC STENTS



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Pigtail



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537350	Advanix Pancreatic Stent Kit ●	No	Pigtail	4 / 1.35	3
M00537360	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	4
M00537370	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	5
M00537380	Advanix Pancreatic Stent Kit ▲	No	Pigtail	4 / 1.35	7
M00537390	Advanix Pancreatic Stent Kit ●	No	Pigtail	5 / 1.67	3
M00537400	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	4
M00537410	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	5
M00537420	Advanix Pancreatic Stent Kit ▲	No	Pigtail	5 / 1.67	7
M00537430	Advanix Pancreatic Stent Kit ●	Yes	Pigtail	5 / 1.67	3
M00537440	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	4
M00537450	Advanix Pancreatic Stent Kit ▲	Yes	Pigtail	5 / 1.67	5



Advanix™ Pancreatic Stent Kits (Stent + RX Pusher) – Straight



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537460	Advanix Pancreatic Stent Kit ● ¹	No	Straight	4 / 1.35	3
M00537470	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	4
M00537480	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	5
M00537490	Advanix Pancreatic Stent Kit ▲ ¹	No	Straight	4 / 1.35	7
M00537500	Advanix Pancreatic Stent Kit ●	No	Straight	5 / 1.67	3
M00537510	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	4
M00537520	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	5
M00537530	Advanix Pancreatic Stent Kit ▲	No	Straight	5 / 1.67	7
M00537540	Advanix Pancreatic Stent Kit ●	Yes	Straight	5 / 1.67	3
M00537550	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	4
M00537560	Advanix Pancreatic Stent Kit ▲	Yes	Straight	5 / 1.67	5

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Stents

PANCREATIC PLASTIC STENTS (CONTINUED)



Advanix™ Pancreatic Single Stents – Pigtail

Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536740	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	2
M00536750	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	3
M00536760	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	4
M00536770	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	5
M00536780	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	6
M00536790	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	7
M00536800	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	8
M00536810	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	9
M00536820	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	10
M00536830	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	11
M00536840	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	12
M00536850	Advanix Pancreatic Stent ■	No	Pigtail	3 / 1.00	13
M00536860	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	2
M00536870	Advanix Pancreatic Stent ●	No	Pigtail	4 / 1.35	3
M00536880	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	4
M00536890	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	5
M00536900	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	6
M00536910	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	7
M00536920	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	8
M00536930	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	9
M00536940	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	10
M00536950	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	11
M00536960	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	12
M00536970	Advanix Pancreatic Stent ▲	No	Pigtail	4 / 1.35	13
M00536980	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	2
M00536990	Advanix Pancreatic Stent ●	No	Pigtail	5 / 1.67	3
M00537000	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	4
M00537010	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	5
M00537020	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	6
M00537030	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	7
M00537040	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	8
M00537050	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	9
M00537060	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	10
M00537070	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	11
M00537080	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	12
M00537090	Advanix Pancreatic Stent ▲	No	Pigtail	5 / 1.67	13

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

1 Stiffer stent material compared to 4F pigtail stents.

2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.




Stents

PANCREATIC PLASTIC STENTS (CONTINUED)

Advanix™ Pancreatic Single Stents – Pigtail (CONTINUED)



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00537100	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	2
M00537110	Advanix Pancreatic Stent ●	Yes	Pigtail	5 / 1.67	3
M00537120	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	4
M00537130	Advanix Pancreatic Stent ▲	Yes	Pigtail	5 / 1.67	5
M00537140	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	6
M00537150	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	7
M00537160	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	8
M00537170	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	9
M00537180	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	10
M00537190	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	11
M00537200	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	12
M00537210	Advanix Pancreatic Stent ▲ ²	Yes	Pigtail	5 / 1.67	13
M00537220	Advanix Pancreatic Stent ■	Yes	Pigtail	7 / 2.33	3
M00537230	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	4
M00537240	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	5
M00537250	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	6
M00537260	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	7
M00537270	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	8
M00537280	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	9
M00537290	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	10
M00537300	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	11
M00537310	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	12
M00537320	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	13
M00537330	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	15
M00537340	Advanix Pancreatic Stent ▲	Yes	Pigtail	7 / 2.33	18



Advanix™ Pancreatic Single Stents – Straight



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536000	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	2
M00536010	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	3
M00536020	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	4
M00536030	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	5
M00536040	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	6
M00536050	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	7
M00536060	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	8
M00536070	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	9

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

1 Stiffer stent material compared to 4F pigtail stents.

2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Stents

PANCREATIC PLASTIC STENTS (CONTINUED)

Advanix™ Pancreatic Single Stents – Straight (CONTINUED)



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536080	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	10
M00536090	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	11
M00536100	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	12
M00536110	Advanix Pancreatic Stent ■	No	Straight	3 / 1.00	13
M00536120	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	2
M00536130	Advanix Pancreatic Stent ● ¹	No	Straight	4 / 1.35	3
M00536140	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	4
M00536150	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	5
M00536160	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	6
M00536170	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	7
M00536180	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	8
M00536190	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	9
M00536200	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	10
M00536210	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	11
M00536220	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	12
M00536230	Advanix Pancreatic Stent ▲ ¹	No	Straight	4 / 1.35	13
M00536240	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	2
M00536250	Advanix Pancreatic Stent ●	No	Straight	5 / 1.67	3
M00536260	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	4
M00536270	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	5
M00536280	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	6
M00536290	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	7
M00536300	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	8
M00536310	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	9
M00536320	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	10
M00536330	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	11
M00536340	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	12
M00536350	Advanix Pancreatic Stent ▲	No	Straight	5 / 1.67	13
M00536360	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	2
M00536370	Advanix Pancreatic Stent ●	Yes	Straight	5 / 1.67	3
M00536380	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	4
M00536390	Advanix Pancreatic Stent ▲	Yes	Straight	5 / 1.67	5
M00536400	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	6
M00536410	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	7
M00536420	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	8
M00536430	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	9
M00536440	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	10

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker

¹ Stiffer stent material compared to 4F pigtail stents.

² Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.



Stents

PANCREATIC PLASTIC STENTS (CONTINUED)

Advanix™ Pancreatic Single Stents – Straight (CONTINUED)



Order Number	Description	Leading Barb	Shape	O.D. (F) / (mm)	Length (cm)
M00536450	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	11
M00536460	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	12
M00536470	Advanix Pancreatic Stent ▲ ²	Yes	Straight	5 / 1.67	13
M00536480	Advanix Pancreatic Stent ■	Yes	Straight	7 / 2.33	3
M00536490	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	4
M00536500	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	5
M00536510	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	6
M00536520	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	7
M00536530	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	8
M00536540	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	9
M00536550	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	10
M00536560	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	11
M00536570	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	12
M00536580	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	13
M00536590	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	15
M00536600	Advanix Pancreatic Stent ▲	Yes	Straight	7 / 2.33	18
M00536610	Advanix Pancreatic Stent ■	Yes	Straight	10 / 3.33	3
M00536620	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	4
M00536630	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	5
M00536640	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	6
M00536650	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	7
M00536660	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	8
M00536670	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	9
M00536680	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	10
M00536690	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	11
M00536700	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	12
M00536710	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	13
M00536720	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	15
M00536730	Advanix Pancreatic Stent ▲	Yes	Straight	10 / 3.33	18

■ Trailing endoscopic marker ● Trailing radiopaque marker ▲ Leading radiopaque marker with trailing endoscopic marker
 1 Stiffer stent material compared to 4F pigtail stents.
 2 Stiffer stent material compared to 5F pigtail and 5F straight stents with No Leading Barb.

Stents

PANCREATIC PLASTIC STENTS (CONTINUED)





NaviFlex™ RX Pancreatic Delivery System & NaviFlex Pushers



Order Number	Description	Diameter (F) / (mm)
M00535000	Pancreatic Stent RX Pusher	3 / 1.00
M00535010	Pancreatic Stent RX Pusher	4 & 5 / 1.35 & 1.67
M00535020	Pancreatic Stent Pusher	3 / 1.00
M00535030	Pancreatic Stent Pusher	4 & 5 / 1.35 & 1.67
M00535060	NaviFlex RX Pancreatic Stent Delivery System	7 / 2.33
M00535080	NaviFlex RX Pancreatic Stent Delivery System	10 / 3.33

Stents

ESOPHAGEAL STENTS



WallFlex™ Fully Covered RMV Esophageal Stents with Permalume™ Covering (including refractory benign indication)



Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (F) / (mm)	System Length (cm)
M00516210†	WallFlex Fully Covered RMV Esophageal Stent	18	25 / 23	10.3	18.5 / 6.2	120
M00516220†	WallFlex Fully Covered RMV Esophageal Stent	18	25 / 23	12.3	18.5 / 6.2	120
M00516230†	WallFlex Fully Covered RMV Esophageal Stent	18	25 / 23	15.3	18.5 / 6.2	120
M00516240†	WallFlex Fully Covered RMV Esophageal Stent	23	28 / 28	10.5	18.5 / 6.2	120
M00516250†	WallFlex Fully Covered RMV Esophageal Stent	23	28 / 28	12.5	18.5 / 6.2	120
M00516260†	WallFlex Fully Covered RMV Esophageal Stent	23	28 / 28	15.5	18.5 / 6.2	120



† Product commercialised only in countries where benign indication is approved.

WallFlex™ Fully Covered RMV Esophageal Stents Longer Loop with Permalume™ Covering (including refractory benign indication)



Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (F) / (mm)	System Length (cm)
M00516270	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	18	25 / 23	10.3	18.5 / 6.2	120
M00516280	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	18	25 / 23	12.3	18.5 / 6.2	120
M00516290	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	18	25 / 23	15.3	18.5 / 6.2	120
M00516300	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	23	28 / 28	10.5	18.5 / 6.2	120
M00516310	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	23	28 / 28	12.5	18.5 / 6.2	120
M00516320	WallFlex Fully Covered RMV Esophageal Stent Longer Loop	23	28 / 28	15.5	18.5 / 6.2	120

Packaged one per box
 Limited availability in all European countries (except UK and France). Please contact your local Territory Manager.



WallFlex™ Partially Covered Esophageal Stents with Permalume™ Covering




Order Number	Description	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (F) / (mm)	System Length (cm)
M00516900	WallFlex Partially Covered Esophageal Stent	18	23 / 23	10.3	7	18.5 / 6.2	120
M00516910	WallFlex Partially Covered Esophageal Stent	18	23 / 23	12.3	9	18.5 / 6.2	120
M00516920	WallFlex Partially Covered Esophageal Stent	18	23 / 23	15.3	12	18.5 / 6.2	120
M00516930	WallFlex Partially Covered Esophageal Stent	23	28 / 28	10.5	7	18.5 / 6.2	120
M00516940	WallFlex Partially Covered Esophageal Stent	23	28 / 28	12.5	9	18.5 / 6.2	120
M00516950	WallFlex Partially Covered Esophageal Stent	23	28 / 28	15.5	12	18.5 / 6.2	120

Packaged one per box
Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

Stents

ESOPHAGEAL STENTS (CONTINUED)



Ultraflex™ Covered Esophageal NG Stents

Proximal Release





Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)
M00513840	Ultraflex Esophageal NG Stent	18	23	10	7
M00513850	Ultraflex Esophageal NG Stent	18	23	12	9
M00513860	Ultraflex Esophageal NG Stent	18	23	15	12
M00514240	Ultraflex Esophageal NG Stent	23	28	10	7
M00514250	Ultraflex Esophageal NG Stent	23	28	12	9

Distal Release

Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)
M00513730	Ultraflex Esophageal NG Stent	18	23	10	7
M00513740	Ultraflex Esophageal NG Stent	18	23	12	9
M00513750	Ultraflex Esophageal NG Stent	18	23	15	12
M00514200	Ultraflex Esophageal NG Stent	23	28	10	7
M00514210	Ultraflex Esophageal NG Stent	23	28	12	9

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621



Ultraflex™ Non-covered Esophageal NG Stents

Proximal Release



Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)
M00513800	Ultraflex Esophageal NG Stent	18	23	7
M00513810	Ultraflex Esophageal NG Stent	18	23	10
M00513820	Ultraflex Esophageal NG Stent	18	23	12
M00513830	Ultraflex Esophageal NG Stent	18	23	15
M00514220	Ultraflex Esophageal NG Stent	23	28	7
M00514230	Ultraflex Esophageal NG Stent	23	28	12

Distal Release



Order Number	Description	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)
M00513700	Ultraflex Esophageal NG Stent	18	23	7
M00513710	Ultraflex Esophageal NG Stent	18	23	10
M00513720	Ultraflex Esophageal NG Stent	18	28	15

Packaged one per box

Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

Stents

ESOPHAGEAL STENTS (CONTINUED)



Polyflex Stent™ Self-expanding Silicone Esophageal Stents




Order Number	Description	Stent Body I.D. (mm)	Stent Flare I.D. (mm)	Stent Length (cm)	System O.D. (mm)
M00514270	Polyflex Esophageal Stent	16	20	9	12
M00514280	Polyflex Esophageal Stent	16	20	12	12
M00514290	Polyflex Esophageal Stent	16	20	15	12
M00514300	Polyflex Esophageal Stent	18	23	9	13
M00514310	Polyflex Esophageal Stent	18	23	12	13
M00514320	Polyflex Esophageal Stent	18	23	15	13
M00514330	Polyflex Esophageal Stent	21	25	9	14
M00514340	Polyflex Esophageal Stent	21	25	12	14
M00514350	Polyflex Esophageal Stent	21	25	15	14

Packaged one per box
Recommended Guidewire: .038in (0.97mm) 260cm Jagwire™ Guidewire, M00556621

Stents

COLONIC AND DUODENAL STENTS



WallFlex™ Colonic Stents





Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565100	WallFlex Colonic Stent	22	27	6	230	10 / 3.33
M00565110	WallFlex Colonic Stent	22	27	9	230	10 / 3.33
M00565120	WallFlex Colonic Stent	22	27	12	230	10 / 3.33
M00565130	WallFlex Colonic Stent	22	27	6	135	10 / 3.33
M00565140	WallFlex Colonic Stent	22	27	9	135	10 / 3.33
M00565150	WallFlex Colonic Stent	22	27	12	135	10 / 3.33
M00565040	WallFlex Colonic Stent	25	30	6	230	10 / 3.33
M00565050	WallFlex Colonic Stent	25	30	9	230	10 / 3.33
M00565060	WallFlex Colonic Stent	25	30	12	230	10 / 3.33
M00565070	WallFlex Colonic Stent	25	30	6	135	10 / 3.33
M00565080	WallFlex Colonic Stent	25	30	9	135	10 / 3.33
M00565090	WallFlex Colonic Stent	25	30	12	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H965180010

Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip



WallFlex™ Duodenal Stents




Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565010	WallFlex Duodenal Stent	22	27	6	230	10 / 3.33
M00565020	WallFlex Duodenal Stent	22	27	9	230	10 / 3.33
M00565030	WallFlex Duodenal Stent	22	27	12	230	10 / 3.33

Packaged one per box

Recommended Guidewires: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H965180010

Stents

COLONIC AND DUODENAL STENTS (CONTINUED)



Wallstent™ Colonic and Duodenal Stents





Order Number	Description	Expanded Stent Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Diameter (F) / (mm)
M00565560	Wallstent Colonic and Duodenal Stent	20	6	230	10 / 3.33
M00565570	Wallstent Colonic and Duodenal Stent	20	9	230	10 / 3.33
M00565620	Wallstent Colonic and Duodenal Stent	20	6	135	10 / 3.33
M00565630	Wallstent Colonic and Duodenal Stent	20	9	135	10 / 3.33
M00565580	Wallstent Colonic and Duodenal Stent	22	6	230	10 / 3.33
M00565590	Wallstent Colonic and Duodenal Stent	22	9	230	10 / 3.33
M00565640	Wallstent Colonic and Duodenal Stent	22	6	135	10 / 3.33
M00565650	Wallstent Colonic and Duodenal Stent	22	9	135	10 / 3.33

Packaged one per box

Recommended Guidewire Systems 230cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 450cm: M00556161 Straight Dream Tip™, M00556171 Angled Dream Tip™ or .035 in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 450cm: M00556601 Straight Tip, M00556611 Angled Tip or .035in (0.89mm) Super Stiff Wire, 500cm: H965180010

Recommended Guidewire Systems 135cm: .035in (0.89mm) Dreamwire™ High Performance Guidewire Stiff Shaft, 260cm: M00556121 Straight Tip, M00556131 Angled Tip or .035in (0.89mm) Jagwire™ Guidewire Stiff Shaft, 260cm: M00557661 Straight Tip



Ultraflex Precision™ Colonic Stents




Order Number	Description	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (cm)	Working Length (cm)	Catheter Tip Diameter (F) / (mm)
M00557330	Ultraflex Precision Colonic Stent	25	30	5.7	100	22 / 7.33
M00557360	Ultraflex Precision Colonic Stent	25	30	8.7	100	22 / 7.33
M00557380	Ultraflex Precision Colonic Stent	25	30	11.7	100	22 / 7.33

Packaged one per box

Recommended Guidewire: .038in (0.97mm) Jagwire™ Guidewire, Order No. M00556621




INDICATED FOR USE IN THE
DASE
PROCEDURE
· CRE™ PRO WIREGUIDED ·
BALLOON DILATOR ·



Dilation

CRE™ CONTROLLED RADIAL EXPANSION BALLOON DILATORS



CRE™ PRO Wireguided Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)




Order Number	Balloon Inflated O.D.		Balloon Length	Inflation Pressure		Minimum Working Channel	Catheter Diameter	Working Length	Units
	(mm)	(F)	(cm)	(ATM)	(kPa)	(mm)*	(F) / (mm)	(cm)	
M00558600	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2	7.5 / 2.50	180	box 1
M00558610	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8/3.2	7.5 / 2.50	180	box 1
M00558620	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2	7.5 / 2.50	180	box 1
M00558630	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2	7.5 / 2.50	180	box 1
M00558640	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7	7.5 / 2.50	180	box 1
M00558650**	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7	7.5 / 2.50	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 Guidewire included and preloaded in balloon catheter.

* Duodenoscope Channel Size

** Esophageal and biliary applications only

This product contains no detectable latex.



CRE™ PRO Wireguided Colonic/Esophageal/Pyloric/Biliary Balloon Dilators (including DASE indication)




Order Number	Balloon Inflated O.D.		Balloon Length	Inflation Pressure		Minimum Working Channel	Catheter Diameter	Working Length	Units
	(mm)	(F)	(cm)	(ATM)	(kPa)	(mm)*	(F) / (mm)	(cm)	
M00558660	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8/3.2	7.5 / 2.50	240	box 1
M00558670	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8/3.2	7.5 / 2.50	240	box 1
M00558680	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8/3.2	7.5 / 2.50	240	box 1
M00558690	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8/3.2	7.5 / 2.50	240	box 1
M00558700	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8/3.7	7.5 / 2.50	240	box 1
M00558710**	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8/3.7	7.5 / 2.50	240	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 Guidewire included and preloaded in balloon catheter.

* Duodenoscope Channel Size


** Esophageal, colonic and biliary applications only

This product contains no detectable latex.



Dilation

CRE™ CONTROLLED RADIAL EXPANSION BALLOON DILATORS (CONTINUED)




CRE™ Fixed Wire Balloon Dilators




Order Number	Inflated O.D. (mm)	Balloon Length (cm)	Inflation Pressure (ATM)	Inflation Pressure (kPa)	Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
M00558330	6-7-8	8	3-6-10	304-608-1013	2.8	6 / 2.00	180	box 1
M00558340	8-9-10	8	3-5.5-9	304-557-912	2.8	6 / 2.00	180	box 1
M00558350	10-11-12	8	3-5-8	304-507-810	2.8	6 / 2.00	180	box 1
M00558360	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 1
M00558361	12-13.5-15	8	3-4.5-8	304-456-810	2.8	6 / 2.00	180	box 5
M00558370	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 1
M00558371	15-16.5-18	8	3-4.5-7	304-456-709	2.8	6 / 2.00	180	box 5
M00558380	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 1
M00558381	18-19-20	8	3-4.5-6	304-456-608	2.8	6 / 2.00	180	box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 This product contains no detectable latex.



Dilation

MAX FORCE™ TTS BALLOON DILATORS



Max Force™ TTS Esophageal Balloon Dilators

Order Number	Inflated O.D. (mm)	Balloon Length (cm)	Inflation Pressure		Minimum Working Channel (mm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Units
			(ATM)	(kPa)				
M00589060	6	6	10	1013	2.8	6 / 2.00	180	box 1
M00589080	8	6	10	1013	2.8	6 / 2.00	180	box 1
M00589100	10	6	7	709	2.8	6 / 2.00	180	box 1
M00589120	12	6	7	709	2.8	6 / 2.00	180	box 1
M00589140	14	6	7	709	2.8	6 / 2.00	180	box 1
M00589150	15	6	5	507	2.8	6 / 2.00	180	box 1
M00589160	16	6	5	507	2.8	6 / 2.00	180	box 1
M00589180	18	6	5	507	2.8	6 / 2.00	180	box 1

Recommended Inflation Device: Alliance™ II Integrated Inflation System Re-usable Handle M00550620, Sterile Syringe M00550601
 This product contains no detectable latex.


Rigiflex™ II Achalasia Balloon Dilators*

Order Number	Inflated O.D. (mm)	Balloon Length (cm)	Catheter Diameter (F) / (mm)	Working Length (cm)	Recommended Guidewire (in) / (mm)**	Units
M00554500	30	10	14 / 4.67	90	.038 / 0.97	box 1
M00554510	35	10	14 / 4.67	90	.038 / 0.97	box 1
M00554520	40	10	14 / 4.67	90	.038 / 0.97	box 1

* Product shipped sterile. Recommended use with Achalasia Pneumatic Hand Pump and Monitor, M00553200


** Guidewire included with Rigiflex™ II Balloon Dilation Catheter

This product contains no detectable latex.



Dilation

INFLATION SYSTEMS



Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance™ II Inflation Handle	box 1
M00550601	⊗ 60ml Syringe/Gauge Assembly	box 5


DILATION ACCESSORIES

Inflation/Deflation Devices and Accessories

Order Number	Description	Units
M00553200	Achalasia Pneumatic Inflator and Monitor	box 1


RX Biliary Balloon Dilators are located in the RX Biliary Section on page 8

Biliary Balloon Dilators are located in the Biliary Section on page 17



Guidewires

HIGH PERFORMANCE



Dreamwire™ High Performance Guidewires



Order Number	Description	O.D. (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556101	Dreamwire Standard	.035 / 0.89	260	Straight	box 2
M00556111	Dreamwire Standard	.035 / 0.89	260	Angled	box 2
M00556121	Dreamwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556131	Dreamwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556141	Dreamwire Standard	.035 / 0.89	450	Straight	box 2
M00556151	Dreamwire Standard	.035 / 0.89	450	Angled	box 2
M00556161	Dreamwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556171	Dreamwire Stiff Shaft	.035 / 0.89	450	Angled	box 2

5cm Jagwire™ Guidewire Tip

10cm Straight Dream Tip™

10cm Angled Dream Tip™



Hydra Jagwire™ High Performance Guidewires




Order Number	Description	O.D. (in) / (mm)	Length (cm)	Dream Tip Style	Units
M00556001	Hydra Jagwire Standard	.035 / 0.89	260	Straight	box 2
M00556011	Hydra Jagwire Standard	.035 / 0.89	260	Angled	box 2
M00556021	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Straight	box 2
M00556031	Hydra Jagwire Stiff Shaft	.035 / 0.89	260	Angled	box 2
M00556041	Hydra Jagwire Standard	.035 / 0.89	450	Straight	box 2
M00556051	Hydra Jagwire Standard	.035 / 0.89	450	Angled	box 2
M00556061	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Straight	box 2
M00556071	Hydra Jagwire Stiff Shaft	.035 / 0.89	450	Angled	box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

Guidewires

HIGH PERFORMANCE (CONTINUED)



Jagwire™ High Performance Guidewires



Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556481	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	260	Straight	box 2
M00556491	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	260	Angled	box 2
M00556461	Jagwire Standard	.025-.018 / 0.64-0.46	260	Straight	box 2
M00557641	Jagwire Standard	.035-.027 / 0.89-0.69	260	Straight	box 2
M00557651	Jagwire Standard	.035-.027 / 0.89-0.69	260	Angled	box 2
M00557661	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Straight	box 2
M00557671	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Angled	box 2
M00556561	Jagwire Standard	.025-.018 / 0.64-0.46	450	Straight	box 2
M00556511	Jagwire Standard	.025-.018 / 0.64-0.46	450	Angled	box 2
M00556521	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	450	Straight	box 2
M00556531	Jagwire Stiff Shaft	.025-.019 / 0.64-0.48	450	Angled	box 2
M00556581	Jagwire Standard	.035-.027 / 0.89-0.69	450	Straight	box 2
M00556591	Jagwire Standard	.035-.027 / 0.89-0.69	450	Angled	box 2
M00556601	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	450	Straight	box 2
M00556611	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	450	Angled	box 2
M00556621	Jagwire Stiff Shaft	.038 / 0.97	260	Straight	box 2

Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.

Jagwire™ High Performance Guidewires – Extendable Jagwire




Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556902	Jagwire Jag Tail	.035 / 0.89	200	Guidewire Extension	box 5
M00556841	Jagwire Standard	.035-.027 / 0.89-0.69	260	Straight	box 2
M00556851	Jagwire Standard	.035-.027 / 0.89-0.69	260	Angled	box 2
M00556861	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Straight	box 2
M00556871	Jagwire Stiff Shaft	.035-.027 / 0.89-0.69	260	Angled	box 2

Jagwire Guidewire may remain in place during sphincterotomy. See "Instructions for Use" for details.


Guidewires

HIGH PERFORMANCE (CONTINUED)



NovaGold™ High Performance Guidewires

Order Number	Description	Guidewire Diameter (in) / (mm)	Length (cm)	Tip Style	Units
M00552001	NovaGold	.018 / 0.46	260	Straight (Shapeable)	box 2
M00552011	NovaGold	.018 / 0.46	480	Straight (Shapeable)	box 2



NaviPro™ 100% Hydrophilic Guidewires

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00556191	NaviPro Standard	.018 / 0.46	260	Angled	box 5
M00556181	NaviPro Standard	.018 / 0.46	260	Straight	box 5
M00556211	NaviPro Standard	.025 / 0.64	260	Angled	box 5
M00556201	NaviPro Standard	.025 / 0.64	260	Straight	box 5
M00556231	NaviPro Standard	.035 / 0.89	260	Angled	box 5
M00556221	NaviPro Standard	.035 / 0.89	260	Straight	box 5
M00556251	NaviPro Stiff Shaft	.035 / 0.89	260	Angled	box 5
M00556241	NaviPro Stiff Shaft	.035 / 0.89	260	Straight	box 5

STANDARD

Amplatz™ Super Stiff Guidewire

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
M00550090	Amplatz Super Stiff Guidewire	.038 / 0.97	260	Straight	box 1

Wallstent™ Super Stiff Guidewire

Order Number	Description	O.D. (in) / (mm)	Length (cm)	Tip Style	Units
H965180011	Wallstent Super Stiff Guidewire	.035 / 0.89	500	Tapered	box 5

RESOLUTION™ CLIP
• PROPHYLACTIC CLIPPING TO REDUCE THE RISK OF DELAYED BLEEDING POST LESION RESECTION •
NEW INDICATION



Contents




Alphabetical Index

Product Code Index

Hemostasis

CLIPPING DEVICES



NEW Resolution 360™ Clip



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00521230	Resolution 360 Clip	235	2.8	11	box 1
M00521231	Resolution 360 Clip	235	2.8	11	box 10
M00521232†	Resolution 360 Clip	235	2.8	11	box 20
M00521233†	Resolution 360 Clip	235	2.8	11	box 40
M00521240	Resolution 360 Clip	155	2.8	11	box 1
M00521241	Resolution 360 Clip	155	2.8	11	box 10
M00521242	Resolution 360 Clip	155	2.8	11	box 20
M00521243	Resolution 360 Clip	155	2.8	11	box 40

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.
 † Please check product availability in your country with your sales representative




Resolution™ Clip



Order Number	Description	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Units
M00522600	Resolution Clip	155	2.8	11	box 1
M00522601	Resolution Clip	155	2.8	11	box 10
M00522602	Resolution Clip	155	2.8	11	box 20
M00522610	Resolution Clip	235	2.8	11	box 1
M00522611	Resolution Clip	235	2.8	11	box 10
M00522612	Resolution Clip	235	2.8	11	box 20

* MR Conditional. Non-clinical testing has demonstrated the Resolution 360™ Clip is Magnetic Resonance (MR) Conditional according to ASTM F2503. Refer to the directions for use for MR conditions of use.

MULTIPLE BAND LIGATORS




Speedband Superview Super 7™ Multiple Band Ligators

Order Number	Description	Minimum Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Units
M00542250	Super 7 Band Ligator	2.8	8.6-11.5	7	box 1
M00542251	Super 7 Band Ligator	2.8	8.6-11.5	7	box 2
M00542253	Super 7 Band Ligator	2.8	8.6-11.5	7	box 4

Irrigation catheter supplied in each kit

Hemostasis


BIPOLAR ELECTROHEMOSTASIS CATHETERS



Injection Gold Probe™ Electrohemostasis/Injection Catheters

Order Number	Description	O.D. (F) / (mm)	Gauge	Needle O.D. (mm)	I.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560150	Injection Gold Probe Catheter	7 / 2.33	25	0.51	0.24	2.8	210	box 1
M00560160	Injection Gold Probe Catheter	10 / 3.33	25	0.51	0.24	3.7	210	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors




Gold Probe™ Electrohemostasis Catheters

Order Number	Description	O.D. (F) / (mm)	Minimum Working Channel (mm)	Working Length (cm)	Units
M00560070	Gold Probe Catheter	7 / 2.33	2.8	300	box 1
M00560071	Gold Probe Catheter	7 / 2.33	2.8	300	box 5
M00560100	Gold Probe Catheter	10 / 3.33	3.7	300	box 1
M00560101	Gold Probe Catheter	10 / 3.33	3.7	300	box 5
M00560220	Gold Probe Catheter	7 / 2.33	2.8	350	box 1

Type Generator/Connector: Please refer to Electrosurgical Bipolar Cable Adaptors

Hemostasis

NEEDLES




Interject™ Sclerotherapy Needles*



Order Number	Description	Sheath Design	Sheath O.D. (mm)	Gauge	Needle O.D. (mm)	I.D. (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518151	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	4	200	box 5
M00518161	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	4	200	box 5
M00518251	Interject Sclerotherapy Needle	Contrast	1.8	23	0.64	0.32	6	200	box 5
M00518261	Interject Sclerotherapy Needle	Contrast	1.8	25	0.51	0.24	6	200	box 5
M00518351	Interject Sclerotherapy Needle	Contrast	2.3	23	0.64	0.32	4	240	box 5
M00518361	Interject Sclerotherapy Needle	Contrast	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology



Interject™ Clear Sclerotherapy Needles*




Order Number	Description	Sheath Design	Sheath O.D. (mm)	Gauge	Needle O.D. (mm)	I.D. (mm)	Maximum Needle Extension Length (mm)	Working Length (cm)	Units
M00518111	Interject Sclerotherapy Needle	Clear	1.8	25	0.51	0.24	4	200	box 5
M00518301	Interject Sclerotherapy Needle	Clear	2.3	23	0.64	0.32	4	240	box 5
M00518311	Interject Sclerotherapy Needle	Clear	2.3	25	0.51	0.24	4	240	box 5

* With Star Catheter Technology

Hemostasis

NEEDLES (CONTINUED)




Optiflo™ Hemostasis Catheters



Order Number	Description	Needle			Sheath O.D. (mm)	Working Length (cm)	Working Channel (mm)	Units
		Gauge	O.D. (mm)	I.D. (mm)				
M00511731	Optiflo Hemostasis Catheter	23	0.64	0.32	2.3	200	2.8	box 5
M00511751	Optiflo Hemostasis Catheter	25	0.51	0.24	2.3	200	2.8	box 5

ELECTROSURGICAL BIPOLAR CABLE ADAPTORS



Electrosurgical Bipolar Cable Adaptors

Order Number	Description	Fix Pin Spacing (mm)	Generator Compatibility*	Units
M00561440	Bipolar Cable Adaptor	18	EMC (Coagulasem Dolley)	box 1
M00561450	Bipolar Cable Adaptor	22	ERBE (VIO International, VIO 100C, ICC Series, VIO 200D), EMC Dolley, Olympus (ESG100, PSD60)**	box 1
M00561460	Bipolar Cable Adaptor	25	Lamidey (Optima 4)	box 1
M00561470	Bipolar Cable Adaptor	28	Valleylab®, Lamidey (Duo), Bovie (Aaron 3250™, Aaron 1250™), Olympus (UES10, UES30, Surgemaster, ESG100™, PSD60™), Bowa, Soring, ConMed™ Sabre Genesys™	box 1
M00561480	Bipolar Cable Adaptor	28 with Third Pin	ConMed® (System 5000), Bovie (Aaron 2250™), Ethicon®	box 1
M00561490	Bipolar Cable Adaptor	32 with Third Pin	ConMed (Sabre 180, Excalibur, System 7550)	box 1
M00561500	Bipolar Cable Adaptor	39	Codman, Malis®	box 1

Bipolar Cable Adaptors are used to connect Gold Probe™ and Injection Gold Probe™ with the above mentioned generators.

* For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative

** Some variations in faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering

ENTIRE POLYPECTOMY SNARE PORTFOLIO
INDICATED FOR BOTH
COLD & HOT
POLYPECTOMY
• BOSTON SCIENTIFIC •



Tissue Resection

EMR DEVICES FOR UPPER GI




Captivator™ EMR



Order Number	Description	Units
M00561600	Captivator EMR for Standard Gastroscopes	box 1
M00561610	Captivator EMR for Large Gastroscopes	box 1

SNARES



Rotatable Snares*




Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	box 10

* With Star Catheter Technology.


Tissue Resection

SNARES (CONTINUED)



Sensation™ Short Throw Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	box 5





Captiflex™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	box 40


Tissue Resection

SNARES (CONTINUED)



Captivator™ Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	box 10




Captivator™ II Snares

Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00561221	10 mm Rounded	Stiff	10	2.4	240	2.8	box 10
M00561222	10 mm Rounded	Stiff	10	2.4	240	2.8	box 20
M00561223	10 mm Rounded	Stiff	10	2.4	240	2.8	box 40
M00561231	15 mm Rounded	Stiff	15	2.4	240	2.8	box 10
M00561232	15 mm Rounded	Stiff	15	2.4	240	2.8	box 20
M00561233	15 mm Rounded	Stiff	15	2.4	240	2.8	box 40
M00561241	20 mm Rounded	Stiff	20	2.4	240	2.8	box 10
M00561242	20 mm Rounded	Stiff	20	2.4	240	2.8	box 20
M00561243	20 mm Rounded	Stiff	20	2.4	240	2.8	box 40
M00561191	25 mm Rounded	Stiff	25	2.4	240	2.8	box 10
M00561192	25 mm Rounded	Stiff	25	2.4	240	2.8	box 20
M00561193	25 mm Rounded	Stiff	25	2.4	240	2.8	box 40
M00561291	33 mm Rounded	Stiff	33	2.4	240	2.8	box 10

Tissue Resection

SNARES (CONTINUED)




Profile™ Pediatric Snares



Order Number	Description	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00562531	Extra Small Oval	Flexible	11	1.9	240	2.0	box 10
M00562551	Small Oval	Flexible	13	1.9	240	2.0	box 10
M00562571	Medium Oval	Flexible	27	1.9	240	2.0	box 10

HOT BIOPSY FORCEPS



Radial Jaw™ 4 Hot Biopsy Forceps

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Units
M00515031	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 5
M00515032	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 20
M00515033	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	2.8	box 40

RETRIEVAL DEVICES


Active Cords

Order Number	Description	Units
M00561250	Active Cord, Fits Olympus® Brand Generator	box 1
M00561270	Active Cord, Fits Endostat™, Endostat™ II, Endostat™ III, Bovie®, ValleyLab®, Aspen Labs and ERBE® Brand Generators	box 1



Biopsy

FORCEPS



Radial Jaw™ 4 Standard Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513383	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 5
M00513381	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 20
M00513382	Radial Jaw 4 Standard Capacity	2.2	160	2.8	Yellow	box 40
M00513393	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 5
M00513391	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 20
M00513392	Radial Jaw 4 Standard Capacity w/Needle	2.2	160	2.8	Yellow	box 40
M00513403	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 5
M00513401	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 20
M00513402	Radial Jaw 4 Standard Capacity	2.2	240	2.8	Orange	box 40
M00513413	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 5
M00513411	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 20
M00513412	Radial Jaw 4 Standard Capacity w/Needle	2.2	240	2.8	Orange	box 40

Radial Jaw™ 4 Large Capacity Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	box 40
M00513313	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 5
M00513311	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 20
M00513312	Radial Jaw 4 Large Capacity w/Needle	2.4	160	2.8	Yellow	box 40
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 5
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 20
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	box 40
M00513333	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 5
M00513331	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 20
M00513332	Radial Jaw 4 Large Capacity w/Needle	2.4	240	2.8	Orange	box 40

Biopsy

FORCEPS (CONTINUED)

Radial Jaw™ 4 Jumbo Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	box 5
M00513353	Radial Jaw 4 Jumbo w/Needle	2.8	160	3.2	Yellow	box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 5
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 20
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	box 40
M00513373	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 5
M00513371	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 20
M00513372	Radial Jaw 4 Jumbo w/Needle	2.8	240	3.2	Orange	box 40

Radial Jaw™ 4 Gastropediatric Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00513443	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 5
M00513441	Radial Jaw 4 Gastropediatric	1.8	160	2.0	Yellow	box 20
M00513453	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 5
M00513451	Radial Jaw 4 Gastropediatric w/Needle	1.8	160	2.0	Yellow	box 20



Multibite™ Multiple Sample Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00510101	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	box 5
M00510121	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	box 5



EndoVive™ Enteral Access Devices

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES

Initial Placement



Standard PEGs

Order Number	Description	Units
M00568381	20Fr Standard PEG Kit, Pull	box 2
M00568391	24Fr Standard PEG Kit, Pull	box 2
M00568401	20Fr Standard PEG Kit, Push	box 2
M00568411	24Fr Standard PEG Kit, Push	box 2



Standard PEGs with ENFit

Order Number	Description	Units
M00509081	EndoVive StdPEG Pull w/ENFit 20Fr	box 2
M00509101	EndoVive StdPEG Pull w/ENFit 24Fr	box 2
M00509091	EndoVive StdPEG Push w/ENFit 20Fr	box 2
M00509111	EndoVive StdPEG Push w/ENFit 24Fr	box 2

PEG Accessories

Order Number	Description	Units
M00580641	Y-Port 20Fr	box 10
M00580651	Y-Port 24Fr	box 10
M00580731	20Fr (6.67mm) Round Bolsters	box 10
M00580781	20Fr (6.67mm) and 24F (8.00mm) C-Clamps	box 10


PEG Accessories with ENFit

Order Number	Description	Units
M00501261	EndoVive Y-Port w/ENFit 20Fr	box 10
M00501271	EndoVive Y-Port w/ENFit 24Fr	box 10

EndoVive™ Enteral Access Devices

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Initial Placement (CONTINUED)



One Step Button Low Profile PEG Kits

Order Number	Description	Length (cm)	Units
M00568450	18Fr (6.00mm) One Step Button Kit, Pull	1.7	box 1
M00568460	18Fr (6.00mm) One Step Button Kit, Pull	2.4	box 1
M00568470	18Fr (6.00mm) One Step Button Kit, Pull	3.4	box 1
M00568480	18Fr (6.00mm) One Step Button Kit, Pull	4.4	box 1
M00568500	24Fr (8.00mm) One Step Button Kit, Pull	1.7	box 1
M00568510	24Fr (8.00mm) One Step Button Kit, Pull	2.4	box 1
M00568520	24Fr (8.00mm) One Step Button Kit, Pull	3.4	box 1
M00568530	24Fr (8.00mm) One Step Button Kit, Pull	4.4	box 1
M00568560	18Fr (6.00mm) One Step Button Kit, Push	2.4	box 1
M00568570	18Fr (6.00mm) One Step Button Kit, Push	3.4	box 1
M00568580	18Fr (6.00mm) One Step Button Kit, Push	4.4	box 1
M00568600	24Fr (8.00mm) One Step Button Kit, Push	1.7	box 1
M00568610	24Fr (8.00mm) One Step Button Kit, Push	2.4	box 1
M00568620	24Fr (8.00mm) One Step Button Kit, Push	3.4	box 1
M00568630	24Fr (8.00mm) One Step Button Kit, Push	4.4	box 1

EndoVive™ Enteral Access Devices


PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Initial Placement (CONTINUED)

One Step Button and Low Profile Button Accessories



Order Number	Description	Length (cm)	Colour Code	Units
M00580401	18Fr (6.00mm) Decompression Tube	1.2	Yellow	box 10
M00580411	18Fr (6.00mm) Decompression Tube	1.7	Yellow	box 10
M00580421	18Fr (6.00mm) Decompression Tube	2.4	Yellow	box 10
M00580431	18Fr (6.00mm) Decompression Tube	3.4	Yellow	box 10
M00580441	18Fr (6.00mm) Decompression Tube	4.4	Yellow	box 10
M00580461	24Fr (8.00mm) Decompression Tube	1.7	Blue	box 10
M00580471	24Fr (8.00mm) Decompression Tube	2.4	Blue	box 10
M00580481	24Fr (8.00mm) Decompression Tube	3.4	Blue	box 10
M00580491	24Fr (8.00mm) Decompression Tube	4.4	Blue	box 10
M00580171	18Fr (6.00mm) Bolus Straight Feeding Adapter		Yellow	box 10
M00580181	18Fr (6.00mm) Continuous Right Angle Feeding Adapter		Yellow	box 10
M00580231	24Fr (8.00mm) Bolus Straight Feeding Adapter		Blue	box 10
M00580241	24Fr (8.00mm) Continuous Right Angle Feeding Adapter		Blue	box 10
M00580271	28Fr (9.33mm) Bolus Straight Feeding Adapter		Green	box 10
M00580281	28Fr (9.33mm) Continuous Right Angle Feeding Adapter		Green	box 10



Through-The-PEG (TTP) Jejunal Feeding Tube




Order Number	Description	Tip Style	Units
M00566320	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	Bent Tip	box 1
M00566330	8.5Fr (2.83mm) TTP Jejunal Feeding Tube	Pigtail Tip	box 1
M00566340	12Fr (4.00mm) TTP Jejunal Feeding Tube	Bent Tip	box 1
M00566350	12Fr (4.00mm) TTP Jejunal Feeding Tube	Pigtail Tip	box 1

EndoVive™ Enteral Access Devices

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)

Replacements




Standard Balloon Replacements

Order Number	Description	Units
M00582040	12Fr (4.00mm) Repl. G-Tube, Straight	box 1
M00582050	14Fr (4.67mm) Repl. G-Tube, Straight	box 1
M00582060	16Fr (5.33mm) Repl. G-Tube, Straight	box 1
M00582070	18Fr (6.00mm) Repl. G-Tube, Straight	box 1
M00582080	20Fr (6.67mm) Repl. G-Tube, Straight	box 1
M00582860	22Fr (7.33mm) Repl. G-Tube, Straight	box 1
M00582090	24Fr (8.00mm) Repl. G-Tube, Straight	box 1
M00582100	28Fr (9.33mm) Repl. G-Tube, Straight	box 1
M00582110	12Fr (4.00mm) Repl. G-Tube, Right Angle	box 1
M00582120	14Fr (4.67mm) Repl. G-Tube, Right Angle	box 1
M00582130	16Fr (5.33mm) Repl. G-Tube, Right Angle	box 1
M00582140	18Fr (6.00mm) Repl. G-Tube, Right Angle	box 1
M00582150	20Fr (6.67mm) Repl. G-Tube, Right Angle	box 1
M00582870	22Fr (7.33mm) Repl. G-Tube, Right Angle	box 1
M00582160	24Fr (8.00mm) Repl. G-Tube, Right Angle	box 1
M00582170	28Fr (9.33mm) Repl. G-Tube, Right Angle	box 1

EndoVive™ Enteral Access Devices

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) DEVICES (CONTINUED)


Replacements



Replacement Gastrostomy Feeding Tube with ENFit



Enfit Kit	Enfit Description	Units
M00509881	12Fr (4.00mm) Repl. G-Tube, Straight	box 5
M00509891	14Fr (4.67mm) Repl. G-Tube, Straight	box 5
M00509901	16Fr (5.33mm) Repl. G-Tube, Straight	box 5
M00509911	18Fr (6.00mm) Repl. G-Tube, Straight	box 5
M00509921	20Fr (6.67mm) Repl. G-Tube, Straight	box 5
M00509931	22Fr (7.33mm) Repl. G-Tube, Straight	box 5
M00509941	24Fr (8.00mm) Repl. G-Tube, Straight	box 5
M00509961	12Fr (4.00mm) Repl. G-Tube, Right Angle	box 5
M00509971	14Fr (4.67mm) Repl. G-Tube, Right Angle	box 5
M00509981	16Fr (5.33mm) Repl. G-Tube, Right Angle	box 5
M00509991	18Fr (6.00mm) Repl. G-Tube, Right Angle	box 5
M00510001	20Fr (6.67mm) Repl. G-Tube, Right Angle	box 5
M00510011	22Fr (7.33mm) Repl. G-Tube, Right Angle	box 5
M00510021	24Fr (8.00mm) Repl. G-Tube, Right Angle	box 5





Low Profile Button Replacements



Order Number	Description	Length (cm)	Colour Code	Units
M00568260	18Fr (6.00mm) Low Profile Button Replacement	1.2	Yellow	box 1
M00568270	18Fr (6.00mm) Low Profile Button Replacement	1.7	Yellow	box 1
M00568280	18Fr (6.00mm) Low Profile Button Replacement	2.4	Yellow	box 1
M00568290	18Fr (6.00mm) Low Profile Button Replacement	3.4	Yellow	box 1
M00568300	18Fr (6.00mm) Low Profile Button Replacement	4.4	Yellow	box 1
M00568320	24Fr (8.00mm) Low Profile Button Replacement	1.7	Blue	box 1
M00568330	24Fr (8.00mm) Low Profile Button Replacement	2.4	Blue	box 1
M00568340	24Fr (8.00mm) Low Profile Button Replacement	3.4	Blue	box 1
M00568350	24Fr (8.00mm) Low Profile Button Replacement	4.4	Blue	box 1
M00568170	28Fr (9.33mm) Low Profile Button Replacement	1.5	Green	box 1
M00568180	28Fr (9.33mm) Low Profile Button Replacement	2.8	Green	box 1
M00568190	28Fr (9.33mm) Low Profile Button Replacement	4.3	Green	box 1


For accessories please refer to the One-Step Button and Low Profile Button Accessories





Pulmonary

ALAIR™ BRONCHIAL THERMOPLASTY (BT) SYSTEM



Alair™ Bronchial Thermoplasty (BT) System

Order Number	Model	Description	Active Electrode Length (mm)	Tip Diameter (mm)	Minimum Working Channel (mm)
M005ATS25020	ATS 2-5	Alair BT Catheter – Global	5.0	~1.5	2.0
M005ATS20000	ATS 200	Alair RF Controller – monopolar (return electrode not included)			

Packaged one per box

Alair™ BT Controller Accessory Kits ATS 201*


Order Number	Model	Description
M005ATS20120	ATS 201	British Type G
M005ATS20130	ATS 201	EU Hybrid Type E/F
M005ATS20140	ATS 201	Australian Type I
M005ATS20150	ATS 201	Danish Type K
M005ATS20160	ATS 201	Brazil IEC 60906-1
M005ATS20180	ATS 201	Swiss Type J

Packaged one per box

* Including Power Cord with Molded Plug, Foot Switch, Operating Manual.

Pulmonary

ENDOSCOPIC BRONCHIAL ULTRASOUND



Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle


Olympus Scope Compatible



Order Number	Description	Minimum Working Channel (mm)	Sheath Diameter (mm)	Units
M00558220	Expect Pulmonary 22ga – Olympus	2.0	1.6	box 1
M00558221	Expect Pulmonary 22ga – Olympus	2.0	1.6	box 5
M00558250	Expect Pulmonary 25ga – Olympus	2.0	1.4	box 1
M00558251	Expect Pulmonary 25ga – Olympus	2.0	1.4	box 5
M00558231	Expect Pulmonary Needle Adaptor – Olympus			box 10

Pulmonary

AIRWAY STENTS




Dynamic™ (Y) Stents



Order Number	Description	Tracheal Width (Internal) (mm)	Bronchial Width (Internal) (mm)	Tracheal Length (mm)	Bronchial Lengths (Right/Left) (mm)
M00570670	Dynamic (Y) Stent	11	8	110	25 / 40
M00570680	Dynamic (Y) Stent	13	10	110	25 / 40
M00570690	Dynamic (Y) Stent	15	12	110	25 / 40

Packaged one per box



Ultraflex™ Covered Tracheobronchial NG Stents – Distal Release

Tracheal Stents



Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Length (mm)
M00564820	Ultraflex Tracheal NG Stent	14	80	65
M00564830	Ultraflex Tracheal NG Stent	16	40	25
M00564840	Ultraflex Tracheal NG Stent	16	60	45
M00564850	Ultraflex Tracheal NG Stent	16	80	65
M00564860	Ultraflex Tracheal NG Stent	18	40	25
M00564870	Ultraflex Tracheal NG Stent	18	60	45
M00564880	Ultraflex Tracheal NG Stent	18	80	65
M00564890	Ultraflex Tracheal NG Stent	20	40	25
M00564900	Ultraflex Tracheal NG Stent	20	60	45
M00564910	Ultraflex Tracheal NG Stent	20	80	65

Bronchial Stents

Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Length (mm)
M00564740	Ultraflex Bronchial NG Stent	8	40	25
M00564750	Ultraflex Bronchial NG Stent	10	30	15
M00564760	Ultraflex Bronchial NG Stent	10	40	25
M00564770	Ultraflex Bronchial NG Stent	12	30	15
M00564780	Ultraflex Bronchial NG Stent	12	40	25
M00564790	Ultraflex Bronchial NG Stent	14	30	15
M00564800	Ultraflex Bronchial NG Stent	14	40	25
M00564810	Ultraflex Bronchial NG Stent	14	60	45

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-22F (5.00-7.33mm)

Pulmonary

AIRWAY STENTS (CONTINUED)



Ultraflex™ Non-covered Tracheobronchial NG Stents – Proximal Release

Tracheal Stents



Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)
M00564710	Ultraflex Tracheal NG Stent	14	40
M00564720	Ultraflex Tracheal NG Stent	14	60


Bronchial Stents

Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)
M00564640	Ultraflex Bronchial NG Stent	8	20
M00564650	Ultraflex Bronchial NG Stent	8	40
M00564660	Ultraflex Bronchial NG Stent	10	20
M00564670	Ultraflex Bronchial NG Stent	10	40
M00564680	Ultraflex Bronchial NG Stent	12	20
M00564690	Ultraflex Bronchial NG Stent	12	40
M00564700	Ultraflex Bronchial NG Stent	14	20

Packaged one per box

Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-22F (5.00-7.33mm)



Ultraflex™ Non-covered Tracheobronchial NG Stents – Distal Release

Tracheal Stents



Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)
M00564530	Ultraflex Tracheal NG Stent	16	40
M00564540	Ultraflex Tracheal NG Stent	16	60
M00564560	Ultraflex Tracheal NG Stent	18	40
M00564570	Ultraflex Tracheal NG Stent	18	60
M00564590	Ultraflex Tracheal NG Stent	20	40
M00564600	Ultraflex Tracheal NG Stent	20	60

Bronchial Stents

Order Number	Description	Maximum Stent O.D. (mm)	Expanded Stent Length (mm)
M00564500	Ultraflex Bronchial NG Stent	10	30
M00564510	Ultraflex Bronchial NG Stent	12	30
M00564520	Ultraflex Bronchial NG Stent	14	30

Packaged one per box


Recommended Guidewire: .035in (0.89mm) Jagwire™ Pulmonary Guidewire, M00515171

Additional Specifications for Ultraflex Tracheobronchial Stents: Delivery System Working Length is 95cm, Maximum Catheter Shaft O.D. is 16F (5.33mm), Compressed Stent on catheter yields Delivery System O.D. range 15-22F (5.00-7.33mm)



Pulmonary

AIRWAY STENTS (CONTINUED)



Polyflex Stent™ Self-expanding Silicone Airway Stents


Order Number	Description	Stent I.D. (mm)	Stent Length (mm)	Delivery System Diameter (mm)
M00570000	Polyflex Airway Stent	8	20	7
M00570010	Polyflex Airway Stent	8	30	7
M00570020	Polyflex Airway Stent	10	20	8
M00570030	Polyflex Airway Stent	10	30	8
M00570040	Polyflex Airway Stent	10	40	8
M00570050	Polyflex Airway Stent	10	50	8
M00570060	Polyflex Airway Stent	12	20	9
M00570070	Polyflex Airway Stent	12	30	9
M00570080	Polyflex Airway Stent	12	40	9
M00570090	Polyflex Airway Stent	12	50	9
M00570100	Polyflex Airway Stent	14	20	9
M00570110	Polyflex Airway Stent	14	30	9
M00570120	Polyflex Airway Stent	14	40	9
M00570130	Polyflex Airway Stent	14	50	9
M00570140	Polyflex Airway Stent	14	60	9
M00570150	Polyflex Airway Stent	16	30	10
M00570160	Polyflex Airway Stent	16	40	10
M00570170	Polyflex Airway Stent	16	50	10
M00570180	Polyflex Airway Stent	16	60	10
M00570190	Polyflex Airway Stent	16	70	10
M00570200	Polyflex Airway Stent	18	30	11
M00570210	Polyflex Airway Stent	18	40	11
M00570220	Polyflex Airway Stent	18	50	11
M00570230	Polyflex Airway Stent	18	60	11
M00570240	Polyflex Airway Stent	18	70	11
M00570250	Polyflex Airway Stent	18	80	11
M00570260	Polyflex Airway Stent	20	40	12
M00570270	Polyflex Airway Stent	20	50	12
M00570280	Polyflex Airway Stent	20	60	12
M00570290	Polyflex Airway Stent	20	70	12
M00570300	Polyflex Airway Stent	20	80	12
M00570310	Polyflex Airway Stent	22	50	13
M00570320	Polyflex Airway Stent	22	60	13
M00570330	Polyflex Airway Stent	22	80	13

Packaged one per box



Pulmonary

BIOPSY FORCEPS




Radial Jaw™ 4 Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515181	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Large Capacity	2.2	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Large Capacity	2.2	100	2.8	Blue	box 20

TRANSBRONCHIAL ASPIRATION NEEDLES



eXcelon™ Transbronchial Aspiration Needle with Syringe



Order Number	Description	Needle			Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Units
		Gauge	O.D. (mm)	I.D. (mm)				
M00564101	eXcelon Transbronchial Aspiration Needle w/Syringe	19	1.07	0.69	15	130	1.8	box 5
M00564111	eXcelon Transbronchial Aspiration Needle w/Syringe	20	0.90	0.58	15	130	1.8	box 5
M00564121	eXcelon Transbronchial Aspiration Needle w/Syringe	21	0.81	0.50	15	130	1.8	box 5

CYTOLOGY BRUSHES


Celebrity™ Endoscopic Cytology Brushes – Bronchoscope



Order Number	Description	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Units
M00516001	Celebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Celebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Celebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Celebrity Cytology Brush	2.0	1.9	150	box 10

Pulmonary

AIRWAY RETRIEVAL BASKETS



ZeroTip™ Airway Retrieval Baskets



Order Number	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Sheath Material
M00513200	ZeroTip Airway Retrieval Basket	0.8	120	12	Polyimide/PTFE
M00513210	ZeroTip Airway Retrieval Basket	1.0	120	16	Polyimide/PTFE

Packaged one per box

DILATION BALLOONS



CRE™ Pulmonary Dilation Balloons



Order Number	Description	Balloon Length (cm)	Initial Balloon O.D. (mm)	Intermediate Balloon O.D. (mm)	Maximum Balloon O.D. (mm)	Working Catheter Length (cm)
M00550300	CRE Balloon	5.5	12	13.5	15	75
M00550310	CRE Balloon	5.5	15	16.5	18	75
M00550320	CRE Balloon	5.5	18	19.0	20	75
M00550330	CRE Balloon	3.0	8	9.0	10	75
M00550340	CRE Balloon	3.0	10	11.0	12	75
M00550350	CRE Balloon	3.0	12	13.5	15	75

Packaged one per box


GUIDEWIRES

Jagwire™ Pulmonary Guidewire



Order Number	Description	O.D. (in) / (mm)	Length (cm)	Units
M00515171	Jagwire	.035 / 0.89	180	box 2

INFLATION SYSTEM



Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance™ II Inflation Handle	box 1
M00550601	⊗ 60ml Syringe/Gauge Assembly	box 5



Find more information on Boston Scientific products
by using your smartphone to scan this code or visit
www.bostonscientific.com/global-endoscopy



For single use only. Do not reuse.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.
ENDO - 431603 - AA Jan 2017 Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

Boston Scientific International SA
Parc Val Saint Quentin - Bâtiment H
78960 Voisins-le-Bretonneux - France
+33 1 39 30 97 00

© 2017 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINEND2396EA

