

Spectra™

Concealable Penile Prosthesis

The Spectra Prosthesis designed for simplified sizing, ease of placement, excellent concealment and rigidity.¹




About the device¹

- ▶ Provides a balance between rigidity and concealment
- ▶ Unique spring and cable design
- ▶ Multiple configurations for a wide range of patient sizes
- ▶ Easy attachment of RTE to cylinder
- ▶ **Diameter:** 9.5 mm, 12 mm and 14 mm
- ▶ **Length:** 12 cm, 16 cm or 20 cm
- ▶ Snap-fit RTEs optimize fit and provide a secure, rigid connection
- ▶ Simplified implantation technique and easy for patients to use
- ▶ Provides the ability to position with minimal spring back
- ▶ **Snap-fit Rear Tip Extenders (RTEs):** add length to the proximal cylinder end of 0.5 to 7.5 cm (excluding 7.0 cm) in 0.5 cm increments for customized fit


Important safety information²

- ▶ Not for patients with urogenital infections or active skin infections in region of surgery
- ▶ Will make latent natural or spontaneous erections, as well as other interventional treatment options impossible
- ▶ Penile shortening, curvature or scarring may result
- ▶ Patients with diabetes, spinal cord injuries or open sores may have an increased risk of infection associated with a prosthesis

Excellent concealment



Simplified sizing



Spectra™ Concealable Penile Prosthesis cylinder length + Snap-fit RTE = Total intracorporal length

Dimensional specifications

Diameter	Length	Product Number
9.5 mm	12 cm	720054-01
9.5 mm	16 cm	720054-02
9.5 mm	20 cm	720054-03
12 mm	12 cm	720074-01
12 mm	16 cm	720074-02
12 mm	20 cm	720074-03
14 mm	12 cm	720056-01
14 mm	16 cm	720056-02
14 mm	20 cm	720056-03

Snap-fit RTEs are used with Spectra cylinders.

Cylinder Diameter	RTE Length(s)	RTE Product Number
9.5 mm	0.5, 1, 1.5, 2, 3 cm (kit)	720170-01
9.5 mm	4, 5, 6 cm (kit)	720171-01
12 mm & 14 mm	0.5 cm	72404320
12 mm & 14 mm	1 cm	72404321
12 mm & 14 mm	2 cm	72404322
12 mm & 14 mm	3 cm	72404323
12 mm & 14 mm	4 cm	72404324
12 mm & 14 mm	5 cm	72404325
12 mm & 14 mm	6 cm	72404326
12 mm & 14 mm	1.5 cm stackable	72404330

The 1.5 cm RTE is the only RTE that can be used in combination with another Snap-fit RTE.

Sizers (Non-sterile)	Product Number
9.5 mm and 12 mm diameters	720153-01

9.5 mm and 12 mm sizers are packaged together.

1. Data on file with Boston Scientific.

2. Spectra Concealable Penile Prosthesis Instructions for Use. American Medical Systems, Inc. 2013.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Caution: The law restricts these devices to sale by or on the order of a physician.

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

The AMS Spectra™ Concealable Penile Prosthesis is intended for use in the treatment of male erectile dysfunction (impotence). Implanting a penile prosthesis will damage or destroy any remaining natural ability to have a natural erection, as well as make other treatment options impossible.

Men with diabetes, spinal cord injuries or skin infections may have an increased risk of infection. Implantation may result in penile shortening, curvature or scarring.

Potential adverse events may include device malfunction/failure leading to additional surgery, device/tissue erosion, infection, and pain/soresness. MH-547821-AA