

# LATITUDE™ NXT

## PATIENT MANAGEMENT SYSTEM



**LATITUDE Wireless Communicator**  
Models 6280 and 6290

### Model 6290 Communicator

#### CRT-D

G124 MOMENTUM™ CRT-D  
G125 MOMENTUM CRT-D  
G126 MOMENTUM CRT-D  
G138 MOMENTUM X4 CRT-D  
G140 INOGEN™ CRT-D  
G141 INOGEN CRT-D  
G146 INOGEN X4 CRT-D  
G148 INOGEN X4 CRT-D  
G150 DYNAGEN™ CRT-D  
G151 DYNAGEN CRT-D  
G154 DYNAGEN CRT-D  
G156 DYNAGEN X4 CRT-D

G158 DYNAGEN X4 CRT-D  
G247 VIGILANT™ X4 CRT-D  
G248 VIGILANT X4 CRT-D  
G447 RESONATE™ X4 CRT-D

#### S-ICD

A209 EMBLEM™ S-ICD  
A219 EMBLEM MRI S-ICD

#### ICD

D010 INOGEN MINI ICD VR  
D011 INOGEN MINI ICD VR  
D012 INOGEN MINI ICD DR  
D013 INOGEN MINI ICD DR  
D020 DYNAGEN MINI ICD VR  
D021 DYNAGEN MINI ICD VR  
D022 DYNAGEN MINI ICD DR  
D023 DYNAGEN MINI ICD DR

D120 MOMENTUM EL ICD VR  
D121 MOMENTUM EL ICD DR  
D140 INOGEN EL ICD VR  
D141 INOGEN EL ICD VR  
D142 INOGEN EL ICD DR  
D143 INOGEN EL ICD DR  
D150 DYNAGEN EL ICD VR  
D151 DYNAGEN EL ICD VR  
D152 DYNAGEN EL ICD DR  
D153 DYNAGEN EL ICD DR  
D232 VIGILANT EL ICD VR  
D233 VIGILANT EL ICD DR  
D400 PERCIVA™ ICD VR  
D401 PERCIVA ICD DR  
D412 PERCIVA ICD VR  
D413 PERCIVA ICD DR  
D432 RESONATE EL ICD VR  
D433 RESONATE EL ICD DR

#### Pacemakers

L100 ESSENTIO™ SR  
L101 ESSENTIO DR  
L110 ESSENTIO SR  
L111 ESSENTIO DR  
L121 ESSENTIO DR EL  
L131 ESSENTIO DR EL  
L300 ACCOLADE™ SR  
L301 ACCOLADE DR  
L310 ACCOLADE SR  
L311 ACCOLADE DR  
L321 ACCOLADE DR EL  
L331 ACCOLADE DR EL

#### CRT-P

U125 VALITUDE™ CRT-P  
U128 VALITUDE X4 CRT-P  
U225 VISIONIST™ CRT-P  
U226 VISIONIST CRT-P  
U228 VISIONIST X4 CRT-P



Available for in-clinic or direct-to-patient distribution via LATITUDE NXT.

### Model 6280 Communicator

#### CRT-D

N050 PUNCTUA™ CRT-D  
N051 PUNCTUA CRT-D  
N118 COGNIS™ CRT-D  
N119 COGNIS CRT-D  
N140 ENERGEN™ CRT-D  
N141 ENERGEN CRT-D  
N160 INCEPTA™ CRT-D  
N161 INCEPTA CRT-D  
N164 INCEPTA CRT-D

#### ICD

E050 PUNCTUA ICD VR  
E051 PUNCTUA ICD VR  
E052 PUNCTUA ICD DR  
E053 PUNCTUA ICD DR  
E102 TELIGEN™ ICD VR  
E110 TELIGEN ICD DR  
E140 ENERGEN ICD VR  
E141 ENERGEN ICD VR  
E142 ENERGEN ICD DR

E143 ENERGEN ICD DR  
E160 INCEPTA ICD VR  
E161 INCEPTA ICD VR  
E162 INCEPTA ICD DR  
E163 INCEPTA ICD DR

#### Pacemakers

K062 ADVANTIO™ SR  
K063 ADVANTIO DR  
K064 ADVANTIO DR EL  
K172 INGENIO™ SR

K173 INGENIO DR  
K174 INGENIO DR EL  
K272 VITALIO™ SR  
K273 VITALIO DR  
K274 VITALIO DR EL

#### CRT-P

V172 INVIVE™ CRT-P  
V173 INVIVE CRT-P  
V272 INTUA™ CRT-P  
V273 INTUA CRT-P



Available for direct-to-patient distribution via LATITUDE NXT.

### Communicator Connectivity Options

#### Internet Adapters

#### USB Cellular Adapters

#### Landline (phone cord included with Communicator)

Refer to the product labeling for specific instructions for use. Rx only.