

Captivator™ COLD Single-use Snare

Boston
Scientific
Advancing science for life™

Round Out Your Treatment Options
with Captivator™ COLD Snare



Optimized for Cold Snaring

[About](#)

[Features](#)

[Cold Snare
Polypectomy](#)

[Ordering](#)

The Captivator COLD Single-use Snare:

1 **Optimized for cold snare polypectomy** and designed for the effective resection of flat lesions and diminutive polyps.

2 **Builds on the strong clinical track record of Captivator II Snares** while offering a design optimized for cold snaring of diminutive polyps.

3 **Is part of a comprehensive detection to resection offering** of devices for biopsy, tissue resection and hemostasis.



Optimized for cold resection.



Round Shape

- Designed for concentric resection and capture of small polyps.



Thin Wire and Robust Sheath



- Designed to enable a precise and clean resection of polyps.



Reinforced Wire Stiffness and Loop Construction

- Designed to provide greater tissue traction and device control for the effective resection of flat lesions.
- Designed to facilitate improved physician control and visualization throughout polyp resection.

Images owned by Boston Scientific.



**Distinctive 1x3
braid pattern wire
construction**

Images are 22x magnification.

About

Features

Cold Snare
Polypectomy

Ordering



Cold Snare Polypectomy



A DISCUSSION ON COLD SNARING: “The Cold Revolution in Polypectomy”

Presented by Douglas Rex, M.D.

Indiana University Health, Indianapolis, Indiana, U.S.

Watch it on www.EndoSuite.com.

About


Features

Cold Snare
Polypectomy

Ordering



Ordering Information



Order Number	Description	Loop Width (mm)	Sheath OD (mm)	Working Length (cm)	Min. Working Channel (mm)	Packaging
M00561101	10mm Rounded – Stiff	10mm	2.4	240	2.8	Box 10
M00561102	10mm Rounded – Stiff	10mm	2.4	240	2.8	Box 20
M00561103	10mm Rounded – Stiff	10mm	2.4	240	2.8	Box 40

[About](#)

[Features](#)

[Cold Snare Polypectomy](#)

[Ordering](#)





Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Images owned by Boston Scientific.

All trademarks are the property of their respective owners.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

©2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-555512-AA June 2018