


WallFlex™ Colonic and Duodenal Soft Stent System with Anchor Lock Delivery System

**Boston
Scientific**
Advancing science for life™

Built on the trusted WallFlex™ platform, the WallFlex Colonic and Duodenal Soft Stent Systems are designed to be a more flexible option to relieve malignant obstructions of the colon and duodenum.




Softer Stent

Designed with a lower axial force which may better conform to torturous anatomy such as the sigmoid, flexures, and duodenal bulb.

Delivery System

New low profile 9Fr delivery system, allowing for placement via a 3.2mm channel scope, which may assist in placement in tighter strictures.



Flexibility

With lower radial and axial forces*, WallFlex Colonic and Duodenal Soft Stents are designed for increased flexibility.

Comprehensive Metal Stent Portfolio

Boston Scientific's robust line of esophageal, biliary and enteral metal stents provides physicians the tools needed to help them effectively manage strictures. BSC's metal stents have been studied in more than 500 publications over the past 25 years, providing physicians with confidence in the therapeutic performance of the portfolio. The addition of WallFlex Enteral Soft adds yet another option designed to meet clinical needs.

*as compared to the original WallFlex Colonic and Duodenal Stent System

WallFlex™ Colonic and Duodenal Soft Stent System with Anchor Lock Delivery System

The WallFlex Colonic Soft Stent is indicated for the palliative treatment of colonic strictures caused by malignant neoplasm and to relieve large bowel obstruction prior to colectomy in patients with malignant strictures. The WallFlex Enteral Duodenal Stent Soft is indicated for palliative treatment of gastroduodenal obstructions produced by malignant neoplasms.

Ordering Information

WallFlex Colonic Soft Stent System with Anchor Lock Delivery System

| Order Number | Diameter (mm) Flare/ Body | Stent Length (cm) | Working Length (cm) | Delivery System Diameter/Working Channel (F/mm) | Guidewire Diameter (inches/mm) |
|--------------|---------------------------|-------------------|---------------------|---|--------------------------------|
| M00555760 | 27 / 22 | 6 | 230 | 9 / 3.2 | 0.035 / 0.89 |
| M00555770 | 27 / 22 | 9 | 230 | 9 / 3.2 | 0.035 / 0.89 |
| M00555780 | 27 / 22 | 12 | 230 | 9 / 3.2 | 0.035 / 0.89 |

WallFlex Duodenal Soft Stent System with Anchor Lock Delivery System

| Order Number | Diameter (mm) Flare/ Body | Stent Length (cm) | Working Length (cm) | Delivery System Diameter/Working Channel (F/mm) | Guidewire Diameter (inches/mm) |
|--------------|---------------------------|-------------------|---------------------|---|--------------------------------|
| M00555860 | 27 / 22 | 6 | 230 | 9 / 3.2 | 0.035 / 0.89 |
| M00555870 | 27 / 22 | 9 | 230 | 9 / 3.2 | 0.035 / 0.89 |
| M00555880 | 27 / 22 | 12 | 230 | 9 / 3.2 | 0.035 / 0.89 |

Recommended guidewires

0.035" WallStent Super Stiff – Order Number H965180011

0.035" Dreamwire Stiff – Order Number M00556161

**Boston
Scientific**
Advancing science for life™

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

All photographs owned by Boston Scientific.
All trademarks are the property of their respective owners.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com/gastro

©2020 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-815206-AA