

Ureteral Catheter Portfolio

Begin your procedures with confidence



Dual Lumen Ureteral Catheter

By helping deliver irrigation fluids to the ureter and drainage of fluids from the urinary tract, our dual-purpose catheter provides access, advancement or exchange of ancillary devices including guidewires.

- ▶ Designed for over-the-wire placement
- ▶ Provides two separate working channels for second wire placement or injection



Axxcess™ Ureteral Catheter

The single lumen catheter with injection hub facilitates injection of contrast or establishment of post-procedure drainage collection.

- ▶ Tapered tip to aid atraumatic access for delivery of contrast or a guidewire
- ▶ Graduated markers to facilitate placement
- ▶ Radiopaque Percuflex™ material to aid in visualization and placement



Flexima™ Ureteral Catheter

The single lumen with injection hub provides access and collection.

- ▶ Various tip designs to facilitate access to the urinary tract
- ▶ Can be placed retrograde
- ▶ Radiopaque Flexima™ material to aid in visualization and placement



Imager™ II Ureteral Catheter

The single lumen, torqueable catheter is designed to facilitate access in PCNL and ureteroscopic procedures.

- ▶ Biocompatible polymer reinforced with a stainless steel braided wire
- ▶ Can be placed retrograde or antegrade

Ordering Information

Dual Lumen Ureteral Catheter

Order Number	UPN	Recommended Guidewire (in)	O.D. Size	Working Length	Injection Lumen	Quantity
405100	M0064051000	0.038	10F	54	0.05	Each

Flexima™ Ureteral Catheter

Order Number	UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
400200	M0064002001	Open End	4F	70	-	Box 20
400201	M0064002011	Open End	5F	70	-	Box 20
400202	M0064002021	Open End	6F	70	-	Box 20
400203	M0064002031	Open End	7F	70	-	Box 20
400204	M0064002041	Open End	8F	70	-	Box 20
400211	M0064002111	Cone Tip	5F	70	8F	Box 20

Axxcess™ Ureteral Catheter

Order Number	UPN	Description	Catheter O.D. Size	Working Length	Quantity
400116	M0064001161	Open End	6F	70	Box 10

Imager™ II Ureteral Catheter

Order Number	UPN	Description	O.D. Size	Working Length	Quantity
400300	M0064003001	C1	5F	65	Box 5
400303	M0064003031	Bern	5F	65	Box 5
400503	M0064005031	Bern	5F	40	Box 5



Please consult your sales representative for more information and ordering details.

To stay up to date on stone management best practices, efficiencies and innovations that have the potential to transform your practice, [visit StoneSmart.com](https://www.stonesmart.com).

Boston Scientific
Advancing science for life™

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.
CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.
Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.
All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.BostonScientific.com

©2022 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-487903-AD AUG 2022