

EUS

**Expect™ Endoscopic Ultrasound
Aspiration Needles**

- Expect Endoscopic Ultrasound Aspiration Needle 43
- Expect Slimline (SL) Endoscopic Ultrasound
Aspiration Needle 43



**Acquire™ Endoscopic Ultrasound
Fine Needle Biopsy (FNB) Device**

- Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB)
Device..... 44

AXIOS™ System

- AXIOS Stent and Electrocautery Enhanced Delivery System ... 44
- AXIOS Stent and Delivery System..... 44

Expect™ Endoscopic Ultrasound Aspiration Needles



Expect Endoscopic Ultrasound Aspiration Needle*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00550000	08714729790051	19 gauge	2.8	1.83	Each
M00550010	08714729790044	22 gauge	2.4	1.65	Each
M00550020	08714729790020	25 gauge	2.4	1.52	Each
M00550040	08714729797463	19 gauge Flex	2.8	1.73	Each
M00550001	08714729790068	19 gauge	2.8	1.83	Box 5
M00550011	08714729790013	22 gauge	2.4	1.65	Box 5
M00550021	08714729790037	25 gauge	2.4	1.52	Box 5
M00550041	08714729797470	19 gauge Flex	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable



Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle*



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555500	08714729848141	19 gauge (1.10mm)	2.8	1.83	Each
M00555510	08714729848158	22 gauge (0.72mm)	2.4	1.65	Each
M00555520	08714729848165	25 gauge (0.52mm)	2.4	1.52	Each
M00555530	08714729848172	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555501	08714729848714	19 gauge (1.10mm)	2.8	1.83	Box 5
M00555511	08714729848721	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555521	08714729848738	25 gauge (0.52mm)	2.4	1.52	Box 5
M00555531	08714729848806	19 gauge Flex (1.14mm)	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device



Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555540	08714729931782	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555560	08714729931805	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555541	08714729931843	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555561	08714729931867	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

AXIOS™ System



AXIOS Stent and Electrocautery Enhanced Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553640	21	10	10	10.8	138	146	Each
M00553650	24	15	10	10.8	138	146	Each

AXIOS Stent and Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553740	21	10	10	10.8	138	146	Each
M00553750	24	15	10	10.8	138	146	Each