


Resection & Retrieval



Resection Devices

EMR for Upper GI

➤ Captivator™ EMR Devices 75

Stiff Snares

➤ Captivator II Single-Use Snares..... 75

➤ Captivator Single-Use Snares..... 76

Speciality Snares

➤ Profile™ Single-Use Snares 77

➤ Single-Use Rotatable Snares..... 77

Snares

➤ Captiflex™ Single-Use Snares 78

➤ Sensation™ Single-Use Short Throw Snares..... 78

Active Cords

➤ Active Cords 79

Retrieval Devices

Retrieval Net

➤ RescueNet™ Retrieval Net..... 80

Retrieval Grasping Forceps

➤ Rescue Rat Tooth Grasping Forceps..... 80

➤ Rescue Rat Tooth Grasping Forceps - Pediatric..... 80

➤ Rescue Alligator Grasping Forceps, Long 80

➤ Rescue Alligator Grasping Forceps..... 80

➤ Rescue Rat Tooth / Alligator Grasping Forceps..... 80

Polyp Traps


➤ TrapEase Single Chamber Poly Trap..... 80

➤ TrapEase Quad Chamber Poly Trap 80



Resection Devices

EMR for Upper GI




Captivator™ EMR Device



Order Number	GTIN	Description	Packaging
M00561600	08714729842675	Captivator EMR Device for Standard Gastrosopes	Box 1
M00561610	08714729842699	Captivator EMR Device for Large Gastrosopes	Box 1
M00561601	08714729842682	Captivator EMR Device for Standard Gastrosopes	Box 5
M00561611	08714729842705	Captivator EMR Device for Large Gastrosopes	Box 5

Stiff Snares



Captivator II Single-Use Snares




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561221	08714729855903	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561222	08714729861300	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561223	08714729861317	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40
M00561231	08714729855910	15mm Rounded – Stiff	15	2.4	2.8	240	Box 10
M00561232	08714729861331	15mm Rounded – Stiff	15	2.4	2.8	240	Box 20
M00561233	08714729861348	15mm Rounded – Stiff	15	2.4	2.8	240	Box 40
M00561241	08714729855927	20mm Rounded – Stiff	20	2.4	2.8	240	Box 10
M00561242	08714729861362	20mm Rounded – Stiff	20	2.4	2.8	240	Box 20
M00561243	08714729861379	20mm Rounded – Stiff	20	2.4	2.8	240	Box 40
M00561191	08714729855934	25mm Rounded – Stiff	25	2.4	2.8	240	Box 10
M00561192	08714729861270	25mm Rounded – Stiff	25	2.4	2.8	240	Box 20
M00561193	08714729861287	25mm Rounded – Stiff	25	2.4	2.8	240	Box 40
M00561291	08714729645801	Extra Large Rounded – Stiff	33	2.4	2.8	240	Box 10



Resection Devices


Stiff Snares (Continued)



Captivator™ Single-Use Snares									
Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging		
M00562301	08714729019312	Small Oval – Stiff	13	2.4	2.8	240	Box 10		
M00562321	08714729019312	Medium Oval – Stiff	27	2.4	2.8	240	Box 10		
M00561311	08714729019251	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 10		
M00562391	08714729019404	Large Oval (thin wire) – Flexible	30	2.4	2.8	240	Box 10		
M00562451	08714729071068	Small Hexagonal – Stiff	13	2.4	2.8	240	Box 10		
M00562341	08714729019350	Medium Hexagonal – Stiff	27	2.4	2.8	240	Box 10		
M00562371	08714729019381	Medium Crescent – Stiff	27	2.4	2.8	240	Box 10		

Resection Devices


Speciality Snares



Profile™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562531	08714729156628	Extra Small Oval – Flexible	11	1.9	2.0	240	Box 10
M00562551	08714729156659	Small Oval – Flexible	13	1.9	2.0	240	Box 10
M00562571	08714729156635	Medium Oval – Flexible	27	1.9	2.0	240	Box 10



Single-Use Rotatable Snares*




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	08714729283980	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	08714729283997	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

*With Star Catheter™ Technology.

Resection Devices


Snares



Captiflex™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562471	08714729019459	Extra Small Oval – Flexible	11	2.4	2.8	240	Box 10
M00562421	08714729019435	Small Oval – Flexible	13	2.4	2.8	240	Box 10
M00562422	08714729501657	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00562401	08714729019411	Medium Oval – Flexible	27	2.4	2.8	240	Box 10
M00562402	08714729501640	Medium Oval – Flexible	27	2.4	2.8	240	Box 40



Sensation™ Single-Use Short Throw Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562691	08714729158158	Small Oval – Flexible	13	2.4	2.8	240	Box 5
M00562692	08714729268826	Small Oval – Flexible	13	2.4	2.8	240	Box 20
M00562693	08714729501671	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00560321	08714729283928	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 5
M00562671	08714729158141	Medium Oval – Flexible	27	2.4	2.8	240	Box 5
M00562672	08714729268819	Medium Oval – Flexible	27	2.4	2.8	240	Box 20
M00562673	08714729501664	Medium Oval – Flexible	27	2.4	2.8	240	Box 40
M00560311	08714729283904	Medium Oval – Medium Stiff	27	2.4	2.8	240	Box 5
M00562651	08714729158110	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 5
M00562652	08714729268802	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 20
M00560331	08714729283935	Medium Crescent – Medium Stiff	27	2.4	2.8	240	Box 5



Section Table of Contents

Prev Page

78

Next Page

Main Table of Contents



Resection Devices


Active Cords

Active Cords			
Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus™ Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each

[Section Table of Contents](#)[Prev Page](#)**79**[Next Page](#)[Main Table of Contents](#)


Retrieval Devices

Retrieval Net



RescueNet™ Retrieval Net							
Order Number	GTIN	Description	Net Size (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGN-538-5	10840253107425	RescueNet Retrieval Net	30	2.5	2.8	230	Box 5

Retrieval Grasping Forceps



Rescue™ Retrieval Single-Use Graspers							
Order Number	GTIN	Description	Jaw Opening Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGR-276-5	00816849017334	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	Box 5
DGR-277-5	00816849017341	Rescue Rat Tooth Grasping Forceps - Pediatric	4.5	1.8	2.0	180	Box 5
DGA-275-5	00816849017310	Rescue Alligator Grasping Forceps, Long	9	2.4	2.8	230	Box 5
DGA-274-5	00816849017303	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	Box 5
DGC-383-5	00816849017327	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	Box 5

Polyp Traps



TrapEase™ Polyp Traps							
Order Number	GTIN	Description					Packaging
STE-297-50	00816849010182	TrapEase Single Chamber Polyp Trap					Box 50
STE-298-25	00816849010199	TrapEase Quad Chamber Polyp Trap					Box 25

